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BIOTRONIK

excellence for life

BIOTRONIK AG
Ackerstrasse 6
CH-8180 Bulach
Switzerland

Tel. +41 44 864 51 11

Fax +41 44 864 50 05

27 March 2020

Manufacturer's Statement

To whom this may concern:

We hereby guarantee that the enclosed document:


Instruction for Use
Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System

is an original document issued by BIOTRONIK AG.

Name: Susanne Minich
Position: Regulatory Affairs Manager

Company: BIOTRONIK AG
Street: Ackerstrasse 6
City / Country: 8180 Bulach - Switzerland

Date and Signature:

Susanne Minich

BIOTRONIK

excellence for life

BIOTRONIK AG
Ackerstrasse 6
CH-8180 Bulach

Official Certification

Seen for authentication of the reverse side signature, acknowledged in our presence by

Ms. Susanne MINICH, born 12th May 1965, Swiss citizen of Wattwil SG, according to her information residing at Grossackerstrasse 81, 8041 Zürich, identified by identity card.

This legalization refers only to the authentication of the signature and not to the contents or validity of the document.

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M. Hofstetter

Instructions For Use

Sirolimus-eluting Resorbable Coronary Magnesium Scaffold System

Compliance Chart

Inflation Pressure (IP)		Scaffold Diameter (mm)(SD)		
	atm	(kPA)	3.0	3.5
	8	(811)	2.85	3.38
	9	(912)	2.92	3.46
NP	10	(1013)	3.00	3.54
	11	(1115)	3.06	3.60
	12	(1216)	3.12	3.66
	13	(1317)	3.18	3.71
	14	(1419)	3.22	3.75
	15	(1520)	3.26	3.79
RBP	16	(1621)	3.29	3.82
	17	(1723)	3.33	3.87
	18	(1824)	3.37	3.90
NP	In vitro tests have shown that the balloons will reach their nominal size at given nominal pressure.			
RBP	In vitro tests have shown that with 95% confidence, 99.9% of the balloons will not burst at or below rated burst pressure. DO NOT exceed RBP.			
Note	Diameters at inflation pressures above RBP are provided for information only. At the pressures above RBP shown in this table, 99% of the balloons with 95% confidence will not burst.			

Maximum Diameter for Post-dilation

Scaffold Diameter (mm) (SD)	Maximum Diameter for Post-dilation (mm)
3.0	3.6
3.5	4.1

Sizes and Associated Total Drug Loads

Length (mm) (L)	Scaffold Diameter (mm) (SD)	
	3.0	3.5
15	166 µg	159 µg
20	226 µg	223 µg
25	287 µg	271 µg

Vessel Sizing

Scaffold Diameter (mm) (SD)	Recommended vessel diameter (mm) (RVD)
3.0	2.7 - 3.2
3.5	3.2 - 3.7

Symbol Legend

Sterilized using ethylene oxide

Keep dry

Do not resterilize

Consult instructions for use

Do not reuse

Keep away from sunlight

Date of manufacture

Manufacturer

Caution

Use by

Batch code

Catalogue number

Upper Temperature limitation

Do not use if package is damaged


MR conditional

Nominal pressure

Rated burst pressure

Patient Implant Card
Inside

Description

Magmaris Sirolimus-eluting Resorbable Coronary Magnesium Scaffold System (hereinafter Magmaris) is intended to improve coronary luminal diameter in patients with symptomatic ischemic heart disease due to stenotic coronary artery **lesions**.

Magmaris is a drug-eluting balloon-expandable resorbable scaffold mounted on a rapid-exchange delivery system.

The scaffold is made from a Magnesium alloy. It has permanent double-eyelet radiopaque markers made from Tantalum at each end. The scaffold itself is radiolucent. Only scaffold markers are visible under fluoroscopy for precise location. The scaffold is covered with a drug-polymer coating consisting of the antiproliferative drug substance Sirolimus formulated in a bioresorbable PLLA (Poly-L-Lactide).

The nominal drug content on each scaffold is 1.4µg Sirolimus per mm² of the scaffold surface area. A matrix of the available product range and nominal total drug loads for each size are provided in the "Sizes and Associated Total Drug Loads" table on page 2.

Magmaris utilizes a rapid-exchange delivery system with a working length of 140 cm. The scaffold is positioned between the distal and proximal radiopaque balloon markers. An inflation/deflation device with 20 ml capacity can be connected to the proximal end of the luer port to inflate/deflate the balloon. The delivery system has a hydrophobic coating on the outer surface of the proximal shaft and a hydrophilic coating on the outer surface of the distal shaft. Shaft markers are located on the hypotube at two locations to indicate catheter length for brachial and femoral access routes.

Magmaris is compatible with most guide wires with a **maximum diameter of 0.014" (0.36 mm)** and guiding catheters with a minimum inner diameter of $\geq 0.070"$ (1.78 mm; 6F) (see "Warnings" section for exceptions).

To indicate when the delivery system tip exits from the guiding catheter, shaft exit markers are located on the hypotube 92 cm (brachial technique) and 102 cm (femoral technique) from the distal end of the delivery system.

How Supplied

Sterile. Non-pyrogenic. Device is sterilized with ethylene oxide. Device is vacuum-packed in Argon gas.

Contents

One (1) Magmaris Drug-eluting Resorbable Metal Scaffold System and one (1) Compliance Chart in a sealed, peel-open pouch.
One (1) Instructions For Use manual including one (1) Patient Implant Card.

Storage

Store in a dark, dry location at max. 25 °C / 77 °F.

Indications

Magmaris is indicated for improving luminal diameter in the treatment of de novo native coronary artery lesions in patients with symptomatic coronary artery disease. The reference vessel diameter should closely match the scaffold nominal diameter (see "Vessel Sizing" table on page 2). The treated lesion length should be less than the nominal scaffold length.

Contraindications

All contraindications for percutaneous transluminal coronary angioplasty (PTCA) are contraindications for this device.

Contraindications for this device and PTCA in general are:

Lesions judged to prevent complete inflation of an angioplasty balloon or proper placement of the scaffold or the scaffold system.

Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.

A known hypersensitivity or allergy to implant materials including Magnesium, Yttrium, Neodymium, Zirconium, Gadolinium and Dysprosium, Tantalum, Poly-L-lactide, Sirolimus or its derivatives.

Furthermore, all procedure-related contraindications as described in the national and international guidelines of the respective medical associations apply.

Warnings

- The use of this product carries the risks associated with coronary artery scaffolding, including scaffold thrombosis, vascular complications and/or bleeding events. Therefore, patients should be carefully selected (see section "Patient **Selection/Individualization of Treatment**") and maintained on clinically adequate post-procedural antiplatelet therapy (acetylsalicylic acid, clopidogrel, ticlopidine or thienopyridine). **The** long-term clinical performance and safety of resorbable drug-eluting scaffold and the long-term risk associated with these implants are unknown.
- It is NOT advised to treat patients with:
 - A lesion with excessive tortuosity proximal to or within the lesion
 - Chronic total occlusions
 - Poor flow distal to the identified lesions
 - Heavily calcified lesions
- Scaffolding across any side branches ≥ 2.0 mm in diameter is NOT advised.
- This device is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of **the** patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. BIOTRONIK will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- DO NOT use device after the "Use by" date indicated on the label.

- DO NOT use if the package is damaged or opened or if any information provided is obscured or damaged.
- DO NOT expose the scaffold to any liquid or humid surface prior to insertion in the body in order to avoid drug coating damage or premature corrosion of the scaffold and drug elution. Exception: flushing of the guide wire lumen just prior to insertion.
- DO NOT attempt to remove or readjust the scaffold on the delivery system as it may damage the scaffold, drug carrier system and/or lead to scaffold embolization. The scaffold cannot be removed and placed on another balloon catheter.
- Magmaris is NOT compatible with:
 - Luge guide wire
 - GuideLiner guide extension catheter
- Manipulate the scaffold system under angiographic guidance when it is in the body.
- Care must be taken to prevent scaffold loss during scaffold placement. To ensure precise location, verify the position of the scaffold using the balloon radiopaque markers before inflating the balloon. After removing the balloon, only the scaffold markers will be visible under fluoroscopy.
- Use only an appropriate balloon inflation medium such as 50:50 mixture of Visipaque or Hexabrix contrast media and saline. NEVER use air or any gaseous medium to inflate the balloon.
- DO NOT exceed the RBP indicated in the "Compliance Chart" table on page 2. Use of a pressure-monitoring device is mandatory to prevent over-pressurization.
- DO NOT exceed the diameter of the vessel proximal and distal to the stenotic lesion by more than 10% when inflating the balloon to avoid causing damage to the vessel.
- If the scaffold needs to be post-dilated, DO NOT exceed the maximum scaffold diameter. The maximum diameter is stated on the label and in the "Maximum Diameter for Post-dilation" table on page 2. ~~Over-dilation may adversely affect the mechanical support provided by the scaffold.~~
- ~~For dissections occurring during or immediately after implantation, a second Magmaris should be used. In case the use of a Magmaris is not possible, non-clinical tests have shown that only a proBIO-coated stent such as BIOTRONIK's ORSIRO Sirolimus-eluting stent or PRO-Kinetic Energy bare metal stent should be considered to avoid premature scaffold degradation. The proBIO coating reduces metal ion diffusion and prevents direct contact of the Magnesium scaffold with stents of dissimilar materials.~~
- Subsequent restenosis may require repeat dilation of the arterial segment containing the scaffold. The long-term outcomes following repeat dilation of endothelialized scaffolds are unknown.
- Percutaneous coronary interventions in the first 120 days may lead to early loss of mechanical support or possible embolization of scaffold segments if not yet covered by new tissue.
 - For balloon angioplasty, DO NOT exceed the maximum scaffold diameter stated on the label and in the "Maximum Diameter for Post-dilation" table on page 2.
 - For stenting, non-clinical tests have shown that only a proBIO-coated stent such as BIOTRONIK's ORSIRO Sirolimus-eluting stent or PRO-Kinetic Energy bare metal stent should be considered to avoid premature scaffold degradation. The proBIO coating reduces metal ion diffusion and prevents direct contact of the Magnesium scaffold with stents of dissimilar materials.
- Non-clinical tests have demonstrated the Magmaris is MR conditional. For the conditions for scanning, refer to the "MRI Safety Information" section.

Drug Interactions

Drugs like Tacrolimus, which act through the same cellular binding protein, may interfere with the efficacy of Sirolimus. Sirolimus is metabolized by CYP3A4. Strong inhibitors of CYP3A4 (e.g. ketoconazole) might cause increased Sirolimus exposure to levels associated with systemic effects, especially if multiple scaffolds are deployed. For a list of drug interaction please refer to the EMA information for the systemic application of Sirolimus:

http://www.ema.europa.eu/docs/en_GB/document_library/EPAR_-_Product_Information/human/000273/WC500046437.pdf

Systemic exposure of Sirolimus should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy. Additional consideration should be given to the potential local drug interactions in the vessel wall when deciding to place a Magmaris device in a patient who is taking a drug with known interactions to Sirolimus or when deciding to initiate therapy with such a drug in a patient that has recently received a Magmaris scaffold.

Precautions

General Precautions

- The patient's exposure to the drug and polymer is directly related to the number of scaffolds and the scaffold length implanted. The nominal total drug load for each size is provided in the "Sizes and Associated Total Drug Loads" table on page 2. See also "Drug Interaction" section.
- Prior to using the product, the scaffold system should be visually examined to verify functionality and ensure that its size is suitable for the specific procedure for which it is to be used.
- Only physicians thoroughly trained and experienced in the performance of PTCA and stent implantation should use this device.
- Use only with water soluble contrast media. Magmaris is compatible with contrast media currently used in clinical practice (ionic and non-ionic). Bench tests were conducted with Visipaque and Hexabrix contrast mediums.
- If the device is withdrawn prior to expansion, DO NOT re-insert as the scaffold and/or the delivery system may have been damaged during the initial attempt to cross the lesion or during withdrawal.

Handling Precautions

- Exercise care during device handling to reduce the possibility of displacement of the scaffold on the balloon and/or accidental breakage, bending, kinking of the scaffold system shaft. Special care must be taken not to handle the scaffold in any way, which may disrupt the drug carrier coating. This is important when/during catheter removal from the packaging, flushing, backloading over the guide wire and advancing it through the hemostatic valve adapter and guiding catheter hub.
- Take care when removing the scaffold system from the spiral packaging as forceful movements may dislocate the protector and the scaffold.
- When removing the protector, always pull at the very distal end of the protector to avoid dislocation of the scaffold. DO NOT touch the part of the protector over the scaffold.
- DO NOT suspend the scaffold within any solution on the sterile field. The scaffold must not be exposed to contact with liquids prior to preparation and delivery as the coating may be susceptible to damage or premature drug elution. Exception: flushing of the guide wire lumen just prior to insertion.

Scaffold Placement Precautions

- Pre-dilation of the lesion is mandatory. DO NOT exceed the reference vessel diameter when selecting a non-compliant balloon for pre-dilation. Ensure that the balloon diameter matches the vessel diameter (see "Vessel Sizing" table on page 2). The residual stenosis before implanting a Magmaris should be approximately 20%. DO NOT continue the procedure until the lumen diameter of the treated section reaches 2 mm or more.
- DO NOT apply negative pressure to the scaffold system at any time prior to placement of the scaffold across the lesion. This may cause loosening and dislodgement of the scaffold.
- Use only guide wires with a diameter of 0.014" (0.36 mm) (refer to the "Warnings" section for exceptions).
- Use guiding catheters and guide extension catheters with a minimum inner diameter of ≥ 0.070 " (1.78 mm; 6F) (refer to the "Warnings" section for exceptions).
- Avoid excessive manipulation of the scaffold during flushing of the guide wire lumen.
- When inserting and positioning the scaffold system, ensure that the hemostatic valve of the guiding catheter is fully open. Manipulation through a partially opened hemostatic valve may damage the scaffold or dislodge it from the centered location of the balloon.
- DO NOT torque the delivery system more than one (1) full turn.
- During the procedure, make sure that the guide wire exit port of the delivery system, 29 cm from the distal tip of the delivery system, remains in the guiding catheter.
- DO NOT apply excessive force while accessing or crossing the lesion. This may damage the scaffold and/or dislodge it from the balloon. If resistance is felt at any time, DO NOT force the passage. If the scaffold system is unable to reach or cross the lesion easily, the procedure should be aborted. In this case, follow the instructions listed in the "Removal of an Unexpanded Scaffold" section.
- DO NOT inflate the balloon if vacuum cannot be held as this indicates a leak in the delivery system. If a vacuum cannot be held, follow the instructions listed in the "Removal of an Unexpanded Scaffold" section.
- The balloon of the delivery system must be inflated up to the nominal pressure (NP) indicated on the label. Maintain inflation pressure for at least 15-30 seconds for full expansion of the scaffold.
- Avoid barotrauma outside the scaffold margins during post-dilation.
- Additional expansion of a deployed scaffold may cause a flow limiting dissection. This may be treated by implantation of a second Magmaris. If this is not possible, non-clinical tests have shown that only a proBIO-coated stent such as BIOTRONIK's ORSIRO Sirolimus-eluting stent or PRO-Kinetic Energy bare metal stent should be considered to avoid premature scaffold degradation. The proBIO coating reduces metal ion diffusion and prevents direct contact of the Magnesium scaffold with stents of dissimilar materials.

- The vascular effect of overlapping Magmaris scaffolds has not been tested and should **not** be employed.
- When treating multiple lesions, the distal lesion should be initially treated, followed by the proximal lesion. Scaffold placement in this order obviates the need to cross the proximal scaffold/stent in placement of the distal scaffold/stent and reduces the chances for dislodging the proximal scaffold/stent.
- The use of brachytherapy treatment, mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters to treat in-scaffold restenosis is not recommended.

Scaffold/Delivery System Removal Precautions

- Scaffold retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications can include bleeding, hematoma or pseudoaneurysm.
- If an unusual resistance is felt at any time during either at lesion access, while crossing the lesion or during the removal of the delivery system (pre- or post-implantation), the scaffold system/delivery system must be removed according to the Instructions listed in the "Removal of an Unexpanded Scaffold" section or the "Removal of the Scaffold System/Delivery System and the Guiding Catheter as a Single Unit" section. Failure to follow correct removal steps for the delivery system or for an unexpanded scaffold system and/or applying excessive force to the system can potentially result in loss or damage to the scaffold and/or delivery system components.

Post-implantation Precautions

- Care must be taken when crossing a newly deployed scaffold with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, a balloon catheter or any other device to avoid disrupting the scaffold placement, wall-apposition, coating or geometry of the scaffold.
- If the patient requires MR imaging, refer to the "MRI Safety Information" section.

Potential Adverse Events When Used As Intended

Possible adverse events associated with PTCA, stent and **scaffold placement:**

- Cardiac events: myocardial infarction or ischemia, abrupt closure of coronary artery, restenosis of treated artery, cardiogenic shock, angina, tamponade, perforation or dissection of coronary artery or aorta, cardiac perforation, emergency cardiac surgery, pericardial effusion, aneurysm formation
- Arrhythmic events: ventricular tachycardia, ventricular fibrillation, atrial fibrillation, bradycardia
- Scaffold system events: failure to deliver scaffold to intended site, scaffold dislodgement from the delivery system, scaffold misplacement, scaffold deformation, scaffold embolization, scaffold thrombosis or occlusion, scaffold fracture, scaffold migration, inadequate apposition or compression of scaffold/s, inflation difficulties, rupture or pinhole of the delivery system balloon, deflation difficulties, withdrawal difficulties, embolization of catheter material
- Respiratory events: acute pulmonary edema, congestive heart failure, respiratory insufficiency or failure
- Vascular events: access site hematoma, hypotension/ hypertension, pseudoaneurysm, arteriovenous fistula formation, retroperitoneal hematoma, vessel dissection or perforation, restenosis, thrombosis or occlusion, compromise of side branch patency, vasospasm, peripheral ischemia, dissection, distal embolization (air, tissue debris, thrombus)
- Neurologic events: permanent (stroke) or reversible (TIA) neurologic event, femoral nerve injury, peripheral nerve injury.
- Bleeding events: access site bleeding or hemorrhage, hemorrhage requiring transfusion or other treatment
- Allergic reactions to contrast media, antiplatelets, anticoagulants, to the drug carrier (PLLA), to Sirolimus and/or Magnesium, Yttrium, Neodymium, Zirconium, Gadolinium and Dysprosium
- Infection and sepsis
- Death

Potential adverse events related to the oral administration of Sirolimus include abnormal liver function tests, anemia, arthralgia, diarrhea, hypercholesterolemia, hypersensitivity (including anaphylactic/anaphylactoid type reactions), hypertriglyceridemia, hypokalemia, infections, interstitial lung disease, leukopenia, lymphoma and other malignancies, as well as thrombocytopenia.

Patient Selection/Individualization of Treatment

Judicious selection of patients according to the indications and contraindications is necessary since the use of this device carries the associated risks of complications listed in the "Potential Adverse Events When Used As Intended" section. The risks and benefit should be considered for each patient before using the device.

Patient selection factors to be assessed should include a judgment regarding risk of long-term antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease and patients with high risk of bleeding in which anticoagulation therapy would be contraindicated.

Antiplatelet/ anticoagulant medication may be used in combination with Magmaris. Physicians should use the information from the current Sirolimus-eluting stent/scaffold literature and national and international guidelines on PCI concerning the selection, dosage, duration and combination of different antithrombotic drugs. Specific needs and the risk profile of individual patients may influence the antiplatelet/anticoagulation regimen to be used.

It is very important that the patient is compliant with the post-procedure antiplatelet/anticoagulant medication recommendation. Premature discontinuation of prescribed antiplatelet/anticoagulant medication could result in a higher risk of thrombosis, myocardial infarction or death. Prior to PCI, if a surgical or dental procedure is anticipated and requires early discontinuation of antiplatelet/anticoagulation therapy, the interventionist and patient should carefully consider whether a drug-eluting scaffold and its associated recommended antiplatelet/anticoagulation therapy is the appropriate PCI choice. Following

PCI, should a surgical or dental procedure be recommended, the risk and benefits of the procedure should be weighed against the possible risk associated with premature **discontinuation of antiplatelet therapy.**

Patients who require premature discontinuation of the antiplatelet/ **anticoagulation therapy** due to significant bleeding, should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet/anticoagulation therapy restarted as soon as **possible at** the discretion of **their** treating physicians. Premorbid conditions that increase the risk of poor initial result or the risk of emergency referral for bypass surgery (diabetes, renal failure and severe obesity) should be reviewed. A review of the vessel location, **reference** vessel size, lesion length and qualitative target lesion characteristics must also be **considered.**

Thrombosis following scaffold implantation is affected by several baseline **angiographic and** procedural factors. These include vessel diameter less than 3.0 mm, intra-procedural thrombus and dissection following scaffold or stent implantation. In patients who have undergone coronary stenting, the persistence of thrombus or **dissection should be** considered a marker for subsequent thrombotic occlusion. Following PCI, the patients should be monitored very carefully during the first month after scaffold implantation.

Use in Special Population

The risks and benefit in special populations should be considered carefully for each patient before use of the device. The safety and effectiveness of the Magmaris has not yet been established in the following patient population:

- Patients with unresolved thrombus at the lesion site.
- Patients with ostial lesions, tandem lesions or lesions located **in the left main coronary artery or bifurcations, in saphenous vein grafts, in previously placed stents or in-stent restenosis.**
- Patients with lesions that have been treated with mechanical atherectomy devices or **laser angioplasty catheters.**
- Patients with tortuous vessels that may impair scaffold placement in the region of the obstruction or proximal to the lesion.
- Patients with diffuse disease, or poor outflow distal to the identified lesion.
- Patients with STEMI and NSTEMI within the previous 72 hours.
- Patients with impaired renal function.
- Patients with brachytherapy treatment of the target lesion.
- Pregnant patients: there are no adequate and well-controlled studies in pregnant women or men intending to father children. Effective contraception should be **initiated** before implanting a Magmaris and for 12 weeks **after implantation.** The **Magmaris** should be used during pregnancy only if the benefit outweighs the risk to **the embryo or fetus.**
- Pediatric patients below the age of 18 years.

Directions for Use

Patient Preparation and Scaffold System Selection

1. Prepare the patient and vascular access site for a PCI procedure according to the institution's standard clinical practice.

Caution: Pre-dilation of the lesion is mandatory. It is advised to use a non-compliant balloon. DO NOT exceed the reference vessel diameter. Ensure that the balloon diameter matches the vessel diameter. The residual stenosis before implanting a Magmaris should be approximately 20%. DO NOT continue the procedure until the lumen diameter of the treated section reaches 2 mm or more.

2. Select the scaffold size to match the diameter of the vessel to achieve a final scaffold diameter to vessel ratio of 1:1 and a full coverage of the lesion over its entire length. Use QCA, IVUS and/or OCT for quantitative lesion evaluation if uncertain about the vessel diameter. Available scaffold diameters are 3.0 mm and 3.5 mm, therefore DO NOT implant into vessel diameters <2.7 mm or >3.7 mm or lengths longer than the available sizes (15, 20 and 25 mm).

Note: It is important to select a scaffold to closely match the target vessel diameter. Refer to the "Vessel Sizing" table on page 2. Quantitative imaging is recommended for the assessment of target vessel diameter at baseline for appropriate scaffold size selection. The scaffold should not be post-dilated more than the maximum diameter stated on the label and in the "Maximum Diameter for Post-dilation" table on page 2. Over-dilation may adversely affect the mechanical support provided by the scaffold.

Note: The treated lesion length should be less than the nominal scaffold length.

Scaffold System Preparation

Warning: DO NOT expose the scaffold to any liquid or humid surface prior to insertion in the body in order to avoid drug coating damage or premature corrosion of the scaffold and drug elution. Exception: flushing of the guide wire lumen just prior to insertion.

Caution: Exercise care during device handling to reduce the possibility of displacement of the scaffold on the balloon and/or accidental breakage, bending, kinking of the scaffold system shaft. Special care must be taken not to handle the scaffold in any way, which may disrupt the drug carrier coating. This is important during catheter removal from the packaging, flushing, backloading over the guide wire and advancing it through the hemostatic valve adapter and guiding catheter hub.

3. Before opening, carefully inspect outer and inner package. DO NOT use if there is any damage to the packaging and/or any information provided upon is obscured.
4. Remove the protection ring containing the scaffold system from the sterile package and place it onto a sterile field.
5. Gently pull out the scaffold system from the protection ring.

Caution: Take care when removing the scaffold system from the protection ring, as forceful movements may dislocate the protector and the scaffold.

6. Carefully remove the scaffold protector.

Caution: When removing the scaffold protector, always pull at the very distal end of the protector to avoid dislocation of the scaffold.

7. Visually check the scaffold crimping for uniformity, no protruding struts, and verify that the scaffold is centered between the proximal and distal balloon markers. DO NOT use if any defects are noted.

Caution: DO NOT apply negative pressure to the scaffold system at any time prior to placement of the scaffold across the lesion. This may cause loosening and dislodgement of the scaffold.

Flushing of the Guide Wire Lumen

8. Connect a syringe containing heparinized normal saline to an appropriately sized flushing needle. Carefully apply the needle to the distal tip of the delivery system and flush the guide wire lumen until fluid exits the guide wire exit port.
9. Remove the syringe and the flushing needle.
10. Leave the prepared scaffold system at ambient pressure.

Insertion and Scaffold Positioning

11. Attach a hemostatic valve to the luer port of the guiding catheter positioned within the vasculature.
12. Position the guide wire under fluoroscopy in accordance with PCI techniques.
13. Backload the proximal end of the guide wire into the distal tip of the delivery system until it exits at the guide wire exit port.
14. Open the hemostatic valve completely.

Caution: When inserting and positioning the scaffold system, ensure that the hemostatic valve of the guiding catheter is fully open. Manipulation through a partially opened hemostatic valve may damage the scaffold or dislodge it from the centered location of the balloon.

15. Carefully insert the scaffold system through the hemostatic valve.
16. Advance the scaffold system through the guiding catheter using fluoroscopic guidance to determine when the delivery system tip approaches the distal tip of the guiding catheter.

Note: The shaft exit markers on the hypotube may be used to approximate when the scaffold system has reached the distal end of the guiding catheter.

17. Carefully advance the scaffold system into the coronary artery over the guide wire while maintaining stable guiding catheter seating and stable guide wire placement across the target lesion.
18. Position the scaffold within the lesion using the balloon radiopaque markers as reference points.

Caution: DO NOT apply excessive force while accessing or crossing the lesion. This may damage the scaffold and/or dislodge it from the balloon. If resistance is felt at any time, DO NOT force the passage. If the scaffold system is unable to reach or cross the lesion easily, the procedure should be aborted. In this case, follow the instructions listed in the "Removal of an Unexpanded Scaffold" section.

19. Verify the scaffold position under angiographic guidance to ensure that the total lesion (proximal to distal margins) has been adequately covered. It is recommended to use an enhanced visualization tool/software to identify the scaffold markers. If the position of the scaffold is not optimal, it should be carefully repositioned or removed (refer to the "Removal of an Unexpanded Scaffold" section). Expansion of the scaffold should not be undertaken if the scaffold is not properly positioned within the target lesion.
20. Sufficiently tighten the hemostatic valve.

Removal of Air from the Delivery System

21. Connect a three-way stopcock to the luer lock of the delivery system.
 22. Prepare and remove air from a 20 ml capacity inflation/deflation device according to manufacturer's recommendations and instructions.
 23. Attach the inflation/deflation device containing 3 ml of balloon Inflation medium to the **stopcock**.
- Warning:** Use only an appropriate balloon inflation medium. In vitro tests were performed with a 50:50 mixture of 50% saline with either Visipaque or Hexabrix. NEVER use air or **any gaseous medium to inflate the balloon**.
24. Open the stopcock **so that an open** fluid path between the catheter and the **Inflation/deflation device is established**.
 25. Pull the plunger of the **inflation/deflation** device and aspirate air from the catheter for at least 35 seconds. Release the plunger to neutral for contrast fill.
- Caution:** DO NOT inflate the balloon if a vacuum cannot be held, as this **Indicates a leak in** the delivery system. If a vacuum cannot be held, follow the Instructions for the **"Removal of an Unexpanded Scaffold"**.
26. Close the stopcock **so that** the fluid path to the catheter is closed and remove all air from the inflation/deflation device through the stopcock.
 27. Repeat steps 24-26 if necessary, to ensure air contained in the balloon and Inflation lumen are removed. Release the inflation/deflation barrel to normal pressure.
 28. Open the stopcock. Set the system on neutral aside for use. The scaffold is now ready to be deployed.

Scaffold Deployment

29. Prior to scaffold expansion, reconfirm the correct position of the scaffold relative to the target lesion. Verify, using angiographic guidance, that the scaffold **has** been positioned correctly relative to the target lesion and that there has been no vessel damage. Inflate the dilatation balloon slowly to expand the scaffold to the diameter **in accordance with** the "Compliance Chart" table on page 2. Maintain inflation pressure for at least 15-30 seconds for full expansion of the scaffold.
- Note:** Diameter values listed in the "Compliance Chart" table on page 2 were assessed in non-clinical tests with no lesion-resistance during inflation. Clinical **conditions may differ**.
- Caution:** Inflate to at least the nominal pressure (NP) indicated on the label and in the "Compliance Chart" table on page 2. DO NOT exceed the rated burst pressure (**RBP**).
30. Use multiple fluoroscopic views to ensure that the balloon is completely expanded along **its length**.
 31. If necessary, the delivery system balloon may be dilated once more in order to achieve optimum seating of the implanted scaffold.
 32. Post-dilate to ensure that the flow has been restored sufficiently and that the scaffold struts are well-apposed and embedded in the vessel wall. Use a non-compliant balloon at high pressure (> 16 atm) at the same nominal size as the used scaffold or up to 0.5 mm larger.
- Caution:** Avoid barotrauma outside the scaffold margins during post-dilation.
- Warning:** DO NOT exceed the maximum scaffold diameter stated on the label and **in the** "Maximum Diameter for Post-dilation" table on page 2. Over-dilation may adversely affect the mechanical support provided by the scaffold.
- Caution:** Care must be taken when crossing a newly deployed scaffold **with an** intravascular ultrasound (IVUS) catheter, a coronary guide wire, a balloon catheter or any **other device to** avoid disrupting the scaffold placement, wall-apposition, coating or geometry of the Magmaris.

Balloon Deflation and Delivery System Removal

33. Deflate the balloon in accordance with standard PCI procedures.
34. Apply negative pressure to the balloon for at least 35 seconds. Deflation of the balloon should be confirmed by absence of contrast within in the balloon.
35. Open the hemostatic valve fully to allow removal of the delivery system.
36. Maintain the position of the guiding catheter and guide wire to prevent them from being drawn into the vessel. Under fluoroscopic control, and maintaining negative pressure on the inflation device, carefully pull the completely deflated delivery system out of the target vessel into the guiding catheter.

Note: If the delivery system cannot be easily withdrawn from the scaffold, slightly advance and retract the delivery system carefully. If resistance is still felt, repeat the operation until it is possible to gently pull the balloon out of the scaffold.

Caution: If resistance is felt whilst pulling the delivery system into the guiding catheter, remove the delivery system and the guiding catheter as a single unit (refer to the "Removal of an Unexpanded Scaffold" section).

37. After removing the delivery system, tighten the hemostatic valve.
38. Inspect the device immediately upon removal from the patient for any signs of breakage or fragmentation.
39. Observation of the patient and angiographic evaluation should be performed periodically in the 15 minutes after implanting the device.
40. After use, dispose the product and packaging in accordance with hospital, administrative and/or local government policy.

Removal of an Unexpanded Scaffold

1. If removal of scaffold system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the delivery system and avoid any acute angle between the floppy part of the delivery system and the guiding catheter.
2. Slowly pull back the scaffold system into the guiding catheter. The entry of the scaffold into the guiding catheter must be performed slowly under fluoroscopic control to avoid dislodgement of the scaffold from its position on the delivery system balloon.

Caution: If resistance is felt whilst pulling the scaffold system into the guiding catheter, remove the scaffold system and the guiding catheter as a single unit (proceed as directed).

3. The lesion must be pre-dilated again or otherwise prepared before a second attempt at scaffolding is undertaken.

Caution: If the product has been withdrawn prior to expansion, DO NOT re-insert it as the scaffold and/or the delivery system may have been damaged during the initial attempt to cross the lesion or during withdrawal.

4. Dispose of the used product in accordance with hospital, administrative and/or local government policy.

Removal of the Scaffold System/Delivery System and the Guiding Catheter as a Single Unit

1. DO NOT retract the scaffold system/delivery system into the guiding catheter. Maintain guide wire placement across the lesion and carefully pull back the scaffold system/delivery system. Position the proximal balloon marker just distal to the tip of the guiding catheter.
2. Advance the guide wire into the artery as distally as safely possible.
3. Tighten the rotating hemostatic valve to secure the scaffold system/delivery system to the guiding catheter.
4. Remove the guiding catheter and the scaffold system/delivery system as a single unit.
5. Dispose of the product in accordance with hospital, administrative and/or local government policy.

MRI Safety Information

Non-clinical tests have demonstrated that the scaffold is MR conditional for:

- Single and overlapped configurations up to total implanted length of 47 mm.
- Overlapped configurations of Magmaris and BIOTRONIK ORSIRO Sirolimus-Eluting Stents or PRO-Kinetic Energy bare metal stents, up to total stented length of 62 mm.

A patient with this scaffold can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T only.
- Maximum spatial gradient magnetic field of 3,000 gauss/cm (30 T/m) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) for 15 minutes of scanning.

Under the scan conditions defined, the scaffold is expected to produce a maximum temperature rise of less than 6.2° C after 15 minutes of continuous scanning.

In non-clinical tests, the image artifact caused by the device extended approximately 4 mm from the scaffold when imaged with a gradient echo pulse sequence or 3 mm when imaged with a spin echo pulse sequence and a 3.0 T MRI system. The artifact may obscure the device lumen.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters. Heating in the presence of local transmit MR coils was not evaluated. It is recommended that local transmit coils not be placed directly over the scaffolds.

Warranty/Liability

The product and each component of its system (hereinafter the product) have been designed, manufactured, tested and packaged with all reasonable care. However, BIOTRONIK has no control over the conditions under which the product is used and a disturbance of the intended function of the product may occur for various reasons. In this respect, the warnings in this product publication/instructions for use are expressly to be considered as an integral part of this disclaimer and provide more detailed information. For this reason, BIOTRONIK disclaims all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability or fitness for a particular purpose of the product. Product descriptions or user guidelines in publications do not constitute any expressed representation or any expressed or implied warranty. BIOTRONIK is not liable for any direct, incidental or consequential damages or medical expenses caused by any use, defect, failure or malfunction of the product whether the claim is based on contract, warranty, tort or otherwise. This does not apply in the case of intention or in the case of gross negligence of legal representatives or executive staff of BIOTRONIK. In commercial transactions relating to merchants, the liability is limited to the compensation of typical damages; compensation for any untypical or incidental damage is excluded. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable in the respective country. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects BIOTRONIK's legitimate interest in limiting its liability or warranty without infringing any mandatory provisions of applicable law. No person has any authority to bind BIOTRONIK to any warranty or liability regarding the product.

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**Addition to the Instruction for Use
Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System**

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The content of the English and Russian document versions is identical.

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Susanne Minich


Official Certification

Seen for authentication of the reverse side signature, affixed in our presence by

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NOTARIAT BÜLACH

Matthias Eisenhut

Matthias Eisenhut, Deputy Notary Public

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S. Overkott

S. Overkott

Addition to the Instruction for Use

**Sirolimus-Eluting Resorbable Coronary Magnesium
Scaffold System**

1. Medical Device Name	3
2. Name and address of the manufacturer, production Site Contact	3
3. Intended Use.....	3
4. Usage Terms, potential user.....	3
5. Scope of use	3
6. Types of contact with human body	3
7. Risks of Medical Device usage	4
8. Technical Description of Medical Device	4
9. Materials, that the medical device is made of	8
10. Contents and Scope of delivery.....	8
11. Criteria of unsuitability of a medical device for use.....	9
12. Packaging Description	9
13. Labelling	9
14. Information on application frequency, maintenance and repairs.....	11
15. Deimplantation and / or replacement	11
16. The list of consumables, the procedure for their use and replacement	11
17. Cleaning, disinfection and sterilization	11
18. Frequency of use	11
19. Shelf life	11
20. Storage and transportation condition.....	11
21. Utilization and disposal	12
22. Environmental protection requirements	12
23. The authorized representative on the territory of Russian Federation.	12
24. List of Applied Standards	13

1. Medical Device Name

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System,
Scope of Delivery:

1. Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold on delivery system, of following sizes: 3,0 mm x 15 mm; 3,0 mm x 20 mm; 3,0 mm x 25 mm; 3,5 mm x 15 mm; 3,5 mm x 20 mm; 3,5 mm x 25 mm - 1 piece.
2. Compliance Chart – 1 piece.
3. Patient Implant Card – 1 piece.
4. Instruction for Use – 1 piece.
5. Addition to Instruction for Use – 1 piece.

hereinafter: «Magmaris», «Stent», «Scaffold».

2. Name and address of the manufacturer, production Site Contact

Name and address of the manufacturer:

BIOTRONIK AG
Ackerstrasse 6, 8180 Bülach, Switzerland
Тел. + 41 (0) 44 864 51 11
Тел/Факс + 41 (0) 44 864 53 53
info.vi@biotronik.com

Production Site Contact:

BIOTRONIK AG
Ackerstrasse 6, 8180 Bülach, Switzerland

CORTRONIK GmbH
Friedrich-Barnewitz-Strasse 4A, Rostock-Warnemünde 18119, Germany

3. Intended Use

The Medical Device is intended to improve coronary luminal diameter in patients with symptomatic ischemic heart disease due to stenotic coronary artery lesions.

4. Usage Terms, potential user

This medical device must be used only by specially trained personnel in specialized medical facilities, when performing surgical operations in stationary conditions.

5. Scope of use

Cardiovascular surgery: during percutaneous transluminal coronary angioplasty / percutaneous coronary intervention.

6. Types of contact with human body

Magmaris Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System has the the following types of contact with human body:

Stent Magmaris	Implantable device with permanent contact (more than 30 days) with internal tissues of the human body
Delivery system of stent Magmaris	A part of the product that has a short-term (less than 24 hours) contact with internal tissues of the human body

7. Risks of Medical Device usage

Risk management methods are applied with regard to the product design, the materials/components and manufacturing processes and the application of the device. Risk management is conducted with reference to international standard EN ISO 14971 Medical devices, application of risk management to medical devices.

Biotronik AG Top Management issues the risk policy which applies to all medical devices manufactured by Biotronik AG. This document defines the criteria for risk acceptability.

Risk analyses are performed to identify and evaluate potential risks that may arise on account of design, manufacturing processes (process FMEA, KKP), and the application of the product. The design FMEA and the application FMEA take account of the results of the usability analysis (HUF). Risk control measures are applied to reduce risks to acceptable levels. The acceptability of residual risks remaining after the implementation of risk control measures in each individual risk analysis is evaluated and documented.

8. Technical Description of Medical Device

Magmaris is a drug-eluting balloon-expandable resorbable scaffold mounted on a fast-exchange 0.014" guidewire compatible delivery system.

It is a drug/device combination product comprising:

A device component: resorbable scaffold pre-mounted on a fast-exchange delivery system with a noncompliant balloon.

A drug component: Drug polymer coating on the scaffold which is a formulation of the active pharmaceutical ingredient (API) Sirolimus in a bioresorbable Poly-L-Lactide (PLLA) polymer coating.

Absorbable Metal Scaffold

Magmaris is a balloon expandable scaffold of magnesium alloy. The magnesium alloy contains >90% magnesium, further additional alloying elements such as Yt-trium, Neodymium, Zirconium, Gadolinium and Dysprosium.

The scaffold design consists of radial revolving loops with connectors in the axial direction, strut thickness is 150 µm, and radiopaque markers are at each end.

The whole scaffold surface is completely (luminal and abluminal) coated with PLLA which contains the active substance Sirolimus.

Table 1: Product Range of Magmaris.

Nº pos.	Designation	number (REF)	nominal Ø at NP (mm)	nominal length (mm)	Usable length of DS (cm)
1	Magmaris 3,0 mm x 15 mm	412526	3,0±0.27	15 ±0,03	140±1
2	Magmaris 3,0 mm x 20 mm	412527	3,0±0.27	20 ±0,03	140±1
3	Magmaris 3,0 mm x 25 mm	412528	3,0±0.27	25 ±0,03	140±1
4	Magmaris 3,5 mm x 15 mm	412529	3,5±0.31	15 ±0,03	140±1
5	Magmaris 3,5 mm x 20 mm	412530	3,5±0.31	20 ±0,03	140±1
6	Magmaris 3,5 mm x 25 mm	412531	3,5±0.31	25 ±0,03	140±1

Table 2: Stent Technical characteristics

Magmaris	Wall thickness (mm)	Outer Ø of scaffold (mm)	inner Ø of inflated scaffold at RBP (mm)	Scaffold length (mm)	Single cell size (edge length of the area, mm)	Total Drug Load (µg)
Magmaris 3,0 mm x 15 mm	0,15±0.05	2.0±0.05	3.29 ± 0.29	15 ±0,03	19.24±0.015	166
Magmaris 3,0 mm x 20 mm	0,15±0.05	2.0±0.05	3.29±0.29	20 ±0,03	19.24±0.015	226
Magmaris 3,0 mm x 25 mm	0,15±0.05	2.0±0.05	3.29±0.29	25 ±0,03	19.24±0.015	287
Magmaris 3,5 mm x 15 mm	0,15±0.05	2.0±0.05	3.82±0.33	15 ±0,03	20.59±0.015	159
Magmaris 3,5 mm x 20 mm	0,15±0.05	2.0±0.05	3.82±0.33	20 ±0,03	20.59±0.015	223
Magmaris 3,5 mm x 25 mm	0,15±0.05	2.0±0.05	3.82±0.33	25 ±0,03	20.59±0.015	271

Scaffolding time ~ 12 weeks

Absorption period of Mg scaffold ~ 12 months

Absorption coating period ~ 36 months

Table 3: Stent Technical characteristics

Magmaris	The force required to inflate scaffold no more than.. (kPA)	The force required to detach scaffold from the unexpanded balloon, not less than... (N)	The force required to deform scaffold, no less than (N/mm)	Min stent bend radius (mm)
Magmaris 3,0 mm x 15 mm	1621	1.5	1.26±0.07	6.00
Magmaris 3,0 mm x 20 mm	1621	1.5	1.26±0.07	6.00
Magmaris 3,0 mm x 25 mm	1621	1.5	1.26±0.07	6.00
Magmaris 3,5 mm x 15 mm	1621	1.5	0.91±0.07	7.25
Magmaris 3,5 mm x 20 mm	1621	1.5	0.91±0.07	7.25
Magmaris 3,5 mm x 25 mm	1621	1.5	0.91±0.07	7.25

Scaffold Coating

The drug-polymer coating consists of the antiproliferative drug Sirolimus (active pharmaceutical ingredient) and a high molecular weight Poly-L-Lactide (PLLA) polymer carrier. The entire scaffold surface is covered by the drugpolymer coating. The composition of the drug polymer coating is shown in Table 4. The drug-polymer coating is specified in Table 5.

Table 4: Drug-polymer coating constituents and quality standards

Function	Materials	CAS Number	Unit amount per Drug Eluting Scaffold
Polymer	PLLA	33135-50-1	$6.6 \pm 0.7 \mu\text{g}/\text{mm}^2$ scaffold surface
Active ingredient	Sirolimus	53123-88-9	$1.4 \pm 0.2 \mu\text{g}/\text{mm}^2$ scaffold surface

Table 5: Drug-polymer coating specifications

[μg]	Acceptance limits
Coating composition	17.5 wt% Sirolimus / 82.5 wt% Poly-L-Lactide L210
Drug load	$1.4 \pm 0.2 \mu\text{g}/\text{mm}^2$ scaffold surface
Layer thickness	$\leq 25 \mu\text{m}$

Drug substance

The Sirolimus drug substance that is incorporated into the scaffold coating is approved in Europe for the prophylaxis of organ rejection in adult patients receiving renal transplants (marketing authorization number EU/1/01/171/007-8) and as an active pharmaceutical ingredient for drug eluting stents. The compound has the molecular formula $\text{C}_{51}\text{H}_{79}\text{NO}_{13}$ and a molecular weight of 914.2 g/mol. Sirolimus blocks the cell cycle progression from phase G1 to S, hence prevents cell proliferation and the proliferation and migration of smooth muscle cells, thereby inhibiting the restenosis process.

Table 6: Product range

	Catalogue Number (Ref. No)	Magmaris scaffold size	Total Drug Load [μg]	Coating mass [μg]
1	412526	Magmaris 3,0 mm x 15 mm	166	946
2	412527	Magmaris 3,0 mm x 20 mm	226	1292
3	412528	Magmaris 3,0 mm x 25 mm	287	1638
4	412529	Magmaris 3,5 mm x 15 mm	159	910
5	412530	Magmaris 3,5 mm x 20 mm	223	1274
6	412531	Magmaris 3,5 mm x 25 mm	271	1546

Drug carrier

The carrier used for the supply and release of Sirolimus, is Poly-L-Lactide (PLLA).

Delivery System

Magmaris utilizes a fast-exchange delivery system. The scaffold is positioned on the balloon between the distal and proximal balloon markers. The distal segment of the delivery system is a dual coaxial shaft composed of the outer inflation/deflation lumen, and the inner guidewire lumen. The guidewire lumen starts at the distal tip and ends at the guidewire exit port.

The proximal segment of the delivery system is composed of a single female luer port (hub) and single lumen stainless steel hypotube that extends to the guidewire exit port. An inflation/deflation device can be connected to the proximal end of the hub to facilitate the inflation/deflation of the balloon. Shaft markers are located on the hypotube at two locations to indicate brachial and femoral access.

The delivery system has a working length of 140 cm. The scaffold system is compatible with most guidewires of maximum 0.014" diameter and with 6F guiding catheters.

Table 7: Delivery System sizes

Usable length (cm)	Internal Ø of delivery systems (mm)	Outer Ø of delivery system (F / mm)	Breaking force, not less than... (N)	Allowable bending radius of the delivery system (mm)	Max. profile/ Ø of delivery system (mm)	Compatible guide wire (mm / ")	Comp. Guide catheter (F/mm)
140 ±1	0.41±0.05	Proximal shaft-2.0 F (0.67mm ±0.02) Distal shaft-2.9 F (0.97mm±0.02)	Tip: 3.0 System: 5.0	Distal outer shaft: 4mm Proximal outer shaft: 10mm	1,5	0,36 mm / 0,014"	6/ 1,78 Inner Ø ≥0,070 "

Product group	The length of the distal shaft [mm]
All sizes	290 mm ±1

Hermeticity of delivery system

The delivery system must not allow leakage of fluid at pressures from 1 up to 16 kPa for 30 sec.

Luer port of delivery system

Type – Luer

The diameter of the small base of the outer cone 3.925 - 3.990 mm

The diameter of the large base of the inner cone is 4,270 - 4,315 mm

Outer cone length - not less than 7,500 mm

The length of the inner cone is not less than 7,500 mm

Taper - 6%

Table 8: Balloon Technical characteristics

Product group (nominal Ø [mm])	Balloon length (mm)	Balloon diameter at NP (mm)	Balloon diameter at RBP (mm)	Inflation time Not more than... (sec)	Deflation time, Not more than... (sec)	Rated balloon wearout (maximum number of inflation-deflation cycles)
Magmaris Ø 3,0 mm, all length	3.0 x 15: 16.25 +0.5mm/-0.25mm	3,00 ±0.24	3,29 ±0.26	15	35	16 atm / 1621 kPa (10 cycles)
	3.0 x 20: 21.75 +0.5mm/-0.25mm			15	35	

	3.0 x 25: 26.75 +0.5mm/- 0.25mm			15	35	
Magmaris Ø 3,5 mm, all length	3.5 x 15: 16.25+0.5 mm / - 0.25mm			15	35	
	3.5 x 20: 21.75+0.5 mm / - 0.25mm	3,54 ±0.28	3,82- 0.3/+0.3 1	15	35	
	3.5 x 25: 26.75+0.5 mm / - 0.25mm			15	35	

9. Materials, that the medical device is made of

Component	Material
Scaffold Magmaris	
Scaffold body	Mg Alloy
Scaffold Marker	Tantalum powder
	Silicone
Drug Coating	Absorbable Poly (L-lactide), Sirolimus, Chloroform
Magmaris Delivery System	
Hub / Kink Protector	Polyamide 12
Hypotube	304 Stainless Steel
Hypotube Shaft Marker	White Printing Ink
Hypotube Tubing	Polyamide 12, Black Pigment 7
Distal tubing	Polyamide 12
Inner Tubing Inner Layer	High density polyethylene
Inner Tubing Middle Layer	Maleic anhydride ethylene copolymer)
Inner Tubing Outer Layer	Polyether Block Amide, Color additive Black Opaque
Balloon	Polyamide 12
Tip	Polyamide 12, Color additive Black
X-Ray Marker	Platinum, Iridium
Hydrophobic Coating (Hypotube)	Lubricant Medical Grade, n-heptane fraction
Hydrophilic Coating (distal shaft)	Base coat Top Coat

10. Contents and Scope of delivery

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System is supplied with:

1. Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold on delivery system, of following sizes: 3,0 mm x 15 mm; 3,0 mm x 20 mm; 3,0 mm x 25 mm; 3,5 mm x 15 mm; 3,5 mm x 20 mm; 3,5 mm x 25 mm - 1 piece.
2. Compliance Chart - 1 piece.
3. Patient Implant Card - 1 piece.
4. Instruction for Use - 1 piece.
5. Addition to Instruction for Use - 1 piece.

11. Criteria of unsuitability of a medical device for use

DO NOT use if the package is damaged or opened or if any information provided is obscured or damaged

12. Packaging Description

Individual / sales packaging

Magmaris is supplied in sterile packaging.

The most distal section of the device (including the Sirolimus-eluting scaffold and folded balloon) is covered by a protective tube (from polyether block amide PEBA) and a protector (from steel 1.4301).

The inner packaging (sizes 240 (± 3) x 240 (± 3) mm) comprises a protective HDPE packaging ring (dispenser) within a pouch of aluminium foil. A Compliance Data Card is supplied inside the foil pouch. After sterilization, air is evacuated out of the pouch, and the pouch is filled with an inert gas (argon) and completely sealed in its final sterile inner packaging (HDPE Tyvek).

The outer packaging (sizes 240 (± 1) x 240 (± 1) x 15 ($\pm 0,5$) mm) comprises one cardboard box (sales package) in which is packed the sterilized inner product packaging and the IFU and the Patient Implant Card (PIC) as one booklet. Product labels are applied on the inner

Transport packaging

Products are packed in the transport packaging. This is a box made of corrugated cardboard. The cardboard box is closed and fixed with adhesive tape.

It is allowed to pack several products in accordance with the overall dimensions of the transport packaging, including an air bubble film, provided that all sets are completely wrapped in a bubble film, which is also placed between cardboard boxes and each side of the transport packaging.

13. Labelling

The following symbols are used on the labeling of the packaging / box to inform about the required information:

- device name,
- Use by date,
- information on sterility and sterilization method,
- information on single usage,
- information on manufacturer,
- other symbols (see the table with Labeling legend).

When delivered to the territory of the Russian Federation for the purpose of circulation, next to the original labeling but without overlapping, the medical device is stickered with a local sticker, which has priority for the Russian Federation.

Legend for symbols used in labeling

Symbols	Description
	Manufacturer
	Catalogue number
	Batch code
	Nominal pressure

Symbols	Description
RBP	Rated burst pressure
STERILE EO	Sterilized using ethylene oxide
	Consult instructions for use
	Caution
	Date of manufacture
	Use by
	Do not reuse
	Upper Temperature limitation
	Do not resterilize
	MR conditional
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	CE Mark 93/42/EEC
	Patient Implant Card Inside
	Mark of conformity to the GOST R system
Nominal Diameter (mm)	Nominal Diameter (mm)
Scaffold Length (mm)	Scaffold Length (mm)
Guide wire diameter max (inch)	Guide wire diameter max (inch)
Guiding Catheter Inner Diameter min (inch)	Guiding Catheter Inner Diameter min (inch)
Maximum Diameter for Post-dilation (Ø max) (mm)	Maximum Diameter for Post-dilation (Ø max) (mm)
Compliance data	Compliance data

Symbols	Description
Pressure	Pressure
Scaffold Inner Diameter (mm)	Scaffold Inner Diameter (mm)
F	French

14. Information on application frequency, maintenance and repairs

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System is a single use device and is not a subject for maintenance and repairs.

15. Deimplantation and / or replacement

Stent Magmaris can not be removed and / or replaced after its implantation.

16. The list of consumables, the procedure for their use and replacement

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System is a single use device and does not contain any consumables.

17. Cleaning, disinfection and sterilization

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System is supplied sterile. Device is sterilized with ethylene oxide.

Resterilization is prohibited.

Cleaning, disinfection and sterilization are not applicable as this device intended for single use only.

18. Frequency of use

Sirolimus-Eluting Resorbable Coronary Magnesium Scaffold System is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

19. Shelf life

Shelf life of the device is 12 (Twelve) months.

20. Storage and transportation condition

Transportation of this medical device can be performed by any type of covered transport in accordance with the rules of transportation of goods. Use a medical device with caution. Keep in original box, in dry, dark, cool place.

During the transportation, storage and handling the following terms should be fulfilled:

Term	Transportation	Storage	Handling
Temperature	Up to $23\pm 1^{\circ}\text{C}$ ($73.4\pm 2^{\circ}\text{F}$)	Up to $25\pm 1^{\circ}\text{C}$ ($77\pm 2^{\circ}\text{F}$)	From 32°C (90°F) to 42°C (108°F)
Humidity	There are no special requirements for air humidity	There are no special requirements for air humidity	Resistant to biological fluid (blood)
Atmospheric pressure	There are no special requirements for atmospheric pressure	There are no special requirements for atmospheric pressure	There are no special requirements for atmospheric pressure

21. Utilization and disposal

After use, the product is a potential biological hazard. The treatment and its disposal must be carried out in accordance with the standards adopted in medical practice, as well as relevant local, state and federal legislative acts, and Instructions **(for the Russian Federation: the disposal of a medical device is carried out in accordance with class B SanPin 2.1.7.2790-10).**

22. Environmental protection requirements

Please accord to the IFU requirements to exclude the negative effect on the human or environment during the handling, shipment and storage of medical device

23. The authorized representative on the territory of Russian Federation

Limited Liability Company "Biotronik" (LLC «Biotronik»),
Nikoloyamskaya street, house 26, building 1A, floor 2, room 6, 109240 Moscow, Russian Federation
Tel.: +7 (495) 789-68-31
Fax: +7 (495) 789-68-31
office@biotronik.ru

24. List of Applied Standards

Standards

Applied Underlying Standard	Title	Compliant with European Standard
EN 556-1:2001 + AC:2006	Sterilization of medical devices. Requirements for medical devices to be designated "STERILE". Requirements for terminally sterilized medical devices	EN 556-1:2001 + AC:2006
ISO 594-1:1986	Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment. General requirements	EN 20594-1:1993 + A1:1997 + AC:1996
ISO 594-2:1998	Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment -- Part 2: Lock fittings	-
EN 980:2008	Symbols for use in the labelling of medical devices	EN 980:2008
EN 1041:2008 + Amd.1:2013	Information supplied by the manufacturer of medical devices	EN 1041:2008+A1:2013
EN 1707:1996	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Lock fittings	EN 1707:1996
ISO 2233:2000	Packaging -- Complete, filled transport packages and unit loads -- Conditioning for testing	EN ISO 2233:2001
ISO 10555-1:2013	Intravascular catheters -- Sterile and single-use catheters -- Part 1: General requirements	EN ISO 10555-1:2013
ISO 10555-4:2013	Intravascular catheters -- Sterile and single-use catheters -- Part 4: Balloon dilatation catheters	EN ISO 10555-4:2013
ISO 10993-1:2009 + Cor.1:2010	Biological evaluation of medical devices. Evaluation and testing within a risk management process	EN ISO 10993-1:2009 + AC:2010
ISO 10993-3:2014	Biological evaluation of medical devices. Tests for genotoxicity, carcinogenicity and reproductive toxicity	EN ISO 10993-3:2014
ISO 10993-4:2002 + Amd.1:2006	Biological evaluation of medical devices. Selection of tests for interactions with blood	EN ISO 10993-4:2009
ISO 10993-5:2009	Biological evaluation of medical devices. Tests for in vitro cytotoxicity	EN ISO 10993-5:2009
ISO 10993-6:2007	Biological evaluation of medical devices. Tests for local effects after implantation	EN ISO 10993-6:2009
ISO 10993-7:2008 + Cor.1:2009	Biological evaluation of medical devices. Ethylene oxide sterilization residuals	EN ISO 10993-7:2008 + AC:2009
ISO 10993-9:2009	Biological evaluation of medical devices. Framework for identification and quantification of potential degradation products	EN ISO 10993-9:2009
ISO 10993-10:2010	Biological evaluation of medical devices. Tests for irritation and skin sensitization	EN ISO 10993-10:2013
ISO 10993-11:2006	Biological evaluation of medical devices. Tests for systemic toxicity	EN ISO 10993-11:2009
ISO 10993-12:2012	Biological evaluation of medical devices. Sample preparation and reference materials	EN ISO 10993-12:2012

ISO 10993-15:2000	Biological evaluation of medical devices. Identification and quantification of degradation products from metals and alloys	EN ISO 10993-15:2009
ISO 10993-18:2005	Biological evaluation of medical devices. Chemical characterization of materials	EN ISO 10993-18:2009
ISO 11135-1:2007	Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices	EN ISO 11135-1:2007
ISO 11138-1:2006	Sterilization of health care products. Biological indicators. General requirements	EN ISO 11138-1:2006
ISO 11138-2:2006	Sterilization of health care products. Biological indicators. Biological indicators for ethylene oxide sterilization processes	EN ISO 11138-2:2009
ISO 11607-1:2006 + Amd.1:2014	Packaging for terminally sterilized medical devices. Requirements for materials, sterile barrier systems and packaging systems	EN ISO 11607-1:2009 + A1:2014
ISO 11607-2:2006 + Amd.1:2014	Packaging for terminally sterilized medical devices. Validation requirements for forming, sealing and assembly processes	EN ISO 11607-2:2006 + A1:2014
ISO 11737-1:2006 + Cor.1:2007	Sterilization of medical devices. Microbiological methods. Determination of a population of microorganisms on products	EN ISO 11737-1:2006 + AC:2009
ISO 11737-2:2009	Sterilization of medical devices. Microbiological methods. Tests of sterility performed in the definition, validation and maintenance of a sterilization process	EN ISO 11737-2:2009
DIN 13273-7:2003	Catheters for medical use - Part 7: Determination of the x-ray attenuation of catheters; Requirements and testing	-
ISO 13485:2003 + Cor.1:2009	Medical devices - Quality management systems - Requirements for regulatory purposes	EN ISO 13485:2012 + AC:2012
ISO 14155:2011 + Cor.1:2011	Clinical investigation of medical devices for human subjects — Good clinical practice	EN ISO 14155:2011
ISO 14630:2012	Non-active surgical implants. General requirements	EN ISO 14630:2012
ISO 14971:2007 + Cor.1:2007	Medical devices. Application of risk management to medical devices	EN ISO 14971:2012
ISO 15223-1:2012	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	EN ISO 15223-1:2012
ISO 25539-2:2012	Cardiovascular implants. Endovascular devices. Vascular stents	EN ISO 25539-2:2012

Laws and Directives

Reference	Title
MDD 93/42/EEC	COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 concerning medical devices (OJ L 169, 12.7.1993, p. 1)
MPD 2001/83/EC	DIRECTIVE 2001/83/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 6 November 2001 on the Community code relating to medicinal products for human use (OJ L 311, 28.11.2001, p. 67)
440/2008 Part B.13/14	Part B: METHODS FOR THE DETERMINATION OF TOXICITY AND OTHER HEALTH EFFECTS -- B.13/14. MUTAGENICITY: REVERSE MUTATION TEST USING BACTERIA
440/2008 Part B.17	Part B: METHODS FOR THE DETERMINATION OF TOXICITY AND OTHER HEALTH EFFECTS -- B.17. MUTAGENICITY -- IN VITRO MAMMALIAN CELL GENE MUTATION TEST
21 CFR 58	GOOD LABORATORY PRACTICE FOR NONCLINICAL LABORATORY STUDIES
PFSB/ELD/OMDE Yakushokuki Notification No.1109-1 (2011)	Guidance on ethylene oxide sterilization residuals on medical devices under the Pharmaceutical Affairs Law after the establishment of applicable Japanese Industrial Standard.
Federal Republic of Germany, ChemG Anhang 1 (zu § 19a Abs. 1) English title: ChemG Appendix 1 (to § 19a Clause 1)	Gesetz zum Schutz vor gefährlichen Stoffen (Chemikaliengesetz - ChemG) - Anhang 1 (zu § 19a Abs. 1) Grundsätze der Guten Laborpraxis (GLP) English title: Chemical act of the Federal Republic of Germany Appendix 1 (to § 19a clause 1) - Principles of Good Laboratory Practice

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