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
TECHNOLAS™
P E R F E C T V I S I O N

ZYOPTIX DIAGNOSTIC WORKSTATION³

ZYOPTIX® Diagnostic Workstation 3

User Manual

Version 2

UM-100006805 GB

SKU: 70005910

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Publishing details

ZYOPTIX® Diagnostic Workstation 3
User Manual
Version 2
[1.7]

Manufacturer

Technolas Perfect Vision GmbH
Messerschmittstr. 1 + 3
80992 München
Germany
www.TechnolasPV.com

For Customer Service or Technical Support:

Contact your local Customer Service or Technical Support.

Copyright

This document is protected by copyright. All rights reserved. Notably, no portion of this publication may be reproduced or made available to third parties, in any form, online or offline, without the written consent of Technolas Perfect Vision GmbH. Further, no portion of this publication may be translated into a language or code that can be utilized for machines, in particular, data processing systems, or the internet. This document may exclusively be used to properly use the products described herein in accordance with the contractual terms for such product.

Copyright © 2014 Technolas Perfect Vision GmbH
Subject to technical changes.

Document No: UM-100006805 GB
Edition: 11/2014

TABLE OF CONTENTS

A General information

A.1	Information about this user manual.....	A - 1
A.2	Explanation of symbols.....	A - 1
A.3	Customer service.....	A - 2
A.4	Product names used in this manual.....	A - 3

1 Introduction

1.1	ZDW 3.....	1 - 1
1.2	Scope of delivery.....	1 - 3
1.3	List of abbreviations.....	1 - 4

2 Safety

2.1	Intended use.....	2 - 1
2.2	Important safety considerations.....	2 - 2
2.3	Customer responsibilities.....	2 - 2
2.4	Responsible organization.....	2 - 3
2.5	Personnel requirements.....	2 - 3
2.6	General hazards.....	2 - 4
2.6.1	Fire hazard.....	2 - 4
2.6.2	Electrical hazard.....	2 - 5
2.6.3	Mechanical hazard.....	2 - 5
2.6.4	Damage to property.....	2 - 5
2.6.5	Warning signals and alarms.....	2 - 6
2.7	Labeling.....	2 - 6
2.7.1	Labels in the work area.....	2 - 6
2.7.2	Labels location on the diagnostic device.....	2 - 7

3 Transportation and storage

3.1	Transportation.....	3 - 1
3.2	Visual inspection after delivery.....	3 - 1

TABLE OF CONTENTS

3.3	Unpacking.....	3 - 2
3.4	Storage before installation.....	3 - 2
4	Installation and commissioning	
5	Service and calibration	
6	Clinical applications and available maps	
6.1	ORBSCAN 3.....	6 - 1
6.1.1	Clinical usage.....	6 - 1
6.1.2	Patient selection criteria.....	6 - 1
6.1.3	ORBSCAN 3 maps.....	6 - 1
6.2	ZYWAVE 3.....	6 - 4
6.2.1	Clinical usage.....	6 - 4
6.2.2	Patient selection criteria.....	6 - 5
6.2.3	ZYWAVE 3 display options of measurement results.....	6 - 5
6.3	Treatment Calculator.....	6 - 12
6.3.1	Clinical usage.....	6 - 12
6.3.2	Available treatment types and functionality.....	6 - 12
6.3.3	Inclusion and exclusion criteria.....	6 - 14
7	Hardware	
7.1	Overview of the main components.....	7 - 1
7.1.1	Center table.....	7 - 2
7.1.2	Lift tables.....	7 - 3
7.1.3	ORBSCAN 3 and ZYWAVE 3 measurement heads.....	7 - 4
7.1.3.1	Positioner with joystick.....	7 - 5
7.1.4	User control screen.....	7 - 6
7.2	Accessories.....	7 - 7
7.3	Operating the hardware.....	7 - 9
7.3.1	Switching on and off.....	7 - 10
7.3.2	Positioning and triggering measurements.....	7 - 10
7.3.3	Adjusting the user control screen settings.....	7 - 11

TABLE OF CONTENTS

7.3.4	Performing system checks.....	7 - 11
7.4	Cleaning instructions.....	7 - 12
8	Software	
8.1	Working with the software.....	8 - 1
8.1.1	Logging in.....	8 - 3
8.1.2	Logging out.....	8 - 4
8.1.3	Changing the user.....	8 - 4
8.2	Calculation settings.....	8 - 4
8.2.1	About the user-defined nomogram.....	8 - 5
8.2.2	Setting a user-defined nomogram.....	8 - 6
8.2.3	Setting default treatment values.....	8 - 7
8.3	Patient administration.....	8 - 8
8.3.1	Selecting patients.....	8 - 8
8.3.2	Searching for patient data.....	8 - 10
8.3.3	Adding patients.....	8 - 15
8.3.4	Editing patient data.....	8 - 18
8.3.5	Deleting patient data.....	8 - 20
8.4	User administration.....	8 - 20
8.4.1	User roles.....	8 - 20
8.4.2	Selecting users.....	8 - 22
8.4.3	Adding users.....	8 - 23
8.4.4	Editing user data.....	8 - 26
8.4.5	Deleting user data.....	8 - 26
8.5	Settings.....	8 - 28
8.5.1	Defining system settings.....	8 - 29
8.5.2	Defining ORBSCAN 3 map settings.....	8 - 31
8.5.3	Performing ORBSCAN 3 calibration checks (alignment PMMA).....	8 - 33
8.5.4	Performing ZYWAVE 3 calibration checks.....	8 - 38
8.5.5	Viewing the product information.....	8 - 44
9	Measurements	
9.1	Preparing for performing measurements.....	9 - 1

TABLE OF CONTENTS

9.1.1	Disinfecting the head and chin rests.....	9 - 1
9.1.2	Selecting or adding a patient.....	9 - 1
9.1.3	Instructing and positioning the patient.....	9 - 2
9.2	Performing ORBSCAN 3 measurements.....	9 - 2
9.2.1	Aligning the patient and starting the ORBSCAN 3 measurement.....	9 - 3
9.2.2	Viewing ORBSCAN 3 measurement results.....	9 - 6
9.3	Performing ZYWAVE 3 measurements.....	9 - 14
9.3.1	Aligning the patient and starting the ZYWAVE 3 measurement.....	9 - 14
9.3.2	Verifying the measurement images.....	9 - 21
9.3.3	Reading the ZYWAVE 3 measurement results.....	9 - 22
9.4	Taking iris images on ZYWAVE 3.....	9 - 23
9.5	Performing the pupilometry on ZYWAVE 3.....	9 - 24
9.6	Working with measurements.....	9 - 27
9.6.1	Selecting and closing measurements.....	9 - 27
9.6.2	Deleting measurements.....	9 - 29
9.6.3	Comparing measurements.....	9 - 29
9.6.3.1	Comparing ORBSCAN 3 measurements	9 - 29
9.6.3.2	Comparing ZYWAVE 3 measurements.....	9 - 31
9.6.4	View data on the USB memory stick.....	9 - 33
9.6.5	Exporting measurements.....	9 - 34
9.6.6	Exporting ZYWAVE 3 pupil images.....	9 - 35
9.6.7	Printing measurement results.....	9 - 37
9.6.7.1	Printing ORBSCAN 3 measurement results.....	9 - 38
9.6.7.2	Printing ZYWAVE 3 display options of measurement results.....	9 - 39
9.6.8	Performing treatment calculations.....	9 - 39
9.6.8.1	Selecting treatments and diagnostic data.....	9 - 40
9.6.8.2	Adding treatments.....	9 - 45
9.6.8.3	Entering and adjusting correction values.....	9 - 53
9.6.8.4	Printing treatment calculations.....	9 - 56
10	Troubleshooting	
10.1	Software messages.....	10 - 2
11	Decommissioning and disposal	

TABLE OF CONTENTS

12 Technical data

12.1	General specifications.....	12 - 1
12.1.1	Dimensions and weights.....	12 - 1
12.1.2	Essential performance.....	12 - 1
12.2	Operating conditions.....	12 - 2
12.3	Transport and storage conditions.....	12 - 2
12.4	Electrical connection.....	12 - 2

13 Index

14 Appendix

,
i
7
3
9
9
0
.5
i3
56

-2

TABLE OF CONTENTS

A General information

To order additional copies of this user manual:

If additional copies of this user manual are required, refer to the SKU number on the cover page and contact your Authorized Technolas Perfect Vision Representative.

A.1 Information about this user manual

This user manual enables safe and efficient handling of the product and its accessories. It is an important part of the product, and must be kept close to the product where it is permanently accessible to the personnel.

Before starting any treatments, the personnel must read the manual thoroughly and understand its contents. Users must comply with all specified safety instructions and operating instructions to ensure safe operation.

In addition, local accident prevention regulations and general safety instructions must be observed in the area where the product is used.

The product images in this user manual are intended to provide a basic understanding of the product and may differ from the actual design you have on site. Screenshots were created using an exemplary patient and exemplary values. Depending on your system's configuration the depicted screenshots may not be identical to the screen shown on your system.

The software user interface may only be available in English. Where applicable, the Appendix section of the user manual contains a list of the English software texts together with the appropriate translation.

A.2 Explanation of symbols

Safety notes

 **DANGER**

Safety notes in this user manual are introduced with symbols and a signal word that describe the severity of the hazard.

Follow these safety notes and proceed with caution to avoid accidents, personal injury, damage to property, and to ensure patient safety.

 **WARNING**

This symbol and the word "DANGER" indicate an imminent hazardous situation that, if not avoided, could cause death or serious injury.

 **CAUTION**

This symbol and the word "WARNING" indicate a possible hazardous situation that, if not avoided, could cause death or serious injury.

NOTICE

This symbol and the word "CAUTION" indicate a possible hazardous situation that, if not avoided, could result in minor or moderate injury.

This symbol and the word "NOTICE" indicate a possible hazardous situation that, if not avoided, could result in damage to property or the environment.

A GENERAL INFORMATION

Hints, recommendations, and examples

This symbol emphasizes useful tips and recommendations as well as information for the efficient and trouble-free use of the product.

Example of ...

This field shows examples of previous instructions or information.

Safety notes in step-by-step instructions

Safety notes are included in step-by-step instructions. This format allows the user to be aware of safety issues while performing tasks. The safety symbols are also part of the step-by-step instructions.

Symbols in this user manual

The following symbols and markings are used in this manual to indicate guidelines, descriptions of results, lists, and references:

Symbol	Indicates
	Step-by-step instructions
	Condition or automatic sequence as a result of action
	References to chapters in this user manual
	Lists or list entries without a certain sequence
[Ctrl]	Keys on the keyboard, mouse buttons, key switches, hardware switches, or hardware buttons
'Start'	Active (clickable) GUI elements
'Patient selection'	Passive elements
c:\Program\...	Text input
'Administration → User Registry'	Menu path

A.3 Customer service

Our customer service division is available to provide technical support and information. If any questions arise that are not answered in this user manual, do not hesitate to contact your Authorized Technolas Perfect Vision Representative for assistance.

A.4 Product names used in this manual

For the sake of better readability and depending on the context, product names may appear in abbreviated form.

Complete product name	Name used in this manual
ZYOPTIX® Diagnostic Workstation 3	ZDW 3*
ORBSCAN® 3 Anterior Segment Analyzer	ORBSCAN 3
Treatment Calculator	Treatment Calculator**
ZYWAVE® 3 Wavefront Analyzer	ZYWAVE 3
PROSCAN	PROSCAN
SUPRACOR™	SUPRACOR
ZYOPTIX® HD	ZYOPTIX HD
ZYOPTIX® ORBSCAN® Ilz Anterior Segment Analyzer	ORBSCAN Ilz
ZYOPTIX® ZYWAVE® II Wavefront Analyzer	ZYWAVE II
TECHNOLAS® TENEO™ 317	TENEO™ 317***
Microsoft® Windows®	Windows

* or generic term 'diagnostic device'

** or generic term 'calculation software'

*** or generic term 'laser system'

or-
not
for

A GENERAL INFORMATION

1 Introduction

This chapter provides general information about the ZDW 3.

1.1 ZDW 3

The diagnostic device combines the two individual devices ZYWAVE® 3 Wavefront Analyzer and ORBSCAN® 3 Anterior Segment Analyzer into one workstation, where the user can operate both of them comfortably without having to change places. ZYWAVE 3 and ORBSCAN 3 each have their own measurement head, while the software for the two devices is installed on the same computer. This enables patient data exchange and easy switching between the two applications.

The software installed on the diagnostic device supports an automatic eye (OD/OS) recognition.

The ZYWAVE 3 and ORBSCAN 3 generated data can be read from the Treatment Calculator software and thus be utilized in the ZYOPTIX treatments at the TECHNOLAS® TENEO™ 317 laser system.

ORBSCAN 3

- Creates data in order to calculate a refractive treatment to be performed on the TENEO™ 317
- GUI supports the display of:
 - quad map view: displays four different maps at the same time
 - difference map view: displays the difference between maps of the same type, but with different measurement date from the same eye (OS or OD)
 - map data and calculation results

ORBSCAN 3 available maps

Map name	Class	Comment
'Anterior map'	Corneal elevation	Both spheres (shapes) are floating. ■ BFS: Best Fit Sphere ■ BFA: Best Fit Aconic ■ Sph./Cyl.: spherocylindrical shape
'Posterior map'	Corneal elevation	
'Pachymetry'	Pachymetry	
'Axial power, keratometric map'	Axial power	
'Axial power, anterior map'	Axial power	
'Axial power, posterior map'	Axial power	
'Axial power, total map'	Axial power	

1 INTRODUCTION

Map name	Class	Comment
'Tangential power, keratometric map'	Tangential power	
'Anterior chamber, apex plane'	Anterior chamber	
'Anterior chamber, epithelium'	Anterior chamber	
'Anterior chamber, endothelium'	Anterior chamber	
'Anterior chamber, normal depth'	Anterior chamber	
'Mean power'	Mean power	

ZYWAVE 3

- Displays images as raw data output:
 - Iris image
 - Centroid image with overlaid centroids connected through horizontal / vertical lines
- Provides iris code with the measurement data for later iris recognition at the TENEO™ 317
- Displays the following measurement display options:
 - Summary
 - 2D Plot
 - Predicted Phoropter Refraction (PPR) vs. Pupil diameter
 - Higher Order Point Spread Function (PSP)
 - Raw data
 - Normal band
 - Difference

External database

- External user, patient and measurement database as well as calibration and calibration check database:
 - Compatible with other devices, such as the TENEO™ 317
 - Provides access to patient data, diagnostic data, treatment data, user data and calibration and calibration check data of the diagnostic device.

1.2 Scope of delivery

Mandatory components

The scope of delivery of the diagnostic device includes the following mandatory components:

- ORBSCAN 3 measurement head with dust cover
- ZYWAVE 3 measurement head with dust cover
- Medical all-in-one PC (referred to as 'user control screen' in the user manual) with power supply
- Keyboard
- Mouse
- Lift tables (identical for ORBSCAN 3 and ZYWAVE 3)
- Center table
- Connection bars (left and right)
- Check tool case including the ORBSCAN 3 and the ZYWAVE 3 calibration check tools
- USB memory stick (approved by Technolas Perfect Vision GmbH)
- User manual ZDW 3 Version 1 or higher.

Optional components

If necessary, the following component can be ordered.

- Database server unit, including power cables and documentation.

1.3 List of abbreviations

Abbreviation	Full name / meaning
D	Diopter
GUI	Graphical User Interface
ID	Identification number; Identifier
IRS	Iris data file
HOA	Higher Order Aberrations
K-value	Keratometric value
HCb	Head control board
LASIK	Laser-Assisted In Situ Keratomileusis
LED	Light-Emitting Diode
OD	Oculus Dexter (right eye)
OS	Oculus Sinister (left eye)
PPR	Predicted Phoropter Refraction
PPR _{3.5mm}	Wavefront-derived refraction at pupil size 3.5 mm
preOp	Pre-operative
postOp	Post-operative
PSF	Point Spread Function
PRK	Photorefractive Keratectomy
PTK	Phototherapeutic Keratectomy
Q-value	Anterior corneal asphericity
®	Registered Trade Mark
SKU	Stock Keeping Unit
™	Trade Mark
TPV	Technolas Perfect Vision

2 Safety

This section provides an overview of the important safety aspects that are necessary to ensure optimal protection of patients and personnel as well as for the safe and trouble-free measurement of patients.

2.1 Intended use

The ZDW 3 has been designed and constructed exclusively for the intended use described in this section.

The medical device 'ZYOPTIX Diagnostic Workstation 3' unites the two individual devices ZYWAVE 3 Wavefront Analyzer and ORBSCAN 3 Anterior Segment Analyzer in an interactive workspace. The recorded measurement data and derived clinical data can be used by the integrated Treatment Calculator to calculate patient-customized treatments by creating a laser ablation pattern for refractive surgery. A shared database is available to store and maintain patient data, examination data and treatment data. The measurements are performed without physical contact to the eye. The device must only be used at adequate facilities by a healthcare professional after training on this device.

The ZYWAVE 3 Wavefront Analyzer measures and displays the ocular aberrations by means of wavefront technology.

The ORBSCAN 3 Anterior Segment Analyzer is designed to scan and display the geometry of the eye's anterior segment and provides a variety of maps of the eye's corneal optical system by means of slit scanning technology and an advanced Placido disc system.

The usage of the diagnostic device is limited to the procedures described in this manual.

Any use that exceeds, or is different from, the specified intended use is considered to be misuse.

Danger of injury due to misuse!

Misuse of the diagnostic device may cause dangerous situations and lead to severe injuries.

- Do not use the device on patients who are not at least 18 years of age.
- Do not use the device with devices and accessories that are not specified in this user manual.
- Do not use the device with software that is not authorized and certified by Technolas Perfect Vision GmbH.

Claims of any type due to damage arising from misuse are excluded.

NOTICE

Damage to property due to inadequate connection of medical and non-medical devices!

Inadequate connection of non-medical devices can cause property damage and lead to severe injuries.

- Always connect non-medical devices to adequate interfaces (e.g. Ethernet).
- Do not connect medical and non-medical devices to multiple plugs.

2.2 Important safety considerations

⚠ WARNING

Danger of injury due to user negligence!

Always double-check the following before measurement:

- Are you measuring the right patient?
- Was the correct eye selected?
- Does the entered data match the patient to be measured?

Failure to verify any of the above-mentioned aspects may result in serious permanent patient injury.

2.3 Customer responsibilities

Customer

The term "customer" refers to a person who operates the product for trade or commercial purposes. A customer can also employ a third party who uses the product, but the customer bears the legal product liability for protecting the user, the personnel, or third parties during the operation of the product.

Responsibilities

The product is used for commercial purposes. Product users are therefore required to observe their obligations with regard to occupational health and safety.

In addition to the safety instructions provided in this user manual, the safety, accident prevention, and environmental regulations applicable at the place of installation must be observed, particularly as follows:

- The customer must ensure that personnel who handle the product have read and understood this user manual.
- The customer must ensure that unauthorized persons have no access to the operating room and the product.
- The customer must ensure that the service and calibration intervals specified in this user manual are observed.
- The customer must provide all users with personal protective equipment, if necessary.
- The customer must ensure that the product is only operated with the installed software.

2.4 Responsible organization

The responsible organisation can be, for example, a public hospital, a private person (ophthalmologist), or a private clinic.

Education and training is included in the term 'use'.

- The responsible organization must ensure that hazards of network coupling of the device must be observed. Such hazards include, for example, virus protection, data security, network protection by firewall systems, and ensuring sufficient bandwidth.
- System modifications and use in combination with other equipment is limited to TPV-approved equipment. Otherwise, the responsible organization is in charge of proper risk management.

NOTICE

Danger of property damage due to use of unauthorized accessories / devices!

Use of unauthorized accessories / devices may result in data loss (danger by viruses or other malicious software).

- Do not use accessories / devices which are not authorized by Technolas Perfect Vision GmbH.

WARNING

Danger of property damage due to unauthorized system modifications!

Unauthorized system modifications may result in data loss or can cause serious injury and substantial damage to property.

- System modifications must be performed by Technolas Perfect Vision GmbH only.

2.5 Personnel requirements

Personnel qualifications

WARNING

Danger of injury if personnel are insufficiently qualified!

If unqualified personnel perform work on the diagnostic device or are in the danger zone of the diagnostic device, hazards may arise that can cause serious injury and substantial damage to property.

- All measurements must be performed by appropriately qualified personnel.

WARNING

Danger due to use by unauthorized persons!

Unauthorized persons who do not meet the requirements described here are not acquainted with the dangers associated with the product.

- Unauthorized persons must not use the product.

Personnel in general

All personnel are expected to perform their work in a reliable manner. Persons with impaired reactions due to the consumption of drugs, alcohol, medication, or any other reason are prohibited.

Required qualifications for diagnostic device personnel:

Authorized Service Technician

Based on their professional, product-related training, professional experience, and knowledge of all relevant regulations, Authorized Service Technicians can service ophthalmologic systems and detect and prevent possible hazards.

User

Users are medical specialists who specialize in the diagnosis, treatment, and surgery of the eye.

Users have been trained by Technolas Perfect Vision staff or personnel authorized by Technolas Perfect Vision in proper handling of the system including special precautionary measures for working with diagnostic and laser-emitting systems. Users know the content of all regulations, guidelines, and standards that are applicable for the safe use of the system and can implement the requirements mentioned in them.

Users have the necessary technical knowledge of the application area of the system for which they are responsible and rigorously follow all hygiene requirements for operating rooms and the use of medical products.

2.6 General hazards

This section specifies several risks that should be avoided based on a risk assessment. In order to reduce health risks and avoid dangerous situations, observe the safety instructions listed in this chapter and the safety instructions contained in other chapters of this user manual.

Failure to follow the directions and safety information given in this user manual may result in severe injuries.

2.6.1 Fire hazard

⚠ WARNING

Danger to life in the event of fire due to highly flammable substances!

Highly flammable substances, liquids, or gases may catch fire and cause serious or fatal injuries.

- **Do not use flammable substances near the product.**
- **Be sure that suitable extinguishing agents are available such as a fire blanket and a fire extinguisher.**

⚠ DANGER**2.6.2 Electrical hazard****Life-threatening hazard due to electric shock!**

There is an imminent life-threatening hazard of electric shock if live parts are touched. Damage to insulation or to specific components can pose a life-threatening hazard.

- Only a qualified electrician should perform work on the electrical equipment.
- Immediately switch off the device and have it repaired if the insulation of the power supply is damaged.
- Keep moisture away from live parts, for example, switch plugs or the computer system. Moisture can cause short circuits.
- Do not use multiple sockets or extension cables since this impacts the function of protective earth and may lead to electrical shock.
- The device must only be connected to mains supply with proper protective earth connection.
- Do not touch the patient at the same time as live parts, for example, USB ports, connectors or back of monitor.

⚠ CAUTION

Danger of injury due to electrical hazards if non-medical devices are placed in the patient vicinity!

- Do not place non-medical devices in the patient vicinity as this may cause dangerous situations and lead to severe injuries.

2.6.3 Mechanical hazard**⚠ CAUTION****Danger of injury from tripping!**

If the power cable is laid in the open, there is a risk that someone can trip over the cable and be injured. The power plug and the unit can also be damaged.

- Lay the power cable in a place where it cannot be tripped over by personnel.
- Always connect the power plug in a safe position where it cannot be damaged or destroyed.

2.6.4 Damage to property**NOTICE****Damage to property due to heavy objects or liquids!**

The product is very sensitive. Heavy objects placed on the product may lead to equipment damage. Liquids placed on the product may spill, which will cause equipment damage.

- Do not put heavy objects on the product.
- Do not put liquids on the product.

NOTICE**Danger of property damage due to malicious software!**

USB memory sticks connected to the product must not be infected with viruses or other malicious software. Malicious software and viruses may cause damage that cannot be repaired.

- Do not use USB memory sticks infected with viruses or other malicious software.
- Before connecting the product to a LAN interface, verify that connected equipment complies with the requirements for the LAN interface.

Visual signals**2.6.5 Warning signals and alarms**

Software and hardware warning messages are displayed in the GUI of the ZDW 3 software, describing the fault that occurs in the diagnostic device.

- For information about the meaning of icons and buttons in the software GUI, refer to *Chapter 8.1 'Working with the software' on page 8 - 1.*
- For a description of possible errors and related causes, and how to rectify them, refer to *Chapter 10 'Troubleshooting' on page 10 - 1.*

2.7 Labeling**⚠ WARNING****Danger of injury due to illegible labels!**

Over time, labels and signs can become dirty or illegible due to wear and tear or fading. If labels become illegible, dangers cannot be recognized, necessary operating instructions cannot be followed, and injuries can occur.

- Always keep the labels in good condition so that they are always legible.
- If a label is damaged or has become illegible, contact an Authorized Service Technician.

2.7.1 Labels in the work area

The following labels must be displayed in the work area.

Symbol	Label	Meaning
	No smoking	Fire hazards exist from flammable or potentially explosive solid, liquid, or gaseous substances.
	Unauthorized access prohibited	Only persons authorized by the customer may enter the operating room.

2.7.2 Labels location on the diagnostic device

Front view

The following figures show the label location on the diagnostic device:


Fig. 1: ZDW 3 front view - labels location

Pos.	Symbol	Label
1		Follow the instructions in the user manual
2	A/N 50003002 ORBSCAN 3 head, assy. S/N: 50003002-0003 27322736	ORBSCAN 3 measurement head: serial number
3	A/N 50004019 ZYWAVE 3 head, assy. S/N: 50004019-0002 27322731	ZYWAVE 3 measurement head: serial number

Side view



Fig. 2: ZDW 3 side view - labels location

Pos.	Symbol	Label
1		Type plate*
2		WEEE symbol**

Pos.	Symbol	Label
3		Laser aperture
4		Laser product - Class 1 with laser beam warning

*Type plate


Fig. 3: Type plate ZDW 3

Label	Symbol	Value
Product type	-	ZYOPTIX® Diagnostic Workstation 3
Manufacturing date		mm/yyyy
Serial number		ZDW 3-XXXX
IEC60601-1	-	Class I
IEC60601-1		Applied part Type BF (chin rest, head rest, patient hand holds)
IEC 60529	-	IP2X
Rated voltage	-	100 - 240 V ~
Power input	-	max. 260 VA
Rated supply frequency	-	50/60 Hz

Label	Symbol	Value
Manufacturer		Technolas Perfect Vision GmbH and address
TPV label number	-	27322649

****Disposal of electronic equipment**

Fig. 4: WEEE symbol

The European Union's WEEE Directive requires manufacturers to label their electronic equipment properly and to provide for recycling at the end of its useful life. Technolas Perfect Vision GmbH fully complies with the WEEE Directive. This symbol on the product indicates that it must be separately collected and recycled in order to protect human health and the environment.

In that regard, Technolas Perfect Vision GmbH will arrange the collection and recycling of end-of-life Technolas Perfect Vision GmbH equipment.

To initiate the collection of end-of-life equipment, please contact an Authorized Service Technician.

3 Transportation and storage

This chapter provides information regarding the transportation, visual inspection and packaging, unpacking, and storage of the diagnostic device before installation.

For information on transport and storage conditions, also refer to Chapter 12.3 'Transport and storage conditions' on page 12 - 2.

3.1 Transportation

The diagnostic device will be delivered to you by a delivery company. Depending on the form of transportation, it may be packed in wooden boxes, or the individual components will be transported with other protective material.

The packaging may be removed only by Authorized Service Technicians.

Do not move or transport the diagnostic device yourself. Any damage caused by unauthorized movement or transport of the diagnostic device will void the warranty.

CAUTION

Hazards due to unauthorized transport!

The product is very sensitive to any changes in temperature, to shocks, or vibrations. Improper transport may result in property damage and impaired treatments, which can cause severe injuries.

The product must only be transported by a delivery company that is authorized by Technolas Perfect Vision GmbH. This requirement applies to delivery, to in-house transports, and to external transports.

3.2 Visual inspection after delivery

The product and the required documents are given to the carrier in perfect condition. The documents state the type and content of the delivered goods. Check the delivery note to ensure the completeness of the delivery. The carrier is responsible for delivering the product promptly and in perfect condition. Should a loss of parts or damage occur during transportation, the delivery company is responsible.

If there are defects or damage that have occurred during transport, perform the following steps:

- Inform the nearest service division immediately. A delay in reporting the damage may result in denial of the damage claim.
- Urge the delivery person to write down the kind of damage on the freight note immediately and sign the freight note as confirmation.

- Under certain circumstances, an expert's opinion on the damage cost or the value of lost goods can only be provided by a representative of the transport insurance company.
- Fill in a form for the insurance company immediately. Payment for lost goods or for repair of damaged goods can only be obtained from the transport insurance company.

If you do not detect any damage of the outside packaging, but the product is damaged on the inside and if you think that this damage has resulted from improper transportation, inform your distributor and the delivery company immediately and insist on an inspection of the product. Keep the packaging material until all questions of the delivery agent and the insurance company have been clarified.

3.3 Unpacking

After delivery, the product is unpacked and installed by an Authorized Service Technician.

WARNING

Hazards due to unauthorized unpacking!

Removal of the transport safety device in the wrong sequence may result in damage.

- **The product must only be unpacked by Authorized Service Technicians.**
- **The product must not be connected before Authorized Service Technicians have made sure that it is in proper working condition before initial operation.**

Technolas Perfect Vision GmbH does not assume any liability for any damage to the product arising from unauthorized unpacking.

The product is packed so that transport damage will not occur under most circumstances. However, should you notice any irregularities or have questions concerning transportation and packaging, contact Technolas Perfect Vision GmbH and the delivery agent as soon as possible.

3.4 Storage before installation

After delivery to your clinic or practice, note the following:

- The individual parts of the shipment must be stored in a secure and dry place before installation and initial use.

4 Installation and commissioning

Initial installation and commissioning is exclusively performed by Authorized Service Technicians.

Unauthorized installation and commissioning will void the warranty.

Danger of injury due to improper calibration!

The diagnostic device may decalibrate during transportation.

- **The diagnostic device must be recalibrated or calibration checked upon installation by an Authorized Service Technician.**

Longer periods out of operation

During longer periods out of operation, the conditions specified under Chapter 3 'Transportation and storage' on page 3 - 1 must be followed.

After longer periods out of operation, the diagnostic device can be started as described in Chapter 8.1 'Working with the software' on page 8 - 1.

Room conditions

For optimum results, make sure that the room where the diagnostic device will be utilized can be darkened. It is recommended that overhead lights be turned off during the measurement since extraneous light illuminating the cornea may interfere with edge detection, which could result in a poor exam. Dim illumination in the room to aid in maneuverability is acceptable.

CAUTION

Danger of injury due to high ambient light level!

The device may be disturbed by high ambient light level. This may lead to erroneous measurements.

- **Make sure that the ambient light level is reduced and that the defined requirements for the operating conditions are met.**

CAUTION

Danger of injury due to high magnetic fields!

The device is disturbed by high magnetic fields (for instance by Magnetic Resonance Imaging, MRI). This may lead to position errors of the motors and to erroneous measurements.

- **Make sure that magnetic fields are reduced as far as possible and that the defined requirements for the operating conditions are met.**

CAUTION

Danger of injury due to high electromagnetic fields!

The device may be disturbed by high electromagnetic fields (X-ray). This may lead to erroneous measurements.

- **Make sure that electromagnetic fields are reduced as far as possible and that the defined requirements for the operating conditions are met.**

5 Service and calibration

All service tasks are solely performed by Authorized Service Technicians.

Service tasks

- Inspection
- Calibration
- Function and (electrical) safety test

For more information about service tasks contact an Authorized Service Technician.

WARNING

Danger due to improperly completed service or calibration!

Service or calibration that is carried out by personnel other than Authorized Service Technicians may result in severe injuries to persons who perform the service, third parties or patients, and may damage the diagnostic device.

- **Service and calibration must only be performed by Authorized Service Technicians.**

NOTICE

Danger of property damage due to improperly performed maintenance!

Improperly performed maintenance may lead to serious injury and significant material damage.

- **Repairs and exchanges of the power-supply cord, fuse, etc. may only be performed by Authorized Service Technicians.**

When to perform service and calibration

The default service and calibration interval for ORBSCAN 3 and ZYWAVE 3 is 6 months.

One month before service and calibration interval will expire, a warning message will be displayed. In this case, contact an Authorized Service Technician as soon as possible for performing service and calibration.

If the service and calibration interval has expired and the device was not calibrated again, it is not possible to perform a new measurement.

Calibration checks

Calibration checks for ORBSCAN 3 and ZYWAVE 3 must be performed by the **user** (§ Chapter 8.5.3 'Performing ORBSCAN 3 calibration checks (alignment PMMA)' on page 8 - 33 and § Chapter 8.5.4 'Performing ZYWAVE 3 calibration checks' on page 8 - 38).

Service life

The expected service life is 5 years.

5 Service and calibration

All service tasks are solely performed by Authorized Service Technicians.

Service tasks

- Inspection
- Calibration
- Function and (electrical) safety test

For more information about service tasks contact an Authorized Service Technician.

WARNING

Danger due to improperly completed service or calibration!

Service or calibration that is carried out by personnel other than Authorized Service Technicians may result in severe injuries to persons who perform the service, third parties or patients, and may damage the diagnostic device.

- Service and calibration must only be performed by Authorized Service Technicians.

NOTICE

Danger of property damage due to improperly performed maintenance!

Improperly performed maintenance may lead to serious injury and significant material damage.

- Repairs and exchanges of the power-supply cord, fuse, etc. may only be performed by Authorized Service Technicians.

When to perform service and calibration

The default service and calibration interval for ORBSCAN 3 and ZYWAVE 3 is 6 months.

One month before service and calibration interval will expire, a warning message will be displayed. In this case, contact an Authorized Service Technician as soon as possible for performing service and calibration.

If the service and calibration interval has expired and the device was not calibrated again, it is not possible to perform a new measurement.

Calibration checks

Calibration checks for ORBSCAN 3 and ZYWAVE 3 must be performed by the **user** (☞ Chapter 8.5.3 'Performing ORBSCAN 3 calibration checks (alignment PMMA)' on page 8 - 33 and ☞ Chapter 8.5.4 'Performing ZYWAVE 3 calibration checks' on page 8 - 38).

Service life

The expected service life is 5 years.

6 Clinical applications and available maps

Use of the device

The use of the device is under the particular responsibility and is at the discretion of the user once the therapeutic benefits and possible complications and side effects have been considered.

6.1 ORBSCAN 3

6.1.1 Clinical usage

ORBSCAN 3

ORBSCAN 3 is a diagnostic device that scans the cornea of the human eye to provide reliable data about the geometrical shape of the outermost surface of the eye (anterior segment of the eye), the thickness of the cornea (pachymetry) and biometric data (anterior chamber depth). The data is obtained by scanning a slit light over the cornea and acquiring the resulting scattering images, combined with the evaluation of a Placido-type pattern to increase accuracy of the anterior surface measurement. Evaluation is done with the help of image processing.

Patient data

The diagnostic device provides all images obtained during data acquisition as raw data for the calculation of maps.

The following images are generated from the raw data and displayed via the diagnostic software:

- Placido disc image
- Alignment image
- Slitscan image: all slitscan images merged within one image
- Eye image averaged

For information about the ORBSCAN 3 images generated from the raw data, refer to Chapter 9.2.2 'Viewing ORBSCAN 3 measurement results' on page 9 - 6.

6.1.2 Patient selection criteria

Measurements performed by ORBSCAN 3 are without physical contact.

The measurements can be performed on patients who are at least 18 years of age.

6.1.3 ORBSCAN 3 maps

The results of ORBSCAN 3 measurements are displayed in maps that visualize the eye's corneal optical system.

From the different map types, four are combined into a **quad map** (Chapter 9.2.2 'Viewing ORBSCAN 3 measurement results' on page 9 - 6), which displays the four maps simultaneously.

Important terms

Keratometry	Related to the curvature of the eye's corneal surface. Measured in diopters (D) or - when expressed as the radius of a circle with the corresponding curvature - in units of mm.
K-reading	(Accurate) Keratometry reading
Pachymetry	Thickness of the eye's cornea
Power	'Power' as used in 'Power map', 'Refractive power', 'Mean power', 'Axial power': referring to the strength of an optical system to refract (change the direction of light beams), or the related curvature of an optical surface in D.
Refraction	A set of three values describing the refractive power – usually of the eye, e.g. as used to describe the power of spectacles. The three parameters are: sphere [D], cylinder (astigmatism) [D], and axis for cylinder [°]. They are typically displayed in the following format: sph. / cyl. / axis, e.g. -4.25/-1.15/170°.

The following ORBSCAN 3 maps are available:

Elevation maps

Elevation maps display the height of the cornea relative to a reference surface. Elevation maps have the appearance of a reversed axial curvature map (as seen with Placido image capturing systems).

Elevation maps	Description
Anterior BFS	Anterior surface corneal topography plotted against the best fit sphere.
Posterior BFS	Posterior surface corneal topography plotted against the best fit sphere. The posterior surface of the cornea will be a smaller sphere size than the anterior surface.
Anterior BFA	Anterior surface corneal topography plotted against an aconic sphere.
Posterior BFA	Posterior surface corneal topography plotted against an aconic sphere. The posterior surface of the cornea will be a smaller sphere size than the anterior surface.
Anterior sph./cyl.	[mm]
Posterior sph./cyl.	[mm]

Mean power map

Mean power maps display the refractive power and the arithmetic average of the principal curvatures for each data point. This map represents the variation of spherical power. This map is useful for determining the exact location and the effect of surface anomalies regardless of where they appear on the cornea since each point is independently defined, i.e. each point is a representation of its local sphere. This map is also useful when viewing extreme surface anomalies such as keratoconus and corneal transplants.

Mean power map Description

Mean power	Uses the standard index of refraction for its calculation.
------------	--

Pachymetry map

Pachymetry maps show the corneal thickness of the measured eye. These maps are useful in determining thin areas of the cornea and the location of the ablation zone where corneal tissue was removed via refractive procedures such as PRK or LASIK.

Axial power maps

Axial power maps represent the traditional method of computing the radius of curvature for the corneal surfaces based on the visual axis of the eye, which is converted to D. These maps are useful in determining separately the contribution of the anterior and posterior surfaces to the overall astigmatism found by axial methods.

Axial power maps Description

Keratometric map	Based on curvature of the surface along planes running from a center point radially to the periphery. The curvature being calculated by measuring perpendicularly from a surface tangent at each point along the plane to the intersection of the central axis. The length of the perpendicular is the radius of curvature which is converted to D using the standard keratometric index. Useful in finding the general curvature of the cornea.
Anterior map	Anterior corneal surface based on the radius of curvature being converted to D using a physiological index to account for each surface interface index of refraction, i.e., air-cornea.
Posterior map	Posterior corneal surface based on the radius of curvature being converted to D using a physiological index to account for each surface interface index of refraction, i.e., cornea-anterior chamber.
Total map	Represents the addition of the anterior-physiological and posterior-physiological maps.

Tangential power map

Tangential power maps compute the radius of curvature for the corneal surfaces and convert this information to D. These maps are useful in determining separately the contribution of the anterior surface to the overall astigmatism found by tangential methods.

Tangential power map**Description**

Keratometric

Anterior surface created from local radii of curvature using the standard keratometric index. Represents data from the Placido rings. Useful in finding small changes from the mid-periphery to the limbus.

Anterior chamber maps

Anterior chamber maps measure the anterior chamber depth from either surface of the cornea to the lens and iris plane.

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Apex Plane

Measured from a plane perpendicular to the axis of fixation on the anterior surface of the cornea to the lens and iris.

Epithelium

Measured from the anterior surface of the cornea to the lens and iris.

Endothelium

Measured from the posterior surface of the cornea to the lens and iris.

Normal Depth

Measures the ray light perpendicular to the anterior surface (a) where it refracts and bends light (b); and the depth between the posterior surface of the cornea to the anterior surface of the lens (c), i.e. $a + b + c = \text{normal depth}$.

6.2 ZYWAVE 3**6.2.1 Clinical usage****ZYWAVE 3**

ZYWAVE 3 is a diagnostic device which measures the corneal aberrations from the deviation in light beams reflected off the eye's retina to define the wavefront deformation as a function of the position on the pupil. Laser light, scattered on the retina, generates the signal data for the calculation of the wavefront deformation. The established measurement data of a patient's eye is saved allowing recalculation of the corneal aberration together with all the personal measurement data and information on the eye.

6.2.2 Patient selection criteria

Measurements performed by ZYWAVE 3 are without physical contact.

The measurements can be performed on patients who are at least 18 years of age.

6.2.3 ZYWAVE 3 display options of measurement results

Important terms

Subjective refraction	Values which are the basis for treatment calculation. Users can manually adjust the subjective refraction values for sphere, cylinder, and axis.
Predicted Phoropter Refraction (PPR)	To translate the wavefront deformation into numbers which are commonly used in the refractive surgery (manifest refraction), a software module analyzes the shape of the wavefront and determines the best spherocylindrical correction to optimize the retinal image, taking into account all higher aberrations. These derived spherocylindrical correction values are called 'Predicted Phoropter Refraction' (PPR).
Zernike RMS	The Zernike polynomials are a set of functions that are orthogonal over the unit circle. Each polynomial represents a particular type of optical aberration. Each wavefront can be described by coefficients which, when taken as a whole, reconstruct the wavefront, while individually they describe the relative amount of each aberration. Root-mean-square (RMS) is the square root of the mean of the squares (quadratic mean) of the Zernike coefficients that give the power of the aberrations for the terms that have been summed. Total RMS is the overall magnitude of all the eye's refractive errors (sphere, cylinder and HOAs).
Higher Order Aberrations (HOAs)	Higher-order aberrations are optical irregularities or imperfections of the eye which cannot be corrected by simple sphere and cylinder corrections. Some of the most important are coma and trefoil (3rd order), spherical aberration and quadrafoil (4th order), and pentafoil (5th order).
Point Spread Function (PSF)	The point spread function or impulse response is the image of a point object formed by the optical system. PSF is an image quality function representing the quality of the images rendered by the optical system.

The following ZYWAVE 3 display options of measurement results are available:

Summary

Fig. 5: 'Summary'

The 'Summary' (Fig. 5) display option outlines the measurement results in form of maps for the wavefront, higher order wavefront and higher order PSF and offers a variety of additional display options which can be viewed by clicking on the respective tab.

Furthermore, the 'Subjective refraction' values, the PPR values and the 'Zernike RMS' are displayed in the lower display area.

For 'Wavefront' and 'Higher order wavefront', you can click 'Color legend'  to change the 'Number of colors' and the 'Color scale type' of the map.

The statistical data on the lower left side of Fig. 5 gives a summary of important parameters.

Location	Description
Examination data	■ Examined eye (OD, OS) ■ Diagnosis (as entered by the user) ■ Dilation status (as entered by the user) ■ Wavefront diameter ($\approx$ diameter of the wavefront information used to plan the surgery)* ■ Exam. diameter (= patient's pupil diameter during wavefront measurement) ■ Database value from the pupilometry (dim or dark light conditions) * For measurement of hyperopic or mixed astigmatism cases the wavefront diameter corresponds to the exam. diameter. For myopic cases and exam. diameters less than 7.7 mm, the wavefront diameter is 10% larger than the exam. diameter. The wavefront diameter is limited to a maximum of 8.5 mm. The refractive indication is determined by the sign of sphere and cylinder of the PPR.
Refractions (sphere, cylinder, axis)	■ Subjective refraction (as entered by the user) ■ PPR values for 3.5 mm ■ Difference values for limits and tolerances based on comparison between 'Subjective refraction' and PPR (3.5 mm)** ■ PPR for exam. diameter ** The difference is based on the minus cylinder convention.
Zernike Root Mean Square	■ Higher Order Zernike Root Mean Square (HORMS) for 5 mm and 6 mm*** ■ HORMS without Zernike coefficient Z400 (spherical aberration) for 5 mm and 6 mm*** ■ Total Zernike RMS for 5 mm and 6 mm*** *** These values are only displayed if the wavefront diameter was larger than 5 mm / 6 mm.

If one of the difference values is given in red, this indicates one of the following limits and tolerances has been exceeded, i.e. the ZYWAVE 3 Predicted Phoropter Refraction at an analysis diameter of 3.5 mm (PPR 3.5 mm) exhibits a difference compared to the entered subjective refraction:

- of greater than 0.75 D in sphere power,
- or a difference of greater than 0.50 D in cylinder power,
- or a difference in cylinder axis of more than 15 degrees when the cylinder is > 0.50 D.

2D Plot

Fig. 6: '2D Plot'

The '2D Plot' (Fig. 6) displays the wavefront map in μm . The values indicated in the left lower corner of the diagram constantly change depending on where you click on the map. They show the cursor position (angle and distance) and the value of the respective spot that the mouse cursor is pointing to.

In the '2D Plot' display option, can be selected

- 'Higher order only',
- 'All', or
- 'None', 'Zernike Polynomials'.

i
i
i

To select the individual 'Zernike Polynomials', click on them.

Clicking on a 'Zernike Polynomial' deactivates the respective Zernike coefficient and changes the corresponding map in such a way that the respective coefficients are no longer included.

The change in the Zernike polynomials display remains active for all wavefront displays of the selected exam, except for the 'Summary' screen. The pyramid icon indicates if and which polynomials have been deactivated. The deactivated polynomials are gray, as shown in the example. In this example all Zernikes are activated (Fig. 6).

PPR vs. Pup. diam.


Fig. 7: 'PPR vs. Pup.diam.'

The 'PPR vs. pup. diam.' (Fig. 7) displays the PPR sphere (blue), cylinder (purple) and axis (green) for different diameters.

To display the corresponding values (sphere/cylinder/axis/pupil diameter) for the diameter that the mouse is pointing to, click with the mouse on any point in the diagram.

Higher order PSF

Fig. 8: 'Higher order PSF'

The 'Higher order PSF' (Fig. 8) displays the simulated PSF (point spread function), but only for the higher order aberrations.

On the left side of the figure, you can select

- 'Higher order', or
- 'None' of the 'Zernike Polynomials'.

To select the individual 'Zernike Polynomials', click on them.

Clicking on the 'Zernike Polynomials' changes the display option in such a way that the respective coefficients are no longer included.

Raw data

Fig. 9: 'Raw data'

The 'Raw data' display option (Fig. 9) enables the user to view the raw data of the selected exam. It shows the three pupil and centroid images that were used for calculating the final wavefront.

The Exam diameter for each picture is given in the table right, next to the images. The table also indicates the three best images out of the five which were taken.

Click an image to enlarge it for a better, more detailed view (↪ Chapter 9.3.2 'Verifying the measurement images' on page 9 - 21).

Normal band

Fig. 10: 'Normal band'

The 'Normal band' (Fig. 10) displays the magnitude of all 2nd to 5th order Zernike coefficients of the measured eye in μm . 'Normal Band' parameters can be displayed for 'PreOp' and 'PostOp' values, the user having the option to select data for 'PreOp' or 'PostOp' eyes and between 5 mm or 6 mm pupil size.

'PreOp' or 'PostOp' is preselected in the 'Normal band' screen, depending on whether you selected 'PreOp' or 'PostOp' in the 'Diagnosis' field. If the status 'PreOp/PostOp' was not automatically detected, no display option is preselected in the 'Normal band' screen but must be selected manually by clicking on the respective check mark.

You can increase or decrease the display by clicking 'Zoom in'  or 'Zoom out' .

Difference


Fig. 11: 'Difference'

The 'Difference' display option (Fig. 11) enables the user to subtract one wavefront examination (or wavefront) from another one. This is a tool for evaluating the development of the patient's eye over time or to compare 'PreOp' and 'PostOp' measurements.

For information about how to compare measurements, refer to  Chapter 9.6.3.1 'Comparing ORBSCAN 3 measurements' on page 9 - 29 and  Chapter 9.6.3.2 'Comparing ZYWAVE 3 measurements' on page 9 - 31.

6.3 Treatment Calculator

6.3.1 Clinical usage

The integrated Treatment Calculator generates an individual treatment data for the patient on the basis of the available diagnostic data from the ORBSCAN 3 and ZYWAVE 3 devices.

The Treatment Calculator uses the available diagnostic data to calculate them. Finally, these calculated ablation patterns are used by TECHNOLAS® TENEO™ 317.

6.3.2 Available treatment types and functionality

The preferred choice of treatment type depends on a variety of factors. Contact your Authorized Technolas Perfect Vision Representative if you need more information on the benefits and features of each treatment type.

The TECHNOLAS® TENEO™ 317 laser system supports various treatment types. They can be classified by the treatment algorithms that are applied for the calculation of the treatment ablation profile.

Depending on the configuration of the laser system, the availability of treatment types may vary. Where applicable, the potential differences in function and display in the laser software are pointed out in the description.

The availability of the SUPRACOR treatment depends on the regulatory approval status in your country. Please contact your Authorized Technolas Perfect Vision Representative for further information.

	PROSCAN	ZYOPTIX HD	SUPRACOR	PTK
Repetition rate	500 Hz	500 Hz	500 Hz	500 Hz
Iris recognition	✓	✓	✓	✓
Z-eye tracking	✓	✓	✓	✓
XY-eye tracking	✓	✓	✓	✓
Myopia	✓	✓	✓	×
Hyperopia	✓	✓	✓	×
Mixed astigmatism	✓	✓	×	×

✓ = available

× = not available

The following Treatment Calculator display options are available:


Fig. 12: Treatment Calculator - display options

PROSCAN treatment

For PROSCAN, the treatment profile calculation is based on the data for treatment refraction and optical zone with an improved asphericity algorithm. It is also possible to customize the treatment by entering the K- and Q-values during the treatment calculation. The K-value corresponds to a keratometric value in diopters fitted to the central area of the cornea to be treated. The Q-value is a measure for the corneal asphericity based on topographical measurement data obtained from ORBSCAN 3.

Iris data for iris recognition at the laser system is transferred as well.

ZYOPTIX HD treatment

ZYOPTIX HD is a wavefront treatment option (all refractive ranges) which considers the compensation of pre-operative higher and lower wavefront irregularities and the reduction of surgically induced spherical aberrations.

The advanced nomogram is automatically included in the treatment option.

SUPRACOR treatment

SUPRACOR is the Excimer-based treatment for presbyopic patients. SUPRACOR consists of a pre-calculated pulse pattern that is mixed into the refractive pulse pattern to combine distance and near vision correction.

Iris data for iris registration at the laser system is transferred as well.

PTK treatment

The PTK mode treats a uniform ablation depth within a predefined optical zone and is used for specific conditions or diseases of the corneal surface.

6.3.3 Inclusion and exclusion criteria

If the ZDW3 diagnostic software is used together with the TECHNOLAS® TENEO™ 317 apply the treatment inclusion and exclusion criteria of the TECHNOLAS® TENEO™ 317 laser system (refer to the TECHNOLAS® TENEO™ 317 user manual for more details).

7 Hardware

7.1 Overview of the main components

This section provides information about the main hardware components (Fig. 13).


Fig. 13: Components overview ZDW 3

- 1 Center table
- 2-3 Lift tables (identical for ORBSCAN 3 and ZYWAVE 3)
- 4-5 Positioners with joysticks (identical for ORBSCAN 3 and ZYWAVE 3)
- 6 ZYWAVE 3 measurement head
- 7 ORBSCAN 3 measurement head
- 8 User control screen

7.1.1 Center table

The center table (Fig. 14) houses the user control screen and serves as a rest table for the keyboard and mouse.


Fig. 14: Center table

Standby switch

The standby switch is located on the center table (Fig. 14/1) and allows to switch on/off the diagnostic device with exception of the user control screen.

7.1.2 Lift tables


Fig. 15: Lift table with up/down buttons

- 1 lift table
- 2 up/down buttons

The lift tables (Fig. 15/1) enable the user to move the ORBSCAN 3 and ZYWAVE 3 measurement heads in z-axis (vertically) with the help of the up/down buttons (Fig. 15/2).

The ORBSCAN 3 and ZYWAVE 3 measurement heads (☞ Chapter 7.1.3 'ORBSCAN 3 and ZYWAVE 3 measurement heads' on page 7 - 4) are each mounted on a lift table which fulfils the same functionalities.

7.1.3 ORBSCAN 3 and ZYWAVE 3 measurement heads

ORBSCAN 3 measurement head

The ORBSCAN 3 measurement head (Fig. 16) consists of the following components:


Fig. 16: ORBSCAN 3 measurement head components - patient's perspective

- 1 Head rest
- 2 Optical components (inside the housing)
- 3 Chin rest adjustment handhold
- 4 Joystick
- 5 Patient hand holds (left and right)
- 6 Chin rest
- 7 Placido ring with aperture

ZYWAVE 3 measurement head

The ZYWAVE 3 measurement head (Fig. 17) consists of the following components:


Fig. 17: ZYWAVE 3 measurement head components - patient's perspective

- 1 ZYWAVE 3 measurement head aperture
- 2 Head rest
- 3 Chin rest
- 4 Patient hand holds (left and right)
- 5 Joystick
- 6 Chin rest adjustment handhold
- 7 Optical components (inside the housing)

7.1.3.1 Positioner with joystick

The ORBSCAN 3 and ZYWAVE 3 measurement heads each have a positioner and a joystick.

The ORBSCAN 3 and ZYWAVE 3 measurement heads are each mounted on a positioner, which provides a mechanical system (joystick [Fig. 18/1] with trigger button [Fig. 18/2]) for movements in the x-, y- and z-axis.

The movements in the x-, y-axis, referred to as alignment, is triggered by a joystick which controls not only the fine adjustments in the x/y-axis, but also the positioning and triggering of measurements by means of a trigger button.

The height of the measurement head is adjusted by turning the joystick clockwise or counterclockwise (adjustment in z-axis).


Fig. 18: Joystick with trigger button

7.1.4 User control screen

The user control screen (Fig. 19) displays all relevant information about patients, measurement data as well as measurement progress, and allows the user to make all necessary data input to operate and control the diagnostic device.


Fig. 19: User control screen (front view)

- 1 [On/Off] switch to switch the user control screen on or off.
- 2 Buttons to adjust settings such as brightness and read light.

For more information about to how to adjust the user control screen settings refer to Chapter 7.3.3 'Adjusting the user control screen settings' on page 7 - 11.

The user control screen is mounted on a pedestal that can be tilted (Fig. 20/1).


Fig. 20: User control screen (back view) - pedestal

- 1 Pedestal

7.2 Accessories

This chapter lists an overview of components, which can be connected to the diagnostic device to enhance functionalities or which are needed for regular maintenance tasks.

Calibration check tools

The diagnostic device is calibrated by our trained service personnel. However, users are required by national standards to perform regular checks of the calibration state. For this purpose, two calibration check tools are included in the delivery (Chapter 1.2 'Scope of delivery' on page 1 - 3).

ORBSCAN 3 Calibration check tool

ZYWAVE 3 Calibration check tool

Human interface devices

User control screen, keyboard and mouse

The diagnostic device is equipped with a user control screen (touch screen), which provides the graphical user interface of the software and which can be used to perform all necessary data input to operate the diagnostic device. However, it is preferred to enter data by an external keyboard and mouse (Fig. 21).


Fig. 21: User control screen, keyboard and mouse

NOTICE

Do not try to operate the touch screen with sharp objects (pens, instrument holders, tools). The touch screen is designed to be used by finger pressure alone and may become damaged.

DANGER

The USB ports are electrical insulated from the ZDW 3 application unit according to IEC 60601-1, 3rd ed. However, the user is responsible for any device plugged into the ports. Defective devices may cause malfunctions and/or pose safety risks to the user.

USB memory stick

NOTICE

Danger of property damage due to use of unauthorized accessories / devices!

Only USB memory sticks delivered with this diagnostic device may be used.

- **Do not use other devices than the one delivered with this diagnostic device.**

Other devices will not work in this diagnostic device.

Only USB memory sticks delivered with this diagnostic device can be used to export measurement data.

The USB memory stick will be inserted into the USB port on the all-in-one PC right side.

Fig. 22: USB memory stick

Network connection

The diagnostic device may be connected to a local area network (LAN) via a standard RJ45 CAT 6 or better cable. The actual configuration of the diagnostic device depends on the network architecture. The installation of network functionality can be executed by trained service personnel only. Contact an Authorized Service Technician for installation.

Database Server Unit

The database is installed on an external Database Server Unit (Fig. 23) which hosts patient data, measuring data and treatment data. The external Database Server Unit is connected to the diagnostic device ('Server mode') and, optionally, to other laser systems.

The external Database Server Unit can also host ZYWAVE II and ORBSCAN IIz measuring data which can be selected and displayed in the ZDW3.

The treatment and measurement data from the diagnostic device are stored on the external Database Server Unit and can be automatically loaded from connected TECHNOLAS® TENEO™ 317 laser systems for treatment of patient's eyes.

If the connection to the external Database Server Unit ('Server mode') is interrupted, the user can switch from the 'Server mode' to the 'Local mode' (connection to the local database on the diagnostic device). In the 'Local mode', the user can manually export or import measurement data using a USB memory stick delivered with the diagnostic device.

The external Database Server Unit is not synchronized with the local database.


Fig. 23: Database Server Unit - front and back view

7.3 Operating the hardware

This section provides information about how to operate the hardware of the diagnostic device.

Switching on**7.3.1 Switching on and off**

1. ➤ Turn the standby switch of the diagnostic device to the [On] position.
2. ➤ Turn on the user control screen.

The diagnostic device is equipped with a user control screen (touch screen), which provides the graphical user interface of the software and which can be used to perform all necessary data input to operate the diagnostic device. However, it is preferred to enter data by an external keyboard and mouse.

⇒ The login page opens.



Fig. 24: Login page

Switching off

1. ➤ Click the **Off** button (1) in the login page to close the ZDW 3 software.
2. ➤ Turn the standby switch of the diagnostic device to the [Off] position.

⇒ The diagnostic device shuts down instantly.

Positioning and triggering the measurement**7.3.2 Positioning and triggering measurements**

For positioning and triggering the measurement, a joystick with trigger button is provided (☞ Chapter 7.1.3.1 'Positioner with joystick' on page 7 - 5).

The height of the measurement head is adjusted by turning the joystick clockwise or counterclockwise (adjustment in z-axis).

For information about the joystick, refer to ☞ Chapter 7.1.3.1 'Positioner with joystick' on page 7 - 5.

For information about aligning the patient and triggering the measurement, refer to:

- Chapter 9.2.1 'Aligning the patient and starting the ORBSCAN 3 measurement' on page 9 - 3
- Chapter 9.3.1 'Aligning the patient and starting the ZYWAVE 3 measurement' on page 9 - 14

7.3.3 Adjusting the user control screen settings

Switching the user control screen on

➔ Switch the user control screen on via the [On] button (Fig. 25/1).

User control screen buttons (front view)


Fig. 25: User control screen (front view)

- 1 On (green light)
- 2 Volume down
- 3 Volume up
- 4 Touch On / Off
- 5 Read light control
- 6 Brightness decrease
- 7 Brightness increase

7.3.4 Performing system checks

Proceed as follows to perform system checks:

- To perform calibration checks, refer to Chapter 8.5.3 'Performing ORBSCAN 3 calibration checks (alignment PMMA)' on page 8 - 33 and to Chapter 8.5.4 'Performing ZYWAVE 3 calibration checks' on page 8 - 38.

7.4 Cleaning instructions

This section provides cleaning instructions for the diagnostic device.

NOTICE

Danger of property damage due to improper cleaning agents!

The product may be damaged if it is cleaned with improper cleaning agents or scratching cleaning materials, such as cloths.

- Do not use solvents or lubricants.
- Do not use ether, acetone, concentrated acids, or bases.
- Do not use abrasive cleaning agents or cleaning agents with ammonia or other possibly caustic products.
- Do not use any steel wool or metal brushes.
- Before using a cleaning agent, test it on an inconspicuous area to check if the cleaning agent damages the material.

NOTICE

Danger of property damage due to aggressive disinfecting agents!

The paint or the surface of the product may be damaged if it is cleaned with improper disinfecting agents.

- Use disinfectants recommended by local institutes for hygiene or microbiology (for example, Meliseptol by Braun).
- You can use any standard neutral to slightly alkaline cleaning agent according to your disinfectant procedures. The maximum concentrations are:
 - For rubbing alcohol (isopropyl alcohol): 70% and 100%
 - For alcohol: 60%
- Use only soft fuzz-free cleaning cloths

When you are not using the diagnostic device, make sure it is powered down and a dust cover is placed over each measurement head to prevent excess dust accumulation.

Periodic cleaning and disinfecting can be done easily by wiping down the surfaces that come in contact with the patient.

What cleaning agents and disinfectants to use

Different cleaning agents must be used on different parts of the diagnostic device. You will find an overview of different cleaning agents in the following table.

Part of the diagnostic device	Cleaning agent
All matte surfaces	All neutral or slightly alkaline cleaning agents
All lacquered surfaces	All neutral or slightly alkaline cleaning agents
Screen of the user control screen	Commercially available cleaning kits for computers
Keyboard and mouse at the user control screen	Dry fuzz-free cleaning cloth
ZYWAVE 3 measurement head aperture	Only soft fuzz-free cleaning cloths Note: Wipe it very softly.
ORBSCAN 3 measurement head aperture with Placido ring	NOTICE Property damage due to cleaning of the ring pattern! Cleaning the ring pattern on the ORBSCAN 3 measurement head aperture will damage the pattern and moisture could enter the device. Do not clean the ring pattern on the ORBSCAN 3 measurement head aperture.

Cleaning the diagnostic device

Precondition:

- The diagnostic device is switched off.
- 1. ➤ Ensure that no moisture enters the device.
- 2. ➤ Dilute the cleaning agent with water.
- 3. ➤ Soak a cleaning cloth in the diluted solution.
- 4. ➤ Wring out the cleaning cloth very well.
- 5. ➤ Wipe the diagnostic device with the cleaning cloth.
- 6. ➤ Wipe the diagnostic device with a dry cloth.

Cleaning the user control screen

Use commercially available cleaning kits for computers to clean the user control screen (screen and keyboard).

1. ➤ Switch off the diagnostic device.
2. ➤ Use a soft dry cleaning cloth to wipe the user control screen and the keyboard.

3. **NOTICE!** Danger of property damage due to aggressive cleaning agents and moisture!

The user control screen is an electronic system and made of plastics.

- Do not use cleaning agents or abrasive powders that dissolve plastics.
- Ensure that no moisture can enter into the user control screen or the screen monitor.
- Ensure that the ventilation slots of the screen monitor remain open.

Follow the instructions on the cleaning kit to clean the user control screen.

8 Software

Overview

This chapter provides information about the software operations available for the diagnostic software.

The diagnostic device comprises the following features:

- Performing of ORBSCAN 3 measurements and viewing and comparing of measurement results
- Performing of ZYWAVE 3 measurements and viewing and comparing of display options of measurement results
- Performing of treatment calculations based on diagnostic data


8.1 Working with the software

You operate the diagnostic software by means of the user control screen (touch screen) and, optionally, using the delivered keyboard and mouse.

Only enabled buttons open the corresponding dialog boxes and only enabled fields can be edited!

GUI element

Function

Indicates possible states of a button: enabled — disabled



The green border allows to click the button.

The green background indicates an active button.

Allows to click the corresponding field and enter the input.

Indicates a confirmation, an information or an error message.

8 SOFTWARE

GUI element	Function
	
	

8.1.1 Logging in

To log in into the diagnostic software, proceed as follows.

Precondition:

- The diagnostic device and the user control screen are switched on.
- The login page is open displaying the installed software version, the login data and the database connection (the 'Server mode' is set as default).

Fig. 26: Database connection

Switching database modes

If the connection to the external database server is interrupted, you can switch from the remote connection (external database) to the local connection (local database on the diagnostic device).

To switch the database modes, click 'Local mode' . The connection to the local database is established.

The local database is not synchronized with the external database server.

Logging in

Fig. 27: Software version and login data

1. Enter your 'Login name' and 'Password' (Fig. 27).
2. Click 'Log in' (Fig. 27).
⇒ The menu page opens (Fig. 28).


Fig. 28: Menu page

- 1 'Calculation settings'
- 2 'Patients'
- 3 'View data on USB stick'
- 4 'Settings'

8.1.2 Logging out

→ In the upper-right corner of every GUI page, login page excluded, click 'Log out' to log out from the diagnostic software.

⇒ The login page opens.

8.1.3 Changing the user

1. → If another user is already logged in, click 'Log out'.

⇒ The login page opens.

2. → Enter another user name and password to change the user.

3. → Click 'Log in'.

8.2 Calculation settings

1. → After your personal login, click 'Calculation settings' in the menu page.

⇒ The 'Calculation settings' dialog box opens displaying:

- Default or user-defined 'Nomogram' values for PROSCAN, ZYOPTIX HD and SUPRACOR treatments (Fig. 29/1),
- 'Nomogram comments' (Fig. 29/2), and
- Default treatment values for 'Standard optical zone diameter [mm]', 'SUPRACOR optical zone diameter [mm]', and 'Flap thickness [μm]' (Fig. 29/3).

Nomogram				Defaults	
	Myopia Base sphere (D)	Hyperopia Base sphere (D)	Mixed astig. Base sphere (D)		
PROSCAN	100	100	100	Standard optical zone diameter [mm]	6.5
ZYOPTIX HD	100	100	100	SUPRACOR optical zone diameter [mm]	6.0
SUPRACOR		100		Flap thickness [μm]	120
SUPRACOR		0.00		Nomogram Comments	
				No modifications	
				✓ Accept	✗ Cancel

1
2
3

Fig. 29: 'Calculation settings' dialog box

2. ➤ Enter the 'Nomogram' and the default treatment values.

3. ➤ If required, click 'Cancel' to reset the changes.

⇒ The menu page opens.

If you want to set calculation settings, repeat step 1 and 2 to continue.

4. ➤ Click 'Accept' to apply the changes.

5. ➤ Click 'Menu' to go back to the menu page.

⇒ The menu page opens.

8.2.1 About the user-defined nomogram

WARNING

Danger of patient injury due to inappropriate nomogram!

The nomogram option is a tool to assist the user in planning a surgical treatment. The user must determine, in each instance, whether the nomogram is appropriate for use in a treatment plan for a particular patient, taking into consideration all available information. The use of incorrect or unevaluated nomogram values may cause serious permanent patient injury. Failure to verify any treatment parameter may result in serious permanent patient injury.

- Ensure you are using correct and verified nomogram values.
- Verify all calculation results before incorporating them into any treatment plan.
- If the nomogram is determined to be inappropriate for incorporation into a treatment plan, consider proceeding to treat the patient without using the nomogram adjustment.

The nomogram option is available as an additional feature for the calculation of PROSCAN treatment, ZYOPTIX HD treatment, and SUPRACOR treatment. A user-defined nomogram may be applied with the goal to improve refractive outcome by adjusting the sphere or the spherical equivalent.

Adjusting calculated refraction

The calculated treatment sphere can be adjusted by a user-defined nomogram with the goal of optimizing the final result. The user-entered nomogram values must be based on the experience and analysis of previous outcomes, regarding the related treatment parameters and the treatment outcome.

The user-defined nomogram is applied to all further treatment calculations, depending on:

- The selected treatment mode,
- The target value,
- The refraction of the PPR for ZYOPTIX HD treatment, or
- The subjective refraction for PROSCAN treatment and SUPRACOR treatment.

For information on how to set an individual nomogram, refer to  Chapter 8.2.2 'Setting a user-defined nomogram' on page 8 - 6.

⚠ WARNING**8.2.2 Setting a user-defined nomogram**

Danger of patient injury due to inadequate Base sphere [%] or Sphere offset [D] values!

Use of an inadequate base sphere [%] or sphere offset [D] values may result in serious permanent patient injury.

- **Always verify the applied nomogram parameters with each treatment calculation.**

Contact an Authorized Service Technician for assistance on initial nomogram settings.

Nomogram adjustment

Unless the user adjusts the nomogram values (§ Chapter 8.2 'Calculation settings' on page 8 - 4), the diagnostic software applies the default nomogram values. The default values for 'Base sphere [%]' are 100%. The 'Sphere offset [D]' for the SUPRACOR treatment is 0.00 D.

In order to bring the actual surgery outcome closer to the targeted outcome, users may modify the default nomogram values according to an analysis of their previous treatments.

For PROSCAN treatment and HD treatment, the user may define a linear adjustment to the laser input sphere value ('Base sphere [%]'). The allowed range for the 'Base sphere [%]' is from 70% to 130%.

For SUPRACOR treatment, it is possible to define an additional nomogram value for 'Sphere offset [D]'. The allowed range for 'Sphere offset [D]' is from -1.50 D to +1.50 D. The allowed range for the 'Base [%]' is from 50% to 200%.

The user-defined nomogram values are applied to adjust the treatment values (correction sphere) based on the subjective refraction values.

Setting a nomogram

For PROSCAN and SUPRACOR treatment, a nomogram for myopia is always applied in minus cylinder convention, for hyperopia in plus cylinder convention. In the HD treatment mode the specific 'Nomogram' is applied.

1. In the 'Nomogram' table, set the 'Base sphere [%]', 'Base [%]', and 'Sphere offset [D]' values as required.

PROSCAN treatment	HD treatment	SUPRACOR treatment
Base sphere for:	Base sphere for:	Base for:
■ Myopia in %	■ Myopia in %	■ Hyperopia in %
■ Hyperopia in %	■ Hyperopia in %	Sphere offset for:
■ Mixed astigmatism in %	■ Mixed astigmatism in %	■ Hyperopia in D

⇒ Invalid nomogram values are highlighted (framed) in red.

2. Enter your 'Nomogram comments' as needed.

The defined nomogram values are saved per user according to the entered user name and password.

3. ➔ Set default treatment values (§ Chapter 8.2.3 'Setting default treatment values' on page 8 - 7) or click 'Accept'.

⇒ The menu page opens.

4. ➔ Click 'Cancel' to cancel your modifications.

⇒ The menu page opens.

8.2.3 Setting default treatment values

The user has the option to set default values for the optical zone diameter and flap thickness. The defined values are automatically applied as **default values** to all further treatment calculations.

WARNING

Danger of patient injury due to inappropriate flap thickness!

Failure to utilize an appropriate flap thickness may result in possible serious permanent patient injury.

- **Carefully verify all adjustments to the flap thickness for each individual case before calculating a treatment file.**

Setting default values

1. ➔ Modify the default value in the 'Standard optical zone diameter [mm]' field.

Default value for the optical zone diameter is: 6.50 mm. The allowed range is from 4.50 mm to 8.50 mm in steps of 0.1 mm.

2. ➔ Modify the default value in the 'SUPRACOR optical zone diameter [mm]' field.

Default value for the optical zone diameter is: 6.00 mm. The allowed range is from 4.5 mm to 8.5 mm in steps of 0.1 mm.

3. ➔ Modify the default value in the 'Flap thickness [µm]' field.

Default value for flap thickness is: 120 µm. The possible range is from 0.00 µm to 300.00 µm in steps of 1 µm.

4. ➔ Click 'Accept'.

⇒ The menu page opens.

5. ➔ Click 'Cancel' to cancel your modifications.

The default values are saved per user according to the entered user name and password.

The defined values are automatically applied as **default values** to all further treatment calculations.

⇒ The menu page opens.

8.3 Patient administration

8.3.1 Selecting patients

To select patients, proceed as follows:

1. ➤ In the menu page, click 'Patients'.
⇒ The 'Select patient' dialog box opens (Fig. 30).


Fig. 30: 'Select patient' dialog box - patient list (example)

i

In the patient list the user can scroll for the last 50 patients which have been processed by the user.

2. ➤ In the 'Select patient' field, enter the first letters of the patient's first or last name, or the patient ID.

⇒ The patient is displayed in the patient list.

Additionally, the number of the patients on the data base and the number of patients have been found according to the entered search criteria is displayed on the top right of the 'Select patient' dialog box.

i

If the patient is not displayed, you need to add the patient first (Chapter 8.3.3 'Adding patients' on page 8 - 15).

3. ➤ To select a patient, double-click the patient's entry.

⇒ The selected patient account opens.


The screenshot shows a patient account interface for 'SMITH John'. It includes fields for patient data (First name, Last name, Date of birth, Patient ID, Gender), subjective refraction (Sphere [D], Cylinder [D], Axis [°], Diagnosis, Dilation), and sections for 'Select/Add measurement' and 'Select/Add treatment'. Numbered callouts indicate: 1. Patient information header; 2. 'Patient data' entry fields; 3. 'Subjective refraction' values; 4. 'Select/Add treatment' area; 5. 'Select/Add measurement' area.

Fig. 31: Patient account (example)

- 1 Patient information
- 2 'Patient data' entry fields
- 3 'Subjective refraction' values
- 4 'Select/Add treatment' area
- 5 'Select/Add measurement' area

If the patient account has already been opened by another user, the patient data is not editable (Fig. 31/2).

4. ➔ To go back to the patient list, click 'Back'.

⇒ The patient list opens.

After selecting a patient, you can:

- Edit patient data & Chapter 8.3.4 'Editing patient data' on page 8 - 18
- Enter 'Subjective refraction' values
- Select or add measurements
 - & Chapter 9.2 'Performing ORBSCAN 3 measurements' on page 9 - 2
 - & Chapter 9.3 'Performing ZYWAVE 3 measurements' on page 9 - 14
 - & Chapter 9.6.1 'Selecting and closing measurements' on page 9 - 27
- Select or add treatments
 - & Chapter 9.6.8.1 'Selecting treatments and diagnostic data' on page 9 - 40 and & Chapter 9.6.8.2 'Adding treatments' on page 9 - 45

8.3.2 Searching for patient data

To perform an advanced search for patient data, proceed as follows:

Precondition:

The 'Select patient' dialog box is open (Fig. 30).

Advanced search

1. In the 'Select patient' dialog box, click 'Advanced search'.

⇒ The 'Advanced search' dialog box opens (Fig. 32).

In the patient list all existing patients on the data base are displayed as default.

Fig. 32: 'Advanced search' dialog box

- 1 Reset search criteria
- 2 Patients total
- 3 Patients found according to the user entered search criteria
- 4 Treatment data
- 5 Measurement data
- 6 Patient data
- 7 Patient list

Now you can perform an advanced search by entering:

- 'Patient data'
- 'Measurement data'
- 'Treatment data'


Patient data

Fig. 33: 'Patient data' area

2. In the 'Patient data' entry fields, enter the first letters of the patient's first or last name, or the patient ID.
3. If necessary, further search criteria can be defined to narrow down the search results.

Search criteria

Born after

Born before

Gender

Subjective refraction

Sphere (from ... to ...)

Cylinder (from ... to ...)

Description

Defines the period in which the date of birth of the patient may lie.

Defines the gender.

Defines the range in which the values may lie.

Action

Enter earliest and latest date of searched date of birth by selecting it in the calendar and clicking 'Set' to confirm.

Select 'Male' or 'Female' or both.

Enter lowest and highest value of searched sphere.

The order of the entered data is not relevant.

Enter lowest and highest value of searched cylinder.

The used cylinder conventions are not relevant.

Refraction types

Emmetropia

Myopia

Hyperopia

Mixed Astigmatism

Enables a search for patients with at least one of the listed refraction types.


Select one or more of the displayed refraction types.

Measurement data

4. To enable the search criteria of this area select 'Measurement data' in the upper right corner (Fig. 34).

Fig. 34: 'Measurement data' area

Search criteria	Description	Action
Device used for examination		
ZYWAVE	Enables a search for patients with measurements performed by ZYWAVE or ORBSCAN.	Select 'ZYWAVE' to search for patients with measurements performed by ZYWAVE.
ORBSCAN		Select 'ORBSCAN' to search for patients with measurements performed by ORBSCAN.
		Select 'ZYWAVE' and 'ORBSCAN' to search for patients with measurements performed by both devices.
OD	Enables a search for patients with measurements of OD or OS or both.	Select 'OD' to search for patients with OD measurements.
OS		Select 'OS' to search for patients with OS measurements.
		Select 'OD' and 'OS' to search for patients with measurements of at least one of both eyes.
ODS	Enables a search for patients with measurements of OD and OS.	Select 'ODS' to search for patients with measurements for both eyes.
Performed after	Enables to define the period in which the measurement may lie.	Define earliest and latest date of searched measurement data by selecting it in the calendar and clicking 'Set' to confirm.
Performed before		
Device used for measurement	Enables a search for measurements performed on a specific device.	Select 'Only measurements performed with following device:'. Click in the entry field to open the drop-down list. Now you can select one or more devices from the drop-down list.

Treatment data

5. ➤ To enable the search criteria of this area select 'Treatment data' in the upper right corner (Fig. 35).

Fig. 35: 'Treatment data' area

Search criteria

Calculated after
Calculated before

Description

Defines the period of the treatment calculation.

Action

Define earliest and latest date of searched treatment calculations by selecting it in the available calendar (for 'Calculated after' and 'Calculated before') and clicking 'Set' to confirm.

Only treatments that have been started

Started after
Started before

Defines the period in which the started treatment may lie.

Select 'Only treatments that have been started' to enable the entry fields 'Started after' and 'Started before'.

Define earliest and latest date of searched treatment by selecting it in the calendar and clicking 'Set' to confirm.

Patients with at least one treatment for

OD
OS

Enables a search for patients with treatments of OD or OS or both.

Select 'OD' to search for patients with OD treatments.

Select 'OS' to search for patients with OS treatments.

Select 'OD' and 'OS' to search for patients with treatments of at least one of both eyes.

ODS

Enables a search for patients with treatments of OD and OS.

Select 'ODS' to search for patients with treatments for both eyes.

⇒ The search results are displayed in the lower area of the screen.


Additionally, the number of the patients on the data base and the number of patients have been found according to the entered search criteria is displayed on the top right of the advanced search dialog box.

6. ➤ To select a patient, double-click the patient's entry.

⇒ The selected patient account opens.

If the patient account has already been opened by another user, the patient data is not editable.

7. ➤ To go back to the patient list, click 'Back'.

⇒ The patient list opens.

After selecting a patient, you can:

- Edit patient data & Chapter 8.3.4 'Editing patient data' on page 8 - 18
- Enter 'Subjective refraction' values
- Select or add measurements
 - & Chapter 9.2 'Performing ORBSCAN 3 measurements' on page 9 - 2
 - & Chapter 9.3 'Performing ZYWAVE 3 measurements' on page 9 - 14
 - & Chapter 9.6.1 'Selecting and closing measurements' on page 9 - 27
- Select or add treatments
 - & Chapter 9.6.8.1 'Selecting treatments and diagnostic data' on page 9 - 40 and & Chapter 9.6.8.2 'Adding treatments' on page 9 - 45

8.3.3 Adding patients

To add patients, proceed as follows:

1. In the menu page, click 'Patients'.
⇒ The patient list opens.


Fig. 36: Patient list (example)

2. In the upper right corner of the patient list, click 'Add patient' (Fig. 36).
⇒ A new patient account opens.

Fig. 37: New patient account (example)

General patient data

3. Enter the first and last name of the patient.
4. Click the 'Date of birth' entry field to open the calendar (Fig. 38).
⇒ A calendar will be displayed.

Fig. 38: Patient data: select date of birth

5. ➤ Select a birth date.
6. ➤ Click 'Set' to confirm or 'Cancel' to delete the selected date of birth.
7. ➤ Enter the 'Patient ID'.

The patient ID can be entered manually in the 'Patient ID' entry field or automatically by clicking the plus button.

Subjective refraction values

8. ➤ Select the patient gender.
9. ➤ Enter the required 'Subjective refraction' values (refer to 'Subjective refraction' values' Table on page 9 - 15).
10. ➤ If you want to restore the original values for the patient data, click 'Undo'.
⇒ A warning message opens (Fig. 39).


Fig. 39: Warning message - restore the original values

- Click 'Yes' to restore the original values for the patient data and then continue to enter patient data or go back to the patient list dialog by clicking the 'Menu' button.
- Click 'No' to continue.

11. ➤ Click 'Accept' to store all entered patient data in the database.
⇒ If the entry is incomplete, a warning message opens informing you which fields are still missing (Fig. 51).


Fig. 40: Missing entries warning message (example)

12. Click 'OK', enter the missing entries and repeat step 11.

After adding the new patient, you can perform a measurement directly by clicking 'ORBSCAN' or 'ZYWAVE'.

- Chapter 9.2 'Performing ORBSCAN 3 measurements' on page 9 - 2
- Chapter 9.3 'Performing ZYWAVE 3 measurements' on page 9 - 14

nd
og

8.3.4 Editing patient data

To edit patient data, proceed as follows:

1. ➤ Select a patient (☞ Chapter 8.3.1 'Selecting patients' on page 8 - 8).
⇒ The patient account opens (Fig. 41) and can be edited.

Fig. 41: Patient account (example)

2. ➤ Edit the required general 'Patient data' (☞ on page 8 - 15).
3. ➤ Additionally, you can edit 'Subjective refraction' values (refer to ☞ "Subjective refraction' values' Table on page 9 - 15).
5. ➤ If you want to restore the original values for the patient data, click 'Undo'.
⇒ A warning message opens (Fig. 39).

Fig. 42: Undo warning message

- Click 'Yes' ☒ Yes to restore the original values for the patient data and continue to enter patient data or to go back to the patient list dialog by clicking the 'Menu' button.
- Click 'No' ☒ No to continue.

6. ➤ Click 'Accept' to store the entries in the patient database.
⇒ If the entry is incomplete, a warning message opens informing you which fields are still missing (Fig. 51).


Fig. 43: Missing entries warning message (example)

7. ➔ Click 'OK', enter the missing entries and repeat step 5.

After editing the patient data, you can

- Select or add treatments (↗ Chapter 9.6.8.1 'Selecting treatments and diagnostic data' on page 9 - 40 and ↗ Chapter 9.6.8.2 'Adding treatments' on page 9 - 45)
- Perform measurements (↗ Chapter 9.2 'Performing ORBSCAN 3 measurements' on page 9 - 2 and ↗ Chapter 9.3 'Performing ZYWAVE 3 measurements' on page 9 - 14)

8.3.5 Deleting patient data

The delete option is only available if the patient you want to delete does not have any measurements or treatments saved. Otherwise, the delete button is not enabled.

To delete patient data, proceed as follows:

1. ➤ In the menu page, click 'Patients'.

⇒ The patient list opens.

2. ➤ In the patient list, select the patient you want to delete (§ Chapter 8.3.1 'Selecting patients' on page 8 - 8).

3. ➤ Click 'Delete'.

⇒ A warning message opens (Fig. 44).




Fig. 44: 'Delete patient' warning message

4. ➤ Click 'Yes' to delete the patient.

⇒ The patient is deleted.

5. ➤ Click 'No' to select another patient.

⇒ The patient list opens.

6. ➤ Repeat from step 2 to step 4.

8.4 User administration

Precondition:

- You are logged in as administrator ('ADMINISTRATOR' role).

8.4.1 User roles

The ZDW 3 GUI provides the following user roles.

- 'USER'
- 'ADMINISTRATOR'

User specific rights

The 'User' is allowed to perform all user related tasks in the diagnostic software that do not involve user administrative tasks. Thus, he can perform

- Patient administration tasks (↗ Chapter 8.3 'Patient administration' on page 8 - 8)
- Calculation settings (↗ Chapter 8.2 'Calculation settings' on page 8 - 4)
- System settings (↗ Chapter 8.5 'Settings' on page 8 - 28)
- ORBSCAN 3 and ZYWAVE 3 measurements (↗ Chapter 9 'Measurements' on page 9 - 1)
- Treatment calculations (↗ Chapter 9.6.8 'Performing treatment calculations' on page 9 - 39)

Administrator specific rights

Additionally to all user related tasks in the diagnostic software, the 'Administrator' (Admin) is allowed to perform user administration tasks.

The following restrictions apply for the administrator role:

- There is **at least one administrator** for the diagnostic software.
- The administrator cannot delete his/her own account.
- The **last** administrator cannot switch to a normal user account.

The personnel requirements are described in ↗ Chapter 2.5 'Personnel requirements' on page 2 - 3.

8.4.2 Selecting users

To select users, proceed as follows:

1. ➤ Open the menu page.
2. ➤ Click 'Settings'.
⇒ The settings page opens (Fig. 45).


Fig. 45: Settings page

3. ➤ Click 'User administration'.
⇒ The 'User accounts' dialog box opens (Fig. 46).



Fig. 46: List of user accounts (example)

4. ➤ Scroll to search for the desired user account.

5. ➔ To select a user, double-click the desired user or click 'Edit'.

⇒ The user data dialog box opens (Fig. 47). Now you can edit the user data.


Fig. 47: User data dialog box

For information about user roles, refer to Chapter 8.4.1 'User roles' on page 8 - 20.

6. ➔ Click 'Accept' to save your entries.

7. ➔ Click 'Cancel' to return to the list of user accounts.


8.4.3 Adding users

To add users, proceed as follows:

1. ➔ Open the menu page.

2. ➔ Click 'Settings'.

[illegible]

4. Click 'Add user'.

1. **Introduction**
 2. **Background**
 3. **Methodology**
 4. **Results**
 5. **Discussion**
 6. **Conclusion**
 7. **References**
 8. **Appendix**
 9. **Index**
 10. **Table of Contents**
 11. **Figure 1**
 12. **Figure 2**
 13. **Figure 3**
 14. **Figure 4**
 15. **Figure 5**
 16. **Figure 6**
 17. **Figure 7**
 18. **Figure 8**
 19. **Figure 9**
 20. **Figure 10**
 21. **Figure 11**
 22. **Figure 12**
 23. **Figure 13**
 24. **Figure 14**
 25. **Figure 15**
 26. **Figure 16**
 27. **Figure 17**
 28. **Figure 18**
 29. **Figure 19**
 30. **Figure 20**
 31. **Figure 21**
 32. **Figure 22**
 33. **Figure 23**
 34. **Figure 24**
 35. **Figure 25**
 36. **Figure 26**
 37. **Figure 27**
 38. **Figure 28**
 39. **Figure 29**
 40. **Figure 30**
 41. **Figure 31**
 42. **Figure 32**
 43. **Figure 33**
 44. **Figure 34**
 45. **Figure 35**
 46. **Figure 36**
 47. **Figure 37**
 48. **Figure 38**
 49. **Figure 39**
 50. **Figure 40**
 51. **Figure 41**
 52. **Figure 42**
 53. **Figure 43**
 54. **Figure 44**
 55. **Figure 45**
 56. **Figure 46**
 57. **Figure 47**
 58. **Figure 48**
 59. **Figure 49**
 60. **Figure 50**
 61. **Figure 51**
 62. **Figure 52**
 63. **Figure 53**
 64. **Figure 54**
 65. **Figure 55**
 66. **Figure 56**
 67. **Figure 57**
 68. **Figure 58**
 69. **Figure 59**
 70. **Figure 60**
 71. **Figure 61**
 72. **Figure 62**
 73. **Figure 63**
 74. **Figure 64**
 75. **Figure 65**
 76. **Figure 66**
 77. **Figure 67**
 78. **Figure 68**
 79. **Figure 69**
 80. **Figure 70**
 81. **Figure 71**
 82. **Figure 72**
 83. **Figure 73**
 84. **Figure 74**
 85. **Figure 75**
 86. **Figure 76**
 87. **Figure 77**
 88. **Figure 78**
 89. **Figure 79**
 90. **Figure 80**
 91. **Figure 81**
 92. **Figure 82**
 93. **Figure 83**
 94. **Figure 84**
 95. **Figure 85**
 96. **Figure 86**
 97. **Figure 87**
 98. **Figure 88**
 99. **Figure 89**
 100. **Figure 90**
 101. **Figure 91**
 102. **Figure 92**
 103. **Figure 93**
 104. **Figure 94**
 105. **Figure 95**
 106. **Figure 96**
 107. **Figure 97**
 108. **Figure 98**
 109. **Figure 99**
 110. **Figure 100**
 111. **Figure 101**
 112. **Figure 102**
 113. **Figure 103**
 114. **Figure 104**
 115. **Figure 105**
 116. **Figure 106**
 117. **Figure 107**
 118. **Figure 108**
 119. **Figure 109**
 120. **Figure 110**
 121. **Figure 111**
 122. **Figure 112**
 123. **Figure 113**
 124. **Figure 114**
 125. **Figure 115**
 126. **Figure 116**
 127. **Figure 117**
 128. **Figure 118**
 129. **Figure 119**
 130. **Figure 120**
 131. **Figure 121**
 132. **Figure 122**
 133. **Figure 123**
 134. **Figure 124**
 135. **Figure 125**
 136. **Figure 126**
 137. **Figure 127**
 138. **Figure 128**
 139. **Figure 129**
 140. **Figure 130**
 141. **Figure 131**
 142. **Figure 132**
 143. **Figure 133**
 144. **Figure 134**
 145. **Figure 135**
 146. **Figure 136**
 147. **Figure 137**
 148. **Figure 138**
 149. **Figure 139**
 150. **Figure 140**
 151. **Figure 141**
 152. **Figure 142**
 153. **Figure 143**
 154. **Figure 144**
 155. **Figure 145**
 156. **Figure 146**
 157. **Figure 147**
 158. **Figure 148**
 159. **Figure 149**
 160. **Figure 150**
 161. **Figure 151**
 162. **Figure 152**
 163. **Figure 153**
 164. **Figure 154**
 165. **Figure 155**
 166. **Figure 156**
 167. **Figure 157**
 168. **Figure 158**
 169. **Figure 159**
 170. **Figure 160**
 171. **Figure 161**
 172. **Figure 162**
 173. **Figure 163**
 174. **Figure 164**
 175. **Figure 165**
 176. **Figure 166**
 177. **Figure 167**
 178. **Figure 168**
 179. **Figure 169**
 180. **Figure 170**
 181. **Figure 171**
 182. **Figure 172**
 183. **Figure 173**
 184. **Figure 174**
 185. **Figure 175**
 186. **Figure 176**
 187. **Figure 177**
 188. **Figure 178**
 189. **Figure 179**
 190. **Figure 180**
 191. **Figure 181**
 192. **Figure 182**
 193. **Figure 183**
 194. **Figure 184**
 195. **Figure 185**
 196. **Figure 186**
 197. **Figure 187**
 198. **Figure 188**
 199. **Figure 189**
 200. **Figure 190**
 201. **Figure 191**
 202. **Figure 192**
 203. **Figure 193**
 204. **Figure 194**
 205. **Figure 195**
 206. **Figure 196**
 207. **Figure 197**
 208. **Figure 198**
 209. **Figure 199**
 210. **Figure 200**
 211. **Figure 201**
 212. **Figure 202**
 213. **Figure 203**
 214. **Figure 204**
 215. **Figure 205**
 216. **Figure 206**
 217. **Figure 207**
 218

Fig. 49: User data dialog box

5. ➤ Enter the required user data ('Personal data' and 'Login data').

The GUI language is user-dependent and can be selected from the language selection list box (Fig. 50/1).

It is possible that the software version installed in your diagnostic device do not support some of the listed GUI languages in this example.

Personal data

Title: _____

First name: John

Last name: Smith

Login data

Login name: Jules

Password: _____

Password confirmation: _____

Language: English (United States) - English (United States) 1

Deutsch - German

English - English

English (United States) - English (United States)

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User

Fig. 50: Selecting language (example)

6. ➤ Select the user role you want to assign and click 'Accept' to save the entries.

⇒ If the entry is incomplete, a warning message opens informing you which fields are still missing (Fig. 51).

Invalid entry!

Invalid entry

1. First name: Field must not be blank.

2. Login name: Field must not be blank.

3. Password missing

4. Password confirmation missing

5. Minimum password length is 4 characters.

OK

Fig. 51: Missing entries warning message (example)

7. ➤ Click 'OK' and enter the missing entries.

8. ➤ Click 'Cancel' to leave the new user account.

⇒ A list of user accounts opens. None of the new user data have been saved.

8.4.4 Editing user data

To edit user data, proceed as follows:


1. ➤ Open the menu page.

2. ➤ Click 'User administration'.

3. ➤ Select a user (☞ Chapter 8.4.2 'Selecting users' on page 8 - 22).

4. ➤ In the user dialog box, edit the required field(s).

5. ➤ Click 'Accept' to save the user data.


⇒ The user data is saved.

6. ➤ Enter the user personal and login data, select the user role you want to assign and click 'Accept' to save the entries.


⇒ If the entry is incomplete, a warning message opens informing you which fields are still missing (Fig. 52).

Fig. 52: Missing entries warning message (example)

7. ➤ Click 'OK' and enter the missing entries.

8. ➤ Click 'Cancel' to leave the new user account.

⇒ The list of user accounts opens. None of the new user data have been saved.

8.4.5 Deleting user data

To delete user data, proceed as follows:

1. ➤ Open the menu page.

2. ➤ Click 'User administration'.

⇒ The list of user accounts opens.

3. ➤ Scroll to search for the desired user account.

4. ➔ Click 'Delete'.

⇒ A warning message opens asking if you really want to delete the selected user (Fig. 53).


Fig. 53: 'Delete user' warning message

5. ➔ Click 'Yes' if you really want to delete the user account.

⇒ The list of user accounts opens. The user account is deleted from the list of user accounts.

6. ➔ Click 'No' if you do not want to delete the user account, or if you want to select another user account.

7. ➔ Continue with step 3 if you want to delete another user account.

8.5 Settings

The 'Settings' allow you to

- View general system settings (↗ Chapter 8.5.1 'Defining system settings' on page 8 - 29),
- Perform user administration tasks (↗ Chapter 8.4 'User administration' on page 8 - 20),
- Define ORBSCAN 3 map settings (↗ Chapter 8.5.2 'Defining ORBSCAN 3 map settings' on page 8 - 31),
- Perform ORBSCAN 3 calibration check (alignment PMMA) (↗ Chapter 8.5.3 'Performing ORBSCAN 3 calibration checks (alignment PMMA)' on page 8 - 33),
- Perform ZYWAVE 3 calibration check (↗ Chapter 8.5.4 'Performing ZYWAVE 3 calibration checks' on page 8 - 38),
- View the product information (↗ Chapter 8.5.5 'Viewing the product information' on page 8 - 44).

Fig. 54: 'Settings'

- 1 'System settings'
- 2 'User administration'
- 3 'ORBSCAN map settings'
- 4 'ORBSCAN calibration check'
- 5 'ZYWAVE calibration check'
- 6 'Information'

8.5.1 Defining system settings

To define system settings, proceed as follows.

1. ➤ Open the menu page.

2. ➤ Click 'Settings'.

⇒ The settings page opens (Fig. 54).

3. ➤ To view the system settings, click 'System settings'.

⇒ A dialog box opens (Fig. 55).

Only the 'Clinic' name can be defined by the user.

General settings

Date format 4/28/2014 12:00:00 AM

Number format 123.45

Clinic


Calibration check

ORBSCAN Last check 4/24/2014

ZYWAVE Last check not available

☒ Accept
 ☐ Cancel

Fig. 55: System settings dialog box

4. ➤ Click 'Accept' to save the entries or changes.

5. ➤ Click 'Cancel' to go back to the system settings.

Elements**Element****Description****'General settings'**

'Date format'

Displays the date format of the selected language.

'Number format'

Displays the number format of the selected language.

'Clinic'

Specifies the name of the institution, to which the device belongs. Mandatory.

'Calibration check'

'ORBSCAN Last check'

Displays the date of the last calibration check.

'ZYWAVE Last check'

8.5.2 Defining ORBSCAN 3 map settings

To define ORBSCAN 3 map settings, proceed as follows.

It is recommended to define the default settings for the ORBSCAN 3 maps in this section because they are valid for all following measurements.

If necessary, the user can then modify some of the ORBSCAN 3 map default settings for a performed measurement (Chapter 9.2.2 'Viewing ORBSCAN 3 measurement results' on page 9 - 6).

1. ➤ Open the menu page.

2. ➤ Click 'Settings'.

⇒ The settings page opens (Fig. 56).

3. ➤ Click 'ORBSCAN map settings'.

⇒ A dialog box opens allowing you to define: 'Overlay settings', 'Pachymetry settings', 'General settings' and 'Quad map settings' (Fig. 56).


Fig. 56: ORBSCAN 3 map settings

For more information about the maps, refer to Chapter 6.1.3 'ORBSCAN 3 maps' on page 6 - 1.

4. ➤ Edit the required field(s).

5. ➤ Click 'Cancel' to return to the settings page.

6. → Click 'Accept' to save the map settings.

⇒ The ORBSCAN 3 map settings are valid for all following measurements.

You can change the map options for a performed measurement too (↪ Chapter 9.2.2 'Viewing ORBSCAN 3 measurement results' on page 9 - 6).

WARNING

Danger of injury due to misuse of acoustic factor!

For pachymetry maps, the acoustic factor allows to adjust ZDW 3 measurements to ultrasound measurements. Misuse of the acoustic factor may lead to incorrect measurements, i.e. a false diagnosis.

- **Do not change the acoustic factor without knowledge of its meaning and consequences in general and the calibration and adjustment situation of the device specifically.**

Elements

Element	Description
'Overlay settings'	
'Circles'	Specifies a series of odd (1, 3, 5, 7) and even (2, 4, 6, 8) circles for displaying the diameter of the eye, in mm.
'Grid (radial)'	Enables a radial grid on all maps. The grid consists of axial lines every 45° which intersect in the eye center and have tick marks at every millimeter shown in black.
'Grid'	Enables a horizontal / vertical grid. The distance between the grid lines is 1 mm.
'Pupil'	Enables the center of the pupil with a white circular shape on all maps.
'Thinnest point'	Enables a yellow dot showing the thinnest area of the cornea on all maps.
'K-lines & K-labels'	Enables the K-lines and K-labels on the keratometric maps. K-lines show the steep and flat axis of astigmatism. K-labels show the steep and flat values of axis of astigmatism.
'Kappa'	Enables to display the angle Kappa at the degree marker on the maps. The angle Kappa is defined as the angle difference between the axis of fixation and the optical axis perpendicular to the lens.
'Numericals'	Enables to show the map data in numeric values in all maps.
'Pachymetry settings'	

Element	Description
'Pachymetry averages'	Enables pachymetry averages (average thickness at center, superior, inferior, temporal and nasal areas of the cornea) on pachymetry maps. Default: off
'Measurement ring diameter [mm]'	Determines the position of the circles in which the pachymetry averages are calculated.
'Measurement sample diameter [mm]'	Determines the size of the circles in which the pachymetry averages are calculated.
'Number of samples (including center)'	Determines the number of circles in which the pachymetry averages are calculated (3, 5, 7, 9, 11 and 13). One circle is always in the center of the map. The others are arranged in a circle around the center.

'General settings'

'Acoustic factor' Gives ultrasound equivalence on pachymetry maps. The possible values are between 0.9 and 1.1.

Note: The acoustic factor depends on the device settings and can be set by the administrator only.

'Quad map settings'

Specifies the default maps displayed in the measurement results. The system's default maps are:

- 'Anterior BFS'
- 'Posterior BFS'
- 'Axial power, keratometric map'
- 'Pachymetry'

Click one map icon  to select another map as your personal default.

8.5.3 Performing ORBSCAN 3 calibration checks (alignment PMMA)

For performing a calibration check—alignment PMMA (Polymethyl methacrylate, material of the calibration check tool), use the ORBSCAN 3 calibration check tool that is part of the scope of delivery.

Calibration reference value

The reference calibration value for performing the calibration check is displayed on the label of the calibration check tool (Fig. 57).

Tolerance range

The calibration check results should deviate max. ± 0.030 mm from the value on the calibration check tool label.

When to perform a calibration check

The user is required to perform a calibration check every 3 months by means of the calibration check tool included in delivery.

i

Make sure to perform calibration checks at regular intervals, as it is possible to perform measurements even if the calibration check results deviate from the reference calibration value by more than the specified tolerance range.

i

However, after an initial setup or severe jamming due to moving the device, the calibration needs to be checked by the user prior to using the ORBSCAN 3 again. If calibration check results deviate from the reference calibration value by more than the specified tolerance range, contact an Authorized Service Technician (Chapter 5 'Service and calibration' on page 5 - 1).

Performing calibration checks


Fig. 57: ORBSCAN 3 calibration check tool


To perform calibration checks (alignment PMMA), proceed as follows.

Precondition:

- The ORBSCAN 3 calibration check tool (Fig. 57) is placed onto the chin rest.
- The calibration interval for ORBSCAN 3 has not expired (Chapter 5 'Service and calibration' on page 5 - 1).
- The settings page is open.

⚠ CAUTION!

If any of the calibration checks is not performed successfully, the system must not be used to avoid serious permanent patient injury.


1. ➔ In 'Settings', click 'ORBSCAN calibration check'.

⇒ The 'Select/Add measurement' dialog box opens (Fig. 58).

Fig. 58: 'Select/Add measurement' dialog box

2. ➔ Click 'PMMA measurement'.

⇒ The alignment PMMA is now ready to start.

i

After about two minutes, a warning message opens asking you if you want to proceed with the alignment.

- Click 'Yes' ☒ Yes if you want to continue with the alignment.
- Click 'No' ☐ No if you want to go back to the 'Select/Add measurement' dialog.

 Start measurement

3. ➤ Align the measurement head as if you were performing a normal ORBSCAN 3 measurement (➤ Chapter 9.2.1 'Aligning the patient and starting the ORBSCAN 3 measurement' on page 9 - 3).

To trigger the measurement, you can also press the joystick button.


Fig. 59: Alignment PMMA

 Start measurement

4. ➔ Click 'Start measurement'.

i

To trigger the measurement, you can also press the joystick button.

⇒ A measurement is performed in the usual way.

The calibration check results (Best fit values) will be displayed directly to the right of 'BFS:' (Fig. 60/1).

i

The 'ORBSCAN 3 calibration check' is saved.


Fig. 60: Calibration check results

5. ➔ Compare the calibration check results (Fig. 60) with the value on the calibration check tool label (Fig. 57).

i

If the calibration check results deviate from the value on the calibration check tool label by more than the specified tolerance range, contact immediately an Authorized Service Technician.

Now you can:

- Print the quad map view of the calibration check results by clicking 'Print'.
- Add comments if you want to enter personal comments for this calibration check by clicking 'Add comments'.
- Export the measurement data on a USB memory stick by clicking 'Export to USB stick'.
- Compare measurements for the selected calibration check data taken at different times by clicking 'Difference'.
- Perform a new ORBSCAN 3 calibration check by clicking 'Another exam'.
- Go back to the settings page by clicking 'Continue'.

8.5.4 Performing ZYWAVE 3 calibration checks

For performing a calibration check, use the calibration check tool that is part of the scope of delivery. The calibration check consists of a trombone check and two checks with two different spheres. The calibration check spheres mimic a human eye with known values.

The ZYWAVE 3 calibration check tool consisting of a holder and two test spheres is placed onto the chin rest.

Test sphere 1 (named in the software 'Sphere #1'): -5

Test sphere 2 (named in the software 'Sphere #2'): +5

Tolerance ranges for the calibration check

Software checks	Tolerance ranges
Trombone	The measured range should maximal deviate 0.3 D from the expected range.
Test spheres	The spherical equivalent of every measurement should deviate max. 0.2 D from the average value as well as from the expected value. The cylinder of every measurement should deviates max. 0.15 D from the average value.

When to perform a calibration check

The user is required to perform a calibration check every 3 months by means of the calibration check tool included in delivery.

Make sure to perform calibration checks at regular intervals.

If the calibration check results deviate from the reference calibration value by more than the specified tolerance range, it is not possible to perform measurements any more. In this case, perform a calibration check again. If the calibration check results still deviate from the reference calibration value by more than the specified tolerance range, contact an Authorized Service Technician.

However, after an initial setup or severe jamming due to moving the device, the calibration needs to be checked by the user prior to using the ZYWAVE 3 again. If calibration check results deviate from the reference calibration value by more than the specified tolerance range, contact an Authorized Service Technician (Chapter 5 'Service and calibration' on page 5 - 1).

To perform calibration checks, proceed as follows:

Precondition:

- The ZYWAVE 3 calibration check tool is placed onto the chin rest.
- The calibration interval for ZYWAVE 3 has not expired (§ Chapter 5 'Service and calibration' on page 5 - 1).
- The settings dialog is open.

Fig. 61: ZYWAVE 3 calibration check tool

⚠ CAUTION!

If any of the calibration checks is not performed successfully, the system must not be used to avoid serious permanent patient injury.

1. ➔ In 'Settings', click 'ZYWAVE calibration check'.

⇒ The ZYWAVE 3 calibration check is now ready to be performed.

During alignment, the following problems may occur:

- a bad lenslet reflex (Fig. 62) or
- an unsharp simulated pupil image (Fig. 63).

f
re
ry
lts
ice
ali-
ra-
eci-
vice

Fig. 62: Alignment: bad lenslet reflex


Fig. 63: Alignment: unsharp pupil image

The following image displays a good alignment (Fig. 64).

Fig. 64: Alignment: good

Trombone check

Trombone check

Fig. 65: 'Trombone check' done and displayed ranges

2. ➔ Click 'Trombone check'.

⇒ A green check mark is displayed if the 'Trombone check' was successful. The 'Measured range' and the 'Expected range' values are displayed in the table on the right side (Fig. 65/1).

i

If the measured range deviates from the expected range by more than the specified tolerance range, the 'Trombone check' will not pass. In this case, contact immediately an Authorized Service Technician.

Calibration check

Start measurement

3. ➔ Align the measurement head in the direction of the first test sphere as if you were performing a normal ZYWAVE 3 measurement (Chapter 9.3 'Performing ZYWAVE 3 measurements' on page 9 - 14).

4. ➔ Click 'Start measurement'.

i

- The software recognize automatically which test sphere is actually being used for the calibration check.
- The calibration check with 'Sphere # 1' and 'Sphere # 2' must be performed three times each.

5. ➔ Repeat step 4 two times again.

- ⇒ A green check mark appears if the calibration check with the first test sphere was successful.

The results of the calibration check is displayed in the 'Result' area of the actually used test sphere.

6. ➔ Align the measurement head in the direction of the second test sphere and repeat from step 4 to step 5.

- ⇒ A dialog with information about the 'Trombone test results' and the 'User calibration tool results' opens when the calibration check of the second test sphere is finished.


Fig. 66: Calibration check info

- 1 Sphere
2 Cylinder values
3 Spherical equivalent values

7. ➔ Click 'OK'.

- ⇒ The green check marks for 'Sphere # 1' and 'Sphere # 2' are displayed together with the calibration check results (Fig. 67/1).


Fig. 67: Calibration check done

8. ➔ If required, click 'Print' to print the calibration check results.

9. ➔ Click 'Accept' to save the performed calibration check.
⇒ The settings page opens.

8.5.5 Viewing the product information

This section provides information about the software version, the product and the manufacturer.


To open the 'Information', proceed as follows.

1. ➤ Open the menu page.

2. ➤ Click 'Settings'.

⇒ The settings dialog box opens.

3. ➤ Click 'Information'.

⇒ The 'Information' dialog box opens (Fig. 68).



Fig. 68: 'Information' dialog box (example)

- 1 Software version
- 2 Product name
- 3 Manufacturer name and address
- 4 ORBSCAN ID
- 5 ZYWAVE ID

4. ➤ To close the 'Information', click 'Close' or 'Back'.

⇒ The settings page opens.

9 Measurements

This section provides information about performing ORBSCAN 3 and ZYWAVE 3 measurements with the ZDW 3 diagnostic software.

CAUTION

Danger of injury due to user inadvertency!

The user is responsible for any outcome when performing a measurement with ZDW 3. The user must supervise the whole measurement process.

Only the user is allowed to perform a measurement in the ZDW 3 software and to interpret and evaluate the measurement results. Failure to observe this responsibility may result in serious permanent patient injury.

The diagnostic software allows performing ORBSCAN 3 and ZYWAVE 3 measurements for the same patient without having to switch between the two applications.

9.1 Preparing for performing measurements

9.1.1 Disinfecting the head and chin rests

For a detailed description about disinfecting the head and the chin rests, refer to  Chapter 7.4 'Cleaning instructions' on page 7 - 12.

9.1.2 Selecting or adding a patient

1.  Click 'Patient' in the menu page.

 The patient list opens.

2.  Select a patient from the patient list ( Chapter 8.3.1 'Selecting patients' on page 8 - 8).

If the patient is not listed, the patient has to be added first  Chapter 8.3.3 'Adding patients' on page 8 - 15.

⚠ CAUTION**9.1.3 Instructing and positioning the patient**

Danger of injury due to incorrect moving of the chin rest!

Moving the chin rest with the patient's head firmly in place may lead to a movement against the muscular tonus.

- **Make sure that the patient is comfortable while moving the chin rest.**
- 1. ➤ Instruct the patient to place the chin on the chin rest, the forehead against the forehead support and to sit upright.
- 2. ➤ Adjust the height of the measurement head to the patient by means of the positioner (joystick) and the patient's head by means of the chin rest adjustment handhold.
- 3. ➤ Ensure that the patient's head is not tilted or turned at an angle to where the bridge of the nose might block the slits scanning the cornea.

9.2 Performing ORBSCAN 3 measurements

This section provides information about performing ORBSCAN 3 measurements.

⚠ CAUTION

Danger of injury due to uncontrolled exposure of an eye in front of this instrument!

The light from this device is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage.

- **Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 35 minutes for the slit lamp illumination during the adjustment of the eye in the image plane.**
- **For uncontrolled exposure of an eye in front of the instrument ocular damage can occur after 260s. Therefore the light is automatically deactivated after 120s when no measurement has been started.**

9.2.1 Aligning the patient and starting the ORBSCAN 3 measurement


Fig. 69: 'Select/Add measurement' dialog box

1. To start the measurement, click 'ORBSCAN' in the 'Select/Add measurement' dialog box.

⇒ The alignment is now ready to start (Fig. 70).

Fig. 70: Alignment ready

⚠ CAUTION! Danger of injury due to missing ability to fixate!
Missing ability to fixate accurately (due to mental impairment, impaired bodily functions, age of the patient) may lead to erroneous measurements.

- Make sure that the patient is able to fixate accurately.

The colors of the  button indicate:

- Gray: no measurement performed
- Green: measurement successful
- Red: measurement not successful

Aligning the patient

2. ➤ Make sure that the patient is comfortable while focusing.

i

After about two minutes, a warning message opens asking you if you want to proceed with the alignment.

- Click 'Yes' ☒ Yes if you want to continue with the alignment.
- Click 'No' ☒ No if you want to go back to the 'Select/Add measurement' dialog.

i

- If the patient has droopy eyelids, there are several possibilities:
 - Instruct the patient to open his eyes wide.
 - Use a cotton bud to open the patient's eyes wide.
 - Or have another person to open the patient's eyes wide, taking care not to press on the eyeball while holding the lid open. (This also works if the patient is unable to keep his eye open during the measurement.)
 - Cover the patient's opposite eye using a tissue or occluder.
 - Instruct the patient to blink a few times, then to open both eyes wide.

3. ➤ Instruct the patient to look at the blue fixation target. If you are taking OD, make sure that the patient is using OD to fixate. The same applies to OS.

■  Shows the patient's **right** eye.

■  Shows the patient's **left** eye.

4. ➤ To focus the pupil of the patient, perform an adjustment twisting the joystick clockwise / counterclockwise, i.e. performing movements in the z-axis (moves the measurement head upwards and downwards) and performing movements in the x-axis (moves measurement head left or right).


5. ➤ To align the two half slits over the central alignment circle, so that the outer radii of the half slits touch each other (Fig. 71/1), perform movements in the y-axis (move the joystick forwards / backwards).

Fig. 71: Alignment of the half slits

6. ➤ Instruct the patient to hold steady, blink once and open both eyes wide and fixate on the light.

Starting the measurement

7. To start the measurement, click 'Start measurement' or press the trigger button of the joystick.

⇒ The measurement results dialog opens displaying the results in form of a quad map (Fig. 72).


Fig. 72: ORBSCAN 3 measurement results

- 1 Selection field for left or right eye and measurement date and time
- 2 Specification for the displayed measured eye with date and time of measurement
- 3 Quad map values
- 4 Quad map view

Now, you can:

- Click 'Add comments' to comment this measurement (this option is only available if no comment is available for this measurement).
- Click 'Export to USB stick' to export measurement results.
- Click 'Difference' to compare two measurements from the same patient if there is another measurement available.
- Click 'Print' to print the measurements results.
- Click 'Delete' to delete this measurement.
The delete option is only available for single measurements within an active measurement session. After closing the measurement session, the delete option is no longer available.
- Click 'Another exam' to perform a new ORBSCAN 3 measurement.
- Click 'Continue' to go back to the patient data dialog.

9.2.2 Viewing ORBSCAN 3 measurement results

Description

This section provides information about viewing ORBSCAN 3 measurement results and changing map settings.

Viewing the quad map

The quad map allows viewing four different maps at the same time.

By default, the quad map (Fig. 73/1) includes an anterior and posterior elevation map, a keratometric map and a pachymetry map.


Fig. 73: Quad map view

In the quad map view, you can

- Zoom the map
- Change colors and the type of the color scale
- Select another map for the quad map view
- Change the overlay, pachymetry and general settings
- Display available OD and OS measurement
- Switch between map and raw data
- Change map settings

i

Precondition: Ensure that the measurement results are open.

Zooming maps

To zoom maps, proceed as follows:

1. Click 'Zoom in'.


Fig. 74: ORBSCAN 3 map enlarged

⇒ The map is displayed enlarged (Fig. 74).

2. Click 'Zoom out' to minimize the map and return to the quad map view.

Changing the color legend

All maps have a color legend that can relate to eye elevation, curvature or pachymetry values. The color legends can be changed while viewing a map by clicking the color legend.

To change the color legend, proceed as follows:

1. Click 'Color legend'.

⇒ The color legend dialog box opens (Fig. 75) displaying the 'Number of colors', 'Color scale type', 'Min./Max. values' and the 'Step size' of the map.

Fig. 75: Color legend dialog box

2. Select the required number of colors and color scale type.

⇒ The number of colors or the color scale type are changed and displayed in the color legend.

Changing map types

To change map types, proceed as follows:

1. Click the map types icon.
⇒ A list with available map types opens.
2. Select the required map type.
⇒ The selected map type is applied. The map header displays the selected map (Fig. 76/1).

Fig. 76: Available ORBSCAN 3 maps for the quad map view

For information about the ORBSCAN 3 maps, refer to Chapter 6.1.3 'ORBSCAN 3 maps' on page 6 - 1.

Parameters and values

The measurement results contain a statistical box that shows patient information including:

- Date and time of the measurement (Fig. 77/1)
- simulated K-readings, 3 and 5 mm zone statistics, best fit values, thinnest point, kappa angle (kappa), kappa intercept, acoustic factor, anterior chamber depth (Ep), white-to-white, and Ø (pupil) (Fig. 77/2).


Fig. 77: Quad map values

Data areas

Element	Description
Measurements drop-down list box	Shows the eye type and the date and time of the measurement. Furthermore, it allows to select between existing measurements of a patient, which have been saved in the patient data.
'OD measurement' / 'OS measurement'	Specifies the date and time of the measurement.
'Sim Ks'	Simulated keratometry readings. Crosses that show local directions of principal curvatures where blue lines mark the flat axis and red lines mark the steep axis.
'Astig.'	Cylinder, [D] Axis, [°]
'Maximum'	Max-K, [D] Axis, [°]

Element	Description
'Minimum'	Min-K, [D] Axis, [°]
'3 mm zone' and '5 mm zone'	Show the average and standard deviation of the mean and astigmatic power and the steep and flat axes of astigmatism.
'Irregularity'	Measure of the variation of mean and astigmatic curvatures of the eye over a specific area (diameter 3 or 5 mm). Irregularities, [D]
'Mean power'	Mean power, [D]
'Astig power'	Astigmatic power, [D]
'Steep axis'	Steep axis, [°]
'Flat axis'	Flat axis, [°]
'Thinnest point'	Displays the thinnest pachymetry in [μm] and [x.x,y.y].
'Kappa'	Displays the angle kappa [°] at the degree marker, that is defined as the angle difference between the axis of fixation and the optical axis.
'Kappa intercept'	Coordinates [x.x, y.y]
'Acoustic factor'	Displays the ultrasound equivalence on pachymetry maps.
'Anterior chamber depth (Ep)'	Detected from the epithelium. [mm]
'White-to-white'	Displays limbus to limbus measurement across the cornea. [mm]
'Ø (pupil)'	As detected during the scanning procedure, [mm].
'Best fit values'	'Anterior map' 'Posterior map'
'Central pow. / Radius'	Curvature of the Aconic surface which best fits the surface. It could be expressed as its inverse, the Radius of curvature. [mm]/[D]
'Astig. D/°'	Difference between K values or Axes of the steepest and flat-test meridians of the aconic surface of reference.
'Asphericity (Q)'	Specifies the asphericity and type (prolate or oblate) of conic section which best fits a corneal meridian. Unless otherwise specified, it refers to the meridian with least curvature (flat-test meridian). Related to Shape Factor (P): Q= P-1.

Element	Description
'Eccentricity (e)'	Describes a conic section and the rate of curvature change away from the apex of the curve.
'Shape factor (P)'	Specifies the asphericity and type (prolate or oblate) of conic section which best fits a corneal meridian. Unless otherwise specified, it refers to the meridian with least curvature (flat-test meridian).
'BFS'	[mm]/[D]
'Apical pow. / Radius'	Curvature of the ellipse which best fits the surface. It could be expressed as its inverse, the Radius of curvature.
'Maximum K.'	[D]
'Minimum K.'	[D]
'Astig. D / °'	Difference between K values or Axes of the steepest and flat-test meridians of a conic section fitting the surface Cylinder, [D] Axis, [°]

Map overlay settings allow to change map parameters, which helps in interpreting the data (☞ 'Changing map settings' on page 9 - 12).

By clicking in a map, a cross marks the position and the concurrent readouts are displayed on the brink of the choosen map (Fig. 78): for example the according radius and meridian, the depth, or the map value (Z-axis).

Displaying OD and OS measurements

To display OD and OS measurements, proceed as follows:

➔ Select the desired 'OD' or 'OS' measurement of the displayed list (Fig. 78/1).

⇒ The displayed quad maps switch from OD to OS (or vice versa).


Fig. 78: Switching between OS and OD

Changing map settings

In the map options dialog box, you can change map parameters, which helps in interpreting the data.

Fig. 79: Open quad map settings

1. ➔ On the left side of the quad map, click  (Fig. 79/1).
⇒ The map options dialog box opens (Fig. 80/1).


Fig. 80: Map options dialog box

2. ➔ Edit the required field(s).
⇒ Your settings are displayed in the maps.

For information about the settings, refer to  Chapter 8.5.2 'Defining ORBSCAN 3 map settings' on page 8 - 31.

3. ➔ Click  to close the map options dialog box.

Switching between map and raw data

To switch between map and raw data, proceed as follows:

The quad map view is set as default for displaying the ORBSCAN 3 measurement results.

The active view is always highlighted in green. ^

Option controls

Element

Description

'Map' / 'Raw'

Switches between the map and a picture of the eye.

Raw

1. Click 'Raw' to switch to the raw data images.

⇒ The raw data images open (Fig. 81). Four raw data images are available:

- Placido disc image: 'Average placido'
- Alignment image: 'Half slit'
- Slitscan image (all slitscan images merged within one image): 'Average slit scan'
- Eye image averaged: 'Whole eye'

Fig. 81: Raw data images

Map

2. Click 'Map' to switch back to the quad map view.

9.3 Performing ZYWAVE 3 measurements

This section provides information about performing ZYWAVE 3 measurements.

9.3.1 Aligning the patient and starting the ZYWAVE 3 measurement

Precondition:

1. Click 'Patients' in the menu page.

⇒ The patient list opens.

2. Select a patient from the patient list (↪ Chapter 8.3.1 'Selecting patients' on page 8 - 8).

If the patient is not listed, the patient has to be added first.

⇒ The patient account opens.

1

Fig. 82: Patient account: subjective refraction

3. → In the 'Subjective refraction' dialog box (Fig. 82), enter the required patient values, specify the diagnosis and the dilation status of the patient eye (☞ "Subjective refraction' values' Table on page 9 - 15).

The 'Diagnosis' and 'Dilation' fields are mandatory for the ZYWAVE 3 measurement (Fig. 82/1). The listed options apply for both eyes (OD and OS). However, the diagnosis as well as the dilation for the partner eye can be changed individually, if required.

Additional to the listed 'Diagnosis' options it is possible to enter manually a personal diagnosis of the patient eye.

'Subjective refraction' values

Element	Description
Ranges for Sphere, Cylinder and Axis	■ 'Sphere [D]' - Min: -25 D - Max: +10 D ■ 'Cylinder [D]' - Min: -10 D - Max: +10 D ■ 'Axis [deg]' - Min: 0° - Max: 180° The default value is 0.00 if no other value is set.
Diagnosis values	■ 'PreOp' ■ 'PostOp' ■ 'One day' ■ 'One week' ■ 'One month' ■ 'Three months' ■ 'Six months' ■ 'One year' ■ 'Cycloplegic' Additional to the listed options it is possible to enter manually a personal diagnosis of the patient eye. The default value does not contain any entry.
Dilation values	■ 'Dilated' ■ 'Undilated'

4. → Click 'Accept' to save the entered values.

5. ➔ To start the measurement, click 'ZYWAVE' in the 'Select/Add measurement' dialog box.

⚠ CAUTION! Danger of injury due to missing ability to fixate!

Missing ability to fixate accurately (due to mental impairment, bodily functions, age of the patient) may lead to erroneous measurements.

- Make sure that the patient is able to fixate accurately.

⇒ The alignment is now ready to be performed.


Fig. 83: Alignment ZYWAVE 3

i

The current pupil size is displayed on the alignment screen.

Centration, focus and fixation

1. ➔ Instruct the patient to only look at the center of the fixation image (green cross). Explain that, although the fixation image will be out of focus (for the scan it is in focus for only fractions of a second), the patient should always concentrate on looking at the center of the image.

i

Inform the patient, that the green cross is a little bit unsharp.

2. ➔ Advise the patient not to look directly at the red laser light dot to avoid accommodation, but slightly above it.
3. ➔ If you are taking OD, make sure that the patient is using OD to fixate. The same applies to OS.
- Shows the patient's **right** eye.
 - Shows the patient's **left** eye.

4. → Center the pupil at the intersection of the crosshairs on the screen by using the joystick:

- Perform movements in the x-axis by shifting the joystick to the side (moves measurement head left or right).
- Perform movements in the z-axis by twisting the joystick clockwise / counterclockwise (moves the measurement head upwards and downwards).
- Perform movements in the y-axis by moving the joystick forwards / backwards (moves measurement head forwards or backwards).

i

- *A yellow cross and a warning symbol will be displayed on the alignment screen if the pupil is not centered enough. In this case, center the crosshairs as far as possible in the middle of the pupil, i.e. on the yellow cross or next to it.*
- *In the picture of the pupil camera, the pupil is overlaid with a circular ring called 'pupil finder' or 'pupil circle'. The major part of the pupil finder should be green.*

5. ➤ Adjust the focus so that the iris structures are in focus and the edge between iris and sclera is a sharp line.

i

Always center the crosshairs in the middle of the pupil, not in the middle of the rings of light points. Do not focus on the reflections of the pupil illumination.

i

The color of the pupil finder is correlated to the acceptance criteria for alignment described in Chapter 9.3.2 'Verifying the measurement images' on page 9 - 21. The pupil finder must be located on the edge between pupil and iris. For this, the complete pupil must be visible. Make sure that it is not covered by the eyelid or eyelashes.

i

If the limbus or pupil have not been detected at all, repeat 'Centration, focus and fixation' on page 9 - 16.

Possible reasons for errors in the pupil detection:

Reflections from the room illumination overlap the iris/pupil.

Eyelids or eyelashes are interfering with the pupil finder.

Reflections from the infrared LEDs are interfering with the pupil finder. In such a case, the most prominent reflections are usually outside the pupil area or on the pupil edge.

Possible reasons for errors in the limbus detection:

Reflections from the room illumination overlap the iris.

Eyelids or eyelashes are covering part of the iris.

Eye is not sufficiently exposed.


Fig. 84: Alignment ZYWAVE 3

- ⇒ The patient eye is good aligned if the major part of the pupil finder is displayed green and the yellow cross becomes green.

Taking iris images

6. →

i

You can take iris images during the measurement session or separately without performing a ZYWAVE 3 measurement.

i

For information about taking iris images separately without performing a ZYWAVE 3 measurement, refer to Chapter 9.4 Taking iris images on ZYWAVE 3' on page 9 - 23.

Performing the pupilometry

7. →

i


You can perform the pupilometry at the begin of a measurement session (the displayed pupilometry values on the online pupilometer are valid for all performed ZYWAVE 3 measurements during this measurement session).

ZYWAVE 3 enables the user to perform a pupilometry prior to the wave-front measurement. The pupil size can be measured over the desired period of time, and a measurement curve is displayed in the online 'Pupilometer' window.

i

For information about performing the pupilometry, refer to Chapter 9.5 'Performing the pupilometry on ZYWAVE 3' on page 9 - 24.

Starting the measurement

8. Click 'Start measurement' to perform the ZYWAVE 3 measurement.

After clicking 'Start measurement' the live image area is bordered green for a few seconds. It indicates that the measurement has been started (Fig. 85).

⇒ The iris images in raw data are displayed allowing you to verify them immediately after the measurement is completed (Fig. 86).


Fig. 85: ZYWAVE measurement started


Fig. 86: ZYWAVE 3 iris images raw data

- 1 **Centroid image with grid lines:** The centroid image shows grid lines that exhibit how the individual points fit together. The analysis should connect each neighboring point with a grid line. Review the centroid images to verify good quality of the centroids (sharpness) as well as good detection of the centroids (e.g. no missing centroids except in the periphery border).
- 2 **Pupil image:** The blue circle shows the found limbus, the green or yellow circles show the edge of the identified pupil. The user must verify that the pupil was correctly identified by the pupil finder. In addition, the focus during the measurement must be checked.
- 3 The table lists the five measurements and highlights the three most repeatable ones. It also shows the Predicted Phoropter Refraction (PPR) that the patient would have if measured with a phoropter**.
- 4 The default setting for the **Wavefront map** displays the higher order wavefront.

** To translate the wavefront deformation into numbers which are commonly used in the refractive surgery (manifest refraction), a software module analyzes the shape of the wavefront and determines the best spherocylindrical correction to optimize the retinal image, taking into account all higher aberrations. These derived spherocylindrical correction values are called "Predicted Phoropter Refraction" (PPR).

Now, you can:

- Verify the measurement images (☞ Chapter 9.3.2 'Verifying the measurement images' on page 9 - 21).
- View the measurement results by clicking 'Results' (☞ Chapter 9.3.3 'Reading the ZYWAVE 3 measurement results' on page 9 - 22).
- Export images by clicking 'Export images to USB stick'.
- Go back to the alignment page by clicking 'Back'.

9.3.2 Verifying the measurement images

Here, you can verify the pupil and centroid images of the new ZYWAVE 3 measurement. To verify measurement images, proceed as follows.

1. ➤ Click one of the pupil pictures (Fig. 86/3) to enlarge and verify it.
⇒ The 'Verification zoom' window opens (Fig. 87).

Fig. 87: Verification zoom window

The color of circle and cross indicates the alignment quality.

- **Green** means that the alignment is inside an acceptable range from the red alignment cross.

The iris is in focus so that in at least three images of the eye the iris structures are sharp and clearly visible, i.e. the pupil finder is placed correctly (green circle) and the limbus has been found (blue circle).

2. ➤ Click 'Video' to start the measurement film.

3. ➤ Select if you want to show the overlays or not by clicking the control switch on or off.

4. ➤ Click 'Close' to close the enlarged image again.

The highlighted number (e.g. picture 1 in Fig. 87) shows which image you selected from the image sequence.


The color of the pupil finder is correlated to the acceptance criteria for alignment described in . If the center of the pupil was inside the acceptance criteria during the measurement, the circle is green, if it was outside the acceptance criteria, the circle is yellow.

The second (blue) circle in the images of the eyes shows the circular limbus finder.

9.3.3 Reading the ZYWAVE 3 measurement results

To read the measurement results, proceed as follows.

1. ➤ After verifying the iris images, click 'Results' to open the display options of measurement results.
⇒ The display options of measurement results are displayed (Fig. 88).
2. ➤ To view the different measurement display options, click the respective tab:

For information about the ZYWAVE 3 measurement display options, refer to & Chapter 6.2.3 'ZYWAVE 3 display options of measurement results' on page 6 - 5.

Fig. 88: ZYWAVE 3 display options of measurement results

- 1 'Summary'
- 2 '2D Plot'
- 3 'PPR vs. Pupil diameter'
- 4 'Higher order PSF'
- 5 'Raw data'
- 6 'Normal band'
- 7 'Difference'

According to the selected tab, you can:

	'Summary'	'2D Plot'	'PPR vs. Pupil diameter'	'Higher order PSF'	'Raw data'	'Normal band'	'Difference'
Comment this measurement by clicking 'Add comments'.	x	—	—	—	—	—	—
This option is available if no comment is available for this measurement.							
Display the available comments for this measurement by clicking 'View comment'.		—	—	—	—	—	—
Export the measurement results by clicking 'Export to USB stick'	x	—	—	—	x	—	—
Print this measurement by clicking 'Print'.	x	x	x	x	x	x	x
Delete this measurement by clicking 'Delete'.	x	x	x	x	x	x	x
The delete option is only available for single measurements within an active measurement session. After closing the measurement session, the delete option is no longer available.							
Go back to the patient data dialog by clicking 'Continue'.	x	x	x	x	x	x	x
Perform a new ZYWAVE 3 measurement by clicking 'Another exam'.	x	x	x	x	x	x	x

9.4 Taking iris images on ZYWAVE 3

You can take iris images during the measurement session or separately without performing a ZYWAVE 3 measurement.

Adjust the room lights prior to measurement in order to achieve an optimum pupil size and verify that the main reflexes of the infrared pupil illumination ring are within the pupil diameter and not interfering with the pupil boundaries.

No bright light sources are allowed.

To take iris images separately without performing a ZYWAVE 3 measurement, proceed as follows.

The diagnostic software enables the user to perform a pupillometry prior to the wavefront measurement.

i

You can perform the pupillometry at the begin of a measurement session (the displayed pupillometry values on the online pupillometer are valid for all performed ZYWAVE 3 measurements during this measurement session).


Fig. 91: Online 'Pupillometer' (1)

Online 'Pupillometer'

Element	Description
Target off (Darkness)	The target illumination is switched off, the green cross is disable for the patient and the red laser light is on.
Target on (Dim light)	The target is illuminated, the green cross is enable for the patient and the red laser light is on.
Blue curve	Represents the current pupil size during the pupillometry.
Current pupil size	Is the current pupil size. This value is automatically renewed approx. 16 times per second. A stabilization of the pupil size can be observed in the online 'Pupillometer', typically after 10-20 seconds. <i>If the pupillometry is reset by means of the joystick, the maximum pupil size will be displayed as the current pupil size.</i>
Maximum pupil size	Is a maximum determinated value and will be displayed inside the online 'Pupillometer' area.

4. ➤ If you want to perform a pupillometry with dark light conditions, click 'Target off'.
- ⇒ The green cross has been switched off (the patient cannot see it anymore). Advise the patient not to look directly at the red laser light dot to avoid accommodation, but slightly above it.

 Stop pupilometer

5. ➤ Click 'Stop pupilometer'.
- ⇒ The maximum pupil size value ('Max. pupil size') for the selected light condition is displayed in the online 'Pupilometer' area directly after 'Target on' (Fig. 92).

i

If the pupillometry is reset by means of the joystick, the maximum pupil size ('Max. pupil size') will be displayed as the 'Current pupil size'.

Now you can:

- Perform a measurement by clicking 'Start measurement'.
- Take iris images by clicking 'Take iris images'.
- Go back to the patient account by clicking 'Back.'

9.6 Working with measurements

This section provides information on how to work with existing measurements.

9.6.1 Selecting and closing measurements

Selecting measurements

To select ZYWAVE 3 and ORBSCAN 3 measurements, proceed as follows.

1. ➤ Select a patient (☞ Chapter 8.3.1 'Selecting patients' on page 8 - 8).
2. ➤ Double-click the required measurement in the 'Select/Add measurement' dialog box to open.

⇒ The measurement results open.

i

If the selected measurement is a ZYWAVE II or an ORBSCAN IIz measurement, following message will be displayed: 'This measurement cannot be displayed. It can only be used for treatment calculation.'

Elements

Description

Indicates ZYWAVE II measurements.

Indicates ORBSCAN IIz measurements.

According to the selected measurement, you can:

ORBSCAN 3 measurements


- Click 'Add comments' to comment this measurement (this option is only available if no comment is available for this measurement). Otherwise, click 'View comments' to display the available comments for this measurement.
- Click 'Export to USB stick' to export measurement results.
- Click 'Difference' to compare two measurements from the same patient if there is another measurement available.
- Click 'Print' to print the measurement results.
- Click 'Back' to go back to the patient data dialog.

ZYWAVE 3 measurements

	'Summary'	'2D Plot'	'PPR vs. Pupil diameter'	'Higher order PSF'	'Raw data'	'Normal band'	'Difference'
Comment this measurement by clicking 'Add comments'.	x	—	—	—	—	—	—
This option is available if no comment is available for this measurement.							
Display the available comments for this measurement by clicking 'View comment'.	x	—	—	—	—	—	—
Export the measurement results by clicking 'Export images to USB stick'	x	—	—	—	x	—	—
Print this measurement by clicking 'Print'.	x	x	x	x	x	x	x
Go back to the patient data dialog by clicking 'Back'.	x	x	x	x	x	x	x

Closing measurements

To close measurements, proceed as follows.

- ➔ In the measurement results page, click 'Patient'.
- ⇒ The patient list opens.

You can also close the measurement by clicking 'Back'. In this case, the patient data page will be displayed.

9.6.2 Deleting measurements

To delete measurements, proceed as follows:

- The measurement results (for ORBSCAN 3) or the display options of measurement results (for ZYWAVE 3) are opened.

The delete option is only available for single measurements within an active measurement session. After closing the measurement session, the delete option is no longer available.

1. Click 'Delete' to delete the measurement.

⇒ A warning message opens (Fig. 93).



Fig. 93: Delete measurement warning message

2. Click 'Yes' to delete the measurement.

⇒ The measurement is deleted.

3. Click 'No' to continue working with the measurement.

9.6.3 Comparing measurements

9.6.3.1 Comparing ORBSCAN 3 measurements


You can compare two measurements and assess the differences between them. The differences are represented in D, mm, or μm , depending on the selected type of map.

A comparison is only possible if at least two different measurements of the same eye are available.

To compare measurements, proceed as follows:

1. Select a measurement (Chapter 9.2.2 'Viewing ORBSCAN 3 measurement results' on page 9 - 6) or start a new measurement (Chapter 9.2 'Performing ORBSCAN 3 measurements' on page 9 - 2).

2. After the measurement opens, click 'Difference'.

⇒ The difference page opens.

By opening of the 'Difference' option, the available measurements will be automatically loaded and can be selected from the 'A' and 'B' drop-down list boxes.

Fig. 94: Compare dialog box

3. ➤ From the 'A' drop-down list box (Fig. 95/1), select the first measurement you want to compare and from the 'B' drop-down list box (Fig. 95/2), select the second measurement you want to compare with the first one.

⇒ The differences between both measurements are displayed on the 'A-B' dialog box (Fig. 95/3).


Fig. 95: Comparing ORBSCAN 3 measurements A-B dialog box

i

You can change the colors and the map types and perform overlay, pachymetry and general settings. If measurements for both eyes are available, you can switch between OD and OS.

4. ➤ Click 'Print' if you want to print the comparison results.

9.6.3.2 Comparing ZYWAVE 3 measurements

You can compare two measurements and assess the differences between them.

The 'Difference' option enables the user to subtract one wavefront examination (or wavefront) from another one. This is a tool for evaluating the development of the patient's eye over time or to compare preOp and postOp measurements.

i

A comparison is only possible if at least two different measurements of the same eye are available.

To compare measurements, proceed as follows:

1. Open the measurement results and click the 'Difference' tab.

⇒ The difference page opens.

By opening of the 'Difference' option, the available measurements will be automatically loaded and can be selected from the 'A' and 'B' drop-down list boxes.

Fig. 96: 'Difference' option

2. ➤ From the 'A' drop-down list box (Fig. 97/1), select the first measurement you want to compare, and from the 'B' drop-down list box (Fig. 97/2), select the second measurement you want to compare with the first one.
 - ⇒ The differences between both measurements are displayed in the 'A-B' dialog box (Fig. 97).

Fig. 97: Comparing ZYWAVE 3 measurements A-B dialog box


You can change the colors and, if measurements for both eyes are available, you can switch between OD and OS.

9.6.4 View data on the USB memory stick

View measurements data on the USB memory stick

To view data on USB memory stick, use the USB memory stick delivered with the diagnostic device.

1. ➤ Insert the USB memory stick into the USB port on the all-in-one PC right side.
2. ➤ Log in into the diagnostic software (Chapter 8.1.1 'Logging in' on page 8 - 3).
3. ➤ Click 'View data on USB stick' to open the patient data on your USB memory stick.
 - ⇒ Now you can view data on USB memory stick.

4. ➔ To view patient data, click the patient name once in the 'Patients' dialog box (Fig. 98).


Fig. 98: 'Measurements' dialog box

- ⇒ If there are any measurements saved for the patient you have selected, these appear in the 'Measurements' dialog box below.

9.6.5 Exporting measurements


You can export ORBSCAN 3 and ZYWAVE 3 measurement data.

Exporting means copying, i.e. the patient data remains on ZDW 3.

To export measurements, proceed as follows:

Precondition:

- The measurement results are opened.
- The 'Export to USB' button is active. This means that a USB memory stick is connected to ZDW 3.

1. Click 'Export to USB' to export this measurement data to your USB memory stick.

⇒ The export options dialog box opens (Fig. 99).



Fig. 99: Export options dialog box

The 'Default' format (binary data) is selected as default.

Do not change this setting. Otherwise a error message will be displayed.


2. Select the desired options and click 'Save'.

⇒ A message opens informing you that the export has been successful (Fig. 100).

Fig. 100: Export successful message

3. Click 'OK'.

⇒ The measurement opens again.

9.6.6 Exporting ZYWAVE 3 pupil images

Exporting ZYWAVE 3 pupil images

To export images, you can use a USB memory stick delivered with the diagnostic device.

1. Log in into the diagnostic software (↪ Chapter 8.1 'Working with the software' on page 8 - 1).

2. Click 'Patients'.


3. Select a patient from the patient list (↪ Chapter 8.3.1 'Selecting patients' on page 8 - 8).

⇒ The patient account opens.

4. ➤ Select a ZYWAVE 3 measurement from the 'Select/Add measurement' dialog box by double-clicking it.
⇒ The measurement opens.
5. ➤ Click the 'Raw data' tab.
⇒ The pupil (upper display part) and centroid (lower display part) images are displayed (Fig. 101).

Fig. 101: Pupil and centroid images

6. ➤ Click 'Export images to USB'.
⇒ The 'Select images' dialog box opens (Fig. 102).


The used images for the wavefront calculation are selected as default.

Fig. 102: 'Select images' dialog box

7. ➤ Select the 'Pupil images' and the 'Centroid images' that you want to export.
8. ➤ Click 'Cancel' to reset the selected images.
⇒ The 'Select images' dialog box closes.

9. Repeat from step 6 to 7 and click 'Accept' to export the selected ZYWAVE 3 pupil images.

⇒ A warning message opens informing you that the export has been successful (Fig. 103).

The images have been saved on your USB memory stick: E:\Technolas\ZyWaveImages.

Fig. 103: Export successful warning message

10. Click 'OK'.

⇒ The measurement opens again.

9.6.7 Printing measurement results

To print the ORBSCAN 3 and ZYWAVE 3 measurement results a network color printer is necessary.

Network printer

The network color printer is an optional component which has to be purchased by the customer in a specialized store.

- The installation of the network color printer must be done by an Authorized Service Technician only.
- The network color printer can be connected to the TPV Gigabit-Ethernet-Switch or directly to the all-in-one PC by an Authorized Service Technician.

Only network color printers supported by HP Universal Print Driver Series can be used with the diagnostic device.

- HP Universal Print Driver Series for PCL6, PCL5, PS are already installed in the delivered diagnostic device.

NOTICE

Danger of property damage due to misplacement of the network color printer!

The network color printer must not be placed in the immediate proximity to patients.

- Place the network color printer at least 1.5 meters away from the diagnostic device.

9.6.7.1 Printing ORBSCAN 3 measurement results

To print ORBSCAN 3 measurement results (e.g. quad map), proceed as follows:

Precondition:

- The required printer is connected to ZDW 3.
- The measurement results are opened.

Map

Raw

Print

1. Click 'Map' to print the quad map and all available values.
2. Click 'Print'.
 - ⇒ The Windows print options open.
3. Select the desired print options in the print options window of your printer, and click 'Print'.
 - ⇒ The ORBSCAN 3 measurement results will be printed.

9.6.7.2 Printing ZYWAVE 3 display options of measurement results

To print ZYWAVE 3 display options of measurement results, proceed as follows:

Precondition:

- The measurement results are opened.
- The required printer is connected to ZDW 3.

1. ➔ Click 'Print'.
 - ⇒ The Windows print options open.
2. ➔ Select the desired print options in the print options window of your printer, and click 'Print'.
 - ⇒ The ZYWAVE 3 measurement results will be printed.

9.6.8 Performing treatment calculations

The following procedures apply to OD and OS.

Preparing treatment calculations

If required, set default treatment values.

For information about setting default treatment values, refer to Chapter 8.2.3 'Setting default treatment values' on page 8 - 7.

9.6.8.1 Selecting treatments and diagnostic data

Selecting treatments

To select treatments proceed as follows.

1. ► Select a patient (Chapter 8.3.1 'Selecting patients' on page 8 - 8).

Fig. 104: 'Select/Add treatment' dialog box

2. ► In the 'Select/Add treatment' dialog box (Fig. 104), double-click the required treatment to open it.

If ZYWAVE II and ORBSCAN IIz measurement data are available in the external Database Server, they can be used in the ZDW3 software for treatment calculations:

Indicates ZYWAVE II measurements.

Indicates ORBSCAN IIz measurements.

- ⇒ A warning message opens if the selected treatment was calculated using different nomogram settings (Fig. 105).


Fig. 105: Different nomogram settings warning message

3. Click 'OK'.

⇒ The selected treatment opens in the treatment calculation page (Fig. 106).

The screenshot shows a treatment calculation interface for a patient named SMITH John. The interface is divided into several sections with numbered callouts:

- 1**: Patient selection icon.
- 2**: Patient name and ID: SMITH John, Patient ID: 967, Patient ID: VBST.
- 3**: Subjective refraction section with fields for Sph (D), Cyl (D), and Axis (°).
- 4**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 5**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 6**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 7**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 8**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 9**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 10**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 11**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 12**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 13**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 14**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 15**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 16**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.
- 17**: IOL power section with fields for IOL power (D), IOL axis (°), and IOL type.

The interface also includes a 'Simulated ablation' section with a 'Comment' field and an 'Accept' button.

Fig. 106: Selected treatment (example): elements

9 MEASUREMENTS

4. ▶ If required, enter and adjust correction values.




For more information about entering and adjusting correction values, refer to Chapter 9.6.8.3 'Entering and adjusting correction values' on page 9 - 53


For more information about the available treatments, refer to the TECHNOLAS® TENEO™ 317 user manual.

Selected treatment: Elements (Fig. 106)

Pos.	Button / Icon	Description
1		Click to open the patient list to search for, select or add a patient.
2	 Subjective refraction values with corresponding PPR values.	Click to enter or adjust the following parameters: subjective sphere, subjective cylinder and subjective axis values fields. The corresponding PPR values are not editable.
3		Click to switch between the plus and minus cylinder convention.
4		Diagnostic parameters for the selected eye: ■ If there are only ZYWAVE measurements available the default 'K-value' is 43,3 and the default 'Q-value' is -0,20. ■ If there are ZYWAVE and /or ORBSCAN measurements available the 'K-value' and 'Q-value' will be calculated on the base of the ORBSCAN measurement. ■ Pachymetry: 'Pachymetry [µm]' initially shows the thinnest pachymetry value obtained from ORBSCAN 3. ■ Flap thickness: The default value for flap thickness is 120 µm.
5	 	'OD' or 'OS' icons indicate the selected eye.
6		List of the imported diagnostic data: ■ 'ZYWAVE' —date of examination at ZYWAVE 3. ■ 'PPR 3.5 mm' —predicted phoropter refraction at pupil size 3.5 mm (for ZYOPTIX HD Treatment). Wavefront-derived refraction from the ZYWAVE 3 measurement. ■ 'Iris data' —available or not. ■ 'Wavefront Ø' —wavefront diameter. ■ 'Z400/HO-RMS' —values from ZYWAVE 3. ■ 'ORBSCAN' —date of examination at ORBSCAN 3.

Pos.	Button / Icon	Description
7	 Diagn. data	Click to select diagnostic data. A list of ORBSCAN 3 and ZYWAVE 3 measurements is displayed.
8	 Iris data	Click to select another iris data if available.
9		Click to discard the calculation and restore the original values.
10	 Log out	Click to log out from the software.
11	 Accept	Click to accept the calculation. If the calculation is not possible, a message informs the user accordingly. It is not possible to select a treatment and the 'Accept' button is not enabled.
12		Detailed treatment maps and 'Comment' field: - Enter a comment for this calculation in the 'Comment' field. - To verify the calculated correction values, use detailed maps.
13	 Print	Click to print a treatment calculation.
14		Click to select a treatment map: ■ Anterior elevation map ■ Radial pachymetry ■ Axial power, keratometric map ■ Simulated ablation ■ Wavefront Full Pupil ■ Wavefront HO Pupil ■ Axial power, estim. postOp keratometric map ■ Anterior elevation, estim. postOp ■ Radial pachymetry, estim. postOp
15		List of treatment values, value and date of the applied nomogram and nomogram comment (not editable).
16		To enter or adjust following parameters if necessary: Target value, optical zone, and calculated correction values (correction sphere, correction cylinder and correction axis).
17	 Back	Click to go back to the patient account.

Selecting diagnostic data

1. If required, click 'Diagn. data' (Fig. 106/7) to select another ORBSCAN 3 and ZYWAVE 3 measurement.
 - ⇒ The 'ORBSCAN 3 measurements' and 'ZYWAVE 3 measurements' dialog boxes opens (Fig. 107/1 and 2).

Fig. 107: ORBSCAN 3 and ZYWAVE 3 measurements dialog boxes

2. Click 'Cancel' to go back to the treatment calculation page and repeat from step 2 to step 5.
3. Select one of the ZYWAVE 3 measurements and, if available, one of the ORBSCAN 3 measurements.
 - ⇒ The 'Accept' button will be activated.

4. Click 'Accept' to open the respective measurement.
 - ⇒ The selected diagnostic data is displayed in the diagnostic data dialog box (Fig. 108) in the treatment calculation page.


Fig. 108: Diagnostic data dialog box

5. ➤ If required, enter and adjust parameters.

For more information about entering and adjusting parameters, refer to
 ➤ Chapter 9.6.8.2 'Adding treatments' on page 9 - 45.

6. ➤ Click 'Undo' to restore the original values.

⇒ A warning message opens asking you if you want to restore the original values (Fig. 109).


Fig. 109: Restore the original values warning message

7. ➤ Click 'Yes' to restore the original values or 'No' to continue.

8. ➤ Enter your comments in the 'Comment' field and click 'Accept'.

⇒ Your comment is displayed in the 'Comment' field.

9.6.8.2 Adding treatments

To add treatments proceed as follows.

1. ➤ Select a patient.
 2. ➤ In the 'Select/Add treatment' dialog box, click 'OD' or 'OS'.
- ⇒ The treatment calculation page opens.

Fig. 110: Treatment calculation page

3. ➔ Select the desired treatment typ for calculation (Fig. 110/1): 'PROSCAN', 'ZYOPTIX HD', 'SUPRACOR' and 'PTK'.

4. → Enter new values or check imported values (Fig. 111).

If required, adjust a value based on your expertise and experience as follows:

Fig. 111: Diagnostic parameters
(example)

Subjective refraction The calculation software uses the initial subjective refraction values stored with the patient data.

The subjective refraction is the basis for treatment calculation. The user can manually adjust the subjective refraction values for sphere, cylinder, and axis.

Depending on the selected treatment, the nomogram modifications directly affect the following correction values:

Subjective sphere:	Correction sphere for all treatments. Correction sphere for ZYOPTIX HD Treatment and SUPRACOR Treatment.
Subjective cylinder:	Correction cylinder for PROSCAN Treatment and SUPRACOR Treatment. For ZYOPTIX HD Treatment: The cylinder value from the ORBSCAN 3 measurement is used.
Subjective axis:	Correction axis for PROSCAN Treatment and SUPRACOR Treatment.

K-value (imported from ORBSCAN 3) Verify the 'K-value' and consider the calculation result. If required, adjust the value.

Q-value (imported from ORBSCAN 3) Verify the 'Q-value' and consider the calculation result. If required, adjust the value.

Pachymetry value 'Pachymetry [μm]' initially shows the thinnest pachymetry value obtained from ORBSCAN 3. If required, the user can enter or modify the value in the 'Pachymetry [μm]' field, based on an alternative measurement technique.

The modification directly affects the calculated residual stromal thickness.

Flap thickness The default value for flap thickness is 120 μm . The user can set an individual default value for flap thickness in the 'calculation settings' page.

The user can also manually enter or modify the flap thickness in the 'Flap thick. [μm]' field.

The residual stroma thickness is adjusted depending on the defined flap thickness.

- ⇒ Calculated correction and treatment values are immediately displayed in the treatment tabs except for PTK (Fig. 112).

Fig. 112: Correction refraction values

Verifying or adjusting calculated parameters

5. → Check the values of the calculated treatment parameters and consider the situation.

If required, adjust a value based on your expertise and experience as follows:

Verify and, if required, adjust the treatment parameters:

Target

Verify the target value and consider the calculation result. If required, adjust the diopter value in the 'Target [D]' field (default value is 0 D or - 0.5 D for SUPRACOR Treatment).

Optical zone

The default value for the optical zone is user-defined. The user can set an individual default value for the optical zone diameter in mm in the 'Calculation settings' page (↪ Chapter 8.2.3 'Setting default treatment values' on page 8 - 7). The user can also manually modify the default value for the optical zone in the 'Opt. zone. [mm]' field.

Residual stromal thickness, treatment zone, count pulses, and the estimated ablation data vary depending on the defined optical zone, nomogram, correction refraction and flap thickness.

SUPRACOR Treatment

6. → For SUPRACOR Treatment, select either 'Mild' or 'Regular' laser pulse patterns.

Verifying or adjusting correction values

⚠ WARNING! Danger of injury due to incorrect diagnostic and measurement parameters!

Omission of verifying diagnostic and measurement parameters may result in serious permanent patient injury.

Always confirm and, if necessary, adjust the following parameters before the confirmation or adjustment of measurement values:

- Subjective refraction values, such as subjective sphere, cylinder, and axis.
- K- and Q-value.

7. ➤ The calculated correction values depend on the user-defined nomogram, the subjective sphere, and the target value (§ Chapter 9.6.8.3 'Entering and adjusting correction values' on page 9 - 53).

For ZYOPTIX HD treatments, the values for 'Corr. sph. [D]', 'Corr. cyl. [D]', and 'Corr. axis [°]' are linked to the PPR from the ZYWAVE 3 wavefront examination. The PPR values are displayed in the diagnostic data pane; the 'Wavefront Ø' is listed in the 'Diagn. data' table.

⇒ Correction and treatment values are immediately calculated and displayed in the respective treatment tab

Target [D]:	Opt. zone [mm]:	
0.00	+6.50	
Corr. sph. [D]:	Corr. cyl. [D]:	Corr. axis [°]:
+2.00	+1.00	0
Nomogram: 100 %		
Treated zone: 10 mm x 10 mm		
Pulse/Duration: 8764		
Max. ablation: 62 µm		
Center ablation: 2 µm		
Residual: 350 µm		

Fig. 113: Correction values (example for PROSCAN Treatment)

Users should be aware of the following points:

- If the entered value is outside the accepted range, the input field is bordered in red.
- If the correction values are modified by the user, the treatment parameters are not recalculated.
- If you have modified the subjective refraction or the target value, the correction values are calculated and updated again based on the current input parameters.

Verify the calculated correction values and consider the calculated treatment values. If necessary, modify the calculated correction values to adjust the targeted outcome as desired:

- Correction sphere
- Correction cylinder
- Correction axis

9 MEASUREMENTS

PTK treatment

8. → For PTK treatment (phototherapeutic keratectomy), the user enters the required ablation and optical zone into the input fields ('Ablation [μm]' and 'Opt. zone [mm]'), (Fig. 114).

PTK

Ablation [μm]: 55 Opt. zone [mm]: 4.00

Nomogram: 0 % Max. ablation: 56 μm
Treated zone: 7.2 mm x 7.2 mm Center ablation: 55 μm
Pulses/Duration: 5509 Residual: 413 μm

4/22/2013 6:17:34 PM PROSCAN modified TT

Fig. 114: PTK treatment tab (example)

Reading iris data

The button 'Iris data' is only enabled if in the database another iris data of the selected eye in undilated status is available.

If necessary, click 'Iris data'.

→ The 'Select alternative reference image' dialog opens.

Select alternate reference image

OD Date: 2/27/2014 10:04:14 AM Dilation: Undilated Exam. diam.: 6.65 mm

OD Date: 2/27/2014 10:17:57 AM Dilation: Undilated Exam. diam.: 6.62 mm

OD Date: 2/27/2014 10:14:14 AM Dilation: Undilated Exam. diam.: 6.65 mm

OD Date: 2/27/2014 10:04:14 AM Dilation: Undilated Exam. diam.: 6.48 mm

OD Date: 2/27/2014 10:04:14 AM Dilation: Undilated Exam. diam.: 6.64 mm

Rotation angle: 0°
Pupil shift total: 131 μm
Pupil shift X: -59 μm
Pupil shift Y: -117 μm

Eye recognized: Use pupil shift compensation: 1

Cancel Use original iris data Accept

Fig. 115: Imported iris data with corresponding pupil pictures

- 10.** ➤ If necessary, select another iris image with a optimum pupil diameter by clicking one of the listed iris images. The pupil images are overlaid with circular rings which indicate the alignment quality.

blue circle	indicates the limbus
yellow circle	indicates the pupil
green circle	indicates a valid measurement
red circle	indicates the center circle and corresponding cross-hairs

- ⇒ The selected reference iris image is displayed in the 'Reference iris' section.

The 'Accept' button is only enabled if in the selected iris image the eye was recognized.

You can click 'Cancel' at any time to go back to the treatment calculation page without saving the performed changes.

- 11.** ➤ If necessary, activate the option 'User pupil shift compensation (PSC)'.

Pupil shift compensation only applies to ZYOPTIX HD treatments.

- 12.** ➤ Click 'Accept' to accept the changes.

⇒ The treatment calculation page opens.

- 13.** ➤ To use original iris data click 'Iris data' again.

⇒ The 'Select alternative reference image' dialog opens.

Click 'Use original iris data' and then 'Accept'.

The treatment calculation page opens.

Now you can:

- Print previously saved treatment calculation by clicking 'Print'.
- Go back to the treatment calculation page by clicking 'Back'.

Verifying the calculated correction values

14. To verify the calculated correction values, click the maps icon for detailed maps (Fig. 116/1).

⇒ A list with the available maps opens. Now, you can select a map.


Fig. 116: Detailed maps icon

15. Optionally, enter a comment on your treatment in the 'Comment' field (Fig. 117).

Fig. 117: 'Comment' field

16. Click 'Undo' to reset your modifications to the original values.

⇒ A warning message opens asking you whether you want to restore the original values or not (Fig. 118).


Fig. 118: Restore the original values warning message

17. When the calculation is successful and all parameters have been specified as needed, accept the appropriate treatment by clicking 'Accept'.

⇒ The calculated treatment is saved and displayed in the 'Select/Add Treatment' dialog box (Fig. 119).


Fig. 119: Calculated and accepted treatments (example)

General remarks

For a detailed list of possible error messages and how to rectify them, refer to Chapter 10 'Troubleshooting' on page 10 - 1.

WARNING

Danger of injury due to incorrect diagnostic and measurement parameters!
The failure to verify diagnostic and treatment parameters may result in serious permanent patient injury.

Always confirm and, if necessary, adjust the following parameters before the confirmation or adjustment of measurement values:

- Subjective refraction values, such as subjective sphere, cylinder, and axis.
- Pachymetry value
- Residual stromal thickness (at least 250 µm)
- K- and Q-values.

If a nomogram has been defined by the user, this is always applied to the corresponding measurement calculation.


Applied values are displayed in the respective measurement tab.

Reading measurements with wavefront size smaller than the default optical zone:


When reading measurements with wavefront size smaller than the default optical zone, the calculation software reduces the (default) optical zone to the wavefront size (but only in the 'ZYOPTIX HD' tab), which is indicated by a warning message. This warning message opens only once and will not be issued again if you change the eyes later on.

- For all measurements recheck and, if applicable, adjust the optical zone.


Invalid input

If an entered parameter or a calculation result is not within the allowed ranges, the input field is bordered in red.

Modified entry field

Modified input is bordered in orange.

Critical residual stroma value (example)

Critical treatment values are highlighted in red.

For a detailed list of the elements on the treatment calculation page, refer to Chapter 9.6.8.1 'Selecting treatments and diagnostic data' on page 9 - 40.

Subjective refraction values with corresponding PPR values

Fig. 120: Subjective refraction values (example)

Click in the available input fields to enter or adjust the following parameters: subjective sphere, subjective cylinder or subjective axis.

The corresponding PPR values are not editable.

Plus or minus cylinder convention

Fig. 121: Plus or minus cylinder convention

Click to switch between plus and minus cylinder convention of subjective refraction. The values for PPR, subjective, and correction refractions are adjusted accordingly.

Switching between the cylinder conventions has no impact on the use and calculation of the nomogram.

If a nomogram has been defined by the user, this is always applied to the corresponding measurement calculation.

Applied values are displayed in the respective measurement tab.

Diagnostic parameters for the selected eye

Fig. 122: Diagnostic parameters for the selected eye (example)

i

- If there are only ZYWAVE measurements available the default 'K-value' is 43,3 and the default 'Q-value' is -0,20.
- If there are ZYWAVE and /or ORBSCAN measurements available the 'K-value' and 'Q-value' will be calculated on the base of the ORBSCAN measurement.
- '**Pachymetry [μm]**': If required, modify the value in the pachymetry field, based on an alternative measurement technique. The modification directly affects the calculated residual stromal thickness.
- '**Flap thick. [μm]**': If required, manually modify the flap thickness in the flap thickness field. The residual stroma thickness is adjusted depending on the defined flap thickness.

i

The default value for flap thickness is 120 μm .

Target value, optical zone and calculated correction values

PROSCAN		
Target [D]:	Opt. zone [mm]:	
0.00	6.50	
Corr. sph. [D]:	Corr. cyl. [D]:	Corr. axis. [°]:
+3.60	-0.70	2

Fig. 123: Target value, optical zone, and calculated correction values (example)

If necessary, click in the available input fields to enter or adjust the following parameters: target value, optical zone and calculated correction values (correction sphere, correction cylinder value and correction axis)

9.6.8.4 Printing treatment calculations

To print treatment calculations, proceed as follows:

Only saved treatment calculations can be printed.

Precondition:

- The treatment calculation is open.
- The required printer is connected to ZDW 3 (☞ Chapter 9.6.7 'Printing measurement results' on page 9 - 37).

1. ➤ Click 'Print'.
 - ⇒ The Windows print options open.
2. ➤ Select the desired print options in the print options window of your printer, and click 'Print'.
 - ⇒ The treatment calculation will be printed.

10 Troubleshooting

The following chapter describes possible causes of faults and the work required to rectify them.

WARNING

Danger of injury due to improper repair!

Improperly executed work in conjunction with troubleshooting may lead to serious injuries and significant property damage.

- Perform only troubleshooting tasks that local personnel can perform according to the troubleshooting chart.
- If you are in doubt about how to perform the troubleshooting tasks that are described in the troubleshooting chart, contact an Authorized Service Technician.
- Troubleshooting tasks that are not described in the troubleshooting chart must be performed only by an Authorized Service Technician.

Standard Windows messages

If any Windows message opens, restart the system. The Windows operating system of ZDW 3 is protected in a way that after restart the initial configuration is restored. Contact an Authorized Service Technician, if errors appear directly at startup.

Troubleshooting charts

-  'User control screen' on page 10 - 1
-  'ZDW 3' on page 10 - 2
-  'ORBSCAN 3' on page 10 - 4
-  'ZYWAVE 3' on page 10 - 6
-  'Treatment calculation' on page 10 - 8

Should the problem persist or should you encounter a problem not listed in the following tables, contact an Authorized Technolas Perfect Vision Representative.

User control screen

Fault description	Cause	Remedy	Personnel
User control screen is blank.	Power cables of ZDW3 are not plugged in, monitor is not switched on or user control screen defect.	■ Switch the device on. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician

10.1 Software messages

ZDW 3

Fault description	Cause	Remedy	Personnel
'A problem has occurred! Switch the device off and switch it on again.' 'If the problem persists, contact an Authorized Service Technician.'	—	■ See message.	User Authorized Service Technician
'Calibration check data is locked. The calibration check data can only be accessed by one user at a time.'	—	■ Perform the calibration check later.	User
'Comment cannot be saved.'	Entered comment is too long.	■ Reduce the number of characters in the comment field.	User
'Date of birth must be in the past.'	—	■ Enter correct date.	User
'Date of birth must not be blank.'	—	■ Enter date of birth.	User
'Data cannot be exported.'	—	■ Check the storage capacity of your USB stick.	User
'Entered password and confirmation password must match.'	—	■ Enter a valid password and identical confirmation password.	User
'Invalid entry!'	—	■ Enter a valid and complete entry.	User
'Invalid date'	—	■ Enter a valid date.	User
'Login name already in use.'	—	■ Enter a different login name.	User
'First name: Field must not be blank.'	—	■ Enter a first name.	User
'Last name: Field must not be blank.'	—	■ Enter a last name.	User
'Measurement cannot be deleted.'	The database connection has been interrupted.	■ Check if the database connection has been re-established. ■ If problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'Measurement cannot be deleted from USB stick.'	The USB stick is probably damaged or not available.	■ Check if the USB stick is inserted in the USB port. ■ If the USB stick is damaged, insert another Technolas Perfect Vision GmbH approved USB stick in the USB port.	User
'Measurement cannot be loaded.'	—	■ Check your USB stick.	User
'Minimum password length is 4 characters.'	—	■ Enter a valid password.	User

Fault description	Cause	Remedy	Personnel
'No printer found!'	—	■ Check whether the printer is switched on. ■ Check whether the correct printer is connected. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'The network connection to the following device has been lost: ORBSCAN/ZYWAVE. Please check if all cables are properly connected and if the device is switched on.' 'Restart the device.'	—	■ See message. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'ORBSCAN/ZYWAVE hardware currently not available. Please verify that the device is switched on and that all cables are properly connected.'	—	■ See message. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'Password confirmation missing.'	—	■ Enter an identical confirmation password.	User
'Password invalid'	—	■ Enter a valid password and identical confirmation password.	User
'Password missing'	—	■ Enter a valid password and identical confirmation password.	User
'Patient data cannot be deleted from USB stick.'	The USB stick is probably damaged or not available.	■ Check if the USB stick is inserted in the USB port. ■ If the USB stick is damaged, insert another Technolas Perfect Vision GmbH approved USB stick in the USB port.	User
'Remaining time for recalibration (days): nn.' 'After the remaining time has expired, it will not be possible to perform measurements. You can only view existing measurements and calculate treatments.'	—	■ Contact an Authorized Service Technician.	Authorized Service Technician
'The following devices need to be recalibrated in one month: ORBSCAN/ZYWAVE.'	—	■ Contact an Authorized Service Technician.	Authorized Service Technician
'The deadline for recalibration has expired. Measurements are not possible until recalibration has been performed.'	—	■ Contact an Authorized Service Technician.	Authorized Service Technician

10 TROUBLESHOOTING

Fault description	Cause	Remedy	Personnel
'The selected treatment cannot be edited.'	The selected treatment has already been performed on the laser system.	—	User
'The USB stick is not available. The measurement cannot be loaded.'	—	■ Insert the USB stick.	User
'The USB stick is not available. The patient data cannot be loaded.'	—	■ Insert the USB stick.	User
'The database is currently not available. Please verify that all cables are properly connected.'	—	■ Contact an Authorized Service Technician.	User
'The remote database is not available. Do you want to use the local database instead?'	—	■ Select the local database if required.	User
'Unknown error.'	—	■ Repeat the measurement. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'Login name or password incorrect.'	—	■ Enter a valid login name and password.	User
'Value is not numeric.'	—	■ Enter numeric values.	User
'Value must be in the range from nn to nn.'	Entered values are outside the defined range.	■ Review the values and make adjustments as appropriate.	User
'Value must be in the range from nn to nn in nn steps.'	Entered values are outside the defined range and steps.	■ Review the values and make adjustments as appropriate.	User

ORBSCAN 3

Fault description	Cause	Remedy	Personnel
'HCB error. Switch the device off and switch it on again.' 'If the problem persists, contact an Authorized Service Technician.'	Hardware error.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'Insufficient image quality due to too much ambient light.'	Too much ambient light.	■ Reduce the ambient light level and make sure that the defined requirements for the operating conditions are met.	User

Fault description	Cause	Remedy	Personnel
'Measurement cannot be evaluated.'	Software problem.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'Measurement failed. Possible causes are: Slits are blocked by the nose. Try rotating the patient's head slightly.' 'Abnormally-shaped cornea. Try centering on the limbus.' 'Reflections from external light sources.' 'The patient blinked or the eyelids were not fully open.' 'The cornea was not well-centered. Repeat the measurement.'	See message.	■ See message.	User
'Measurement failed. Possible causes are: 'Slits are blocked by the nose. Try rotating the patient's head slightly.' 'Abnormally-shaped cornea. Try centering on the limbus.' 'The patient blinked or the eyelids were not fully open.' 'The cornea was not well-centered.' 'Not enough image information available to perform proper eye registration.' 'Repeat the measurement.'	See message.	■ See message.	User
'Measurement failed. Repeat the measurement. If the problem persists, contact an Authorized Service Technician.' '(error code nnnn)'	Depends on the error code number.	■ Repeat the measurement. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'Slits cannot be detected. The slits may be shifted, tilted, incomplete, blurry, etc.'	See message.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'The ORBSCAN hardware has reported an error. Please restart the device and try again. If the problem persists, contact an Authorized Service Technician.'	Hardware error.	■ Restart the device and try again. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'The HCB is responsive, but the camera is not. Restart the device. If the problem persists, contact an Authorized Service Technician.'	—	■ Restart the device and try again. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'The alignment has automatically terminated earlier than expected. If the problem persists, restart the device.'	—	■ If the problem persists, restart the device.	User

Fault description	Cause	Remedy	Personnel
'The measurement cannot be compared to another measurement because required data is missing. The measurement may have been imported from an OTE file.'	ORBSCAN 3 measurement is compared with ORBSCAN IIZ measurement.	■ Comparison is not possible.	User
'The measurement failed. Do you want to repeat the measurement?'	—	■ If desired, repeat the measurement.	User
'During measurement significant eye movement was detected. Possible causes are: 'The patient moved excessively or blinked. Repeat the measurement.' 'Abnormal corneal shape. Try centering on the limbus.' 'Not enough image information available to perform proper eye registration. Repeat the measurement.'	See message.	■ Repeat the measurement. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician

ZYWAVE 3

Fault description	Cause	Remedy	Personnel
'Although the wavefront diameter of this measurement is sufficient for a ZYOPTIX HD treatment, at least one of the three used lenslet images indicates an eye lid coverage beyond the recommended limit.' 'Actual coverage: nn mm.' 'Recommended limit: nn mm.' 'Using this measurement in a ZYOPTIX HD treatment is not recommended.'	See message.	■ Perform a new measurement.	User
'Taking iris images failed.'	—	■ Align and focus patient's eye.	User
'Enter a diagnosis for at least one eye!'	—	■ See message.	User
'Enter the dilation status for at least one eye!'	—	■ See message.	User
'For this sphere a sufficient number of measurements is already available.'	This sphere has been already measured three times.	■ Align the measurement head in the direction of the other sphere.	User
'The measurement cannot be evaluated.'	Calibration check calculation failed.	■ Repeat calibration check. ■ If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician

Fault description	Cause	Remedy	Personnel
'The estimated spherical equivalence of the test sphere is nn D. Please use the test spheres with -5 D or +5 D.'	Test sphere not recognized.	- Repeat calibration check. - If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'The results of the preliminary and the final measurements are significantly different. Please repeat the measurement.'	—	See message.	User
'The latest calibration which is available for the connected ZYWAVE is incomplete. Recalibrate the device and perform all steps. The following data is missing: nn'	ZYWAVE calibration incomplete data.	Contact an Authorized Service Technician.	Authorized Service Technician
'The ZYWAVE cannot be identified. Make sure that the device is properly connected.'	ZYWAVE 3 cameras will not be identified.	Contact an Authorized Service Technician.	Authorized Service Technician
'The ZYWAVE hardware does not respond. Please restart the device. If the problem persists, contact an Authorized Service Technician.'	—	- Switch off the standby switch for a few seconds and restart the device. - If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician
'Initializing the ZYWAVE hardware for a measurement failed. If the problem persists, restart the device. If the problem still persists, contact an Authorized Service Technician.'	—	See message.	User Authorized Service Technician
'The previous measurement failed. The hardware is not ready for a new measurement. Please restart the device.'	—	See message.	User
'The measurement failed because of an insufficient number of processable images.'	—	Repeat the measurement.	User
'Measurement failed.'	—	Repeat the measurement.	User
'Unable to determine the optimum measurement position.'	—	Repeat the measurement.	User
'Unable to determine the optimum measurement position because of insufficient number of images suitable for evaluation.'	—	Repeat the measurement.	User
'The ZYWAVE hardware does not respond within the expected time frame. If the problem persists, check the network connection and restart the device.'	—	See message.	User

Fault description	Cause	Remedy	Personnel
'The measurement failed because the pupil size was too instable. Repeat the measurement.'	—	■ See message.	User
'Target position failure. If the problem persists, restart the system by powering off. Then switch off the power supply for a few seconds and restart the device.'	—	■ See message.	User
'The latest calibration check for this ZYWAVE (serial number: nn) has failed. If the failure is due to an imprecise alignment during the calibration check, redo the check. If the calibration check continues to fail, contact an Authorized Service Technician and have the device recalibrated.'	—	See message.	User Authorized Service Technician
'The difference between subjective refraction and PPR is so large that both refractions must be considered suspicious. This measurement may not be suitable as a basis for a wavefront-guided treatment.' 'Make sure that the subjective refraction is correct and repeat the measurement.'	—	■ See message.	User
'The optimum scanning position calculated during prescan was suspicious. This measurement may not be suitable for a wavefront-guided treatment.' 'Make sure that the subjective refraction is correct and repeat the measurement.'	—	■ See message.	User
'Failure to save the calibration check result to the external database.'	The database connection has been interrupted.	■ Check if the external database connection has been re-established. If the problem persists, contact an Authorized Service Technician.	User Authorized Service Technician

Treatment calculation

Fault description	Cause	Remedy	Personnel
'Ablation depth exceeds safety range.'	The treatment calculation leads to an ablation depth of greater than the safe range.	■ As appropriate, choose different treatment parameters or choose a different treatment mode.	User
'Calculated residual stromal thickness less than 250 µm.'	The treatment calculation results in a residual stromal thickness of less than 250 µm.	■ Review the selected treatment parameters, such as flap thickness or optical zone, and make adjustments as appropriate. ■ Do not treat patients with residual stroma < 250 µm.	User

Fault description	Cause	Remedy	Personnel
'Calculation error'	—	■ Review the values and make adjustments as appropriate.	User
'Calculation error. Calculated values are outside the supported range.'	—	■ Review the values and make adjustments as appropriate.	User
'Correction refraction outside the supported range.'	The calculated correction refraction values are out of supported range.	■ Review the subjective refraction and treatment values and make adjustments as appropriate.	User
'CRC error in data transferred to ZTC_lib517.dll.'	Hardware error.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'CRC error in results from ZTC_lib517.dll (ablation data).'	Hardware error.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'CRC error in results from ZTC_lib517.dll (laser pulse list).'	Hardware error.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'CRC error in results from ZTC_lib517.dll (localResult).'	Hardware error.	■ Contact an Authorized Service Technician.	Authorized Service Technician
'Entered flap thickness less than 100 µm!'	The entered flap thickness is below the allowed value.	■ Review the selected flap thickness and make adjustments as appropriate, unless you are certain you want to continue with the original flap thickness.	User
'Flap thickness out of range.'	The entered flap thickness is not within the allowed values.	■ Review the selected flap thickness and make adjustments as appropriate, unless you are certain you want to continue with the original flap thickness.	User
'Input parameter exceeds safety range.'	The entered parameter exceeds safety range.	■ Review the entered value and make adjustments as appropriate.	User
'Input parameter out of range.'	The entered parameter is out of the defined range.	■ Review the entered value and make adjustments as appropriate.	User

Fault description	Cause	Remedy	Personnel
'No wavefront data available (ZYWAVE diagnosis data required).'	—	■ Import ZYWAVE diagnosis data. ■ Perform a wavefront measurement at the ZYWAVE.	User
'Not enough centroids in ZYWAVE measurement.'	The wavefront measurement is based on too few centroids, which may be an indication of a low quality centroid image resulting in unreliable wavefront information.	■ Consider using a different ATE file, if available. ■ Alternatively, repeat the wavefront measurement at the ZYWAVE.	User
'Not enough laser pulses.'	The number of laser pulses required for the ablation profile is out of the supported range.	■ Evaluate the situation carefully. ■ As appropriate, adjust the correction and/or optical zone.	User
'Offset between PPR and correction refraction too large.'	The difference between correction and PPR refraction is larger than the supported range.	■ Review the subjective refraction and correction and make adjustments as appropriate.	User
'Optical zone cannot be larger than wavefront diameter.'	The entered optical zone is larger than the imported wavefront diameter.	■ Review the optical zone and make adjustments as appropriate, unless you are certain you want to continue with the original wavefront diameter.	User
'Pachymetry out of range.'	The entered pachymetry is invalid.	■ Review the selected pachymetry and make adjustments as appropriate.	User
'Please confirm treatment data! Different cylinder conventions used for PPR and corr. refraction.'	—	■ Review the used cylinder conventions and make adjustments as appropriate.	User
'Please confirm treatment data! Different cylinder conventions used for PPR and subj. refraction.'	—	■ Review the used cylinder conventions and make adjustments as appropriate.	User
'Please confirm treatment data! Different cylinder conventions used for subj. and corr. refraction.'	—	■ Review the used cylinder conventions and make adjustments as appropriate.	User
'PPR and correction refraction are of different type.'	The refraction types (myopia, hyperopia or mixed astigmatism) of PPR refraction and correction are not the same.	■ Review the correction and make adjustments as appropriate.	User

Fault description	Cause	Remedy	Per- sonnel
'PTK parameters invalid.'	—	■ Review the values and make adjustments as appropriate.	User
'Subjective and correction refraction are of different type.'	The refraction types (myopia, hyperopia or mixed astigmatism) of subjective refraction and correction are not the same.	■ Review the correction and make adjustments as appropriate.	User
'Subjective refraction invalid.'	The entered subjective refraction values are invalid.	■ Review the subjective refraction values and make adjustments as appropriate.	User
'Target invalid.'	The entered parameter is out of the defined range or not numerical.	■ Review the entered value and make adjustments as appropriate.	User
'The advanced nomogram result is out of range (PPR sphere minus calculated correction sphere must be within range -3 to +2 D).'	—	■ Evaluate the situation. ■ Verify the subjective refraction.	User
'The current settings cannot be saved. Do you want to cancel the changes?'	—	■ See message.	User
'The original treatments were calculated using different nomogram settings. 'New calculations will apply the current nomogram settings!'	The nomogram settings were adjusted after a treatment had been calculated or different nomogram values are used from different installations of the device. If the user attempts to recalculate this treatment, the diagnostic system alerts the user that the actual nomogram settings will be applied to the recalculation.	■ Verify that the current nomogram settings are appropriate.	User
'The difference between subjective axis and PPR axis is more than 15 degrees!'	Depending on the case, this warning message may appear alone or in combination.	■ Evaluate the situation. ■ Review the subjective refraction and make adjustments as appropriate.	User
'The difference between PPR axis and correction axis is more than 15 degrees!'	Depending on the case, this warning message may appear alone or in combination.	■ Evaluate the situation. ■ Review the correction values and make adjustments as appropriate.	User

Fault description	Cause	Remedy	Personnel
'The difference between subjective axis and correction axis is more than 15 degrees!'	Depending on the case, this warning message may appear alone or in combination.	■ Evaluate the situation. ■ Review the subjective refraction and make adjustments as appropriate.	User
'The difference between subjective cylinder and PPR cylinder is more than 0.5 D!'	Depending on the case, this warning message may appear alone or in combination.	■ Evaluate the situation. ■ Review the subjective refraction and make adjustments as appropriate.	User
'The difference between subjective sphere and PPR sphere is more than 0.75 D!'	Depending on the case, this warning message may appear alone or in combination.	■ Evaluate the situation. ■ Review the subjective refraction and make adjustments as appropriate.	User
'There is no clinical experience with Q-values higher than +0.5!'	The Q-value has been adjusted to a value above +0.5.	■ Review the treatment values and consider the situation carefully before proceeding.	User
'The patient you intend to treat with SUPRACOR is younger than forty years. Please check if the patient meets the inclusion criteria.'	Entered date of birth indicates that the patient is younger than 40 years.	■ If the patient is indeed younger than 40 years, refer to the patient selection criteria in this user manual and verify carefully whether the patient is eligible for SUPRACOR.	User
'Treatment not supported.'	—	■ Verify if you have selected the suitable patient for calculation.	User
'Treatment not supported in SUPRACOR.'	SUPRACOR only supports hyperopic and myopic presbyopia treatments.	■ Verify if you have selected the suitable patient for calculation.	User
'Wavefront size too small.'	Wavefront size smaller than the defined values: 4.5 mm for hyperopic and 5 mm for other wavefronts.	■ Import different ZYWAVE diagnosis data. ■ Perform a new wavefront measurement at the ZYWAVE.	User
'Wavefront size smaller than default optical zone. Optical zone has been reduced to wavefront size!'	The calculation software automatically reduces the (default) optical zone to the wavefront size.	■ For SUPRACOR or ASPHERIC treatments recheck and, if applicable, adjust the optical zone. ■ For PERSONALIZED treatments it is not possible to calculate an enlarged optical zone.	User

11 Decommissioning and disposal

WARNING

Risk of injury in conjunction with improper decommissioning and disposal!

- **Decommissioning and disposal is performed only by Authorized Service Technicians.**

Technolas Perfect Vision GmbH will arrange for the collection and recycling of (a) end-of-life Technolas Perfect Vision GmbH equipment and (b) end-of-life equipment from other manufacturers when a Technolas Perfect Vision GmbH product is purchased to replace the older equipment and has identical functionality. To initiate the collection of end-of-life equipment, contact an Authorized Service Technician.

SKU: 70005910 * Version 2 * UM-100006805 GB

12 Technical data

12.1 General specifications

12.1.1 Dimensions and weights

Diagnostic device dimensions (complete system)

Data	Value	Unit
Weight (without database server)	170	kg
Dimensions	218 x 88 x 147	cm

Measurement head

Data	Value	Unit
Weight	25	kg
Dimensions	53 x 45 x 54	cm
Movement ranges (measurement head)	x ± 40	mm
	y ± 37	mm
	z ± 15	mm
Chin rest (stroke: up/down)	≥ 55	mm

Center table

(with PC and distribution box)

Data	Value	Unit
Dimensions	60 x 66 x 122	cm
Weight	32	kg

Lift tables

Data	Value	Unit
Dimensions	61 x 56 x 94	cm
Weight	32	kg

12.1.2 Essential performance

Measurement accuracy in defined ranges as well as matching the patient and measurement data are considered as essential performance of the ZDW 3.

12.2 Operating conditions

Ambient operating conditions

Data	Value	Unit
Ambient temperature	10 to 40	°C
Relative humidity, maximum (non condensing)	30 to 50	%
Operating Altitude [maximum height above sea level (SL)]	3000	m

Operating time

Data	Value	Unit
Operating life	5	years

12.3 Transport and storage conditions

- Store the diagnostic device in a dry and dust-free location.
- Protect the diagnostic device from direct sunlight.
- Avoid mechanical vibration.

Data	Value	Unit
Ambient temperature	-10 to +55	°C
Relative humidity, maximum (non condensing)	10 to 90	%
Altitude [maximum height above sea level (SL)]	3000	m

12.4 Electrical connection

Data	Value	Unit
Electrical safety class (IEC 60601-1)	I	
International Protection Rating (IEC 60529)	IP2X	
Rated voltage	100 - 240	V ~
Power input	max. 260	W
Rated supply frequency	50 / 60	Hz

13 Index

1, 2, 3 ...

3 mm zone.....	9 - 6
5 mm zone.....	9 - 6

A

Abbreviations.....	1 - 4
Accessories.....	7 - 7
Acoustic factor.....	9 - 6
Acronyms.....	1 - 4
Admin	
rights.....	8 - 21
Anterior.....	6 - 3
Anterior BFA.....	6 - 2
Anterior BFS.....	6 - 2
Anterior chamber depth.....	9 - 6
Anterior chamber map	
Apex Plane.....	6 - 4
Endothelium.....	6 - 4
Epithelium.....	6 - 4
Normal Depth.....	6 - 4
Anterior chamber maps.....	6 - 4
Apex Plane.....	6 - 4
Application.....	8 - 1
Astig.....	9 - 6
Astig power.....	9 - 6
Authorized Service Technician.....	2 - 3
Authorized Technolas Perfect Vision Represen- tative.....	2 - 3
Axial power map.....	6 - 3
Anterior.....	6 - 3
Keratometric.....	6 - 3
Posterior.....	6 - 3
Total.....	6 - 3

Æ

Ø (pupil).....	9 - 6
----------------	-------

B

Brands.....	A - 3
-------------	-------

C

Calibration.....	5 - 1
Calibration check	
performing	8 - 33, 8 - 38
Calibration check tool.....	7 - 7
Changing	
user.....	8 - 4
Changing colors.....	9 - 7
Changing map settings.....	9 - 12
Changing map types.....	9 - 8
Circles.....	9 - 12
Closing	
software.....	8 - 4
Color scale type	
Automatic.....	9 - 6
Custom.....	9 - 6
Default.....	9 - 6
Normal Band.....	9 - 6
Comments.....	9 - 6
Component.....	7 - 1
Customer's	
obligations.....	2 - 2
responsibility.....	2 - 2
Customer Service.....	A - 2
Cylinder convention.....	9 - 54

D

Default	
treatment value.....	8 - 7
Displaying OD and OS measurements.....	9 - 11
Disposal of electronic equipment.....	2 - 10
Draw circles at.....	9 - 12

E

Editing	
nomogram.....	8 - 6
Electrical connection.....	12 - 2
Electronic equipment, disposal.....	2 - 10

Elevation map		
Anterior BFA.....	6 - 2	
Anterior Sph/Cyl.....	6 - 2	
Posterior BFA.....	6 - 2	
Posterior BFS.....	6 - 2	
Posterior Sph/Cyl.....	6 - 2	
Elevation maps.....	6 - 2	
Anterior BFS.....	6 - 2	
EMC Declaration.....	14 - 7	
Endothelium.....	6 - 4	
Epithelium.....	6 - 4	
Explanation of symbols.....	A - 1	
F		
Flap thickness.....	8 - 7	
Flat axis.....	9 - 6	
G		
General settings.....	9 - 12	
Acoustic factor.....	9 - 12	
Radius of curvature.....	9 - 12	
Grid.....	9 - 12	
H		
Hardware		
operating.....	7 - 9	
Hardware component.....	7 - 1	
Hardware messages		
visual.....	2 - 6	
Human interface devices.....	7 - 8	
I		
Information page.....	8 - 44	
Intended use.....	2 - 1	
K		
Kappa.....	9 - 6	
Keratometric.....	6 - 3, 6 - 4	
Keratometry.....	6 - 2	
K-Labels.....	9 - 12	
K-Lines.....	9 - 12	
K-Lines & K-Labels.....	9 - 12	
K-Reading.....	6 - 2	
L		
Labels.....	2 - 6	
Logging in		
software.....	8 - 1	
Logging off		
software.....	8 - 4	
M		
Map.....	9 - 6	
Maps		
Anterior chamber.....	6 - 4	
Axial power.....	6 - 3	
Elevation.....	6 - 2	
Mean power.....	6 - 3	
Pachymetry.....	6 - 3	
Tangential power.....	6 - 4	
Max.....	9 - 6	
Mean power.....	9 - 6	
Mean power map.....	6 - 3	
Measurement		
closing.....	9 - 27, 9 - 28	
comparing.....	9 - 29, 9 - 31	
deleting.....	9 - 29	
exporting	9 - 34	
performing	9 - 1, 9 - 2, 9 - 14	
printing	9 - 38, 9 - 39, 9 - 56	
viewing	9 - 6	
Measurement results page.....	9 - 6	
Min.....	9 - 6	
N		
Network connection.....	7 - 8	
Nomogram		
adjustment.....	8 - 6	
parameters.....	8 - 5	
user-defined.....	8 - 5, 8 - 6	
Normal Depth.....	6 - 4	
Number of colors.....	9 - 6	
O		
O.....	9 - 6	
Obligations, customer.....	2 - 2	
Operating conditions.....	12 - 2	

Operating the hardware.....	7 - 9	Product names.....	A - 3
Operating time.....	12 - 2	Proscan treatment.....	6 - 13
Optical zone diameter.....	8 - 7	PTK treatment.....	6 - 13
Optional component.....	1 - 3	Pupil.....	9 - 12
Options		Pupil diameter.....	9 - 6
default treatment values.....	8 - 7	R	
user-defined nomogram.....	8 - 6	r.....	9 - 6
ORBSCAN 3		Radius of curvature.....	9 - 12
software messages.....	10 - 4	Raw.....	9 - 6
Overlay settings.....	9 - 12	Raw data.....	9 - 6
Circles.....	9 - 12	Refraction.....	6 - 2
Draw circles at.....	9 - 12	adjusting.....	8 - 5
Grid.....	9 - 12	S	
K-Labels.....	9 - 12	Safety	
K-Lines.....	9 - 12	considerations.....	2 - 2
K-Lines & K-Labels.....	9 - 12	Hazards.....	2 - 4
Thinnest point.....	9 - 12	signals.....	2 - 6
P		Safety instructions	
Pachymetry.....	6 - 2	explanation.....	A - 1
Pachymetry averages.....	9 - 12	Service.....	5 - 1
Pachymetry map.....	6 - 3	Setting	
Pachymetry settings.....	9 - 12	default treatment values.....	8 - 7
Draw circles at.....	9 - 12	flap thickness.....	8 - 7
Pachymetry averages.....	9 - 12	optical zone diameter.....	8 - 7
Pages		user-defined nomogram.....	8 - 6
Information.....	8 - 44	Shutting down.....	8 - 4
Measurement results.....	9 - 6	Signs.....	2 - 6
Parameter		Sim-K's.....	9 - 6
nomogram value.....	8 - 6	Software.....	8 - 1
Patient		logging in.....	8 - 1
adding.....	8 - 15	logging out.....	8 - 4
deleting.....	8 - 20	Software messages	
editing.....	8 - 18	ORBSCAN 3.....	10 - 4
selecting.....	8 - 8, 8 - 10	treatment calculation.....	10 - 8
Patient administration.....	8 - 8	troubleshooting chart.....	10 - 1
Personnel.....	2 - 3	visual.....	2 - 6
Posterior.....	6 - 3	Standard optical zone diameter.....	8 - 7
Posterior BFA.....	6 - 2	Stickers.....	2 - 6
Posterior BFS.....	6 - 2	Storage conditions.....	12 - 2
Power.....	6 - 2	SUPRACOR optical zone diameter.....	8 - 7

SUPRACOR treatment.....	6 - 13	Troubleshooting chart	
Switching between map and raw data.....	9 - 13	software messages.....	10 - 1
Symbols		user control screen.....	10 - 1
in this user manual.....	A - 1	U	
T		USB memory stick.....	7 - 8
Tangential power map.....	6 - 4	Use.....	2 - 1
Keratometric.....	6 - 4	User	
Technical support.....	A - 2	changing.....	8 - 4
Thinnest point.....	9 - 6, 9 - 12	profiles.....	8 - 20
Total.....	6 - 3	rights.....	8 - 20
Transport.....	3 - 1	User (Ophthalmologist).....	2 - 3
visual inspection.....	3 - 1	User administration.....	8 - 20
Transport conditions.....	12 - 2	User-defined nomogram.....	8 - 5, 8 - 6
Treatment		W	
functionality.....	6 - 12	WEEE directive.....	2 - 10
overview.....	6 - 12	White-to-white.....	9 - 6
Proscan.....	6 - 13	Z	
SUPRACOR.....	6 - 13	Z.....	9 - 6
ZYOPTIX HD.....	6 - 13	Zooming maps.....	9 - 7
Treatment calculation		ZYOPTIX HD treatment.....	6 - 13
software messages.....	10 - 8		
Treatment modes			
PTK.....	6 - 13		

Appendix

The Appendix section contains the following documents:

- Declaration of Conformity as per Directive 93/42/EEC Annex II
- Declaration of Conformity as per Directive 2011/65/EU
- EMC Declaration
- User information TPV Medical Device Network
- National safety regulations*)

**) File the applicable national safety requirements of your country in the respective national language in the Appendix section of this user manual.*

**A Declaration of Conformity as per Directive 93/42/EEC
Annex II**

EU DECLARATION OF CONFORMITY

Technolas Perfect Vision GmbH
Messerschmittstraße 1 + 3
80992 München
Germany

hereby declares that below described medical devices for ophthalmological applications meet the requirements of the


Directive 2011/65/EU

The object of the declaration described below is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

Name of Product: **Zyoptix Diagnostic Workstation 3**
SKU: **50007000**

This declaration of conformity is issued under the sole responsibility of Technolas Perfect Vision GmbH.

Munich, August 11, 2014

Rupert VEITH
Plant Manager/ VP Operations

C EMC Declaration

1. Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The ZDW 3 is intended for use in the electromagnetic environment specified below. The customer or user of the ZDW 3 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The ZDW 3 uses RF energy only for its internal function. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The ZDW 3 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC61000-3-3	Complies	

2. Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The ZDW 3 is intended for use in the electromagnetic environment specified below. The customer or user of the ZDW 3 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ZDW 3 requires continued operation during power mains interruptions, it is recommended that the ZDW 3 be powered from an uninterruptible power supply or a battery.
Magnetic field at Power frequency (50/60Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the AC mains voltage prior to the application of the test level.

3. Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The ZDW 3 is intended for use in the electromagnetic environment specified below. The customer or user of the ZDW 3 should assure that it is used in such an environment.

Portable and mobile RF communication equipment should be used no closer to any part of the ZDW 3, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff}	Recommended separation distance $d = 1.2\sqrt{P}$ (Conducted)
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz (Radiated) $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz (Radiated)

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.

In the vicinity of devices bearing the label shown below interference is possible.

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ZDW 3 is used exceeds the applicable RF compliance level above, the ZDW 3 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ZDW 3.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

4. Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the ZDW 3

The ZDW 3 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the ZDW 3 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF-communications equipment (transmitters) and the ZDW 3 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter meters (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

D User information TPV Medical Device Network

SKU: 70005910 * Version 2 * UM-100006805 GB

User information TPV Medical Device Network

Version 4

Copyright © 2014 Technolas Perfect Vision GmbH

1	Introduction	
1.1	Overview on the data flow.....	5
2	Network solutions	
2.1	Direct connection.....	6
2.2	Patched connections to the clinic network.....	7
2.3	Patched connections of the clinic network and VLAN switch of the clinic.....	8
3	IP addresses	
4	Requirements and possible risks	

1 Introduction

This user information explains:

- how to integrate ZYOPTIX® Diagnostic Workstation 3 (ZDW 3), TECHNOLAS® TENEO™ 317 (TENEO 317), ZYOPTIX® Diagnostic Workstation (ZDW), and the 'Database Server Unit' (DSU) into an IT infrastructure network in your practice or clinic
- and provides the information necessary for running and servicing the IT infrastructure network.

The document specifies the requirements stated in IEC 80001-1:2010; German version EN 80001-1:2011 Application of risk management for IT-networks incorporating medical devices Part 1: Roles, responsibilities and activities along with the potential risks resulting from non-conformance with those requirements.

Abbreviations

DSU	Database Server Unit
IT	Information Technology
IP	Internet Protocol
TENEO	TECHNOLAS TENEO 317
TPV	Technolas Perfect Vision
VLAN	Virtual Local Area Network
WAN	Wide Area Network
WLAN	Wireless Local Area Network
ZDW	ZYOPTIX Diagnostic Workstation
ZDW 3	ZYOPTIX Diagnostic Workstation 3

1.1 Overview on the data flow

Fig. 1 illustrates the data flow between TPV medical devices and DSU.


Fig. 1: Data flow between TPV medical devices and DSU

2 Network solutions

The following sections give an overview about the three possible solutions for the network configuration. The solutions below may be combined in one network, as far as the solution fulfills the requirements according to Chapter 4 'Requirements and possible risks' on page 11. If for any reason none of the solution is possible, contact an Authorized Service Technician.

2.1 Direct connection

The medical devices are connected directly using the 8-port ethernet switch provided by TPV. All devices are completely separated from the clinic network. No network infrastructure of the clinic is used. The network switch is placed near the server.

Advantages

- No preparation by clinic necessary.

Disadvantages

- No long distances between devices possible.
- Drilling between the rooms is necessary.
- Loose cables may cause tripping hazard.

Network solutions


Fig. 2: Direct connections

Direct connection ZDW 3 to printer:

It is also possible to directly connect the printer to the dedicated LAN port of the ZDW 3.

2.2 Patched connections to the clinic network

The medical devices are connected to network outlets of the clinic. The combining element is an 8-port ethernet switch provided by TPV that is placed near the 'Database Server Unit'. The clinic patches the connections by hardware.

Advantages

- No loose cables.
- No drilling between the rooms is necessary.
- Long distances between devices possible.

Disadvantages

- Preparation by clinic necessary.
- Clinic IT personnel must check clinic HW connections on regular intervals to still fulfill the requirements according to § Chapter 4 'Requirements and possible risks' on page 11.

Patched connections to the clinic network: Example TENE0 317


Fig. 3: Using clinic outlets and TPV switch

A clinic hardware connection is also possible for all other TPV devices, such as, ZDW 3 and ZDW, separately as well as simultaneously.

2.3 Patched connections of the clinic network and VLAN switch of the clinic

The medical devices are connected to network outlets of the clinic. There is no network switch provided by TPV. The combining element is a configured VLAN switch of the clinic (in the server room). All medical devices must be located in one VLAN.

If you choose a solution that includes VPN, be sure to fulfill all requirements according to Chapter 4 'Requirements and possible risks' on page 11.

Advantages

- No loose cables.
- No drilling between the rooms is necessary.
- Long distances between devices possible, even between separate buildings and over VPN.

Disadvantages

- Hardware preparation by clinic necessary.
- Configuration of the VLAN switch by clinic necessary.
- Clinic IT personnel must check clinic HW connections and Configuration of VLAN switch on regular intervals to still fulfill the requirements according to Chapter 4 'Requirements and possible risks' on page 11.

Patched connections of the clinic network and VLAN switch of the clinic


Fig. 4: Using a clinic VLAN network

3 IP addresses

Fig. 5 gives an overview of the IP addresses of each connected device for the standard configuration. The IP addresses are predefined by TPV and must not be changed.


Fig. 5: IP addresses standard configuration

When several medical devices of the same type are in one medical device network, the last number of the IP address must be increased by one. For example, if two TECHNOLAS TENE0 317 systems are within the medical network, they will have the following IP configuration:

First TECHNOLAS TENE0 317	IP = 10.10.254. 137 / 255.255.255.0
Second TECHNOLAS TENE0 317	IP = 10.10.254. 138 / 255.255.255.0
First ZYOPTIX Diagnostic Workstation 3	IP = 10.10.254. 150 / 255.255.255.0
Second ZYOPTIX Diagnostic Workstation 3	IP = 10.10.254. 151 / 255.255.255.0

4 Requirements and possible risks

The responsible organization must ensure that the IT infrastructure in which the medical devices will be integrated meets the requirements listed below. If the IT infrastructure does not yet comply with all conditions, prepare the infrastructure as needed.

Some of the requirements are only relevant for certain network solutions.

Location

- The 'Database Server Unit' must be located 1.5 m outside the patient environment.
- The network printer must be located 1.5 m outside the patient environment. This is also applicable for a printer directly connected to the ZDW 3, this means, it is not allowed to position the printer in the center table of the ZDW 3.

DSU and network printer do only provide protection against electrical hazards and EMC according to the applicable industry standard for devices of its kind.

Hardware requirements

Requirement

RJ45 network cables complying with CAT 5 standard.

If ethernet switches are used, bandwidth of the switches must be at least 100 Mbit.

Note: Recommended is Gigabit – Ethernet (GigE, 1000 Mbit).

Risk if requirement is not met

No or slow connection to DSU.

No or slow connection to DSU.

In use cases where DB communication is involved the GUI reacts slowly.

Network conditions

Requirement

The TPV medical devices must only be located in a private network or VLAN network separated from the clinical network.

Only equipment provided by TPV in connection with the DSU is allowed in the private network or VLAN network (except for managed switches). Older types of equipment are not allowed in the medical device network.

Risk if requirement is not met

Potential virus infections, hacker attacks, denial of service attacks.

Note: Medical devices do not receive regular updates of the operating system since such changes require new reva-
lidation of the medical device.

Potential virus infections or hacker attacks.

Note: Medical devices do not receive regular updates of the operating system since such changes require new reva-
lidation of the medical device.

Requirement

In addition to TPV equipment the usage of managed switches is allowed.

The VLAN must be configured on the ports where TPV medical devices are connected. (Either use a separate switch and directly connect the medical devices to it, or use your existing managed switches and configure a VLAN on them.)

VLAN may use existing patch panel technology of the clinic.

The medical devices must not have internet access.

It must be possible to use IP addresses / IP network masks that result in a separated IP network.

Risk if requirement is not met

If programmable routers instead of managed switches are used: Low risk of virus infection.

Wrongly configured VLAN ports will probably result in no connection between devices.

Wrongly patched VLAN ports will probably result in no connection between devices.

Potential virus infections, hacker attacks, denial of service attacks.

Note: Medical devices do not receive regular updates of the operating system since such changes require new revalidation of the medical device.

If there is no separated IP network, some devices may gain internet access.

VPN requirements

Requirement

If VPN gateways are used, only special dedicated hardware devices are allowed. No software solution is allowed on both sides of the VPN connection.

Strong encryption must be used. For example, IP-Sec. PPTP is not allowed.

The bandwidth for the VPN tunnel must be at least 1 MBit/s (recommended 10 MBit/s).

Only routing between the connected VLANs is allowed (no connection to the clinic network or the internet allowed).

The latest vendor-provided firmware must be used.

Risk if requirement is not met

If unsafe IT equipment is used, virus infections or hacker attacks are possible.

If no or unsafe encryption is used, data theft or man-in-the-middle attacks are possible.

■ No or slow connection to DSU.

■ In use cases where communication is involved the GUI reacts slowly.

Potential risk of virus infections, hacker attacks, denial of service attacks.

Obsolete firmware increases the risk of security leaks.

Further requirements

Requirement

Wireless technologies, such as WLAN are currently not allowed in the medical device network.

Strong passwords for all devices.

The medical device network must be administered and serviced according to *DIN EN 80001-1:2011*.

All network configurations including (but not limited to) VPN configurations, managed switch configurations, patch panel changes, and so on, must be under control and regular inspection by clinic IT.

Risk if requirement is not met

- Data security may be compromised if insufficient encryption is used.
- Slow connection and quality degradation.

Weak passwords increase the risk of unauthorized access to the system.

Risks from network configurations, configuration changes and their change control may not be handled adequately.

Changes in network settings or configurations may lead to failure of one or more of the requirements in above chapters.

24. 6. 2016

[Signature]

URNr. H 2060 / 2016

I hereby certify that the signature above is true, acknowledged in my presence, of

Mr. Rupert **Veith**,
born on March 24th, 1970,
business address: Messerschmittstr. 1+3, 80992 Munich, Germany,
personally known to me,

acting on behalf of


Technolas Perfect Vision GmbH

with the seat in Munich.

(address: Messerschmittstr. 1+3, 80992 Munich, Germany)

After careful examination of the commercial register of the local court of Munich made on June 24th, 2016, I hereby certify that the above-mentioned company is registered under section B, No. 178055 and that Mr. Rupert **Veith** holder of a special statutory authority ("Prokura") is authorized to represent it solely.

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