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1. Spicer EK, Horton R, Bloem L, et al. Isolation of cDNA clones coding for human tissue factor: Primary structure of the protein and cDNA. Proc Natl Acad Sci USA. 1987; 84:5148-52.
2. Quick AJ. Hemorrhagic diseases and thrombosis. Philadelphia: Lea and Febiger; 1966.
3. Bader R, Mannucci PM, Tripodi A, et al. Multicentric evaluation of a new PT reagent based on recombinant human tissue factor and synthetic phospholipids. Thromb Haemost. 1994;71:292-9.
4. Poller L. The Prothrombin Time. WHO/LAB/98.3. 1998.
5. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays. Approved Guideline - Fifth Edition. CLSI document H21-A5 [ISBN 1-56238-657-3]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2008.
6. Llamas P, Santos AB, Outeiriño J, et al. Diagnostic utility of comparing fibrinogen Clauss and prothrombin time derived method. Thromb Res. 2004; 114: 73-4.
7. De Cristofaro R, Landolfi R. Measurement of plasma fibrinogen concentration by the prothrombintime-derived method: applicability and limitations. Blood Coagul Fibrinolysis. 1998; 9: 251-9.
8. Wagner C, Dati F. Fibrinogen. In: Thomas L, ed. Clinical Laboratory Diagnostics. Frankfurt: TH-Books Verlagsgesellschaft, 1998:609-12.
9. Tripodi A, Chantarangkul V, Clerici M, et al. Laboratory control of oral anticoagulant treatment by the INR system in patients with the antiphospholipid syndrome and lupus anticoagulant. Results of a collaborative study involving nine commercial thromoplastins. Br J Haematol. 2001;115:672-8.
10. Tobu M, Iqbal O, Messmore HL, et al. Influence of different anticoagulant agents on fibrinopeptide A generation. Clin Appl Thromb Hemost. 2003;9:273-92.
11. Fenyvesi T, Joerg I, Harenberg J. Influence of Lepirudin, Argatroban, and Melagatran on prothrombin time and additional effect of oral anticoagulation. Clin Chem. 2002;48:1791-4.
12. Van den Besselaar AM, Hoekstra MM, Witteveen E, et al. Influence of Blood Collection Systems on the Prothrombin Time and International Sensitivity Index Determined With Human and Rabbit Thromboplastin Reagents. Am J Clin Pathol. 2007; 127: 724-9.
13. Ansell J, Hirsh J, Hylek E, et al. Pharmacology and management of the vitamin K antagonists. Chest. 2008; 133:160S-198S.

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© 2010 Siemens Healthcare Diagnostics Products GmbH.

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1. Cooper J, Douglas AS. Fibrinogen level as a predictor of mortality in survivors of myocardial infarction. Fibrinolysis 1991; 5: 105-8.
2. Cook NS, Ubben D. Fibrinogen as a major risk factor in cardiovascular disease. Trends Pharmacol Sci. 1990; 11: 444-51.
3. Technical Manual of the American Association of Blood Banks 9th ed. Arlington: AABB; 1985: 417.
4. Clauss A. Gerinnungsphysiologische Schnellmethode zur Bestimmung des Fibrinogens. Acta Haematol. 1957; 17: 237-46.
5. Borgström S. On the prothrombin index in acute affections of the pancreas. Acta Chir Scand. 1945; 90: 419-30.
6. Jacobsson K. Studies on the determination of fibrinogen in human blood plasma. Scand J Clin Lab Invest. 1955; 7 (Suppl 14): 1-54.
7. Clinical and Laboratory Standards Institute (CLSI). Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline - Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400. Wayne, PA 19087-1898 USA, 2008.
8. Morse EE, Panek S, Menga R. Automated fibrinogen determination. Amer J Clin Pathol. 1971; 55: 671-6.
9. CLSI. Procedure for the determination of fibrinogen in plasma; Approved Guideline. 2nd Edition. CLSI Publication H30-A2. Villanova, PA. 2001.

10. Wagner C, Dati F. Fibrinogen. In: Thomas L, ed. Clinical Laboratory Diagnostics. Frankfurt/Main, Germany: TH-Books Verlagsgesellschaft mbH 1998: 609-12.
11. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory*; Approved Guideline – Third Edition. CLSI document C28-A3 (ISBN 1-56238-682-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2008.
12. Okuno T, Selenko V. Plasma fibrinogen determination by automated thrombin time. *Amer J Med Tech.* 1972; 38: 196-201.

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© 2008 Siemens Healthcare Diagnostics Products GmbH.

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SIEMENS

LA Control High

Lupus Anticoagulant Quality Control Plasma, High Level
for use with LA 1 Screening Reagent and LA 2 Confirmation Reagent (DRVVT)

Intended use

LA Control High is a high positive control plasma for use in Lupus Anticoagulant clotting test assays, specifically LA 1 Screening Reagent (DRVVT) and LA 2 Confirmation Reagent (DRVVT) from Siemens Healthcare Diagnostics on automated coagulation analyzers.

Summary and explanation

"Lupus Anticoagulant" (LA) was a term used in the early 1970s to describe a condition whereby non-specific circulatory coagulation inhibitors were found to be present particularly in cases of systemic lupus erythematosus (SLE)¹. Currently LA is recognized as being associated with a variety of thrombotic syndromes².

Reagents

For *in-vitro* diagnostic use.

Caution: Treat as potentially infectious.

This product has been collected from lupus anticoagulant patients under FDA approved conditions. All human plasmas prepared for lyophilization have been found negative when tested for antibodies to HIV-1/2 and HCV. Non-reactive for hepatitis B surface antigen (HBsAg) by FDA required tests.

Pack contents:

LA Control High, No. OQWD 11

6 vials each for 1 mL

6 vials lyophilized platelet poor human plasma from LA positive and normal donors, buffers and stabilizers, sodium azide.

When disposing of sodium azide, always flush with large volumes of water to avoid the possibility of an explosive residue forming in metal plumbing.

Additional materials required but not provided

Purified water, USP or equivalent

Micropipet, 1.0 mL

Storage Conditions

Store at +2 - +8 °C

The lyophilized plasma should remain stable at least until the expiry date stated on the vial when stored at +2 - +8 °C. Following reconstitution, plasma should be used within 4 hours at +20 °C (open vial) or 8 hours at +2 - +8 °C (closed vial). Reconstituted plasma may be frozen for 1 week at -20 °C.

Test procedure

NB: all standard precautionary practices (hygiene, spill procedures, disposal of biological materials, etc) must be observed.

1. Reconstitute with the volume of purified water stated on the vial label.
2. Mix gently, do not shake.
3. Allow to stand for at least 15 minutes before use.
4. Prepare Lupus Anticoagulant testing reagents according to the manufacturers specifications.
5. Test undiluted QC plasma at the same time as the patient samples.
6. Record results of clotting times and/or ratios on Quality Control charts as required.

Performance characteristics

LA Control High should give a ratio between 1.7 and 2.4 using a BCT* instrument with LA 1 Screening Reagent and LA 2 Confirmation Reagent. Other reagent / instrument combinations may yield different results.

A new range based on mean $\pm$ 2 SD should be established with each lot of controls and reagents in each test system.

Within-run and total precision were calculated for the LA Control High by testing twice daily over twenty days using an automated coagulation analyzer. Within-run and total precision for both volumes of control were less than 3% with LA 1 Screening Reagent and LA 2 Confirmation Reagent.

Bibliography

1. Feinstein D.I., Rapaport, S.P. Acquired inhibitors of blood coagulation. Progress in hemostasis and thrombosis, Vol 1. Ed. Spaet T., New York, Grune and Stratton, pp 75-95, 1972
2. Schleider M.A. et al. A clinical study of the lupus anticoagulant. Blood 48:499-509, 1976

* BCT is a trademark of Siemens Healthcare Diagnostics.

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

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
SIEMENS

LA Control Low

Lupus Anticoagulant Quality Control Plasma, Low Level

for use with LA 1 Screening Reagent and LA 2 Confirmation Reagent (DRVVT)

Intended use

LA Control Low is a low positive control plasma for use in Lupus Anticoagulant clotting test assays, specifically LA 1 Screening Reagent (DRVVT) and LA 2 Confirmation Reagent (DRVVT) from Siemens Healthcare Diagnostics on automated coagulation analyzers.

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Reagents

For *in-vitro* diagnostic use.

Caution! Treat as potentially infectious!

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Performance characteristics

LA Control Low should give a ratio between 1.35 and 1.65 using a BCT* instrument with LA 1 Screening Reagent and LA 2 Confirmation Reagent. Other reagent / instrument combinations may yield different results.

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Within-run and total precision were calculated for the LA Control Low by testing twice daily over twenty days using an automated coagulation analyzer. Within-run and total precision for both volumes of control were less than 3 % with LA 1 Screening Reagent and LA 2 Confirmation Reagent.

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Dade® DIMERTEST® Latex Assay

INTENDED USE

The Dade® DIMERTEST® Latex Assay is intended for the rapid qualitative or semiquantitative evaluation of circulating derivatives of cross-linked fibrin degradation products (XL-FDP) in human plasma.

SUMMARY AND EXPLANATION OF THE TEST

During blood coagulation, fibrinogen is converted to fibrin by the activation of thrombin. The resulting fibrin monomers polymerize to form a soluble gel of noncross-linked fibrin. This fibrin gel is then converted to cross-linked fibrin by thrombin activated Factor XIII to form an insoluble fibrin clot. Production of plasmin, the major clot-lysing enzyme, is triggered when a fibrin clot is formed. Fibrinogen and fibrin are both cleaved by the fibrinolytic enzyme plasmin to yield degradation products, but only degradation products from cross-linked fibrin contain D-dimer^{1,2}. Therefore, cross-linked fibrin degradation products (XL-FDP) are a specific marker of fibrinolysis.

TEST PRINCIPLE

Dade® DIMERTEST® Latex is a rapid agglutination assay utilizing latex beads coupled with a highly specific D-dimer monoclonal antibody. XL-FDP present in a plasma sample bind to the coated latex beads, which results in visible agglutination occurring when the concentration of D-dimer is above the threshold of detection of the assay.

REAGENTS

Composition

1. DIMERTEST® Latex Beads: a 0.83 % suspension of latex particles coated with murine anti-D-dimer monoclonal antibody, 10 mg/mL BSA and 0.1 % sodium azide.
2. DIMERTEST® Positive Control: a solution containing purified human D-dimer fragment, 5 mg/mL BSA and 0.1 % sodium azide.
3. DIMERTEST® Negative Control: a buffer solution containing 5 mg/mL BSA and 0.1 % sodium azide.
4. DIMERTEST® Buffer: a 10 mM phosphate buffer solution containing 0.1 % sodium azide.

Warnings and Precautions

1. For In-Vitro Diagnostic Use Only.
2. Harmful if swallowed (R 22). If swallowed, seek medical advice immediately and show this container or label (S 46).
3. CAUTION: All reagents in the Dade® DIMERTEST® Latex Assay contain sodium azide (0.1 %) as preservative. Do not ingest or allow to contact skin or mucous membranes. Sodium azide may form explosive azides in metal plumbing. Use proper disposal procedures.
4. CAUTION! POTENTIAL BIOHAZARD
The Positive Control in the Dade® DIMERTEST® Latex Assay contains components of human origin. Each individual blood donation intended for the production of this reagent is tested for Hepatitis B surface antigen, anti-HCV, anti-HIV1 and anti-HIV2. Only donations with negative findings are employed. As complete absence of infectious agents can never be assured, all materials derived from human blood should be treated as potentially infectious and handled with due care following the precautions recommended for biohazardous material.

Storage and Stability

Storage: Store at 2 to 8 °C. DO NOT FREEZE.
Stability: Refer to outer package and vial labels for expiration date.

Indication of Reagent Deterioration

Reagent deterioration is indicated by failure of the Latex Beads to agglutinate with the Positive Control, agglutination with the Negative Control, or evidence of microbial contamination.

SPECIMEN COLLECTION AND PREPARATION

Plasma prepared from whole blood anticoagulated with sodium citrate is recommended. The use of EDTA and heparin will result in an increased level of false positive reactions. After separation of the plasma by centrifugation (1500 x g for 15 minutes), specimens may be tested directly for the presence of XL-FDP. Defibrination of the plasma is not recommended.

Plasma storage/stability: -20 °C: 2 weeks

Thaw frozen specimens rapidly at 37 °C and centrifuge before testing. Refer to CLSI/NCCLS publication H21-A5 for further instructions on specimen collection, handling and storage⁴.

PROCEDURE

Materials Provided

Cat No. 84233-60	Dade® Dimerest® Latex Assay 60 tests
Latex Beads	1 x 2.0 mL White Cap
Negative Control	1 x 0.6 mL Black Cap
Positive Control	1 x 0.6 mL Yellow Cap
Buffer	1 x 20.0 mL
Test Cards	x 10:8-well, for agglutination reaction
Stirrers	x 60: for mixing

Instructions for Use

Cat No. 84233-61L	Dade® Dimerest® Latex Vial 2.0 mL
Latex Beads	1 x 2.0 mL White cap

Instructions for Use

Materials Required but not Provided
Precision pipettes and tips - 20 µL and 100 µL
Plastic test tubes and rack
Stopwatch or timing device
Disposable gloves
Tissue (for wiping dropper bottle tips)

Important!

- Equilibrate reagents to room temperature (20 to 25 °C) before use.
- Latex Beads should be mixed by inversion immediately prior to use.
- Prior to each use, the dropper bottle tips must be wiped dry with a tissue.
- Dropper bottles must be held vertically when dispensing drops of reagent.

For correct drop delivery...

Wipe tip dry prior to use

Qualitative Method

1. Mark (or make note of) positions on the test slide for specimens and, as needed, for positive and negative controls.
2. Hold the Latex Beads dropper bottle vertically and place one drop of the reagent within a well on a Test Card.
3. AVOID touching the surface of the Test Card.
4. Accurately pipette 20 µL of undiluted plasma or add one drop of control solution inside the same well next to the drop of Latex Beads.
5. Mix the Latex Beads and sample with a stirrer until the Latex is uniformly distributed.
6. Rock the Test Card gently by hand for exactly 3 minutes.
7. At exactly 3 minutes, check for agglutination under a strong light source.

Note: If test reading is delayed beyond 3 minutes, the latex suspension may dry out giving a false agglutination pattern. If this is suspected, the specimen must be re-tested.

Semiquantitative Method

1. Prepare serial dilutions of the test plasma with Buffer as follows:
1:2 dilution 100 µL plasma plus 100 µL Buffer solution
1:4 dilution 100 µL 1:2 dilution plus 100 µL Buffer solution
1:8 dilution 100 µL 1:4 dilution plus 100 µL Buffer solution
2. Test each dilution as described in the qualitative method.

Quality Control

It is recommended that both Positive and Negative Controls be included in each batch of tests to ensure proper functioning of the system. Control solutions should be tested by the same procedures as patient samples. For the qualitative screening procedure, a positive result (agglutination) can be obtained by substituting DIMERTEST® Positive Control for the plasma in steps 3-6. Conversely a negative result (no agglutination) can be obtained by substituting DIMERTEST® Negative Control for the plasma in steps 3-6. When performing the semiquantitative procedure, it is recommended that testing of serial dilutions of the Positive Control be conducted. DIMERTEST® Latex Positive Control consists of a solution of human D-dimer at a level of approximately 0.80 mg/L (800 ng/mL).

RESULTS

A. Qualitative Assay

For the qualitative assay protocol, the following pattern of results should be obtained:

Undiluted Plasma	D-dimer (XL-FDP) concentration
Negative	Less than 0.20 mg/L (200 ng/mL)
Positive	Greater than 0.20 mg/L (200 ng/mL)

Note: All values in mg/L (ng/mL) are approximate

B. Semiquantitative Assay

Approximate levels of XL-FDP, containing the D-dimer domain, for specimen dilutions are shown in Table 1. As with all semiquantitative tests, some variability in dose-response can be expected.

Table 1

Approximate Range D-dimer (XL-FDP) mg/L (ng/mL)	Sample Dilution			
	Undiluted	1:2	1:4	1:8
< 0.20 (< 200)	-	-	-	-
0.20 - 0.40 (200 - 400)	+	-	-	-
0.40 - 0.80 (400 - 800)	+	+	-	-
0.80 - 1.60 (800 - 1600)	+	+	+	-
1.60 - 3.20* (1600 - 3200*)	+	+	+	+

*+ = agglutination, - = no agglutination
Levels of XL-FDP greater than 3.2 mg/L (3200 ng/mL) can be estimated by further dilutions beyond 1:8.

LIMITATIONS OF THE PROCEDURE

Clinical diagnosis should not be based on the result of the Dade® DIMERTEST® Latex Assay alone. Clinical signs and other relevant test information should be included in the diagnostic decision.

EXPECTED VALUES

A positive result, indicating active fibrinolysis, should be obtained with the Dade® DIMERTEST® Latex Assay when XL-FDP (D-dimer) levels are at or greater than approximately 0.20 mg/L (200 ng/mL). Plasma specimens from normal subjects are expected to give negative results because their plasma XL-FDP concentrations are typically less than 0.20 mg/L (200 ng/mL). Due to many variables that may affect results, each laboratory should establish its own normal range.

Elevated levels of XL-FDP (containing the D-dimer domain) have been demonstrated in patients by a combination of immunoprecipitation and gel electrophoresis techniques^{1,2,3,6}. Monoclonal antibodies allow the specific detection of the D-dimer domain⁷. Monoclonal antibody based D-dimer assays are of diagnostic value in disseminated intravascular coagulation (DIC) and acute vascular diseases, including pulmonary embolism (PE) and deep venous thrombosis (DVT), conditions that are difficult to detect reliably by clinical examination⁸.

The amount of XL-FDP detected in a specimen will depend on several interrelated factors in vivo, such as the severity of the thrombotic episode, the rate of cross-linked fibrin formation, and the time elapsed after the thrombotic event until blood is drawn from the patient.

Elevated levels of XL-FDP as an indication of reactive fibrinolysis have also been reported in surgery, trauma, sickle cell disease, liver disease, severe infection, sepsis, inflammation, and malignancy^{10,11}. D-dimer levels also rise during normal pregnancy but very high levels are associated with complications¹².

The Dade® DIMERTEST® Latex Assay does not cross-react with fibrinogen, factor XIIIa cross-linked fibrinogen¹³, or fibrinogen degradation products¹⁴.

SPECIFIC PERFORMANCE CHARACTERISTICS

Plasma from one hundred and seventy (170) apparently healthy, voluntary blood donors was tested using the Dade® DIMERTEST® Latex Assay. A negative result was obtained for

ne hundred and sixty-two (162) of the samples. This equates to a specificity of 95.3% (162/170).

One hundred and forty-five (145) plasma samples from patients judged to be suffering from, or having a high probability for thrombotic episode, were tested by the Dade[®] DIMERTEST[®] Latex Assay and another agglutination reference method. The correlation coefficient was $r = 0.94$ and the regression equation was $y = 1.19x$. Intra-assay (within run) reproducibility was determined for 10 replicates of 3 plasma samples that contained different levels of XL-FDP. The results were equivalent for all replicates. Inter-assay (run-to-run) reproducibility was determined using 10 plasma samples with XL-FDP titers ranging from 1 to 16. In 10 runs, the replicates of these specimens did not vary by more than one titer.

In an anticoagulant study of 50 parallel citrated, EDTA and heparin plasma samples, the correlation between the titers obtained with the Dade[®] DIMERTEST[®] Latex Assay and the expected titers (based on ELISA XL-FDP values) was $r = 0.91$ for citrated samples, $r = 0.73$ for EDTA samples and $r = 0.78$ for heparin samples. Citrate is the anticoagulant of choice. In a study of samples from patients with rheumatoid arthritis, 17 were found to agglutinate with the Dade[®] DIMERTEST[®] Latex Assay. In all 17 samples, the agglutination could be inhibited by the addition of the D-dimer specific monoclonal antibody DD3B6/22, but not with a non-specific monoclonal antibody of the same subgroup, IgG₃. This suggests that the Dade[®] DIMERTEST[®] Latex Assay is insensitive to rheumatoid factor disturbances. No assay interference was demonstrated with the Dade[®] DIMERTEST[®] Latex Assay with spiked specimens containing potential interferents at the following concentrations:

Bilirubin	0.2 mg/mL	Hemoglobin	5.0 mg/mL
Lipids (triglycerides)	30 mg/mL	Protein (gamma globulin)	0.06 g/mL

Bibliography

- Gaffney, P.J. Distinction between Fibrinogen and Fibrin Degradation Products in Plasma. Clin. Chim. Acta. 65 (1): 109-115; 1975.
- Lane, D.A. et al. Characterisation of Serum Fibrinogen and Fibrin Fragments Produced During Disseminated Intravascular Coagulation. Br. J. Haematol. 40 (4): 609-615; 1978.
- Whitaker, A.N. et al. Identification of D-dimer-E complex in Disseminated Intravascular Coagulation. Thromb. Res. 18 (3-4): 453-459; 1980.
- CLSI/NCCLS Publication H21-A5 - Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline - Fifth Edition; 2008
- Graeff, H. et al. Detection and Relevance of Crosslinked Fibrin Derivatives in Blood. Semin. Thromb. Hemost. 8 (1): 57-68; 1982.
- Yoshioka, K. et al. Distinction between Fibrinogen and Fibrin Products Produced during Disseminated Intravascular Coagulation in Childhood. Eur. J. Pediatr. 138 (1): 46-48; 1982.
- Ryatt, D.B. et al. An Immunoassay for Human D-dimer using Monoclonal Antibodies. Thromb. Res. 31 (6): 767-778; 1983.
- Elms, M.J. et al. Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma using Monoclonal Antibody-Coated Latex Particles. Am. J. Clin. Pathol. 85 (3): 360-364; 1986.
- Whitaker, A.N. et al. Measurement of Cross-Linked Fibrin Derivatives in Plasma: an Immunoassay using Monoclonal Antibodies. J. Clin. Pathol. 37 (8): 882-887; 1984.
- Hunt, F.A. et al. Serum Crosslinked Fibrin (xDP) and Fibrinogen/Fibrin Degradation Products (FDP) in Disorders Associated with Activation of the Coagulation or Fibrinolytic Systems. Br. J. Haematol. 60 (4): 715-722; 1985.
- Smith, R.T. et al. Fibrin Degradation Products in the Postoperative Period. Evaluation of a New Latex Agglutination Method. Am. J. Clin. Pathol. 60 (5): 644-647; 1973.
- Nolan, T.E. et al. Maternal Plasma D-dimer Levels in Normal and Complicated Pregnancies. Obstet. Gynecol. 81 (2): 235-238; 1993.
- Greenberg, C.S. et al. Measurement of Plasma Fibrin D-dimer Levels with the Use of Monoclonal Antibody Coupled to Latex Beads. Am. J. Clin. Pathol. 87 (1): 94-100; 1987.

DIMEREST[®] is a trademark of Sekisui Diagnostics, LLC.

Dade[®] is a trademark of Siemens Healthcare Diagnostics. Manufactured for Siemens Healthcare Diagnostics Products GmbH by Sekisui Diagnostics, LLC, Stamford, CT, USA.

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

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1. Penner, J.A. and Hiss, R.G. Heparin therapy in the 70's. *Univ. Mich. Med. Cent.* 40:57, 1974.
2. Vessler, S. Prevention of venous thromboembolism by low-dose heparin. *Mod. Con. Cardiovascular Dis.*, 45:105, 1976.
3. Triplett, D.A. Laboratory Evaluation of Coagulation. *American Soc. of Clinical Pathologists*, Chicago, 271-313, 1982.
4. Linhardt, R., Grant, A., Cooney, C.L. and Langer, R. Differential anticoagulant activity of heparin fragments prepared using microbial heparinase. *J. Biol. Chem.* 257:7310-7313, 1982.
5. Hull, E.D. and Kingdon, H.S. Use of heparinase to eliminate heparin inhibition in routine coagulation assays. *J. Lab. Clin. Med.* 79:1027-1034, 1972.
6. Baugh, R.F., Deemar, K.A. and Zimmermann, J.J. Heparinase in the Activated Clotting Time Assay: Monitoring heparin independent alterations in coagulation function. *Anesthesia Analgesia*, 72:201-205, 1992.
7. Zimmermann, J.J., Oddie, K., Langer, R. and Cooney, C.L. The release of heparinase from the periplasmic space of *Flavobacterium heparinum* by three-step osmotic shock. *App. Biochem. and Biotech.*, 30:137-148, 1991.
8. Cumming, A.M., Jones, G.R., Worsley, R.T. and Cundell, R.B. *In vitro* neutralization of heparin in plasma prior to the activated partial thromboplastin time test: an assessment of four heparin antagonists and two anion exchange resins. *Thromb. Res.*, 41:43-66, 1986.
9. Thompson, A.R. and Counts, R.B. Removal of heparin and protamine from plasma. *J. Lab. Clin. Med.*, 88:922-929, 1976.
10. Wenz, B. and Burns, E. Rapid removal of heparin from plasma by affinity filtration. *Am. J. of Clin. Path.*, 96:385-389, 1991.
11. Sandset, P.M., Abildgaard, U. and Larsen, M.D. Heparin induces release of extrinsic pathway inhibitor. *Thromb. Res.*, 50:803-818, 1988.
12. Lindahl, A.K., Abildgaard, U. and Staalesen, R. The anticoagulant effect in heparinized blood and plasma resulting from interactions with Extrinsic Pathway Inhibitor. *Throm. Res.*, 64:155-168, 1991.
13. Abildgaard, U., Lindahl, A.K. and Sandset, P.M., Heparin requires both Antithrombin and Extrinsic Pathway Inhibitor for its anticoagulant effect in human blood. *Haemostasis*, 21:254-267, 1991.

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All rights reserved.
Siemens Inc.