0D 21 DC

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of XIAMEN SCIENOX BIOLOGICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Huang Xiuling

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
(Date: Sep. 25, 2020)

证明书查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Xiamen Sciendox Biological Technology Co., Ltd		
Technical Document - Feces Sample liquefaction processor Sciendox 1A		Page 1 of 24
Prepared	Lingling Zheng	General Manager Assistant
Verified	Lingling Zheng	General Manager Assistant
Approved	Qinghua Liao	General Manager

OPERATION DOCUMENTATION

FECES SAMPLE LIQUEFACTION PROCESSOR SCIENDOX 1A

Registration Certificate № XXXX fromXXXXXX

Manufacturer:

Xiamen Sciendox Biological Technology Co., Ltd, People's Republic of China

Address: Room 402, Building 3, №2054 of Wengjiao West Road, Haicang District 361006 Xiamen, Fujian, People's Republic of China

Phone : 0592—5752336 / 5752337 Fax : 0592-5561477

Email : xindao2011@126.com

Web : <http://www.sciendox.com/>

Manufacturing site

Xiamen Sciendox Biological Technology Co., Ltd, China

Room 402, Building 3, №2054 of Wengjiao West Road, Haicang District 361006 Xiamen, Fujian, People's Republic of China

**AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION
/IMPORTER:**

REAL CARE LLC

ADDRESS: 111524, Moscow, Perovskaya street, house 1, building 10, floor 2, room VII room 3

Phone: +7(495)234-57-31

e-mail: lab@int-med.ru

web:<http://www.int-med.ru>

INSTRUCTION VERSION:

Version 1/6.10.2017

Contents

Chapter 1: General instructions.....	5
Chapter 2: Basic Security Rules.....	5
2.1. WARNING AND ATTENTION.....	5
2.1.1. WARNING	5
2.1.2. GROUND CONNECT MARK	5
2.1.3. ATTENTION MARK.....	5
2.2. SECURITY GUIDANCE.....	5
2.2.1. OPENING AND REPAIR	5
2.2.2. BIOLOGICAL SAFETY	5
Chapter 3: Designation. Potential users. Description. Contraindications. Side effects.....	6
3.1. Designation	6
3.2. Potential users	6
3.3. Brief procedure description	6
3.4. Description of target analyte, with indication of qualitative, semi-quantitative or quantitative type	6
3.5. The specific disorder, condition, or risk factor for which MD is intended to detect, identify, or differentiate.....	6
3.6. Type of sample analyzed	6
3.7. Contraindications and side effects.....	6
Chapter 4: Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device.....	6
Chapter 5: Operating conditions	7
Chapter 6: Package contents	7
Chapter 7: Medical device construction and description of sample process flow.	8
7.1. MEDICAL DEVICE CONSTRUCTION.....	8
7.2. MAIN CONTROL PANELS.....	8
7.2.1. ADD FLUID CONTROL INTERFACE	8
7.2.2. MIXING, FILTER CONTROL INTERFACE	9
7.3. CONSUMABLES	10
7.3.1. SAMPLE CONTAINER, APPLIED SOLUTION TUBE (PLASTIC), SAMPLING SPOON.	10
7.3.2. PLASTIC TUBE RACK.....	11
7.3.3. BIO-EKON SAMPLE CONTAINER RACK	11
7.3.4. OTHER COMPONENTS AND REAGENTS	11
7.4. SAMPLE PROCESS FLOW.....	12
7.4.1. SAMPLE PROCESS FLOW DESCRIPTION.....	12
7.4.2. ADD FLUID.....	12
7.4.3. MIXING	13

7.4.4. FILTERING	13
7.4.5. Sample preparation completion.	14
Chapter 8: Installation and connection.....	14
8.1 UNPACKING.....	14
8.2 INSTALLATION	14
8.3 ELECTRICAL CONNECTION.....	14
SALINE FEED COLLECTION	15
Chapter 9: Device operation.....	15
9.1 PREPARATION FOR WORK.....	15
9.2 COLLECTION OF BIOMATERIAL AND DETERMINATION OF THE PROCESSING PARAMETERS	15
9.3 SAMPLE PREPARATION	16
Chapter 10: Daily maintenance	19
10.1 SALINE REPLACEMENT	19
10.2 TUBE MAINTAIN (CLEAN, DRAIN)	19
10.3 PERISTALTIC PUMP TUBE REPLACEMENT	19
10.4 FUSES REPLACEMENT	19
10.5 PUNCTURE NEEDLE.....	19
10.6 MAINTENANCE PERIOD REFERENCE TABLE.....	19
Chapter 11: TECHNICAL SERVICE.....	20
Chapter 12: CURRENT REPAIR	20
Chapter 13: STORAGE AND TRANSPORTATION. SERVICE LIFE.....	20
Chapter 14: WARRANTY OBLIGATIONS	20
Chapter 15: DATA TO DESTRUCTION AND UTILIZATION OF MEDICAL DEVICE.....	21
Chapter 16: Requirements to the Electromagnetic Compatibility.....	21
Chapter 17: Symbols	24

Chapter 1: General instructions

Current operation instruction describes characteristics operating procedures of Feces Sample Liquefaction Processor Sciendox 1A (hereinafter processor Sciendox 1A).

Before device operation please read the Operation Instruction.

Due to continuous product improvement the design of the product can be changed without degrading characteristics.

For clinical laboratory diagnostic.

Chapter 2: Basic Security Rules

2.1. WARNING AND ATTENTION

2.1.1. WARNING

Before Sciendox1A operation you must observe the current operation instruction.

2.1.2. GROUND CONNECT MARK

Must ground connect mark.

Must ground connect through socket.

According to medical equipment safety our device Classify to 1 category B type, common equipment. Device power must have safe ground terminal.

When use in humid environment we must install leakage protector.

2.1.3. ATTENTION MARK

This mark in the operating instruction below highlights important recommendations and cautions.

2.2. SECURITY GUIDANCE

2.2.1. OPENING AND REPAIR

Please do not open the device to observe and modify without permit. In case of the following you should unplug the power supply and ask professional person:

- Power supply cable, signal cable, or socket broken;
- Device in water;
- Device to be dropped and broken, or framework broken;
- Operate according to instruction book, but still cannot work;
- Device character have big change, after exclude sample, reagent, and operate factor. But still cannot work. Then need maintain.

2.2.2. BIOLOGICAL SAFETY

Feces sample maybe a potential pollution source,

The device parts, consumables and waste fluid which contact feces sample are should be disposed according to local law.

Chapter 3: Designation. Potential users. Description. Contraindications. Side effects.

3.1. Designation

Feces sample liquefaction processor Sciendox 1A is intended for automatic preparation of feces samples.

3.2. Potential users

Feces sample liquefaction processor Sciendox 1A may be used in clinical diagnostic laboratories and other medical institutions by experienced medical staff.

Before device operation please read the Operation Instruction.
Work should be carried out in protective gloves.

3.3. Brief procedure description

Sample preparation is carried out by soaking the sample in a physiological solution, mixing and filtering the sample with a two-stage filter. The apparatus brings the sample to the state of the suspension, filtering the conglomerates and passing the elements necessary for the study. Sciendox 1A allows to standardize sample preparation, producing it at a high quality level, not available with a manual method.

For the device operation is necessary: sample collection container, sampling spoon, applied solution tube, saline.

3.4. Description of target analyte, with indication of qualitative, semi-quantitative or quantitative type

Not applicable

3.5. The specific disorder, condition, or risk factor for which MD is intended to detect, identify, or differentiate.

Not applicable

3.6. Type of sample analyzed

Stool specimen

3.7. Contraindications and side effects

There are no contraindications and side effects at this time.

This device does not content medicines, materials of animal and (or) human origin.

The list of standards used by manufacturer is available upon request.

Chapter 4: Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device

Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device are represented in the Table 1.

Table1. Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device

	Value
Productive capacity, samples per hour	90 - 200 samples per hour
Supply voltage	230V
AC frequency	50Hz
Power consumption, W	Not more than 65VA
Dimensions LxWxH, mm	375x285x670
Weight, kg	21.4
Fuse current	2A
Acceleration of centrifugation	not less than 29340 m / s^2 (2992 g)
Setting ranges	Stirring intensity: 0 to 99 times

Chapter 5: Operating conditions

The device is used only indoor.
 The temperature should be between 15°C- 30°C;
 Relative humidity –no more than 70%;
 Air pressure 86.0 kPa to 106.0kPa;
 Avoid moist.
 Please do not move after installed.

Chapter 6: Package contents


1. 1A host 1 piece
2. Bio-Ekon sample container rack 1 piece
3. Plastic tube rack 1 piece
4. Sample container 10 pieces “Start kit”
5. Applied solution tube (plastic) 10 pieces “Start kit”
6. Sampling spoon 10 pieces “Start kit”
7. Bio-Ekon Diluent bottle 1500mL 1 piece
8. Pump pipe (L:150mm, $\phi 3*5\text{mm}$) 2 pieces
9. Apparatus power cord 1 piece
10. Fuse 2 pieces
11. Hex Wrench set 1 suit, including:

- Wrench 1.5 mm
 - Wrench 2 mm
 - Wrench 2.5 mm
 - Wrench 3 mm
 - Wrench 4 mm
 - Wrench 5 mm
 - Wrench 6 mm
 - Wrench 7 mm
12. M5 Phillips screwdriver 1 piece
 13. M3 Phillips screwdriver 1 piece
 14. Operation instruction 1 piece

Chapter 7: Medical device construction and description of sample process flow.

7.1. MEDICAL DEVICE CONSTRUCTION


The device is made of light colored plastics. Functionally the device is divided into two parts. In the left part the sample is diluted. In the right side - mix and filtration. (see Picture1.)

Picture1. Device's main view and functions.

7.2. MAIN CONTROL PANELS


7.2.1. ADD FLUID CONTROL INTERFACE

Picture2. Add fluid control interface

1. **Power indicator light:** After power on ,this light will always bright
2. **Alarm indicator light:** when alarm happen ,this light will twinkling
3. **Digital add fluid indicator:** this indicator shows current fluid volume or settled fluid volume
4. **Value increase key:** Under setting condition, press this key will increase value.
5. **Value reduce key:** Under setting condition, press this key will reduce value.
6. **Setting key:** press this key(6)for parameters change. When the left part of the indicator(3) is twinkling, please press keys(4 or 5) to set volume of the added fluid and press setting key (6) again for confirmation.
7. **Stop key:** under working condition press this key, device will force to stop working.
8. **Start key:** press this key and device start work

7.2.2. MIXING, FILTER CONTROL INTERFACE

Picture 3. Mixing, filter control interface

1. **Power indicator light:** After power on, this light will always bright.
2. **Alarm indicator light:** when alarm happen, this light will twinkling.
3. **Digital parameter indicator:** this indicator shows current saved mixing strength and filter time. The first bit and the second bit indicator show mixing strength "T" shows

mix strength times, the third bit and fourth bit show filter time "S" means filter time.

4. **Value increase key:** Under setting condition, press this key will increase value.

5. **Value reduce key:** Under setting condition, press this key will reduce value.

6. **Setting key:** press this key(6)for parameters change. When the left part of the indicator(3)is twinkling, please press keys(4 or 5) to set sample mixing intensity and press setting key (6) again for confirmation.


When the right part of the indicator(3) is twinkling, please press keys(4 or 5) to set filtration time and press setting key (6) again for confirmation.

7. **Stop key:** under working condition press this key, device will force to stop working.

8. **Start key:** press this key and device start work

7.3. CONSUMABLES

7.3.1. SAMPLE CONTAINER, APPLIED SOLUTION TUBE (PLASTIC), SAMPLING SPOON.

Picture 4. Consumables

- Special items for using only with Processor Sciendox1A.
- Sample container has two-stage integrated filter with holes 2 mm and 250 μ m - 300 μ m.
- Disposable products.
- Do not use if damaged or contaminated.
- Forbid to use other brands of consumables.
- Only buy from designated agent.

Sample container

It's indicated for feces sample collection and processing

Applied solution tube (plastic)

It is a test tube where the filtrate goes

Sampling spoon

It is intended for feces collection.

- If you press the button on the end of the spoon, the front end is automatically reset.
- Convenient for collection of a certain (fixed) amount of sample



7.3.2 PLASTIC TUBE RACK


Special item to use with Processor Sciendox 1A for placement and temporary storage of applied solution tubes for filtrate.

7.3.3 BIO-EKON SAMPLE CONTAINER RACK

Special item to use with Processor Sciendox 1A for placement and temporary storage of sample containers.

7.3.4 OTHER COMPONENTS AND REAGENTS


Bio-Ekon Diluent bottle 1500mL 1 piece 	It is a plastic container with cap. Saline comes from the container into the device through a special tube. Saline: Used to dilute and soak the sample.
Pump pipe 2 pieces 	It's necessary to connect the peristaltic pump
Apparatus power cord 1 piece 	It is an electricity cord for connecting to the power supply
Fuses 2 pieces 	Performs a protective function, opens the circuit in which it is included, interrupting the current if it exceeds a specified value for a certain time.
Hex Wrench 1 suit, including Wrench 1.5 mm Wrench 2 mm Wrench 2.5 mm Wrench 3 mm Wrench 4 mm Wrench 5 mm Wrench 6 mm Wrench 7 mm 	It's necessary for equipment installation

M5 Phillips screwdriver 1 piece 	It's necessary for equipment installation
M3 Phillips screwdriver 1 piece 	It's necessary for equipment installation

7.4. SAMPLE PROCESS FLOW

7.4.1. SAMPLE PROCESS FLOW DESCRIPTION


Sample process flow process includes the addition of a diluent (saline), soaking, mixing, filtration.

Picture5. Sample container in the position of dilution.

7.4.2. ADD FLUID


Place the sample container at the dilution position (see pic. 5.).The device punctures the green part of the container with a needle. Then the device adds saline (see pic. 6) in the amount installed on the control panel (see Fig. 2) to the container through a needle, the standard volume is 4 ml.


Picture 6. Add fluid

7.4.3. MIXING

When the device finish adding fluid put the container to mixing position (pic.7). After that, it is necessary to pull out the tongue for the filtrate tube, insert the tube into the special opening and close the tongue to the stop (pic. 8). Start the mixing and filtering process from the control panel (pic. 3).

Picture7. Sample container in the position of mixing and filtration.


Picture 8. The tube in the position of filtering

The device mixes the sample according to the set parameters. The sample is mixed at high speed for transformation into a homogeneous suspension. The standard number of mixing times is 4.

7.4.4. FILTERING

The finished suspension is fed into the applied solution tube for filtrate under pressure, through a piston mechanism (pic. 9), the standard filtration time is 4 seconds.

Picture 9. Sample filtering

7.4.5. Sample preparation completion.

After feeding the sample to the tube for the filtrate, the sample preparation is completed.


Chapter 8: Installation and connection

8.1 UNPACKING

Carefully remove the device and its components from the transport box. Check the integrity of the content. Remove the protective pads. If the device or its components are damaged, immediately contact the supplier.

8.2 INSTALLATION

Install the device as shown in picture 10 (1 m of free space on all sides and 10 cm from the wall)

Picture10. Installation scheme

Do not plug in during the installation.

8.3 ELECTRICAL CONNECTION

Use the original cable for electrical connection.

Input power supply: 230V, 50Hz. Use the fuse for 2A.

Input power must have ground connection!

Connect power after device installation and confirmed safety condition!

SALINE FEED COLLECTION

Settle diluent bottle for saline and attach the PE tube to the diluent bottle.

All liquids must be filtered by $<50\mu\text{m}$ filter screen to avoid clogging in the inner tubes.

8.4 TUBES CLEANING

Before using the new device, you must drain and clean the tubes (see Section 10.2 of this manual).

After cleaning the tubes and setting up the device, you can begin preparing the sample for analysis.

Chapter 9: Device operation

9.1 PREPARATION FOR WORK

Before work:

1. Put the empty sample container in the position of dilution;
2. Check the integrity of the sample container and the saline quantity. During sample preparation process regularly check the liquid amount in the container and add if necessary;
3. Put the switch on.

9.2 COLLECTION OF BIOMATERIAL AND DETERMINATION OF THE PROCESSING PARAMETERS

Collection of biomaterial:

Use only new sample container to collect the biomaterial.

Follow the instruction printed on the sample container label.

Biomaterial does not contain external particles to avoid the sample damage.






To achieve the correct concentration of the final sample, it is necessary to take the sample quantity as recommended in Table 2. The description of the types of feces in the Bristol stool chart is given in Table 3.

Table 2. Recommendations for biomaterial's collection.

Feces consistency	Type according to Bristol stool chart	Quantity of spoons
Solid stool	1,2	1
Normal stool	3,4,5	1
Liquid stool	6,7	2-4

If present, take part with blood, mucus, etc., and do not mix the sample with water, urine, or other liquids.

Table 3. Bristol stool chart

Type	Description	Picture
Type 1	Separate hard lumps like nuts (hard to pass)	
Type 2	Sausage-shaped, but lumpy	
Type 3	Like a sausage but with cracks on its surface	
Type 4	Like a sausage or snake, smooth and soft	
Type 5	Soft blobs with clear cut edges	
Type 6	Fluffy pieces with ragged edges, a mushy stool	
Type 7	Watery, no solid pieces	

Determination of sample preparation parameters:

To achieve the correct concentration of the final sample, it is necessary to determine the amount of saline added, the number of mixing cycles and the filtration time as recommended in Table 4.

Table 4 Recommendations for sample preparation

Stool consistency	Type	Saline volume, ml	The number of mixing cycles	Filtration time, c
Solid stool	1,2	4	6-8	4
Normal stool	3,4	4	4	4
Liquid stool	5,6,7	1-2	1-3	4

Note: Feces volume + saline volume ≤ 4 ml.

9.3 SAMPLE PREPARATION

Saline adding.

- The lid of the sample container should be closed during the entire period of the sample preparation.
 - The amount of saline can adjust between 0ml to 5 ml depending on the requirements for analysis.
- See the instructions for setting the standard amount of saline -4 ml.

Table 5. Step-by-step adjustment of the amount of saline per 4 ml.







№	Operate instruction	Operate key	Indicate condition
1	Connect power plug to socket.	—	Show current liquid volume parameter 00.00 mL
2	Press SET key	SET	The first number twinkle 00.00 mL
3	Adjust value to zero by use value increase and reduce key.	+ —	Press this key ,will increase Press this key ,will reduce 00.00 mL
4	Press SET key	SET	The second number twinkle 00.00 mL
5	Adjust value to 4	+ —	Press this key ,will increase Press this key will deduce 04.00 mL
6	Press SET key	SET	Four number twinkle mean storage setting and show current liquid 4ML 04.00 mL
	Press START key.	START	Dilution is started.

Mixing and filtration

After saline adding the diluted solution mixed and filtered according to setting parameters. See the instructions for setting the standard parameters of intensity (8) and filtration time (4 sec).

Table 6. Step-by-step adjustment of mixing intensity (8) and mixing time (for 4 sec).

№	Operate instruction	Operate key	Indicate condition
1	Connect power plug to socket.	—	Show current mixing strength parameter 00.00 s
2	Press SET key	SET	The first number twinkle 00.00 s

№	Operate instruction	Operate key	Indicate condition
3	Use value increase/reduce key adjust "mixing time" (T) to zero	+ —	 S(+)increase value,(-) reduce value
4	Press SET key	SET	 S the second value twinkle
5	use value increase/reduce key make "mixing time"(T) second value adjust to 8	+ —	 S(+) increase value (-)reduce value.
6	Press SET key	SET	Save add liquid volume and show current mixing intensity. The third value twinkle, setting filter time
7	Use value increase/reduce key adjust "filter time" (S) first value to zero	+ —	 S(+) increase value (-)reduce value.
8	Press SET key	SET	the fourth value twinkle
9	Use value increase /reduce key adjust "fliter time"(S) second value to 4	+ —	 S(+) increase value (-)reduce value
10	Press SET key	SET	 S Filter time value is set.
11	Press START key.	START	The mixing and filtration is started

Can choose mixing intensity from 0 to 99 times.


The filtration time can be chosen from 0 to 99 sec.

Setting one of the parameters to zero means close this function.

EXAMPLES:


Example 1: Mixing and filtration.

If set mixing to 8 times and filter time as 4 sec

Indicators will show:  S


Example 2: Only mixing.

If set mixing to 8 times and filter time as 0 sec

Indicators will show:  S

Example 3: Only filtration

If set mixing to 0 times and filter time as 4 sec

Indicators will show: 

- For convenience, number the applied solution tubes for the filtrate immediately after removing it from the device.
- After setting the parameters, keep the appliance door closed all the time during the sample preparation.

Chapter 10: Daily maintenance

10.1 SALINE REPLACEMENT

Daily saline replacement is recommended. Otherwise the bottle for saline can become unsuitable for use.

All liquids must be cleaned with a filter with a purity of <50µm. in order to avoid clogging of the inner tubes.

The saline bottle needs to be replaced if damaged or deteriorated.

10.2 TUBE MAINTAIN (CLEAN, DRAIN)

When replace saline or eliminate the bubbles in the tubes, drain the tubes. Make sure that there is enough liquid in the saline bottle, and an empty sample container with the open lid is set to the adding fluid position (see figure 5). Set the diluents volume as 5 ml then set start on the adding fluid control interface (see fig. 2). Repeat the operation until the tubes is filled with fluid.

10.3 PERISTALTIC PUMP TUBE REPLACEMENT

Check the pump tubes regularly, replace if damaged or deteriorated. The peristaltic pump replacement needs to be made every 3-6 months. Only qualified person can replace it. Please change liquid bottle once it is broken or aging.

10.4 FUSES REPLACEMENT

Unplug the device, replace the fuse. Technical characteristics: 230V, 50Hz, fuse for 2A

10.5 PUNCTURE NEEDLE

The needle is sharp –be careful when handling.

The needle replacement needs to make every 1-2 years only by qualified person.

10.6 MAINTENANCE PERIOD REFERENCE TABLE

Table 7. Maintenance period reference table

№	Item	Period	Instuction
1	Diluent check	Everyday	Check the liquid volume, add if necessary
2	Peristaltic pump tube check and change	Once a month	Check the tube, replace if damaged or deteriorated
3	Puncture needle change	Irregular	Check the needle, replace if damaged or deteriorated
4	Inner tubes change	Irregular	Check the inner tubes, replace if damaged or deteriorated

Correct operation and regular maintenance can prolong the service life of the instrument.

CHAPTER 11: TECHNICAL SERVICE

Only authorized person should service the system or its components. For maintenance and repair, contact the authorized representative of the manufacturer.

CHAPTER 12: CURRENT REPAIR

Only authorized personnel should service the system or its components. For maintenance and repair, contact the authorized representative of the manufacturer. Please do not open the device to observe and modify without permit. In case of the following you should unplug the power supply and ask professional person:

Power supply cable, signal cable, or socket broken;

Device in water;

Device to be dropped and broken, or framework broken;

Operate according to instruction book, but still cannot work;

Device character have big change, after exclude sample, reagent, and operate factor. But still cannot work. Then need maintain.

Only authorized personnel should service the system or its components. For maintenance and repair, contact the authorized representative of the manufacturer (see at the end of this document).

CHAPTER 13: STORAGE AND TRANSPORTATION. SERVICE LIFE.

Transportation and storage environment

Temperature	From 0 °C to 55 °C
Relative Humidity	No more than 85 %
Atmospheric Pressure	50,0 to 106,0 kPa,

The product is transported by all modes of transport in covered vehicles in accordance with the requirements applicable to each type of transport.

The expected service life for Feces sample liquefaction processor Sciendox 1A is 7 years.

Guaranteed usage period 12 months.

Guaranteed storage period 12 months.

Concretization of the professional level of potential users.

Medical laboratory technician (laboratory assistant) or other specialist.
Operators must be trained.

CHAPTER 14: WARRANTY OBLIGATIONS

All questions related to the operation and quality of the product can be addressed to the authorized representative of the manufacturer in the territory of the Russian Federation.

The manufacturer guarantees the conformity of the product to the established requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation.

The manufacturer is not liable for any damage or damage caused by using the product until proven guilty.

The manufacturer is not responsible for any damage or damage caused by using the product after

the expiry date printed on the package.

For questions about the quality of a medical product, manufactured by Xiamen Sciendox Biological Technology Co., Ltd, China, please contact:

REAL CARE LLC

111524, Moscow, Perovskaya street, house 1, building 10, floor 2, room VII room 3

Telephones: +7(495)234-57-31, +7(495)780-74-88, +7926-613-04-00

CHAPTER 15: DATA TO DESTRUCTION AND UTILIZATION OF MEDICAL DEVICE

In case of clinic, warehouse, etc. medical device is utilized/destroyed in accordance with medical waste circulation rules of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation»

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

CHAPTER 16: Requirements to the Electromagnetic Compatibility

The Sciendox 1A stool sample preparation apparatus must be connected and operated in 230 V and 50Hz AC power grids. The power consumption of the device is not more than 65 VA.

The device Sciendox 1A is designed for use in the electromagnetic environment described below. The customer or user of the device must ensure that the Sciendox 1A stool sample preparation Apparatus is used in such the same conditions described below

Radiation testing	Compliance	Electromagnetic environment-manual
Radiations CISPR 11 (industrial interference From industrial, scientific, medical and household high frequency devices)	Class A	The apparatus is intended for use in all non-residential premises, and may be used in residential and other premises directly connected to the public low voltage electrical network supplying buildings used for residential purposes, if the following precautions are observed:  Attention: This equipment (apparatus) is intended for use by medical professionals only. This equipment (apparatus) may cause radio frequency interference or disrupt adjacent equipment. You may need to take mitigation measures, such as reorientation in space, moving the machine, or shielding.
The radiation of the harmonic components IEC 61000-3-2	complies with requirements	
Voltage fluctuations/flickering radiation IEC 61000-3-3	complies with requirements	

The device Sciendox 1A is designed for use in the electromagnetic environment described below. The customer or user of the device must ensure that the Sciendox 1A stool sample preparation Apparatus is used in such the same conditions described below			
Tests for interference immunity	IEC 60601 test level	Level of compliance	Electromagnetic environment-rules and regulations
Electrostatic discharge, IEC 61000-4-2	8 kV contact Air discharge at 2, 4, 8, 15 kV	8 kV contact Air discharge at 2, 4, 8, 15 kV	The floor should be wooden, concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients /bundles, IEC 61000-4-4	2 kV for power supply lines 1 kV for input and output lines	2 kV for power supply lines «Inapplicable»	The quality of the power supply must meet standard commercial or hospital conditions.
IEC 61000-4-5 current surge	1 kV - differential mode • 2 kV in - phase voltage with respect to ground • 2 kV signal (input/output) relative to ground.	1 kV - differential mode • 2 kV in - phase voltage with respect to ground • «Inapplicable»	The quality of the power supply must meet standard commercial or hospital conditions.

The voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	0% UT; 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single-phase at 0° 0% UT; 250/300 cycles	0% UT; 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single-phase at 0° 0% UT; 250/300 cycles	The quality of the power supply must meet standard commercial or hospital conditions. If the user of the device needs to work with it continuously during a power interruption, it is recommended to use an uninterruptible power supply or a battery.
Industrial frequency magnetic field (50/60 Hz) iec 61000-4-8	30 A/m	30 A/m	Magnetic fields of industrial frequency should be typical for standard commercial enterprises or medical institutions.

The device Sciendox 1A is designed for use in the electromagnetic environment described below. The customer or user of the device must ensure that the Sciendox 1A stool sample preparation Apparatus is used in such the same conditions described below			
Noise immunity test	Test level IEC 60601	Level of compliance	Electromagnetic environment-manual
Conductive RF interference, IEC 61000-4-6	3 V, 6 V	3 rms, 6 watts (Rms is the RMS value of the voltage) 3 V/m	Portable and mobile radio frequency communication equipment shall be used at a distance not less than the recommended spatial separation value from any part of the device and shall be calculated according to the formula for the corresponding frequency of the transmitter. The recommended spatial
Radiated RF interference, 61000-4-3	3 V/m		

	3 V / m at 150 kHz-80 MHz; 6 V / m at 150 kHz-80 MHz and 80% AM in 1 kHz 10 V / m at 80 MHz-2.7 GHz	3 V / m at 150 kHz-80 MHz; 6 V / m at 150 kHz-80 MHz and 80% AM in 1 kHz 10 V / m at 80 MHz-2.7 GHz	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{12}{V_2} \right] \sqrt{P}$ $d = \left[\frac{12}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{23}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ spacing: here P is the maximum permissible output power of the transmitter in watts (W) according to the manufacturer, d is the recommended spatial spacing in meters (m). The field strength of stationary RF transmitters, determined by electromagnetic examination of the object, must be below the level of compliance in each frequency range. Interference may occur in the vicinity of equipment  marked with the following symbol:
--	---	---	--

The Sciendox 1A stool sample preparation apparatus is designed for use in an electromagnetic environment with controlled radiated radio interference. The customer or user of the machine can help prevent electromagnetic interference by observing the minimum distance between portable and mobile RF communication devices (transmitters) and the machine, as recommended below according to the maximum output power of the communication equipment.

Maximum rated output power of transmitter Watt	Space interval depending on the frequency of the transmitter (m)			
	150 w -80 MHz beyond the industrial, scientific and medical range $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	150 w -80 MHz beyond the industrial, scientific and medical range $d = \left[\frac{12}{V_2} \right] \sqrt{P}$	80 MHz - 800 MHz $d = \left[\frac{12}{E_1} \right] \sqrt{P}$	800 MHz - 2,5 GHz $d = \left[\frac{23}{E_1} \right] \sqrt{P}$
0,01	0,12	0,2	0,4	1
0,1	0,37	0,64	1,3	2,6
1	1,17	2	4	8
10	3,7	6,4	13	26
100	11,7	20	40	80

For transmitters with a maximum output power not listed above, the recommended spatial spacing d in meters (m) can be estimated using the formula applicable to the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer's instructions.

NOTE 1: at 80 - 800 MHz, spatial spacing is applied for a higher frequency range.

NOTE 2: These principles may not apply in all situations. The propagation of electromagnetic waves is affected by their absorption and reflection from structures, objects and people.

Note

* Information on potential electromagnetic or other interference complies with IEC 60601-1-2: fourth edition (2014) Medical electrical equipment, part 1-2: Auxiliary standard: electromagnetic compatibility. Requirements and tests.

* The radiation characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11, class A). If used in a residential environment (which normally requires CISPR 11, class B), this equipment may not provide adequate protection for radio frequency services. You may be required to take action to mitigate impacts, such as relocation or repositioning of equipment.

Attention!

* Do not use this equipment in conjunction with other equipment or in combination with other equipment, as this may result in improper operation. If such use is necessary, this equipment and other equipment should be checked to ensure that they operate normally.

• Use of components, sensors, and cables other than those specified or provided by the manufacturer of this equipment may result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment, as well as improper operation.

* Portable RF communication equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the machine, including the cables specified by the manufacturer.

CHAPTER 17: Symbols

MARKING	DEFINITION
	Must ground connect mark: Must ground connect though socket
	Serial number
	Manufacturer
	Authorized representative in EC
	Refer to the instructions for use.
	Warning! Refer to the instructions for use.
	MD for in vitro diagnostics
	Caution! Need special disposal!
	Biohazard Sign

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