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English text

Instruction for Use
Single-use wand FLOW 50

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The components are designed to be operated as a single unit

DEVICE DESCRIPTION

The FLOW 50 Wand (Wand) is a bipolar, RF electrosurgical device used for specific indications in orthopedic and arthroscopic surgeries. The Wand utilizes the AMBIENT feature which provides accurate (± 3 °C) real-time temperature monitoring of the intraarticular irrigating fluid. The Wand is capable of operating in five Modes (Table 1).

Table 1: Modes and Levels for the FLOW 50 Wands

Action	Mode	Level	Mode Level Indicator
Ablation	Lo	Lo Minus (-)	
		Lo (+)	
		Lo Plus (+)	
Ablation	Med	Med Minus (-)	
		Med (+) *	
		Med Plus (+)	
Ablation	Hi	Hi Minus (-)	
		Hi (+)	
		Hi Plus (+)	
Coagulation	Coag	Coag (+) *	
		Coag Plus (+)	
Vacuum	Vac	Not Adjustable	N/A

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Table 2: Output Voltage and Power for the WEREWOLF COBLATION System with the FLOW 50 Wands

Mode/Level	Output Voltage (Vrms)	350 Ohm Load Output Power (Watts)
Lo -	257	189
Lo •	279	222
Lo +	289	239
Med -	262	196
Med •	279	222
Med +	301	259
Hi -	301	259
Hi •	314	281
Hi +	323	298
Mode/Level	Output Voltage (Vrms)	50 Ohm Load Output Power (Watts)
Coag •	65	85
Coag +	85	145

INDICATIONS FOR USE

The WEREWOLF COBLATION System, comprised of the FLOW 50 Wand and the WEREWOLF Controller is indicated for the resection, ablation, and coagulation of soft tissues and hemostasis of blood vessels during arthroscopic and orthopedic procedures. ablation, and coagulation of soft tissue and hemostasis of blood vessels in the following arthroscopic and orthopedic procedures:

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Joint	Ablation/Debridement
All Joints (Hip, Knee, Shoulder, Wrist, Ankle, Elbow)	- Articular Cartilage - Bursectomy - Chondroplasty - Fascia - Ligament - Scar Tissue - Soft Tissue - Synovectomy - Tendon
Knee	- ACL/PCL - Notchplasty
Shoulder	- Acromioplasty - Subacromial Decompression
Joint	Excision/Resection
All Joints (Hip, Knee, Shoulder, Wrist, Ankle, Elbow)	- Articular Labrum - Capsule - Cysts - Ligament - Loose Bodies - Plica Removal - Scar Tissue - Soft Tissue - Synovial Membrane - Tendon
Hip	Acetabular Labrum
Knee	- Capsular Release - Cartilage Flaps - Discoid Meniscus - Lateral Release - Meniscal Cystectomy - Meniscectomy - Villusectomy
Shoulder	- Frozen Shoulder Release - Glenoid Labrum
Wrist	Triangular Fibrocartilage (TFCC)

CONTRAINDICATIONS

The Wand is contraindicated in:

- Any non-arthroscopic or non-orthopedic surgical procedure.
- Any arthroscopic or orthopedic procedures where conductive irrigating fluid (e.g., saline or Ringer's Lactate) is not used.
- Patients with cardiac pacemakers or other electronic device implants.
- Patients for whom arthroscopic or orthopedic procedures are contraindicated for any reason.

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- The Wand should never be used as the primary outflow pathway for irrigation. Instead, it should be

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used as a secondary outflow pathway (in conjunction with a primary pathway) to prevent thermal injury to the user and patient.

• Ensure adequate circulation of irrigating fluid to prevent heating, loss of joint distension, elevated fluid temperatures that may cause thermal injury to the user or patient, and/or damage to the electrode tip. Do not use pre-warmed irrigating fluid with the Wand as it may result in tissue damage or thermal injury.

• Ensure that the Wand is completely surrounded by irrigating fluid when activated to prevent tissue damage or thermal injury.

• The Wand fluid outflow varies within and across Levels and Modes. Periodically check the Wand tubing for obstructions in flow. If the Wand's outflow is obstructed, do not attempt to clear the blockage. Instead, a new Wand must be used.

• The Wand suction line should not directly contact the patient as it may cause a burn. Drapes between the line and skin are recommended.

• Damaged shaft insulation of the Wand could lead to unwanted electrical conduction resulting in patient burns. Therefore, the Wand must be periodically checked to ensure device integrity. If the shaft insulation is damaged, use a new Wand.

• When not using the Wand, place it in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact of the Wand tip may cause thermal injuries to the user or patient.

• In Lo Mode, the system may have minimal or negligible outflow, which may result in heating of the Wand's shaft. Therefore, the Wand should be activated in short bursts avoiding continuous activation which may cause thermal injury to the user and patient.

General

• The System should be inspected prior to, and periodically during, surgery for any damage as this could adversely impact performance, safety and clinical outcome. Do not use if damaged or malfunctioning. Back-up devices should always be available to avoid any delay prior to or during surgery.

• If the Wand's finger switches become unresponsive, reposition fingers on the switch and depress again. If unsuccessful, use the foot control or open another Wand.

• The Wand is designed for ablation and/or coagulation only, and not for mechanical displacement of tissue through applied force. Do not use the Wand as a lever to enlarge a surgical site or gain access to tissue as this may result in a bent or detached electrode, damage to device, and/or a cracked spacer leading to patient injury.

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Damage includes bending or detachment of the distal end of the device.

Do not contact metal objects while activating the Wand as damage to the tip and/or shaft insulation may occur resulting in device malfunction or patient injury.

- The rate and depth of tissue ablation is affected by the set point selected, the amount of pressure on the tissue, the integrity of the electrode, and the speed at which the Wand is passed over the target tissue.

- Check the Wand tip periodically to ensure the electrode is intact and the Wand is functioning as intended. If the electrode is not intact, discontinue use and check the joint cavity by appropriate means. Excessive electrode wear includes, but is not limited to erosion, thinning, or detachment of the metal at the distal end of the device.

- Electrode failure and reduced life expectancy may occur if the recommended ablation time for each Mode is exceeded.

- When inserting and removing the FLOW 50 Wand, use caution and avoid contacting the distal portion of cannula with the shaft or electrode face as this may lead to potential electrode/shaft damage and/or patient injury.

- Ensure that the Wand suction tubing is properly loaded into the Integrated Fluid Outflow Regulator, primed, and irrigating solution is flowing through the suction tubing prior to treating tissue.

Failure to prime the suction tubing may result in tissue damage and/or decreased Wand performance.

- Do not touch the electrode at the Wand tip to exposed skin while power is being applied as injury may result.

- Do not contact the non-targeted tissue with the activated Wand tip as injury may result.

- Do not activate the controls until the Wand is in close proximity with the targeted tissue as injury may result.

- Do not withdraw the Wand while power is being applied as injury may result.

- Use only conductive irrigating fluid (e.g. saline, Ringer's Lactate).

Do not use non-conductive irrigating fluid (e.g. sterile water, air, gas, glycine) as this will compromise the performance of the Wand.

- Do not clean, resterilize, or reuse the Wand as this may damage or compromise the performance resulting in malfunction, failure, and/or patient injury and may expose the patient to the risk of infectious diseases.

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• Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire. Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as methanol or alcohol), as electrosurgical device has the potential for providing a source of ignition.

• Do not use flammable agents for cleaning and/or disinfection of the Controller as an electrosurgical device has the potential for providing a source of ignition.

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Precautions

• Lo Mode is recommended for use in debridement of articular cartilage.

• Vac Mode should not be used on tissue. Vac Mode provides high suction and RF energy which is intended to be used only to remove fine, less dense, free floating debris in order to improve surgical site visibility. Do not attempt removal of more coarse, dense debris as obstruction of the Wand outflow may occur.

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• Evaluate patients for predisposing medical problems that may be aggravated by the stress of surgery.

• Excessive wear of or damage to the electrode from vigorous use against very dense tissue or bony surfaces may lead to compromised performance resulting in patient injury.

• The Wand is not designed for prolonged ablation. Electrodes wear under use and the rate of wear is dependent on variables that cannot be replicated in all cases. Factors of electrode wear include, but are not limited to, rate of ablation, operating Mode and Level, and prolonged use against dense tissue or bony surfaces.

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• Do not insert other cables between the Wand cable and the Controller.

• As with other electrosurgical units, the electrode and cables can provide paths for high frequency current. Position the cable to avoid contact with patient or other leads. Other electrical equipment may experience interference when positioned near

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the system.

• High frequency (HF) electrosurgical equipment such as the WEREWOLF controller may adversely affect the operation of other electronic equipment.

• Monitoring electrodes (e.g., electrocardiography (ECG)) should be positioned as far as possible from the device tip when HF surgical equipment and physiological monitoring equipment are used simultaneously on the same patient. Monitoring needle electrodes (e.g., electromyography (EMG)) are not recommended.

• Monitoring equipment incorporating high frequency current limiting devices are recommended.

POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. Potential complications accompanying surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis); nerve injury or palsy; treatment or implant failure; and secondary surgical intervention to address complications associated with surgery or treatment.

DIRECTIONS FOR USE

Note: The FLOW 50 Wand is to be used only with the WEREWOLF controller.

Note: Prior to initial use, the User Manual for the WEREWOLF controller and the Instructions for Use for the Wand should be carefully reviewed. Particular attention should be directed to the indications, warnings, and precautions.

OPERATOR TRAINING REQUIREMENTS

• Use of the Wand requires an appropriate level of surgical skill and experience. Review the Instructions for Use and the User Manual, which is provided with the Controller, prior to clinical use.

• Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/ assistance with similar surgical procedures is recommended.

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• Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using the Wand.

• The surgeon should advise the patient of known risks and complications associated with surgical procedure and use of this product.

SYSTEM USE

Before Surgery

1. Prepare the patient pre-operatively according to standard procedures.

2. Switch on the Controller. Follow the instructions in the User Manual.

3. Connect a new sterile Wand to the Controller. The Controller will prompt you with an animation.

4. Open Controller door. The Controller will prompt you with an animation.

5. Insert Wand tubing, matching the red portion of tubing to the red portion of the Controller and the grey portion of tubing to the grey portion of the Controller. The Controller will prompt you with an animation.

NOTE: If the Wand tubing is not properly installed, the Controller will not allow activation of RF energy.

6. Close the door. The Controller will prompt you with an animation.

7. Attach the tubing connector to a waste collection receptacle.

8. If desired, connect the foot control to the Controller. Follow the instructions provided in the User Manual.

9. Once the Wand is connected, the Controller automatically displays the default Mode and Level. Priming Step: Enter the joint space ensuring the Wand is not activated. Confirm that the tip is fully submerged in conductive irrigating fluid and not contacting tissue. Press and immediately release the Ablation or Coagulation foot pedal or Wand finger switch (any audible tones and visible icons on the controller will stop once the pedal/switch is released). This action will start the Integrated Fluid Outflow Regulator, which may continue to run even after the pedal/button is released.

NOTE: There will be no audible tones and the visible ablate icons on the Controller will not be illuminated when the Wand is deactivated during the priming step.

CAUTION: Continuing to hold the pedal/button during priming will continuously activate the ablate/coagulate functions.

Contact with tissues should be avoided to minimize the risk of inadvertent damage.

10. Confirm irrigating solution is flowing through the suction tubing between the Wand and the Integrated Fluid Outflow Regulator.

11. Once the Wand is primed, and fluid is in the suction tubing, the Wand is ready to be used.

CAUTION: The Wand should be activated away from the tissue and moved into close proximity for ablation or coagulation.

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During Surgery

1. Insert the Wand into the treatment area through a prepared entry point. Ensure the Wand tip is surrounded by irrigating fluid when in use

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CAUTION: When inserting and removing the Wand, use caution and avoid contacting the distal portion of cannula with the shaft or electrode face as this may lead to potential electrode/shaft damage and/or patient injury.

2. To activate an Ablation Mode using the integrated finger switch, press and hold the YELLOW button. The Ablation Mode and Level

will be displayed on the Controller and a tone will be heard. The Wand should be activated prior to contacting the tissue and then moved into light contact with the targeted tissue. While activated, the Wand should be kept in constant motion. Refer to the User Manual for the Controller on how to activate using the foot control.

NOTE: Lo Mode is recommended for the debridement of articular cartilage. For treating articular cartilage, the Wand should be positioned on its side (45 to 90 degrees) so that only the edge of the electrode contacts the targeted tissue. This will help to control ablation depth and enhance visibility.

NOTE: The Wand is designed to be used in each Mode and electrode wear rates will vary. For example, on meniscus tissue, the FLOW 50 Wand has been shown to maintain its integrity for the duration of the following ablation times:

Lo Mode: 10 minutes

Med Mode: 4 minutes

Hi Mode: 2 minutes

Vac Mode: 4 minutes, 30 seconds

CAUTION: Check the Wand tip periodically to ensure the electrode is intact and the Wand is functioning as intended. If the electrode is not intact, discontinue use and check the joint cavity by appropriate means.

CAUTION: Electrode failure and reduced life expectancy may occur if the recommended ablation time for each Mode is exceeded.

3. The Ablation Mode may be changed by momentarily pressing the BLACK Mode button on the Wand handle or foot control, or by selecting the Mode on the Controller touch screen.

4. **Vac Mode** provides high suction and RF energy which is intended to be used only to remove finer, less dense, free floating debris in order to improve surgical site visibility. The **Vac Mode** is activated by holding down the BLACK Mode button on the Wand handle or foot control. Once the BLACK Mode button

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is released, the system will revert to the most recent Ablation Mode and Level.

NOTE: The BLACK Mode button on the Wand handle and foot control can change Ablation Modes and enable the **Vac Mode**.

Levels within each Mode can only be adjusted using the touch screen on the Controller.

CAUTION: The Wand should not be in contact with tissue when **Vac Mode** is activated.

CAUTION: Do not attempt removal of more coagulated debris as obstruction of the Wand outlet may occur.

5. The Wand incorporates the Ambient feature that monitors the real-time temperature of the intra-articular irrigating fluid (between 20–60 °C, ± 3 °C). The Controller automatically detects this function and enables the temperature display on the touch screen. During use, if the measured temperature exceeds this set point, an audible and visual alarm is activated to alert the operator. In addition, the operator can select a temperature alarm set point and alarm volume or this feature can be disabled by touching the thermometer on the Controller touch screen.

NOTE: The user may elect to continue to use the product after the alarm has sounded. Use and performance of the Wand is not affected when the alarm has sounded.

NOTE: This system also contains a feature to prevent elevated Wand tubing temperatures. If elevated temperatures are detected, the RF energy will cycle on and off until lower temperatures are achieved.

To activate the Coag Mode using the integrated finger switch, press and hold the BLUE button. The Coag Mode and Level will be displayed on the Controller and a tone will be heard. The Controller automatically displays the default Coag Level. Refer to the User Manual for the Controller on how to activate using the foot control.

6. Once RF activation has been off for more than 3 minutes the suction will stop. To restart suction, activate RF.

After Surgery

1. Remove the Wand from the joint space.
2. Disconnect the Wand from Controller by unscrewing the Connector and pulling on the connector body.
3. Open the Controller door and remove Wand tubing.

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4. Dispose of the Wand (see Disposal Section below).

SAFETY FEATURES

The safety features of the WEREWOLF controller are outlined in the User Manual which is provided with the Controller.

HOW SUPPLIED

• The Wand is supplied sterile for SINGLE USE ONLY. The Wand is sterile only if the packaging is not opened, damaged or broken.

If the packaging is damaged, do not use the Wand.

• Do not use the Wand after the "use-by" date printed on the label. Wand performance and patient safety may be compromised if used after the expiration date.

• The Wand is sterilized using ethylene oxide.

• The Wand is not made with natural rubber latex.

• Wand materials do not contain DEHP.

• Rated voltage for the FLOW 50 Wand is 340 Vrms.

DISPOSAL

The device must be disposed of in accordance with applicable government regulations, and the Medical Waste Management Plan for the facility.

MAINTENANCE AND TROUBLESHOOTING

Refer to the User Manual for the WEREWOLF controller for operating and troubleshooting guidelines for the Controller.

CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/ OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

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Product Complaints and Return Authorization

All questions, complaints, and return authorization requests should be directed to Customer Service or an authorized representative.

ArthroCare Corporation


7000 West William Cannon
Austin, TX 78735 U.S.A.
(800) 343-5717
CUSTOMER SERVICE
www.arthrocare.com


Smith & Nephew Orthopaedics GmbH

Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

SYMBOLS KEY	
	Catalogue number
	Lot number
	Use-by Date
	Contents
	Consult Instructions for Use
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Keep dry
	Keep Away From Sunlight
	Caution

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STERILE EO	Sterilized Using ethylene oxide
	Manufacturer
EC REP	Authorized representative in the European Community
Rx only	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
	CE mark and identification number of Notified Body. This product meets the essential requirements of Medical Device Directive (93/42/EEC).
Icons Used on the Wand/ Controller	
 ABLATE	Ablation
Lo ~	Lo Mode Ablation
Med ≈	Med Mode Ablation
Hi ≡	Hi Mode Ablation
 COAG	Coagulation

Med 	Режим абляции Средний (Med)	Vac 	Vac Mode
Hi 	Режим абляции Высокий (Hi)	™Trademark of Smith & Nephew. This product may be covered by one or more US Patents. See smithnephew.com/patents for details. ©Smith & Nephew 2019	
 COAG	Коагуляция		
Vac 	Режим Отсос (Vac)		

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SYSTEM DESCRIPTION

The WEREWOLF COBLATION System is an electrosurgical system consisting of the:

- Bipolar radiofrequency (RF) Controller;
- Disposable, sterile Wand;
- Reusable, non-sterile foot control; and
- Reusable, non-sterile power cord.

The components are designed to be operated as a single unit.

DEVICE DESCRIPTION

The FLOW 90 Wand (Wand) is a bipolar, RF electrosurgical device used for specific indications in orthopedic and arthroscopic surgeries.

The Wand features AMBIENT™ technology which provides accurate ($\pm 3^\circ\text{C}$ between $20-60^\circ\text{C}$) real-time temperature monitoring of the intra-articular irrigating fluid. The Wand is designed for use up to 7 minutes of cumulative active ablation and is to be used only with the WEREWOLF™ Controller (Controller). The Wand is capable of operating in five Modes (Table 1).

Table 1: Modes and Levels for the FLOW 90 Wands

Action	Mode	Level	Mode Level Indicator
Ablation	Lo	Lo Minus (-)	
		Lo (*)	
		Lo Plus (+)	
Ablation	Med	Med Minus (-)	
		Med (*)	
		Med Plus (+)	
Ablation	Hi	Hi Minus (-)	
		Hi (*) *	
		Hi Plus (+)	
Coagulation	COAG	Coag (*) *	
		Coag Plus (+)	
Vacuum	Vac	Not Adjustable	N/A

*Default Modes and Levels for FLOW 90

Note: Lo Mode is recommended for the debridement of articular cartilage.

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Table 2: Output Voltage for the WEREWOLF COBLATION System with the FLOW 90 Wands

Mode/Level	Output Voltage (Vrms)
Lo -	271
Lo •	277
Lo +	284
Med -	287
Med •	295
Med +	300
Hi -	305
Hi •	325
Hi +	328
Mode/Level	Output Voltage (Vrms)
Coag •	85
Coag +	105

INDICATIONS FOR USE

The FLOW 90 Wand, used with the WEREWOLF COBLATION System, is indicated for the resection, ablation, and coagulation of soft tissues and hemostasis of blood vessels in the following arthroscopic and orthopedic procedures:

Joint	Ablation/Debridement	Excision/Resection
All Joints (Hip, Knee, Shoulder, Wrist, Ankle, Elbow)	- Articular Cartilage - Bursectomy - Chondroplasty - Fascia - Ligament - Scar Tissue - Soft Tissue - Synovectomy - Tendon	- Articular Labrum - Capsule - Cysts - Ligament - Loose Bodies - Plica Removal - Scar Tissue - Soft Tissue - Synovial Membrane - Tendon
Hip		Acetabular Labrum
Knee	- ACL/PCL - Notchplasty	- Capsular Release - Cartilage Flaps - Discoid Meniscus - Lateral Release - Meniscal Cystectomy - Meniscectomy - Villousctomy
Shoulder	- Acromioplasty - Subacromial Decompression	- Frozen Shoulder Release - Glenoid Labrum
Wrist		Triangular Fibrocartilage

CONTRAINDICATIONS

The Wand is contraindicated in:

- Any non-arthroscopic or non-orthopedic surgical procedure.

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- Any arthroscopic or orthopedic procedures where conductive irrigating fluid (e.g., saline or Ringer's Lactate) is not used.
- Patients with cardiac pacemakers or other implanted electronic devices.
- Patients for whom arthroscopic or orthopedic procedures are contraindicated for any reason.

WARNINGS and PRECAUTIONS

⚠ Warnings

Potential User/Patient Thermal Injury

- Do not use pre-warmed irrigating fluid with this Wand.
- The Wand suction is not intended to replace primary outflow for irrigating fluid. Wand suction should be used in conjunction with another source of outflow.
- Ensure adequate flow and circulation of irrigating fluid to prevent excessive heating of the fluid. Periodically check the Wand tubing for obstructions or restrictions in flow. If the Wand outflow is obstructed, do not attempt to clear the blockage.
- Use a new Wand.
- Damaged shaft insulation of the Wand could lead to unwanted electrical conduction. Periodically check the Wand to ensure device integrity. If the shaft insulation is damaged, use a new Wand.
- Do not contact the non-targeted tissue with the activated Wand tip. Do not activate the controls until the Wand is in close proximity with the targeted tissue. Do not withdraw an activated Wand.
- When not using the Wand, place it in a clean, dry, highly visible area not in contact with the user or patient. Inadvertent contact of an active or inactive Wand tip may cause user or patient thermal injury.
- In Lo Mode, the system may have minimal to negligible outflow, which may result in heating of the Wand shaft. Therefore, avoid continuous activation of the Wand which may cause thermal injury to the user and patient.

General

- The System should be inspected prior to, and periodically during surgery for any damage, as this could adversely impact performance, safety, and clinical outcome. Do not use if damaged or malfunctioning. Back-up devices should always be available to avoid any delay prior to or during surgery.

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• Use caution when inserting and removing the wand with a cannula. Contacting the Wand electrode and/or insulation against edges of the cannula may lead to Wand damage.

• Do not contact metal objects while activating the Wand as damage to the Wand insulation and/or tip may occur resulting in device malfunction or patient injury.

• Do not bend, modify or alter the Wand in any way prior to or during surgery as this may damage the device and lead to failure and/or user or patient injury.

• Use only conductive irrigating fluid (e.g. saline, Ringer's Lactate). Do not use non-conductive irrigating fluid (e.g. sterile water, air, gas, glycine) as this will compromise the performance of the Wand.

• The Wand is not designed for mechanical displacement of tissue through applied force. The Wand is indicated for the ablation and coagulation. Do not use the Wand as a lever to enlarge a surgical site or gain access to tissue as this may result in a bent or detached electrode, damage to the device, and/or a cracked spacer leading to patient injury.

• Do not allow patient contact with grounded metal objects as injury may result.

• Electrode failure and reduced device performance may occur if the recommended ablation time is exceeded. Check the Wand tip periodically to ensure electrode integrity. If a portion of the electrode detaches from the Wand, discontinue use and check the joint cavity by appropriate means.

• Do not clean, resterilize, or reuse the Wand as this may damage or compromise the performance resulting in malfunction, failure, and/or patient injury and may expose the patient to the risk of infectious diseases.

Flammability/Electrical Shock

• Observe fire precautions at all times. Sparking and heating associated with electrosurgery may be an ignition source.

• Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or heated have the potential for providing a source of ignition.

• Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen), or in close proximity to volatile solvents (such as methanol or alcohol) as the electrosurgical device has the potential for providing a source of ignition.

• Do not use flammable agents for cleaning of the Controller as the electrosurgical device has the potential for providing a source of ignition.

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• Do not allow fluid to contact the end of the cable connector which mates to the Controller as electrical shock may occur.

Precautions

- Wand is only to be used with compatible arthroscopic instruments and equipment.
- If the Wand finger switches become unresponsive, reposition fingers on the switch and depress again. If unsuccessful, use the foot control or use a new Wand.
- Neuromuscular stimulation is an expected consequence of electrosurgical procedures. Ensure that adequate precautions are taken.
- Evaluate patients for predisposing medical problems that may be aggravated by the stress of surgery.
- The rate and depth of tissue ablation is affected by the setting, the amount of pressure applied to the tissue, the size of the electrode, and the speed at which the Wand is passed over the target tissue.
- Use the lowest Controller setting necessary to achieve the desired tissue effect.
- Lo Mode is recommended for the debridement of articular cartilage.
- The Wand is not designed for prolonged ablation. The Wand electrode wears with use. The rate of wear is dependent on a combination of variables that cannot be replicated for every situation. Factors of electrode wear include, but are not limited to, rate of ablation, operating Mode and Level, and prolonged use against dense tissue or bony surfaces.
- Do not insert other cables between the Wand cable and the Controller.
- As with other electrosurgical units, the electrode and cables can provide paths for high frequency (HF) current. Position the cable to avoid contact with patient or other leads. Other electrical equipment may experience interference when positioned near the System.
- HF electrosurgical equipment such as the WEREWOLF Controller may adversely affect the operation of other electronic equipment.
- Monitoring electrodes (e.g., electrocardiography (ECG)) should be positioned as far as possible from the device tip when HF surgical equipment and physiological monitoring equipment are used simultaneously on the same patient. Monitoring needle electrodes (e.g., electromyography (EMG)) are not recommended.
- Monitoring equipment incorporating HF current-limiting devices are recommended.

POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. Potential complications accompanying

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surgery with a device may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis); nerve injury or palsy; treatment or implant failure; and secondary surgical intervention to address complications associated with surgery or treatment.

DIRECTIONS FOR USE

Note: The FLOW 90 Wand is to be used only with the WEREWOLF Controller.

Note: The User Manual for the WEREWOLF Controller and the Instructions for Use for the Wand should be carefully reviewed prior to initial use. Particular attention should be directed to the indications, warnings, and precautions.

USER RECOMMENDATIONS

- Use of the Wand requires an appropriate level of surgical skill and experience. Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using the Wand.

System Preparation

1. Prepare the patient pre-operatively according to standard procedures.
2. Install the WEREWOLF Controller according to the instructions in the User Manual. Press the Main Power switch (I/O) on the front of the Controller to activate the Controller.
3. Connect a new sterile Wand connector to the Controller. The Controller will prompt you with animation.
4. Open the Controller door. The Controller will prompt you with animation.
5. Insert the Wand tubing. Match the connector colors on the tubing with the colors of the fitting on the Controller. The Controller will prompt you with animation.

Note: If the Wand tubing is not properly installed the Controller will not allow activation of R energy.

6. Close the door. The Controller will prompt you with animation.
7. Attach the tubing connector at the end of the tubing to a waste collection receptacle.
8. If desired, connect the foot control to the Controller. Follow the instructions provided in the User Manual.

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9. When the Wand is installed in the Controller, the Controller sets the Wand Coag Mode and Ablation Mode and Level to default.

Select the appropriate settings on the Controller touch screen following the instructions provided in the User Manual.

Wand Use

Note: The electrode at the distal tip of the Wand is applied to tissue to effect ablation and coagulation. The Wand is bipolar. A neutral electrode is not used with this Wand.

CAUTION: Use the lowest Controller setting necessary to achieve the desired tissue effect.

CAUTION: The Wand is not designed for mechanical displacement of tissue through applied force. Do not use the Wand as a lever to enlarge a surgical site or gain access to tissue as this may result in a bent or detached electrode, damage to the device, and/or a cracked spacer leading to patient injury.

1. Insert the Wand into the treatment area through a prepared entry point. Ensure the Wand tip is surrounded by irrigating fluid when in use.

CAUTION: Use caution when inserting and removing the Wand with a cannula. Contacting the Wand electrode and/or insulation against edges of the cannula may lead to Wand damage.

Ablation

2. To activate an Ablation Mode using the integrated finger switch, press and hold the YELLOW button. The Ablation Mode and Level will be displayed on the Controller and a tone will be heard.

The Wand should be activated prior to contacting the tissue then moved into light contact with the targeted tissue.

While activated, the Wand should be kept in constant motion.

Refer to the User Manual for the Controller for how to activate using the foot control.

Note: Lo Mode is recommended for the debridement of articular cartilage. For treating articular cartilage, the Wand should be positioned on its side (45 to 90 degrees) so that only the edge of the electrode contacts the targeted tissue. This will help to control ablation depth and enhance visibility.

Note: During periods of continuous activation, temperatures up to 48 °C on the Wand shaft and up to 60 °C on the Wand suction tubing may be reached.

3. To change the Ablation Mode, momentarily press the BLACK button on the Wand handle or foot control, or select the Mode on the Controller touch screen. Levels within each Ablation Mode

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can be adjusted only by using the Controller touch screen.

CAUTION: Electrode failure and reduced device performance may occur if the recommended ablation time is exceeded.

Check the Wand tip periodically to ensure electrode integrity.

If a portion of the electrode detaches from the Wand, discontinue use and check the joint cavity by appropriate means.

Note: Reduced device performance may be caused by obstructed or restricted flow. Periodically check the Wand tubing for obstructions or restrictions in flow. If the Wand outflow is obstructed, do not attempt to clear the blockage and use a new Wand.

Coagulation

4. To activate the Coag Mode using the integrated finger switch, press and hold the BLUE button. The Coag Mode and Level will be displayed on the Controller and a tone will be heard.

Refer to the User Manual for the Controller on how to activate using the foot control.

Vacuum

Note: Vac Mode provides RF energy and high suction intended to only be used to remove fine, free floating debris to improve surgical site visibility. Attempting to remove large, dense debris may obstruct the Wand outflow.

5. To activate the Vac Mode, press and hold the BLACK button on the Wand handle or foot control. When the BLACK button is released, the System will return to the most recent Ablation setting.

Note: Suction through the Wand will continue for 3 minutes after the last activation of any Mode. After 3 minutes, suction through the Wand will stop. To restart suction, activate any Mode.

AMBIENT

6. AMBIENT technology provides accurate ($\pm 3^\circ\text{C}$ between $20-60^\circ\text{C}$) real-time temperature monitoring of the intra-articular irrigating fluid. The user can select a temperature alarm set point and alarm volume, or this feature can be disabled by touching the thermometer on the Controller touch screen. If the measured temperature exceeds this set point, an audible and visual alarm is activated to alert the user. Performance of the Wand is not affected when the alarm is activated. **Note:** The Wand contains a feature to prevent elevated Wand tubing temperature. If an elevated temperature is detected, the energy to the Wand will cycle on and off until the Wand tubing temperature is reduced.

System Shut Down

1. Remove the Wand from the treatment area.

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2. Press the Main Power switch (I/O) on the front of the Controller to deactivate the Controller. The lights on the Controller will go off within five seconds.

3. Disconnect the Wand connector from the Controller. Do not separate the cable from the Wand handle or connector.

4. Open the Controller door and remove the Wand tubing.

5. Dispose of the Wand (See Disposal Section below).

SAFETY FEATURES

The safety features of the WEREWOLF Controller are outlined in the User Manual provided with the Controller

HOW SUPPLIED

- The Wand is supplied sterile for SINGLE USE ONLY. The Wand is sterile only if the packaging is not opened or damaged. If the packaging is damaged, do not use the Wand.

- Do not use the Wand after the "use-by" date printed on the label. Wand performance and patient safety may be compromised if used after the expiration date.

- The Wand is sterilized using ethylene oxide.

- Wand materials are not made with DEHP or natural rubber latex.

- Rated voltage for the FLOW 90 Wand is 340 Vrms.

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MAINTENANCE AND TROUBLESHOOTING

Refer to the User Manual for the WEREWOLF Controller for operating and troubleshooting guidelines for the Controller.

CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

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All questions, complaints, and return authorization requests should be directed to

ArthroCare Corporation
7000 West William Cannon Austin,
TX 78735 U.S.A. (CША)
(800) 343-5717
ОБСЛУЖИВАНИЕ КЛИЕНТОВ
www.arthrocare.com

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen Germany
ec.rep@smith-nephew.com

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Customer Service or an authorized representative



ArthroCare Corporation
7000 West William Cannon
Austin, TX 78735 U.S.A.
(800) 343-5717
CUSTOMER SERVICE
www.arthrocare.com

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

SYMBOLS KEY

	Catalogue number
	Lot number
	Use-by Date
	Consult Instructions for Use
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Keep dry
	Keep Away From Sunlight

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™Trademark of Smith & Nephew.

These products may be covered by one or more U.S. patents.

See smith-nephew.com/patents for details.

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The FLOW 50 Wand (Wand) is a bipolar, RF electrosurgical device used for specific indications in orthopedic and arthroscopic surgeries. The Wand utilizes the AMBIENT feature which provides accurate ($\pm 3^{\circ}\text{C}$) real-time temperature monitoring of the intraarticular irrigating fluid. The Wand is capable of operating in five Modes.

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The Wand features AMBIENT™ technology which provides accurate ($\pm 3^{\circ}\text{C}$ between $20-60^{\circ}\text{C}$) real-time temperature monitoring of the intra-articular irrigating fluid. The Wand is designed for use up to 7 minutes of cumulative active ablation and is to be used only with the WEREWOLF™ Controller (Controller). The Wand is capable of operating in five Modes.

The Smith & Nephew FLOW 50 and FLOW 90 Wands are a bipolar, high frequency electrosurgical devices used for specific indications in orthopedic and arthroscopic surgeries.

The FLOW 50 and FLOW 90 consists of an active electrode, a ceramic spacer, a stainless steel shaft, internal suction tubing, a handle with integrated finger controls, electrical cable, and suction tube set. The electrical cable and suction tube set connect to the WEREWOLF Controller. The proximal end of the suction tube set connects to a male fitting of standard hospital vacuum systems. The wand is controlled using the three finger switches, integrated into the handle that enable the Ablation Mode switching (Lo, Med, Hi) as well as activation of the Wand (Vac, Coag, or Ablate).

The technological characteristics of the FLOW 90 Wand are the same as the predicate FLOW 50 Wand. FLOW 50 Wand and FLOW 90 Wand are identical in intended use, fundamental scientific technology, and principle of operation. The wands vary by the electrode configuration, the tubing/shaft angle near the wand tip, and the patient-contacting materials.

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

Table 2: Modes and Levels for the FLOW 50 and FLOW 90 Wands:

Action	Mode	Level
Ablation	Lo 	Lo Minus (-)
		Lo (•)
		Lo Plus (+)
Ablation	Med 	Med Minus (-)
		Med (•) *
		Med Plus (+)
Ablation	Hi 	Hi Minus (-)
		Hi (•)
		Hi Plus (+)
Coagulation	 COAG	Coag (•) *
		Coag Plus (+)
Vacuum	Vac 	Not Adjustable

*Default Modes and Levels for FLOW 50 and FLOW 90

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The following labeling of FLOW 50 /FLOW 90 Wands is applied on the primary packaging template in Russia:

- LOT number/ BATCH code;
 - Name;
 - Quantity;
 - Shelf life (sign "Use by..."); catalogue number; following the packaging;
 - Information about manufacture;
 - Information about authorized representative;
 - Sterilized method: EO. Sterile. Non-toxic.
 - Number of license;
 - Indication, method of application, contraindications – following the IFU;
 - Warnings (following the IFU);
 - Storage conditions;
 - Information on compliance with the technical regulations of the Customs Union;
 - IFU : <http://www.smith-nephew.com/about-us/where-we-operate/europe/russia/ifus/>;
- Information about the name of the medical product is provided in accordance with the registration certificate.

Labeling in Russian is given above put so that not to block marking of the manufacture on consumer packing, applied directly on consumer packing each medical product in the package.

The following symbols may appear on the packaging and in the documentation of sterile wands FLOW 50 /FLOW 90.

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Interpretation of symbols used on the package:

	Catalogue number
	Lot number
	Use-by Date
	Contents
	Quantity
	Consult Instructions for Use
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Keep dry
	Keep Away From Sunlight
	Caution
	Sterilized Using ethylene oxide

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Environmental conditions for surgical intervention

Temperature	10 °C–40 °C
Humidity	10 %–75 % (without condensation)
Pressure	70 kPa–106 kPa

Operation conditions (for use in the patient's body)

Temperature	32 °C–42 °C
Humidity	До 100%
Pressure	70 kPa–106 kPa

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Recommended transportation and storage ranges are presented below.

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The device is subject to the Directive 2002/96 / EC (European Community) for waste electrical and electronic equipment (WEEE). In accordance with the registration status of this directive, this device cannot be disposed of with ordinary electrical and electronic equipment.

After use, the product may be potentially hazardous. Dispose of the product in accordance with local laws and regulations.

The manufacturer provides product safety information upon request.

The class of medical waste products that had contact with the internal environment of the body – B.

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Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. DO NOT REUSE. **THIS WARRANTY IS IN LIEU OF**

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(7000 West William Cannon Drive Austin TX 78735-8531, USA)

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