State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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1. Country (Pays): **United States of America**

This public document (Le présent acte public) **Biological Indicators for Sterilization Quality Control**

2. has been signed by (a été signé par) **Ryan J. Bartz**

3. acting in the capacity of (agissant en qualité de) **Notary Public, State of Minnesota**

4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Ryan J. Bartz, Notary Public, State of Minnesota

Certified Attesté

5. at (à) **St. Paul, Minnesota**

6. Date (Date) **11/15/2019**

7. by (par) **Secretary of State, State of Minnesota**

8. File No (Sous n°) **1118405400107**

9. Seal / Stamp (Sceau / Timbre):

10. Signature (Signature):

Steve Simon

Steve Simon

Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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
To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that attached IFU for:

**Biological Indicators for Sterilization Quality Control:
Attest Rapid Readout Biological Indicators designed for rapid and reliable monitoring
of the 270°F (132°C), gravity steam sterilization process.**

in English language is truthful and accurate


Regina Feferman
International Regulatory Affairs Associate
3M Health Care

11/14/19


Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 14th day of November 2019.

Ryan Bartz, Notary Public
My Commission Expires: January 31, 2024

Instruction for Use

Biological Indicators for Sterilization Quality Control:

Attest Rapid Readout Biological Indicators designed for rapid and reliable monitoring of the 270°F (132°C), gravity steam sterilization process.

Product Description

The Attest Rapid Readout Biological Indicators for rapid and reliable monitoring of the 270°F (132°C), gravity steam sterilization process (3M Attest 1291 Rapid Readout Biological Indicator (blue cap)) is a dual readout biological indicator system specifically designed for rapid and reliable monitoring of the 270°F (132°C), gravity steam sterilization process when used in conjunction with either the 3M Attest 290 Auto-reader or the 3M Attest Auto-reader 390. The sterilizer must have a 1.5 minute come-up time to kill the 3M Attest 1291 Rapid Readout Biological Indicator in 3 minutes. If the come-up time is < 1.5 minutes, extend the cycle time to 4 minutes to achieve consistent kill.

The 3M Attest 1291 Rapid Readout Biological Indicator detects the presence of *Geobacillus stearothermophilus* by detecting the activity of alpha-glucosidase, an enzyme present within the organism. The presence of the enzyme is detected by reading fluorescence produced by the enzymatic breakdown of a non-fluorescent substrate. This creates a fluorescence change, which is detected by the auto-reader. A fluorescence change indicates a steam sterilization process failure. The 3M Attest 1291 Rapid Readout Biological Indicator also detects the presence of *G. stearothermophilus* organisms by a visual color change reaction. Biochemical activity of the *G. stearothermophilus* organism produces acid by-products that cause the media to change color from purple to yellow. A visual pH color change also indicates a steam sterilization process failure. Due to the high sensitivity of the 1-hour fluorescent results, however, there is no advantage to incubating the Attest 1291 Rapid Readout Biological Indicator beyond 1 hour.

Intended Use

Sterilization quality control

Scope of use: Sterilization

Contraindications

None.

Potential consumer

Only for professional use

Possible side effects

No known possible side effects.

Warnings

There is a glass ampoule inside the plastic vial of the biological indicator. To avoid the risk of serious injury from flying debris due to a ruptured biological indicator:

- Allow the biological indicator to cool for the recommended time period before crushing. Crushing or excessive handling of the biological indicator before cooling may cause the glass ampoule to burst.
- Wear safety glasses and heat-resistant gloves when removing the biological indicator from the sterilizer.
- Wear safety glasses when crushing the biological indicator.
- Handle the biological indicator by the cap when crushing and tapping.
- Do not use your fingers to crush the glass ampoule.
- Do not roll the biological indicator between fingers to wet the spore strip.

Precautions

Do not use the 3M Attest 1291 Rapid Readout Biological Indicator to monitor sterilization cycles which it is not designed to challenge:

1. 132°C (270°F) or 121°C (250°F) vacuum-assisted steam sterilization cycles.

2. 121°C (250°F) gravity steam sterilization cycles.
3. Dry heat, chemical vapor, ethylene oxide or other low temperature sterilization processes.

Monitoring Frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that every steam sterilization load be monitored with an appropriate biological indicator.

Directions for Use

1. Identify the Attest rapid readout biological indicator by writing the sterilizer and load number and processing date on the indicator label. Do not place another label or indicator tape on the vial or on the cap.

2. Place the Attest rapid readout biological indicator in an appropriate process challenge device according to recommended practices. Do not place the Attest rapid readout biological indicator in direct contact with a chemical indicator. Fluorescent residue could transfer to the rapid readout biological indicator and affect the result.

Appropriate process challenge devices for specific loads containing:

- A. Unwrapped metal instruments or hard goods run at 132°C (270°F) for 3 minutes in a gravity displacement cycle:

- Place an Attest rapid readout biological indicator in an unwrapped instrument tray with a representative number and type of instruments normally processed. AAMI suggests placing a biological indicator in an empty tray.

- B. Unwrapped metal instruments or hard goods with porous items run at 132°C (270°F) for 10 minutes in a gravity displacement cycle:

- Place an Attest rapid readout biological indicator in an unwrapped instrument tray with a representative number and type of instruments normally processed. AAMI suggests placing a biological indicator in an empty tray but include porous items if applicable.

- C. Wrapped metal instruments or hard goods run at 132°C (270°F) for 10 minutes in a gravity displacement cycle:

- Place an Attest rapid readout biological indicator in a wrapped instrument tray with a representative number and type of instruments normally processed. AAMI suggests placing a biological indicator in an empty tray, include porous items if applicable.

- D. Container systems run at 132°C (270°F) in a gravity displacement cycle:

- Place an Attest rapid readout biological indicator in the area determined by product testing to provide the greatest challenge to the sterilization process.

- E. Patient Care trays run at 132°C (270°F) for 3 minutes a gravity displacement cycle:

- Place an Attest rapid readout biological indicator in a tray or container of medical devices that is going directly to patient use to ensure that the medical devices are properly sterilized. Follow aseptic techniques when retrieving the Attest rapid readout biological indicator from the tray or container for incubation.

3. Place the process challenge device in the most challenging area for the sterilant. This is typically on the bottom shelf, near the door and over the drain.

4. Process the load according to recommended practices.

5. After completion of the cycle, fully open the sterilizer door for a minimum of 5 minutes prior to removing the Attest rapid readout biological indicator.

6. When the Attest rapid readout biological indicator is not contained in a test pack or other heat absorbing packaging material, remove the Attest rapid readout biological indicator from the sterilizer and allow to cool for an additional 10 minutes prior to crushing.

7. When the Attest rapid readout biological indicator is contained in a test pack or other heat absorbing packaging material, the test pack or any other heat absorbing packaging material should be removed from the sterilizer and opened up for 5 minutes to dissipate heat prior to removing the biological indicator. Then allow the biological indicator to cool outside the test pack for an additional 10 minutes prior to crushing.

8. Check the throughput chemical indicator (CI) on the label of the Attest rapid readout biological indicator. A color change from rose to brown/black confirms that the biological indicator has been exposed to the steam sterilization process. This CI color change does not indicate that the sterilization process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process monitoring controls and investigate placement of the biological indicator in the sterilizer.

9. While wearing safety glasses, press the cap down. Crush the glass ampoule of the biological indicator in the crusher well of the auto-reader. Hold the biological indicator by the cap and tap on a hard surface until media wets strip at bottom of the vial, and then place the biological indicator in an auto-reader incubation well configured to incubate Attest 1291 Rapid Readout Biological Indicators. See Attest Auto-reader Operator's Manual for further details.

10. Each day that a processed Attest 1291 Rapid Readout Biological Indicator is incubated, crush, tap and incubate at least one non-processed Attest 1291 Rapid Readout Biological Indicator to use as a positive control. Write a "C" (for "control") and the date on the label. The positive control should be from the same manufacturing date and lot number as the processed biological indicator. Incubating a positive control helps ensure:

- correct incubation temperatures are met,
- viability of spores have not been altered due to improper storage temperature, humidity or proximity to chemicals,
- capability of media to promote rapid growth, and
- proper functioning of Attest Auto-reader.

11. Incubation and Reading:

Incubate the positive control and processed Attest 1291 Rapid Readout Biological Indicators for 1 hour at $60 \pm 2^{\circ}\text{C}$ ($140 \pm 3^{\circ}\text{F}$) in a 3M Attest 290 Auto-reader or a 3M Attest Auto-reader 390. See the applicable auto-reader Operator's Manual for the proper use of this equipment. The auto-readers automatically take readings and may indicate a positive result in less than 1 hour. The final fluorescent negative biological indicator reading is made at 1 hour. After the final reading is obtained the processed biological indicator may be discarded.

The processed biological indicator and the positive control may also be further incubated at 60°C for a visual pH color change result. Examine the biological indicator for early detection of positive results (i.e., media turns yellow) at convenient time intervals such as 6, 12 and 18 hours. The final negative reading (media remains purple) for a visual pH color change is made at 24 hours. The positive control should show a yellow color change of the growth media within 24 hours.

Interpretation of Results:

Fluorescent Result

The positive (unprocessed) control biological indicator must provide a positive result [red light on the 3M Attest 290 Auto-reader or plus symbol (+) on the LCD display of the 3M Attest Auto-reader 390]. If the positive biological indicator control reads negative [green light on the 3M Attest 290 Auto-reader or minus symbol (-) on the LCD display of the 3M Attest Auto-reader 390], refer to the applicable auto-reader Operator's Manual Troubleshooting section. Retest the auto-reader with a new positive control biological indicator. The processed biological indicator results are not valid until the positive biological indicator control reads positive.

With a processed biological indicator, a positive [red light on the 3M Attest 290 Auto-reader or plus symbol (+) on the LCD display of the 3M Attest Auto-reader 390] means a sterilization process failure has occurred. A negative [green light on the 3M Attest 290 Auto-reader or minus symbol (-) on the LCD display of the 3M Attest Auto-reader 390] after 1 hour of incubation indicates an acceptable sterilization process.

pH Color Change Result (Optional)

The appearance of a yellow color in the processed indicator demonstrates bacterial growth and a sterilization process failure. No color change (i.e. media remains purple) indicates an adequate sterilization process. A final negative result is made after 24 hours of incubation. The positive control indicator should show a color change from purple to yellow for the processed indicator results to be valid.

12. Immediately act on any positive biological indicator results. Always retest the sterilizer and do not use sterilizer for processing loads until three consecutive biological indicator results are negative.

Validation of Reduced Incubation Time (Readout Reliability Data)

The 3-hour and 48-hour incubation times have been correlated with a 7-day incubation period. Sterilized indicators were examined daily for detection of a visual pH color change. The 1-hour fluorescence change reading and the 24-hour visual pH color change readings were compared to the 7-day visual pH color change readings to determine the readout reliability of the indicator. Readout reliability of the Attest™ 1291 rapid readout biological indicator was determined using the sensitivity calculation described below:

Sensitivity =

$$\frac{(\text{Number of Growth Positives after 168 hours}) - (\text{Number of False Negatives}) \times 100}{\text{Number of Growth Positives after 168 Hours}}$$

Attest 1291 Rapid Readout Biological Indicators

132°C (270°F) Gravity Displacement Steam Sterilization Process Validation of Reduced Incubation Time - Readout Reliability Summary

Sterilization Process	Incubation Temperature	# Tested	# Growth Positives at 168_hours	Growth		Fluorescence	
				#_False Negatives at 24 hours	Sensitivity at 24 hours	# False Negatives at 60 minutes	Sensitivity at 60 minutes
132°C (270°F) Gravity Displacement	60°C (140°F)	880	533	13	97.6%	0	100%

These data demonstrate that ≥97% of the 7-day (i.e. 168 hours) growth visual positives were detected by fluorescence within 1 hour of incubation and by the visual pH color within 24 hours of incubation. The 3M™Attest™1291 Rapid Readout Biological Indicator therefore meets readout reliability of ≥97% for the 1-hour fluorescence results and the 24-hour visual color change results.

Storage and Shelf Life

- Best stored under normal room conditions: 15-30°C (59-86°F) and 35-60% relative humidity.
- Do not store these biological indicators near sterilants or other chemicals.
- Attest 1291 rapid readout biological indicators have a 2-year shelf life from the date of manufacture. The expiration date is indicated on the BI and packaging by the 4-digit year and 2-digit month of expiration (e.g. 2017-11).
- All of the information to the right of the lot-in-a-box symbol indicates the lot number (e.g.,  2017-11AD).

Transportation: During transportation, the products can be exposed to short-term temperature changes, which do not affect the safety and products specifications.

Usage: at room temperature 15-30 °C. Values of humidity and pressure are not regulated.

Disposal

Dispose of used Attest biological indicators according to your healthcare facility's policy. You may wish to autoclave any positive biological indicators at 250°F/121°C for at least 15 minutes, or at 270°F/132°C for 10 minutes in a gravity displacement steam sterilizer, or at 270°F/132°C for 4 minutes in a vacuum assisted steam sterilizer.

For processed BIs observed to have a broken media ampoule after the sterilization process, or if the ink of the CI appears smeared or runny (if included) leave the BI and CI in the sterilization pouch and place the pouch in packaging compatible with steam sterilization and steam sterilize per the parameters above.

Packaging/contents should be disposed in accordance with the local/regional/national/ international

regulations.

Medical products should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class B).

Maintenance and repair data

Not applicable. Single use only.

Sterilization

These products are not supplied sterile.

Data on cells, tissues and animal organs contained in medical product

Drugs: Medical product does not contain and does not deliver drugs

Human blood derivatives: Medical product does not contain and does not deliver human blood derivatives.

Materials of animal origin: materials of human origin are not intentionally added in the product.

Warranties

The products are evaluated at coated spore jumbo stage (after spores coated to strips that will be placed in the indicators as spore strips). The data generated from jumbo evaluation are used for the quality certificates.

Content of packaging

50 BIs/box, 4 boxes/case

Labeling

Products are supplied in a package with labelling containing the following information: product name and reference number; lot code and shelf life, name and address of the manufacturer; contents, storage/transportation conditions; trademarks information; and symbols – Attention! Refer to Operating Documents, incubation time.

Technical characteristics

Glass ampoule: Diameter 6.4 ± 0.2 mm, Height $39.6 + 1.0 / - 0.8$ mm

Plastic sleeve: Diameter 8.5 ± 0.05 mm, Height 50.8 ± 0.5 mm

Plastic cap: Length 16.6 ± 0.13 mm, Diameter 10.7 ± 0.13 mm

Legal manufacturer

3M Company, 3M Health Care, USA

3M Center, 2510 Conway Avenue, Bldg. 275-5W-06, Saint Paul, Minnesota, 55144, USA

Authorized representative on the territory of the Russian Federation:

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10. Signature (Signature):

Steve Simon

Steve Simon
Secretary of State
State of Minnesota

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
Biological Indicators for Sterilization Quality Control:
Attest Biological Indicators for ethylene oxide (E0) sterilization.
Attest Biological Indicators for ethylene oxide (E0) sterilization in Test Pack

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

Attest Biological Indicators for ethylene oxide (EO) sterilization.

Attest Biological Indicators for ethylene oxide (EO) sterilization in Test Pack

Product Description

The 3M Attest 1264 Biological Indicator (green cap) is designed for monitoring of ethylene oxide (EO) sterilization processes. The presence of *Bacillus atrophaeus* spores is detected by a visual color change (media turns yellow). The yellow color change indicates a sterilization process failure. The final readout of a negative result (media remains green) is made after 48 hours of incubation.

The Attest Biological Indicator for EO in Test Pack (1278) is specifically designed to routinely challenge ethylene oxide (EO) sterilization processes. This convenient disposable pack is designed to be as resistant as the routine EO biological indicator test pack recommended by AAMI.

Each pack contains a 3M Attest 1264 Biological Indicator (green cap). The chemical indicator on the Attest biological indicator label and the chemical indicator dot on the pack turn from red to green when EO processed. Attest biological indicator controls are provided.

The presence of *Bacillus atrophaeus* spores is detected by a visual color change (media turns yellow). The yellow color change indicates a sterilization process failure. The final readout of a negative result (media remains green) is made after 48 hours of incubation.

Intended Use

Sterilization quality control

Scope of use: Sterilization

Contraindications

None.

Potential consumer

Only for professional use

Possible side effects

No known possible side effects.

Warnings

- Butyl, neoprene, or nitrile gloves should be worn when removing the biological indicator from the load prior to the completion of aeration.
- There is a glass ampule inside the plastic vial of the biological indicator.
 - o Crushing or excessive handling of the biological indicator before cooling may cause the glass ampule to burst.
 - o Wear safety glasses when crushing the biological indicator.
 - o Handle the biological indicator by the cap when crushing.
 - o Do not use your fingers to crush the glass ampule.

Precautions

Do not use the Attest 1264 biological indicator to monitor sterilization cycles which it is not designed to challenge:

1. Steam sterilization cycles.
2. Dry heat, chemical vapor, or other low temperature sterilization processes.

Precautions Test pack

To ensure the test pack delivers the intended challenge:

- DO NOT open pack prior to sterilization;

- DO NOT reuse test pack.

Do not use the Attest 1278 EO Pack to monitor sterilization cycles which it is not designed to challenge:

1. Steam sterilization cycles;
2. Dry heat, chemical vapor, or other low temperature sterilization processes.

Monitoring Frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that an Attest biological indicator be used to monitor every ethylene oxide sterilization load.

Directions for Use

1. Identify the Attest biological indicator by writing the sterilizer, load number, and processing date on the indicator label.
2. Place the Attest 1264 biological indicator in a Process Challenge Device (BI PCD) as recommended by professional association guidelines or national standards for hospital practice (e.g., AAMI routine BI syringe test pack).
3. Place the BI PCD in a full load in the most challenging area for the sterilant. This is generally in the center of the load.
4. Process the load as usual.
5. After sterilization two options exist:
 - a. Aeration of BI PCD Prior to Incubation
Aerate the other packs from the load for the same amount of time as the BI PCD. At the end of the aeration cycle, remove the biological and chemical indicators from the BI PCD. Dispose of the remaining BI PCD components according to the healthcare facility's policy.
 - b. Removal of Biological Indicator from BI PCD Prior to Aeration
Open the sterilizer door according to the manufacturer's instructions; transfer the load to the aerator and retrieve the BI PCD. (If a sterilizer with aeration capability is used, follow the manufacturer's instructions for opening the door and remove the BI PCD prior to the aeration cycle.) Remove the biological and chemical indicator from the BI PCD. Return the remainder of the BI PCD to the load for aeration according to the healthcare facility's policy. When it is necessary to handle items that are not fully aerated, butyl, neoprene, or nitrile gloves should be worn. The breathing zone of personnel should be monitored to verify the safety of the practices followed.
6. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to an EO sterilization process. This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
7. Incubate the biological indicator at $37 \pm 2^{\circ}\text{C}$ ($99 \pm 3^{\circ}\text{F}$).

Attest Incubator 240 volt (International Usage)

Model 129 (28 indicators)

- a. Position indicator in metal block (see Figure 1). Place bottom of the indicator into the incubator's metal heating block so that the indicator is at an angle of approximately 45° .


Figure 1

- b. Push the indicator straight back (see Figure 2). This crushes the media ampule and activates the indicator. Be sure that the cap will remain above the metal block when the indicator is pushed back.


Figure 2

- c. Push the activated indicator down to seat it in the metal heating block. (See Figure 3). Be sure that the cap remains above the metal block when seated in the incubator.


Figure 3

8. Incubate at least one unprocessed Attest biological indicator (positive control) each day a processed indicator is incubated. The positive control indicator should be from the same manufacturing date and lot number as the processed indicator in the incubator.
9. Write a "C" and a date on the label of the positive control indicator. Crush and incubate the control at $37 \pm 2^\circ\text{C}$ ($99 \pm 3^\circ\text{F}$).

The purpose of the positive control is to ensure:

- correct incubation conditions
- viability of indicators (incorrect storage conditions could adversely affect even those indicators which are within their stated shelf life)
- capability of media to promote rapid growth.

10. Incubate processed and control biological indicators for 48 hours at $37 \pm 2^\circ\text{C}$ ($99 \pm 3^\circ\text{F}$).

Incubation Times:

Early Detection	12 hours
	18 hours
	24 hours

Final Detection 48 hours

11. The appearance of a yellow color in the processed indicator demonstrates bacterial growth and a sterilization process failure. No color change indicates an adequate sterilization process. A final

negative result is made after 48 hours of incubation. The positive control indicator should show a yellow color change for the processed indicator results to be valid.

12. Record the processed and control biological indicator results. Act on any positive test as soon as the first evidence of growth is noted. Always retest the sterilizer and do not use the sterilizer until the biological indicator test results are negative.

Directions for Use (Test Pack)

1. Use an Attest EO pack to test loads containing wrapped telescopic instruments, plastic, rubber, metal instruments, and equipment.

2. Place the EO pack in a full load in the most challenging area for the sterilant. This is generally in the center of the load.

3. Process the load as usual.

4. After sterilization two options exist:

- A. Aeration of EO pack prior to incubation

Aerate the other packs from the load for the same amount of time as the test pack. At the end of the aeration cycle, check to see that the chemical indicator dot on the test pack turned from red to green which indicates the pack has been exposed to the EO

sterilization process. If the chemical indicator dot is unchanged, check the sterilization process. Remove the biological indicator from the test pack. Dispose of the remaining test pack components according to the healthcare facility's policy.

- B. Removal of biological indicators from EO pack prior to aeration

Open the sterilizer door according to the manufacturer's instructions; transfer the load to the aerator and retrieve the test pack. If a sterilizer with aeration capability is used, follow the manufacturer's instructions for opening the door and remove the test pack prior to the aeration cycle. Check to see that the chemical indicator dot on the test pack turned from red to green which indicates the

pack has been exposed to the EO sterilization process. If the chemical indicator dot is unchanged, check the sterilization process.

Remove the biological indicator from the test pack. Return the remainder of the test pack to the load for aeration according to the healthcare facility's policy. When it is necessary to handle items that are not fully aerated, butyl, neoprene, or nitrile gloves should be worn. The breathing zone of personnel should be monitored to verify the safety of the practices followed.

5. Check the chemical indicator on the label of the biological indicator. A color change from red to green indicates the pack has been exposed to the EO sterilization process. The color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator on the label is unchanged, check the sterilization process.

6. Identify the Attest biological indicator by writing the sterilizer, load number, and processing date on the indicator label.

7. Discard the test pack. Using a test pack more than once will invalidate subsequent test results.

8. Incubate the biological indicator at $37 \pm 2^{\circ}\text{C}$ ($99 \pm 3^{\circ}\text{F}$).

Attest Incubator 220/240 volt

(International Usage)

Model 129 (28 indicators)

- A. Position indicator in metal block (see Figure 1). Place bottom of the indicator into the incubator's metal heating block so that the indicator is at an angle of approximately 45° .

- B. Push the indicator straight back (see Figure 2). This crushes the media ampule and activates the indicator. Be sure that the cap will remain above the metal block when the indicator is pushed back.

- C. Push the activated indicator down to seat it in the metal heating block (see Figure 3). Be sure that the cap remains above the metal block when seated in the incubator.


Figure 1


Figure 2


Figure 3

9. Incubate at least one non-processed Attest biological indicator (positive control) each day a processed indicator is incubated. The positive control indicator should be from the same lot code as the processed indicator.

10. Write a "C" and a date on the label of the positive control indicator. Crush and incubate the control at $37 \pm 2^\circ\text{C}$ ($99 \pm 3^\circ\text{F}$).

The purpose of the positive control is to ensure:

- correct incubation conditions
- viability of indicators (incorrect storage conditions could adversely affect even those indicators which are within their stated shelf life)
- capability of media to promote rapid growth

11. Incubate the processed and control biological indicators for 48 hours at $37 \pm 2^\circ\text{C}$ ($99 \pm 3^\circ\text{F}$) for a visual color

change readout. Examine the indicators for early detection of positive results (media turns yellow) at convenient time intervals such as 18 and 24 hours. The final negative reading (media remains green) is made after 48 hours of incubation. The positive control should show a yellow color change in the media. If the positive control remains green, check the incubator temperature and storage conditions. Retest the incubator with a new positive control. The processed indicator results are not valid until the positive control shows a yellow color change.

With a processed indicator, a visual color change to yellow indicates bacterial growth and a sterilization process failure. No color change (media remains green) after 48 hours of incubation indicates an adequate sterilization process.

12. Record the processed and control biological indicator results. Act on any positive test as soon as the first evidence of growth is noted. Always retest the sterilizer and do not use the sterilizer until the biological indicator test results are negative.

Storage/ Shelf Life

- Best stored under normal room conditions: $15-30^\circ\text{C}$ ($59-86^\circ\text{F}$), 35-60% relative humidity.
- Do not store these biological indicators near sterilants or other chemicals.
- Attest biological indicators have a 24 month shelf life. The expiration date is indicated on the BI and packaging by an hourglass symbol followed by the 4-digit year and 2-digit month of expiration (e.g. ⌚ 2017-11).
- **Transportation:** During transportation, the products can be exposed to short-term temperature changes, which do not affect the safety and products specifications. Transportation of medical device is carried out by all type of transport.
- **Usage:** at room temperature $15-30^\circ\text{C}$. Values of humidity and pressure are not regulated.

Disposal

Dispose of used Attest biological indicators according to your healthcare facility's policy. You may wish to autoclave any positive biological indicators at $250^\circ\text{F}/121^\circ\text{C}$ for at least 15 minutes, or at $270^\circ\text{F}/132^\circ\text{C}$ for

10 minutes in a gravity displacement steam sterilizer, or at 270°F/132°C for 4 minutes in a vacuum assisted steam sterilizer.

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Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations.

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Content of packaging

100 BIs/box, 4 boxes/case

Test pack: 25 packs and 25 controls/box, 2 boxes/case

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
**Biological Indicators for Sterilization Quality Control:
Attest Rapid Readout Biological Indicator for Vaporized Hydrogen Peroxide
Sterilization**

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

Attest Rapid Readout Biological Indicator for Vaporized Hydrogen Peroxide Sterilization

Product Description

The Attest Rapid Readout Biological Indicator for Vaporized Hydrogen Peroxide Sterilization (3M Attest Rapid Readout Biological Indicator 1295 (pink cap, referred to hereinafter as the 1295 BI)) is a self-contained biological indicator specifically designed for rapid and reliable routine monitoring of vaporized hydrogen peroxide sterilization processes when used in conjunction with the 3M Attest Auto-reader 490H or the 3M Attest Auto-reader 490 having software version 4.0.0 or greater (hereinafter referred to as the 490H Auto-reader or 490 Autoreader having software version 4.0.0 or greater). The 1295 BI is a single-use device.

A schematic illustrating the design of the 1295 BI is provided in Figure 1. The self-contained design includes a carrier with spores of *Geobacillus stearothermophilus* and a media ampoule containing bacteriological growth medium which meets the requirements for growth promoting ability specified in ANSI/AAMI/ISO 11138- 1:2006/(R)2010. The spore carrier and media ampoule are contained in a plastic vial topped with a pink cap. A chemical process indicator printed with stripes which change in color from blue towards pink upon exposure to vaporized hydrogen peroxide is located on the top of the cap.

The 1295 BI utilizes the α -glucosidase enzyme system, which is generated naturally within growing cells of *Geobacillus stearothermophilus*. The α -glucosidase in its active state is detected by measuring the fluorescence produced by the enzymatic hydrolysis of a non-fluorescent substrate, 4-methylumbelliferyl- α -D-glucoside (MUG). The resultant fluorescent by-product, 4-methylumbelliferone (MU), is detected in the Auto-reader. The presence of fluorescence within the specified incubation time for the 1295 BI in the Auto-reader indicates a sterilization process failure.

The 1295 BI can also indicate the presence of *G. stearothermophilus* organisms by a visual pH color change reaction. Biochemical activity of the *G. stearothermophilus* organism produces metabolic byproducts that cause the media to change color from purple to yellow which also indicates a sterilization process failure. Use of this indication method is optional and is typically restricted to special studies.


Figure 1: Design of the 3M Attest Rapid Readout Biological Indicator

Readout Time

The rapid readout result has been correlated with a 7-day visual pH color change result following the FDA's Reduced Incubation Time protocol. The time to result is determined by the software version programmed on the Auto-reader.

24-minute Fluorescent Result

1295 BIs incubated in a 490H or 490 Auto-reader having software version 4.0.0 or greater have a 24 minute reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

4-hour Fluorescent Result

1295 BIs incubated in a 490H Auto-reader having software versions less than 4.0.0 have a 4-hour reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

Due to the high reliability of the fluorescent result there is no advantage to incubating 1295 BIs after the fluorescent result has been determined by the 490H Auto-reader or 490 Auto-reader having software version 4.0.0 or greater and recorded.

Intended Use

Sterilization quality control

Scope of use: Sterilization

Indications for Use

Use the 3M Attest Rapid Readout Biological Indicator 1295 in conjunction with the 3M Attest Auto-reader 490H or 490 Auto-reader having software version 4.0.0 or greater as a standard method of routine monitoring of vaporized hydrogen peroxide sterilization processes in the following systems:

AMSCO V-PRO 1 Low Temperature Sterilization System (Lumen cycle), AMSCO V-PRO 1 Plus Low Temperature Sterilization System (Lumen and Non Lumen cycles), AMSCO V-PRO maX Low Temperature Sterilization System (Lumen, Non Lumen, and Flexible cycles), AMSCO V-PRO 60 Low Temperature Sterilization System (Lumen, Non Lumen and Flexible cycles) and in STERRAD 100S, STERRAD NX (Standard and Advanced cycles), STERRAD 100NX (Standard, Flex, Express and Duo cycles) systems, STERRAD NX with AllClear Technology (Standard and Advanced cycles) and STERRAD 100NX with AllClear Technology (Standard, Flex, Express and Duo cycles).

Contraindications

None.

Potential consumer

Only for professional use

Possible side effects

No known possible side effects.

Warnings

There is a glass ampoule inside the plastic vial of the biological indicator (BI). To avoid the risk of serious injury from peroxide burns:

- Wear safety glasses and gloves when removing the 1295 BI from the sterilizer.
- Wear safety glasses and gloves when activating the 1295 BI.
- Handle the 1295 BI by the cap when crushing or flicking.

Residual hydrogen peroxide may be trapped within the 1295 BI if the media ampoule is damaged during the sterilization process. If a broken ampoule is observed after processing, avoid direct contact with the 1295 BI as it may result in hydrogen peroxide burns. Follow the disposal instructions provided at the end of this document.

Precautions

1. Do not use the 1295 BI to monitor sterilization cycles which it is not designed to challenge:

- a. Steam sterilization cycles;
- b. Dry heat sterilization cycles; or
- c. Ethylene oxide sterilization processes.

2. To reduce the risk associated with incorrect results:

- Before sterilization, inspect 1295 BI to verify media ampoule is intact and process indicator stripes are blue. Do not use any 1295 BIs which have a broken media ampoule or process indicator stripes which are not blue.
- Do not place tape or labels on 1295 BI prior to sterilization or incubation in the Auto-reader.
- Activate and incubate the 1295 BI within 1 hour after the completion of the sterilization cycle.
- Do not incubate a 1295 BI if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new biological indicator.
- After 1295 BI activation, ensure media has flowed to the spore growth chamber.

3. To ensure the product functions as intended throughout the labeled shelf life, store 1295 BIs in the resealable foil pouch until use.


Monitoring Frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that every sterilization load be monitored with a biological indicator.

Directions for Use

1. Remove 1295 BI from sealed foil pouch, then reseal foil pouch if other 1295 BIs remain in foil pouch. Do not place any labels or indicator tape on the vial or on the cap.
2. Place the 1295 BI in a sterilization pouch indicated for use in vaporized hydrogen peroxide sterilization processes. Seal the sterilization pouch.
3. Place the pouched BI in the sterilizer chamber, with the white side of the pouch facing up and the clear plastic side facing down. When there is adequate space in the loaded sterilizer chamber, place the pouched BI directly on the sterilizer chamber rack or shelf. The sterilizer manufacturer should be consulted to identify the area of the chamber to place the BI.
4. Process the load according to recommended practices.
5. After completion of the cycle, don safety glasses and gloves and remove the pouched BI from the sterilizer. Inspect the 1295 BI to verify the media ampoule is intact. If a chemical indicator was included in the pouch with 1295 BI, inspect the CI to assure the ink of the CI is not smeared or runny. If the BI media ampoule is intact and the ink of the CI (if included) appears typical remove the BI from the sterilization pouch and proceed to Step 6. If the media ampoule is broken or if the ink of the CI appears smeared or runny, leave them in the sterilization pouch and follow the disposal instructions. Retest the sterilizer using a new 1295 BI and CI (if included).
6. Check the process indicator on the top of the cap of the 1295 BI. A color change of the stripes from blue towards pink confirms that the 1295 BI has been exposed to the vaporized hydrogen peroxide sterilization process. This color change does not indicate that the process was sufficient to achieve sterility. If the process indicator is unchanged, check the sterilizer physical monitors.
7. Identify the 1295 BI by writing the load number, sterilizer, and processing date on the indicator label.
8. Activate the 1295 BI.

While wearing safety glasses and gloves, place the 1295 BI in an Attest Biological Indicator Activator. Close and squeeze the activator to close the 1295 BI cap and crush the media ampoule (see pictures at right).

Immediately remove the BI and flick it (see of the vial. If the media hasn't filled the growth chamber, hold the BI by the cap and flick it until media fills the growth chamber. Place the activated 1295 BI in a 490H Auto-reader incubation well which is color-coded pink or 490 Auto-reader having software version 4.0.0 or greater incubation well which is color-coded black and wait for the result. See the Autoreader Operator's Manual for further information related to its use.

NOTE: Activate and incubate the 1295 BI within 1 hour of the completion of the sterilization cycle.

9. Each day that a processed 1295 BI is incubated, activate and incubate at least one nonprocessed 1295 BI to use as a positive control. Follow the activation instructions provided in Step 8 above. Write a "C" (for "control") and the date on the BI label. The positive control should be from the same lot code as the processed biological indicator. The positive control BI helps confirm:

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10. Incubation and Reading:

Incubate the positive control and processed 1295 BIs at $60 \pm 2^{\circ}\text{C}$ in a 490H Auto-reader or 490 Auto-reader having software version 4.0.0 or greater. See the Auto-reader Operator's Manual for the proper use of this equipment. The Auto-reader will indicate a positive result as soon as it is obtained. The final negative 1295 BI reading is made at:

- 24 minutes in 490H or 490 Auto-readers having software version 4.0.0 or greater
- 4 hours in 490H Auto-readers having a software version less than 4.0.0.

After the results are displayed and recorded, the 1295 BIs may be discarded.

Fluorescent Result

The positive control (unprocessed) 1295 BI must provide a positive fluorescent result (+ symbol on the 490H Auto-reader or 490 Auto-reader having software version 4.0.0 or greater LCD display). Processed 1295 BI results are not valid until the positive control reads fluorescent positive. If the positive control reads negative (- symbol on the LCD display), check the Auto-reader Operator's Manual Troubleshooting Guide. Retest the 490H Auto-reader or 490 Auto-reader having software version 4.0.0 or greater with a new positive control. With processed 1295 BIs, a final negative reading (- symbol on the LCD display) indicates an acceptable sterilization process. A positive (+ symbol on the LCD display) result for a processed 1295 BI indicates a sterilization process failure. Act immediately on any positive results for processed BIs. Determine the cause of the positive BI following facility policies and procedures. 3M recommends retesting the sterilizer according to your facility policy before processing additional loads.

Optional Visual pH Color Change Result

The 1295 BI is normally discarded after the fluorescent result has been recorded. If, however, special studies are desired, 1295 BIs may be further incubated for a visual pH color change result. In the case of the positive control BI, a yellow color change of the growth media will appear within 24 hours. Any observation of a yellow color within the vial indicates a positive result. In the case of a processed 1295 BI, a media color change from purple to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media remains purple, can be assessed at 7 days. To avoid media dry-out, it is recommended that the 1295 BI be transferred to a humidified incubator operating at 60°C after the fluorescent result has been recorded.

Storage/ Shelf Life

- Best stored under normal room conditions: $15\text{--}30^{\circ}\text{C}$ ($59\text{--}86^{\circ}\text{F}$), 35-60% relative humidity.
- Do not store these biological indicators near sterilants or other chemicals.
- Attest 1295 biological indicators have a 24 month shelf life. The expiration date is indicated on the BI and packaging by an hourglass symbol followed by the 4-digit year and 2-digit month of expiration.
- **Transportation:** During transportation, the products can be exposed to short-term temperature changes, which do not affect the safety and products specifications. Transportation of medical device is carried out by all type of transport.
- **Usage:** at room temperature $15\text{--}30^{\circ}\text{C}$. Values of humidity and pressure are not regulated.

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Technical characteristics

Glass ampoule: Diameter 6.4±0.2 mm, Height 39.6+1.0/-0.8 mm

Plastic sleeve: Diameter 12.0±0.1 mm, Height 60±0.5 mm

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**Biological Indicators for Sterilization Quality Control:
Attest Rapid Readout Biological Indicators for ethylene oxide sterilization**

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

Biological Indicators for Sterilization Quality Control:
Attest Rapid Readout Biological Indicators for ethylene oxide sterilization

Product Description

The 3M Attest 1294 Rapid Readout Biological Indicator (green cap) is specifically designed for the rapid, reliable monitoring of ethylene oxide (EO) sterilization processes when used in conjunction with the 3M™ Attest Auto-readers 290G or 390G. The 3M Attest rapid readout biological indicator detects the activity of a naturally occurring *Bacillus atrophaeus* enzyme, beta-glucosidase, which is one of the enzymes involved in spore outgrowth and normal vegetative cell function. The Auto-reader detects the activity of a *B. atrophaeus* enzyme by reading a fluorescent product produced by the enzymatic breakdown of a substrate in the media. Detection of the enzyme indicates a sterilization process failure. The final fluorescent readout of a negative indicator result is made after 4 hours of incubation when sterilization conditions at gas inject are $\geq 35\%$ RH. The 4-hour fluorescent readout has been validated and is routinely verified at standard ethylene oxide test conditions ($54^{\circ}\text{C} \pm 1^{\circ}\text{C}$, 600 ± 30 mg/l ethylene oxide, and $60\% \pm 10\%$ relative humidity).

Intended Use

Sterilization quality control

Scope of use: Sterilization

Contraindications

None.

Potential consumer

Only for professional use

Possible side effects

No known possible side effects.

Directions for Use

1. Identify the Attest rapid readout biological indicator by writing the sterilizer number, load number, and processing date on the indicator label. Do not place another label or indicator tape on the biological indicator vial or on the cap.

Place an Attest rapid readout biological indicator in an appropriate test package or tray according to recommended practices. Avoid placing the biological indicator in direct contact with a chemical indicator which could transfer fluorescent residue to the biological indicator. Appropriate test packages or trays for loads containing:

Wrapped telescopic instruments, plastic, rubber, metal instruments, and equipment

• Attest 1294 rapid readout biological indicator in an AAMI routine test pack*.

***AAMI Routine Test Pack**

• Place one Attest rapid readout biological indicator into the barrel of a plastic syringe with the cap of the biological indicator facing the syringe tip and the bottom facing the plunger.

• Remove the tip guard of the syringe.

• Place the syringe between the folds of one cotton surgical towel.

• Place a chemical indicator beside the syringe.

• Place the surgical towel in one peel pouch or wrapper.

Container systems

• Attest 1294 rapid readout biological indicator in the area(s) determined by product testing to be the most resistant.

3. Place the test package or tray in a full load in the most challenging area for the sterilant. This is generally in the center of the load.

4. Process the load as usual.

5. After sterilization, two options exist:

4. Aeration of Test Pack Prior to BI Incubation: Aerate the BI test pack for the same amount of time as the other packs from the load. At the end of the aeration cycle, remove the biological and chemical indicators from the test pack. Dispose of the remaining test pack components according to the healthcare facility's policy.

Removal of Biological Indicator from Test Pack Prior to Aeration: Open the sterilizer door according to the manufacturer's instructions, transfer the load to the aerator and retrieve the test pack. (If a sterilizer with aeration capability is used, follow the manufacturer's instructions for opening the door and remove the test pack prior to the aeration cycle). Remove the biological and chemical indicators from the test pack. Return the remainder of the test pack to the load for aeration according to the healthcare facility's policies.

6. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to the EO sterilization process. This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

7. While wearing safety glasses, press cap down. Crush the glass ampoule of the biological indicator in the crusher well of the auto-reader. Hold the biological indicator by the cap and tap on a hard surface until media wets strip at bottom of the vial, and then place the biological indicator in an auto-reader incubation well. See Attest Auto-reader Operator's Manual for further details.

8. Crush, tap and incubate at least one non-sterilized Attest rapid readout biological indicator (positive control) each day a sterilized indicator is incubated. Write a "C" and a date on the label. The positive control should be from the same lot number as the sterilized indicator. Incubating a positive control helps ensure:

- correct incubation temperatures are met;
- viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals;
- capability of media to promote rapid growth; and
- proper functioning of Auto-reader components.

9. Incubate the positive control and processed indicator for 4 hours at $37 \pm 2^\circ\text{C}$ ($99 \pm 3^\circ\text{F}$) and read the results in the 3M Attest 290G Auto-reader or 3M Attest Auto-reader 390G. See the Auto-reader Operator's Manual for the proper use of this equipment. The auto-readers automatically take readings and may indicate a positive result in less than 4 hours. The final fluorescent negative biological indicator reading is made at 4 hours. If sterilizer chamber relative humidity is 35%, record the final reading and discard the processed indicator. If sterilizer chamber relative humidity is unknown or $< 35\%$, continue to incubate the biological indicator for a pH color change result as low relative humidity may increase the fluorescent readout time beyond 4 hours. A humidified incubator will be required to avoid media dry out.

Interpretation of Results:

Fluorescent Results:

The positive biological indicator control must provide a positive result [red light on the 3M Attest 290G Auto-reader or plus symbol (+) on the LCD display of the 3M Attest Auto-reader 390G]. If the positive biological indicator control reads negative [green light on the 3M Attest 290G Auto-reader or minus symbol (-) on the LCD display of the 3M Attest Auto-reader 390G], refer to the applicable 3M Attest Auto-reader Operator's Manual Troubleshooting section. Retest the auto-reader with a new positive biological indicator control. The sterilized biological indicator results are not valid until the positive biological indicator control reads positive.

With a processed biological indicator, a positive [red light on the 3M Attest 290G Auto-reader or plus symbol (+) on the LCD display of the 3M Attest Auto-reader 390G] means a sterilization process failure has occurred. A negative [green light on the 3M Attest 290G Auto-reader or minus symbol (-) on the LCD display of the 3M Attest Auto-reader 390G] after 4 hours of incubation and chamber relative humidity 35% indicate an acceptable sterilization process.

pH Color Change Results

The appearance of a yellow color in the processed indicator demonstrates bacterial growth and a sterilization process failure. No color change indicates an adequate sterilization process. A final negative result is made after 7 days of incubation. The positive control indicator should show a color change from blue-green to yellow for the processed indicator results to be valid.

10. Immediately act on any positive biological indicator results. Always retest the sterilizer and do not use sterilizer for processing loads until the biological indicator results are negative.

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The 7-day results may be used to validate the 4-hour readings or in cases of unknown or low (<35%) chamber relative humidity. A humidified incubator will be required to avoid media dry out before 7 days. Indicators are examined daily for detection of a visual color change. After exposure to standard EO test conditions, the 4-hour fluorescent results are compared to the 7-day visual readings to determine the readout reliability of the indicator. Readout reliability of the Attest™ 1294 rapid readout biological indicator is determined using the sensitivity calculation in the following table.

Attest 1294 Rapid Readout Biological Indicators

Joslyn EO BIER Vessel, Oxyfume® 12 Gas

600 mg EO/l, 37°C, 60% RH and 600 mg EO/l, 54°C, 60% RH

Validation of Reduced Incubation Time (Readout Reliability Summary)

Sterilization Process	Incubation Temperature	Growth		Fluorescence	
		# Tested	# Growth Positives 168_hrs.	# False Negatives 4 hrs.	Sensitivity
37°C	37°C (99°F)	1110	752	0	100%
54°C	37°C (99°F)	1300	842	0	100%

$$\text{Sensitivity} = \frac{(\text{Number of Growth Positives after 168 hours}) - (\text{Number of False Negatives})}{\text{Number of Growth Positives after 168 Hours}} \times 100$$

The table shows that all the 7- day (i.e., 168 hours) growth visual positives were detected by fluorescence within 4 hours of incubation. This means no false negative indicators were detected after 4 hours of incubation. Based on the claimed 97% readout reliability of the 4-hour biological indicator results, there is no advantage to incubating the Attest 1294 rapid readout biological indicator beyond 4 hours when sterilizer chamber relative humidity is greater than or equal to 35%.

Storage and Shelf Life

- Best stored under normal room conditions: 15-30°C (59-86°F).
- Do not store these biological indicators near sterilants or other chemicals.
- Attest 1294 rapid readout biological indicators have a 2-year shelf life from the date of manufacture. The expiration date is indicated on the BI and packaging by the 4-digit year and 2-digit month of expiration (e.g. 2017-07).
- All of the information to the right of the lot-in-a-box symbol indicates the lot number (e.g., 2017-07AD).
- Transportation: During transportation, the products can be exposed to short-term temperature changes, which do not affect the safety and products specifications. Transportation of medical device is carried out by all type of transport.
- Usage: at room temperature 15-30°C. Values of humidity and pressure are not regulated.

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Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations. Medical products should be disposed by the local public law organizations in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for Medical Waste Treatment" waste Class B).

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Apostille

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1. Country (Pays): United States of America

This public document (Le présent acte public) Biological Indicators for Sterilization Quality Control

2. has been signed by (a été signé par) Ryan J. Bartz

3. acting in the capacity of (agissant en qualité de) Notary Public, State of Minnesota

**4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Ryan J. Bartz, Notary Public, State of Minnesota**

Certified Attesté

5. at (à) St. Paul, Minnesota

6. Date (Date) 11/15/2019

7. by (par) Secretary of State, State of Minnesota

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
Attest Super Rapid Readout Biological Indicators for rapid and reliable monitoring of the 270°F (132°C) and 275°F (135°C) dynamic-air-removal steam sterilization cycles

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

Attest Super Rapid Readout Biological Indicators for rapid and reliable monitoring of the 270°F (132°C) and 275°F (135°C) dynamic-air-removal steam sterilization cycles

Product Description

Attest Super Rapid Readout Biological Indicators for rapid and reliable monitoring of the 270°F (132°C) and 275°F (135°C) dynamic-air-removal steam sterilization cycles (The 3M Attest Super Rapid Readout Biological Indicator 1492V (brown cap, referred to hereinafter as the 1492V BI)) is a self-contained biological indicator specifically designed for rapid and reliable qualification testing and routine monitoring of 270°F (132°C) and 275°F (135°C) dynamic-air-removal (pre-vacuum and steam-flush pressure-pulse (SFPP)) steam sterilization processes when used in conjunction with a 3M Attest Auto-reader 490 (hereinafter referred to as the 490 Auto-reader) or a 3M Attest Auto-reader 490H having software version 4.0.0 or greater (hereinafter referred to as the 490H Auto-reader). The 1492V BI is a single-use device.

A schematic illustrating the design of the 1492V BI is provided in Figure 1.


Figure 1.

The self-contained design includes a carrier with spores of *Geobacillus stearothermophilus* and a media ampoule containing bacteriological growth medium which meets the requirements for growth promoting ability specified in ISO 11138-1:2017. The spore carrier and media ampoule are contained in a plastic vial topped with a brown cap. A chemical process indicator which changes from pink to light brown or darker upon exposure to steam is located on the top of the cap.

The 1492V BI utilizes the α -glucosidase enzyme system, which is generated naturally within growing cells of *Geobacillus stearothermophilus*. The α -glucosidase in its active state is detected by measuring the fluorescence produced by the enzymatic hydrolysis of a non-fluorescent substrate, 4-methylumbelliferyl- α -D-glucoside (MUG). The resultant fluorescent by-product, 4-methylumbelliferone (MU), is detected in the Auto-reader. The presence of fluorescence within the specified incubation time of the 1492V BI in the 490 or 490H Auto-reader indicates a steam sterilization process failure.

The 1492V BI can also indicate the presence of *G. stearothermophilus* organisms by a visual pH color change reaction. Biochemical activity of the *G. stearothermophilus* organism produces metabolic by-products that cause the media to change color from purple to yellow which also indicates a steam sterilization process failure. Use of this indication method is optional and is typically restricted to special

studies.

Readout Times

The super rapid readout and the optional 48-hour visual pH color change incubation times have been correlated with a 7-day incubation period following the FDA's Reduced Incubation Time protocol. Processed indicators were examined at 48 hours and 7 days for detection of a visual pH color change. The fluorescent result and the 48-hour visual pH color change readings were compared to the 7-day visual pH color change readings to determine the readout time of the indicator. The time to a fluorescent result is determined by the software version programmed on the Auto-reader.

24-minute Fluorescent Result

1492V BIs incubated in a 490 or 490H Auto-reader having software version 4.0.0 or greater have a 24 minute reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

1-hour Fluorescent Result

1492V BIs incubated in a 490 Auto-reader having a software version less than 4.0.0 have a 1 hour reduced incubation time result that correlates to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

48-hour Visual pH Color Change Result

1492V BIs incubated in either a 490 or 490H Auto-reader have 48 hour reduced incubation time results that correlate to the 7 day (168 hours) visual readout result $\geq 97\%$ of the time.

Due to the high reliability of the fluorescent result there is no advantage to incubating 1492V BIs after the fluorescent result has been determined by the 490 or 490H Auto-reader and recorded.

Note: To determine the software version on a 490 or 490H Auto-reader, depress the Information Button once. The Auto-reader will display two lines of text on the LCD display. The top row indicates (left to right):

- Serial number
- Software version (V X.Y.Z)
- Current date and time.

Intended Use

Sterilization quality control

Scope of use: Sterilization

Indications for Use

Use the 3M Attest Super Rapid Readout Biological Indicator 1492V in conjunction with the 3M Attest Auto-reader 490 or Attest Auto-reader 490H having software version 4.0.0 or greater to qualify or monitor 270°F (132°C) to 275°F (135°C) dynamic-air-removal steam sterilization cycles.

Contraindications

None.

Potential consumer

Only for professional use

Possible side effects

No known possible side effects.

Warnings

There is a glass ampoule inside the plastic vial of the biological indicator (BI). To avoid the risk of serious injury from flying debris due to a ruptured BI:

- Allow the BI to cool for the recommended time period before activating. Activating or excessive handling of the BI before cooling may cause the glass ampoule to burst.
- Wear safety glasses and gloves when removing the BI from the sterilizer.
- Wear safety glasses when activating the BI.
- Handle the BI by the cap when crushing or flicking.
- Do not use your fingers to crush the glass ampoule.

Precautions



1. DO NOT use the 1492V BI to monitor sterilization cycles which it is not designed to challenge:
 - a. Gravity-displacement steam sterilization cycles;
 - b. 250°F (121°C) dynamic-air-removal steam sterilization cycles;
 - c. Dry heat, chemical vapor, ethylene oxide or other low temperature sterilization processes.
2. To reduce the risk associated with incorrect results:
 - Do not place tape or labels on 1492V BI prior to sterilization or incubation.
 - Do not incubate a 1492V BI if, after processing and before BI activation, it is observed to have a broken media ampoule. Retest the sterilizer with a new biological indicator.
 - After processing, allow BI to cool for 10 minutes before incubation.
 - After BI activation, ensure media has flowed to the spore growth chamber.

Monitoring Frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, 3M recommends that every steam sterilization load be monitored with a biological indicator

Directions for Use

1. Identify the 1492V BI by writing the load number, sterilizer, and processing date on the indicator label. Do not place another label or indicator tape on the vial or on the cap.
2. Place the 1492V BI in a representative tray configuration or Process Challenge Device (PCD) as recommended by professional association guidelines or national standards for healthcare facility practice. Do not place the 1492V BI in direct contact with a chemical indicator as residue could transfer to the biological indicator and affect the result.
3. Place the PCD in the most challenging area of the sterilizer. This is typically on the bottom shelf, over the drain, however, the sterilizer manufacturer should be consulted to identify the area of the chamber least favorable to sterilization.
4. Process the load according to recommended practices.
5. After completion of the cycle take the PCD out of the sterilizer, and remove the 1492V BI.
6. Allow the 1492V BI to cool for 10 minutes prior to activation.
7. Check the process indicator on the top of the cap of the 1492V BI. A color change from pink to light brown or darker confirms that the 1492V BI has been exposed to the steam process. This color change does not indicate that the steam process was sufficient to achieve sterility. If the process indicator is unchanged, check the sterilizer physical monitors.
8. To activate the biological indicator, place it in an Attest Biological Indicator Activator. Close and squeeze the activator to close the 1492V BI cap and crush the media ampoule. Then remove the BI and flick it (see pictures). Visually verify that media has flowed into the growth chamber at the bottom of the vial. If the media hasn't filled the growth chamber, hold the BI by the cap and flick it until media fills the growth chamber. Place the activated 1492V BI:
 - in any well of a 490 or 490H Auto-reader having software version 4.0.0 or greater
 - or
 - in an incubation well of a 490 Auto-reader having software version less than 4.0.0 which is color-coded brown (i.e., configured to incubate 1492V BIs) and wait for the result. See the Auto-reader Operator's Manual for further information related to its use.

9. Each day that a processed 1492V BI is incubated, activate and incubate at least one non-processed 1492V BI to use as a positive control. Follow the activation instructions provided in Step 8 above. Write a "C" (for "control") and the date on the BI label. The positive control should be from the same lot code as the processed biological indicator. The positive control BI helps confirm:

- correct incubation temperatures are met,
- viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals,
- capability of media to promote rapid growth, and
- proper functioning of the Auto-reader.

10. Incubation and Reading:

Incubate the positive control and steam processed 1492V BIs in a 490 Auto-reader or a 490H Auto-reader having software version 4.0.0 or greater. See the Auto-reader Operator's Manual for the proper use of this equipment. The Auto-reader will indicate a positive result as soon as it is obtained. The final negative 1492V BI reading is made at:

- 24 minutes in 490 or 490H Auto-readers having software version 4.0.0 or greater
- 1 hour in 490 Auto-readers having software version less than 4.0.0.

After the results are displayed and recorded, the 1492V BIs may be discarded.

Interpretation of Results

Fluorescent Results

The positive control (unprocessed) 1492V BI must provide a positive fluorescent result (+ on the 490 or 490H Auto-reader LCD display). Processed 1492V BI results are not valid until the positive control reads fluorescent positive. If the positive control reads negative (- on the LCD display), check the Auto-reader Operator's Manual Troubleshooting Guide. Retest the Auto-reader with a new positive control.

With processed 1492V BIs, a positive (+ on the LCD display) result indicates a sterilization process failure. A final negative processed 1492V BI reading (- on the LCD display) after the specified incubation time indicates an acceptable sterilization process.

Act immediately on any positive results for processed BIs. Determine the cause of the positive BI following facility policies and procedures. Always retest the sterilizer and do not use sterilizer for processing loads until three consecutive BI results are negative.

Optional Visual pH Color Change Result

The 1492V BI is normally discarded after the fluorescent result has been recorded.

If, however, special studies are desired, 1492V BIs may be further incubated for a visual pH color change result. After activation and during incubation, the white Nonwoven Material will absorb the bromocresol purple indicator, the pH-sensitive indicator dye in the growth media, and appear blue. In the case of the positive control BI a yellow color change of the growth media and/or Nonwoven Material will appear within 48 hours. Any observation of a yellow color within the vial indicates a positive result.

In the case of a processed 1492V BI, a media and/or Nonwoven Material color change from purple to yellow indicates a sterilization process failure. A negative pH color change result, i.e., media and Nonwoven Material remain purple/blue, can be assessed at 48 hours.

Storage/ Shelf Life

- Best stored under normal room conditions: 15-30°C (59-86°F), 35-60% relative humidity.
- Do not store these biological indicators near sterilants or other chemicals.
- Attest 1492V biological indicators have a 24 month shelf life. The expiration date is indicated on the BI and packaging by an hourglass symbol followed by the 4-digit year and 2-digit month of expiration.
- **Transportation:** During transportation, the products can be exposed to short-term temperature changes, which do not affect the safety and products specifications. Transportation of medical device is carried out by all type of transport.
- **Usage:** at room temperature 15-30°C. Values of humidity and pressure are not regulated.

Disposal

Dispose of used biological indicators according to your health care facility policy. You may wish to steam sterilize any positive biological indicators at 132°C (270°F) for 4 minutes or at 275°F (135°C) for 3 minutes in a dynamic-air-removal steam sterilizer prior to disposal.

For processed BIs observed to have a broken media ampoule after the sterilization process, or if the ink of the CI appears smeared or runny (if included) leave the BI and CI in the sterilization pouch and place the pouch in packaging compatible with steam sterilization and steam sterilize per the parameters above. Packaging/contents should be disposed in accordance with the local/regional/national/ international regulations.

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Technical characteristics

Glass ampoule: Diameter 6.4±0.2 mm, Height 39.6±1.0/-0.8 mm

Plastic sleeve: Diameter 12.0 ± 0.1 mm, Height 60 ± 0.5 mm

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