15

ANGIODYNAMICS

NanoKnife

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Manufacturer of Electrode probes for use with the NanoKnife System:

AngioDynamics, Inc.
603 Queensbury Ave.
Queensbury, NY 12804
USA

Use Conditions:

The probe is intended for use with an AngioDynamics NanoKnife generator only. The physician must read the NanoKnife generator's User Manual thoroughly before operating the unit.

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Classifications:

Council Directive 93/42/EEC Concerning Medical Devices
Class IIb

Physical Specifications (Without Packaging)

Dimensions (Length):

The Single Electrode Probe - 15 cm.

The Single Electrode Probe - 25 cm.

The Single Electrode Activation Probe - 15 cm.

The Single Electrode Activation Probe - 25 cm.

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Indications for Use

Electrode probe for use with the NanoKnife System is a medical device for electrosurgical manipulations in General practical surgery practice.

Cleaning and disinfection

This device is sterilized by ethylene oxide and intended for single patient use.

Do not use if package is received open or damaged. Do not re-sterilize. Reprocessing may compromise the integrity of the device and / or lead to device failure.

Utilization

The utilization of the Electrode probes for use with the NanoKnife System must be in accordance with all applicable federal, national and local laws and regulations.

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NanoKnife

Instructions for use.

Electrode probes for use with the NanoKnife System

angiodynamics

603 Queensbury Avenue
Queensbury
NY 12804, USA

Stephen A. Trowbridge
Stephen A. Trowbridge
Sr. VP and General Counsel
AngioDynamics, Inc.
06 DEC 2013

AngioDynamics and NanoKnife are registered trademarks of AngioDynamics, Inc. U.S. and Foreign Patents Pending.

© AngioDynamics, Inc. 2010 all rights reserved

**CAUTION**

Federal (USA) law restricts this device to sale by or use under the order of a physician.

This device is sterilized by ethylene oxide and intended for single patient use. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Do not use if package is received open or damaged. Do not re-sterilize. Reprocessing may compromise the integrity of the device and / or lead to device failure.

The probe is intended for use with an AngioDynamics NanoKnife generator only. The physician must read the generator's User Manual thoroughly before operating the unit.

Avoid short circuiting the electrodes when delivering pulses.

INTENDED USE

The Single Electrode Probe and Single Electrode Activation Probe are intended for use with the NanoKnife System.

POTENTIAL ADVERSE EFFECTS

Adverse effects that may be associated with the use of the NanoKnife System include, but are not limited to:

- Arrhythmia
- Hemorrhage
- Infection
- Damage to critical anatomical structure (nerve, vessel, duct)*
- Pneumothorax
- Unintended mechanical perforation
- Bradycardia
- Vagal stimulation, asystole
- Muscle contraction

*Adverse effect may occur but has not been reported

WARNINGS

Do NOT use a device with damaged insulation.

Do NOT attach anything to the device unless it is supplied by AngioDynamics and indicated for use with this device. Attachments may damage the insulation and contribute to patient injury.

MR Unsafe: This device should not be used in conjunction with MRI image guidance as the device has not been tested for MRI compatibility.

PRODUCT DESCRIPTION

The Single Electrode Probe and Single Electrode Activation Probe are probes with 19 gauge trocar needles available in 15 cm and 25 cm lengths. Including the insulation, diameter is 18 gauge.


Figure 1: Probe configurations

The probes are comprised of:

1. Active electrode, length adjustable in 0.5 cm increments from 0 – 4 cm via the thumb slide.
2. Insulation sleeve.
3. 19 gauge needle with depth markers and echogenic needle tip.
4. 10 foot connection cable.
5. Sterile labels for identifying the generator port number the probe will be connected to (supplied only with the Single Electrode Activation Probe).


Figure 2: Probe components

The Single Electrode Activation Probe is blue, while the standard Single Electrode Probe is white. Both are available in 15 cm and 25 cm lengths. The blue Single Electrode Activation Probe is required to activate the generator which will allow other standard Single Electrode Probes to function with it. Note: only one Single Electrode Activation Probe is required to activate the generator but a minimum of two probes (white or blue) are required for an ablation to be executed. The Single Electrode Activation Probe may be attached to any output port of the generator. Depending on the size of the soft tissue area to be ablated, a maximum of six probes can be used at a time. Probes may be re-positioned after each ablation to cover a larger area, as guided by the generator software.

INSTRUCTIONS FOR USE

1. Refer to the generator's user manual and become familiar with operation of the generator.
2. Turn on the generator.
3. Inspect all of the devices and packaging for damage prior to use. Do NOT use any devices that are damaged or if the sterile barrier is breached.
4. Using sterile technique, open the package and carefully remove the device. Remove and dispose of the protective shipping sheath covering the needle.
5. Sterile labels are included with the Single Electrode Activation Probe that may be used to identify the connector and corresponding handles that will be attached to each port of the generator during the procedure.

Figure 3: Labels for connector and handle

6. Attach the cable connectors to the front panel of the generator. One of the Single Electrode Activation Probes will be required to activate the generator. Note: If excessive force is required to connect the cable to the generator, check to see if the pin in the cable connector is bent.
7. Program the generator with the parameters required to achieve the desired ablation zone. Refer to the generator manual for specific programming instructions.
 - a. Typical results (based on in vitro porcine data) for two probes spaced 1 cm apart with the active electrode length set at 2 cm.

Generator Settings:

Probe +	Probe -	Voltage	Pulse Length	N. Pulses
1	2	2000 V	100 μ s	90


Figure 4: Typical ablation results based on in vitro porcine data.

- b. Typical results (based on in vitro porcine data) for two probes spaced 1.5 cm apart with the active electrode length set at 2 cm.

Generator Settings:

Probe +	Probe -	Voltage	Pulse Length	N. Pulses
1	2	2500 V	100 μ s	90


Figure 5: Typical ablation results based on in vitro porcine data.

8. Adjust the exposure of the electrode using the thumbslide, pressing down on the proximal raised end to unlock. Exposure may be adjusted from 0 to 4.0 cm in 0.5 cm increments.


Figure 6: Electrode exposure adjustment.

9. Using imaging guidance, the probes should be placed in the targeted tissue area as parallel as possible. Markings on the needle can assist in placement of the device. The markings extend up the needle in 1.0 cm increments.


Figure 7: Depth markings on needle.

10. Verify the position of the probes with imaging.
Note: Needles should be verified to be as close to parallel as possible with imaging.
11. Verify the length of the exposed electrode indicated by the thumbslide on the handle (see Figure 6).
12. Once all desired probes are placed, deliver the programmed pulse sequence. Refer to generator manual for pulse delivery instructions.
13. After the ablation is complete, review the graphs on the Pulse Generation screen to ensure that the pulses were delivered.

SYMBOL DEFINITIONS

Do Not Reuse		Keep Away From Sunlight	
Sterilized using Ethylene Oxide		Do Not Use If Package is Damaged	
Do Not Restertize		Catalog Number	
Caution		Batch Code	
Keep Dry		Use By	
Upper Temperature Limit			

Stephen A. Trowbridge
Stephen A. Trowbridge
Sr. VP and General Counsel
AngioDynamics, Inc.
c/o DEC 2013

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