
"APPROVE"
ETHICON, LLC
Highway 183 Km 8.3,
San Lorenzo, Puerto Rico 00754, USA


Stephanie Saati
Manager, RA Global Strategy Execution
Stephanie Saati

(signature)

« 09 » Oct 2020 г.

Subscribed and sworn to before me
this 09th day of October 2020

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INSTRUCTION FOR USE

TEMPORARY CARDIAC PACING WIRE IN
CONFIGURATIONS: TPW10, TPW30, TPW92

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INSTRUCTION FOR USE

TEMPORARY CARDIAC PACING WIRE IN CONFIGURATIONS: TPW10, TPW30, TPW92


1. Temporary Cardiac Pacing Wire TPW10: TPW10 Wire – 12 pcs, IFU;
2. Temporary Cardiac Pacing Wire TPW30: TPW30 Wire – 12 pcs, IFU;
3. Temporary Cardiac Pacing Wire TPW92: TPW92 Wire – 24 pcs, IFU.

(Hereinafter may be referred to as - Temporary Cardiac Pacing Wire, TPW, the device)

DESCRIPTION

This Temporary Cardiac Pacing Wire is a sterile single-use product. Temporary cardiac pacing wires consist of the following configurations and optional accessories:

1. The monopolar leads consist of one insulated multifilament stainless steel conductor coated with colored polyethylene.
2. The intracorporeal end of the wire has a section of exposed, uninsulated wire electrode which terminates with an attached stainless steel needle.
3. The exposed, uninsulated section of the intracorporeal end of the wire is either straight or has one or more multiple pre-formed curves (pre-formed wave).
4. The extracorporeal end of wire has a straight needle with breakaway tip attached.
5. The wire length is 60 cm.
6. Wires range in multiple diameters from 0 to 2-0, depending on the product code.

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The lead consists of:

- Item 1 – Insulated multifilament conductor coated with colored polyethylene
- Item 2 – Uninsulated, multifilament wire (intracorporeal) electrode
- Item 3 – Optional wave pre-formed into the uninsulated, intracorporeal wire end
- Item 4 – Intracorporeal curved needle

Item 5 – Extracorporeal straight breakaway needle

PRINCIPLES OF OPERATION/PERFORMANCE

The Temporary Cardiac Pacing Wire provides a conductive connection between an external pacemaker and the myocardium.

Electrical parameters for Temporary Cardiac Pacing Wires:

- Lead conductor resistance ($9.1 \pm 2.1 \Omega$)

INDICATIONS

Temporary Cardiac Pacing Wire is intended for use in temporary epicardial cardiac pacing or monitoring and should be removed after temporary pacing has been discontinued.

CONTRAINDICATIONS

Use of temporary cardiac pacing wires is contraindicated for permanent cardiac pacing or for monitoring.

WARNINGS AND PRECAUTIONS

Before using the Temporary Cardiac Pacing Wire, it should first be established that this device is compatible with all products from other manufacturers that may be required for use. This particularly applies when checking the safety of the connection between pacing wire and pacemaker or pacemaker cable.

Care should be taken regarding the following:

- Avoid disruption of the polyethylene coating or damage to the conductive wire, i.e., by twisting, manipulation with instruments, or by tying knots.

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- Avoid any possible unintentional contact between the Temporary Cardiac Pacing Wire and any equipment.
- Particular care should be taken to ensure that no components which might touch each other and conduct electricity are present at the junction (e.g., projecting needle tips).
- If a connection cable (not supplied) is being used between the pacing wire and the pacemaker and temporary disconnection is necessary, be careful to disconnect the connecting cable in such a way as to prevent damage to the pacing wire itself or accidental dislodgment of the pacing wire from the myocardium.
- Ensure that all the electrical junctions between the pulse generator and temporary cardiac pacing wire(s) are intact to provide optimum performance, following relevant instructions for use where applicable.
- An implanted temporary cardiac pacing wire provides a low-resistance electrical connection to the myocardium. Particular care is therefore required when other devices liable to generate electrical energy are used concurrently (e.g., defibrillators, electrosurgical devices, and other devices that generate electromagnetic frequency). Follow the relevant instructions for use.
- The pacing wires must be placed in the myocardium and trimmed in such a way that the risk of injury when placing and extracting the pacing wires is minimized.
- Reuse of this device (or portions of this device) may create a risk of product degradation and cross-contamination, which may lead to infection or transmission of bloodborne pathogens to patients and users.

ADVERSE REACTIONS

None Known.

POTENTIAL COMPLICATIONS OF TEMPORARY CARDIAC ELECTRICAL PACING

Potential complications of temporary cardiac electrical pacing include:

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- Avoid any possible unintentional contact between the Temporary Cardiac Pacing Wire and any equipment.

- Particular care should be taken to ensure that no components which might touch each other and conduct electricity are present at the junction (e.g., projecting needle tips).

- If a connection cable (not supplied) is being used between the pacing wire and the pacemaker and temporary disconnection is necessary, be careful to disconnect the connecting cable in such a way as to prevent damage to the pacing wire itself or accidental dislodgment of the pacing wire from the myocardium.

- Ensure that all the electrical junctions between the pulse generator and temporary cardiac pacing wire(s) are intact to provide optimum performance, following relevant instructions for use where applicable.

- An implanted temporary cardiac pacing wire provides a low-resistance electrical connection to the myocardium. Particular care is therefore required when other devices liable to generate electrical energy are used concurrently (e.g., defibrillators, electrosurgical devices, and other devices that generate electromagnetic frequency). Follow the relevant instructions for use.

- The pacing wires must be placed in the myocardium and trimmed in such a way that the risk of injury when placing and extracting the pacing wires is minimized.

- Reuse of this device (or portions of this device) may create a risk of product degradation and cross-contamination, which may lead to infection or transmission of bloodborne pathogens to patients and users.

ADVERSE REACTIONS

None Known.

POTENTIAL COMPLICATIONS OF TEMPORARY CARDIAC ELECTRICAL PACING

Potential complications of temporary cardiac electrical pacing include:

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- Impairment or loss of sensing due to lead dislodgment, or fracture, changes in the myocardium, or failure of the pulse generator.
- Increased pacing threshold or loss of pacing due to lead dislodgment or fracture, changes in the myocardium, or failure of the pulse generator.
- Arrhythmias due to myocardial irritability.
- Infections, which can be systemic, local, or myocardial.
- Bleeding or myocardial damage.

STERILITY

The product is sterilized by irradiation. Do not use if the package is damaged or opened. Discard opened, unused product. Do not resterilize.

USAGE CONDITIONS

The medical device is intended for use in the standard hospital operating room environment.

Medical Device has the implantable part. The implantable part is intended to be used in the body's internal environment.

MRI

MR Conditional

This Temporary Cardiac Pacing Wire is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 1,400 Gauss/cm (extrapolated).
- Transmit/receive RF head coil, only.
- Maximum MR system reported, transmit/receive RF head coil specific absorption rate (SAR) of 3.2 W/kg for 15 minutes of scanning (i.e. per pulse sequence).
- Normal Operating Mode of operation for the MR system.

In non-clinical testing and numerical modeling, the Temporary Cardiac Pacing Wire produced a temperature rise of less than 0.5°C at a maximum MR system calculated transmit/receive RF head coil SAR of 3.2-W/kg in 1.5-Tesla/64-MHz and 3-Tesla/128-MHz MR systems.

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Additional MRI Safety Information

- Do not perform an MRI procedure if the Temporary Cardiac Pacing Wire is damaged or otherwise not functioning properly.
- The Temporary Cardiac Pacing wire must be disconnected from the pulse generator prior to entry into the MR system room.
- The use of a transmit/receive RF body coil to perform an MRI examination in a patient with the Temporary Cardiac Pacing Wire may cause patient injury due to excessive MRI-induced heating.

CONTENTS

The Temporary Pacing Wire with attached needles is packaged in a paper folder that is sealed into a Tyvek sterile barrier pouch. Twelve pouched TPWs are then placed into a sales carton tray with an IFU. The tray is then placed into a sleeve and then wrapped with a dust cover.

APPLICATION

Removing Pacing Wires from Package

Temporary Cardiac Pacing Wire is offered in various horizontal and vertical configurations with 12 or more pouches per box. After removing the plastic outer wrap, open the box and select one Temporary Cardiac Pacing Wire(s) pouch. Over the sterile field, carefully peel back the clear polymeric portion of the pouch and withdraw the paper folder. Open the paper folder by lifting the tab(s). Then remove the device(s) to the sterile field. Avoid disruption of the polyethylene coating or damage to the conductive wire, i.e., by twisting, or by manipulation with instruments.

Intracorporeal Placement (Applicable to All Pacing Wires)

The intracorporeal end of the Temporary Cardiac Pacing Wire is attached to the myocardium as indicated by surgical circumstances, clinical judgment, and the preference of the operator. The temporary cardiac pacing wires should be placed in a position which allows uncomplicated removal and minimizes the risk of trauma when treatment is complete. The following is a suggested method for placement of the intracorporeal end(s) of the Temporary Cardiac Pacing Wire:

For all product types (non-color coded and color coded), surgeon clinical judgment determines the location and technique employed for wire placement. When applicable,

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the light and dark blue wire coatings may be helpful, after implantation, in differentiating the positive and negative wires or in differentiating atrial and ventricular wire placement.

VENTRICULAR/ATRIAL USE. Utilize the curved needle on the intracorporeal end of the Temporary Cardiac Pacing Wire to place the wire into the myocardium. Care should be taken in the placement of the wire to minimize the risk of unintended trauma to tissue upon removal (see Figure 1 for ventricular example). If a needle is not utilized on the intracorporeal end of the pacing wire or to minimize trauma to the tissue and facilitate easy removal, the pacing wire can be brought into contact with the myocardium using a small loop of suture passed through the surface of the myocardium. Use scissors to cut excess exposed wire and the attached needle when present (see Fig. 2 for ventricular example).

OPTIONAL ATRIAL USE. If the surgeon chooses, the following technique can be used to minimize trauma to the atrial tissue, provide close contact with the atrial tissue and facilitate easy removal: place the uninsulated wire in a fold created on the epicardial surface of the atrial tissue using a suture loop. Use scissors to cut excess exposed wire and the attached needle when present (see Fig. 2).

PRE-FORMED WAVE USE. For wires that include the optional pre-formed wave in the exposed wire of the intracorporeal end, the wave should be gently pulled into and left embedded in the myocardium. Use scissors to cut excess exposed wire and the attached needle when present (see Fig. 2).

Extracorporeal Placement (Applicable to All Pacing Wires)

The extracorporeal end(s) of the Temporary Cardiac Pacing Wire is managed as indicated by surgical circumstances, clinical judgment, and the preference of the operator. The following is a suggested method for placement of the extracorporeal end(s) of the Temporary Cardiac Pacing Wire:

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The extracorporeal end of the Temporary Cardiac Pacing Wire is brought through the chest wall using the attached needle.

For products with **Straight Breakaway extracorporeal needles**, there is no exposed wire at the point of swage and the straight needle is scored, eliminating the need to cut the needle off.

The needle is broken off at the score mark (see Fig. 3) to create a connector pin (stainless steel electrode).

The connector pin is attached to the connecting cable of the pacemaker or, as determined by clinical need, directly to the temporary pacing unit or monitoring unit per the manufacturer's instructions (see Fig. 4).

Note: If a pacemaker will be attached later, clinical judgment should be exercised to determine the optimal means of isolating the electrodes in such a way that no exposed metal remains.

IMPORTANT: Before using the cardiac pacing wire, it is necessary to check that the Temporary Cardiac Pacing Wire electrode and the pacemaker cable or pacemaker device are fully compatible. The electrode, once the needle is broken-off, is 20 - 38 mm length by 0.7 - 1.0 mm diameter, depending on the product code. Particular care should be given to ensuring that no components which could conduct electricity might touch each other at the connection point (e.g., projecting needle tips).

Removal of Temporary Cardiac Pacing Wire

The Temporary Cardiac Pacing Wire is removed as indicated by surgical circumstances, clinical judgment, and the preference of the operator. The following is a suggested method for removal of the intracorporeal end(s) of the Temporary Cardiac Pacing Wire:

When temporary pacing has been discontinued, the lead can be removed with gentle traction.

Care should be taken to prevent traction on the lead that may cause inadvertent premature removal.

TRANSPORTATION

Transportation conditions: Do not exceed -60°C or 60°C

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STORAGE

No special storage conditions required. Do not use after the expire date. Recommended storage conditions: Do not exceed -60 °C or 30 °C.

EXPIRATION DATE

Shelf life of the product is 5 years.
Warranty period of storage is 5 years.

WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements.

There are no warranties that extend beyond the description on the face hereof. Suitability of use this medical device for any particular surgical procedures should be determined by the user in conformance with the 'manufacturer's instructions for use'

DISPOSAL

Discard open packaging and unused devices.

Discard used needles in "sharps" containers.

Dispose the product and its residual elements or waste items in accordance with local laws and regulations.

After being used the products are potentially infectious waste (class B medical devices under medical wastes classification).

SYMBOLS USED ON LABELING

	See instructions for use
	Do not reuse
	MR Conditional
	Number of units

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LOT	Batch number
STERILE R	Sterile unless package is damaged or opened. Method of sterilization: Irradiation
	Manufacturer
	Use by – year and month
EC REP	Authorized Representative in the European Community
CE	CE-Mark and Identification Number of Notified Body. The product meets the essential requirements of Medical Device Directive 93/42/EEC.
REF	Catalogue number
	Do not use if packaging is damaged
	Temperature range
	See instructions for use
	Non-pyrogenic
PG	Symbol of conformity

 Ethicon, LLC, Highway 183 Km 8.3, San Lorenzo, Puerto Rico 00754, USA

AUTHORIZED REPRESENTATIVE OF MANUFACTURER

"Johnson & Johnson" LLC, 121614, Russia, Moscow, st. Krylatskaya, building 17/2,

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SPECIFICATIONS

Product Code	Description	Parts of electrode	Type
TPW10	1 pc of TPW consisting of 2 atraumatic needles (SS 420) on each end of the wire (SS 316L) with blue PE insulation	Wire	Monopolar, 2-0 (3 Ph. Eur.) / 2-0 (3.0 metric), length 60 cm
		Curved Needle	RB-1: (length 17mm, 1/2c, Round bodied/ Taperpoint)
		Straight Needle	SKS (length 60mm, Cutting/ Cutting Edge Reverse, breakaway)
TPW30	1 pc of TPW consisting of 2 atraumatic needles (SS 420) on each end of the wire (SS 316L) with blue PE insulation	Wire	Monopolar, 2-0 (3 Ph. Eur.) / 2-0 (3.0 metric), length 60 cm
		Curved Needle	SH: (length 26mm, 1/2c, Round bodied/ Taperpoint)
		Straight Needle	SKS-3: (length 90 mm, Cutting / Cutting Edge Reverse, breakaway)
TPW92:	2 pcs of TPW. Each TPW consist of: wire (SS 316L) with two atraumatic needles (SS 420). One wire has dark blue PE insulation. Another wire has light blue PE insulation	Wire	Monopolar, 2-0 (3 Ph. Eur.) / 2-0 (3.0 metric), length 60 cm
		Curved Needles	SH: (length 26mm, 1/2c, Round bodied/ Taperpoint)
		Straight Needles	SKS-3: (length 90 mm, Cutting/ Cutting Edge Reverse, breakaway)

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TPW10 characteristics	Range
TPW quantity in the primary package	1
TPW quantity in the secondary package/Delivery set	12
Wire Gauge Size	2-0 (3 Ph. Eur.) / 2-0 (3.0 metric)
Wire nominal length	60 cm
Total Length	63...66 cm (24.75 – 26 in)
Coated Diameter	0.559...0.610 mm (22 – 24 mils)
Uncoated Diameter	0.310...0.338 mm (12.2 – 13.3 mils) (2-0 U.S.P.)
Lead conductor resistance	(9.1 ± 2.1) Ω
Wire material	316L SS
Wire insulation	Polyethylene
Wire insulation color	Blue
Needle material	420 SS
Needle coating	Silicone
Number of Atraumatic Needles per Wire	Double Armed
RB-1 Needle Shape	Curve
Radius of RB-1 curvature	5.41±0.13 mm
SKS Needle Shape	Straight
RB-1 Needle Length	17 mm*
SKS Needle Length	60 mm*
RB-1 Needle Diameter	0,6858 mm (27 mils)*
SKS Needle Diameter	1,016 mm (40 mils)*
RB-1 Diameter in the binding zone	No more than 1.15 of needle diameter in the beginning of the binding zone
SKS Diameter in the binding zone	No more than 1.15 of needle diameter in the beginning of the binding zone
RB-1 Needle Pull-Off	Minimum Average 10.7 N Minimum Individual 4.45 N
SKS Needle Pull-Off	Minimum Average 10.7 N Minimum Individual 4.45 N
RB-1 and SKS Needle hardness due to Vickers*	HV 4900-6475 N/mm ²
Roughness parameters of RB-1 and SKS needles surface, no more than	on the shank and point – Ra 0.32 μm; in the binding zone – Rz 20 μm

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Sharpness of the stitching part of the RB-1 and SKS needles, no more than	0,025 mm
Breakaway of SKS Needle	Straight needles intended for direct insertion into external pacemaker must break at the score marks
Connections of electrodes with Cardiac pacemakers	Non-specific mechanical tightening connector
Leakage current	≤ 2.0 mA
TPW Weight	(0.78±0.025) g
Lead pacing impedance	From 21,34 to 28,34 Ω
Lead sensing impedance	From 80,26 to 129,31 Ω

TPW 30 characteristics	Range
TPW quantity in the primary package	1
TPW quantity in the secondary package/Delivery set	12
Wire Gauge Size	2-0 (3 Ph. Eur.) / 2-0 (3.0 metric)
Wire nominal length	60 cm
Total Length	63...66 cm (24.75 – 26 in)
Coated Diameter	0.559...0.610 mm (22 – 24 mils)
Uncoated Diameter	0.310...0.338 mm (12.2 – 13.3 mils) (2-0 U.S.P.)
Lead conductor resistance	(9.1 ± 2.1) Ω
Wire material	316L SS
Wire insulation	Polyethylene
Wire insulation color	Blue
Needle material	420 SS
Needle coating	Silicone
Number of Atraumatic Needles per Wire	Double Armed
SH Needle Shape	Curve
Curve of SH Needle	8,76±0.13 mm
SKS-3 Needle Shape	Straight
SH Needle Length	26 mm*

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SKS-3 Needle Length	90 mm*
SH Needle Diameter	0,6858 mm (27 mils)*
SKS-3 Needle Diameter	1,016 mm (40 mils)*
SH Diameter in the binding zone	No more than 1,15 of needle diameter in the beginning of the binding zone
SKS-3 Diameter in the binding zone	No more than 1,15 of needle diameter in the beginning of the binding zone
SH Needle Pull-Off	Minimum Average 10.85 N (2.44lbs) Minimum Individual 4.45 N (1.00lbs)
SKS-3 Needle Pull-Off	Minimum Average 10.7 N Minimum Individual 4.45 N
SH and SKS-3 Needle hardness due to Vickers	HV 4900-6475 N/mm²
Roughness parameters of SH and SKS-3 needles surface, no more than	on the shank and point – Ra 0.32 µm; in the binding zone – Rz 20 µm
Sharpness of the stitching part of the SH and SKS-3 needles, no more than	0,025 mm
Breakaway of SKS-3 Needle	Straight needles intended for direct insertion into external pacemaker must break at the score marks
Connections of electrodes with Cardiac pacemakers	Non-specific mechanical tightening connector
Leakage current	≤ 2.0 mA
TPW Weight	(1.035±0.025) g
Lead pacing impedance	From 21,34 to 28,34 Ω
Lead sensing impedance	From 80,26 to 129,31 Ω

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TPW 92 characteristics	Range
TPW quantity in the primary package	2
TPW quantity in the secondary package/Delivery set	24
Wire Gauge Size	2-0 (3 Ph. Eur.) / 2-0 (3.0 metric)
Wire nominal length	60 cm
Total Length	63...66 cm (24.75 – 26 in)
Coated Diameter	0.559...0.610 mm (22 – 24 mils)
Uncoated Diameter	0.310...0.338 mm (12.2 – 13.3 mils) (2-0 U.S.P.)
Lead conductor resistance	$(9.1 \pm 2.1) \Omega$
Wire material	316L SS
Wire insulation	Polyethylene
Wire insulation color	Dark blue and Light Blue
Needle material	420 SS
Needle coating	Silicone
Number of Atraumatic Needles per Wire	Double Armed
SH Needle Shape	Curve
Curve of SH Needle	8.76±0.13 mm
SKS-3 Needle Shape	Straight
SH Needle Length	26 mm*
SKS-3 Needle Length	90 mm*
SH Needle Diameter	0,6858 mm (27 mils)
SKS-3 Needle Diameter	1,016 mm (40 mils)
SH Diameter in the binding zone*	no more than 1,15 of needle diameter in the beginning of the binding zone
SKS-3 Diameter in the binding zone*	no more than 1,15 of needle diameter in the beginning of the binding zone
SH Needle Pull-Off	Minimum Average 10.85 N (2.44lbs) Minimum Individual 4.45 N (1.00lbs)
SKS-3 Needle Pull-Off	Minimum Average 10.7 N Minimum Individual 4.45 N
SH and SKS-3 Needle hardness due to Vickers	HV 4900-6475 N/mm ²
Roughness parameters of SH and SKS-3 needles surface, no more than	on the shank and point – Ra 0.32 μm ; in the binding zone – Rz 20 μm

Sharpness of the stitching part of the SH and SKS-3 needles, no more than	0,025 mm
Breakaway of SKS-3 Needle	Straight needles intended for direct insertion into external pacemaker must break at the score marks
Connections of electrodes with Cardiac pacemakers	Non-specific mechanical tightening connector
Leakage current	≤ 2.0 mA
TPW Weight	(1.035 $\pm$ 0.025) g
Pre-formed wave length	13.33...13.97 mm (0.525"...0.550")
The distance from the needle to the pre-formed wave	≥ 23.1 mm (0.91")
Pre-formed wave height	9.9...13.5 mm (0.39"...0.53")
Lead pacing impedance	From 21,34 to 28,34 Ω
Lead sensing impedance	From 80,26 to 129,31 Ω

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