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Instruction for use

Nasal drops sterile MARIMER

«APPROVED»

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Position
Responsible Pharmacist/Quality Director

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*I hereby certify accuracy, correctness and
reliability of this document
translated into Russian*

CONTENTS

NAME OF MEDICAL DEVICE	3
INFORMATION ABOUT MANUFACTURER.....	3
DEVELOPER / MANUFACTURER.....	3
MANUFACTURING SITE OF MEDICAL DEVICE.....	3
INTENDED USE OF MEDICAL DEVICE	3
INDICATIONS FOR USE.....	3
CONTRAINDICATIONS.....	3
SIDE EFFECTS	3
CONDITIONS OF USE	4
DESCRIPTION, COMPOSITION AND SUPPLY PACKAGE OF MEDICAL DEVICE.....	4
POSODOLOGY AND METHOD OF ADMINISTRATION.....	4
TECHNICAL DESCRIPTION OF MEDICAL DEVICE.....	5
PRECAUTIONS FOR USE	5
STERILIZATION METHODS	5
OPERATIONAL, STORAGE AND TRANSPORTATION CONDITIONS.....	5
SHELF-LIFE	5
ENVIRONMENTAL REQUIREMENTS	5
LABELLING.....	6
PROCEDURE FOR DISPOSAL OR DESTRUCTION OF MEDICAL DEVICE.....	6
LIST OF INTERNATIONAL STANDARDS.....	6
LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN	7
INFORMATION ABOUT DRUG PRODUCTS, MEDICINAL PRODUCTS AND PHARMACEUTICAL SUBSTANCES	7
WARRANTY	7
RECLAMATION.....	7

NAME OF MEDICAL DEVICE

Nasal drops sterile MARIMER, composed of:

- Nasal drops in a single-dose 5 mL dropper vial – 10 pcs.
- Package leaflet – 1 pce.

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Tel./Fax +33 2 31 47 15 15/+33 2 31 43 95 00

MANUFACTURING SITE OF MEDICAL DEVICE

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INTENDED USE OF MEDICAL DEVICE

Nasal drops sterile MARIMER is intended for irrigation, rinsing of nasal passages and sinuses, for the complex treatment and prophylaxis of ENT diseases, viral infections (ARVI, flu), allergic and vasomotor rhinitis

INDICATIONS FOR USE

Nasal drops sterile MARIMER are used for the newborns, children and adults:

- for daily cleaning of the nasal cavities, including for the maintenance of the protective properties of the mucous membrane under conditions of increased dryness or air pollution (air conditioning, central heating, etc.);
- in case of nasal congestion to encourage elimination of mucus;
- for the complex treatment and prophylaxis of ENT diseases, viral infections (ARVI, influenza), allergic and vasomotor rhinitis.

CONTRAINDICATIONS

Device has no documented contraindications.

SIDE EFFECTS

Side effects when using the device in accordance with the instructions for use have not been identified.

CONDITIONS OF USE

Used in medical facilities, as well as at home. Refers to personal use devices. Carrying out procedures with this device at home does not require special training and special skills.

DESCRIPTION, COMPOSITION AND SUPPLY PACKAGE OF MEDICAL DEVICE

MARIMER nasal drops sterile is a sterile isotonic seawater solution in a single-use dropper vial. The solution is administered self-sustained by compressing the dropper vial. The trace elements in the composition of sterile isotonic seawater solution improve the function of the ciliary epithelium, which increases the resistance of the nasal mucosa to pathogens and viruses.

MARIMER nasal drops sterile are packaged in a 5 ml single-use dropper vials made of low-density polyethylene and then placed in amount of 10 pcs. in a tamper-evident carton along with a package leaflet.

Composition per 100 mL: seawater – 28 mL, purified water – up to 100 mL.
pH 6,5-8,5
Osmolality 270 - 310 mOsm/kg

POSODOLOGY AND METHOD OF ADMINISTRATION

For daily cleansing of the nasal cavities: rinsing of each nasal passage 1 to 4 times a day.

For complex treatment, prophylaxis of diseases and nasal congestion: rinsing of each nasal passage 4 to 6 times a day before additional local treatment.

Single-use dropper vial provides easy and safe use of the product.

Single-use dropper vial is opened by twisting off the nozzle.

For instillation and rinsing:

- Primarily, each nasal passage must be cleaned by blowing. In infants, it is possible to use an aspirator, for example, a nasal pediatric aspirator MARIMER (MA № ФСЗ 2011/09100).
- Infants are advised to rinse while lying down, children - sitting or standing, tilting the head to one side (it is recommended to take into account the general condition of the child and his age characteristics).
- Carefully insert the tip into the nasal passage and squeeze the dropper vial slightly.
- Repeat the procedure with a different nasal passage.
- Straighten the head, wipe the nose with a tissue or blow the nose.

TECHNICAL DESCRIPTION OF MEDICAL DEVICE

Nasal drops sterile MARIMER are manufactured in compliance with international and national standards.

Dropper vial of 5 ml ($\pm 5\%$) is made of material: low density polyethylene.

Used materials comply with the requirements of GOST R 52770 and the standards of the GOST ISO 10993 series by the biological effect on the patient. The product is biologically safe.

The product is sterile. Sterilization is carried out using aseptic production methods according to GOST ISO 13408-2-2007 (sterilizing filtration).

Weight of the product in a dropper vial is 6.45-7.90 g.

PRECAUTIONS FOR USE

- Store at the temperature from 15 to 30 °C.
- Keep out of reach of children!
- For infants – instill the solution with a minimum pressure to avoid any risk of affection of the middle ear.
- Dropper vial with sterile isotonic seawater solution is single-use and individual.
- Do not use opened or damaged dropper vial with sterile isotonic seawater solution due to a risk of rapid bacterial contamination.

STERILIZATION METHODS

Sterile: sterilization using the sterilizing filtration method.

OPERATIONAL, STORAGE AND TRANSPORTATION CONDITIONS

This medical device is transported by means of all types of transport in roofed vehicles in accordance with the requirements and rules for the transport of goods operating on vehicles of each type.

Transportation conditions: At a temperature from 15 to 30 °C and a relative humidity of NMT 80%.

During operation the device is resistant against influence of climatic aspects: temperatures from 15 to 30 °C and a relative humidity of NMT 80%.

Store at the temperature from 15 to 30 °C.

SHELF-LIFE





3 years from the date of manufacture. Do not use the device after the expiration date indicated on the package.

ENVIRONMENTAL REQUIREMENTS

This medical device does not adversely affect humans or the environment when used, transported or stored.

LABELLING

Graphic symbols on the labelling of medical product

Graphic symbols	Description
	Manufacturer
	Read the package leaflet carefully before use
	Sterilization using aseptic techniques
	Do not reuse
	"Throw in the bin" - the device and/or its package are disposed with household waste at home.
	Temperature range
	Complies with European Directive 93/42/EEC
	Catalogue number
	Batch number
	Expiry date
	Manufacturing date

PROCEDURE FOR DISPOSAL OR DESTRUCTION OF MEDICAL DEVICE

Special precautions for the destruction of an unused/used device are not required. The treatment and its disposal must be carried out in accordance with the standards accepted in medical practice, as well as relevant local, state and federal laws and instructions. In medical institutions, device and/or its package are disposed of as Class A medical waste in accordance with SanPin 2.1.7.2790-10. At home, the device and/or its packaging are disposed of with household waste.

LIST OF INTERNATIONAL STANDARDS

Nasal drops sterile MARIMER complies with the following standards:

Standard	Name
EN ISO 13485	Standard: Medical devices - QUALITY MANAGEMENT SYSTEMS - Requirements for regulatory purposes.
EN ISO 14971	Medical devices - Application of risk management to medical devices.
EN ISO 10993-1	Biological evaluation of medical devices - Part 1 Evaluation and testing within a risk management process.

EN ISO 10993-5	Biological evaluation of medical devices - Part 5 Tests for in vitro cytotoxicity.
EN ISO 10993-10	Biological evaluation of medical devices - Part 10 Tests for irritation and skin sensitization.
EN 556-2	Sterilization of medical devices – Requirements for medical devices to be designated "STERILE" – Part 2 Requirements for aseptically processed medical devices.
EN 1041	Information supplied by the manufacturer of medical device.
ISO 15223-1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1 General requirements.
EN ISO 13408-1	Aseptic processing of health care products — Part 1 General requirements.
EN ISO 13408-2	Aseptic processing of health care products — Part 2 Filtration.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Does not contain materials of animal and (or) human origin.

INFORMATION ABOUT DRUG PRODUCTS, MEDICINAL PRODUCTS AND PHARMACEUTICAL SUBSTANCES

Does not contain drug products, medicinal products or pharmaceutical substances.

WARRANTY

This device is designed, manufactured, tested and packaged in compliance with all relevant requirements. The manufacturer guarantees the quality of the medical device until the expiration date, provided that its package is not damaged, and storage conditions are met.

RECLAMATION

Reclamations to be send to the address of the authorized representative of the manufacturer in the Russian Federation:

Limited liability company «Abbott Laboratories» («Abbott Laboratories» LLC)
125171, Moscow, Leningradskoe shosse, 16A, build. 1
Tel./fax: +7(495)-258-42-80/81

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