Neuville-sur-Oise, March 10th, 2020

RE: Product Registration File for Russia

To Whom It May Concern,

Please find enclosed documents compiled for product registration of RevohemTM Arachidonic acid in Russia.

Product:

Product name	Ref. No.	Packaging size
Revohem TM Arachidonic acid	BV413997	3 x 0.625 mL

Documents enclosed:

Contents	Total N° of Page(s)
• Instruction for use (English for IFU in Russia)	6

Signature:

Name: Claire DUNOIS
Position: Clinical study and Validation Director
HYPHEN BioMed

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Paris, le

18 mai 2020

Name of the *in vitro* diagnostics medical device

Revohem Arachidonic acid, Reagent for the platelet aggregation testing

Intended use

For *in vitro* diagnostics
Measurement of platelet aggregation

Indication for use

Revohem Arachidonic acid is intended to be used for measurement and quantification of the maximal aggregation rate (%) of platelet-rich plasma (PRP) to assess the platelet function using the automated coagulation analyzers. Revohem Arachidonic acid is used for healthcare professional in aid to diagnosis for platelet function disorders or to assess responsiveness to antiplatelet drugs.

Intended user

The users are qualified medical professionals in a clinical diagnostic laboratory environment

Medical Device Use Risks

The class of potential risk of using a medical device in accordance with the nomenclature classification of medical devices is 2a.

There is no direct potential risk to the patient or person who is performing the analysis because the reagent is not in direct contact with the human body and is not introduced into the body. However, indirect risks may be the result of hazardous events associated with the erroneous result possibly yielded by the medical device.

Risk analysis has been performed in accordance with ISO 14971.

Risk assessment included: risk analysis, the device properties affecting safety, potential hazards, potential clinical harm and associated serious consequences, analysis of the clinical risk/ benefit ratio. Based on the analysis, all identified and expected hazards and associated risks were identified and reduced to acceptable levels.

Test Method Principles

When arachidonic acid is added to the platelet rich plasma (PRP) from a healthy subject, it is converted to thromboxane A₂ (TXA₂) by the action of cyclooxygenase and TX synthase. TXA₂ is a potent platelet activator and is required for platelet aggregation following certain types of stimulus.

After light transmission of PRP and platelet poor plasma (PPP) is set as 0 % and 100 % respectively, arachidonic acid is added to the PRP.

TXA₂ is generated and then platelets start to aggregate causing a change in light transmission which is measured at a wavelength of 660 nm.

The maximal aggregation (%) is calculated by the maximal amount of change in light transmission caused by platelet

aggregation induced by arachidonic acid.

Medical Device Description

Components

- Revohem Arachidonic acid reagent - 3 vials, 12 mmol / L arachidonic acid after reconstitution,
- Instruction for use.

Packaging

Revohem Arachidonic acid is filled in the clear glass transparent vial with white screw-on caps. The reagents (3 vials) are packed in the outer packs made of bleached board and the outer packs are taped over with adhesive labels.

Reagent	Inner package	Outer package	
		Number in package, vials.	Size (length x width x height), mm
Arachidonic acid reagent	glass vials with white PP screwing caps	3	80 x 36 x 58

Labelling

Labelling of the inner packaging of the components contains:

- Manufacturer's name;
- Full and abbreviated name of the product;
- Volume of the vial (after reconstitution);
- Batch number;
- Expiry date;
- Storage conditions;
- Symbol of: « *In vitro* diagnostic medical device(IVD)»;
- Toxicity symbol;
- CE mark;
- «SYSMEX» graphic symbol
- «Reconstitute with» symbol;
- Identification Number.

The outer package labeling contains:

- Name and address of manufacturer;
- Full and abbreviated name of the reagent;
- Volume of the vial (After reconstitution);
- Components (number of vials) of reagents
- Catalogue number;
- Batch number;
- Expiry date;
- Storage conditions;
- Symbol of: «Medical device for *in vitro* diagnostics (IVD)»;
- Toxicity symbol;
- Hazard statement;
- Symbol of «Consult Instructions for Use»;
- CE mark;
- «SYSMEX» graphic symbol
- Barcode;
- «Reconstitute with» symbol;
- Identification number;

- Statement: «Sysmex is registered with the United States Patent and Trademark Office».
Labeling is in several languages.

Additional labelling for the Russian Federation in Russian includes the following:

- Name of the medical device;
- Components (Number of vials);
- Registration certificate (MA) number;
- Information on the representative agent in the territory of the Russian Federation.
- Symbol of "Do not reuse"

Contraindications

When used by specially trained and qualified personnel and taking into account the intended use, no contraindications have been identified.

Possible side effects

When used by the qualified and trained personnel and taking into account the appropriate safety precautions there have been no side effects detected.

Warnings and Precautions

1. This reagent is for in-vitro diagnostic use only. Do not use for any other purpose. If the reagent is used in a way other than that described, performance cannot be guaranteed.
2. When handling blood samples, take suitable measures to avoid the potential risk of infection.
3. In case of reagent contact with the mouth or eyes, rinse thoroughly with water and seek medical advice.
4. Reagent vials are fragile and may possibly crack or break whilst being handled. Please take care when handling, especially when opening.
5. Use the reagent within the expiration date displayed on the label.
6. Do not mix old reagent with fresh reagent.
7. Air bubbles generated in the reagent might cause abnormal results.
8. Handle the reagent with care to prevent contamination with microorganisms and dust.
9. Do not use the reagent vials for other purposes.
10. Reading the instruction for use (IFU) is mandatory.
11. Only use citrated blood samples.
12. Vials are closed under vacuum. Carefully remove the stopper in order to avoid any loss of powder when either opening the vials or reconstituting the reagent.
13. Arachidonic acid is susceptible to oxidation and turns yellow as it oxidizes.
14. Oxidized arachidonic acid deteriorates and it may affect reagent performance.
15. To avoid oxidation by exposure to air for a long time, please reconstitute immediately after opening and keep with its vial cap.
16. The hazard and toxic information written on the label is classified by Regulation (EC) No 1272/2008. Please see the Safety Data Sheet for the country-specific classification.

Warning Phrases

- H319 Causes serious eye irritation.
- P264 Wash skin thoroughly after handling.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P305 If in eyes: wash with water thoroughly for several minutes. Remove the contact lenses if any. Keep washing with water.
- P338 If irritation persists: get medical advice/aid.
- P313 Harmful if swallowed
- H302 Causes skin irritation.
- H315 May cause respiratory tract irritation
- H318 Causes serious harm to eyes

Component	Substance under Class 1272/2008 (CLP)	Concentration
Butylated hydroxyanisole	H302, H315, H319, H335, H351, H318	<1%
Sodium metabisulphite		<3%

Additional Required Equipment and Materials

- 0.9% NaCl solution;
- Distilled water;
- Medical examination gloves;
- Sysmex automated coagulation analyser with the platelet aggregation measurement;
- Sample cup, 4 mL
- Sample tube with stirring bar;
- Sample cup 4 mL holder;
- SB set tool;
- Plastic or siliconized glass tubes with sodium citrate 3.2%;
- Vacuum systems for blood sampling with sodium citrate 3.2%;
- Venipuncture needles 21G or bigger;
- Laboratory centrifuge (allowing centrifugal acceleration of 200g and 1500g);
- Variable volume pipettors;
- Single-use tips.

Reagent preparation

Remove the reagents from the refrigerator and allow them to equilibrate to ambient temperature. Reconstitute the arachidonic acid reagent with 0.625 µl of distilled water and mix thoroughly (using a vortex mixer). Leave the vial for 30 minutes (with occasional mixing) to obtain a homogeneous suspension (the concentration of reconstituted arachidonic acid is 12 mmol / L). Powder dissolution time should not exceed 30 min. Dilute the reconstituted arachidonic acid 1.5 fold with a 0.9% NaCl solution and dispense it into a sample cup for use (working concentration of arachidonic acid is 8.0 mmol / L). The diluted arachidonic acid solution must be immediately loaded on the analyser after dilution and be used within 8 hours.

Platelet aggregation testing is strongly influenced by sample handling. Local guidelines should be followed with regard to sample collection.

1. Blood collection

- (1) Blood sampling should be by clean venepuncture with minimal venous stasis, using needles no smaller than 21 gauge. Discard any samples with in vitro haemolysis.
- (2) Use only plastic or siliconized glass collection tubes, with citrate as the anticoagulant (9 parts blood: 1 part sodium citrate solution).
- (3) After blood collection, mix the tube(s) gently.

2. Centrifugation

- (1) Centrifuge the blood at 200 g for 10 minutes at room temperature.
- (2) Carefully transfer the PRP with a plastic pipette into a plastic tube and avoid bubbles. Avoid the transfer of other cell components.
- (3) Store the PRP in a closed tube at room temperature.
- (4) Centrifuge the primary tube once more at 1500 g for 15 minutes.
- (5) Transfer the resulting PPP to a plastic tube and store it at room temperature until use.

3. Perform and complete the platelet aggregation measurements within 4 hours after blood sampling.

Examination procedure

This kit is adapted only for Sysmex automated coagulation analysers with the platelet aggregation analysis function enabled.

- 1. The analyser adds 20 µL of 0.9% NaCl solution to 140 µL of PPP and the transmitted light intensity at a wavelength of 660 nm is set as the baseline (100%).
- 2. 20 µL of arachidonic acid is added to 140 µL of PRP (Final concentration of arachidonic acid: 1.0 mmol/L) and the change in transmission at 660 nm is measured. The analyser calculates the maximal aggregation rate from the transmitted light intensity.

Transportation and storage conditions

The medical device is transported in covered vehicles in accordance with the rules for the carriage of goods in force for each type of transport at ambient temperature. The medical device should be stored at a temperature of + 2 ° C to + 8 ° C until the expiry date indicated on the label. The shelf life of an unopened reagent is 18 months after the manufacture date.

Storage and shelf life after the first opening

Table 1. Storage temperature and stability

Storage Temperature	Stability
From 18 to 25 °C	24 hours
From 2 to 8 °C	7 days
From -25 to -15 °C	2 months

The reconstituted arachidonic acid can be frozen to be used only once.

Control method

Control materials and calibrators are not available. It is recommended to validate each new lot of reagent using samples collected by normal subject.

Interpretation of results

- 1. The results are evaluated by examining the aggregation curve and the % maximum aggregation. These will vary depending on the instrument type, and each laboratory should determine their own normal values.
- 2. Results should be interpreted according to the patient's clinical condition, platelet count, potential influences of medication, life style, nutrition, and pre-analytical conditions.
- 3. Abnormal curves should be confirmed by re-testing.

Performance characteristics

Table 2. Example of the maximal platelet aggregation: normal and abnormal (%)

	Normal	Abnormal
Maximal aggregation (%)	90.2	2.3

Abnormal: A patient with Glanzmann disease


Figure 1. Example of normal and abnormal platelet aggregation curves

Biological reference intervals

Biological reference intervals should be established by each laboratory in the actual working conditions.

Limitations of the examination procedure

The use of this reagent is validated on the Sysmex automated coagulation analysers with the function to measure platelet aggregation. HYPHEN BioMed cannot take responsibility for patient results obtained by the use of HYPHEN BioMed reagents on unauthorised analysers, therefore it is the responsibility of the user to validate the performance of the test when modifications to these instructions or use of the reagent were done. No significant interference has been observed on CS analyzer for bilirubin concentration up to 28 mg/dL, hemoglobin up to 200 mg/dL and for triglycerides up to 286 mg/dL. If the number of platelets is less than 150 × 10⁹/L or greater than 480 × 10⁹/L, the test results might be affected by it.

Analytical Properties

Highly platelet aggregation sample: PRP from healthy, non-bleeding problem individual

Lowly platelet aggregation sample: PRP with antiplatelet drug

1. Sensitivity

Test methods:

Measure a highly platelet aggregation sample and a lowly platelet aggregation sample for 3 times each, and calculate the difference of the maximal aggregation rates (%) between a highly platelet aggregation sample and a lowly platelet aggregation sample.

Results:

- a) The maximum aggregation rate is 50% or more for a highly platelet aggregation sample.
- b) The difference of the maximum aggregation rates between the measured results of a highly platelet aggregation and lowly platelet aggregation sample is 15% or more.

2. Accuracy

Test methods:

Measure a highly platelet aggregation sample and a lowly platelet aggregation sample for 3 times each.

Results:

The maximum aggregation rate is 50% or more in highly platelet aggregation sample, and below 50% in lowly platelet aggregation sample.

3. Repeatability

Test methods:

Measure a highly platelet aggregation sample and a lowly platelet aggregation sample for 5 times each, and calculate the coefficient of variation (CV%) of the maximal aggregation rate (%).

Results:

The coefficient of variation (CV%) is 20% or less in highly platelet aggregation sample.
The maximal coefficient of variation (CV) is 30% or less, or standard deviation (SD) is within 9% in highly platelet aggregation sample.

4. Lot-to-lot (Between-lot imprecision)

Test methods:

A highly platelet aggregation sample is tested with 3 different batches of reagent, 5 times for each batch. Calculate the mean value of 15 measurements and the coefficient of variation (CV%) of the maximal aggregation rate (%).

Results:

Repeated tests of highly platelet aggregation sample with 3 different lots of reagents, the coefficient of variation of the obtained results (CV%) is 20% or less.

Table 3. Repeatability (Within-run) An example for information

	Mean maximal aggregation (%)	Within run repeatability (CV,%)	Within-run SD
Normal sample	83.4	3.1	2.6
Abnormal sample	6.6	16.1	1.1

The results have been obtained from 30 repeated

measurements.

Manufacturer’s Warranty

The manufacturer guarantees that the reagent complies with the requirements of the applicable national and international standards. The manufacturer is solely responsible for the performance of the medical device only if it is used in accordance with the instructions, and its storage and transportation are carried out in accordance with the specified conditions and rules. The manufacturer guarantees the safety of the medical devices, the absence of unacceptable risks to human health, as well as to the environment, if the reagents are used in accordance with the intended purpose and conditions specified by the manufacturer.

Reagents are to be used before the expiry date printed on the package.

This limited warranty supersedes any other express or implied warranties, including any implied warranties of merchantability for commercial purpose or for any other use. The manufacturer is not responsible for any accidental, indirect or intentional losses.

Single Use Device

The medical device is single use only.

Cleaning, Disinfection, Sterilization

Single use device, delivered non-sterile, non-sterilizable.

Environment protection requirements in use of the medical device

When used as intended and subject to the precautions recommended by the manufacturer, the reagent is not a source of harmful substances that pollute the environment in excess of the established standards. Prevent the spread, leakage and runoff of material, if possible; do not allow the reagent to enter sewage.

Disposal Procedures

Disposal of waste should be disposed of in accordance with national regulations and in consultation with environmental authorities. The waste class is assigned in accordance with national regulations. Dispose of the hazardous or bio-infected materials in accordance with your facility rules and with the federal, and local regulations. The remaining samples and waste obtained during the various stages of processing the sample, as well as all components of the medical device that come into contact with biological samples, should be disposed of by using disinfectants in accordance with the national requirements for the safe handling of pathogenic microorganisms and parasitic pathogens, and with the requirements for medical waste management (in the Russian Federation - as class B waste (epidemically hazardous) in accordance with SanPiN 2.1.7.2790-10)

Cardboard boxes are disposed of as Class A waste, expired reagents and empty vials as Class G waste. Packages must be completely emptied.

Waste disposal and destruction shall be performed by properly licensed organizations at specialized sites in accordance with the requirements provided for by federal laws and in accordance with mandatory environmental

safety requirements, using methods agreed upon with local authorities responsible for the sanitary and epidemiological welfare of the population.

References

(1) Born G. V.: Aggregation of blood platelets by adenosine diphosphate and its reversal.: Nature 1962, **194**, 927–9.

(2) Cattaneo M., Cerletti C., Harrison P. et al.: Recommendations for the Standardization of Light Transmission Aggregometry: A Consensus of the Working Party from the Platelet Physiology Subcommittee of SSC / ISTH.: J Thromb Haemost 2013, **11**, 1183–9.

(3) Lawrie A. S., Kobayashi K., Lane P. J. et al.: The automation of routine light transmission platelet aggregation. : Int J Lab Hematol 2014, 36 (4), 431–8.

(4) CLSI. Platelet Function Testing by Aggregometry; Approved Guideline. CLSI document H58-A (ISBN 1-56238-683-2). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA 2012.

manufacturer

 **HYPHEN BioMed**
155 rue d’Eragny, 95000 Neuville-sur-Oise,
France

Authorized Representative of Manufacturer in the Russian Federation

Sysmex RUS Limited Liability Company (Sysmex RUS LLC)
Registered Address: 1 Magistralniy tucic, 11, building 10, floor 2,
space V, rooms 1-23, 25-39, Moscow, 123290
Postal Address: Magistralniy tucic, 11, building 10, office 1020,
Moscow, 123290
Telephone: +7(495)781-67-72

Information on Product
(Revohem Arachidonic acid regent)

Revohem Arachidonic acid Reagent

RCNS	0.625 mL
------	----------

 × 3

List of the applicable International Documents/Standards

The medical device complies with the requirements of the following standards and directives:

Standard	Title of standard
ISO 9001	Quality Management Systems. Requirements
ISO 13485	Medical devices — Quality management systems — Requirements for regulatory purposes
ISO 14971	Medical devices — Application of risk management to medical devices
EN 13641	Elimination or reduction of risk of

Standard	Title of standard
	infection related to in vitro diagnostic reagents
November 9 2004 , decree « NOR: SOCT0412152 A »	Chemical Products Safety Data Sheet – Part 1: content: procedure for sections and transfers
EN 18113-1	In vitro diagnostic medical devices — Information supplied by the manufacturer (labelling) — Part 1: Terms, definitions and general requirements
EN 18113-2	In vitro diagnostic medical devices — Information supplied by the manufacturer (labelling) — Part 2: In vitro diagnostic reagents for professional use
EN 15223-1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements
EN 13612	Laboratory assessment of the technical properties of the medical devices
EN 13640 EN ISO 23640	In vitro diagnostic medical devices. Evaluation of stability of in vitro diagnostic reagents
EN ISO 17511	In vitro diagnostic medical devices — Measurement of quantities in biological samples — Metrological traceability of values assigned to calibrators and control materials
Directive 98/79 CE	Directive 98/79/EC by European Parliament and Council dated 27.10.1998 on medical devices for <i>in vitro</i> diagnostics.
Regulations (EC) #1069/2009	Health rules as regards animal by-products and derived products not intended for human consumption
EN 62366	Application of usability engineering to medical devices
Regulations (EC) 1907/2006 REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH)
Regulations (EC) CLP: 1272/2008	Classification, labelling and packaging of substances and mixtures

	Manufacturer
	Use by
	Batch Code
	In vitro diagnostic medical device
	Catalogue number
	Warming
	CE mark
	Temperature Range
	Consult Instructions for Use
	Reconstitute with
	Do not reuse

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