
Instructions for Use
Merit PAK
Dated May 26, 2020
Consisting of 29 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 26th day of May 2020, before me, Jennifer Fordham, a notary public, personally appeared Cory Marsh, Manager, Regulatory Operations proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public

My Commission Expires: June 27, 2021

Cory Marsh
Manager, Regulatory Operations
Merit Medical Systems, Inc.
26 May 2020

**Instructions for use for
MEDICAL DEVICE**

**« System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened,
PAK Hydrophyllic »**

**Manufactured by
Merit Medical Systems, Inc., the USA**

**Instructions for use for
MEDICAL DEVICE**


« System for percutaneous mini-access in kits Merit PAK, models: PAK»

PRODUCT DESCRIPTION

The Merit PAK: model PAK utilizes a small coaxial introducer with dilator and guide wire for placement of larger diameter guide wires into the vasculature system when a small needle stick is preferred.

The Merit PAK: model PAK consists of a combination of the following components, see unit label for specific product components. These components may be packaged in a single pouch or may be packaged separately.

- Coaxial introducer/dilator 4F(1.3mm), 5F (1.7 mm) - 1 pcs.
- Guide Wire 0.018" (0.46mm) - 1 pcs.
- Haemostasis valve with 3-way stopcock - 1 pcs.
- Echo Enhanced Introducer Needle 21G (40 mm or 70 mm length) -1 pcs. (if necessary)
- Instruction for use -1 pcs

models: PAK, PAK Stiffened

1. Co-axial introducer/dilator 4F/5F
2. Guide Wire 0.018" (0.46mm)
3. Echo Enhanced Introducer Needle 21G (40 or 70 mm length)
4. Haemostasis valve with 3-way stopcock

CATEGORIZATION OF DEVICE(S)

Devices risk class is 2a

The categorization of medical device "System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophyllic" see below.

DEVICE	CATEGORIZATION
System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophyllic (introducers, dilators,	Externally Communicating Device with Circulating Blood Contact for a Limited Duration (< 24 hours). Products are intended for invasive use (use temperature from 32 ° C to 44 ° C, humidity 100% (direct contact with

<u>DEVICE</u>	<u>CATEGORIZATION</u>
guidewires, needles)	blood)
Hemostasis valve	Externally Communicating Device with Indirect Blood Contact for a Limited (< 24 hours) Duration.

INTENDED USE

System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophilic is designed to perform the "micropuncture" technique, which helps in facilitating percutaneous placement of medical devices of larger diameter (catheters, guide wires, balloons, stents) into the patient's body.

WARNINGS

Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. Do not withdraw with force. Do not advance the guide wire if resistance is met.

CAUTIONS

- Read instructions prior to use.
- This device is intended for single use only. Do not reuse or resterilize.
- This device is sterile if package is unopened or undamaged.
- This device is non-pyrogenic.
- This device should be used by clinicians with adequate training in the use of the device.
- Utilize appropriate anticoagulant therapy for patient during procedure.
- Prior to use, ensure that the sheath and dilator are the appropriate size

CONTRAINDICATIONS

There are no known contraindications.

POSSIBLE COMPLICATIONS

The possible complications are (but not limited) air embolism, thrombus formation, infection, pseudo aneurysm formation, hematomas, guide wire embolization, bleeding, perforation or laceration of the vessel wall and standard risks connected with percutaneous diagnostic and (or) interventional procedures.

INFORMATION ABOUT POTENTIAL USERS AND USE CONDITIONS

To be use only in medical facilities by medical staff (interventional surgeons, vascular surgeons, specialists in diagnostic and surgical departments and centers) who has been specially trained how to work with this type of medical device.

Application conditions

Room temperature range	from + 15°C to + 30°C
Room humidity	< 80% (without condensation)
Atmospheric pressure	84.0-106.7 kPa (630-800mm Hg)
Room elevation	Max 1500 m a.s.l.

INSTRUCTIONS FOR USE

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
2. Insert the 21 gauge introducer needle using standard technique.
3. Carefully advance the flexible end of the 0.018" (0.46mm) guide wire through needle. Advance the guide wire as far as appropriate. Verify correct positioning. WARNING: Do not advance guide wire if resistance is met.
4. Remove needle while maintaining the 0.018" (0.46mm) guide wire in position. WARNING: Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. To avoid guide wire damage during manipulation, remove the introducer needle and proceed to step 5.
5. Insert the coaxial introducer/dilator pair over the 0.018" (0.46mm) guide wire and advance to the desired position.
6. Remove dilator and 0.018" (0.46mm) guide wire leaving the introducer in position. NOTE: Place a finger over the hub of the introducer to minimize blood loss and risk of air aspiration.
7. Insert 0.035" (0.89mm) or 0.038" (0.97mm) guide wire through introducer.
8. Remove introducer, leaving 0.035" (0.89mm) or 0.038" (0.97mm) guide wire in place.

MAIN PRODUCT SPECIFICATION

CHARACTERISTIC	SPECIFICATION
Coaxial couple (introducer) , hemostasis valve	
Introducer, Inner Diameter (ID)	4F(1,3 mm): 0.50 mm 5F(1,7 mm): 0.50 mm
Outer Diameter (OD)	4F 1,3 mm 5F 1,7 mm
Introducer, effective length	4F (1,3 mm): 10 cm 5F (1,7 mm): 10 cm
Dilator, Outer Diameter (OD)	4F (1,3 mm) 5F (1,7 mm)
Dilator, Inner Diameter (ID)	4F(1,3 mm): 0.45 mm 5F(1,7 mm): 0.46 mm
Dilator, effective length	4F,5F-12.5 cm
Dilator, taper length	4F (1,3 mm): 5,0 mm 5F (1,7 mm): 5,0 mm
Dilator,weight	≤ 5 g
Introducer, weight	≤ 5 g
Hemostasis valve	≤ 15 g
Hub Tensile Test	Force to remove hub shall be greater than 2.25 lbs (10N)
Freedom from leakage from introducer	No leakage at test pressure of 300 kPa
Burst strength (Introducer and introducer to hub bond)	Minimum burst strength - 4F /5F sizes ≥ 3.37 lbs (14.99N)
Burst strength (introducer shaft to the marker tip segment bond)	Minimum burst strength - 4F/ 5F sizes ≥ 3.37 lbs (14.99N)

CHARACTERISTIC	SPECIFICATION
Introducer Hub/Cap Snap Fit Tensile force to remove cap from hub	Minimum Force > 9.0 lbs. (40.03N)
Sidearm Pull Test	Minimum Peak Break Force $\geq$ 3.37 lbs (14.99N)
Sidearm Leak Pressure Test	Sidearm shall not leak or detach at $\leq$ 90 psi (6.328 кг/см ²)
Sidearm flow rate	> 600 ml/min (200mm Hg with water)
Hub to dilator bond strength	Minimum Burst strength: 4F/ 5F sizes $\geq$ 2.25 lbs
Dilator to Introducer Hub Snap Fit	4 F Minimum Force: 4.75 (21,13H) 5 F Minimum Force: 4.5 (20,02 H)
Bond Strength – Dilator (Hub to Rod)	Minimum Burst Strength 4F/5F > 2.25 lbs (10,01 H)
Dilator to Introducer Hub Snap Fit	Minimum Force: 4.5 lbs (20,02)
GuideWire 0,018" (0.46 mm)	
Outside Diameter	0.46 mm (0,018")
Product Lengths	40 cm
Minimum Breaking Strength (Pull Strength)	10 N / 2.25 lbs minimum
The strength of the connection rod and spiral (tensile strength)	$\geq$ 0.5 lbs (2.2 H). Minimum tensile force
Wight (in a case)	$\leq$ 5 g
Echo Enhanced Introducer Needle 21G (0.8 mm)	
The outer diameter of the needle	0,80 mm
Inner diameter	0,547 mm
Length of needle tube:	70 mm 40 mm
Tensile strength at break structure "needle-hub" (mounting strength)	No less than 44N
Bending force	15N (maximum deviation 0.50 mm)
Type of sharpening	trihedral
Allowable angle deviations needle without compromising the integrity of	short cut at an angle of $17 \pm 2^\circ$
needle roughness	No more than 0,4 mcm
The radius of the working portion of blunting needle	No more than 0,03 mm
Weight	$\leq$ 2 g

MATERIAL SPECIFICATIONS

Component	Materials
Coaxial introducer/dilator Model PAK	
Dilatator 4F/5F (tube)	HDPE Barium sulphate
Dilatator 4F/5F (hub)	HDPE
Dilatator 4F/5F (nut)	ABS Colorants: 4F Red /5F Light Gray

Introducer 4F/5F (tube)	Pebax Barium sulphate Titanium Dioxide
Introducer 4F/5F (hub)	Polyamide Titanium Dioxide Ink for Printing: 4F Red /5F Light Gray
Protection Tube	LDPE
Haemostasis valve with 3-way stopcock	
Valve (body)	Polycarbonate
Valve (sheath)	Medical silicone
Valve (cup)	Polycarbonate Colorant: Blue
Stopcock Body	Polycarbonate
Stopcock Handle	High viscosity Acetal homopolymer Colorant: Blue
Stopcock Lubricant	Polydimethylsiloxane lubricant
Echo Enhanced Introducer Needle 21G (0.8 mm) length 40 mm, 70 mm.	
Needle Tube:	Stainless steel
Needle Base	Polycarbonate Printing Ink Color: black
Protective tube:	LDPE
Nithinol Guidewire 0,018" (0,46 mm)	
Nithinol Guidewire 0,018" (0,46 mm)	Nithinol /Pd
Package materials	
Package	Tyvek® 1059B / LLPDE

STERILIZATION METHOD

Ethylene Oxide (EO) sterilization. This product does not provide for a user re-sterilization.

DESCRIPTION OF PACKAGING WITH INDICATION OF NUMBER OF PRODUCTS AND ACCESSORIES IN PACKING



System for percutaneous mini-access in kits Merit PAK, models PAK, PAK Stiffened includes coaxial pair (introducer / dilator) 4 F or 5 F, 0.018 "(0.46 mm) guidewire, 21G (0.8mm) introducer needle, hemostatic valve with 3-way stopcock. Kit components are packaged in individual sterile disposable packaging of Tyvek® 1059B. The kit contents may vary depending on the physician's need in a certain component and the clinical case.

Five pcs of individual-packed kits are placed in a transport box of 200 gr. Natural Kraft E-flute carton.

INFORMATION ON LABELING OF PRODUCTS

Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

Description of the marking symbol is represented in the table.

	Manufacturer
	Authorized representative in the European Community
	Manufacturing date
	Do not reuse
	Catalogue #
	Batch Code
	Sterilized by ethylene oxide
	Caution: United States federal law restricts this device to sale by or approximately a physician
	Use by date
	Attention, Refer to the instructions for use
	Do not use, if package is damaged
	Non-pyrogenic

TRANSPORTATION, STORAGE AND HANDLING

The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures. Rotate inventory so that device is used prior to the Use By date on the package label. Products released from the transport packaging should be kept at an ambient condition:

Storage conditions:

Temperature +5°C to +27 °C (environment)
Relative humidity 10% to 90 % (non-condensing)
Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Transport conditions:

Temperature -35°C to +55 °C (environment)
Relative humidity 10% to 90 % (non-condensing)
Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Usage environment:

Temperature range for use - the body temperature of 32 °C to 44 °C
Relative humidity: 100% (blood contact)
Atmospheric pressure 86-106 kPa.

THE SHELF LIFE

Of sterile devices – 36 months (3 years)

DISPOSAL

To avoid health and environmental risks during disposal of medical devices they shall be disposed of in compliance with the hospital-approved protocol.

Medical Device " System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophyllic " according to the Sanitary Rules and Regulations 2.1.7.2790-10 is related to Class B contaminated and potentially pathological wastes (epidemiologically hazardous wastes).

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/ or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

LIMITED WARRANTY AND DISCLAIMER

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

LIMITATION OF LIABILITY

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER

PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

COMPLAINT

The organization, which takes claims from consumers about the product quality on the Russian Federation territory:

Limited Liability Company "Merit Technologies"

(LLC "Merit Technologies"),

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

Tel / fax: +7 (495) 221-89-02

Instructions for use for

MEDICAL DEVICE


« System for percutaneous mini-access in kits Merit PAK, models: PAK Stiffened»

PRODUCT DESCRIPTION

The Merit PAK: model PAK Stiffened utilizes a small coaxial introducer with dilator and guide wire for placement of larger diameter guide wires into the vasculature system when a small needle stick is preferred.

The Merit PAK: model PAK Stiffened consists of a combination of the following components, see unit label for specific product components. These components may be packaged in a single pouch or may be packaged separately.

- Coaxial introducer/dilator 4F(1.3mm), 5F (1.7 mm) - 1 pcs.
- Guide Wire 0.018" (0.46mm) - 1 pcs.
- Haemostasis valve with 3-way stopcock - 1 pcs.
- Echo Enhanced Introducer Needle 21G (40 mm or 70 mm length) -1 pcs. (if necessary)
- Instruction for use -1 pcs

models:PAK, PAK Stiffened

1. Co-axial introducer/dilator 4F/5F
2. Guide Wire 0.018" (0.46mm)
3. Echo Enhanced Introducer Needle 21G (40 or 70 mm length)
4. Haemostasis valve with 3-way stopcock

CATEGORIZATION OF DEVICE(S)

Devices risk class is 2a

The categorization of medical device "System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophillic" see below.

DEVICE	CATEGORIZATION
System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophillic (introducers, dilators, guidewires, needles)	Externally Communicating Device with Circulating Blood Contact for a Limited Duration (< 24 hours). Products are intended for invasive use (use temperature from 32 ° C to 44 ° C, humidity 100% (direct contact with blood)
Hemostasis valve	Externally Communicating Device with Indirect Blood

<u>DEVICE</u>	<u>CATEGORIZATION</u>
	Contact for a Limited (< 24 hours) Duration.

INTENDED USE

System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophilic is designed to perform the "micropuncture" technique, which helps in facilitating percutaneous placement of medical devices of larger diameter (catheters, guide wires, balloons, stents) into the patient's body.

WARNINGS

Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. Do not withdraw with force. Do not advance the guide wire if resistance is met.

CAUTIONS

- Read instructions prior to use.
- This device is intended for single use only. Do not reuse or resterilize.
- This device is sterile if package is unopened or undamaged.
- This device is non-pyrogenic.
- This device should be used by clinicians with adequate training in the use of the device.
- Utilize appropriate anticoagulant therapy for patient during procedure.
- Prior to use, ensure that the sheath and dilator are the appropriate size for the access vessel and devices to be use.

CONTRAINDICATIONS

There are no known contraindications.

POSSIBLE COMPLICATIONS

The possible complications are (but not limited) air embolism, thrombus formation, infection, pseudo aneurysm formation, hematomas, guide wire embolization, bleeding, perforation or laceration of the vessel wall and standard risks connected with percutaneous diagnostic and (or) interventional procedures.

INFORMATION ABOUT POTENTIAL USERS AND USE CONDITIONS

To be use only in medical facilities by medical staff (interventional surgeons, vascular surgeons, specialists in diagnostic and surgical departments and centers) who has been specially trained how to work with this type of medical device.

Application conditions

Room temperature range	from + 15°C to + 30°C
Room humidity	< 80% (without condensation)
Atmospheric pressure	84.0-106.7 kPa (630-800mm Hg)
Room elevation	Max 1500 m a.s.l.

INSTRUCTIONS FOR USE

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
2. Insert the 21 gauge introducer needle using standard technique.
3. Carefully advance the flexible end of the 0.018" (0.46mm) guide wire through needle. Advance the guide wire as far as appropriate. Verify correct positioning. WARNING: Do not advance guide wire if resistance is met.
4. Remove needle while maintaining the 0.018" (0.46mm) guide wire in position. WARNING: Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. To avoid guide wire damage during manipulation, remove the introducer needle and proceed to step 5.
5. Insert the coaxial introducer/dilator pair over the 0.018" (0.46mm) guide wire and advance to the desired position.
6. Remove dilator and 0.018" (0.46mm) guide wire leaving the introducer in position. NOTE: Place a finger over the hub of the introducer to minimize blood loss and risk of air aspiration.
7. Insert 0.035" (0.89mm) or 0.038" (0.97mm) guide wire through introducer.
8. Remove introducer, leaving 0.035" (0.89mm) or 0.038" (0.97mm) guide wire in place.

MAIN PRODUCT SPECIFICATION

CHARACTERISTIC	SPECIFICATION
Coaxial couple (introducer) , hemostasis valve	
Introducer, Inner Diameter (ID)	4F(1,3 mm): 0.50 mm 5F(1,7 mm): 0.50 mm
Outer Diameter (OD)	4F 1,3 mm 5F 1,7 mm
Introducer, effective length	4F (1,3 mm): 10 cm 5F (1,7 mm): 10 cm
Dilator, Outer Diameter (OD)	4F (1,3 mm) 5F (1,7 mm)
Dilator, Inner Diameter (ID)	4F(1,3 mm): 0.45 mm 5F(1,7 mm): 0.46 mm
Dilator, effective length	4F,5F-12.5 cm
Dilator, taper length	4F (1,3 mm): 5,0 mm 5F (1,7 mm): 5,0 mm
Dilator,weight	≤ 5 g
Introducer, weight	≤ 5 g
Hemostasis valve, weight	≤ 15 g
Hub Tensile Test	Force to remove hub shall be greater than 2.25 lbs (10N)
Freedom from leakage from introducer	No leakage at test pressure of 300 kPa
Burst strength (Introducer and introducer to hub bond)	Minimum burst strength - 4F /5F sizes ≥ 3.37 lbs (14.99N)
Burst strength (introducer shaft to the marker tip segment bond)	Minimum burst strength - 4F/ 5F sizes ≥ 3.37 lbs (14.99N)

CHARACTERISTIC	SPECIFICATION
Introducer Hub/Cap Snap Fit Tensile force to remove cap from hub	Minimum Force > 9.0 lbs. (40.03N)
Sidearm Pull Test	Minimum Peak Break Force ≥ 3.37 lbs (14.99N)
Sidearm Leak Pressure Test	Sidearm shall not leak or detach at ≤ 90 psi (6.328 кг/см ²)
Sidearm flow rate	> 600 ml/min / 200mm Hg with water
Hub to dilator bond strength	Minimum Burst strength: 4F/ 5F sizes ≥ 2.25 lbs
Dilator to Introducer Hub Snap Fit	4 F Minimum Force: 4.75 (21,13H) 5 F Minimum Force: 4.5 (20,02 H)
Bond Strength – Dilator (Hub to Rod)	Minimum Burst Strength 4F/5F > 2.25 lbs (10,01 H)
Dilator to Introducer Hub Snap Fit	Minimum Force: 4.5 lbs (20,02)
GuideWire 0,018" (0.46 mm)	
Outside Diameter	0.46 mm (0,018")
Product Length	40 cm
Minimum Breaking Strength (Pull Strength)	10 N / 2.25 lbs minimum
The strength of the connection rod and spiral (tensile strength)	≥ 0.5 lbs (2.2 H). Minimum tensile force
Wight (in a case)	≤ 5 g
Echo Enhanced Introducer Needle 21G (0.8 mm)	
The outer diameter of the needle	0,80 mm
Inner diameter	0,547 mm
Length of needle tube:	70 mm 40 mm
Tensile strength at break structure "needle-hub" (mounting strength)	No less than 44N
Bending force	15N (maximum deviation 0.50 mm)
Type of sharpening	trihedral
Allowable angle deviations needle without compromising the integrity of	short cut at an angle of $17 \pm 2^\circ$
needle roughness	No more than 0,4 mcm
The radius of the working portion of blunting needle	No more than 0,03 mm
Weight	≤ 2 g

MATERIAL SPECIFICATIONS

Component	Materials
Coaxial introducer/dilator Model PAK	
Dilator 4F/5F (tube)	Stainless steel Polyamide Barium Sulphate Colorant blue
Dilator 4F/5F (hub)	HDPE

Dilatator 4F/5F (nut)	ABS Colorants: 4F Red /5F Light Gray
Introducer 4F/5F (tube)	Pebax Barium sulphate Titanium Dyoxide
Introducer 4F/5F (hub)	Polyamide Titanium Dyoxide Ink for Printing: 4F Red /5F Light Gray
Protection Tube	LDPE
Haemostasis valve with 3-way stopcock	
Valve (body)	Polycarbonate
Valve (sheath)	Medical silicone
Valve (cup)	Polycarbonate Colorant: Blue
Stopcock Body	Polycarbonate
Stopcock Handle	High viscosity Acetal homopolymer Colorant: Blue
Stopcock Lubricant	Polydimethylsiloxane lubricant
Echo Enhanced Introducer Needle 21G (0.8 mm) length 40 mm, 70 mm.	
Needle Tube:	Stainless steel
Needle Base	Polycarbonate Printing Ink Color: black
Protective tube:	LDPE
Nithinol Guidewire 0,018" (0,46 mm)	
Nithinol Guidewire 0,018" (0,46 mm)	Nithinol /Pd
Package materials	
Package	Tyvek® 1059B / LLPDE

STERILIZATION METHOD

Ethylene Oxide (EO) sterilization. This product does not provide for a user re-sterilization.

DESCRIPTION OF PACKAGING WITH INDICATION OF NUMBER OF PRODUCTS AND ACCESSORIES IN PACKING

System for percutaneous mini-access in kits Merit PAK, models PAK, PAK Stiffened includes coaxial pair (introducer / dilator) 4 F or 5 F, 0.018 "(0.46 mm) guidewire, 21G (0.8mm) introducer needle, hemostatic valve with 3-way stopcock. Kit components are packaged in individual sterile disposable packaging of Tyvek® 1059B. The kit contents may vary depending on the physician's need in a certain component and the clinical case.

Five pcs of individual-packed kits are placed in a transport box of 200 gr. Natural Kraft E-flute carton .

INFORMATION ON LABELING OF PRODUCTS

Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

Description of the marking symbol is represented in the table.

Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

Description of the marking symbol is represented in the table.

	Manufacturer
	Authorized representative in the European Community
	Manufacturing date
	Do not reuse
	Catalogue #
	Batch Code
	Sterilized by ethylene oxide
	Caution: United States federal law restricts this device to sale by or approximately a physician
	Use by date
	Attention, Refer to the instructions for use
	Do not use, if package is damaged
	Non-pyrogenic

TRANSPORTATION, STORAGE AND HANDLING

The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures. Rotate inventory so that device is used prior to the Use By date on the package label. Products released from the transport packaging should be kept at an ambient condition:

Storage conditions:

Temperature +5°C to +27 °C (environment)

Relative humidity 10% to 90 % (non-condensing)

Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Transport conditions:

Temperature -35°C to +55 °C (environment)

Relative humidity 10% to 90 % (non-condensing)

Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Usage environment:

Temperature range for use - the body temperature of 32 °C to 44 °C

Relative humidity: 100% (blood contact)

Atmospheric pressure 86-106 kPa.

THE SHELF LIFE

Of sterile devices – 36 months (3 years)

DISPOSAL

To avoid health and environmental risks during disposal of medical devices they shall be disposed of in compliance with the hospital-approved protocol.

Medical Device " System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophyllic " according to the Sanitary Rules and Regulations 2.1.7.2790-10 is related to Class B contaminated and potentially pathological wastes (epidemiologically hazardous wastes).

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/ or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

LIMITED WARRANTY AND DISCLAIMER

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

LIMITATION OF LIABILITY

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

COMPLAINT

The organization, which takes claims from consumers about the product quality on the Russian Federation territory:

Limited Liability Company "Merit Technologies"
(LLC "Merit Technologies"),

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6
Tel / fax: +7 (495) 221-89-02

**Instructions for use for
MEDICAL DEVICE**


«System for percutaneous mini-access in kits Merit PAK, models: PAK Hydrophillic»

PRODUCT DESCRIPTION

The Merit PAK: model PAK Hydrophillic utilizes a small introducer with 3-way stopcock, dilator and guide wire for placement of larger diameter guide wires into the vasculature system when a small needle stick is preferred.

The Merit PAK: model PAK Hydrophillic consists of a combination of the following components, see unit label for specific product components. These components may be packaged in a single pouch or may be packaged separately.

- Hydrophylic introducer 4F (1.3 mm) with 3-way stopcock - 1 pcs.
- Dilator 4F (1.3 mm) – 1 pcs.
- Guide Wire 0.018" (0.46mm) - 1 pcs.
- Echo Enhanced Introducer Needle 21G (0.8 x 40 mm) -1 pcs. (if necessary)
- Echo Enhanced Introducer Needle 21G (0.8 x 70 mm) -1 pcs. (if necessary)
- Instruction for use – 1 pcs

model: PAK Hydrophillic

2. Guide Wire 0.018" (0.46mm)
3. Echo Enhanced Introducer Needle 21G (40 or 70 mm length)
7. Dilator
8. Hydrophilic Sheath introducer 4F with 3-way stopcock

CATEGORIZATION OF DEVICE(S)

Devices risk class is 2a

The categorization of medical device "System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophillic" see below.

DEVICE	CATEGORIZATION
System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophillic (introducers, dilators, guidewires, needles)	Externally Communicating Device with Circulating Blood Contact for a Limited Duration (< 24 hours). Products are intended for invasive use (use temperature from 32 ° C to 44 ° C, humidity 100% (direct contact with blood)
Hemostasis valve	Externally Communicating Device with Indirect Blood

<u>DEVICE</u>	<u>CATEGORIZATION</u>
	Contact for a Limited (< 24 hours) Duration.

INTENDED USE

System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophilic is designed to perform the "micropuncture" technique, which helps in facilitating percutaneous placement of medical devices of larger diameter (catheters, guide wires, balloons, stents) into the patient's body.

WARNINGS

Withdrawal, pull back, or manipulation of the guide wire distal tip through the needle tip may result in breakage or embolization. Do not withdraw with force. Do not advance the guide wire if resistance is met.

CAUTIONS

- Read instructions prior to use.
- This device is intended for single use only. Do not reuse or resterilize.
- This device is sterile if package is unopened or undamaged.
- This device is non-pyrogenic.
- This device should be used by clinicians with adequate training in the use of the device.
- Utilize appropriate anticoagulant therapy for patient during procedure.
- Prior to use, ensure that the sheath and dilator are the appropriate size

CONTRAINDICATIONS

Radial access is contraindicated if there is an abnormal Allen's test, radial pulse, or insufficient dual arterial supply

POSSIBLE COMPLICATIONS

The possible complications are (but not limited) air embolism, thrombus formation, infection, pseudo aneurysm formation, hematomas, guide wire embolization, bleeding, perforation or laceration of the vessel wall and standard risks connected with percutaneous diagnostic and (or) interventional procedures.

INFORMATION ABOUT POTENTIAL USERS AND USE CONDITIONS

To be use only in medical facilities by medical staff (interventional surgeons, vascular surgeons, specialists in diagnostic and surgical departments and centers) who has been specially trained how to work with this type of medical device.

Application conditions

Room temperature range	from + 15°C to + 30°C
Room humidity	< 80% (without condensation)
Atmospheric pressure	84.0-106.7 kPa (630-800mm Hg)
Room elevation	Max 1500 m a.s.l.

INSTRUCTIONS FOR USE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the device. The techniques and procedures described do not represent all medically acceptable protocols, nor are they intended as a substitute for the clinician's experience and judgment in treating any specific patient.

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
2. Remove the PAK Hydrophilic components from package using proper aseptic technique.
3. Flush all components with heparinized saline or suitable isotonic solution. Be sure to wet the outer surface of the sheath introducer to activate the hydrophilic coating. The sheath should not be used in a dry state.

Warning: After flushing side port, turn stopcock to off position to maintain flush in side port and prevent bleed back upon insertion into the vessel.

Warning: Do not wipe outer surface of the sheath introducer with dry gauze.

4. Insert vessel dilator into PAK Hydrophilic Sheath through hemostasis valve and snap into place.

Warning: Dilator must be securely snapped into place to avoid damage to the vessel.

5. Insert appropriate access needle into vessel.
 - a. If a metal access needle is used, while holding the access needle, place the flexible end or J end of the guide wire through access needle into vessel.
 - b. If an access needle with inner metal needle and outer plastic cannula is used, after appropriate access is obtained, remove the inner metal needle. While holding the plastic cannula portion of the access needle, place the flexible end or J end of the guide wire through the plastic cannula into the vessel.

Note – Refer to product labeling for appropriate guide wire compatibility with the system components.

Warning: Do not reinsert the inner metal needle into the plastic cannula at any time.

Warning: Do not advance the guide wire if resistance is met. Determine the cause of resistance before proceeding.

6. Hold guide wire in place while removing access needle. Apply manual pressure above puncture site during needle removal and until the introducer/ dilator assembly is placed.

Warning: If a needle with a metal cannula is used, do not withdraw the guidewire after it has been inserted because it may damage the guide wire.

7. Insert the introducer/dilator assembly over the guide wire into the vessel. Using a rotating motion, advance the introducer/dilator assembly through the tissue into the vessel.

Warning: Ensure that the surface of the sheath is wet prior to insertion; the sheath should not be used in a dry state.

Warning: During insertion, hold assembly near distal tip while passing over the guide wire and into the vessel to avoid buckling.

8. After introducer/dilator assembly has been placed into vessel, detach the dilator from the introducer by bending the dilator hub down slightly (this will unsnap the dilator hub from the introducer cap). While holding the sheath, carefully remove the dilator and guide wire together, leaving the sheath introducer in the vessel.

9. Aspirate from the side port extension to remove any potential air or debris. After aspiration, flush the side port with a suitable solution.

Warning: Stopcock handle must be turned to the off position (toward the sheath hub) to prevent inadvertent blood loss.

10. Use caution when inserting and removing selected device(s) (wires, catheters, etc.) into PAK Hydrophilic Sheath.

Note: Hold the sheath in place when inserting, positioning, or removing the devices. Always exchange or remove devices slowly through the sheath.

11. REMOVAL: The sheath should be removed within 24 hours. Compression on the vessel, above the puncture site, should be started as the sheath is slowly removed. Non-occlusive compression should be used to achieve hemostasis once the sheath is removed.

Note: Collected fibrin at the tip of the sheath may be aspirated via the side arm tubing prior to removal of the sheath.

12. Discard the sheath appropriately.

MAIN PRODUCT SPECIFICATION

CHARACTERISTIC	SPECIFICATION
Coaxial couple (introducer) , hemostasis valve	
Introducer, Inner Diameter (ID)	4F(1,3 mm): 0.50 mm
Outer Diameter (OD)	4F 1,3 mm
Introducer, effective length	4F (1,3 mm): 7 cm
Dilator, Outer Diameter (OD)	4F (1,3 mm)
Dilator, Inner Diameter (ID)	4F(1,3 mm): 0.45 mm
Dilator, effective length	4F-12.5 cm
Dilator, taper length	4F (1,3 mm): 5,0 mm
Dilator, weight	≤ 5 g
Hydrophilic introducer 4F(1.3 mm) with 3-way stopcock (PAK Hydrophilic only)	≤ 15 g
Hub Tensile Test	Force to remove hub shall be greater than 2.25 lbs (10N)
Freedom from leakage from introducer	No leakage at test pressure of 300 kPa
Burst strength (Introducer and introducer to hub bond)	Minimum burst strength - 4F /5F sizes ≥ 3.37 lbs (14.99N)
Burst strength (introducer shaft to the marker tip segment bond)	Minimum burst strength - 4F/ 5F sizes ≥ 3.37 lbs (14.99N)
Introducer Hub/Cap Snap Fit Tensile force to remove cap from hub	Minimum Force > 9.0 lbs. (40.03N)
Dilator drag force, valve (PAK Hydrophilic only)	Merit Medical Specification: ≤ 2,0 oz (0,057 kg)
Introducer tip dilator drag (PAK Hydrophilic only)	Merit Medical Specification: ≤ 30 oz (0,850 kg)
Suture ring tear force (PAK Hydrophilic only)	Merit Medical Specification: ≥ 3,37 lb (14,99 N)

CHARACTERISTIC	SPECIFICATION
Introducer wet modulus (PAK Hydrophillic only)	Merit Medical Specification: 4F: 100,000-152,984 psi (7 031 - 10 756 kg/cm ²)
Introducer kink extension (PAK Hydrophillic only)	Merit Medical Specification: 4F: ≥ 14.23 mm
Introducer kink recovery angle (PAK Hydrophillic only)	Merit Medical Specification: at least 90 degrees
Dilator dry flex modulus (PAK Hydrophillic only)	Merit Medical Specification: 4F: 6.7 ± 0.50 gs (0,066 $\pm$ 0,005 H)
Coating lubricity (PAK Hydrophillic only)	Merit Medical Specification: ≤ 38.2 gf (0,375 H) average with 500 g grip force (4,903 H)
Coating coverage/adherence (PAK Hydrophillic only)	Merit Medical Specification: 3 uncoated areas ≥ 5 mm ² or 1 uncoated area ≥ 15 mm ²
Introducer coated length (PAK Hydrophillic only)	Merit Medical Specification: Distance between the end of strain relief and coating $\leq 0.25''$ (0.64 cm)
Introducer tip insertion force (PAK Hydrophillic only)	4F: $\mu \leq 1,489.9$ g (with a STDEV ≤ 334 g)(3,275 H))
Sidearm Pull Test	Minimum Peak Break Force ≥ 3.37 lbs (14.99N)
Sidearm Leak Pressure Test	Sidearm shall not leak or detach at ≤ 90 psi (6.328 kg/cm ²)
Sidearm flow rate	> 600 ml/min at 200mm Hg with water
Hub to dilator bond strength	Minimum Burst strength: 4F/ 5F sizes ≥ 2.25 lbs
Dilator to Introducer Hub Snap Fit	4 F Minimum Force: 4.75 (21,13H) 5 F Minimum Force: 4.5 (20,02 H)
Bond Strength – Dilator (Hub to Rod)	Minimum Burst Strength 4F/5F > 2.25 lbs (10,01 H)
Dilator to Introducer Hub Snap Fit	Minimum Force: 4.5 lbs (20,02)
GuideWire 0,018" (0.46 mm)	
Outside Diameter	0.46 mm (0,018")
Product Lengths	40 cm
Minimum Breaking Strength (Pull Strength)	10 N / 2.25 lbs minimum
The strength of the connection rod and spiral (tensile strength)	≥ 0.5 lbs (2.2 H). Minimum tensile force
Wight (in a case)	≤ 5 g
Echo Enhanced Introducer Needle 21G (0.8 mm)	
The outer diameter of the needle	0,80 mm
Inner diameter	0,547 mm
Length of needle tube:	70 mm 40 mm
Tensile strength at break structure "needle-hub" (mounting strength)	No less than 44N
Bending force	15N (maximum deviation 0.50 mm)
Type of sharpening	tri-hedral
Allowable angle deviations needle without compromising the integrity of	short cut at an angle of $17 \pm 2^\circ$

CHARACTERISTIC	SPECIFICATION
needle roughness	No more than 0,4 mcm
The radius of the working portion of blunting needle	No more than 0,03 mm
Weight	≤ 2 g

MATERIAL SPECIFICATIONS

Component	Materials
Coaxial introducer/dilator	
Dilatator (hub)	HDPE Colorant: Red
Dilatator (tube)	Polypropylene Barium Sulphate Colorant: RED
Introducer (hub) 4F	ETFE
Introducer (Tube) 4F	ETFE <i>Hydrophilling coating</i>
Suture Ring	Petrothene Colorant: RED
Strain relief	Thermoplastic elastomer
Valve cup	Polycarbonate Color Red
Valve (body)	Elastosil
Lubricant	Medical silicone
Sidearm tube	Carbothane Cyanoacrylic adhesive
3-way stopcock	<i>Stopcock body</i> Polycarbonate <i>Stopcock arm</i> Acetal Colorant RED
3-way stopcock	ABS
Echo Enhanced Introducer Needle 21G (0.8 mm) length 40 mm, 70 mm.	
Needle Tube:	Stainless steel
Needle Base	Polycarbonate Printing Ink Color: black
Protective tube:	LDPE
Nithinol Guidewire 0,018" (0,46 mm)	
Nithinol Guidewire 0,018" (0,46 mm)	Nithinol /Pd
Package materials	
Package	Tyvek® 1073 / EVA/Surlin/EVA

STERILIZATION METHOD

Ethylene Oxide (EO) sterilization. This product does not provide for a user re-sterilization.

DESCRIPTION OF PACKAGING WITH INDICATION OF NUMBER OF PRODUCTS AND ACCESSORIES IN PACKING

The PAK Hydrophilic kit includes a 4F hydrophilic introducer with a 3-way stopcock, a 0.018" (0.46mm) guidewire, a 21 G introducer needle and a dilator. Kit components are packaged in



individual sterile disposable packaging of high density polyethylene (HDPE) Tyvek® 1073B. The kit contents may vary depending on the physician's need in a certain component and the clinical case.

Five pcs of individual-packed kits are placed in a transport box of 200 gr. Natural Kraft E-flute carton

INFORMATION ON LABELING OF PRODUCTS

Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

Description of the marking symbol is represented in the table.

	Manufacturer
	Authorized representative in the European Community
	Manufacturing date
	Do not reuse
	Catalogue #
	Batch Code
	Sterilized by ethylene oxide
	Caution: United States federal law restricts this device to sale by or approximately a physician
	Use by date
	Attention, Refer to the instructions for use
	Do not use, if package is damaged
	Non-pyrogenic

TRANSPORTATION, STORAGE AND HANDLING

The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures.

Rotate inventory so that device is used prior to the Use By date on the package label. Products released from the transport packaging should be kept at an ambient condition:

Storage conditions:

Temperature +5°C to +27 °C (environment)
Relative humidity 10% to 90 % (non-condensing)
Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Transport conditions:

Temperature -35°C to +55 °C (environment)
Relative humidity 10% to 90 % (non-condensing)
Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Usage environment:

Temperature range for use - the body temperature of 32 °C to 44 °C
Relative humidity: 100% (blood contact)
Atmospheric pressure 86-106 kPa.

THE SHELF LIFE

Of sterile devices – 36 months (3 years)

DISPOSAL

To avoid health and environmental risks during disposal of medical devices they shall be disposed of in compliance with the hospital-approved protocol.

Medical Device " System for percutaneous mini-access in kits Merit PAK, models: PAK, PAK Stiffened, PAK Hydrophillic " according to the Sanitary Rules and Regulations 2.1.7.2790-10 is related to Class B contaminated and potentially pathological wastes (epidemiologically hazardous wastes).

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/ or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

LIMITED WARRANTY AND DISCLAIMER

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with

any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

LIMITATION OF LIABILITY

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

COMPLAINT

The organization, which takes claims from consumers about the product quality on the Russian Federation territory:

Limited Liability Company "Merit Technologies"
(LLC "Merit Technologies"),

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6
Tel / fax: +7 (495) 221-89-02

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