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NOTARIAL CERTIFICATE

LAW FIRM IDEA & NOTARY OFFICE INC.

(Cheongjin-dong, Samsong Bldg.) 701, 702,
34, Jong-ro 3-gil, Jongno-gu, Seoul, Republic of Korea

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DECLARATION

We (Declarer) : Sewon Medix Inc.

do hereby solemnly and sincerely declare.


1. That the attached document(s) :

Instructions for use for dental implant system

■ is (are) true and correct.

DEC 16, 2020

Jin Soo Lee/ CEO

Sewon Medix Inc.

369-16, Deokpo-dong, Sasang-gu, Busan, Korea

Instructions for Use for Dental Implant System

English Version 4.1 / Date : November 2020

Dental Implant system consist from:

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- мини (mini): 4x7.5mm, 4x9.5mm, 4x11.5mm, 4.5x7.5mm, 4.5x9.5mm, 4.5x11.5mm ;
- стандарт (regular): 4x8.4mm, 4x10.4mm, 4x12.4mm, 4.5x8.4mm, 4.5x10.4mm, 4.5x12.4mm, 5x8.4mm, 5x10.4mm, 5x12.4mm, 5.5x8.4mm, 5.5x10.4mm, 5.5x12.4mm, 6x8.4mm, 6x10.4mm, 6x12.4mm, 6.5x8.4mm, 6.5x10.4mm, 6.5x12.4mm ; 9x8.4mm, 9x10.4mm, 9x12.4mm ;
- стандарт S с закруглением (S Curve Regular): 4x6.5mm, 4x7.5mm, 4x8.5mm, 4x9.5mm, 4x10.5mm, 4x11.5mm, 4.5x6.5mm, 4.5x7.5mm, 4.5x8.5mm, 4.5x9.5mm, 4.5x10.5mm, 4.5x11.5mm ;
4x7.4mm, 4x8.4mm, 4x9.4mm, 4x10.4mm, 4x11.4mm, 4x12.4mm, 4.5x7.4mm, 4.5x8.4mm, 4.5x9.4mm, 4.5x10.4mm, 4.5x11.4mm, 4.5x12.4mm, 5x7.4mm, 5x8.4mm, 5x9.4mm, 5x10.4mm, 5x11.4mm, 5x12.4mm, 5.5x7.4mm, 5.5x8.4mm, 5.5x9.4mm, 5.5x10.4mm, 5.5x11.4mm, 5.5x12.4mm, 6x7.4mm, 6x8.4mm, 6x9.4mm, 6x10.4mm, 6x11.4mm, 6x12.4mm, 6.5x7.4mm, 6.5x8.4mm, 6.5x9.4mm, 6.5x10.4mm, 6.5x11.4mm, 6.5x12.4mm, 9x8.4mm, 9x9.4mm, 9x10.4mm, 9x11.4mm, 9x12.4mm ;

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- мини тип шестигранник (mini Hex type) (D x G/H x H) : 4.5x1x5.5mm, 4.5x1.5x5.5mm, 4.5x2x5.5mm, 4.5x2.5x5.5mm, 4.5x3x5.5mm, 4.5x3.5x5.5mm, 4.5x4x5.5mm, 4.5x4.5x5.5mm, 4.5x5x5.5mm, 4.5x1x7mm, 4.5x1.5x7mm, 4.5x2x7mm, 4.5x2.5x7mm, 4.5x3x7mm, 4.5x3.5x7mm, 4.5x4x7mm, 4.5x4.5x7mm, 4.5x5x7mm ;
- стандарт тип шестигранник (regular Hex type) (D x G/H x H) : 4.5x1x4mm, 4.5x1.5x4mm, 4.5x2x4mm, 4.5x2.5x4mm, 4.5x3x4mm, 4.5x3.5x4mm, 4.5x4x4mm, 4.5x4.5x4mm, 4.5x5x4mm, 4.5x1x5.5mm, 4.5x1.5x5.5mm, 4.5x2x5.5mm, 4.5x2.5x5.5mm, 4.5x3x5.5mm, 4.5x3.5x5.5mm, 4.5x4x5.5mm, 4.5x4.5x5.5mm, 4.5x5x5.5mm, 4.5x1x7mm, 4.5x1.5x7mm, 4.5x2x7mm, 4.5x2.5x7mm, 4.5x3x7mm, 4.5x3.5x7mm, 4.5x4x7mm, 4.5x4.5x7mm, 4.5x5x7mm, 5x1x4mm, 5x1.5x4mm, 5x2x4mm, 5x2.5x4mm, 5x3x4mm, 5x3.5x4mm, 5x4x4mm, 5x4.5x4mm, 5x5x4mm, 5x1x5.5mm, 5x1.5x5.5mm, 5x2x5.5mm, 5x2.5x5.5mm, 5x3x5.5mm, 5x3.5x5.5mm, 5x4x5.5mm, 5x4.5x5.5mm, 5x5x5.5mm, 5x1x7mm, 5x1.5x7mm, 5x2x7mm, 5x2.5x7mm, 5x3x7mm, 5x3.5x7mm, 5x4x7mm, 5x4.5x7mm, 5x5x7mm ;

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

5.5x1x4mm, 5.5x1.5x4mm, 5.5x2x4mm, 5.5x2.5x4mm, 5.5x3x4mm, 5.5x3.5x4mm, 5.5x4x4mm, 5.5x4.5x4mm, 5.5x5x4mm, 5.5x1x5.5mm, 5.5x1.5x5.5mm, 5.5x2x5.5mm, 5.5x2.5x5.5mm, 5.5x3x5.5mm, 5.5x3.5x5.5mm, 5.5x4x5.5mm, 5.5x4.5x5.5mm, 5.5x5x5.5mm, 5.5x1x7mm, 5.5x1.5x7mm, 5.5x2x7mm, 5.5x2.5x7mm, 5.5x3x7mm, 5.5x3.5x7mm, 5.5x4x7mm, 5.5x4.5x7mm, 5.5x5x7mm, 6x1x4mm, 6x1.5x4mm, 6x2x4mm, 6x2.5x4mm, 6x3x4mm, 6x3.5x4mm, 6x4x4mm, 6x4.5x4mm, 6x5x4mm, 6x1x5.5mm, 6x1.5x5.5mm, 6x2x5.5mm, 6x2.5x5.5mm, 6x3x5.5mm, 6x3.5x5.5mm, 6x4x5.5mm, 6x4.5x5.5mm, 6x5x5.5mm, 6x1x7mm, 6x1.5x7mm, 6x2x7mm, 6x2.5x7mm, 6x3x7mm, 6x3.5x7mm, 6x4x7mm, 6x4.5x7mm, 6x5x7mm, 6.5x1x4mm, 6.5x1.5x4mm, 6.5x2x4mm, 6.5x2.5x4mm, 6.5x3x4mm, 6.5x3.5x4mm, 6.5x4x4mm, 6.5x4.5x4mm, 6.5x5x4mm, 6.5x1x5.5mm, 6.5x1.5x5.5mm, 6.5x2x5.5mm, 6.5x2.5x5.5mm, 6.5x3x5.5mm, 6.5x3.5x5.5mm, 6.5x4x5.5mm, 6.5x4.5x5.5mm, 6.5x5x5.5mm, 6.5x1x7mm, 6.5x1.5x7mm, 6.5x2x7mm, 6.5x2.5x7mm, 6.5x3x7mm, 6.5x3.5x7mm, 6.5x4x7mm, 6.5x4.5x7mm, 6.5x5x7mm ;

- Стандарт тип шестигранное основание (regular Hex-Guide type) (D x G/H x H) : 4.5x1x4mm, 4.5x1.5x4mm, 4.5x2x4mm, 4.5x2.5x4mm, 4.5x3x4mm, 4.5x3.5x4mm, 4.5x4x4mm, 4.5x4.5x4mm, 4.5x5x4mm, 4.5x1x5.5mm, 4.5x1.5x5.5mm, 4.5x2x5.5mm, 4.5x2.5x5.5mm, 4.5x3x5.5mm, 4.5x3.5x5.5mm, 4.5x4x5.5mm, 4.5x4.5x5.5mm, 4.5x5x5.5mm, 4.5x1x7mm, 4.5x1.5x7mm, 4.5x2x7mm, 4.5x2.5x7mm, 4.5x3x7mm, 4.5x3.5x7mm, 4.5x4x7mm, 4.5x4.5x7mm, 4.5x5x7mm, 5x1x4mm, 5x1.5x4mm, 5x2x4mm, 5x2.5x4mm, 5x3x4mm, 5x3.5x4mm, 5x4x4mm, 5x4.5x4mm, 5x5x4mm, 5x1x5.5mm, 5x1.5x5.5mm, 5x2x5.5mm, 5x2.5x5.5mm, 5x3x5.5mm, 5x3.5x5.5mm, 5x4x5.5mm, 5x4.5x5.5mm, 5x5x5.5mm, 5x1x7mm, 5x1.5x7mm, 5x2x7mm, 5x2.5x7mm, 5x3x7mm, 5x3.5x7mm, 5x4x7mm, 5x4.5x7mm, 5x5x7mm, 5.5x1x4mm, 5.5x1.5x4mm, 5.5x2x4mm, 5.5x2.5x4mm, 5.5x3x4mm, 5.5x3.5x4mm, 5.5x4x4mm, 5.5x4.5x4mm, 5.5x5x4mm, 5.5x1x5.5mm, 5.5x1.5x5.5mm, 5.5x2x5.5mm, 5.5x2.5x5.5mm, 5.5x3x5.5mm, 5.5x3.5x5.5mm, 5.5x4x5.5mm, 5.5x4.5x5.5mm, 5.5x5x5.5mm, 5.5x1x7mm, 5.5x1.5x7mm, 5.5x2x7mm, 5.5x2.5x7mm, 5.5x3x7mm, 5.5x3.5x7mm, 5.5x4x7mm, 5.5x4.5x7mm, 5.5x5x7mm, 6x1x4mm, 6x1.5x4mm, 6x2x4mm, 6x2.5x4mm, 6x3x4mm, 6x3.5x4mm, 6x4x4mm, 6x4.5x4mm, 6x5x4mm, 6x1x5.5mm, 6x1.5x5.5mm, 6x2x5.5mm, 6x2.5x5.5mm, 6x3x5.5mm, 6x3.5x5.5mm, 6x4x5.5mm, 6x4.5x5.5mm, 6x5x5.5mm, 6x1x7mm, 6x1.5x7mm, 6x2x7mm, 6x2.5x7mm, 6x3x7mm, 6x3.5x7mm, 6x4x7mm, 6x4.5x7mm, 6x5x7mm, 6.5x1x4mm, 6.5x1.5x4mm, 6.5x2x4mm, 6.5x2.5x4mm, 6.5x3x4mm, 6.5x3.5x4mm, 6.5x4x4mm, 6.5x4.5x4mm, 6.5x5x4mm, 6.5x1x5.5mm, 6.5x1.5x5.5mm, 6.5x2x5.5mm, 6.5x2.5x5.5mm, 6.5x3x5.5mm, 6.5x3.5x5.5mm, 6.5x4x5.5mm, 6.5x4.5x5.5mm, 6.5x5x5.5mm, 6.5x1x7mm, 6.5x1.5x7mm, 6.5x2x7mm, 6.5x2.5x7mm, 6.5x3x7mm, 6.5x3.5x7mm, 6.5x4x7mm, 6.5x4.5x7mm, 6.5x5x7mm ;

- мини тип не шестигранник (mini Non-Hex type) (D x G/H x H) : 4.5x1x5.5mm, 4.5x1.5x5.5mm, 4.5x2x5.5mm, 4.5x2.5x5.5mm, 4.5x3x5.5mm, 4.5x3.5x5.5mm, 4.5x4x5.5mm, 4.5x4.5x5.5mm, 4.5x5x5.5mm, 4.5x1x7mm, 4.5x1.5x7mm, 4.5x2x7mm, 4.5x2.5x7mm, 4.5x3x7mm, 4.5x3.5x7mm, 4.5x4x7mm, 4.5x4.5x7mm, 4.5x5x7mm ;

- стандарт тип не шесигранник (regular Non-Hex type) (D x G/H x H) : 4.5x1x4mm, 4.5x1.5x4mm, 4.5x2x4mm, 4.5x2.5x4mm, 4.5x3x4mm, 4.5x3.5x4mm, 4.5x4x4mm, 4.5x4.5x4mm, 4.5x5x4mm, 4.5x1x5.5mm, 4.5x1.5x5.5mm, 4.5x2x5.5mm, 4.5x2.5x5.5mm, 4.5x3x5.5mm, 4.5x3.5x5.5mm, 4.5x4x5.5mm, 4.5x4.5x5.5mm, 4.5x5x5.5mm, 4.5x1x7mm, 4.5x1.5x7mm, 4.5x2x7mm, 4.5x2.5x7mm, 4.5x3x7mm, 4.5x3.5x7mm, 4.5x4x7mm, 4.5x4.5x7mm, 4.5x5x7mm, 5x1x4mm, 5x1.5x4mm, 5x2x4mm, 5x2.5x4mm, 5x3x4mm, 5x3.5x4mm, 5x4x4mm, 5x4.5x4mm, 5x5x4mm, 5x1x5.5mm, 5x1.5x5.5mm, 5x2x5.5mm, 5x2.5x5.5mm, 5x3x5.5mm, 5x3.5x5.5mm, 5x4x5.5mm, 5x4.5x5.5mm, 5x5x5.5mm, 5x1x7mm, 5x1.5x7mm, 5x2x7mm, 5x2.5x7mm, 5x3x7mm, 5x3.5x7mm, 5x4x7mm, 5x4.5x7mm, 5x5x7mm, 5.5x1x4mm, 5.5x1.5x4mm, 5.5x2x4mm, 5.5x2.5x4mm, 5.5x3x4mm, 5.5x3.5x4mm, 5.5x4x4mm, 5.5x4.5x4mm, 5.5x5x4mm, 5.5x1x5.5mm, 5.5x1.5x5.5mm, 5.5x2x5.5mm, 5.5x2.5x5.5mm, 5.5x3x5.5mm, 5.5x3.5x5.5mm, 5.5x4x5.5mm, 5.5x4.5x5.5mm, 5.5x5x5.5mm, 5.5x1x7mm, 5.5x1.5x7mm, 5.5x2x7mm, 5.5x2.5x7mm, 5.5x3x7mm, 5.5x3.5x7mm, 5.5x4x7mm, 5.5x4.5x7mm, 5.5x5x7mm, 6x1x4mm, 6x1.5x4mm, 6x2x4mm, 6x2.5x4mm, 6x3x4mm, 6x3.5x4mm, 6x4x4mm, 6x4.5x4mm, 6x5x4mm, 6x1x5.5mm, 6x1.5x5.5mm, 6x2x5.5mm, 6x2.5x5.5mm, 6x3x5.5mm, 6x3.5x5.5mm, 6x4x5.5mm, 6x4.5x5.5mm, 6x5x5.5mm, 6x1x7mm, 6x1.5x7mm,

6x2x7mm, 6x2.5x7mm, 6x3x7mm, 6x3.5x7mm, 6x4x7mm, 6x4.5x7mm, 6x5x7mm, 6.5x1x4mm, 6.5x1.5x4mm, 6.5x2x4mm, 6.5x2.5x4mm, 6.5x3x4mm, 6.5x3.5x4mm, 6.5x4x4mm, 6.5x4.5x4mm, 6.5x5x4mm, 6.5x1x5.5mm, 6.5x1.5x5.5mm, 6.5x2x5.5mm, 6.5x2.5x5.5mm, 6.5x3x5.5mm, 6.5x3.5x5.5mm, 6.5x4x5.5mm, 6.5x4.5x5.5mm, 6.5x5x5.5mm, 6.5x1x7mm, 6.5x1.5x7mm, 6.5x2x7mm, 6.5x2.5x7mm, 6.5x3x7mm, 6.5x3.5x7mm, 6.5x4x7mm, 6.5x4.5x7mm, 6.5x5x7mm ;

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Information regarding the use of Fixture, Prosthetic components, Instruments

Purpose:

Dental Implant System designed for use in dental implant surgery. To provide a single or multi-element restoration design.

• Description of the system

All SewonMedix instruments are intended to be used for dental implant procedures and optimized for the SewonMedix implant system. They should not be used along with other similar products.

General Description

IH2 Fixture is dental implant made of titanium metal intended to be surgically placed in the bone of the upper or lower jaw arches. Fixture is made of pure titanium metal and supplied sterile. The surfaces are SLA (Sandblasted, Large grit and Acid etched) and RBM (Resorbable Blast Media) treated.

IH Prosthetic System is device made of titanium, titanium alloy, CCM alloy and POM intended for use as an aid in prosthetic restoration. It consists of Abutments (Healing Abutment, Solid Abutment, Cement Abutment, Angled Abutment, Temporary Abutment, Multi-unit Abutment and FreeMilling Abutment), Solid Protect Cap, Abutment Screws, Healing Cap, Cylinders (CCM Cylinder, Temporary Cylinder, Plastic Cylinder) and Cylinder Screws. Its surfaces are partially Tin coated or uncoated. IH Prosthetic System is supplied non-sterile and should be sterilized before use. It has a various diameters, gingival heights, heights and angles. IH Prosthetic System is device made of titanium, titanium alloy, CCM alloy and POM. Angled Abutment, Temporary Abutment, Healing Cap, Temporary Cylinder and FreeMilling Abutment are made of titanium. Healing Abutment, Solid Abutment, Cement Abutment, Multi-unit Abutment,

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

Abutment Screw and Cylinder Screw are made of titanium alloy. CCM Cylinder is made of Cobalt-Chrome-Molybdenum alloy. Solid Protect Cap and Plastic Cylinder are made of POM.

• Description of the system

SEWONMEDIX products are dental implants and materials are as follows.

Classification	Material
Fixture	Titanium
Abutment	Titanium, CCM(Cr-Co-Mo)
Prosthetic material and surgical instruments	Titanium, Stainless steel, Polymer

See our product catalog/home page (www.sewonmedix.com) for detailed information on each product.

For product code, specification, date of manufacture, and expiry date, see the label attached to the product.

Indication for Use

The SEWONMEDIX IMPLANT SYSTEM is designed for dental implant surgery; it is placed on the maxillary or mandibular alveolar bone through a surgical operation, and after osseointegration with the alveolar bone, it can replace a lost tooth by connecting the abutment post. The SEWONMEDIX IMPLANT SYSTEM is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including: cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. The abutment is intended for use with a dental implant fixture to provide support for prosthetic restorations such as crown, bridge, or overdenture.

1. Selection criteria/indications

Analysis of local and systemic contraindications, "Normal" wound healing capacity, effective oral hygiene/healthy remaining teeth, growth of maxilla and mandible complete, good general medical condition, adequate supply of healthy jaw bone.

2. Local examination

Anatomy of the alveolar ridge, intermaxillary relationships, deep bite, quality and thickness of the mucosa, study models and bite registration on the articulator, radiographic findings.

• Use of the SEWONMEDIX system

1. Preparations before use

- ① The dentist should be fully aware of the product and the procedure and should explain the limits of implants to the patient prior to use. And the patient should clearly understand the implant's functionality and its aesthetic limits.
- ② Selecting and fixing a proper implant is an important factor affecting the lifespan of the implant, so it is required to strictly follow indications, contraindications, precautions and warnings.
- ③ Thoroughly examine the condition of the patient before procedure and prepare a procedural planning for the patient. Failure to do it may lead to implant loss.
- ④ Read X-ray pictures well and make the right choice of the product for the patient.

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

- ⑤ Check the packaged condition of the product and expiry date and also check the product for damage.
- ⑥ The product is supplied "sterile". Do not use the product if the packaging has been previously opened or broken.
- ⑦ Contamination of a surgical instrument may lead to implant loss, so check if surgical instruments are kept in a clean condition.

2. Directions for use

Please use the product according to general implant procedure guides.

<Fixture>

- ① Take an X-ray picture and choose an implant site.
- ② Make an incision in the gum with a scalpel at the chosen site.
- ③ If the bone is exposed, position the implant at the implant site using a guide drill.
- ④ Use an initial drill to drill an initial hole to fit the length of the fixture at the site. The depth of the drilled hole can be measured by a depth gauge.
- ⑤ After the Initial Drill is used, select a middle drill and a final drill based on the diameter of the fixture to conduct drilling.
- ⑥ The dentist should select the final drill to fit the diameter of the fixture to be implanted based on the patient's bone quality, and should conduct drilling with proper rotational speed maintained.
- ⑦ Check if the implant inside the ampoule is properly connected with the Machine Fixture Driver suitable for each product. After ensuring the fit, remove the implant carefully from the ampoule and implant it in a selected position by rotating in a clockwise direction. At this time, please make sure that implantation torque above 50 Ncm should not be applied.
- ⑧ After implantation is completed, connect the cover screw to the implant with less than 15Ncm torque and then suture the surgical site to protect the implant.
- ⑨ After the implantation, the patient should avoid drinking alcohol and smoking and should maintain good oral hygiene in order to prevent contamination and infection.
- ⑩ Allow three to six month healing time depending on the patient, and then the dental practitioner should check for osseointegration and decide to conduct the prosthetic procedure.

<Abutment>

- ① Take an impression after soft tissue is healed. (It is possible immediately after the healing abutment removal in the case of One Stage Surgery.)
- ② Make the right type and size of abutment by comparing the radiograph with the taken impression and by considering adjoining teeth.
- ③ Non-sterile products should be sterilized in an autoclave at 132 °C for at least 15 minutes, and then allow a drying time of more than 15 minutes.
- ④ Connect the chosen abutment to the fixture in the mouth.
- ⑤ Take an impression at either fixture level or at abutment level in order to manufacture the final prosthesis.

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

- ⑥ In the case of cement-retained prosthesis, use a hard finishing plaster to do a wax-up on the abutment, and manufacture the prosthesis.
- ⑦ Attach the final prosthesis to the implant or the abutment in the mouth and connect it to the fixture and then check the occlusion between the upper and lower jaws.
- ⑧ Recommended tightening torque in tightening abutment screws is as follows.

Classification	Recommended torque
Cylinder Screw	20N·cm
Abutment Screw (Mini)	20N·cm
Abutment Screw (Regular)	30N·cm
Temporary Abutment	20N·cm

Modification or Customization of Abutment

The Abutment may be modified if required due to the patient’s anatomy. And IH Temporary Abutment and IH FreeMilling Abutment are designed especially for patient-specific customization.

IH Temporary Abutment is temporarily used until final prosthesis is made to maintain esthetic appreciation. The minimum post height of the IH Temporary Abutment is 4mm. The post diameter, wall thickness of the IH Temporary Abutment cannot be modified. There is no angular correction for the IH Temporary Abutment.

The minimum post height of the IH Free Milling Abutment is 4mm and the minimum post diameter is 0.4mm. The wall thickness can be modified up to 0.2mm. There is no angular correction for the IH FreeMilling Abutment.


IH Temporary Abutment and IH FreeMilling Abutment are modified by hand-milling.

Multi-unit loaded Restorations

Non-Hex types of all IH Abutments are intended for multi-unit loaded restorations. Solid Abutment, Cement Abutment (Non-Hex type), IH Angled Abutment(Non-Hex type), IH Temporary Abutment(Non-Hex type), Multi-unit Abutment and FreeMilling Abutment(Non-Hex type) are intended for multi-unit restorations.

IH Multi-unit Abutment can be used for single-unit loaded restoration when it is combined with Multi-unit Cylinder. The minimum post height of the IH Multi-unit Abutment combined with Multi-unit Cylinder is to be at least 4mm. (No angular correction)

• Precautions for use

1. Contraindications

Placement of dental implant may be contraindicated based on patient's medical condition.

The dental implant should not be used for the following patients:

- ① Patients with a history of heart failure or artery hardening;
- ② Uncooperative patients or patients with mental or physical disorder which may result in unstable implant, failure to fix the implant or complications in management after operation;
- ③ Patients with the problem in bone remodeling, microcirculation or blood clotting function;
- ④ Pregnant patients; and
- ⑤ Patients with high blood pressure or diabetes mellitus.
- ⑥ Patients with bone disease (osteoporosis, osteomalacia) or Bone metabolism disorders.

2. Procedural precautions of surgery

- ① When the dental practitioner places the final prosthesis in the patient's mouth, he or she must check osseointegration of the implanted fixture by using radiographs, the stability testing machine and percussion test.
- ② Torque within the specified torque range should be applied in tightening abutment screws. If a screw is tightened beyond the torque range, it may cause damage to the driver and the function part and also cause screw fracture. Also, tightening torque less than the torque range may lead to loosening. Thus, recommended tightening torque described in Directions for Use section must be followed.
- ③ In connecting the screw to the fixture, special attention should be paid not to apply excessive force until tightening of the first screw thread is confirmed when the screw is exactly in line with tightening direction. If the tightening of the first screw thread is not smooth, it may lead to damage to the fixture and the screw thread of the abutment, causing product damage or tightening failure.
- ④ Prosthesis should be designed in such a way that masticatory force is not applied excessively leaning to one side. If not, fracture of prosthesis due to stress concentration may be caused.

Instruments:

Purpose

All SewonMedix instruments are intended to be used for the surgical procedure of the SewonMedix implant system

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

procedure. The Drill is used to secure a space for implant placement. Cortical Drill is used to remove cortical bone to prevent over-torque in placing the fixture in hard bone. The Driver is used for fixture placement and prosthesis setting. Also, other instruments are used as aid tools in the surgical procedure.

• Using procedure

Use the instruments as described in general implant surgical procedure. For detailed information, refer to the product catalog or our company's homepage (www.sewonmedix.com).

Contraindications

If bone resorption is underway to such a degree that a proper width and height of alveolar bone cannot be secured around the implant after implant placement, the SewonMedix implant system should not be used. Lack of osseointegration or implant failure can occur for patients with insufficient bone, poor bone quality, excessive smoking habit, blood disorder, or uncontrolled diabetes.

3. Warnings

- ① Improper technique in either implant placement or restoration can result in implant failure and substantial loss of surrounding bone. So, implants should be used by trained and experienced dentists.
- ② Implants must not be used again and they should be applied based on their original purposes.
- ③ If the implant is damaged or improperly handled, it should be removed.
- ④ If the implant is improperly chosen or if the implant site is improper and fixed condition is unstable, the lifespan of the implant may be reduced.
- ⑤ Defected products should be returned.
- ⑥ Special attention should be paid in handling implants to prevent implant damage or deformation.

For safe and effective use of the SewonMedix implant system, we strongly recommend that dental implant procedures be performed by a person who has completed professional training courses because such procedures are highly specialized and complicated. Improper patient selection and technique may result in implant failure or bone loss. The SewonMedix implant system should be used, following only the described instructions for use. Any arbitrary change in using methods is not accepted for dental implants. Do not use an electro-surgical instrument or laser around the fixture or abutment because an electro shock may occur. Implant mobility, bone resorption and chronic infection mean implant failure. If an implant is infected by a patient's body fluid for one cause or another, the implant cannot be used for another patient.

4. Procedural Precautions of Restoration

Allow sufficient healing period for osseointegration to occur after implantation. The dental practitioner should closely monitor patients and decide to move to the next step. Excessive pressure including chewing pressure should not be applied to the dental fixture during the healing period.

5. Potential unfavorable consequences

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

Potential unfavorable consequences associated with the use of dental implants may include:

- ① Dehiscence of bone, delayed healing, swelling;
- ② Hematoma, infection, and inflammation; and
- ③ General hypersensitivity etc.

6. Side Effect

- ① Implant surgical procedure may accompany with risk such as swelling of certain parts, tear and opening of tissues, temporary heightened sensations, edema, embolism, bleeding etc.
- ② Loss of sensation on lower lip, surgical side effects on jaw during mandibular surgery, and nose and surrounding tissue during maxillary surgery may result from the operations. Most of symptoms are expected to be temporary though permanent paralysis may be resulted.
- ③ Necrosis of mucosa of gingiva (gum tissue), and infections on cell tissue may occur. However, these are common with localized treatment.

7. Adverse Effect

- ① Low stability during early stage, failure of osseointegration, and inappropriate prosthodontics may occur even after the surgery.
- ② Low quantity and quality of reaming bone, poor oral hygienic by patients, lack of cooperation by patients, and metabolic diseases such as diabetes may cause the reason for failure in some patients.

• Handling & Storage

The product is intended only for single use. It must not be re-used.

Keep the product in a dry place at room temperature. Keep away from direct sunlight.

*Cleaning & Sterilization Information

- IH Prosthetic System

The prosthetic components supplied non-sterile should be performed following user sterilization method provided by SewonMedix Inc. Remove the prosthetic components from packaging prior to sterilization. We recommend the product to be sterilized with moist heat sterilization using the gravity displacement autoclave for 15 minutes at temperatures 132 °C (269.6 °F) and to be dried for 15 minutes. FDA-cleared accessories (e.g., pouch, wrap, tray) are to be used during the sterilization of the prosthetic components. The recommended sterilization parameters listed above have been validated to a SAL of 10⁻⁶ in accordance with a FDA-recognized standards such as ISO 17665 using bioburden sterilization method. After using the surgical instruments, select a suitable method of cleaning that removes all visible contamination from the product and distilled water. All non-sterile components of this system **MUST** be sterilized prior to use.

Information regarding the use of **OVERDENTURE COMPONENTS**

Information regarding the use accessories of the attachment for dentures.

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

• Device description and indications for use

The LT/BT Plastic Sockets are retentive elastic components for the construction of over-dentures attachments.

The ANALOG, IMPRESSION, BLOCK OUT SPACER line is composed by prefabricated accessories for the correct use of the attachment for dentures. The Ti-CAP is prefabricated for the construction of removable prostheses.

• Contraindications

At the present state of knowledge there are no contraindications to the use of the LT/BT Plastic socket, Ti-CAP, ANALOGS, IMPRESSIONS, BLOCK OUT SPACER; at the date of emission of these instructions for use no undesirable side effects have been reported.

• Caution

The choice of the attachment is demanded to the dentist or to the dental technician according to the prosthetic project.

• Material

All the materials used for the LT/BT PLASTIC SOCKET, Ti-CAP, ANALOG, IMPRESSION COPPING, BLOCK OUT SPACER are of quality and meet the international standards of reference. The patient has to inform the dentist about potential allergies. If allergic reactions occur search the source and remove it.

Classification	Material
Plastic socket	Polyamide, POM
Ti-Cap	Titanium
Prosthetic material and surgical instruments	Stainless steel, POM, Titanium

• Warnings and single use

The LT/BT PLASTIC SOCKET, Ti-CAP, ANALOG, IMPRESSION COPPING, BLOCK OUT SPACER are single-use products; their reuse undermines the functionality and may affect the patient's health.

Improper use or reuse of LT / BT PLASTIC SOCKET, Ti-CAP, ANALOG, IMPRESSION COPPING, PROTECTIVE DISK are prohibited.

• Washing and sterilization

Dentist is required to pack and store the prefabricated metal parts in cleared FDA 510(k) pouches (paper/plastic film) suitable for sterilizing and have to sterilize them before use with the following parameters: steam sterilization with prevacuum cycles at 135°C-3minutes, drying time 16minutes, pressure 2.06 bar. It is recommended to sanitize LT/BT PLASTIC SOCKET, Ti-CAP, ANALOG, IMPRESSION COPPING, BLOCK OUT SPACER with benzalkonium chloride solutions before use. The dentist is responsible of the sterilization.

• Maintenance and periodic care

Dentists have the responsibility to keep the proper functionality and retention of the attachments assuring the safety of

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

the patient by constant maintenance. In order to maintain the high quality standard offered by the . LT/BT PLASTIC SOCKET, Ti-CAP, ANALOG, IMPRESSION COPPING, BLOCK OUT SPACER we recommend the substitution of the retentive elastic components yearly.

• Guidelines for the patients

Patients are recommended to follow the indications provided by the dentist, to attend periodical controls and perform daily accurate hygiene.

Safety and compatibility has not been evaluated about Ti-Cap in MR environments. Heating, migration, or image artifacts have not been tested in the MR environment. Safety has not known in MR environment.

Checking up the patient with this device may result in injury to the patient.

• Technical specs

Allow to manage the degree of retention of the prosthesis. Available in various retentive forces distinguishable by different colors.

LT Plastic Socket Pin type		LT Plastic Socket No-Pin type		BT Plastic Socket(Ø2.2 Ball type)	
Plastic socket	Retention force	Plastic socket	Retention force	Plastic socket	Retention force
Violet cap	1800gr./17.65N	Violet cap	1500gr./14.71N	Clear cap	Normal size:1300gr.
Clear cap	1500gr./14.71N	Clear cap	1000gr./9.8N	Pink cap	Normal size:900gr.
Pink cap	900gr./8.83N	Pink cap	800gr./7.84N	Yellow cap	Normal size:500gr.
Yellow cap	700gr./6.86N	Yellow cap	600gr./5.88N	Green cap	Normal size:350gr.
Black cap	No retention(only for laboratory use)	Black cap	No retention(only for laboratory use)		

- ANALOG:

stainless steel device to reproduce the position and the shape of the attachment into the laboratory master model.

- BLOCK OUT SPACERS:

DISPOSABLE plastic devices. used to protect the gum and to eliminate the undercuts during the fixation, with self-curing resin, of the metallic Ti-Cap + Plastic socket in the patient's mouth.

- Ti-Cap

prefabricated metal device designed to house the PLASTIC SOCKET; available in titanium or anodized titanium. The outer profile is designed to be cured directly in the prosthesis or to be cemented or welded into a cast reinforcement. The cementation must be done with anaerobic self-curing composite cement for metal to metal bonding.

• Instructions for use

- ① The LT/BT PLASTIC SOCKET has to be inserted inside a prefabricated Ti-CAP embodied the resin denture.
- ② To insert or remove the PLASTIC SOCKET inside their Ti-CAP without any damage, use the appropriate PLASTIC SOCKET TOOL.
- ③ Before fixing the Ti-CAP inside the denture resin, put the BLOCK OUT SPACER around the neck of the attachment, connect the Ti-CAP with the PLASTIC SOCKET on it and create enough room inside the denture; in correspondence of the Ti-CAP and check if there is enough space for the Ti-CAP; once done, fix it with self curing resin.

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax +82 51 3031714

- ④ To avoid the loss of the performances of the PLASTIC SOCKET replace them every 12 months.
- ⑤ Any use of the "PLASTIC SOCKET" which does not follow the present instructions is considered improper.

- Analog

place the correct ANALOG into the corresponding IMPRESSION COPING, make sure the two elements are correctly coupled and proceed with the production of the master model according to the technique chosen in the laboratory.

- Ti-Cap

Prefabricated devices that can be welded or cemented into cast reinforcements or cured directly in the denture resin.

•Precautions for use for instruments

- The hand-piece should be moved vertically up and down during bone drilling.(pumping).
- Saline solution must be applied during bone drilling in order to reduce frictional heat generated from bone drilling.
- Laser marking is a dimension indication line, and the bottom of laser marking is the standard. For laser marking of all drills, the 11.5mm and 13mm levels are marked in bold line.
- The drill and the driver are recommended to be used for 50 times respectively. The maximum allowable torque is 20Ncm for mini size abutment and 30Ncm for regular size abutment.
- If the drill or the driver is used exceeding the recommend number of use, implant failure may occur due to loss of product function caused by abrasion and fracture of the driver, or due to bone heating caused by declined cutting force of drill. Please, follow the recommendations.
- The dental practitioner should properly determine the bone condition of the patient and apply appropriate cutting speed.
- The length of drill stopper means the remained actual length of twist drill when the drill stopper is installed on drill.
- Small instruments should be fixed with a wire before use because they may be moved into the throat during the procedure.
- Please, immediately detach and clean all the instruments used for the procedure. Store them at the room temperature after drying them.
- Do not keep instruments in a place where the risk of infection exists.
- This product is a medical device and must be sterilized in an autoclave before use (at the temperature of 132°C, for 15 minutes, and drying for 30 minutes).
- Do not use hydrogen peroxide, which is generally used as a disinfecting and cleaning agent, since TiN coating, laser marking, and color coding may be damaged or discolored.
- The warranty period of this product is one year from the date of purchase.

•Methods of cleaning and sterilization

Instruments must be cleaned and sterilized before re-use. Please, observe the following instructions.

<Cleaning>

- Remove the contaminants (skin tissue, blood, bone chips etc.) remaining on the instrument used for the procedure.
- Submerge the instrument for 10 minutes in alkaline cleaning solution or other similar cleaning solution.
- Then, remove all visible contaminants with a soft brush.

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713

Fax. +82 51 3031714

- Use distilled water to completely wash the remaining cleaning solution on the surface of the instrument.
- Place the instrument into the cleaning container in an ultrasonic cleaner and fill it with water.
- Perform ultrasonic cleaning at the room temperature for 10 minutes.
- Wipe the washed instrument with a non-abrasive cloth and then keep it in a dry place.

<Sterilization>

- Sterilize non-sterile instruments in an autoclave before use.
- Sterilize the instruments in an autoclave (for 15 minutes at 132 °C) and then dry for 30 minutes.
- Store the sterilized instruments in a clean environment between 10°C and 30°C.

• **Temperature of use, transport and storage**

Do not damage the packaging during the transportation; store in a clean, dry place, in the original packaging away from sunlight and heat sources.

Temperature condition of storage : 1° ~ 30°

Expiration date for ;

- Implants : 5 year
- Prosthetic components : No expiration date
- Instruemnsts : No expiration date

• **Disposal**

Follow the local laws for the disposal of the medical devices. It is suggested to wash and sterilize the products before which disposal for all products. After the use the device has to be disposed as a contaminated healthcare waste.

• **Warranty**

All SEWONMEDIX Inc. products are guaranteed to be free from defects in material and workmanship at the time of shipping. All of our products are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals, for any failure due to misuse or application which is not in accordance with the designer's intentions. Products will lose their right to carry the CE sign in case they will be changed or repaired at another place than or a company authorized by SEWONMEDIX Inc. to do the change or repair.

Head office: #29, Sasang-ro, 375beon-gil, Sasang-gu, Busan, 617-814, Korea

Telephone : 82. 51. 303. 1713 Fax : 82. 51. 303. 1714

SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713


Fax. +82 51 3031714

[Pictograph for Fixture and Abutment]

	Do not reuse		CAUTION
	Batch Code		Consult Instructions for Use
	Manufacturer		EC Representative
	Do not resterilize		Temperature Limitation
	Sterilized using irradiation		Use-by date
	Catalogue number		Date of manufacture
	Do not use if package is damaged		Non-sterile
	For USA Only: Federal law restricts this device to sale by or the order of a dentist		

- The non-sterile symbol is applied to only Abutment, which is end user sterilization.
- The sterilized using irradiation symbol is applied for only Fixture.

[Pictograph] for instruments

	Batch Code		CAUTION
	Manufacturer		Consult Instructions for Use
	Temperature Limitation		EC Representative
	Catalogue number		Do not use if package is damaged
	Non-sterile		


SEWONMEDIX Inc.

29, Sasang-ro 375beon-gil, Sasang-Gu, Busan, Republic of Korea

Tel. +82 51 3031713


Fax. +82 51 3031714

[Pictograph] for components

	Batch Code		CAUTION
	Manufacturer		Consult Instructions for Use
	Temperature Limitation		EC Representative
	Catalogue number		Do not use if package is damaged
	Non-sterile		Do not reuse

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appeared before me and admitted
said principal`s subscription to
the attached DECLARATION. ----

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This is hereby attested on this
17th day of Dec. 2020 at this office.

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(Cheongjin-dong, Samsong Bldg.) 701, 702,
34, Jong-ro 3-gil, Jongno-gu, Seoul,
Republic of Korea

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.243.

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