ACKNOWLEDGMENT

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California
County of Contra Costa

On 5 June 2015 before me, Margarita Maria Cazares, Notary Public
(insert name and title of the officer)

personally appeared Jennifer Yonemura
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are
subscribed to the within instrument and acknowledged to me that he/she/they executed the same in
his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the
person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing
paragraph is true and correct.

WITNESS my hand and official seal.

Signature


Margarita Maria Cazares

(Seal)

«I CERTIFY ACCURACY, CORRECTNESS AND
RELIABILITY OF THIS DOCUMENT TEXT»

Jennifer Yonemura
Regulatory Affairs Director
Cerus Corporation, USA

Concord, CA USA
05 June 2015

Instruction for use

INTERCEPT Platelet Processing Set with Dual Storage Containers

ver. 2015-06

INTERCEPT Platelet Processing Set with Dual Storage Containers

INDICATION FOR USE

INTERCEPT Illuminator CERUS (device INTERCEPT Illuminator).

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria, and parasites as well as contaminating donor leukocytes in platelet components. Platelets processed using the INTERCEPT Set (INTERCEPT Platelets), are indicated for support of patients requiring platelet transfusions, according to clinical practice guidelines. INTERCEPT Platelets may be stored up to 7 days from time of collection, at 20-24°C with continuous agitation. INTERCEPT Platelets stored up to 7 days have been shown to adequately prevent and control bleeding. Any extension of platelet storage time should be evaluated according to applicable local policies and regulations.

Platelet additive solutions approved for use with INTERCEPT: InterSol, SSP+

SCOPE OF USE

Transfusion medicine.

COMPOSITION

Each INTERCEPT set includes one 17.5mL 3mM amotosalen hydrochloride solution (formula : amotosalen hydrochloride 101mg, sodium chloride 924mg, water for injection q.s.100ml), one illumination container, one container with compound adsorption device, two INTERCEPT platelet storage containers.

6 INTERCEPT Sets are wrapped in an aluminium foil package with tamper-evident. INTERCEPT Sets are sterilized by a combination of steam and radiation.

CONTRAINDICATIONS

Use of INTERCEPT Platelets is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

CONDITIONS OF USE

For healthcare professionals only. The correct application of INTERCEPT set requires the appropriate conditions of specialized blood services or transfusion departments and healthcare professionals properly trained.

PRECAUTIONS

Do not use if: aluminium foil package with tamper-evident has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container.

INTERCEPT Set is single use only.

All the following conditions must be met for pathogen inactivation:

- Platelets must be prepared with additive solution for platelet re-suspension in the volume range specified in Table 1.
- Platelet content, volume, and red blood cell count must be within the ranges specified in Table 1.
- Platelets mixed with amotosalen must be exposed to UVA light dose from the INTERCEPT Illuminator. No other source of UVA light may be used.
- Platelets collected on Day 0 must be exposed to UVA light by end of Day 1.
- After illumination, platelets must be agitated in CAD container in accordance with the durations specified in Table 2.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, platelet product must be discarded.

Warnings: Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device should be used according to manufacturer's instructions for use.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses, bacteria and parasites, there is no pathogen inactivation process that has been shown to eliminate all pathogens.

Neonatal patients who require platelet transfusion during phototherapy for treatment of hyperbilirubinemia should not be treated with phototherapy devices that emit light less than 425nm to avoid the theoretical potentiation of erythema resulting from interaction between UVA light and psoralen.

POSSIBLE ADVERSE EFFECTS

The adverse platelet transfusion reactions may be the consequence of immune reactions.

- Febrile non-hemolytic transfusion reactions (FNHTRs)
- Hemolytic Transfusion Reactions
- Allergic and anaphylactic reactions. Allergic reactions are very rare and highly variable in severity. Manifestations may occur as isolated pruritus and urticaria as the only dermal manifestations. Systemic reactions may include bronchoconstriction, hypotensive reactions and shock.
- Other Immunologically Mediated Reactions Caused by Transfusion of Platelets

The healthcare personnel should be aware of potential adverse effects. A suspicion of an adverse reaction should be reported and be appropriately treated.

EQUIPMENT AND MATERIALS

Materials Needed: One (1) INTERCEPT Platelet Processing Set with Dual Storage Containers.

Equipment Needed: INTERCEPT Illuminator device, Sterile Connecting Device (SCD), Tube Sealer, Flatbed Agitator.

Temperatures should be controlled to meet applicable regulations for platelet processing.

Used and unused INTERCEPT sets should be discarded like any used blood containers, as biohazardous waste.

INSTRUCTIONS FOR USE

A- Preparation of Platelets

Platelets should be prepared with additive solution. Platelets are suspended in approximately 35% plasma and 65% additive solution (see Table 1). Platelet products within the following ranges have been shown to be acceptable for use with this processing set.

Table 1

Platelet Content	Volume	Plasma content	Additive Solution Content	RBC Count
2.5 – 7.0 x 10 ¹¹	300-420ml	32-47%	53-68%	<4 x 10 ⁶ RBC/mL
7.1 – 8.0 x 10 ¹¹	375-420ml			

B- Amotosalen Addition to Platelets

1. Remove set from package. Unwrap only illumination container from organizer.
2. Weld tubing from platelet container to amotosalen container tubing using Sterile Connecting Device.

3. Label INTERCEPT set containers using appropriate donation identification. See warning in Section C. While labeling storage containers, separate them to ensure they do not adhere to one another.
 4. Hang platelets and break both cannulae on amotosalen container.
 5. Allow platelets and amotosalen solution to flow into illumination container, marked by number "1".
 6. Ensure that platelets are completely transferred to illumination container by expressing air from platelet and amotosalen container into illumination container.
 7. When platelets are in illumination container, mix thoroughly by gentle agitation to ensure complete mixing of amotosalen and platelets.
 8. Express air from platelets into amotosalen container.
 9. Express a small amount of platelet and amotosalen mixture into tubing, filling about 4cm of tubing.
 10. Seal tubing between illumination container and amotosalen container, so that tubing is no longer than approximately 4cm from illumination container inlet port.
- Warning:** During illumination, tubing must be held within large compartment of illumination tray.
11. Remove and discard empty platelet and amotosalen containers.

C- Illumination

Illuminate platelets.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use.

Warning: Platelets in illumination container must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with platelets. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

D- Agitation with CAD

Warning: Do not fold or bend Compound Adsorption Device.

1. Unwrap set from organizer.
2. Hang platelets, break cannula and allow platelets to flow into Compound Adsorption Device container, marked by number "2".
3. Express air from Compound Adsorption Device container into illumination container.
4. Seal tubing close to inlet port of Compound Adsorption Device container.
5. Remove and discard empty illumination container.
6. Place Compound Adsorption Device container on flatbed agitator for the duration specified in Table 2.

Table 2

Suspension Medium	Compound Adsorption Device Agitation Duration
Additive Solution (53–68%)	6 to 16 hours

E- Transfer to Storage Containers

1. Remove platelets from agitator and hang platelets.
2. Close clamp on platelet sampling pouch.
3. Break cannula and allow platelets to flow into each of the two storage containers, marked by number "3."
4. Express air from storage containers into CAD container and clamp above the Y fitting.
5. Redistribute volume between the storage containers as needed. Note: Steps 4 and 5 may be performed in the reverse order.
6. Ensure appropriate donor identification is attached to each filled platelet storage container.
7. Seal tubing to each storage container close to outlet of the Y fitting.

8. Remove and discard empty CAD container. Place platelets on platelet agitator. The INTERCEPT Platelet process is now complete.

F- Sampling Platelet Product (optional)

1. Mix INTERCEPT platelets by gently agitating platelet storage container.
2. Open clamp to platelet sampling pouch and squeeze several times.
3. Allow sampling pouch to fill with platelets. Seal tubing.
4. Remove sample pouch.
5. Transfer sample to appropriate laboratory tube immediately.

CAUTION:

Diethylhexylphthalate (DEHP) is known to be released from polyvinyl chloride (PVC) medical devices; increased leaching can occur with extended storage or increased surface area contact. The INTERCEPT sets only have tubing components, container ports and, if included, an in-line filter that contain PVC; all containers and other parts are PVC-free. During use of INTERCEPT set, blood components are in contact with PVC for a brief period of time (approx. <15 minutes). Based on limited surface area contact and minimal contact time, DEHP levels in blood components after use of INTERCEPT set are estimated to be well below those resulting from other medical applications containing PVC tubing (e.g. hemodialysis, intravenous fluid administration, extracorporeal membrane oxygenation and cardiopulmonary bypass procedures). The risks associated with DEHP released to the blood components must be weighed against the benefits of therapeutic transfusion and inactivation of harmful viruses, bacteria and other pathogens.

TRANSPORTATION AND STORAGE CONDITIONS

Do not store above 25°C.

Do not vent. Do not freeze.

Protect the pack and tubing from sharp objects.

Unused INTERCEPT sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 8 hours.

Keep set in light-protective (aluminium) package until time of use. Protect from direct sunlight and strong UVA light source.

Transportation of the device should be in accordance with the regulatory requirements.

Device should be stored in accordance with the regulatory requirements.

SHELF LIFE

2 years.

Do not use after the expiry date printed on the aluminium package.

DISPOSAL

Used and unused or expired INTERCEPT Processing sets should be disposed like any used blood containers, as biohazardous waste.

ASSURANCE

This device was designed, developed, manufactured, examined and packed in accordance with the requirements of the applicable national and international standards and regulations.

Assurance of "Cerus Corporation" applies to the device which is compliance with the instructions such as not expired, unopened, undamaged and correctly kept.

Explanation of symbols

Do not use if the product, its sterile barrier system or its packaging is damaged or shows any sign of deterioration.

Do not freeze. Do not store above 25°C.

DS

Dual Storage Containers

Platelet Storage Container

Refer to Instructions for Use

Sterile fluid path. Sterilized by a combination of steam and radiation.

Non pyrogenic fluid path.

Protect from direct sunlight and strong UVA light source.

Do not vent.

Contains or presence of phthalates: Di-(2-ethylhexyl)phthalate (DEHP).

Do not re-use.

Apply label here.

CE – mark. Product conforms with the essential requirements in the MDD: 93/42/EEC

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InterSol – a trademark of Fenwal, Inc.
SSP+ – a trademark of MacoPharma