All-Russian

Classification of Products code 2


32.50.50.190

**OPERATING MANUAL
for the medical product**

"Full-Thickness Resection device of upper GI system set",
Manufactured by Ovesco Endoscopy AG, Germany

Meets standard international and Russian requirements

December 02, 2020

Christoph Damm, Dipl.-Ing. (FH)

Head of Quality Management & Regulatory Affairs

Ovesco Endoscopy AG

Friedrich-Miescher-Strasse 9
72076 Tübingen
Germany

ooo|ovesco

Ovesco Endoscopy AG
Friedrich-Miescher-Straße 9
D-72076 Tübingen
www.ovesco.com

Table of contents

1. Product indication
2. Supply kit
3. General application and operation details
4. Indications for application
5. Contraindications
6. Possible complications
7. Precaution measures
8. Product description
9. Operation conditions and rules
10. Product materials
11. Product technical specifications
12. Marking
13. Production diagram
14. Package
15. Transportation
16. Storage
17. Technical maintenance
18. Minor repair
19. Disposal details
20. List of the national and international standards
21. International certificates and conformity certificates data sheet
22. Warranty certificate

1. PRODUCT PURPOSE

It is intended for application during the procedure of full-thickness resection of mucosa of upper GI.

2. SUPPLY KIT

Full-thickness resection device system set of upper GI composed of:	1 pcs
1. FTRD system	1 pcs
1.1 FTRD clip	1 pcs
1.2 The thread	RZN 2017/5662 dated 01.04.2020
1.3 Application cap FTRD	1 pcs
1.4 The ring on the application cap FTRD	RZN 2017/5662 dated 01.04.2020
1.5 hand wheel FTRD	RZN 2017/5662 dated 01.04.2020
1.6 PE endoscopic sleeve	RZN 2017/5662 dated 01.04.2020
1.7 Plastic tape	RZN 2017/5662 dated 01.04.2020
1.8 FTRD snare	RZN 2017/5662 dated 01.04.2020
1.9 thread retriever	RZN 2017/5662 dated 01.04.2020
2. FTRD coagulation probe	RZN 2017/5662 dated 01.04.2020
3. FTRD grasper forceps	RZN 2017/5662 dated 01.04.2020
4. Ballon	1 pcs (upon request)
5. Guidewire	1 pcs (upon request)
6. Coagulation dissector	1 pcs (upon request))
7. IFU	1 pcs

(hereinafter – “FTRD set”).

3. GENERAL APPLICATION AND OPERATION DETAILS

Product application conditions – surgery and endoscopy departments of the medical treatment facilities.

FTRD set is a single-use product.

FTRD set is a sterile product (EO sterilization).

Depending on the potential application risk, FTRD set belongs to the 2b class (according to the Order of the Ministry of Public Health of the Russian Federation №4Н dd. 06.06.2012).

The way of protection according to GOST 14254-96 – IPX0.

Type of the medical devices according to classification of the medical device is 328560. Isolation strength should meets the requirement of i.8.5.2.1 GOST 60601-1-2010 for the devices that work from the power 220V, 50 Hz

4. INDICATIONS FOR APPLICATION

FTRD set is applied for application during the procedure of full-thickness resection of mucosa.

5. CONTRAINDICATIONS

FTRD set shall not be applied in case any contraindications exist for treatment by flexible endoscopy methods.

- FTRD set must not be used near a stenosis.
 - FTRD set must not be used in case of adhesions of the colorectum with adjacent organs.
- FTRD set must not be used if it is clear that the tissue to be resected cannot be completely excised and the resection would cut through diseased tissue.
- FTRD SET must not be used in conjunction with cardiac stimulators and any other active implants .

6. POSSIBLE COMPLICATIONS

- Haemorrhages in the area of the procedure
- Damage to/perforation of the colon or rectum wall
- Damage to adjacent organs
- The FTRD set may not provide sufficient haemostasis

- The bowel lumen may narrow in the area of intervention.
- The application site of the FTRD System may become infected

7. PRECAUTION MEASURES

WARNINGS

Check components for damages and other defects prior to application of FTRD set.

Defective parts can cause incorrect device operation or patient's injuries.

Incorrect application can cause injuries for the patient: for example, a clip can be released or can fall from the application cap accidentally.

The electrical burst power must exceed or equal the maximum peak output voltage with which the instrument is operated in combination with an

appropriate HF generator at an appropriate operating mode/setting.

These instructions for use must be carefully read before use. The user should have fully understood how the instrument operates and how it should be

handled.

All components of the FTRD set are for single-use only. Materials used in the product are not intended for the processing and re-sterilization. Device components are designed for single-use for one patient and are not intended for multiple use.

Processing and reusing of the single-use devices can cause harm to the patient (inadequate functioning, insufficient sterility, infection, corrosion and other quality depreciation types). Prior to application, please, ensure that the package is secure and check the expiry date. It is prohibited to guarantee products sterility in case the package is damaged or the shelf life has expired.

MRT-information

FTRD clips are MRT-compatible. Patients treated by this device can be safety scanned immediately after placement under the following conditions:

- Static magnetic field ≤ 3 tesla;
- Magnetic field spatial gradient ≤ 720 gauss/cm.

Heating connected with MRT

FTRD clips produced heating for less than $1,7^{\circ}\text{C}$ at the MRT system maximum, considering that average specific absorption rate (SAR) of the whole body was equal to 3 W/kg for 15 min in 3 Tesla (3-Tesla/128MHz, Radiating, Software G3.0-052B, General Electric Healthcare, Milwaukee, Wisconsin) of MRT system.

In the result, MRT-related test for FTRD clips heating under 3 Tesla with the use of a transmit-receive radio-frequency coil on the MRT system provided the result of the average specific absorption rate (SAR) equal to 3 W/kg (for example, associated with the value of 2.8 W/kg measured in the calorimeter), and demonstrated the maximum heating level which became apparent in connection with these specific conditions at the level of $\leq +1,7^{\circ}\text{C}$.

Information about distortions

Quality of the MRT image can be put under doubt in case the zone of concern is located directly or closely to the point of FTRD clips location. Due to this, MRT image parameters optimization is required for compensation of this device presence.

Impulse mode	T1-SE	T1-SE	GRE	GRE
Size of the	543mm ²	112mm ²	984mm ²	402mm ²

signal-free

zone

Planar

Parallel

Perpendicular


Parallel

Perpendicular

orientation

8. PRODUCT DESCRIPTION

FTRD set:

1. FTRD hand wheel (RZN 2017/5662 dated 01.04.2020)

By turning the FTRD hand wheel, the thread is tensioned, the white ring at the distal cap end is pulled and the FTRD clip is applied. To simplify the procedural step, it is indicated by a 1.

2. FTRD applicator cap

The FTRD applicator cap is mounted on the tip of the endoscope. The FTRD applicator cap is fitted with a clip and a HF snare.


An applicator ring sits on the FTRD applicator cap below the clip. The thread is secured to the FTRD applicator cap and the applicator ring, and is used to release the clip.

The FTRD applicator cap is designed for use with standard flexible endoscopes with an outer diameter of 11.5–13.2 mm. The thread is secured to the applicator ring and cap and must not be removed.

3. FTRD snare (RZN 2017/5662 dated 01.04.2020)


The FTRD snare is a monopolar high-frequency instrument (HF). The FTRD snare is made of stainless steel and is integrated into the distal end of the FTRD applicator cap. The flexible shaft is guided externally along the endoscope. The FTRD snare consists of:

- Sliding handle (A) to close the snare
- HF connection, 3-mm plug (B)

- Flexible shaft (C)
- Snare safety lock to prevent the snare from being accidentally released (D). To simplify the procedural step, this is indicated by a 2.
- Snare that is integrated in the FTRD applicator cap (E).

The electrical burst power must exceed or equal the maximum peak output voltage with which the instrument is operated in combination with an appropriate HF generator at an appropriate operating mode/setting.


4. FTRD clip

The FTRD clip consists of a super-elastic shape memory alloy (Nitinol®) and is designed to grasp large volumes of tissue. Thanks to its material properties, following application into the tissue, the clip returns to its original, slack shape and therefore compresses the tissue between the clip branches.

The FTRD clip material – Nitinol® – is MR-conditional.

In comparison with the OTSC clip (14t), the FTRD clip has additional teeth around the edge to optimise tissue fixation in the border area.

FTRD clips are placed under mechanical pressure on the application cap.

Compression zone width – clipping width (a)

The table below shows clips and application cap sizes

Clipping width [a]	Maximum outer diameter of the application cap [b]
11 mm	23 mm

Size of the clipped tissue varies depending on the tissue size drawn into the application cap. Therefore, clipped tissue zone size can exceed the clipping width.


5. PE Endoscope sleeve (RZN 2017/5662 dated 01.04.2020)

Once the FTRD applicator cap is fitted, the endoscope cover, which is already positioned on the flexible shaft of the FTRD snare, is secured to the FTRD applicator cap using the adhesive tape supplied and is unrolled to the proximal end of the endoscope. The cover is designed to protect the tissue when inserting the FTRD System into the body, and to better guide the snare externally along the endoscope insertion tube.

6. Thread retriever (RZN 2017/5662 dated 01.04.2020)

The thread retriever supplied is an aid used in the assembly of the FTRD System. It is used to pull the thread through the working channel.

The FTRD hand wheel is inserted into the rubber-rim valve on the endoscope's working channel and secured to the handle using the Velcro strap (4).


The two cable ties that are fixed on the snare shaft are removed. The endoscope sleeve is slid over the endoscope, starting from the distal side of the FTRD applicator cap (5). To secure the FTRD applicator cap on the distal end of the endoscope, the thread that is attached to the FTRD applicator cap is pulled through the endoscope's working channel using the FTRD thread retriever (6-A) or a similar instrument (6-B). The knot on the thread is clamped in the designated slot on the hand wheel (7). The thread is wound up by turning the hand wheel clockwise (7).


The FTRD applicator cap must be positioned at the tip of the endoscope so that the thread passes through the working channel vertically and does not conceal any of the endoscope functions (8). Ensure that the thread is not wound around the snare shaft. The distal end of the endoscope must abut against the stopper in the FTRD applicator cap (18). The FTRD snare shaft is positioned in the recess on the silicone sleeve (9). The figure above (8) shows the correct position of the FTRD applicator cap on the endoscope tip. The thread is further wound up using the FTRD hand wheel (7) until slight tension is felt. The endoscope sleeve is slid to the FTRD applicator cap and the cable tie on the side that is proximal to the cap is removed.

The end of the sleeve must be secured to the endoscope approximately 2–3 mm above the FTRD applicator cap using the adhesive tape supplied. The adhesive tape must not be touching the FTRD applicator cap – it should be applied just above the cap (10).


The endoscope sleeve is pulled back to the proximal end of the endoscope using the small tube on which it is positioned. In doing so, it unrolls. The remaining cable tie and small tube are removed (11).

For safety reasons, the instrument's HF connection to the HF generator should not be established until shortly before the clip is applied. (Observe the instructions for use of the HF generator.)

7. FTRD grasper (RZN 2017/5662 dated 01.04.2020)

Application aid "FTRD® Grasper" is an instrument for grasping and manipulation tissue using flexible endoscopes with a working channel diameter of minimum 2.8 mm, for use alone or minimum 3.2 mm together with the FTRD® System.

DESCRIPTION OF COMPONENTS

FTRD Grasper consists of:

- handle (A)
- sliding handle (B)
- flexible shaft, length 220 cm (C)
- jaw (D)

FTRD Grasper is a part of FTRD set and it's a sterile device.


8. FTRD marking probe (RZN 2017/5662 dated 01.04.2020)

The FTRD Marking Probe is a monopolar high-frequency instrument (HF) for marking tissue with coagulation points in flexible endoscopy. The minimum diameter of the working channel is 2.8 mm.

FTRD Marking Probe is a part of FTRD set and it's a sterile device.


9. Dilatation balloon

Dilatation balloon is intended for expand of the esophagus for providing of FTRD distal cap before using of FTRD set.


The length of the balloon is 1800 mm, diameter of expanded balloon is 20 mm, length of the balloon is 60 mm.

10. Guidewire

Nitinol guidewire intended for moving the instrument during the full-thickness resection. The length of the guidewire is 3700 mm, diameter of guidewire is 0,89 mm.

11. Coagulation dissector

Coagulation dissector intended for hemostasis during full-thickness resection procedure.

Length of the dissector: 1650 mm

Diameter of the dissector: 2,5 mm

Length of the working part: 13 mm

Length of the cups: 6,5 mm

Angle of the cups: 90 degree

9. OPERATION CONDITIONS AND RULES

9.1. Preparation for application

Prior to the product application it is required to check the components for cracks and other defects. Replace all defective parts into new one. Defective parts can cause failures or cause injuries to the patient.

Incorrect application can also cause injury to the patient or doctor, since a clip can accidentally turn around or jump off the cap.

Pay attention to the package damages and shelf life expiry date. In case the package is damaged or the shelf life has already expired, product sterility cannot be guaranteed.

All OVESCO products are designed only for the single-use application. Materials do not withstand processing and re-sterilization. Components are intended for application for one patient and do not withstand multiple use.

Recycling and repeated application of the single-use tools can cause harm to the patient (improper operation, insufficient sterility, contamination, corrosion and other destruction types).

WARNING/PRECAUTION MEASURES

- The doctor shall ensure that FTRD set are applied in compliance with the operation manual.
- FTRD clips are mounted on the FTRD Application cap with a high mechanical tension. Due to this, the risk of the injury in the result of improper treatment exists both for the doctor and the patient. Potential problems also include unintentional clip release or release in the non-target point.
- Ovesco products are compatible with each other. Prior to Ovesco items application together with the products of another manufacturers, the doctor shall ensure that the products are compatible in terms of the target application scope.

9.2 Before the application

- Use the product only as intended
- The user must ensure that the size of the lesion to be resected does not exceed the capacity of the FTRD applicator cap.

- The FTRD System is designed for use with standard flexible endoscopes with a working channel diameter of at least 3.2 mm.
- If an instrument is to be inserted through the same working channel as the FTRD System thread, ensure that the thread remains slightly tensioned. Otherwise, it may wind itself around the instrument.
- It is recommended that the lesion to be resected is marked in advance to provide better visualisation (via coagulation points using the FTRD Marking Probe).
- The FTRD snare must be connected to the monopolar outlet of a HF generator using a suitable monopolar cable and must only be used with monopolar coagulation and cutting current.
- The FTRD snare can be connected to standard HF generators. Do not use the instrument on patients with cardio stimulators or other active implants.

9.3 Application

Safety measures to avoid HF burn:

- Always use a gloves before using
- Before activating the HF generator, ensure that the working end of the product does not touch any electrically conductive accessory.
- Always keep the working end of the snare or coagulation probe in the user's field of vision during HF activation.
- Never put down the product on, or directly next to, the patient
- Switch off the HF generator's automatic activation mode when accessories are used endoscopically

Preparations:

It is recommended that the target tissue is marked with the FTRD Marking Probe, before using the FTRD System. To do it, please use a coagulation probe that is a part of FTRD set, insert the instrument to the working channel and make a coagulation of a tissue.

If the patient has a narrow esophagus you have to use a dilatation balloon. Provide the balloon through the working channel of the scope and inflate the balloon by inflation device (not included). Making a dilatation and take out the balloon. If the patient has a hemorrhage during the procedure, you should use a coagulation dissector. Provide the balloon through the working channel of the scope and connect the instrument to the electrosurgical unit. Using default settings make an endoscopic hemostasis.

After the dilatation, use a guidewire to provide the FTRD set through the guidewire for good positioning of the set.

The procedure:


After the tissue was marked by coagulation probe, the system should be mount according to item 8. During the mount, the scope should be on the vision filed, at the same time the compression is not allowed.

Before insertion, the distal end of the scope should be coated with lubricant.

Once the application site has been reached, the FTRD snare is connected to the HF generator via the connection cable.

Grasp the tissue:

The required resection site is aimed at with the FTRD applicator cap. The target tissue is grasped with the FTRD grasper. The tissue is pulled into the FTRD applicator cap by pulling the forceps into instrumental channel of the scope until the markings are within the cap.

In doing so, the user must ensure that the size of the lesion and the size of the FTRD applicator cap match and that the lesion fits wholly within the FTRD applicator cap.

Due to the increased risk of sucking in adjacent tissue or organs at the same time, the tissue must not be sucked into the FTRD applicator cap.

Depending on the anatomy, it may be more difficult to collect the tissue due to fixation.

Clip application:

By turning the FTRD hand wheel clockwise, the white ring on the applicator cap is pulled forwards and the clip is released.

In doing so, it is important that the tissue is kept under slight tension using the FTRD grasper, for example, until resecting. As a result, the tissue and the distal end of the endoscope will be held in place.


When using an instrument to manipulate tissue through the working channel, ensure that the instrument tip is sufficiently drawn into the applicator cap and remains in this position while the clip is triggered. Otherwise, there is a risk that the clip may be applied to the instrument shaft or parts of the branches, and the instrument may become fixed to the tissue as a result.

Make sure that the clip has been properly applied before resectioning begins. If the clip has been successfully applied, the white ring will have moved forwards to the edge of the cap. The white ring must therefore remain visible and be observed in the endoscopic image. Only when the white ring has moved fully forwards to the edge of the cap can you be sure that the clip has been applied. If the clip has not been properly applied, this may result in perforation when resectioning is performed.

Resection:

After the clip has been applied, the safety lock on the FTRD snare is removed. By pulling the sliding handle backwards, the snare is firmly positioned around the target tissue. HF current is then applied to the FTRD snare.

The tissue is cut as a result of tightening the snare by pulling it back fully into the shaft of the FTRD applicator cap while applying HF current. During the cutting procedure, the snare must be pulled continually. The HF generator must be set so that a sufficient tissue-cutting effect is produced. After cutting the grasped tissue, the power supply is disconnected. By removing the endoscope from the body, the resected tissue can be retrieved from the FTRD applicator cap.

Before resecting, it is important to ensure that the clip is positioned and the target tissue is securely within the FTRD applicator cap above the clip. In addition, the tissue must be securely grasped with the grasping forceps. If the tissue is not securely fixed in the FTRD applicator cap, this may result in perforation next to the clip caused by the FTRD applicator cap sliding unintentionally.

Inspection:

After applying the FTRD System, carefully inspect the application site with the endoscope. When doing so, ensure especially that the FTRD clip is correctly positioned; check for possible haemorrhages, perforations, unacceptable lumen restrictions (stenosis) or other abnormalities. These may result in further endoscopic or surgical procedures being necessary. Check the resected tissue for extraluminal structures that were also resected. If this is the case, it may result in subsequent extraluminal perforations or haemorrhages.

Disposal:

To dismantle the FTRD System, the thread is unwound from the hand wheel. The FTRD applicator cap, which is connected to the thread and the FTRD snare, is carefully removed from the endoscope, along with the endoscope sleeve.

After use, the FTRD applicator cap, FTRD snare, endoscope sleeve and FTRD hand wheel are disposed of as waste according to rules of utilisation.

10. PRODUCT MATERIALS

I. FTRD set

No.	System component name	The system	Material/trademark	Standard
1	FTRD clip	FTRD system	Niti-S, Nitinol	ASTM F2063
2	The thread	FTRD system	PE Dyneema Purity UG dtex 110 TS100	ISO 10993 – 3,4,5,6,10,11,12
3	Application cap	FTRD system	SAN LURAN 388S	EN 71-9
4	The ring on the application cap	FTRD system	Hostaform MT12U01	USP Class VI
5	Hand wheel	FTRD system	ABS-plastic (Novodur HD-15) and (Maxithen 5BA6157) Stainless steel 1.4301 Silicon rubber	USP Class VI EN ISO 7153 AISI 303 USP Class VI FDA 177.2600
6	Thread retriever	FTRD system	Shaft: Vestamid L, Polyamid 12 Hook: Ultramid, Polyamid 6	FDA 21 CFR 177.1500
7	FTRD snare	FTRD system	Terlux MABS 2812HD ABS-plastic (Novodur HD-15) and (Maxithen 5BA6157), PET, PA, Stainless steel 1.4301/1.4305/1.4310	USP Class VI
8	PE endoscopic sleeve	FTRD	Low Density	USP Class VI

		system	Polyethylen (LDPE)	FDA 177.2600
9	Plastic tape	FTRD system	Soft plastic PVC	-
10	FTRD coagulation probe	FTRD coagulation probe	PP Stamylan P 16M10 Maxithen HP 5BA6187 Stainless steel 1.4305 HDPE B5429	ISO 1133, ISO 527-1/-2, ISO 75-1/2 EN 10088-1, DIN 17449, AISI 316, UNS S31600, ISO 10993 – 3,4,5,6,10,11,12
11	FTRD Grasper forceps	FTRD Grasper forceps	Stainless steel 1.4301 ABS-plastic PTFE	EN ISO 7153, AISI 303 USP Class VI
12	Dilatation balloon	Dilatation balloon	ABS-plastic (Novodur HD-15) and (Maxithen 5BA6157) Low Density Polyethylen (LDPE) Polyamid PA	USP Class VI FDA 177.2600
13	Guidewire	Guidewire	Niti-S, Nitinol	ASTM F2063
14	Coagulation dissector	Coagulation dissector	Stainless steel 1.4301 ABS-plastic PTFE	EN ISO 7153, AISI 303 USP Class VI

11. PRODUCT TECHNICAL SPECIFICATIONS

- 11.1 Length of the thread of the FTRD system – 2200 ± 75 mm
- 11.2 Maximum outer diameter of the application cap of the FTRD system – 21 ± 2 mm
- 11.3 Length of the distal cap of FTRD system – 37 ± 2 mm
- 11.4. Length of nitinol clip – 23 ± 2 mm
- 11.5 Clipping width (width of the clip's compression zone) of the FTRD system – 13 ± 1 mm
- 11.6 Length of the teeth of nitinol clip – $3 \pm 0,1$ mm
- 11.7 Length of FTRD snare – 2200 ± 75 mm
- 11.8 Diameter of FTRD snare – $2,3 \pm 0,1$ mm
- 11.9 FTRD snare has a HF connector to connect the snare an active cord, the diameter of HF connector – $8,4 \pm 0,1$ mm
- 11.10 Length of thread retriever – $2200 \pm 2,3$ mm
- 11.11 Diameter of thread retriever – $2,3 \pm 0,1$ mm
- 11.12 FTRD snare and FTRD coagulation probes maximum electrical burst power is 50 watts.
- 11.13 Length of FTRD coagulation probe – 2200 ± 75 mm
- 11.14 Diameter of FTRD coagulation probes – $2,2 \pm 0,1$ mm
- 11.15 Length of FTRD grasper forceps – 2200 ± 75 mm
- 11.16 Diameter of FTRD grasper forceps – $2,3 \pm 0,1$ mm
- 11.17 Maximum angle of branches opening of FTRD grasper forceps– $130 \pm 5^\circ$
- 11.18 Diameter of the distal part of the scope – $11,5 - 13,2$ mm
- 11.19 Weight of the FTRD system set should be:
- FRTD clip – 40 ± 10 g
 - Thread – 10 ± 5 g
 - Distal cap FTRD - 80 ± 10 g
 - The ring on the distal cap FTRD - 20 ± 5 g
 - Hand wheel FTRD - 100 ± 10 g
 - Thread retriever - 50 ± 10 g
 - FTRD snare – 200 ± 50 g
 - FTRD coagulation probe 200 ± 50 g
 - FTRD grasper forceps 300 ± 50 g

- PE endoscopic sleeve 100 ± 50 g
- plastic tape 100 ± 50 g
- guidewire 150 ± 50 g
- dilatation balloon 150 ± 50 g
- coagulation dissector 200 ± 50 r

11.20 The length of the dilatation balloon is 1800 ± 50 mm, diameter of the expanded balloon is 20 ± 2 mm, length of the balloon is 60 ± 5 mm

11.21 Length of the guidewire is 3700 ± 75 mm, diameter of the guidewire is $0,89 \pm 0,05$ mm.

11.22 Length of the dissector is 1650 ± 50 mm, diameter of the dissector is $2,5 \pm 0,2$ mm, length of the working part is $13 \pm 0,5$ mm, length of the cup is $6,5 \pm 0,5$ mm, angle of the cup is 90 degree.

11.23 Cracks, pockets, dints, scratches, chipped places, burrs, burns, slag, particles of the rubbing and polishing materials are not acceptable.

11.24 Surface roughness parameters shall be as follows:

- for clips not more than $0.1 \mu\text{m}$
- for other devices – from 0.1 to $0.32 \mu\text{m}$

11.25 Tools shall be manufactured from corrosion resistant materials. The items shall be corrosion resistant in the operation, transportation and storage conditions.

11.26 The tools shall withstand temperature and air humidity impact in the process of transportation and storage for the conditions under group 5 (OZh4) under GOST 15150.

11.27 Set of the accessories should meet the requirements of the resistance to climatic factors of the environment according to UHL 4.2 under GOST 15150.

11.28 The tools in transport packaging should be sustainable in the process of transportation under GOST R 50444.

11.29 Surface roughness parameters shall be as follows:

- for clips not more than $0.1 \mu\text{m}$

11.30 Filament transporter must smoothly pass through the working channel of the tool without jamming.

11.31 The main requirements to electro surgery generators, that provides operability of the device:

- Electro surgery generators cut power should be at least 40 watts

- Electro surgery generators coagulation power should be at least 20 watts
- Electro surgery generators must be operable with the network 220

All electro surgery generators meets these requirements are acceptable to work with the medical device.

Recommended settings of electro surgery generators:

In coagulation mode with FTRD coagulation probe – power 20 watts.

In cutting mode with FTRD snare – power 40 watts.

11.32 The diameter of the thread of the system: $0,4 \pm 0,1$ mm

11.33 The dimensions of the ring on the distal FTRD cap:

- Outer diameter – $20,85 \pm 0,05$ mm
- internal diameter – $17,3 \pm 0,05$ mm

11.34 The dimensions of the hand wheel:

- Length – $66,3 \pm 0,5$ mm
- Width – $42 \pm 0,5$ mm
- Height – $35 \pm 0,5$ mm

11.35 The dimensions of the PE endoscopic sleeve:

- Length – 1800 ± 50 mm
- Width – 40 ± 5 mm

11.36 Sterilisation - ethylene oxide (EO)

11.37 The condition of the using of medical device:

- Temperature – $10 - 35$ °C
- Humidity – 40-100 %

12 MARKING

12.1. Product marking shall comply with the user's manual and requirements of GOST R 50444-92.

12.2. Main information on the package shall be as follows:

- name and/or trademark of the manufacturer,
- manufacturer's address and telephone,
- Description of the package content,
- inscriptions "Sterile", "Non-toxic".
- Batch (series) number,
- shelf life

Additional information application is also possible:

"Warning! Read the manual"

LOT

Batch number

REF

Reference number

Manufacture date

Ethylene-oxide sterilized (EO).

Best before "year-month"

Single-use only

Don't sterilize


Do not use in case of damaged package

Latex-free

13 PRODUCTION DIAGRAM

14. PACKAGE

14.1. Consumer package design shall provide tools hermetic tightness and protection from external mechanical damages.

14.2. FTRD set is packed into blister package placed into cardboard box. The original package provides best security during transportation.

14.3 The product can't be identified by an external scanner.

14.4 The material of packaging that is subject to final sterilization: high density polyethylene (HDPE), high strength paper (Tyvek).

14.5 Overall dimensions of the package that is subject to final sterilization:

- endoscopic system of application of clips FTRD with clips – 385 ± 5 mm x 211 ± 5 mm x 45 ± 5 mm;
- probe for coagulation FTRD - 200 ± 5 mm x 300 ± 5 mm;
- forceps for gripping FTRD - 200 ± 5 mm x 350 ± 5 mm;
- dilatation balloon - 385 ± 5 mm x 211 ± 5 mm;

- guidewire - 385 ± 5 mm x 211 ± 5 mm;

- coagulation dissector 390 ± 5 mm x 260 ± 5 mm;

14.6 The material of packaging that is subject to final sterilization should not be leachable, toxic and should not have a smell.

14.7 The material of packaging that is subject to final sterilization should not have holes, cracks, wrinkle or local thickenings and/or thinning.

14.8 Physical and chemical characteristics of the packaging that is subject to final sterilization must meet the requirements for medical devices, and must be consistent with packaging and final sterilization processes.

14.9 Components of packages should not react, contaminate, switch to products and should not adversely affect the products.

14.10 Process indicators that are printed on packages must meet the requirements ISO 11140-1.

14.11 Opening of packages by delamination should be continuous and homogeneous, without disruption of the material which may worsen its aseptic properties.

14.12 The packing density should be in the range 71 g/m^2 – 80 g/m^2 .

15. TRANSPORTATION

The product is placed into a corresponding package together with the operation documentation. Transportation by all roofed transport means is allowed at the temperature mode for the group 5 conditions (OZh4) under GOST 15150.

16. STORAGE

Keep the product in a dry place protected from adverse impact of the temperature, humidity, ventilation, sunlight, dust and air containing salt, etc.

After transportation pay attention to its consequences and possible damages.

Avoid places of chemicals storage of gas formation.

Protect from moisture penetration.

The product withstands the temperature and air humidity in the process of transportation and storage for 5th group conditions (OZh4) under GOST 15150.

Product shelf life – 3 years.

17. TECHNICAL MAINTENANCE

Single-use product, is not subject to maintenance.

18. MINOR REPAIR

The FTRD set is a single-use item and is not subject to repair.

In case of the package integrity damage or product damage, please, address manufacturer's representative:

"Endo Stars", LLC, h. 27 Bolshaya Monetnaya Street, lit. A, room 12H, St. Petersburg 197101, tel./fax: 8-800-555-56-57

19. DISPOSAL DETAILS

The product does not contain any precious metals.

The product withdrawn from the operation, is subject to the disposal.

In accordance with the Decree of the Chief state sanitary officer of the Russian dd. 9th December, 2010 N 163 "Concerning the approval of the sanitary regulations and standards", this product belongs to the class B.

Class B wastes are subject to mandatory disinfection/neutralization.

For collection of the class B wastes it is necessary to use single-use puncture-proof water-proof reservoirs (containers) of yellow color or with yellow marking. The container shall have a tightly closing cover, excluding the possibility of inadvertent opening.

20. LIST OF THE NATIONAL AND RUSSIAN STANDARDS

The list of the applied RF national standards and international standards for the purpose of safety, effectiveness and quality provision, as well as standards with the control methods: GOST R 50444-92, Directive 93/42/EEC, EN ISO 13485:2012+AC:2012, GOST ISO 11607 - 2011

Product manufacture is performed in accordance with the Directive 93/42/EEC dd. 08.11.2015, appendix 2 for the class 2B items.

21. INTERNATIONAL CERTIFICATES DATA SHEET

- Certificate EN ISO 13485:2016 dd. 20.05.2019
- EC Certificate for the Quality Assurance System according to Directive 93/42/EEC dd. 20.05.2019
- Certificate of the company's registration dd. 09.07.2019

WARRANTY CERTIFICATE

Note of Sale

Product category _____

Manufacturer _____

Seller Stamp

Serial/identification number _____

Sale date _____

Selling company _____

The product is tested and no defects found.

Warranty period – 1 year.

For warranty questions, please, refer to:

Exclusive representative on the Russian Federation territory Endo Stars LLC, legal address: h. 27 Bolshaya Monetnaya Street, lit. A, room 12H, St. Petersburg 197101, tel./fax: 8-800-555-56-57

ATTENTION! Make sure that the Warranty Certificate is filled in correctly. The Serial/identification number must correspond exactly to the one in the Certificate. All fields in the Certificate must be filled in.

In case of inaccurate or incomplete filling of the Certificate, the warranty for the product is regarded as invalid.

Manufacturer's warranty is provided in accordance with the Civil Code of the Russian Federation on the usual terms.

The product should not bear the marks of mechanical or thermal action, as well as of various chemical agents. The product use is allowed only in accordance with the operating manual.

WARRANTY CENTER NOTE

DATE	DESCRIPTION	WARRANTY CENTER
------	-------------	-----------------

Notary Daniel Westa

Schaffhausenstrasse 77 • D-72072 Tübingen
Phone +49 (0)7071 / 13870-0 • Fax +49 (0)7071 13870-99



mail@notar-westa.de
www.notar-westa.de

notarised copy

Herewith I certify that this copy is an entire copy of the original document
which was presented to me today.

city of Tübingen, December 2nd, 2020
federal state of Baden-Wuerttemberg / Germany

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