State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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1. **Country (Pays): United States of America**
This public document (Le présent acte public) Instructions for Use
2. **has been signed by (a été signé par) Ryan J Bartz**
3. **acting in the capacity of (agissant en qualité de) Notary Public, State of Minnesota**
4. **bears the seal / stamp of (est revêtu du sceau / timbre de)**
Ryan J Bartz, Notary Public, State of Minnesota

Certified Attesté

5. **at (à) St. Paul, Minnesota**
6. **Date (Date) 11/25/2020**
7. **by (par) Secretary of State, State of Minnesota**
8. **File No (Sous n°) 1195334000091**
9. **Seal / Stamp (Sceau / Timbre):**
10. **Signature (Signature):**

Steve Simon

Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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Truthful and Accurate Statement

This is to certify that the attached document is the Instructions for Use for
Unitek™ Gemini Clear Orthodontic Ceramic Brackets.

Kristin Totushek, Regulatory Affairs Manager
3M Unitek Corporation

18 Nov 2020

(Date)

State of Minnesota)
County of Ramsey) ^{SS}

Subscribed and sworn to before me this 18 day of November 2020.

Ryan J. Bartz, Notary Public

Unitek™ Gemini Clear orthodontic ceramic brackets in variants:

1. Unitek™ Gemini Clear orthodontic ceramic brackets:

1.1. Unitek™ Gemini Clear orthodontic ceramic brackets in kits of 10 pcs., identical: 117-301 (Upper Left Central), 117-302 (Upper Right Central), 117-303 (Upper Left Lateral), 117-304 (Upper Right Lateral), 117-305 (Upper Left Cuspid, with Hook), 117-306 (Upper Right Cuspid, with Hook), 117-307 (Upper Left Bicuspid, with Hook), 117-308 (Lower Anterior), 117-309 (Lower Left Cuspid, with Hook), 117-310 (Lower Right Cuspid, with Hook), 117-311 (Lower Left 1st Bicuspid, with Hook), 117-312 (Lower Right 1st Bicuspid, with Hook), 117-313 (Lower Left 2nd Bicuspid, with Hook), 117-314 (Lower Right 2nd Bicuspid, with Hook), 117-316 (Upper Right Bicuspid, with Hook);

1.2. Unitek™ Gemini Clear orthodontic ceramic brackets in kits of 12 pcs. (U/L; 3x3):

- bracket kit 117-104 consist of: 117-301 (Upper Left Central) - 1 pc., 117-302 (Upper Right Central) - 1 pc., 117-303 (Upper Left Lateral) - 1 pc., 117-304 (Upper Right Lateral) - 1 pc., 117-305 (Upper Left Cuspid, with Hook) - 1 pc., 117-306 (Upper Right Cuspid, with Hook) - 1 pc., 117-308 (Lower Anterior) - 4 pcs., 117-309 (Lower Left Cuspid, with Hook) - 1 pc., 117-310 (Lower Right Cuspid, with Hook) - 1 pc.;

1.3. Unitek™ Gemini Clear orthodontic ceramic brackets, in kits of 20 pcs. (U/L; 5x5):

- bracket kit 117-100 consist of: 117-301 (Upper Left Central) - 1 pc., 117-302 (Upper Right Central) - 1 pc., 117-303 (Upper Left Lateral) - 1 pc., 117-304 (Upper Right Lateral) - 1 pc., 117-305 (Upper Left Cuspid, with Hook) - 1 pc., 117-306 (Upper Right Cuspid, with Hook) - 1 pc., 117-307 (Upper Left Bicuspid, with Hook) - 2 pcs., 117-316 (Upper Right Bicuspid, with Hook) - 2 pcs., 117-308 (Lower Anterior) - 4 pcs., 117-309 (Lower Left Cuspid, with Hook) - 1 pc., 117-310 (Lower Right Cuspid, with Hook) - 1 pc., 117-311 (Lower Left 1st Bicuspid, with Hook) - 1 pc., 117-312 (Lower Right 1st Bicuspid, with Hook) - 1 pc., 117-313 (Lower Left 2nd Bicuspid, with Hook) - 1 pc., 117-314 (Lower Right 2nd Bicuspid, with Hook) - 1 pc.

2. Unitek™ Gemini Clear orthodontic ceramic brackets, with APC™ II adhesive:

2.1. Unitek™ Gemini Clear orthodontic ceramic brackets, with APC™ II adhesive in kits of 5 pcs., identical: 3117-301 (Upper Left Central), 3117-302 (Upper Right Central), 3117-303 (Upper Left Lateral), 3117-304 (Upper Right Lateral), 3117-305 (Upper Left Cuspid, with Hook), 3117-306 (Upper Right Cuspid, with Hook), 3117-307 (Upper Left Bicuspid, with Hook), 3117-308 (Lower Anterior), 3117-309 (Lower Left Cuspid, with Hook), 3117-310 (Lower Right Cuspid, with Hook), 3117-311 (Lower Left 1st Bicuspid, with Hook), 3117-312 (Lower Right 1st Bicuspid, with Hook), 3117-313 (Lower Left 2nd Bicuspid, with Hook), 3117-314 (Lower Right 2nd Bicuspid, with Hook), 3117-316 (Upper Right Bicuspid, with Hook);

2.2. Unitek™ Gemini Clear orthodontic ceramic brackets, with APC™ II adhesive in kits of 12 pcs. (U/L; 3x3):

- bracket kit 3117-104 consist of: 3117-301 (Upper Left Central) - 1 pc., 3117-302 (Upper Right Central) - 1 pc., 3117-303 (Upper Left Lateral) - 1 pc., 3117-304 (Upper Right Lateral) - 1 pc., 3117-305 (Upper Left Cuspid, with Hook) - 1 pc., 3117-306 (Upper Right Cuspid, with Hook) - 1 pc., 3117-308 (Lower Anterior) - 4 pcs., 3117-309 (Lower Left Cuspid, with Hook) - 1 pc., 3117-310 (Lower Right Cuspid, with Hook) - 1 pc.;

2.3. Unitek™ Gemini Clear orthodontic ceramic brackets, with APC™ II adhesive, in kits of 20

pcs. (U/L; 5x5):

- bracket kit 3117-100 consist of: 3117-301 (Upper Left Central) - 1 pc., 3117-302 (Upper Right Central) - 1 pc., 3117-303 (Upper Left Lateral) - 1 pc., 3117-304 (Upper Right Lateral) - 1 pc., 3117-305 (Upper Left Cuspid, with Hook) - 1 pc., 3117-306 (Upper Right Cuspid, with Hook) - 1 pc., 3117-307 (Upper Left Bicuspid, with Hook) - 2 pcs., 3117-316 (Upper Right Bicuspid, with Hook) - 2 pcs., 3117-308 (Lower Anterior) - 4 pcs., 3117-309 (Lower Left Cuspid, with Hook) - 1 pc., 3117-310 (Lower Right Cuspid, with Hook) - 1 pc., 3117-311 (Lower Left 1st Bicuspid, with Hook) - 1 pc., 3117-312 (Lower Right 1st Bicuspid, with Hook) - 1 pc., 3117-313 (Lower Left 2nd Bicuspid, with Hook) - 1 pc., 3117-314 (Lower Right 2nd Bicuspid, with Hook) - 1 pc.

Package contents:

Unitek™ Gemini Clear orthodontic ceramic brackets:

- Orthodontic ceramic brackets according to the variants;
- Identification card (for sets with different brackets - 117-104, 117-100) - 1 pc.

The following medical devices name can be used in this Instructions for use: Unitek™ Gemini Clear, Unitek™ Gemini Clear precoated with APC™ II, brackets, medical device, product.

REF 117-XXX, 3117-XXX

Intended Use

3M Unitek™ products are intended for use during orthodontic treatment of malocclusions (imperfect positioning of the teeth when the jaws are closed).

Potential Users

An Orthodontist is a dental specialist who is concerned with the study and treatment of malocclusions (imperfect positioning of the teeth when the jaws are closed), which may be a result of tooth irregularity, disproportionate jaw relationships or both.

Indications for use

This product is intended for use during orthodontic treatment of malocclusions (imperfect positioning of the teeth when the jaws are closed).

Contraindications

None.

Potential Adverse Events (Potential Complications)

Enamel Damage
Allergic Reaction
Patient Injury

Contamination

Accessories or products to be used in combination with device

Instruments

The following Hand Instruments are recommended to be used with the Unitek™ Gemini Clear:

Unitek™ Debonding Instrument (REF 900-850) (RC No. FSZ 2011/10227).

It is also possible to use other instruments manufactured by 3M Unitek Corporation.

These items are sold separately.

Adhesives

For Brackets without APC™ II adhesive, it is recommended to use orthodontic adhesives manufactured by 3M Unitek Corporation.

Archwires

It is recommended to use 3M Unitek Corporation archwires of the following cross-section (size) with Unitek™ Gemini Clear Brackets:

Wire archwires:

Slot 0.018:

Round - diameter up to 0.49 mm.

Square - side up to 0.49 mm.

Rectangular - sides up to 0.49 and 1.85 mm, respectively.

Slot 0.022:

Round - up to 0.59 mm in diameter.

Square - up to 0.59 mm.

Rectangular - sides up to 0.59 and 1.85 mm, respectively.

Curing lamps

For Unitek™ Gemini Clear Brackets with APC™ II adhesive, it is recommended to use the Ortholux Luminous Curing Light (RC No. FSZ 2012/13215).

Other Orthodontic Products

For Unitek™ Gemini Clear Brackets, it is recommended to use ligatures, elastic bands, intraoral and extraoral ligatures manufactured by 3M Unitek Corporation.

Compatibility with products from other companies

3M Unitek Corporation has not conducted such research.

Warnings

- Brackets are for single use only. Use of recycled brackets may compromise material integrity, bond strength, patient safety and overall product performance.

- Due to the hardness of ceramic brackets, bonding brackets in occlusion should be avoided to prevent wearing of enamel surfaces during all phases of treatment.
- Deep bite cases should be opened prior to bonding with Unitek™ Gemini Clear Ceramic Brackets to prevent wearing of enamel surfaces. Unitek™ Elastomeric Ligatures with Guard (REF 406-430, 406-429) are designed to help protect opposing tooth surfaces from occlusal interference and from potential enamel wear. Bite plates and other intrusion mechanics may also be necessary to prevent bracket to tooth contact. Indications for using Unitek™ Elastomeric Ligatures with Guard modules include deep bite cases, protection of cusp tips during retraction and situations in which tooth to bracket contact is possible during finishing.
- Instruct patients not to chew or bite on hard substances such as hard candy, ice, carrots, etc. Careful and thorough patient instruction is key to avoiding appliance or enamel damage.
- Bonding of ceramic brackets to compromised teeth (i.e., with large restorations, peg laterals or preexisting conditions) can increase the risk of tooth damage.
- Bonding to porcelain crowns or facings may cause chipping or breakage of the crown or facing during treatment or debonding.
- It is recommended to debond Unitek™ Gemini Clear Ceramic Brackets using the Unitek™ Debonding Instrument, REF 900-850.
- It is recommended that extra precautions are taken when debonding ceramic brackets to ensure that bracket parts are held securely.
- If a bracket fractures during treatment or debonding, use a diamond bur to carefully remove the ceramic parts. Failure to follow the correct debonding procedure may lead to tooth damage.

Additionally, for Unitek™ Gemini Clear coated with adhesive APC™ II.

Warning: This product contains acrylate monomers. Acrylate monomers are known to produce allergic skin reactions in certain sensitive individuals. May cause eye and skin irritation.

Precautions: Avoid eye and skin contact. Wear gloves when handling this material.

First Aid:

Eye contact: Immediately flush with plenty of water. See a physician if irritation persists.

Skin contact: Wash affected area with soap and water. See a physician if irritation persists.

Precaution: Porcelain Crowns or Facings: Care must be taken when bonding to porcelain crowns or facings as debonding may cause chipping, delamination or breakage of the crown. Do not bond to porcelain crowns that have thin facings or that appear to be compromised. Prepare the porcelain crown to be bonded by using a porcelain primer such as 3M ESPE RelyX Ceramic Primer. Follow instructions included with the porcelain primer.

Bracket Identification and Bracket Positioning

Unitek™ Gemini Clear brackets are identifiable by a color coded recessed distal gingival dot.

Bracket Bonding

Remove adhesive flash at bonding to minimize staining of the adhesive and reduce the risk of bracket breakage during the debonding procedure. Unitek™ Gemini Clear brackets feature mechanical retention on the bonding base. No special primers or pretreatments are necessary.

Unitek™ Gemini Clear brackets can be used with traditional direct or indirect bonding methods. No change in bonding technique is necessary. For bonding procedures please follow the adhesive

manufacturer's recommendations. When placing Unitek™ Gemini Clear brackets, it is suggested that the bracket be placed in a sliding motion, occlusal to gingival, forcing excess adhesive to the incisal edge of the bracket for easier clean-up. Care must be taken when cleaning up adhesive flash so as not to disturb the final positioning before adhesive curing.

Bracket Bonding for Unitek™ Gemini Clear coated with adhesive APC™ II.

A. Tooth Preparation

1. Isolate teeth using the Dry Field System or a combination of retractors, dri-angles and cotton rolls.
2. Prophyl teeth with an oil-free pumice or paste. Figure 1. Rinse with water. Isolate teeth.
3. Air dry thoroughly using oil and moisture-free air source. Figure 2.

B. Acid Etching

1. If using Transbond Plus Self Etching Primer, proceed to Option 3: Priming Surfaces with Transbond Plus Self Etching Primer steps below. Otherwise, continue with phosphoric acid etching.
2. Apply Unitek™ Etching Gel (REF 712-039 or 712-044) to tooth surfaces following the instructions provided with the etching gel delivery system. If using other gel etching systems, consult the manufacturer's instructions for proper technique and recommended etching times.
3. Rinse with water.
4. Air dry thoroughly using oil and moisture-free air source.

C. Priming Teeth

Option 1: Priming Surfaces with Transbond XT Primer

For detailed instructions on Transbond XT Primer, please refer to REF 011-519.

Option 2: Priming Surfaces with Transbond MIP Moisture Insensitive Primer

For detailed instructions on Transbond MIP Moisture Insensitive Primer, please refer to REF 011-563.

Option 3: Priming Surfaces with Transbond Plus Self Etching Primer

For detailed instructions on Transbond Plus Self Etching Primer, please refer to REF 011-581.

D. Bracket Placement and Curing

1. Appliance Removal – Determine the tooth to be bonded. Select the adhesive coated appliance and peel back the tab lid of the individual appliance blister.

Note: To protect adhesive from premature curing, DO NOT open blister packs until ready for appliance placement. Blisters cannot be re-sealed.

2. To remove the adhesive coated appliance from the blister, grasp the mesial/distal sides of the appliance with a bracket placement instrument and pull quickly in an upward and straight motion. Figure 4.

Note: The appliance has been pre-oriented in the blister to facilitate placement. Figure 5.

3. Lightly place appliance onto tooth surface, adjust to final position and then press firmly.

Note: If final positioning will be delayed, cover patient's mouth with a dark colored mask or other article to avoid premature curing by ambient light.

4. Gently remove excess adhesive around appliance base without disturbing the appliance.


5. Hold the curing light perpendicular to ceramic appliance surface.

Precaution: Please follow manufacturer's directions for the handling, proper use and recommendations for eye protection when using a curing light.

Refer to the chart to determine the curing conditions for optimal bond strength. If your curing light does not appear on this chart, please refer to your light manufacturer's instructions for curing conditions.

Ortholux LED Curing Light (App. 1000 mW/cm ²) (LED)	Ortholux Luminous Curing Light (App. 1600 mW/cm ²) (LED)
5 seconds through the bracket	3 seconds through the bracket

6. Archwires may be placed immediately after curing the last appliance.

Debonding Instrument

Use the recommended Unitek™ Debonding Instrument (REF 900-850) when debonding.

Important:

Clean the jaws of the debonding instrument after debonding each bracket to ensure even contact and force distribution.

Debonding Procedure

1. Remove adhesive flash around base of bracket to be debonded. Note: Failure to remove flash around bracket base, especially on the mesial-distal sides, may result in incomplete debonding.
2. Place the jaws of the Unitek™ Debonding Instrument on the mesial and distal sides of the bracket.
3. Gently squeeze instrument handles until the Unitek™ Gemini Clear bracket collapses. Gently rock the bracket in the mesial or distal direction to completely separate the bracket from the enamel, if necessary.

Note: Maintain the hold on the bracket to keep the debonded bracket parts in the instrument tips.

Bracket Rebonding Procedure

In the case of a spontaneous bond failure, it may be necessary to rebond a bracket. The following steps are recommended:

1. Carefully inspect the bracket for any damage. Brackets that have fractured through the vertical debonding slot cannot be rebonded and must be replaced. If a tie-wing is cracked, replace the bracket.
2. Remove any excess adhesive. Extra care must be taken to prevent chipping or breaking of the bracket. Use a hand scaler to remove any excess flash from around the edges of the bracket. **Do not use a bur. Do not attempt to scrape adhesive from the base of the bracket or attempt to micro-etch the adhesive as this may damage the bracket's bonding surface.**
3. If the bracket has been contaminated (e.g., moisture), rinse the bracket in isopropyl alcohol and allow to dry.
4. Prepare the tooth surface and bond the bracket using the procedure as described by the adhesive manufacturer.

Torquing Precautions

Unitek™ Gemini Clear ceramic brackets are capable of withstanding all normal torque requirements. However, care should be taken when making large torquing activations since large corrections with full size stainless steel wires may result in bracket failure and should be avoided.

Warranty and Limited Remedy

3M Unitek™ warrants that all products are free from defects in materials and manufacture. 3M Unitek™'s sole obligation and buyer's sole remedy in the event of any claimed defect shall be, at 3M Unitek™'s option, repair of the product, replacement of the product, or refund of the purchase price. Written notice of claimed defect must be received by 3M Unitek™ at the address below within reasonable time after discovery not to exceed one year from the date of delivery.

Waiver of Implied Warranties

ALL IMPLIED WARRANTIES, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR USE ARE DISCLAIMED.

Limitation of Liability

Except where prohibited by law, 3M Unitek™ will not be liable for any loss or damages arising from any product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Materials of animal and (or) human origin

Medical devices do not contain or deliver human blood derivatives.
No materials of animal origin are intentionally being added to the products.

Maintenance and repair

Not applicable

Cleaning and disinfection

Not applicable

Sterilization

Not applicable. Not supplied in a sterile state. In the orthodontic industry, it is common for orthodontic brackets to not be supplied in a sterile state.

While the orthodontic brackets are produced in a clean but not sterile environment, there is no need for cleaning or disinfection prior to use. Brackets are labeled for single use only:

Orthodontic brackets are not sterile devices; therefore, there is no need for testing of micro-organisms.

Operating, storage and transportation conditions

Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

For Unitek™ Gemini Clear coated with adhesive APC™ II

1. Do not expose materials to elevated temperatures or intense light. Bring to room temperature prior to use.
2. Do not store material in proximity to eugenol-containing products as this may inhibit proper curing of the adhesive.
3. This system is designed to be used at room temperature (68-77°F, 20-25°C), the humidity and pressure are not regulated.

Do NOT freeze.

Storage between 35-80°F, 2-27°C, the humidity and pressure are not regulated.

4. Shelf-life at room temperature is per the expiration date on the individual blister pack lid. Rotate inventory to optimize shelf life.

5. Do NOT expose uncured APC™ II appliances to ambient light for extended periods of time. When not in use, store blisters in original packaging or dispensing unit at room temperature.

Do not expose materials to intense light as this may cause the adhesive to cure.

No special storage requirement for uncoated brackets. Product is designed to be used at 20-25°C and stored at room temperature (5-40°C), the humidity and pressure are not regulated.

Disposal

Dispose of contents/container in accordance with applicable local/regional/national/international regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" (to the waste class B).

Technical Characteristics

		Gemini Clear			
		.018 Slot		.022 Slot	
Bracket Dimension Definition	Units	Min	Max	Min	Max
Bracket Depth	mm	1,63	2,44	1,63	2,44
Bracket Height	mm	3,30	4,45	3,30	4,45
Bracket Length	mm	3,48	4,39	3,48	4,39
Hook Length	mm	0,76	1,02	0,76	1,02
Angulation (ISO 27020 3.2)	Degrees	-4	9	-4	12
Angle of Torque (ISO 27020 3.1)	Degrees	-18	18	-23	18
Slot Height (ISO 27020 3.14)	mm	0,46	0,49	0,55	0,59
Slot Depth (ISO 27020 3.13)	mm	1,14	1,85	1,14	1,85
Bracket In - Out (ISO 27020 3.10)	mm	0,71	1,60	0,71	1,60
Slot Length (ISO 27020 3.15)	mm	2,51	3,68	2,51	3,68

Weight, g	
Brackets without adhesive (1 pc.)	Brackets with adhesive (1 pc.)
0,08±10%	0,1±10%

Information about medical device manufacturer:

Manufacturer


3M Unitek Corporation
2724 South Peck Road, Monrovia, California, 91016, USA

Autorised representative of the manufacturer in Russian Federation

Russia, 108811, Moscow, Poselenie Moskovskiy, Kievskoe shosse, 22 km, 6, bld.1
Tel. +7 (495) 784 7474

International Symbols

Symbol	Description
	Batch code
	Part Number, Catalog Number
	Keep away from sunlight
	Temperature limitation
	Use by date
	Caution, consult accompanying documents or Attention; see instructions for use

	Do not re-use, single use or use only once
	CE Mark, conformity to European Directives with Notified Body Number
	Authorized representative in the European Community
	Manufacturer
	Peel back to open.
	Do not reseal.
	Palmer notation – teeth numbering system.
	Palmer notation – teeth numbering system.

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