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Erdemoglu and Mungan T. Significance of detecting insulinlike growth factor binding protein-1 in cervicovaginal secretions: comparison with nitrazine test and amniotic fluid volume assessment. Acta Obstet Gynecol Scand (2004) 83:622-626.

Jain K and Morris PG. A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. J Obstet Gynaecol (1998) 18:33-36.

Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. J Obstet Gynecol Res (1998) 24:411-417.

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Rutanen E-M et al. Evaluation of a rapid strip test for insulinlike growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clin Chim Acta (1996) 253:91-101.

Rutanen E-M, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the RO Mcheck Membrane Immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta (1993) 214:73-81.

See Annex 1: Analytical studies and Clinical evidence

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Oy

Me

Tee

(

(sig


04.

INSTRUCTIONS FOR USE

Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test)

“APPROVED”

Oy Medix Biochemica Ab, Finland
Member of the Board of Directors
Teemu Alahuhta

(signature)


04.05.2020

This is to certify that TEEMU ALAHUHTA
is/are authorized to solely/jointly sign for the
company OY MEDIX BIOCHEMICA AB

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1. Name of medical product and Contents of delivery:

Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test)
(hereinafter "Test" or "Actim PROM Test")

Package options:

1. Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 1 test, comprised of:

- one sterile polyester swab for specimen collection;
- one dipstick;
- one tube of Specimen Extraction Solution (0.5 ml);
- Instructions for Use.

2. Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 10 tests, comprised of:

- 10 sterile polyester swabs for specimen collection;
- 10 dipsticks;
- 10 tubes of Specimen Extraction Solution (0.5 ml);
- Instructions for Use.

3. Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 20 tests, comprised of:

- 20 sterile polyester swabs for specimen collection;
- 20 dipsticks;
- 20 tubes of Specimen Extraction Solution (0.5 ml);
- Instructions for Use.

2. Developer, manufacturer, site of manufacture**Medical product manufacturer:**

Oy Medix Biochemica Ab
Klovinpellontie 3, FI-02180 Espoo, Finland,
Tel.: +358 9 547 680, fax: +35895053441, medix@medixbiochemica.com

Medical product developer:

Oy Medix Biochemica Ab
Klovinpellontie 3, FI-02180 Espoo, Finland,
Tel.: +358 9 547 680, fax: +35895053441, medix@medixbiochemica.com

Site of manufacture:

Oy Medix Biochemica Ab
Noljakantie 13, 80130 Joensuu, Finland,
medix@medixbiochemica.com,
Tel.: +358 9 547 680, fax: +35895053441

3. Authorized representative in the Russian Federation:

LLC Medicana Pharm,

3, 1st Pavlovsky Side-Str., Moscow, 115093

Tel./fax: 981-01-03, 937-21-26, 980-78-53,

E-mail – info@medicanapharm.ru

4. Medical product description

Purpose	Actim PROM Test is a visually interpreted, qualitative immunochromatographic dipstick test for detection of amniotic fluid in vaginal secretions during pregnancy. Actim PROM Test detects IGFBP-1, which is a major protein in amniotic fluid and a marker of the amniotic fluid in a vaginal sample. The test is intended for professional use to help diagnose the premature rupture of fetal membranes (PROM) in pregnant women.
Domain of usage	Clinical laboratory diagnostics
Intended user	Medical worker
Indications	The premature rupture of fetal membranes (PROM)
Contra indications	No
Sterility	Product components are supplied non-sterile, not subject to sterilization by user, except for the sterile polyester swab for specimen collection, sterilized by ethylene oxide.
Content of medicines	Does not contain medicines
Number of uses	Single-use product
Maintenance and repairs	Not applicable to a single-use product
Potential risk class	1

5. Product appearance and design:


Fig. 1. Test components (top to bottom):

- tube of Specimen Extraction Solution (0.5 ml);
- dipstick;
- sterile polyester swab for specimen collection.


Fig. 2. Dipstick design

(1) – Dip area; (2) – Result area; (3) – Test line; (4) – Control line

6. Product description and principle

Actim PROM Test is a visually interpreted, qualitative immunochromatographic dipstick test for detection of amniotic fluid in vaginal secretions during pregnancy. The Test detects IGFBP-1, which is a major protein in amniotic fluid and a marker of the amniotic fluid in a vaginal sample. The test is intended for professional use to help diagnose the rupture of fetal membranes (ROM) in pregnant women. The product is a test system for rapid in vitro diagnostics of ROM in pregnant women. The test is intended for professional use and can be applied both in clinical laboratory conditions and at any other place where a patient is located. It is a single use product.

The concentration of IGFBP-1 (insulin-like growth factor binding protein-1) in amniotic fluid is 100 to 1,000 times higher than in maternal serum. IGFBP-1 is not usually present in the vagina, but after rupture of fetal membranes, amniotic fluid with a high concentration of IGFBP-1 mixes with vaginal secretions. In Actim PROM Test, a specimen of vaginal secretion is taken with a sterile polyester swab and the specimen is extracted into Specimen Extraction Solution. The presence of IGFBP-1 in the solution is detected using a

dipstick. The -makes it possible to detect enhanced concentrations of IGFBP-1 (above 25 µg/l in the extracted sample).

The test is based on immunochromatography. It involves two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick (designated as (1) on fig. 2) is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains IGFBP-1 it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if IGFBP-1 is bound to them, they bind to the catching antibody. A blue line (test line - designated as (3) on fig. 2) will appear in the result area if the concentration of IGFBP-1 in the sample exceeds the detection limit of the test. A second blue line, the control line (designated as (4) on fig. 2), confirms correct performance of the test.

7. Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- If rupture of fetal membranes has occurred but the leakage of amniotic fluid has ceased more than 12 hours before the specimen is taken, IGFBP-1 may have been degraded by proteases in the vagina and the test may give a negative result.
- A negative test result is an indication on the present condition and cannot be used to predict the forthcoming.

8. Precautions

- For external use only
- Dispose of after use. A test cannot be used more than once.
- Do not use after shelf life expiry
- Do not use if a pouch is pierced or is not sealed properly.
- Keep away from children.

NOTES!

- No quantitative interpretation should be made based on the test results.
- The test requires about 150 µl of extracted sample to ensure proper performance of the test.
- Care must be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminum foil pouch or the seals of the pouch are not intact.
- Use only the swab provided with the kit.
- Improper sampling may lead to false negative result.
- When dipping, be careful to hold the dipstick in position (with the dip area in the sample extract) until the sample liquid front reaches the result area.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.
- At five minutes the appearance of any faint-to-dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 5 minutes.
- If only the control line is visible, the result should be interpreted as negative only after 5 minutes have elapsed.

- A positive Actim PROM Test result, although detecting the presence of amniotic fluid in the sample, does not locate the site of the rupture.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

9. Specimen type and collection

The specimen is vaginal secretion that is extracted into the Specimen Extraction Solution provided. A sample is obtained using a sterile polyester swab (provided in the kit). The sample should be collected prior to performing digital examination and/or transvaginal ultrasound. Take care not to touch anything with the swab before taking the sample. Separate the labia and carefully insert the tip of the swab into the vagina toward the posterior fornix until resistance is met. Alternatively the sample can be taken from the posterior fornix during a sterile speculum examination. The swab should be left in the vagina for 10-15 seconds to allow it to absorb the vaginal secretion (stage 5 on fig. 3).

Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10-15 seconds (stage 6 on fig. 3). Press the swab against the wall of the Specimen Extraction Solution tube to remove any remaining liquid from the swab. Discard the swab.

Specimens should be tested as soon as possible after extraction but in any case no more than 4 hours after specimen collection and extraction. If a specimen cannot be tested within this time it should be frozen. After thawing, the specimens should be mixed and tested as described below.

10. Test procedure and interpretation of the results

- If stored refrigerated, allow the aluminum foil pouch and the Extraction Solution tube to reach room temperature.
- Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick.
- The dipstick must be used shortly after its removal from the foil pouch.
- Place the yellow dip area into the extracted sample (stage 7 on fig. 3) and hold it there until you see the liquid front enter the result area (stage 8 on fig. 3). Remove the dipstick from the solution and place it in a horizontal position.
- The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at 5 minutes. Do not pay attention to any lines appearing later than 5 minutes.
- If two blue lines, the test line and the control line, appear, the test result is positive. If one blue line, the control line, appears, the test result is negative (stage 9 on fig. 3).
- If the control line or any lines do not appear at 5 minutes, the test is invalid (stage 10 on fig. 3). In this case, the study should be repeated using a new test.


Fig. 3. Testing outline

11. Product specifications:

Swab for specimen collection

Specification	Value
Swab length, mm	150±2
Swab weight, g	1.4±0.1

Tube of Specimen Extraction Solution

Specification	Value
Tube dimensions (height x diameter), mm	(65 x 11) ±1
Tube lid dimensions (height x diameter), mm	(9 x 14) ±1
Tube weight with the lid and solution, g	2.8±0.1
Solution volume in the tube, ml	0.50±0.05
Solution pH	6.0±0.1
Solution appearance	Colorless transparent liquid free from foreign matters, odorless

Dipstick

Specification	Value
Dipstick length, mm	120±1
Dipstick width, mm	6±1
Dipstick thickness, mm	0,8 – 2,0
Dipstick weight, g	0.7±0.1

Analytical sensitivity

The analytical sensitivity (detection limit) of the Actim PROM Test was identified by evaluating different concentrations of IGFBP-1 in extracted samples on three different lots of the Actim PROM Test. Two different operators each interpreted ten devices run at each concentration under various lighting conditions for a total of 60 determinations per level. The Actim PROM Test limit of detection is approximately 25 µg/l in extracted sample (calibration in accordance with Behring Institut IGFBP-1 (PP12) requirements, preparation of Batch 307/323). The measuring range of the Actim PROM is approximately 25-500,000 µg/l in extracted sample.

Analytical specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP-1 proteins at concentrations ranging from 10-5,000 µg/l of each protein in extracted sample using three different lots of the Actim PROM Test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim PROM Test is specific to human IGFBP-1.

Inter- and intra-assay reproducibility

A panel of specimens consisting of samples of different IGFBP-1 concentration levels was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using three different lots of the Actim PROM Test. Repeatable results were obtained. For reproducibility a study of Actim PROM Test was conducted at three separate sites. Test operators (n=9) tested each level on five different days using one lot of the Actim PROM Test. A total of 360 tests were performed (120 per site) with a total of 45 tests per sample type. The overall reproducibility of the Actim PROM Test is 97 % (350/360) with no significant differences within runs (replicates tested by one operator), between runs (five different days), between sites (three sites) or between operators (nine operators).

Interfering substances

The following substances, conditions and micro-organisms were tested with Actim PROM Test and were found not to affect Actim PROM Test performance

TESTED SUBSTANCE/CONDITION/ MICRO-ORGANISM	INTERFERENCE (NO/YES)
Whole blood	No
Semen	No
Urine	No
pH 3.5-8.5	No
Ultrasound Transmission Gel	No
Pevaryl (active ingredient: econatzol.nitras)	No
Gyno-Trosyd (tioconazol)	No
Flagyl (metronidazole)	No
Canesten (clotrimazol)	No
Personal Lubricant	No

Baby Oil	No
Baby Powder	No
Feminine Deodorant Suppositories	No
Vaginal Gel	No
Shower and bath products	No
Candida albicans	No
Gardnerella vaginalis	No
Neisseria gonorrhea	No
Chlamydia trachomatis	No
HSV-1	No
HSV-2	No

Clinical (diagnostic) sensitivity and clinical (diagnostic) specificity

The Actim PROM Test has been evaluated in several clinical studies as an instrument for detection of premature rupture of fetal membranes.

The studies showed that sensitivity and specificity of the test are very high, namely, sensitivity is 86-100%, and specificity is 74 – 100%.

12. Source materials

Dipstick:

Component	Material
Plastic frame	APET, polyester
Sample pad / Conjugate pad	Hydrophilic porous polyethylene, containing about 1-3 µg of murine monoclonal antibodies to IGFBP-1 (insulin-like growth factor binding protein-1), put on latex.
Membrane	Nitrocellulose ethereous membrane put on a polyester frame containing about 0,5-1,5 µg of murine monoclonal antibodies to IGFBP-1 (insulin-like growth factor binding protein-1)
Absorbent pad	Cellulose
Capping tapes	Polyester, polyethylene

Antibodies used in the dipstick:

Antibody	Function	Amount of component put on the dipstick
Monoclonal antihuman murine antibody IGFBP-1	Test line	0,5-1,5 µg
Polyclonal rabbit antimurine antibody IgG	Control line	1-3 µg
Monoclonal antihuman murine antibody IGFBP-1	Conjugated with latex particles and spread over the sample / conjugate	1-3 µg

	pad	
--	-----	--

Specimen Extraction Solution:

Component	Material
Plastic tube	Polypropylene
Lid	High density polypropylene
Specimen extraction solution (pH 6.0)	Sodium phosphate solution (Na ₂ PO ₄ x H ₂ O), 3.5 mg
	Sodium chloride (NaCl), 1.5 mg
	Ethylenediamine tetraacetic acid (EDTA), 1.9 mg
	Additive: Proclin 300, 0.5 mg
	Polysorbate: Tween-20, 0.15 ml
	Bovine serum albumin, 5 mg
	Aprotinin 25 KIU

Sterile swab for specimen collection:

Made by Puritan Medical Products Company LLC, Maine, U.S.

Component	Material
Grip	Polystyrol (made by Puritan Medical Products Company LLC, Maine, U.S.)
Body	Polyester (made by Puritan Medical Products Company LLC, Maine, U.S.)

Packaging materials

Title	Material
Individual aluminum dipstick pouch	SC 50 g/sq.m PE/LD 12 g/sq.m ALU 9 µm ION 23 g/sq.m
Inner pack	PET 12 + PE 50
Outer container	Cardboard 410 g/sq.m

Dryer

Title	Material
Dryer	Silica gel 0.5 g

List of materials of human and/or animal origin

Materials of human origin: not used

Materials of animal origin:

Source materials of animal origin used for product manufacture are murine monoclonal antibodies, polyclonal rabbit antibodies, and bovine albumin.

Bovine albumin used for product manufacture is obtained from sources certified as free from BSE (bovine spongiform encephalopathy).

Murine monoclonal antibodies and polyclonal rabbit antibodies are procured from a safe source, the material corresponds to requirements.

Raw materials are processed in line with proper laboratory practices.

13. Labeling

13.1 Labeling on the outer container:

Label on every box (outer container) contains the following information:

- name, address and website of the manufacturer;
- address of the site of manufacture;
- product name;
- number of tests for which the package is intended;
- catalogue number (REF) depending on product version:
REF 30830ETAC for Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 1 test,
or
REF 30831ETAC for Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 10 tests,
or
REF 30832ETAC for Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 20 tests;
- batch code (LOT);
- barcode;
- "In Vitro Diagnostic Medical Device" symbol;
- "Consult instructions for use" symbol;
- symbol of compliance with primary EU directive requirements;
- "Single use" symbol;
- storage conditions – "Temperature limitation" symbol;
- shelf life - "Use By" symbol;

In addition, a label for circulation on the territory of the Russian Federation is attached to the outer container. This label contains the following:

- full product name;
- number of products in the package;
- list of components;
- transportation conditions;
- safety precautions;
- storage conditions;
- name of the manufacturer;
- name and telephones of the distributor in Russia;
- number and date of the Registration Certificate;

Fig. 4

the source, the

membranes

membranes

membranes

uter


Fig. 4. Outer package as in the case of Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 10 tests


Fig. 5. Label layout for circulation on the territory of the Russian Federation on the outer package of the product as in the case of Premature Rupture of Fetal Membranes (Actim PROM Test)

13.2 Labeling on the inner pack:

The inner pack of rapid tests carries the following information:

- name, address and website of the manufacturer;
- address of the site of manufacture;
- product name;
- number of tests for which the pack is intended;
- catalogue number (REF);
- batch code (LOT);
- contents description;
- diagrammatic instructions for use;
- "In Vitro Diagnostic Medical Device" symbol;
- "Consult instructions for use" symbol;
- symbol of compliance with primary EU directive requirements;
- "Single use" symbol;
- storage conditions – "Temperature limitation" symbol;
- shelf life - "Use By" symbol;

In addition, a label for circulation on the territory of the Russian Federation is attached to the inner pack. This label contains the following:

- full product name;
- list of components;
- transportation conditions;
- safety precautions;
- storage conditions;
- name of the manufacturer;
- name and telephones of the distributor in Russia;
- number and date of the Registration Certificate;

e of the

pack.


Fig. 6. Labeling on the inner pack as in the case of Premature Rupture of Fetal Membranes (Actim PROM Test)

13.3 Labeling on the individual dipstick pouch:

The individual dipstick pouch carries the following information:

- name of the manufacturer;

- abbreviated product name (Actim PROM);
- product name;
- catalogue number (REF);
- batch code (LOT);
- “In Vitro Diagnostic Medical Device” symbol;
- symbol of compliance with primary EU directive requirements;
- storage conditions – “Temperature limitation” symbol;
- shelf life - “Use By” symbol;

In addition, a label for circulation on the territory of the Russian Federation is attached to the individual dipstick pouch. This label contains the following:

- product name (dipstick – 1 pce)


Fig. 7. Individual dipstick pouch

13.4 Labeling on the individual pouch of the sterile polyester swab for specimen collection:

The individual pouch of the sterile polyester swab for specimen collection carries the following information:

- name, address and website of the manufacturer;
- product name;
- catalogue number (REF);
- batch code (LOT);
- shelf life - “Use By” symbol;
- “Single use” symbol;
- “In Vitro Diagnostic Medical Device” symbol;
- “Do not use if package is damaged” symbol;
- “Sterilization by ethylene oxide” symbol;
- barcode;
- country of manufacture;
- symbol of compliance with primary EU directive requirements;

In addition, a label for circulation on the territory of the Russian Federation is attached to the individual

pouch of the specimen collection swab. This label contains the following:

- product name (sterile polyester swab for specimen collection – 1 pce)


Fig. 8. Individual pouch of the sterile polyester swab for specimen collection

13.5 Labeling on the tube of Specimen Extraction Solution (0.5 ml)

The tube of Specimen Extraction Solution (0.5 ml) carries the following information:

- name of the manufacturer;
- name of tests for which the tube is intended (Actim Partus / PROM Test)
- catalogue number (REF);
- batch code (LOT);
- shelf life - "Use By" symbol;
- "In Vitro Diagnostic Medical Device" symbol;
- storage conditions – "Temperature limitation" symbol;
- symbol of compliance with primary EU directive requirements;
- tube contents volume;

In addition, a label for circulation on the territory of the Russian Federation is attached to the tube of Specimen Extraction Solution. This label contains the following:

- product name (tube of Specimen Extraction Solution (0.5 ml) – 1 pce)


Fig. 9. Tube of Specimen Extraction Solution (0.5 ml)

13.6 Explanation of label symbols

Symbol	Definition
	Use By...
	Contains sufficient for <n> tests
	In Vitro Diagnostic Medical Device
	Single use
	Batch code
	Catalogue number
	Manufacturer

	Consult instructions for use
	Compliance with primary EU directive requirements
	Temperature limitation
	Do not use if package is damaged
	Sterilization by ethylene oxide

14. Packaging

14.1 Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 1 test:

Individual dipstick package: one dipstick with the dryer (silica gel) is placed in a retort pouch of aluminum film (fig. 7).

Individual package of the sterile polyester swab for specimen collection: one sterile swab is placed in a paper retort pouch (fig. 8);

Specimen Extraction Solution (0.5 ml) is put in a lidded tube (fig. 9);

An inner pack (fig. 6) contains:

- one sterile polyester swab for specimen collection;
- one dipstick;
- one tube of Specimen Extraction Solution (0.5 ml).

An outer package (cardboard box) (fig. 5) contains:

- one inner pack with products;
- Instructions for Use.

Fig. 10. Outer package as in the case of Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 1 test

14.2 Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 10 tests:

Individual dipstick package: one dipstick with the dryer (silica gel) is placed in a retort pouch of aluminum film (fig. 7).

Individual package of the sterile polyester swab for specimen collection: one sterile swab is placed in a paper retort pouch (fig. 8);

Specimen Extraction Solution (0.5 ml) is put in a lidded tube (fig. 9);

An inner pack (fig. 6) contains:

- one sterile polyester swab for specimen collection;
- one dipstick;
- one tube of Specimen Extraction Solution (0.5 ml).

An outer package (cardboard box) (fig. 5) contains:

- 10 inner packs with products;
- Instructions for Use.

14.3. Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 20 tests:

Individual dipstick package: one dipstick with the dryer (silica gel) is placed in a retort pouch of aluminum film (fig. 7).

Individual package of the sterile polyester swab for specimen collection: one sterile swab is placed in a paper retort pouch (fig. 8);

Specimen Extraction Solution (0.5 ml) is put in a lidded tube (fig. 9);

An inner pack (fig. 6) contains:

- one sterile polyester swab for specimen collection;
- one dipstick;
- one tube of Specimen Extraction Solution (0.5 ml).

An outer package (cardboard box) (fig. 5) contains:

- 20 inner packs with products;
- Instructions for Use.


Fig. 11. Outer package as in the case of Rapid Test for In Vitro Diagnostics of Premature Rupture of Fetal Membranes (Actim PROM Test) for 20 tests

Package
Maximum

Specifica

Outer co

(length x

Outer co

(length x

Outer co

(length x

Inner pa

Individu

Individu

15

Store th
unopene
pack is
compon
If storag
Use the
Shelf lif

10

The pro
external
to the s
Transpo
Handlin
handlin
If trans

1

Use the
Use the

1

In line
For the
medica
items
Regula

The m

n PROM

ed in a

Package dimensions:

Maximum dimensional inaccuracies are $\pm 5\%$

Specifications	Value
Outer container of Actim PROM Test for 1 test (length x width x height), mm	112x16x35
Outer container of Actim PROM Test for 10 tests (length x width x height), mm	236 x 110 x 70
Outer container of Actim PROM Test for 20 tests (length x width x height), mm	235 x 121 x 70
Inner pack (length x width), mm	233 x 99
Individual swab pouch (length x width), mm	193 x 40
Individual dipstick pouch (length x width), mm	170 x 55

15. Storage conditions

Store the product at the temperatures of $+2...+25\text{ }^{\circ}\text{C}$ and relative air humidity of up to 80% in the unopened manufacturer's original container until the "Use By" date indicated on the package. If the test pack is opened, each component must be used until the "Use By" date indicated on the specific component.

If storage conditions are not observed, then the medical product is unfit for use.

Use the product immediately after opening the foil pack.

Shelf life: 24 months.

16. Transportation

The product should be transported in the manufacturer's original container that protects the product from external impact using all types of transport in closed vehicles in line with transportation rules applicable to the specific type of transport.

Transportation conditions: temperature of $+2...+25\text{ }^{\circ}\text{C}$ and relative air humidity of up to 80%

Handling instructions on the container and transport pack must always be observed in the course of handling operations and transportation.

If transportation conditions are not observed, then the medical product is unfit for use.

17. Usage conditions

Use the product indoors at the temperature of $+18...+25\text{ }^{\circ}\text{C}$ and relative air humidity of up to 80%.

Use the dipstick immediately after opening the foil pouch.

18. Disposal

In line with clinic protocols and local legislation requirements.

For the purposes of disposal, based on the morphological composition, the product can be classified as medical waste of epidemiological hazard (Infected and potentially infected waste. Materials and tools, items contaminated with blood and/or other biological liquids) – Class B according to Sanitary Regulations and Norms 2.1.7.2790-10.

of Fetal

19. Warranties

The manufacturer guarantees the compliance of the medical product performance to that stated in the

operational documentation, subject to the use in accordance with the instructions and intended use, provided for by the manufacturer, by the indicated "use by" date.

For all questions, technical advice and support contact the authorized representative of the manufacturer in the territory of the Russian Federation:

LLC Medicana Pharm,
3, 1st Pavlovsky Side-Str., Moscow, 115093
Tel./fax: 981-01-03, 937-21-26, 980-78-53,
E-mail – info@medicanapharm.ru

20. Applicable standards

IVDD 98/79/EC	In Vitro Diagnostic Medical Devices Directive
ISO 15223-1:2016	Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied Part 1: general requirements.
ISO 15223-2:2010	Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied Part 2: symbol development, selection and validation.
EN ISO 13485:2016 AC2016, AC2018	Medical devices -- Quality management systems -- Requirements for regulatory purposes
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices
EN ISO 18113-1:2012	In vitro diagnostic medical devices -- Information supplied by the manufacturer (labelling) -- Part 1: Terms, definitions and general requirements
EN ISO 18113-2:2012	In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling) – Part 4: In vitro diagnostic reagents for professional use
EN ISO 23640:2015	In vitro diagnostic medical devices. Evaluation of stability of in vitro diagnostic reagents
EN 13612: 2002 / AC: 2002	Performance evaluation of in vitro diagnostic medical devices

21. Risk management

All risks detected in the course of risk assessment were reduced to an acceptable level.

22. Information on development and manufacture

The quality management system that supports development and manufacture is certified in accordance with EN ISO 13485:2016.

, provided

facturer in

Literature:

Erdemoglu and Mungan T. Significance of detecting insulinlike growth factor binding protein-1 in cervicovaginal secretions: comparison with nitrazine test and amniotic fluid volume assessment. *Acta Obstet Gynecol Scand* (2004) 83:622-626.

Jain K and Morris PG. A clinical study to evaluate the usefulness of the MAST test in diagnosing prelabour rupture of membranes. *J Obstet Gynaecol* (1998) 18:33-36.

Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. *J Obstet Gynecol Res* (1998) 24:411-417.

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Rutanen E-M et al. Evaluation of a rapid strip test for insulinlike growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. *Clin Chim Acta* (1996) 253:91-101.

Rutanen E-M, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the RO Mcheck Membrane Immunoassay in the diagnosis of ruptured fetal membranes. *Clin Chim Acta* (1993) 214:73-81.

See Annex 1: Analytical studies and Clinical evidence

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