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Hedström J et al. Urine trypsinogen-2 as a marker of acute pancreatitis. *Clin Chem* (1996) 42:685-690.

Jang T et al. Point-of-care Urine Trypsinogen Testing for the Diagnosis of Pancreatitis. *Acad. Emerg. Med.* (2007) 14:29-34.

Kylänpää-Bäck M-L et al. Reliable screening for acute pancreatitis with rapid urine trypsinogen-2 test strip. *Br J Surg* (2000), 87:49-52.

Kylänpää-Bäck M-L et al. Comparison of urine trypsinogen-2 test strip with serum lipase in the diagnosis of acute pancreatitis. *Hepato-Gastroenterology* (2002) 49:1130-1134.

Kemppainen E et al. Rapid measurement of urinary trypsinogen-2 as a screening test for acute pancreatitis. *N Engl J Med* (1997) 336: 1788-1793.


“
O
M
T
(
(s
04

INSTRUCTIONS FOR USE

Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis)

“APPROVED”

Oy Medix Biochemica Ab, Finland
Member of the Board of Directors
Teemu Alahuhta

(signature)

04.05.2020


This is to certify that _____
is/are authorized to ~~solely~~ ^{jointly} sign for the

company OY MEDIX BIOCHEMICA AB

Espoo, Finland

04. 05. 2020

Ex officio:

SATU SYDÄNTALVI
ylitarkastaja, julkinen notaari
överinspektör, notarius publicus
Senior Officer, Notary Public

1. Name of medical product:

Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) (hereinafter "Test" or "Actim Pancreatitis")

Package options:

1. Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 10 tests, comprised of:

- 10 dipsticks;
- Instructions for Use.

2. Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests, comprised of:

- 20 dipsticks;
- Instructions for Use.

2. Developer, manufacturer, site of manufacture**Medical product manufacturer:**

Oy Medix Biochemica Ab
Klovinpellontie 3, FI-02180 Espoo, Finland,
Tel.: +358 9 547 680, fax: +35895053441, medix@medixbiochemica.com

Medical product developer:

Oy Medix Biochemica Ab
Klovinpellontie 3, FI-02180 Espoo, Finland,
Tel.: +358 9 547 680, fax: +35895053441, medix@medixbiochemica.com

Site of manufacture:

Oy Medix Biochemica Ab
Noljakantie 13, 80130 Joensuu, Finland,
medix@medixbiochemica.com,
Tel.: +358 9 547 680, fax: +35895053441

3. Authorized representative in the Russian Federation:


LLC Medicana Pharm,
3, 1st Pavlovsky Side-Str., Moscow, 115093
Tel./fax: 981-01-03, 937-21-26, 980-78-53,
E-mail – info@medicanapharm.ru

4. Medical product description

Purpose	Actim Pancreatitis test is a visually interpreted, qualitative, immunochromatographic dipstick test for use when acute pancreatitis is suspected. The test is based on detection of elevated levels of tripsynogen-2 in urine.
Domain of usage	Clinical laboratory diagnostics

Intended user	Intended for a professional user: clinical laboratory diagnostician or laboratory technician (medical technologist)
Indications	Suspected acute pancreatitis
Contra indications	No
Sterility	The product is supplied non-sterile, not subject to sterilization by user.
Content of medicines	Does not contain medicines
Number of uses	Single-use product
Maintenance and repairs	Not applicable to a single-use product
Potential risk class	1

5. Product appearance and design:

(1) – Dip area; (2) – Result area; (3) – Test line; (4) – Control line

Fig. 1. Dipstick appearance and design

6. Product description and principle

Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) is a qualitative test for rapid *in vitro* diagnostics of acute pancreatitis on the basis of urine samples. The product is intended for professional use and can be applied both in clinical laboratory conditions and at any other place where a patient is located.

Levels of trypsinogen-2 have been shown to increase in the urine of patients with acute pancreatitis.

The Test makes it possible to detect enhanced concentration (more than 50 µg/l) of trypsinogen-2 in urine.

The Test is based on immunochromatography. It involves two types of monoclonal antibodies to human trypsinogen-2 spread on the dipstick. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch labelled particles and indicate a positive result.

When the dip area of the dipstick (designated as (1) on fig. 1) is placed in a urine sample, the dipstick absorbs liquid, which starts to flow up the dipstick under the impact of capillary forces. If the sample contains trypsinogen-2 it binds to the antibody attached to the latex particles. The particles are carried by the liquid flow and, if trypsinogen-2 is bound to them, they bind to the catching antibody. A blue line (test line – designated as (3) on fig. 1) will appear in the result area if the concentration of trypsinogen-2 in the sample exceeds the cut-off value for the test. A second blue line (control line – designated as (4) on fig. 1) confirms correct performance of the test.

7. Limitations of the test

The test is intended for professional *in vitro* diagnostic use only.

8. Precautions

- For external use only.
- Dispose of after use. A test cannot be used more than once.
- Do not use after shelf life expiry
- Do not use if a pouch is pierced or is not sealed properly.
- Keep away from children.

Attention!

- No quantitative interpretation should be made based on the test results.
- Care must be taken when placing the dipstick in sample. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing
- Do not use the dipstick if its aluminum foil pouch is not intact.
- When dipping, be careful to hold the dipstick in position (with the dip area in the urine sample) until the sample liquid front reaches the result area.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- Excessive amount of biotin in urine sample may interfere for the formation of control line.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.
- At 5 minutes, the appearance of any faint-to dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 5 minutes.
- If only the control line is visible, the result should be interpreted as negative only after 5 minutes have elapsed.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.


9. Specimen type

The specimen is urine. Urine samples should be tested as soon as possible after collection. Samples may be stored in a refrigerator at +2...+8 °C for 24 hours. If the sample needs to be stored for longer, it should be frozen. The sample should be brought to room temperature before testing.

10. Test procedure and interpretation of the results

- If stored refrigerated, allow the aluminum foil pouch to reach room temperature.
- Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick.
- If necessary, identifying marks may be written on the upper pink part of the dipstick.
- The dipstick must be used shortly after its removal from the foil pouch.
- Place the yellow dip area into the sample and hold it there until you see the liquid front enter the result area.
- Remove the dipstick from the sample and place it horizontal.
- The result can be interpreted as positive as soon as two blue lines become visible in the result area. If only one blue line, control line, becomes visible at 5 minutes, the result can be interpreted as negative.

- Do not pay attention to any lines appearing later than 5 minutes.
- If two blue lines, the test line and the control line, appear, the test result is positive. If one blue line, control line, appears, the test result is negative.
- If the control line does not appear, the test is invalid. In this case, the study should be repeated using a new test.

Place the yellow dip area into the sample

Control the liquid front in the result area

Assess the result at 5 minutes.

The test is invalid if the control line does not appear or no lines appear at 5 minutes

11. Product specifications:

Specification	Value
Dipstick length, mm	120±1
Dipstick width, mm	6±1
Dipstick thickness, mm	0.8 – 2,0
Dipstick weight, g	0.7±0.1

Analytical sensitivity

The analytical sensitivity (detection limit) of the Actim Pancreatitis test was evaluated using samples with different concentrations of human trypsinogen-2 on three different lots of Actim Pancreatitis test. The detection limit of the Actim Pancreatitis test is approximately 50 µg/l of human trypsinogen-2 in sample and the results remain positive at least to 100 000 µg/l.

Repeatability. Inter- and intra-assay reproducibility

Intra- (repeatability) and inter- (reproducibility) assay precision were evaluated using samples containing 0-100 000 µg/l of human trypsinogen-2 in urine. For repeatability the samples were tested on the same day using ten replicates at the concentrations of 0 – 100 µg/l and three replicates at the concentrations of 30 000 and 100,000 µg/l. The testing was done on three different lots of Actim Pancreatitis test, the total amount being 198 tests in one day. For reproducibility the samples were tested on the five consecutive workdays using two replicates per sample on three different lots of Actim Pancreatitis test, the total amount being 240 tests in five days. Repeatability and reproducibility were at 100% based on results of those evaluations.

High dose hook effect:

Urine samples with a high trypsinogen-2 concentration were tested using two series of PANCREATITIS dipsticks. Trypsinogen-2 concentration in urine samples was also measured using the Trypsinogen-2 IEMA test of Medix Biochemica Oy.

PANCREATITIS dipsticks of all series produced a correct result event at very high trypsinogen-2 concentrations in urine: up to 348,390 µg/l.

Interfering substances

The following interfering substances were tested with the Actim Pancreatitis test. No interference of the tested substances was observed with the performance of the Actim Pancreatitis test

INTERFERING SUBSTANCES OR CONDITIONS	INTERFERENCE (NO/YES)
Glucose	No
Albumin	No
Hemoglobin	No
Acetylsalicylic acid	No
Ascorbic acid	No
Salicylic acid	No
Bilirubin	No
Biotin	No
pH 5.0-9.0	No

Cross reactions:

Trypsin-2 enters into cross reaction with the pancreatitis test at the concentration of 2,000 µg/l.

Clinical (diagnostic) sensitivity and clinical (diagnostic) specificity

The Actim Pancreatitis test was reviewed in the clinical research process as a tool for acute pancreatitis diagnostics.

The clinical group consisted of 500 randomly selected patients (306 males and 194 females), who were admitted to the intensive care unit of Central Hospital under the Helsinki University or to the intensive care unit of the Helsinki City Hospital.

The criteria for diagnosing acute pancreatitis were distinctive clinical signs and very high concentrations of amylase (over 900 u/l in serum, over 6,000 u/l in urine), clinical manifestations coupled with enhanced concentration of amylase (over 300 u/l in serum or over 2,000 u/l in urine) and results of computer tomography or ultrasound research indicating acute pancreatitis, or a distinctive clinical picture and positive computer tomography research contrasting to those of patients who did not have enhanced amylase concentrations at the time of diagnosing.

The research showed that the test for trypsinogen-2 detection in urine with the help of dipsticks helped diagnose acute pancreatitis more precisely than measurement of amylase content in serum or urine.

The research brought the following results: diagnostic sensitivity of 94%, diagnostic specificity of 95%.

Clinical research

The test was estimated in the course of clinical research.

12. Source materials

Dipstick:

Component	Material
Plastic frame	APET, polyester
Sample pad / Conjugate pad	Hydrophilic porous polyethylene, containing about 1-3 µg of murine monoclonal

	antibodies to human Tripsynogen-2, put on latex.
Membrane	Nitrocellulose ethereous membrane put on a polyester frame containing about 1-2 µg of murine monoclonal antibodies to human Tripsynogen-2
Absorbent pad	Cellulose
Capping tapes	Polyester, polyethylene

Antibodies used in the dipstick:

Antibody	Function	Amount of component put on the dipstick
Monoclonal antihuman murine antibody Trypsinogen-2	Membrane	1-2 µg
Monoclonal antihuman murine antibody Trypsinogen-2	Conjugated with latex particles and distributed over the sample / conjugate pad	1-3 µg
Monoclonal antihuman murine antibody (TSH)	Control line	1-2 µg

Packaging materials

Title	Material
Individual aluminum dipstick pouch	SC 50 g/sq.m PE/LD 12 g/sq.m ALU 9 µm ION 23 g/sq.m
Outer container	Cardboard TVK 310 g/sq.m

Dryer

Title	Material
Dryer	Silica gel

List of materials of human and/or animal origin

Materials of human origin: not used

Materials of animal origin:

Source materials of animal origin used for product manufacture are murine monoclonal antibodies. Murine monoclonal antibodies are procured from a safe source, the material corresponds to requirements.

Raw materials are processed in line with proper laboratory practices.

13. Labeling

13.1 Labeling on the outer container:

- name, address and website of the manufacturer;
- address of the site of manufacture;

- product name;
- number of tests for which the package is intended;
- catalogue number (REF) depending on product version;
REF 32731ETAC for Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 10 tests,
or
REF 32732ETAC for Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests;
- batch code (LOT);
- barcode;
- "In Vitro Diagnostic Medical Device" symbol;
- "Consult instructions for use" symbol;
- symbol of compliance with primary EU directive requirements;
- "Single use" symbol;
- storage conditions – "Temperature limitation" symbol;
- shelf life - "Use By" symbol;

In addition, a label for circulation on the territory of the Russian Federation is attached to the outer container. This label contains the following:

- full product name;
- number of products in the package;
- list of components;
- transportation conditions;
- safety precautions;
- storage conditions;
- name of the manufacturer;
- name and telephones of the distributor in Russia;
- number and date of the Registration Certificate;

Fig. 2. Outer package as in the case of Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests

Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests:

- Dipsticks - 20 pcs.;
 - Instructions for use.
 Number: 20 rapid tests;
 To be stored at: +2 °C.... +25 °C
 To be transported at: +2 °C.... +25 °C, humidity <85%
 Manufacturer: Oy Medix Biochemica Ab
 Klovipellontie 3, FI-02180 Espoo, Finland
 Distributor in Russia: LLC Mediana Pharm
 3, block 1, office 2, 1st Pavlovsky Side-Str., Moscow, 115093
 Tel./fax: 981-01-03, 937-21-26, 980-78-53
 Registration Certificate: No. FSZ 2010/07124 dd.

For external use only
Do not use after shelf life expiry
Keep away from children

Fig. 3. Label layout for circulation on the territory of the Russian Federation as in the case of Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests

13.2 Label on individual dipstick pouch:

The individual dipstick pouch carries the following information:

- name of the manufacturer;
- product name (Actim Pancreatitis);
- catalogue number (REF);
- batch code (LOT);
- “In Vitro Diagnostic Medical Device” symbol;
- symbol of compliance with primary EU directive requirements;
- storage conditions – “Temperature limitation” symbol;
- shelf life - “Use By” symbol;

In addition, a label for circulation on the territory of the Russian Federation is attached to the individual dipstick pouch. This label contains the following:

- product name (dipstick – 1 pce)


Fig. 4. Individual dipstick pouch


Fig. 5. Dipstick without an individual pouch.

13.3 Explanation of label symbols

Symbol	Definition
	Use By...

	Contains sufficient for <n> tests
	In Vitro Diagnostic Medical Device
	Single use
	Batch code
	Catalogue number
	Manufacturer
	Consult instructions for use
	Compliance with primary EU directive requirements
	Temperature limitation

14. Packaging

14.1 Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 10 tests:

Individual dipstick package: one dipstick with the dryer (silica gel) is placed in a retort pouch of aluminum film (fig. 4).

An outer package (cardboard box) contains:

- Dipsticks in individual packages – 10 pieces;
- Instructions for use – 1 piece.

14.2 Rapid Test for In Vitro Diagnostics of Acute Pancreatitis (Actim Pancreatitis) for 20 tests:

Individual dipstick package: one dipstick with the dryer (silica gel) is placed in a retort pouch of aluminum film (fig. 4).

An outer package (cardboard box) contains:

- Dipsticks in individual pouches – 20 pieces;
- Instructions for use – 1 piece.

14.3 Package dimensions:

Maximum dimensional inaccuracies are $\pm 5\%$

Outer container of Actim Pancreatitis for 10 tests (length x width x height), mm	180 x 60 x 60
Outer container of Actim Pancreatitis for 20 tests (length x width x height), mm	180 x 60 x 60
Individual dipstick pouch (length x width), mm	170 x 55

15. Storage conditions

Store the product at the temperatures of +2...+25 °C and relative air humidity of up to 80% in the unopened manufacturer's original container until the "Use By" date indicated on the package.

If storage conditions are not observed, then the medical product is unfit for use.

Use the dipstick immediately after opening the individual pouch.

Shelf life: 24 months.

16. Transportation

The product should be transported in the manufacturer's original container that protects the product from external impact using all types of transport in closed vehicles in line with transportation rules applicable to the specific type of transport.

Transportation conditions: temperature of +2...+25 °C and relative air humidity of up to 80%

Handling instructions on the container and transport pack must always be observed in the course of handling operations and transportation.

If transportation conditions are not observed, then the medical product is unfit for use.

17. Usage conditions

Use the product indoors at the temperature of +18...+25 °C and relative air humidity of up to 80%

18. Disposal

In line with clinic protocols and local legislation requirements.

For the purposes of disposal, based on the morphological composition, the product can be classified as medical waste of epidemiological hazard (Infected and potentially infected waste. Materials and tools, items contaminated with blood and/or other biological liquids) – Class B according to Sanitary Regulations and Norms 2.1.7.2790-10.

19. Warranties

The manufacturer guarantees the compliance of the medical product performance to that stated in the operational documentation, subject to the use in accordance with the instructions and intended use, provided for by the manufacturer, by the indicated "use by" date.

For all questions, technical advice and support contact the authorized representative of the manufacturer in the territory of the Russian Federation:

LLC Medicana Pharm,

3, 1st Pavlovsky Side-Str., Moscow, 115093

Tel./fax: 981-01-03, 937-21-26, 980-78-53,

E-mail – info@medicanapharm.ru

20. Applicable standards

IVDD 98/79/EC	In Vitro Diagnostic Medical Devices Directive
EN 980:2008	Symbols to be used with medical device labels
EN ISO 13485:2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes
EN ISO 14971:2013	Medical devices - Application of risk management to medical devices
EN ISO 18113-1:2011	In vitro diagnostic medical devices -- Information supplied by the manufacturer (labelling) -- Part 1: Terms, definitions and general requirements
EN ISO 18113-2:2011	In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling) -- Part 4: In vitro diagnostic reagents for professional use
EN 13641:2002	Elimination or reduction of risk of infection related to in vitro diagnostic reagents
EN 13640:2002	Stability testing of in vitro diagnostic reagents
EN 13612:2002	Performance evaluation of in vitro diagnostic medical devices

21. Risk management

All risks detected in the course of risk assessment were reduced to an acceptable level.

22. Information on development and manufacture

The quality management system that supports development and manufacture is certified in accordance with EN ISO 13485:2016.

Literature:

Chen et al. Rapid Urinary Trypsinogen-2 Test Strip in the Diagnosis of Acute Pancreatitis. *Pancreas* (2005) 30:243-247.

Delcenserie et al. Diagnostic de la pancreatite aigue par detection du trypsinogene II urinaire (pancreatitis test). Communication orale: Club Francais du pancreas (1999).

Hedström J et al. Urine trypsinogen-2 as a marker of acute pancreatitis. *Clin Chem* (1996) 42:685-690.

Jang T et al. Point-of-care Urine Trypsinogen Testing for the Diagnosis of Pancreatitis. *Acad. Emerg. Med.* (2007) 14:29-34.

Kylänpää-Bäck M-L et al. Reliable screening for acute pancreatitis with rapid urine trypsinogen-2 test strip. *Br J Surg* (2000), 87:49-52.

Kylänpää-Bäck M-L et al. Comparison of urine trypsinogen-2 test strip with serum lipase in the diagnosis of acute pancreatitis. *Hepato-Gastroenterology* (2002) 49:1130-1134.

Kemppainen E et al. Rapid measurement of urinary trypsinogen-2 as a screening test for acute pancreatitis. *N Engl J Med* (1997) 336: 1788-1793.

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