State of Minnesota

SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Document Regarding 3M Red Dot Prewired ECG Monitoring Electrodes

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Red Dot™ Neonatal/Infant Prewired Monitoring Electrodes in English is truthful and accurate.

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
3M Red Dot Prewired ECG Monitoring Electrodes

Product Description:

The information provided in this insert applies to the 3M Red Dot 2269T Neonatal and Infant Monitoring Electrodes. The 3M Red Dot 2269T monitoring electrodes are disposable with a silver/silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 3 days wear, are radiolucent and repositionable.

Intended use

The Red Dot Monitoring ECG Electrodes are intended for use in ECG Monitoring.

They can be used in all ECG monitoring applications where standard infant/neonatal ECG monitoring electrodes are used. These electrodes can be used for short term and long term monitoring.

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Caution:

The 3M Red Dot 2269T Neonatal and Infant Monitoring Electrode was determined to be MR-Conditional via non-clinical testing. For information on the use of a 3M Red Dot 2269T Neonatal and Infant monitoring electrode during an MRI procedure please see below.

- U.S.A. federal law restricts this device to sale by or on the order of the physician.


Precautions:

- To prevent dry-out, electrodes should be used immediately after opening the package.
- Do not use an electrode if the conductive adhesive is dry and no longer sticky.
- Do not prepare neonate and infant skin with alcohol or an abrader to minimize skin irritation.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Assess electrode sites periodically.

Instructions For Use:

Step 1:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring of neonates and infants.
- Clean the electrode site with a mild soap and water solution, or water only. The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Leave the skin moist, but not wet. Do not prepare neonates' or infants' skin with alcohol or abrasion.

Step 2:

- Remove electrodes from sealed package. Separate and straighten lead wires.
- Remove electrode from the liner.

Step 3:

- Apply to prepared skin site and smooth the entire electrode down with gentle pressure in a circular motion around the outer edge first and then the inner element.
- Repeat application for remaining electrodes.

Step 4:

- Connect electrode lead wire to patient cable only.

Note:

- This electrode can be repositioned upon initial application. Select a new skin site for reapplication.

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2269T MRI Safety Information - MR Conditional

Non-clinical testing demonstrated that the Neonatal Monitoring Electrode with Pre-Attached Lead Wire, 2269T, is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)(extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined for the 3M Red Dot 2269T electrode is expected to produce a

- Lead positioned in a straight configuration (i.e., no loops or cross over points)

Artifact Information

In non-clinical testing, the image artifact caused by the 3M Red Dot 2269T electrode extends approximately 5-mm from the Electrode when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

IMPORTANT NOTE: To ensure patient safety and proper recording of the electrocardiogram during an MRI procedure, use only electrocardiograph (EKG/ECG) electrodes tested specifically with regard to the MRI environment and recommended by the manufacturer of the MR Conditional monitoring equipment. EKG electrodes should only be used according to the specific guidelines provided by the manufacturers of the MR Conditional, monitoring devices. To prevent injuries (e.g., burn), for a patient undergoing an MRI procedure, EKG electrodes must be removed if they are not connected to acceptable MR Conditional monitoring equipment.

Shelf Life:

Shelf life is 2 years. For shelf life, refer to the expiration date that is printed on each package.

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3M, Red Dot and Trace Prep are trademarks of 3M. Used under license in Canada.

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APOSTILLE

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1. Country: United States of America

This Public Document: Document Regarding 3M Red Dot ECG Monitoring Electrodes with Foam tape and Sticky Gel

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Red Dot™ ECG Monitoring Electrode with Foam Tape and Sticky Gel 2560 series in English is truthful and accurate.

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3M Red Dot ECG Monitoring Electrodes with Foam tape and Sticky Gel

Product Description:

The information provided in this insert applies to the 3M™ Red Dot™ Monitoring Electrodes with Foam tape and Sticky Gel 2560 series with snap connector and sticky gel.

The 3M™ Red Dot™ Monitoring Electrodes are disposable with a silver/ silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 5 days wear.

Intended use

3M™ Red Dot™ Monitoring Electrodes with Foam tape and Sticky Gel are intended for use in ECG Monitoring.

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Warning:

The 3M™ Red Dot™ Monitoring Electrodes with with a metal snap connector have not been tested for use during Magnetic Resonance Imaging (MRI) procedures. The use of conductive patient-connected devices and patient lead wires/electrodes in MRI procedures may result in serious patient burns. If a physician orders either of these devices to be used during MRI procedures, the user should not allow loops in patient lead wires or allow the lead wires to come in contact with the patient's skin.

Caution:

U.S.A. federal law restricts this device to sale by or on the order of the physician.

Precautions:

- To prevent dryout, electrodes should be used within 45 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time. ○ Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin abrasion should be part of the skin preparation. Use 3M Red Dot™ Trace Prep as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the edge of the electrode.

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APOSTILLE

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1. Country: United States of America

This Public Document: Document Regarding 3M Red Dot Repositionable Monitoring ECG Electrodes

2. has been signed by Irma Vargas

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Irma Vargas, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: November 20, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9808230-9

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

3M Health Care Business

3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M Red Dot Repositionable Monitoring Electrode, 2600 series, with snap connector and conductive adhesive in English is truthful and accurate

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3M Red Dot Repositionable Monitoring ECG Electrodes

Product Description

The information provided in this insert applies to the 3M Red Dot Repositionable Monitoring Electrode, 2600 series, with snap connector and conductive adhesive. The 3M Red Dot Repositionable Monitoring Electrodes are disposable with a silver/silver-chloride sensing element and are designed for ECG monitoring. The 3M Red Dot Repositionable Monitoring Electrodes are radiolucent and may be safely worn up to 3 days.

Intended use

The Red Dot Repositionable Monitoring ECG Electrodes are intended for use in ECG Monitoring.

This device is not offered as a sterile product, nor is it intended for sterilization before use. These devices should only be used by trained medical professionals.

Contraindications: Contraindications are not available

Caution

The 3M Red Dot Repositionable Monitoring Electrode 2600 series were determined to be MR-Conditional via non-clinical testing. For information on the use of a 3M Red Dot Repositionable monitoring electrode during an MRI procedure please see below.

U.S.A. Federal Law restricts this device to sale by or on the order of the physician.

Precautions

- To prevent dry-out, electrodes should be used immediately after opening the package.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
 - During surgical procedures, place ECG electrodes as far as possible from the electrosurgical site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin. Application sites with heavy or curly hair should be clipped.
- For best trace quality, skin impedance reduction should be part of the skin preparation. Use 3M Red Dot Trace Prep 2236 as part of the skin preparation.
- Attach the lead-wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.

- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a circular motion.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the tab or the edge of the electrode.

Note

- Duration of wear up to 3 days may vary depending on skin condition and environmental factors.
- This electrode can be repositioned upon initial application. Select a new skin site for reapplication.

Contents: 5 electrodes/strip; 1 strip/bag; 200 bags/case

Technical specification

Electrode length 1.413 inches (35.9 mm) +/- 10%.

Electrode width 1.412 inches (35.9 mm) +/- 10%.

Diameter of electrode stud (snap) 0.375 inches (9.5 mm) +/- 10%.

Medical devices do not contain or deliver medicinal substances.

Medical devices do not contain or deliver human blood derivatives.

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2600 Series MRI Information - *MR Conditional*

3M Red Dot Repositionable Monitoring Electrode 2600 series was determined to be **MR-Conditional**.

Non-clinical testing demonstrated that the 3M Red Dot repositionable monitoring electrodes are MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

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MRI-Related Heating

In non-clinical testing, the 3M™ Red Dot™ repositionable monitoring electrode produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change +1.5°C

Therefore, the MRI-related heating experiments for the 3M Red Dot repositionable monitoring electrode at 3- Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/ kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.5°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the 3M™ Red Dot™ repositionable monitoring electrode. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 2-mm relative to the size and shape of the 3M™ Red Dot™ repositionable monitoring electrode.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
2670 Signal Void Size	1,096 mm ²	88 mm ²	1,401 mm ²	92 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

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Thank you for choosing 3M Red Dot Electrodes. If you are not completely satisfied or if you have questions or comments, in the U.S.A., please call the 3M™ Health Care Customer Helpline at 1-800-228-3957, in Canada 1-800-364-3577. Outside of the U.S.A. and Canada, please contact your 3M Distributor or your 3M Sales Representative.

Explanation of Symbols

	This product and package do not contain natural rubber latex.
	Caution, see instructions for use
	Do not reuse
	Use by date
	Batch code
	Manufacturer
	CE mark
	EC Rep

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APOSTILLE

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1. Country: **United States of America**

This Public Document: Document Regarding 3M Red Dot Monitoring ECG Electrodes

2. has been signed by **Irma Vargas**

3. acting in the capacity of **Notary Public, State of Minnesota**

4. bears the seal/stamp of **Irma Vargas, Notary Public, State of Minnesota**

CERTIFIED

5. at: **St. Paul, Minnesota**

6. Date: **November 20, 2017**

7. by: **Secretary of State, State of Minnesota**

8. File No: **9808230-8**

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

3M Health Care Business

3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110

3M

To whom it may concern


In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Red Dot™ Monitoring Electrodes with snap connector in English is truthful and accurate.

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3M Red Dot Monitoring ECG Electrodes

Intended Use and Product Description: The information provided in this insert applies to the 3M Red Dot Monitoring Electrodes 2200 series with snap connectors (2248, 2271). The 3M Red Dot monitoring electrodes are disposable with a silver/silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 3 days wear.

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Warning:

This device has not been tested for use during Magnetic Resonance Imaging (MRI) procedures. The use of conductive patient-connected devices and patient lead wires/electrodes in MRI procedures may result in serious patient burns. If a physician orders this device to be used during MRI procedures, the user should not allow loops in patient lead wires or allow the lead wires to come in contact with the patient's skin.

Caution: U.S.A. federal law restricts this device to sale by or on the order of the physician.

Precautions:

- To prevent dryout, electrodes should be used within 30 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.

- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin. Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin abrasion should be part of the skin preparation. Use 3M Red Dot Trace Prep or the electrode's built-in abrader as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Remove the paper liner from the tab if applicable.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the tab or the edge of the electrode.

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SECRETARY OF STATE

APOSTILLE

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1. Country: **United States of America****This Public Document: Document Regarding 3M Red Dot Radiolucent Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series**2. has been signed by **Irma Vargas**3. acting in the capacity of **Notary Public, State of Minnesota**4. bears the seal/stamp of **Irma Vargas, Notary Public, State of Minnesota**

CERTIFIED

5. at: **St. Paul, Minnesota**6. Date: **November 20, 2017**7. by: **Secretary of State, State of Minnesota**8. File No: **9808230-7**

9. Seal/Stamp:

10. Signature:

A handwritten signature in black ink, reading "Steve Simon", written over a horizontal line.

Steve Simon
Secretary of State

3M Health Care Business

3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110


3M

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M Red Dot 2570 series (2570, 2570-3, and 2570-5) Radiolucent Monitoring Electrodes with snap connector and sticky gel in English is truthful and accurate.


Regina Feferman
International Regulatory Affairs Associate
3M Health Care

Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 17th day of November 2017.

Irma Vargas, Notary Public
My Commission Expires: January 31, 2020

3M Red Dot Radiolucent Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series

Product Description:

The information provided in this insert applies to the 3M Red Dot 2570 series (2570, 2570-3, and 2570-5) Radiolucent Monitoring Electrodes with snap connector and sticky gel. The 3M Red Dot Monitoring Electrodes are disposable with a silver/silver chloride sensing element and are designed for ECG monitoring. These electrodes have been tested for up to 5 days wear and are radiolucent.

Intended use

3M Red Dot Radiolucent Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series are intended for use in ECG Monitoring.

This device is not offered as a sterile product, nor it is intended for sterilization before use. These devices should only be used by trained medical professionals.

Contraindications: Contraindications are not available.

Warning:

Caution:

• **The 3M Radiolucent Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series were determined to be MR-Conditional via non-clinical testing.** For information on the use of a 3M Red Dot 2570 series Monitoring Electrode during an MRI procedure please see below.

• U.S.A. federal law restricts this device to sale by or on the order of the physician.

Precautions:

- To prevent dryout, electrodes should be used within 45 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical site to minimize RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.

Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin abrasion should be part of the skin preparation. Use 3M Red Dot Trace Prep or the electrode's built-in abrader as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.

- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- To minimize skin trauma during electrode removal, pull up slowly and gently from the edge of the electrode.

Contents

Catalog Number	Contents
2570	50/bag - 20 bag /case
2570-3	3/bag - 200 bag/case
22570-5	5/bag - 200 bag/case

Technical specification

Electrode length 1.56 inches (39.6 mm) +/- 10%.

Electrode width 1.38 inches (35 mm) +/- 10%.

Diameter of electrode stud (snap) 0.375 inches (9.5 mm) +/- 10%.

Medical devices does not contain or deliver medicinal substances.

Medical devices does not contain or deliver human blood derivatives.

No materials of animal origin are intentionally being added to the products.

2570 MRI Information - MR Conditional

Red Dot Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series

were determined to be MR-Conditional via non-clinical testing.

A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the Red Dot Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, hDX, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: *Highest temperature change* +1.5°C

Therefore, the MRI-related heating experiments for the 2570 Series Cardiac Monitoring Electrode at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9-W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.5°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Red Dot Monitoring Electrodes with Foam Tape and Stick Gel 2570 Series. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 2-mm relative to the size and shape of the Red Dot Monitoring Electrodes with Foam

Tape and Stick Gel 2570 Series.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
2570 Signal Void Size	1,483 mm ²	130 mm ²	1,822 mm ²	273 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Important Note: This ECG Electrode that underwent MRI testing was not connected to an electrocardiograph (ECG) lead and/or an MR conditional monitoring system. Therefore, the safety of using this ECG Electrode attached to an ECG lead and/or MR conditional monitoring system is unknown. Furthermore, an MR conditional monitoring system typically requires use only with ECG electrodes specifically designed and tested with that specific equipment. It is unknown if this ECG Electrode is acceptable for use with an MR conditional ECG monitoring system designed to be used to monitor a patient undergoing an MRI procedure.

Shelf Life:

Shelf life is 2 years.

For shelf life, refer to the expiration date that is printed on each package.

Storage and Transportation conditions:

Normal room conditions: 15-30°C (59-86°F), 35-60% relative humidity.

During transportation, the products may experience short term temperature, humidity and pressure variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Working conditions:

Normal room temperature of 15-30°C (59-86°F), The values of humidity and pressure are not regulated

Disposal

Dispose of content/container in accordance with the local/regional/national/international regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B).

Maintenance and repair

This medical device isn't subject to maintenance and repair

Thank you for choosing 3M Red Dot Electrodes. If you are not completely satisfied or if you have questions or comments, in the U.S.A., please call the 3M Health Care Customer Helpline at 1-800-228-3957, in Canada 1-800-364-3577. Outside of the U.S.A. and Canada, please contact your 3M Distributor or your 3M Sales Representative.

Explanation of Symbols

	This product and package do not contain natural rubber latex.
	Caution, see instructions for use

	Do not reuse
	Use by date
	Batch code
	Manufacturer
	CE mark
	EC Rep

Warranty

3M warrants this product will be free from defects in material and manufacture. Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Made in Canada

 3M Company, 3M Health Care,
3M Center, 2510 Conway Avenue, Bldg. 275-5W-06, St. Paul, Minnesota, 55144, USA
1 800 364-3577
www.3M.com

Manufacturing site:

3M CANADA CO.
400 Rt. 100, Morden Manitoba, Canada R6M 1Z9

Authorized representative on the territory of the Russian Federation:

ZAO "3M Russia"
Russia 121614, Moscow, Krylatskaya str. bldg.17/3
Tel.: +7 (495) 784 7474 / Fax: +7 (495) 784 7475
e-mail 3mra@mmm.com

3M, Red Dot and Trace Prep are trademarks of 3M. Used under license in Canada.
Issue Date: 2017-10

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