INSTRUCTION FOR USE

The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma specimen

LOT: NCP20070003, NCP20070004

Test-cassette – 25; buffer, LOT 20070525 – 1; droppers – 25; instruction for use – 1

NAME OF DIAGNOSTIC DEVICE

The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma specimen, LOT: NCP20070003, NCP20070004

Test-cassette – 25; buffer – 1; droppers – 25; instruction for use – 1

Short title: "2019-nCoV IgG/IgM Rapid Test Cassette (whole blood, serum or plasma specimen)", diagnostic device, kit.

INTENDED USE/PURPOSE

The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019-nCoV in whole blood, serum or plasma specimen - is a qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019-nCoV in whole blood (finger and venous), serum or plasma specimen.

The diagnostic device is intended for single use

For clinical laboratory diagnostics in vitro

FUNCTIONAL PURPOSE

It is used as an auxiliary tool in the diagnosis of primary and secondary coronavirus infection

SCOPE

Clinical laboratory diagnostics

POTENTIAL CUSTOMER

Clinical and diagnostic laboratories, laboratory divisions of medical and preventive institutions

PROFESSIONAL LEVEL OF POTENTIAL CUSTOMER

Only for professional in vitro diagnostics by properly qualified medical staff. Medical staff are obliged to work with gloves on

PROPOSED CLINICAL USE

▣ Suspected infection with SARS-CoV-2 coronavirus (earlier name 2019-nCoV), when common symptoms occur and contact with an infected person

▣ To detect COVID-19 virus infection in an asymptomatic course of the disease

COUNTERINDICATIONS

▣ Expiration date

▣ Loss of sealing

▣ Improper storage and transportation conditions

▣ If used in accordance with the instructions for use there are no contraindications

POTENTIAL SIDE EFFECTS

There are no side effects associated with the use of a diagnostic device

GENERAL PROPERTIES

The kit is meant to be used for the study of 25 samples and is intended for use in healthcare institutions

Time for analysis : 10 minutes

STERILIZATION

Non-sterile

The supplied medical device is not intended for user sterilization or other applicable microbiological control

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Description:

Test cassette is a flat white quadrat plate with an immunochromatographic test strip inside. The test cassette is packed in a sealed foil moisture-proof pouch, which contains a package with adsorbent granules. A bottle of buffer solution contains 3 ml of buffer solution. Pipettes for transferring a single-use sample are packed in a zip-bag. Test cassettes, a bottle of buffer solution, pipettes and instructions for use are packed in a consumer cardboard box. One set is designed to perform 25 rapid analyses.

Note:

- One sealed aluminum foil pouch includes 1 test cassette and a desiccant;
- The product contains no materials that come into direct or indirect contact with the body of the patient and the personnel using the product, when meeting the requirements of the operating documentation (instructions for use);
- The composition of the product no medicinal products and pharmaceutical substances

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection containers
- Scarifiers
- Centrifuge (for plasma only)
- Lancets (for fingerstick whole blood only)
- Timer
- Capillary tubes (for fingerstick whole blood only)
- Pipette
- Disposable gloves
- Alcohol pad

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible by SARS-CoV-2 (2019-nCoV). Currently, the patients infected by the novel coronavirus SARS-CoV-2 (2019-nCoV) are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

TEST PRINCIPAL

This test consists of two components, an IgG component and an IgM component. In the IgG component, anti-human IgG is coated in IgG test line region. During testing, the specimen reacts with 2019-nCoV antigen-coated particles in the test cassette. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the anti-human IgG in IgG test line region, if the specimen contains IgG antibodies to 2019-nCoV. A colored line will appear in IgG test line region as a result of this. Similarly, anti-human IgM is coated in IgM test line region and if specimen contains IgM antibodies to 2019-nCoV, the conjugate-specimen complex reacts with anti-human IgM. A colored line appears in IgM test line region as a result.

Therefore, if the specimen contains 2019-nCoV IgG antibodies, a colored line will appear in IgG test line region. If the specimen contains 2019-nCoV IgM antibodies, a colored line will appear in IgM test line region. If the specimen does not contain 2019-nCoV antibodies, no colored line will appear in either of the test line regions, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

PRECAUTIONS

1. This instructions must be read completely before performing the test. Failure to follow directions in instructions may yield inaccurate test results.
2. For professional in vitro diagnostic use only. Do not use after expiration date.
3. Do not eat, drink or smoke in the area where the specimens or kits are handled.
4. Do not use test if sealed foil pouch is damaged.
5. Handle all specimens as if they contain infectious agents.
6. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.

when specimens are assayed.

8. Please ensure that an appropriate amount of samples are used for testing. Too much or too little sample size may lead to deviation of results.

9. The used test should be discarded according to local regulations.

10. Humidity and temperature can adversely affect results.

SPECIMEN COLLECTION AND PREPARATION

The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma specimen can be performed using whole blood, serum and plasma specimen.

Both Fingerstick Whole Blood and Venipuncture Whole Blood can be used

- To collect Fingerstick Whole Blood specimens:
- Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger
- Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using a capillary tube:
- Touch the end of the capillary tube to the blood until filled to approximately 20 μ L. Avoid air bubbles.

Serum

- Separate the serum or plasma from blood as soon as possible to avoid hemolysis. Only clear, non-hemolyzed specimens can be used.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 7 days, for long term storage, serum/plasma specimens should be kept below -20°C.

Whole blood collected by venipuncture

- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.
- EDTA K2, Heparin sodium, Citrate sodium and Oxalate potassium can be used as anti-coagulants.

DIRECTIONS FOR USE

Bring tests, specimens, buffer, and/or controls to room temperature (15-30°C) before use

1. Bring the sealed foil pouch to room temperature before opening. Remove the test cassette from the sealed foil pouch and use it within one hour.
2. Place the test cassette on a clean and level surface.

For Serum or Plasma Specimens:

- To use a dropper: Hold the dropper vertically, draw the specimen up to the Fill Line (approximately 10 μ L), and transfer the specimen to the sample well of the test cassette, then add 2 drops of buffer (approximately 80 μ L) and start the timer. Avoid trapping air bubbles in the specimen well.
- To use a pipette: To transfer 10 μ L of specimen to the sample well of the test cassette, then add 2 drops of buffer (approximately 80 μ L) and start the timer.


For Venipuncture Whole Blood specimen:

- To use a dropper: Hold the dropper vertically, draw the specimen about 1 cm above the fill line and transfer 1 full drop (approx. 20 μ L) of specimen to the sample well(S). Then add 2 drops of buffer (approximately 80 μ L) and start the timer.
- To use a pipette: To transfer 20 μ L of whole blood to the specimen well(S), then add 2 drops of buffer (approximately 80 μ L), and start the timer.

For Fingerstick Whole Blood specimen:

- To use a dropper: Hold the dropper vertically, draw the specimen about 1 cm above the fill line and transfer 1 full drop (approx. 20 μ L) of specimen to the sample well(S). Then add 2 drops of buffer (approximately 80 μ L) and start the timer.
- To use a capillary tube: Fill the capillary tube and transfer approximately 20 μ L of fingerstick whole blood specimen to the specimen well (S) of test cassette, then add 2 drops of buffer (approximately 80 μ L) and start the timer. See illustration below.
- Wait for the colored line(s) to appear. The test result should be read at 10 minutes. Do not interpret the result after 20 minutes.

Note: It is suggested not to use the buffer, beyond 6 months after opening the vial care sites.

INTERPRETATION OF RESULTS

IgG POSITIVE:* Two colored lines appear. One colored line should always appear in the control line region (C) and another line should be in the IgG line region.

IgM POSITIVE:* Two colored lines appear. One colored line should always appear in the control line region (C) and another line should be in the IgM line region.

IgG and IgM POSITIVE:* Three colored lines appear. One colored line should always appear in the control line region (C) and two test lines should be in the IgG line region and IgM line region.

*NOTE: The intensity of the color in the test line regions may vary depending on the concentration of 2019-nCoV antibodies present in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE: One colored line appears in the control line region (C). No line appears in the IgG region and IgM region.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1.The test Procedure and the Interpretation of test Result must be followed closely when testing for the presence of SARS-CoV-2 virus specific antibodies in the serum, plasma or whole blood specimen from individual subjects. For optimal test performance, proper sample collection is critical. Failure to follow the procedure may give inaccurate results.

2.The test result will only indicate the presence of IgG and IgM antibodies to SARS-CoV-2 in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infections.

3.Positive results can be associated with both the current strain of coronavirus and other previously identified types of coronavirus, such as HKU1, NL63, OC43, 229E.

4.This kit of reagents is screening and is recommended for laboratory diagnostics of COVID-19. The results obtained with the test should be considered with other clinical findings from other laboratory tests and evaluations.

5.The Kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma specimen is for in vitro diagnostic use only. This test should be used for detection of IgG and IgM antibody to SARS- CoV-2 in whole finger and venous blood, serum or plasma specimens as an aid in the diagnosis of patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation and the results of other laboratory tests. Neither the quantitative value nor the rate of increase in the concentration of IgG or IgM antibodies to 2019-nCoV can be determined by this qualitative test.

6.If the test result is negative or non-reactive and clinical symptoms persist, It is recommended to resample the patient a few days later or test with a molecular diagnostic device to rule out infection in these individuals.

7.The hematocrit level of the whole blood can affect the test results. Hematocrit level needs to be between 25% and 65% for accurate results.

8.The test will show negative results under the following conditions: The titer of the novel coronavirus antibodies in the sample is lower than the minimum detection limit of the test, or the antibody utilized in the test, or the novel coronavirus antibody has not appeared at the time of

donation or to inform about the status of COVID-19.

10.Not for checking blood donations.

11.In the early infection, anti-SARS-COV-2 antibodies concentrations may be below detectable level. Therefore it is not recommended to use the test in early diagnosis of COVID-19.

12.The continued presence or absence of antibodies cannot be used to determine the success or failure of therapy.

13.Results from immunosuppressed patients should be interpreted with caution.

14.At this time, it is unknown how long IgM or IgG antibodies may persist following infection.

PERFORMANCE CHARACTERISTICS

Performance of IgG Specimen

Relative sensitivity: 100% (95%CI*: 86.0%~100%); Relative specificity: 98.0% (95%CI*: 89.4%~99.9%);

Accuracy: 98.6% (95%CI*: 92.3%~99.96%). *Confidence Intervals

Performance of IgM specimen

Relative sensitivity: 85.0% (95%CI*: 62.1%~96.8%); Relative specificity: 96.0% (95%CI*: 86.3%~99.5%);

Accuracy: 92.9% (95%CI*: 84.1%~97.6%). *Confidence Intervals

CROSS REACTIVITY

There was no cross-reaction with Anti-influenza A virus, anti-influenza B virus, anti-Adenovirus, Anti-Measles, HAMA, RF, non-specific IgG, non-specific IgM, anti-EV71, anti-Parainfluenza virus, anti-RSV, HBsAg, anti-Syphilis, anti-H. Pylori, anti-HIV and anti-HCV positive specimens at 10 and 20 minutes.

INTERFERING SUBSTANCES

100mg/dL Triglyceride, 1000mg/dL hemoglobin, 20mg/dL Ascorbic Acid, 15mmol/L Total cholesterol, 60mg/dL Bilirubin conjugated, 60mg/dL Bilirubin unconjugated, 2000mg/dL human Albumin, 1mg/ml anti-E.coli, 1mg/ml human IgG antibody, 1mg/ml human IgM antibody, 40mg/L Ribavirin, 30mg/L Oseltamivir, 200mg/L Levofloxacin, 100mg/L Azitromycin and 250mg/L Cefuroxime Axitilhave indicated no interference to the performance of The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma.

EXPIRATION DATE

The kit is stable for 24 months from the date of manufacture.

The expired kit is not applicable.

STORAGE, STABILITY, TRANSPORTATION

Store the test at 2-30°C.

DO NOT FREEZE!

Humidity: below 70 %

Keep away from sunrays

The kit is stable for 24 months from the date of manufacture. (not confirmed by real-time research). The buffer is stable for 180 days after the first opening.

Use the test cassette for no more than one hour after opening the sealed foil pouches at 15-30°C.

Increased humidity and temperature may adversely affect the test result.

The test cassette from a damaged sealed foil pouch is not suitable for tests.

Keep the kit out of the reach of children!

It should be transported in the original packaging that protects the product from external influences, on all types of covered transport and in covered vehicles in accordance with the transport regulations applicable to this type of transport at a temperature of +2 to +30°C for the entire shelf life.

During loading and unloading operations and transportation, you must strictly comply with the requirements of the handling marks applied on the packaging and transport containers.

It is not recommended to use the buffer 6 months after opening the bottle.

requirements of international standards Directive 98/79/EC, ISO 13485.

On issues related to quality and handling The kit of "2019-nCoV IgG/IgM Rapid Test Cassette" for qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to 2019 in whole blood, serum or plasma specimen please contact: Closed Joint Stock Company «Scientific Production Association «MEDCAR», Moscow region, Ramenskoye, st. Moskvoretskaya, d. 62, room. 206

Tel: +7 (495)540-45-96. Email: med@medcar.ru. <https://medcar.ru/>

The manufacturer of diagnostic device

Hangzhou AllTest Biotech Co., Ltd. 550, Yinhai Street, Hangzhou Economic and Technologic Development Area,- 310018, ZJ, P.R. China.

Information about authorized representative in the territory of the Russian Federation

Closed Joint Stock Company «Scientific Production Association «MEDCAR»

Moscow region, Ramenskoye, st. Moskvoretskaya, d. 62, room. 206. Tel: +7 (495)540-45-96. Email: med@medcar.ru. <https://medcar.ru/>


RECOVERY

Disposal must be made in accordance with local waste management regulations. Contaminated packaging must be disposed of in the same manner as the product.

Non-Contaminated packaging materials may be recycled. Contact your local service providers for further informations.

BIBLIOGRAPHY

1. Westgard JO, Barry PL,Hunt MR, Groth T. A multi-rule Shewhart for quality control in clinical chemistry, Clinical Chemistry 1981;27:493-501

Hangzhou AllTest Biotech Co., Ltd.
550, Yinhai Street, Hangzhou Economic and Technologic Development Area,- 310018, ZJ, P.R. China.

Information about authorized representative in the territory of the Russian Federation Closed Joint Stock Company «Scientific Production Association «MEDCAR» Moscow region, Ramenskoye, st. Moskvoretskaya, d. 62, room. 206 Tel: +7 (495)540-45-96 Email: med@medcar.ru <https://medcar.ru/>

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