APPROVED BY

Sophie REHAULT
Quality and Regulatory Affairs Director

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November 05th 2019

INSTRUCTIONS FOR USE

CATCH+ Device for Thrombectomy of Cerebral Arteries

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AN. **REHAULT**
certificat exclusivement la certifier la signature

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These Instructions for Use have been prepared to provide information for registering a medical device in the Russian Federation. The information is provided as required in accordance with the Rules on Registration of Medical Devices (Decree No. 1416 of the Government of the Russian Federation dated December 27, 2012). The contents of the Instructions for Use comply with the requirements of Order No. 11Н of the Ministry of Healthcare of the Russian Federation "On Approval of Requirements to the Contents of Technical and Operational Documentation of the Manufacturer (Producer) of the Medical Device" dated January 19, 2017.

The device is developed, manufactured and sold in accordance with the requirements of DIN EN ISO 13485:2016 + AC : 2017-07, EN ISO 13485:2016 + AC:2016, ISO13485:2016; the device is CE marked (complies with the requirements of Annex II, excluding Section 4, Annex II.4 of Directive 93/42 / EEC).

1. Device name:

CATCH+ Device for Thrombectomy of Cerebral Arteries:

Design options:

- CATCHMINI 3x5;
- CATCHMINI 3x10;
- CATCHMINI 3x15;
- CATCHMINI 3x20;
- CATCH 4x20;
- CATCHMAXI 6x20;
- CATCHMAXI 6x30;
- CATCHMAXI 6x40;
- Instructions for Use (supplied with each device).

All devices are supplied as follows:

CATCH+ Device for Thrombectomy of Cerebral Arteries in the design option – 1 pc, in a sterile package. The sterile package and instructions for use are packed in individual cardboard packaging.

Device classification in accordance with the Nomenclature Classifier of Medical Devices approved by Order No. 4Н of the Ministry of Healthcare of the Russian Federation dated June 6, 2012 – 333980;

OKPD2 – 32.50.50.190 – Instruments and devices used for medical purposes, not included in other groups;

Class depending on the potential risk degree of medical application in accordance with the Nomenclature Classifier of Medical Devices approved by Order No. 4Н of the Ministry of Healthcare of the Russian Federation dated June 6, 2012 – 3;

Type of contact with the human body – short-term (less than 24 hours) contact with circulating blood.

2. Full name of the medical device manufacturer

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3. *Intended use of the device and principles of operation*

3.1 Intended use of the device: The device is intended for use to restore blood flow in the neurovasculature by removing clots in patients experiencing acute ischemic stroke due to a large cerebral artery occlusion.

3.2 Principle of operation of the device:

The CATCH+ Device for Thrombectomy of Cerebral Arteries (hereinafter referred to as the Device, the CATCH+ Device) is a self-expandable system mounted on its own pusher and designed for catching and removing clots in patients. The device is intended for use to restore blood flow in the neurovasculature by removing clots in patients experiencing acute ischemic stroke due to a large cerebral artery occlusion. The Stentriever design makes it possible to catch a clot within cells and restore blood flow. The device is placed and deployed inside an occluded vessel within the thrombus formation and is then withdrawn in its deployed state. The mechanical thrombectomy procedure does not take more than 1 hour.

The appropriate pusher should be used to position the CATCH+ Device so that the effective length portion of the device would extend past each side of the thrombus in the vessel when fully deployed. If necessary, and if the device is not damaged, the same CATCH+ Device can be used for up to two flow restoration procedures in the same patient after cleaning the device with saline solution. The CATCH+ Device for Thrombectomy of Cerebral Arteries should not be used for more than two flow restoration procedures and more than three restoration attempts in the same vessel.

The Device and its pusher (also fitted with a radiopaque marker) are placed inside a PTFE (polytetrafluoroethylene) insertion tube and packaged in a peelable pouch and a cardboard box, together with the instructions for use.

The CATCH+ Device is intended for single patient use only and is sterilized with ethylene oxide. In fact, these products are surgically invasive devices intended for transient use (less than 60 minutes) in contact with the central circulatory system (arteriae cerebrales).

3.3 Information on potential consumers

The CATCH+ Device for Thrombectomy of Cerebral Arteries should only be used by medical specialists trained in the field of interventional neuroradiology and ischemic stroke treatment. *Conditions of use:* Under fluoroscopic control in medical treatment facilities and healthcare centers, only when prescribed by a medical specialist.

4. *Functional characteristics of the device*

The CATCH+ Device for Thrombectomy of Cerebral Arteries is a self-expandable device designed for capturing and removing clots in patients and mounted on its own pusher designed for delivery of the capture device to the site of the clot formation. The Stentriever design makes it possible to catch a clot within cells and restore blood flow. The CATCH+ Device for Thrombectomy of Cerebral Arteries is placed and deployed inside an occluded vessel within the thrombus formation and is then withdrawn in its deployed state.

With the pusher, the CATCH+ Device for Thrombectomy of Cerebral Arteries is positioned in the vessel so that the effective length portion of the device would extend past each side of the thrombus when fully deployed.

Note: All devices are visually identical and differ only in size.

5. Risks of using the medical device, contraindications, expected and predictable side effects:

5.1 Indications:

The device is intended for use to restore blood flow in the neurovasculature by removing clots in patients experiencing acute ischemic stroke due to a large cerebral artery occlusion. This treatment option is possible for patients who fail intravenous tissue plasminogen activator (IV t-PA) therapy or ineligible for it.

The CATCH+ Device should only be used by physicians trained in interventional neuroradiology and treatment of ischemic stroke. The CATCH+ Device should be used under fluoroscopic control and with the administration of appropriate anticoagulants.

5.2 Contraindications:

The CATCH+ Device for Thrombectomy of Cerebral Arteries is contraindicated for use in the following cases:

- Patients with known hypersensitivity to nickel-titanium alloy.
- Patients with stenosis located near the thrombus site that may prevent safe withdrawal of the CATCH+ Device for Thrombectomy of Cerebral Arteries.
- Patients with angiographic evidence of carotid dissection.
- These devices are contraindicated for use in newborns, premature babies and infants.

5.3 Potential complications

The complications that may occur as a result of the mechanical thrombectomy procedure include all possible complications that may occur during or after the intervention. Most adverse effects are attributed to clinical progression that may occur in the patient, or to the risks associated with any endovascular brain surgery:

- Haematoma and bleeding at the puncture site;
- Perforation, rupture, dissection or other arterial injuries;
- Vessel spasms;
- Vessel dissection;
- Difficulty when removing the device;
- Subarachnoid hemorrhage;
- Cramp of joints;
- Mental change;
- Neurological deterioration, including stroke and death;
- Ischemia;
- Infectious disease;
- Air embolism;
- Intracranial hemorrhage;
- Vascular occlusion;
- Pseudoaneurysm formation, including at the femoral artery access site;
- Postoperative hemorrhage;
- Distal embolization, including embolization of the previously uninvolved area;
- Adverse reaction to antiplatelet/anticoagulant agents or contrast media;
- Deformation, collapse, fracture or malfunction of the device(s);

- Thrombosis (acute and subacute);
- Arteriovenous fistula;
- Hemorrhagic conversion.

5.4 Precautions

The specific clinical risks associated with the use of the CATCH+ Device for mechanical thrombectomy in patients with acute ischemic stroke are as follows:

- Device toxicity;
- Failure to deploy the device and remove the clot;
- Traumatic nature of the device resulting in arterial dissection;
- Device collapse with material dispersion inside arteries;
- Inadequate sliding within a microcatheter resulting in the system blockage;
- Radiopaque marker detachment and migration through the vasculature;
- Failure to perform thrombectomy and thrombus fragmentation;
- Occlusion of previously unaffected arteries;
- Device-related adverse events.

6. Technical specifications of the device

6.1. Description of basic functional characteristics:

The CATCH+ Device for Thrombectomy of Cerebral Arteries is a self-expandable device designed for capturing and removing clots in patients and mounted on its own pusher designed for delivery of the capture device to the site of the clot formation. The Stentriever design makes it possible to catch a clot within cells and restore blood flow. The CATCH+ Device for Thrombectomy of Cerebral Arteries is placed and deployed inside an occluded vessel within the thrombus formation and is then withdrawn in its deployed state.

With the pusher, the CATCH+ Device for Thrombectomy of Cerebral Arteries is positioned in the vessel so that the effective length portion of the device would extend past each side of the thrombus when fully deployed.

Note: All devices are visually identical and differ only in size.


Figure 1.1 Appearance of the CATCH+ Device for Thrombectomy of Cerebral Arteries


Figure 1.2 Appearance of the CATCH+ Device for Thrombectomy of Cerebral Arteries in an introducer

Note: an introducer is not a component of the medical device and is used solely to protect the medical device from creases, bends and damage.


Figure 1.3 Appearance of the CATCH+ Device for Thrombectomy of Cerebral Arteries, capture device component


Figure 1.4 Appearance of the CATCH+ Device for Thrombectomy of Cerebral Arteries, proximal part of the pusher


Figure 1.5 Drawing of the general view of the CATCH+ Device for Thrombectomy of Cerebral Arteries with indication of components

Figure 2. Appearance of the capture device, location of the markers

Note: «ORX» - radiopaque marker

The CATCH+ Device for Thrombectomy of Cerebral Arteries includes a full range of thrombectomy devices for acute ischemic stroke treatment.

Table 1. Technical specifications of the CATCH+ Device for Thrombectomy of Cerebral Arteries.

Specification	CATCH MINI 3x5	CATCH MINI 3x10	CATCH MINI 3x15	CATCH MINI 3x20	CATCH 4x20	CATCH MAXI 6x20	CATCH MAXI 6x30	CATCH MAXI 6x40
Total length, mm	2013	2018	2024	2028	2031	2031	2043	2051
Pusher length, mm	2000±1.5	2000±1.5	2000±1.5	2000±1.5	2000±1.5	2000±1.5	2000±1.5	2000±1.5
Effective length of the capture device, mm	5	10	15	20	20	20	30	40
Total length of the capture device, mm	13.67	18.67	23.67	28.67	31.65	31.65	42.85	51.10
Undeployed capture device diameter, mm	0.41	0.41	0.41	0.41	0.52	0.60	0.60	0.60
Deployed capture device diameter (inside the vessel), mm	4	4	4	4	5	6.5	6.5	6.5
Cell diameter in deployed state, mm (max)	3.38	3.38	3.38	3.38	3.40	4.04	4.04	4.04
Proximal diameter of the pusher, mm	Ø 0.38 ± 0.005	Ø 0.38 ± 0.005	Ø 0.38 ± 0.005	Ø 0.38 ± 0.005	0.47± 0.005	0.47± 0.005	0.47± 0.005	0.47± 0.005
Distal diameter of the pusher, mm	Ø 0.18 ± 0.01							
Weight, g	0.01	0.01	0.02	0.02	0.02	0.03	0.03	0.03
Compatible guide catheter diameter, mm	1.67	1.67	1.67	1.67	1.67	2	2	2
Number of radiopaque markers at the distal end of the Catch+	3	3	3	3	3	3	3	3
Number of radiopaque markers at the pusher	1	1	1	1	1	1	1	1
Marker dimensions at the distal end of the Catch+, mm	0.15x0.27							
Overall dimensions of marker at the pusher, mm	0.29x0.39x0.50							
Min. tear force, N	> 6 N							
Compatible delivery microcatheter	VASCO+10 or	VASCO+10 or	VASCO+10 or	VASCO+10 or	VASCO+18	VASCO+21	VASCO+21	VASCO+21 or

	any microcatheter with a min. inner diameter of 0.017" (0.41 mm)	any microcatheter with a min. inner diameter of 0.017" (0.41 mm)	any microcatheter with a min. inner diameter of 0.017" (0.41 mm)	any microcatheter with a min. inner diameter of 0.017" (0.41 mm)	or any microcatheter with a min. inner diameter of 0.021" (0.52 mm)	or any microcatheter with a min. inner diameter of 0.024" (0.60 mm)	or any microcatheter with a min. inner diameter of 0.024" (0.60 mm)	any microcatheter with a min. inner diameter of 0.024" (0.60 mm)
Intended use, max. vessel diameter, mm	3	3	3	3	4	6	6	6
Internal diameter of the introducer, mm	218							
External diameter of the introducer, mm	220							
Introducer length, mm	2250							
External diameter of the proximal part of the pusher, mm	14 mm in the widest part (see Fig. 6)							
Internal diameter of the proximal part of the pusher, mm	1.4 mm in the narrowest part (see Fig. 6)							
Radial force, N	0.15							
Capture device wire diameter, mm	0.09 – 0.17 mm							

Note: Unspecified tolerances are $\pm 1\%$. These design options of the CATCH+ Device for Thrombectomy of Cerebral Arteries have identical structures and differ only in size.

6.2 General requirements applicable to all devices:

Device appearance: The device surface has no foreign particles, bubbles or stains. The surface is smooth and free from technological and external defects and provides the minimum level of vascular injuries during use. The surface type is atraumatic, the presence of 3 distal radiopaque markers, without rupture of pushers. Fixation of the ring on the pusher.

Resistance to biological fluids: devices are resistant to biological fluids and body tissue liquids that they come into contact with during operation: resistant to aggressive biological fluids (blood).

Controllability: The same form before and after use.

Repeated insertion / withdrawal into a compatible microcatheter: 100% withdrawal. Good repeated insertion.

Deployment: Deploys fully. Uniform and symmetric deployment.


Compressive strength: During compression the diameter remains unchanged. After the expiration date the characteristics do not change.

Radial force: After the expiration date the characteristics do not change.

6.3 Description of subassemblies and parts:

The CATCH+ Device for Thrombectomy of Cerebral Arteries consists of a capture device with 3 radiopaque markers at the distal part of the device and its pusher (also with a radiopaque marker). The device is placed inside an introducer designed to protect the device from breakage and damage.

Figure 3. Drawing of the MD "CATCH+ Device for Thrombectomy of Cerebral Arteries" (models: CATCHMINI 3X5, CATCHMINI 3X10, CATCHMINI 3X15, CATCHMINI 3X20) with indication of the device components

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INDICE **E**

Figure 5. Drawing of the MD “CATCH+ Device for Thrombectomy of Cerebral Arteries” (model: CATCH 4x20) with indication of the device components



Figure 6. Drawing of the capture device of the MD "CATCH+ Device for Thrombectomy of Cerebral Arteries" (model: CATCH 4x20)

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Rep.: Nb: Matière: NITi

Dessinateur: K.Yappare
Date: 19/07/2018

Observation: Utilisation du tube NITi Ref: TN5.00X4.79-105

Vérificateur: E. Largen
Date: 19/07/2018

Ech: 5:1

Stent CATCH+4x20

Approbateur: R.Barthelemy
Date: 19/07/2018

Date d'Application: 19/07/2018

A4

10, rue Croix Vigneron 95160 Montmorency

Ref:

D810.01

INDICE G

Figure 7. Drawing of the MD "CATCH+ Device for Thrombectomy of Cerebral Arteries" (models: CATCHMAXI6x20, CATCH MAXI6x30, CATCH MAXI6x40) with indication of the device components



Figure 8. Drawing of the capture device of the MD "CATCH+ Device for Thrombectomy of Cerebral Arteries" (models: CATCHMAXI6x20, CATCH MAXI6x30, CATCH MAXI6x40)

Ref stent	Nb de cellule sur l'axe	L1	L2	L3
STCATCHMAXI6X20	4	31,85	16,8	11,2
STCATCHMAXI6X30	6	42,85	28	22,4
STCATCHMAXI6X40	8	51,10	39,25	33,6

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
Ref.: D810.09

INDICE E

Figure 9. Drawing of the component «Proximal part of the pusher»

7. Description of accessories and medical devices used with the subject device:

The CATCH+ Device for Thrombectomy of Cerebral Arteries should only be introduced using the microcatheter indicated on the label. CATCHMINI is compatible with a catheter having a min. inner diameter of 0.41 mm (VASCO+10). CATCH 4×20 is compatible with a catheter having a min. inner diameter of 0.52 mm (VASCO+18). CATCHMAXI is compatible with a catheter having a min. inner diameter of 0.60 mm (VASCO+21). All CATCH+ models have identical packages and sterilized with ethylene oxide gas.

Please refer to the instructions supplied with all interventional devices and materials to be used in combination with the CATCH+ Device for Thrombectomy of Cerebral Arteries for their intended use, contraindications and potential complications. Other accessories compatible with the device which are NOT supplied together with the CATCH+ should be selected based on experience and preferences of the physician:

- 5-6F balloon guide catheter (1.67-2 mm)
- Microcatheter (use only the microcatheter indicated on the label)
- 0.014" (0.35 mm) or 0.012" (0.30 mm) pusher
- Aspiration device
- Saline / heparin-saline continuous flushing set (syringe from 20 cc to 50 cc)
- Two rotating hemostatic valves (RHV)
- Infusion stand
- Femoral arterial lock

8. Information on the use of a human medicinal product, animal or human origin materials in the medical device:

Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not contain any medicinal products, pharmaceutical substances, animal or human origin materials.

9. Information on the procedure of installation, assembly, setting, and calibration required to commission the medical device:

The CATCH+ Device for Thrombectomy of Cerebral Arteries (hereinafter referred to as the Device, the CATCH+ Device) is a self-expandable system mounted on its own pusher and designed for catching and removing clots in patients. The device is intended for use to restore blood flow in the neurovasculature by removing clots in patients experiencing acute ischemic stroke due to a large cerebral artery occlusion. The Stentriever design makes it possible to catch a clot within cells and restore blood flow. The device is placed and deployed inside an occluded vessel within the thrombus formation and is then withdrawn in its deployed state. The mechanical thrombectomy procedure does not take more than 1 hour.

The appropriate pusher should be used to position the CATCH+ Device so that the effective length portion of the device would extend past each side of the thrombus in the vessel when fully deployed. If necessary, and if the device is not damaged, the same CATCH+ Device can be used for up to two flow restoration procedures in the same patient after cleaning the device with saline solution. The CATCH+ Device for Thrombectomy of Cerebral Arteries should not be used for more than two flow restoration procedures and more than three restoration attempts in the same vessel.

The Device and its pusher (also fitted with a radiopaque marker) are placed inside a PTFE (polytetrafluoroethylene) insertion tube and packaged in a peelable pouch and a cardboard box, together

with the instructions for use.

The CATCH+ Device is intended for single patient use only and is sterilized with ethylene oxide. In fact, these products are surgically invasive devices intended for transient use (less than 60 minutes) in contact with the central circulatory system (arteriae cerebrales).

9.1 Method of application:

Preparation and application procedure

Preparation

- Administer anticoagulant and antiplatelet agents in accordance with standard clinical practice guidelines.
- Use X-ray angiography to determine the location and size of the area which is subject to thrombectomy.
- Select the CATCH+ Device of the appropriate size.
- To achieve optimal performance of the CATCH+ Device and reduce the risk of thromboembolic complications, maintain continuous flushing between
 - a) the femoral arterial sheath and the guide catheter,
 - b) the microcatheter and the guide catheter,
 - c) the microcatheter and the guide catheter and the CATCH+ Device;
- Check all connections to make sure that no air enters the guide catheter or the microcatheter during continuous flushing.
- Position a suitable balloon guide catheter as close to the thrombus site as possible using a standard method. The balloon guide catheter should be appropriately sized to remove a clot if so required in the following steps.
- Connect a RHV to the guide catheter fitting, and then connect a tube to the continuous flushing system line.
- Select a microcatheter suitable for advancing the CATCH+ Device.
- Connect the second RHV to the microcatheter fitting, and then connect a tube to the continuous flushing system line.
- Set the flushing system line rate in accordance with standard clinical practice guidelines.
- Use an appropriate pusher to advance the microcatheter until the end of the microcatheter is positioned distal to the thrombus so that the usable length portion of the CATCH+ Device would extend past each side of the thrombus in the vessel when fully deployed. Tighten the RHV around the microcatheter.
- Delivery of the CATCH+ Device for Thrombectomy of Cerebral Arteries:
 - Partially insert the distal end of a pusher sheath into the RHV connected to the microcatheter. Tighten the RHV and make sure that fluid flows out of the proximal end of the pusher sheath.
 - Loosen the RHV and advance the pusher sheath until it is firmly seated in the microcatheter hub. Tighten the RHV around the pusher sheath to prevent reverse blood flow, but not so tight as to damage the CATCH+ Device during its introduction into the microcatheter. Make sure that there are no air bubbles trapped in the system.
 - Transfer the CATCH+ Device into the microcatheter by advancing the pusher in a smooth and continuous manner. As soon as the flexible portion of the pusher enters the microcatheter shaft, loosen the RHV and remove the pusher sheath over the proximal end of the pusher. Once completed, tighten the RHV around the pusher.

Leaving the pusher sheath in place may interrupt normal supply of flushing solution and allow reverse blood flow into the microcatheter.

- Visually verify that flushing solution is normally supplied. Once confirmed, loosen the RHV to advance the pusher.
- Carefully advance the CATCH+ Device under fluoroscopy until its distal markers align with the end of the microcatheter.

The CATCH+ Device for Thrombectomy of Cerebral Arteries should be positioned so that the usable length portion of the device would extend past each side of the thrombus in the vessel when fully deployed.

Note: *If excessive resistance force is observed during delivery of the CATCH+ Device for Thrombectomy of Cerebral Arteries, stop the delivery procedure and determine the cause of such resistance. Advancement of the CATCH+ Device for Thrombectomy of Cerebral Arteries against resistance may result in damage to the device and/or injury to the patient.*

Deployment of the CATCH+ Device for Thrombectomy of Cerebral Arteries:

- Loosen the RHV around the microcatheter. To deploy the CATCH+ Device for Thrombectomy of Cerebral Arteries, secure the pusher to maintain the device position and carefully withdraw the microcatheter in a proximal direction.
- Extend the microcatheter until it approaches the proximal marker of the CATCH+ Device. Tighten the RHV to prevent any movement of the pusher. The usable length of the deployed CATCH+ device should extend past each side of the thrombus.
- Tighten the RHV around the microcatheter. Use angiography to assess the thrombectomy status of the treated vessel.
- Adjunctive drug therapy can be used in accordance with standard clinical practice guidelines (thrombolytics, etc.).

Withdrawal of the CATCH+ Device for Thrombectomy of Cerebral Arteries:

- Inflate the guide catheter balloon to occlude the vessel as labeled on the balloon guide catheter.
- Make sure that the proximal marker on the CATCH+ Device is positioned inside the microcatheter. Loosen the RHV around the microcatheter to permit movement without loss of sealing.
- To remove a thrombus, slowly withdraw the microcatheter and the CATCH+ Device as a unit to the guide catheter tip while applying syringe aspiration to the guide catheter. Never advance the deployed CATCH+ Device in a distal direction.

Note: *Make sure that the microcatheter covers the proximal marker.*

Withdraw the CATCH+ Device and the microcatheter from the guide catheter while applying vigorous syringe aspiration to the guide catheter. Continue aspirating the guide catheter until the CATCH+ Device and the microcatheter are withdrawn from the guide catheter.

Note: *If withdrawal through the guide catheter is difficult, deflate the balloon and then withdraw simultaneously with the CATCH+ Device for Thrombectomy of Cerebral Arteries.*

Withdraw the guide catheter, microcatheter and CATCH+ Device as a whole through the sheath while maintaining aspiration. Remove the sheath if necessary.

Warning: *If excessive resistance force is observed during withdrawal of the CATCH+ Device, stop the procedure and determine the cause of such resistance.*

Do not perform more than three subsequent retrieval attempts in the same vessel using the CATCH+

Device.

- Open the guide catheter RHV so that the microcatheter and the CATCH+ Device could be removed without resistance. Be careful to avoid any interaction with the intervention site and to prevent air from entering the system.
- Aspirate the guide catheter to make it free from any thrombus material.
- Deflate the guide catheter balloon.
- ***If additional flow restoration attempts are required:***
- with a new CATCH+ Device for Thrombectomy of Cerebral Arteries, then repeat the steps described above starting with the Preparation section;
- with the same CATCH+ Device for Thrombectomy of Cerebral Arteries, then:
 - a. Clean the device with saline solution. Do not use solvents and do not autoclave.
 - b. Carefully inspect the device for damage. If there is any damage, discard the device and use a new CATCH+ Device for subsequent flow restoration attempts following the steps described above starting with the Preparation section. Use of the damaged device may result in additional damage to the device or injury to the patient.

Warning: Do not use each CATCH+ Device for Thrombectomy of Cerebral Arteries for more than two blood flow restoration procedures.

Repositioning of the CATCH+ Device for Thrombectomy of Cerebral Arteries:

If repositioning (e.g., redeployment) of the CATCH+ Device for Thrombectomy of Cerebral Arteries is necessary, proceed as follows:

- Loosen the RHV around the microcatheter and the pusher. Hold the pusher firmly in its position under fluoroscopy to prevent the CATCH+ Device from moving.
- Carefully reposition the CATCH+ Device for Thrombectomy of Cerebral Arteries by advancing the microcatheter over the CATCH+ device until the distal markers of the CATCH+ device appear at the end of the microcatheter.

Note: If significant resistance is observed during the repositioning process, immediately stop and proceed to the Thrombectomy Device Withdrawal section above.

Precautions for use

- Do not use if package is opened or damaged. The devices are sterile if package is undamaged.
- The CATCH+ Device for Thrombectomy of Cerebral Arteries is intended for single use only. Do not reuse. Reuse of the device may result in a high risk of contamination to the patient and performance degradation.
- Do not resterilize the devices.
- Store in a dry and dark place at room temperature.
- Do not use past the expiration date.
- These devices should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology.
- The catheter and CATCH+ Device should be used together with fluoroscopic control and appropriate anticoagulants.
- Follow the instructions for use provided with the corresponding systems and injectants.
- Be careful not to touch the basket during various operations, such as removal from package or insertion

into the guide catheter, as this can cause its displacement or damage.

- Do not advance the CATCH+ Device through the stent, since there is a significant risk of stent displacement.
- Perfuse the guide catheter throughout the entire procedure.

10. Requirements to the premises intended for installation of the medical device:

The CATCH+ Device for Thrombectomy of Cerebral Arteries should be used under the following conditions: Ambient temperature from 0 °C to +40 °C, relative humidity from 0 % to 90 % (non-condensing), and atmospheric pressure from 700 hPa to 1,060 hPa. All basic sterility requirements to the premises and medical personnel clothes should be met to avoid the risk of contamination to the patient.

11. Information for checking correctness of installation of the medical device and possibility of safe operation/use, including:

a) Scope and frequency of maintenance, including cleaning and disinfection of the medical device: The CATCH+ Device for Thrombectomy of Cerebral Arteries is not subject to maintenance, cleaning or disinfection.

b) The list of information, keys and access passwords required for assembly and maintenance of the medical device provided by the manufacturer: The CATCH+ Device for Thrombectomy of Cerebral Arteries does not contain any software, information, keys or access passwords required for assembly and is not subject to maintenance.

c) List of consumables, as well as procedure for use and replacement thereof: The CATCH+ Device for Thrombectomy of Cerebral Arteries does not contain any consumables.

d) Need in calibrating the medical device to ensure appropriate and safe operation throughout the service life: Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not require calibration to ensure appropriate and safe operation throughout the service life.

e) Methods of mitigating risks related to installation, calibration or maintenance of the medical device: To mitigate the risks related to installation of the CATCH+ Device for Thrombectomy of Cerebral Arteries, read the Instructions for Use carefully. For additional information please refer to Section 9 "Information on the Procedure for Installation, Assembly, Setting, and Calibration Required to Commission the Medical Device".

The device is intended for single use and does not require maintenance. Calibration is also not required, all the necessary settings are preset by the manufacturer.

f) Information on assembly, adjustment, setting, calibration and other actions required to commission the medical device and for correct use (operation) thereof: Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not require assembly, adjustment, setting, calibration and other actions to commission the medical device and for correct use (operation) thereof.

g) Information on the list of main specifications for the use of the medical device, transportation and storage conditions:

Transportation conditions: Ambient temperature from 0°C to +40°C, relative humidity from 0 % to 90 % (non-condensing), and atmospheric pressure from 700 hPa to 1,060 hPa, keep dry. The devices should be transported by air, rail, water or road transport in covered vehicles in accordance with the rules of transportation applicable to this type of transport.

Storage conditions: Store the device in its original manufacturer's package in a dry and clean place at controlled temperature. The shelf life is indicated on each device label. Ambient temperature from 0°C to 40°C, relative humidity from 0 % to 90 % (non-condensing), and atmospheric pressure from 700 hPa to 1,060 hPa.

Operation conditions: Ambient temperature from 0°C to +40°C, relative humidity from 0 % to 90 % (non-condensing), and atmospheric pressure from 700 hPa to 1,060 hPa. The device is operated in the human body at a temperature from + 32°C to + 42°C. The product is resistant to body fluids (see paragraph 6.2.).

h) List of the national standards used by the manufacturer of the medical device:

The medical device "CATCH+ Device for Thrombectomy of Cerebral Arteries" is a disposable sterile product that meets the requirements of standards and control methods envisaged by the standards of the Russian Federation:

GOST R 50444-92 "Devices, apparatuses and equipment for medical use. General specifications"

GOST EN 556-1-2011 "Sterilization of medical devices. Requirements for medical products in the category "sterile". Part 1. Requirements for medical devices subject to final sterilization"

GOST ISO 11607-1-2018 "Packaging for medical products subject to final sterilization. Part 1. Requirements for materials, barrier systems for sterilization and packaging systems"

GOST ISO 11607-2-2018 "Packaging for medical products subject to final sterilization. Part 2. Validation requirements for formation, sealing and assembly processes"

GOST R IEC 62366-2013 "Medical products. Design of medical devices with regard to operational usability"

GOST R ISO 15223-1-2014 "Medical products. Symbols used for marking on medical products, labels and accompanying documentation. Part 1. Basic requirements"

Requirements for the Instructions for Use and Technical Documentation

GOST R ISO 10993-2-2009 "Medical products. Evaluation of biological effects of medical devices. Part 2. Requirements for the treatment of animals"

GOST ISO 10993-1-2011 "Medical products. Evaluation of biological effects of medical devices. Part 1. Evaluation and Research"

GOST ISO 10993-4-2011 "Medical products. Evaluation of biological effects of medical devices. Part 4. Research of products associated with blood"

GOST ISO 10993-5-2011 "Medical products. Evaluation of biological effects of medical devices. Part 5. Cytotoxicity studies: in vitro methods"

GOST ISO 10993-7-2016 "Medical products. Evaluation of biological effects of medical devices. Part 7. Residual content of ethylene oxide after sterilization"

GOST ISO 10993-10-2011 "Medical products. Evaluation of biological effects of medical devices. Part 10. Research of irritating and sensitizing action"

GOST ISO 10993-11-2011 "Medical products. Evaluation of biological effects of medical devices. Part 11. Studies of general toxic effects"

GOST ISO 10993-12-2015 "Medical products. Evaluation of biological effects of medical devices. Part 12.

Preparation of samples and control samples"

GOST ISO 10993-13-2016 "Medical products. Evaluation of biological effects of medical devices. Part 13. Identification and quantification of degradation products of polymeric medical devices"

GOST ISO 10993-15-2011 "Medical products. Evaluation of biological effects of medical devices. Part 15. Identification and quantification of degradation products of metals and alloys"

GOST R 52770-2016 "Medical products. Safety requirements. Methods of Sanitary and Chemical and Toxicological Testing"

GOST 31214-2016 "Medical products. Requirements for samples and documentation submitted for toxicological, sanitary and chemical tests, tests for sterility and pyrogenicity"

GOST 31209-2003 "Containers for blood and its components. Requirements for chemical and biological safety and test methods "p. 5.3

State Pharmacopoeia XIII edition, OFS 1.2.4.0005.15 Pyrogenicity.

State Pharmacopoeia XIII edition, OFS 1.2.4.0003.15 Sterility.

12. Information on the sterile condition of the medical device, sterilization method and procedure in case of breach of sterile packaging: Each CATCH+ Device for Thrombectomy of Cerebral Arteries is supplied as a disposable sterile product. The CATCH+ Device for Thrombectomy of Cerebral Arteries is placed into a terminally sterilized package and is sterilized with ethylene oxide gas. The device should not be used and should be disposed of immediately if the package sterility is compromised. Resterilization of the device may result in a high risk of contamination to the patient and performance degradation.

13. Information on processing of the medical device for reuse: The device is supplied sterile and should not be subject to processing, disinfection or resterilization. After use, the device should not be reused and should be disposed of immediately. Reuse of the device may result in a high risk of contamination to the patient and performance degradation.

14. Information required to identify medical devices to obtain a safe combination and information on the known restrictions of concurrent use of medical devices:

Section 7 "Description of Accessories and Medical Devices Used with the Subject Device" contains the information on the devices which may be used with the CATCH+ Device for Thrombectomy of Cerebral Arteries. Please refer to the instructions for use for each BALT medical device for all the necessary information on the use and known restrictions.

15. Information on the nature, type, intensity and distribution of radiation during the use of the medical device:

Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not require RF energy and therefore does not create electromagnetic fields and interference induced thereby.

16. Information on the precautions taken if

A) The medical device is defective, malfunctions or features deviations of operation when these issues may affect safety of the medical device, including visual defects:

Please refer to Section 9 for the information on identifying potential failures, malfunctions or deviations in operation that may affect safety of the medical device, including those identified by means of visual inspection for the possibility of their elimination.

B) The medical device is exposed to external factors related to the use of the medical device with other

medical devices:

Refer to the instructions for use of the devices used with the CATCH+ Device for Thrombectomy of Cerebral Arteries (Section 7 of the Instructions for Use)

C) A risk of electromagnetic interference created by the medical device affecting other medical equipment: Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not require RF energy and therefore does not create electromagnetic fields and interference induced thereby.

17. Warnings and precautions taken by the consumers when using the medical device containing a human medicinal product, carcinogenic, mutagenic or toxic materials, possible extraction or leaching whereof results in sensitization, allergic reaction or negatively affects fertility: Not applicable. The CATCH+ Device for Thrombectomy of Cerebral Arteries does not contain any medicinal products, pharmaceutical substances, animal or human origin materials, carcinogenic, mutagenic or toxic materials, possible extraction or leaching whereof results in sensitization, allergic reaction or negatively affects fertility.

18. Warnings and precautions taken by the consumer when disposing of the medical device, accessories and consumables used therewith:

Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). Package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

If the product has been damaged before use, or the product has expired, it must be disposed of in the same manner and conditions.

19. Information on the circumstances requiring consultation of the consumer with a health care professional:


The device should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology. Conditions of use: Medical treatment facilities and healthcare centers, only when prescribed by a medical specialist. In view of the foregoing, it follows that in any case, the consumer should consult with a health care professional to determine the necessity to perform the procedure, and as a consequence, the use of the CATCH+ Device for Thrombectomy of Cerebral Arteries prior to the cerebral embolization.

20. Information on the first issue or the latest revision of the operating documentation:

Version No.	Revision details
Instructions for Use (version 1)	March 2015
Revision (version 2)	September 2017
Revision (version 3)	May 2018

21. Layout of the terminally sterilized package label and the individual cardboard box label for the CATCH+ Device for Thrombectomy of Cerebral Arteries for distribution in the Russian Federation:

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
EAN13:


3 700481 315538

Note: The label layout is applicable to all design options of the MD CATCH+ Device for Thrombectomy of Cerebral Arteries.

The label contains the following information:

- Logo of the manufacturer;
- Device name and model;
- Expiration date;
- EAN barcode;
- Quantity;
- Schematic drawing of the device, including explanatory text:
 - 3xORX (location of 3 distal radiopaque markers);
- Length of capture device (i.e. 31 mm);
- Information: Nominal (max.) vessel diameter;
- Information: Capture device in deployed state;
- Information: Device length;
- Information: Min. inner diameter of the microcatheter;
- Information: Compatible delivery microcatheter;
- Warning: "Do not overtorque the CATCH+ Device for Thrombectomy of Cerebral Arteries"
- Information on the manufacturer;


As well as the following conventional symbols

	Catalogue number
	Lot number
	Use by
	Date of manufacture
	CE mark
	Sterilized with ethylene oxide
	Caution! Consult instructions for use
	Do not reuse
	Do not use if package is damaged
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	European Article Number

	Manufacturer
	Do not re sterilize
	Fragile
	QR code
	Apyrogenic
НЕ ТОКСИЧНО	-
РУ № _____	Registration Certificate Number

Note: The markings applied to the pouch and the box are identical.

Appearance of the box:

Cardboard packages with devices are placed into a corrugated cardboard shipper box.
One shipper box contains max. 20 cardboard packages.

The shipper box label contains the following information:

- Device name and model;

- Information on the manufacturer;
- Quantity;
- Lot number;
- Batch number;
- Date of manufacture;
- Expiration date;
- "Keep away from sunlight" symbol;
- "Keep dry" symbol.

22. **Procedure and conditions of disposal or destruction of the medical device:** Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). Package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

23. **Guarantees**

Shelf life – 5 years. Sterilization validity – 5 years.

Address for complaints:

Center for Perinatal Medicine, Limited Liability Company

Legal/actual address: 20/4 Tovarishchesky per., 109004, Moscow

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