COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, MICHELLE MIRZAEI, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Onyx Liquid Embolic JFu (Onyx 3+L)
Description of documents

Is a true and exact copy of the correct and complete original document.

Michelle Mirzaei
Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California

County of ORANGE

Subscribed and sworn to (or affirmed) before me on this 9th day of September,
20 19 by MICHELLE MIRZAEI _____, proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

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This page is intentionally left blank.

TABLE OF CONTENTS

Onyx™ Liquid Embolic System

English.....	4
Français.....	5
Deutsch.....	7
Italiano.....	9
Español.....	11
Svenska.....	13
Nederlands.....	15
Português.....	17
Suomi.....	19
Dansk.....	21
Ελληνικά.....	23
Čeština.....	25
Magyar.....	27
Русский.....	29
Polski.....	31
Türkçe.....	33
Norsk.....	35
Slovenčina.....	37
Română.....	39
Български.....	41
한국어.....	43
عربي.....	47
Symbol Glossary.....	49

English Instructions for Use

EN

Onyx™ Liquid Embolic System

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

DESCRIPTION

Onyx™ LES is a non-adhesive liquid embolic agent comprised of EVOH (ethylene vinyl alcohol) copolymer dissolved in DMSO (dimethyl sulfoxide), and suspended micronized tantalum powder to provide contrast for visualization under fluoroscopy. The Onyx™ Liquid Embolic System (LES™) consists of a 1.5 ml vial of Onyx™ LES, a 1.5 ml vial of DMSO, two 1 ml Onyx™ LES delivery syringes and one 1 ml DMSO syringe. A DMSO compatible delivery micro catheter that is indicated for use in the neuro vasculature is used to access the embolization site. The Onyx™ LES Syringe-Catheter Interface adapter, which is sold separately, is an accessory to the Onyx™ LES system that provides an interface between the Onyx™ LES delivery syringe and the delivery catheter.

PRINCIPLE OF OPERATION

Onyx™ LES is delivered through a micro catheter into the AVM under fluoroscopic control. The DMSO solvent dissolves into the blood and interstitial fluids, causing the EVOH copolymer and suspended tantalum to precipitate in situ into a spongy, coherent embolus. Onyx™ LES immediately forms a skin as the polymeric embolus solidifies from the outside to the inside, while traveling more distally in the vessel. Since Onyx™ LES is non-adhesive, the micro catheter can be left in place while slow, controlled injections are performed. Post embolization angiography can be conducted with the delivery micro catheter in place, enabling the physician to make additional injections through the same micro catheter, if necessary.

INDICATIONS FOR USE

Embolization of lesions in the peripheral and neurovasculature, including arteriovenous malformations and hypervascular tumors.

HOW SUPPLIED

Onyx™ LES is available in three product formulations, Onyx™-18 LES (6% EVOH), Onyx™-20 LES (6.5% EVOH) and Onyx™-34 LES (8% EVOH):

- Onyx™-18 LES and Onyx™-20 LES: Recommended when feeding pedicle injections will be conducted close to the nidus.
- Onyx™-34 LES: Recommended for embolizing higher flow and larger fistulous components.

CONTRAINDICATIONS

Onyx™ LES is not indicated for use with premature infants (<1,500 g) or individuals with significant liver function impairment.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Hematoma.
- Arterial thrombosis.
- Ischemic events due to embolic migration, vasospasm, thrombosis.
- Catheter entrapment.
- Catheter rupture.

- Device migration and cast movement.
- Hemorrhagic accidents: vascular rupture/perforation.
- Hemodynamic changes induced by the embolization may result in hemorrhagic complications.
- Hemorrhagic complications related to attempts to remove entrapped catheter.
- These ischemic or hemorrhagic complications may result in various functional neurological deficits and possibly death.

WARNINGS

Due to the possibility of electrical arcing with the tantalum metal in the Onyx™ LES material, use of monopolar electrocautery devices for surgical resection of bAVMs or arteriovenous fistula embolized with Onyx™ LES should be avoided. Bipolar devices should be used with caution.

This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

- Use only the Onyx™ LES syringes to inject DMSO and Onyx™ LES. Other syringes may not be compatible with DMSO.
- Use only ev3 DMSO compatible micro catheters. Other micro catheters may not be compatible with DMSO and their use can result in thromboembolic events due to catheter degradation.
- Always consider the potential for Onyx™ LES to interact with other embolic agents, e.g., cyanoacrylates, coated coils, particles and/or embolic spheres.
- Performing embolization to occlude blood vessels is a high-risk procedure. The procedure should be carried out by a specialist with the appropriate neuro or peripheral interventional training, and a thorough knowledge of the pathology to be treated, angiographic techniques, and super-selective embolization.
- The patient's anatomic anatomy must be amenable to micro catheter tip placement at a location that is distal to arterial branches that potentially may feed a cranial nerve.
- Adjunctive coil use should be considered if angiography shows that venous drainage of the AVM appears almost simultaneously with arterial opacification.
- Premature solidification of Onyx™ LES may occur if micro catheter luer contacts saline, blood or contrast of any amount.
- Failure to continuously mix Onyx™ LES for the required time may result in inadequate suspension of the tantalum, resulting in inadequate fluoroscopic visualization during delivery. Inject Onyx™ LES immediately after mixing. If Onyx™ LES injection is delayed, tantalum settling can occur within the syringe resulting in poor visualization of Onyx™ LES during injection.
- Using the Onyx™ LES Syringe-Catheter Interface Adapter will reduce labeled deadspace of the micro catheter. Failure to make the appropriate allowances may result in unintended embolization.
- Reference the Onyx Syringe-Catheter Interface Adapter Instructions for Use for a list of compatible micro catheters and deadspace information.

When injecting Onyx™ LES, fluoroscopic visualization should reveal Onyx™ LES traveling through the catheter lumen. It is recommended to obtain fluoroscopic imaging prior to reaching the minimum catheter + adapter deadspace in order to visualize the embolic material before it exits the tip of the catheter.

- Only use thumb pressure to inject Onyx™ LES. Using palm of hand to advance plunger may result in catheter rupture due to over pressurization in the event of catheter occlusion.
- Inject Onyx™ LES and DMSO at a slow, steady rate, not to exceed 0.3 ml/minute. Animal studies have shown that a rapid injection of DMSO into the vasculature may lead to vasospasm and/or angioneurosis.

- Do not allow more than 1 cm of Onyx™ LES to reflux back over catheter tip. Angioarchitecture, vasospasm, excessive Onyx reflux, or prolonged injection time may result in difficult catheter removal and catheter entrapment. Excessive force to remove an entrapped catheter may cause serious intracranial hemorrhage. The long term effects of an entrapped catheter that is left in a patient is unknown, but could potentially include clot formation, infection or catheter migration.
- DO NOT interrupt Onyx™ LES injection for longer than two minutes prior to re-injection. Solidification of Onyx™ LES may occur at the catheter tip resulting in catheter occlusion, and use of excessive pressure to clear the catheter may result in catheter rupture.
- Adequate fluoroscopic visualization must be maintained during Onyx™ LES delivery or non-target vessel embolization may result. If visualization is lost at any time during the embolization procedure, HALT Onyx™ LES delivery until adequate visualization is re-established. Testing has shown that over-pressurization and rupture can occur if 0.05 ml of Onyx™ LES is injected and is not visualized exiting the catheter tip.
- Increased resistance to Onyx™ LES injection – Stop injection. Do not attempt to clear or overcome resistance by applying increased injection pressure. If this occurs, determine the cause of resistance (for example, Onyx™ LES occlusion in catheter lumen), and replace catheter. Use of excessive pressure may result in catheter rupture and embolization of unintended areas.
- After using a micro catheter with Onyx™ LES, do not attempt to clear the micro catheter or to inject any material through it. Attempts to clear catheter may lead to embolus or embolization of unintended area.
- If Onyx™ LES escapes outside the vascular space, as might occur if the vessel wall is compromised, a subacute inflammatory response to the material may occur.

PRECAUTIONS

Inspect product packaging prior to use. Do not use if sterile barrier is opened or damaged.

Use prior to expiration date.

Review both the catheter and interface adapter Instructions for Use prior to using Onyx™ LES.

Verify that the catheters and accessories used in direct contact with the Onyx™ LES polymer are clean and compatible with the material and do not trigger precipitation or degrade with contact. Refer to the respective Warnings and Directions for Use sections.

Failure to wait a few seconds to retrieve micro catheter after Onyx™ LES injection may result in fragmentation of Onyx™ LES into non-target vessels.

Difficult catheter removal or catheter entrapment may be caused by one or more of the following factors:

- Long catheterization time
- Angio-architecture, very distal arteriovenous malformation fed by afferent, lengthened, small, or tortuous pedicles
- Vasospasm
- Reflux
- Injection time

To reduce the risk of catheter entrapment, carefully select catheter placement and manage reflux to minimize the factors listed above.

Onyx™ LES reflux along the distal tip of the micro catheter: Apart from the risk of ischemic complications due to unintended embolization, significant reflux may result in entrapment of the micro catheter causing difficult removal. The amount of reflux that can be accepted must always be compared to the angio-architecture of

the malformation to minimize risk of unintended embolization or difficult catheter removal. In general, minimize the reflux to less than 1 cm along the distal tip of the micro catheter. All other factors may affect this limit.

- Should catheter removal become difficult, the following technique may assist in catheter retrieval:
 - Carefully pull the catheter to assess any resistance to removal.
 - If resistance is felt, remove any "slack" in the catheter.
 - Gently apply traction to the catheter (approximately 3-4 cm of stretch to the catheter).
 - Hold this traction for a few seconds and release. Assess traction on vasculature to minimize risk of hemorrhage.
 - This process can be repeated intermittently until catheter is retrieved.
 - When using ev3 micro catheters, do not apply more than 20 cm of traction to catheter, to minimize risk of catheter separation.
- For entrapped catheters:
 - Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system, rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.
 - This is accomplished by stretching the catheter and cutting the shaft near the entry point of vascular access allowing the catheter to remain in the artery.
 - If catheter breaks during removal, distal migration or coiling of the catheter may occur.
 - Same day surgical resection should be considered to minimize risk of thrombosis.

STORAGE

Store Onyx™ LES and DMSO between -20° and 55°C. Prior to use, maintain product temperature between 19° and 24°C. If product freezes due to exposure to colder temperatures, thaw at room temperature before use.

DIRECTIONS FOR USE

1. Shake Onyx™ LES at least 20 minutes on an Onyx™ LES mixer* at a setting of 8. Continue mixing until ready to inject Onyx™ LES per step 7.
 2. Confirm micro catheter placement with injection of contrast agent per institutional procedure.
 3. Flush contrast from micro catheter hub with 10 ml of saline. Leave the syringe connected.
 4. Filling catheter deadspace: aspirate approximately 0.8 ml of ev3 sterile DMSO into an Onyx™ LES 1 ml syringe. Inject the DMSO into the delivery micro catheter in sufficient volume to fill the catheter deadspace. Refer to delivery catheter labeling for deadspace volume.
 5. When filling the syringe with Onyx™ LES, aspirate approximately 1.1 ml of Onyx™ LES so that the piston head is beyond the 1 ml graduation line on the syringe barrel.
- NOTE:** This extra volume of Onyx™ LES will be used to fill (prime) the interface adapter.
6. Attach the interface adapter to the syringe.
 7. Inject Onyx™ LES through the interface adapter until the distal end of the piston head is even with the 1 ml graduation line on the syringe barrel.
 8. Any excess Onyx™ LES on the tip of the interface adapter may be wiped off with a clean, dry particle-free cloth.
 9. Remove the DMSO syringe from the catheter, overfill and wash the luer hub with the balance of the DMSO.

10. Immediately connect the Onyx™ LES syringe firmly to the micro catheter hub, making sure there is no air in the hub during the connection.
11. Inject Onyx™ LES to displace DMSO. Based on clinical practice, it is recommended that Onyx™ LES be injected at a steady rate of 0.16 ml/minute (0.25 ml/90 seconds). Do not exceed 0.3 ml/minute.

WARNINGS

Failure to continuously mix Onyx™ LES for the required time (20 minutes) may result in inadequate fluoroscopic visualization during delivery.

Inject Onyx™ LES immediately after mixing. Tantalum settling can occur within the syringe resulting in poor visualization of Onyx™ LES during injection.

Use thumb pressure only to inject Onyx™ LES. Using palm of hand to advance plunger may result in catheter rupture due to over-pressurization in the event of catheter occlusion.

Stop injection if Onyx™ LES is not visualized exiting catheter tip. If the catheter becomes occluded, over-pressurization can occur. During Onyx™ LES injection, continuously verify that Onyx™ LES is exiting the catheter tip. Testing has shown that over-pressurization and rupture can occur if 0.05 ml of Onyx™ LES is injected and is not visualized exiting the catheter tip.

Stop injection if increased syringe resistance is felt. Increased injection force may be due to catheter occlusion. Continued injection into an occluded catheter will result in catheter rupture.

Do not interrupt Onyx™ LES injection for longer than two minutes prior to re-injection. Solidification of Onyx™ LES may occur at the catheter tip resulting in catheter occlusion, and use of excessive pressure to clear the catheter, may result in catheter rupture.

Stop injection if abnormal dilation of the micro catheter diameter is observed. The increase in micro catheter diameter can signify a blocked catheter. The increase in micro catheter diameter can best be determined by double incidence fluoroscopy at low magnification or with "road mapping" procedure, examining the longest possible section of the micro catheter.

Adequate fluoroscopic visualization must be maintained during Onyx™ LES delivery or non-target vessel embolization may result. If visualization is lost at any time during the embolization procedure, HALT Onyx™ LES delivery until adequate visualization is re-established.

12. If a second Onyx™ LES syringe is required for an embolization, do not remove the adapter from the catheter. When the new Onyx™ LES syringe is ready, simply remove the empty syringe from the adapter's proximal end and connect the new syringe to the adapter.
13. Upon completion of Onyx™ LES injection, wait a few seconds, slightly aspirate syringe, and then gently pull the catheter to separate the Onyx™ LES cast.

The following is an optional catheter retrieval technique:

- Upon completion of Onyx™ LES injection, wait a few seconds and then slightly aspirate.
- Remove all slack from the catheter by putting a few centimeters of traction on the catheter to create a slight tension in the catheter system.
- Firmly hold the catheter and then pull it using a quick wrist snap motion (from left to right) 10 – 15 cm to remove the catheter from the Onyx™ LES cast.

NOTE: It is not necessary to pull the catheter any more than 20 cm.

NOTE: This technique has not been qualified for the Rebar catheters.

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* Scientific Industries Genie 2, Model No(s): 120V SI-0240, 240V SI-0251, Vial Attachment No. DA-0570-010

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Michelle Mirzaei

Signature of document custodian (affiant)

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
State of California

County of ORANGE

Subscribed and sworn to (or affirmed) before me on this 9th day of September,
20 19 by MICHELLE MIRZAEI, proved

to me on the basis of satisfactory evidence to be the person~~(s)~~ who appeared before me.

Steve Oberlander
Signature (Seal)

INSTRUCTIONS FOR USE

Onyx™ Syringe-Catheter Interface Adapter

70739-001 Rev. 02/13

Mull Hall 09sep19

TABLE OF CONTENTS

Onyx™ Syringe-Catheter Interface Adapter

English	3
Français	4
Deutsch	6
Italiano	8
Español	10
Svenska	12
Nederlands	14
Português	16
Suomi	18
Dansk	19
Ελληνικά	21
Česky	23
Magyar	24
Русский	26
Polski	28
Türkçe	30
Norsk	32
Slovenčina	33
Română	35
Български	37
한국어	39
عربي	42
Symbol Glossary	43

English Instructions for Use

EN

Onyx™ Syringe-Catheter Interface Adapter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

DESCRIPTION

This device is an Onyx™ LES and DMSO compatible adapter used to provide an interface between an MTI 1 ml syringe and the 1.5F UltraFlow HPC™ micro catheter or Marathon™ micro catheter during an Onyx™ LES AVM embolization. The proximal end of the adapter incorporates a standard ISO female luer design to facilitate connection to the syringe. The distal end is designed specifically to fit the hub of the 1.5F UltraFlow HPC™ micro catheter or Marathon™ micro catheter.

PRINCIPLE OF OPERATION

This device, attached to an MTI 1 ml syringe filled with Onyx™ LES, will effectively reduce the total volume of DMSO in the hub of the 1.5F UltraFlow HPC™ micro catheter or Marathon™ micro catheter during connection, thereby minimizing the total amount of DMSO that may be injected. The adapter reduces the labeled dead space of the UltraFlow™ micro catheter or Marathon™ micro catheter by approximately 0.04 ml (see dead space table below). In addition, because the adapter eliminates mixing of Onyx™ LES and DMSO in the hub of the catheter during connection, the adapter helps maintain a distinct interface between concentrated Onyx™ LES and DMSO during the injection.

Dead Space Conversion Table

Catheter	Minimum Catheter Dead Space	
	Catheter	Catheter + Onyx™ Adapter
1.5F UltraFlow HPC™ (Model: No. 105-S065/66)	0.26 ml	0.22 ml
Marathon™ (Models No. 105-S055)	0.23 ml	0.19 ml

INDICATIONS FOR USE

The Micro Therapeutics' Onyx™ Syringe-Catheter Interface is intended for use as an accessory to the Onyx™ Liquid Embolic System (LES) 1.5F UltraFlow HPC micro catheter and Marathon™ micro catheter during Onyx™ LES injection.

CONTRAINDICATIONS

- Not for use with Onyx™ HD 500 LES.

- Not for use with any embolic material other than Onyx™ LES for the treatment of AVMs.
- Not for use with Rebar™ micro catheter or any catheter other than the 1.5F UltraFlow HPC™ micro catheter and Marathon™ micro catheter.
- Not for use during DMSO injection.
- Not for use during angiography or saline injection.
- Not for use with guidewires or introducers.

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- Onyx™ LES leakage at adapter connections.
- Premature Onyx™ LES solidification in the adapter.
- Increased delivery pressure due to damage or obstruction in the adapter.
- Catheter overpressurization and/or catheter rupture.
- Unintended embolization.

WARNINGS

- Use only with the MTI 1 ml Onyx™ LES (white) syringe to inject Onyx™ LES. Other syringes may not be compatible with DMSO.
- Use only with AVM Onyx™ LES formulations. Do not use with Onyx™ HD 500 LES.
- Use only with the MTI 1.5F UltraFlow HPC™ micro catheter or Marathon™ micro catheter. Other micro catheters may not be compatible with DMSO and are not compatible with the interface adapter.
- After attaching the adapter to the syringe, failure to fill (prime) the adapter with Onyx™ LES prior to connecting it to the catheter may result in introduction of air into the catheter.
- Premature solidification of Onyx™ LES may occur if the adapter contacts saline, blood or contrast of any amount.
- The adapter will reduce the UltraFlow™ micro catheter or Marathon™ micro catheter labeled dead space by approximately 0.04 ml (See the Dead Space Conversion Table in the Principle of Operations section). Failure to make the appropriate allowance may result in unintended embolization.
- Increase resistance to Onyx™ LES injection – Stop injection. Do not attempt to clear or overcome resistance by applying increased injection pressure. If this occurs, determine the cause of resistance (for example, Onyx™ LES occlusion in interface adapter), and replace adapter. Use of excessive pressure may result in catheter rupture and/or embolization of unintended areas.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Review both the Onyx™ LES and micro catheter instructions for use prior to using this device.
- Inspect product packaging prior to use. Do not use if sterile barrier is opened or damaged.
- Use prior to expiration date.
- Verify that the adapter is clean and free from any contamination including saline, blood and/or contrast media prior to connecting it to an Onyx™ LES filled syringe.
- Review the dead space table in the Principle of Operations section in order to make allowance for the catheter dead space reduction provided by the adapter prior to beginning an Onyx™ LES injection.

DIRECTIONS FOR USE

1. Prepare the Onyx™ LES according to the Onyx™ LES Instructions for Use.
2. When filling the syringe with Onyx™ LES, aspirate approximately 1.1 ml of Onyx™ LES so that the piston head is beyond the 1 ml graduation line on the syringe barrel.
Note: This extra volume of Onyx™ LES will be used to fill (prime) the interface adapter.
3. Attach the interface adapter to the syringe.
4. Inject Onyx™ LES through the interface adapter until the distal end of the piston head is even with the 1 ml graduation line on the syringe barrel.
5. Any excess Onyx™ LES on the tip of the interface adapter may be wiped off with a clean, dry particle free cloth.
6. Ensure that the infusion catheter has been previously filled (prime) with DMSO according to the Onyx™ LES Instructions for Use.
7. Remove the DMSO syringe from the catheter, overfill and wash the luer hub with the balance of the DMSO.
8. Immediately connect the Onyx™ LES syringe firmly to the micro catheter hub, making sure there is no air in the hub during the connection.
9. Inject Onyx™ LES to displace DMSO. Based on clinical practice, it is recommended that Onyx™ LES be injected at a steady rate of 0.16 ml/minute (0.25 ml/90 seconds). Do not exceed 0.3 ml/minute.
10. Proceed with the Onyx™ LES injection according to the Onyx™ LES Instructions for Use.
11. If a second Onyx™ LES syringe is required for an embolization, do not remove the adapter from the catheter. When the new Onyx™ LES syringe is ready, simply remove the empty syringe from the adapter's proximal end and connect the new syringe to the adapter.

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CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

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On August 7, 2020

Date

before me, Hayley George - Notary Public

Here Insert Name and Title of the Officer

personally appeared Michelle Nushine Mirzaee

Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature

Hayley George
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Appendix for the Instruction

Document Date: _____ Number of Pages: 6

Signer(s) Other Than Named Above: None

Capacity(ies) Claimed by Signer(s)

Signer's Name: Michelle Nushine Mirzaee

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© Scientific Industries GmbH, Model N. 120V SI 0240, 240V SI 0251, Vial Attachment No. DA-0570-010