Instruction for use

MEDICAL PRODUCTS

Dental Implants and Accessories

Manufacturer: **Green Implant System Technologie GmbH & Co. KG**,
Petershäuser Weg 162 D - 32425 Minden Germany, Tel.: +49 (0) 571-3921640,
Fax: +49 (0) 571-3921640, E-mail: info@greenimplant.de

Postanschrift:
Petershäuserweg 162, 32425 Minden
Tel: +49 (0) 571 3921640
Fax: +49 (0) 571 3921639

Bankverbindung:
Stadtsparkasse Porta Westfalica
BLZ: 49051990
Konto-Nr. 0000019166

Steuer-Nr.
335/5725/5258
BIC: WELADED1PWF
IBAN: DE 15 4905 1990 0000 0191 66

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BIC: WELADED1PWF

IBAN: DE 15 4905 1990 0000 0191 66

17. Wide straight titanium abutment;
18. Plastic temporary straight anatomical abutment;
19. Plastic straight abutment;
20. Narrow straight titanium abutment;
21. Zirconium Abutment;
22. Straight titanium abutments for turning;
23. Angled Titanium Abutment 15 °& 25°;
24. Anatomic titanium abutment angled 15 °& 25; °
25. Angular aesthetic anatomical titanium abutment 15 °& 25 °;
26. Angular Titanium Abutment .Slim Platform 15°& 25°;
27. Plastic temporary angular anatomical abutment 15 °& 25 °;
28. Plastic angled abutment 15 °&25 °;
29. Angular milling titanium abutment for 15 °& 25 °
30. Wide angular titanium abutment for milling 15 °& 25 °
31. Angular Zirconium Abutment 15 °& 25 °;
32. Locator abutment;
33. Abutment with ball attachment;
34. Esthetic screw abutment -connector;
35. Screw abutment - connector;
36. Esthetic abutment screw;
37. Plastic Abutment for casting;
38. Screw titanium abutment for Immediate Loading;
39. Angular titanium temporary abutment for Immediate Loading of 15 °& 25 °;
40. Fixation Screw for aesthetic screwed abutment;
41. Surgical screw cap for the implant;
42. Screw surgical Short
43. Surgical screw.
44. Long surgical screw.
45. Short screw to secure the abutment;
46. Long screw to secure the abutment;
47. Transfer implant for open tray;
48. Transfer implant for closed tray;
49. Transfer Abutment to open tray;
50. Transfer Abutment for closed tray;
51. Transfer implant clip for closed tray;
52. Transfer implant Slim Platform
53. Implant analog

Postanschrift:

Petershagerweg 162, 32425 Minden

Tel: +49 (0) 571 3921640

Fax: +49 (0) 571 3921639

Bankverbindung:

Stadtsparkasse Porta Westfalica

BLZ: 49051990

Konto-Nr. 0000019166

Steuer-Nr:

335/5725/5258

BIC: WELADED1PWF

IBAN: DE 15 4905 1990 0000 0191 66

54. Abutment analog
55. Titanium angled adapter 17 ° & 30 °;
56. Protective plastic cap;
57. Prosthetic accessories for overdenture restoration
58. Removable clip attachment system Removable Teeth.
59. Nylon cap.
60. Titanium cap.
61. Casting titanium abutment.
62. Molded plastic simulated abutment.
63. The plastic simulated cap for aesthetic abutment/connector.
64. The Plastic simulated cap for abutment/connector.
65. Simulated Plastic cap for the implant
66. The plastic simulated cap for temporary crowns.
67. Silicone insert (cap) is soft, standard and ekstra soft;
68. Pin parallelism
69. Kit for implant prosthetics: abutment screw abutment transfer, simulated plastic cap, plastic cap for a temporary crown, the connector for the transfer, the connector for the plastic cap temporary crown, similar abutment or implant.
70. Kit for implant prosthetics: screw abutment, transfer, simulated plastic cap, plastic cap for a temporary crown, the connector for the transfer, the connector for the plastic cap temporary crown, or implant abutment analog, titanium screw, nylon cap, protective disc.
71. Kit for implant prosthetics: attachment abutment, transfer, simulated plastic cap, plastic cap for a temporary crown, the connector for the transfer, the connector for the plastic cap temporary crown, analog attachment abutment or implant titanium screw nylon cap, protective disc.

Postanschrift:

Weslöhnerweg 162, 32425 Minden

Telefon: +49 (0) 571 3921640

Telefax: +49 (0) 571 3921639

Bankverbindung:

Stadtsparkasse Porta Westfalica

BLZ: 49051990

Konto-Nr. 0000019166

Steuer-Nr.

335/5725/5258

BIC: WELADED1PWF

IBAN: DE 15 4905 1990 0000 0191 66

Instructions for use

Product Description

Implant System Technologie GmbH & Co. KG" developed products : dental implants , which are different : diameter , length and size of the hexagon 2.1mm or 2.43mm , abutments , superstructure and tools. made of titanium alloy, of medical G- 5

surface of the implant has a double treatment (sandblasting and acid treatment) .

While the lower surface treatment unit implant processed (sandblasting) and larger coarse grains of aluminum oxide to produce a rougher surface texture and deep . The top module is processed (sandblasting) fine not coarse grains of aluminum oxide , in order not to breach the form and function of the blades pitch mikrozuby .

1. Normative documents regulating the quality and safety

1. Certificate of compliance with ISO 9001:2008;
2. Certificate of compliance with ISO 13485:2003, EN ISO 13485:2012;
3. EC certificate of conformity to international standards.
4. Declaration of conformity with international standards.

2. Intended use :

Implants "Implant System Technologie GmbH & Co. KG" is intended (s) for osseointegration in the jawbone and further used as a support for the restoration of the tooth (s) different types (removable and fixed) dental prosthetic designs on the upper and lower jaw. Used in all types of bone (D1-D4). Almost solves almost all clinical cases . The best results when implanted in all divisions of the upper and lower jaw.

4. Statement

Dental implants "Implant System Technologie GmbH & Co. KG" are used in full or partial edentulous mandible or maxilla . Dental implants are the backbone for functional and aesthetic restoration of missing teeth in upper or lower jaw . Recovery may include the installation of crowns, bridges , full and partial dentures , which are connected by means of an implant abutments .

Dental implants 3.0mm diameter and 3.3 mm have a lower mechanical strength than the implants with a diameter of 3.75 mm and 4.2 mm , etc. Such implants should only be used according to indications , implying a minimum load . They can not be used to restore a single tooth , especially in the sides of the upper and lower jaw.

5. Technical characteristics implants:

"Green Implant System Technologie GmbH & Co." Includes eleven types of implants that are adapted to different types of implant in bone tissue at various kinds of its capacity .

1. Green Classic Titanium Implant: Screw , samonorezayuschy , antiratotsionny conical implant hexagon 2.43mmd , variable double-helical thread 2x0.6mm (1.2mm) and a variable four spiral mikrozuby 4x0.3 mm (1.2mm) . Internal thread implant 1-72 " under the abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameter, mm : 3.3, 3.75 , 4.2, 5.0, 6.0
Length, mm : 6.8, 10 , 11.5, 13, 16 .

2 . Green Universal Titanium Implant - Screw , samonorezayuschy , anti-rotation tapered implant with internal hexagon and 2.1mmd 2.43mmd variable double stranded thread 2x1.1mm (2.2mm) and a variable chetyrehspiralnoy mikrorezboy 4x0.4 mm (1.6mm) . Internal thread implant 1-72 " under the abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameter, mm : 3.0 , 3.3, 3.75 , 4.2, 5.0, 6.0
Length, mm : 8, 10 , 11.5, 13, 16 , 18

3 . Green Respect Titanium Implant Screw , samonorezayuschy , antiratotsionny tapered implant with internal hexagon and 2.1mmd 2.43mmd variable double stranded thread 2x1.1mm (2.2mm) and a variable chetyrehspiralnoy mikrorezboy 4x0.55 mm (2.2mm) . Coronary narrowed portion of the implant in the form of inverted cone diameter or 3.0mmd 3.8mmd . Internal thread implant 1-72 " under the abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameters in mm: 3.75 , 4.2, 5.0, 6.0
Length in mm: 10 , 11.5, 13, 16 .

4 . Green Classic Plus Titanium Implant Modular , screw , samonorezayuschy , antiratotsionny tapered implant with internal hexagon and 2.1mmd 2.43mmd variable double stranded thread 2x0.6mm (1.2mm) and a variable chetyrehspiralnoy mikrorezboy 4x0.4 mm (1.6mm) . Internal thread implant 1-72 " under the abutment screw . The implant consists of an upper and lower base removable modules. Removable module is connected to the base by means shestigranika 2.43md , inside the castle shestigranika and abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameters in mm: 3.75 , 4.2, 5.0, 6.0
Length, 8 , 10, 11.5 , 13, 16 .

5 . Green Universal Plus Titanium Implant Modular , screw , samonorezayuschy , antiratotsionny tapered implant with internal hexagon and 2.1mmd 2.43mmd variable double stranded thread 2x1.1mm (2.2mm) and a variable chetyrehspiralnoy mikrorezboy 4x0.4 mm (1.6mm) . Internal thread implant 1-72 " under the abutment screw . Modular implant consists of : lower and upper basic plug-in modules . Removable module is connected to the base by means shestigranika 2.43md , inside the castle shestigranika and abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameters in mm: 3.75 , 4.2, 5.0, 6.0
Length, 8 , 10, 11.5 , 13, 16 .

6. Green Respect Plus Titanium Implant Modular , screw , samonorezayuschy , antiratotsionny tapered implant with internal hexagon and 2.1mm 2.43mm variable double stranded thread 2x1.1mm (2.2mm) and a variable chetyrehspiralnoy mikrorezboy 4x0.55 mm (2.2mm) . Internal thread implant 1-72 " under the abutment screw . The implant consists of an upper and lower base removable modules. Removable module is connected to the base by means shestigranika 2.43mm , inside the castle shestigranika and abutment screw . Coronary narrowed portion of the implant in the form of inverted cone diameter of 3.0mm 3.8mm . Double surface treatment: sand , aluminum oxide and etching.
Diameters in mm: 3.75 , 4.2, 5.0, 6.0
Length, 8 , 10, 11.5 , 13, 16 .

7. Green Prestige Implant Kit $\varnothing$ 3.75, 4.2, 5.0, 6.0mm, length 10 , 11.5, 13, 16 mm

Kit for dental implantation , comprising:

- Green Universal Titanium Implant

- titanium screw to close the internal threads of the implant with a two-stage implantation

Titanium screw - healing abutment (gingiva shaper) is installed on the implant in the case of implantation in a single stage or step disclosure implants normal platform 4,5 ppm, the height -3mm

- Direct the locking titanium abutment.

- Plastic carrier abutment + implant is used to hold the implant + abutment in a sterile package and for twisting and untwisting of the implant into the jawbone at the surgical stage. The lower part of the plastic carrier is used for making a temporary plastic crown.

8. Green Ideal Implant Kit $\varnothing$ 4.2, 5.0, 6.0mm, length 10 , 11.5, 13, 16 mm

Kit for dental implantation , comprising:

- Green Respect Titanium Implant

- Titanium screw to close the internal threads of the implant with a two-stage implantation

Titanium screw - healing abutment (gingiva shaper) is installed on the implant in the case of implantation in a single stage or step disclosure implants normal platform 4,5 ppm, the height -3mm .

- Direct the locking titanium abutment

- Plastic carrier abutment + implant is used to hold the implant + abutment in a sterile package and for twisting and untwisting of the implant into the jawbone at the surgical stage. The lower part of the plastic carrier is used for the manufacture of a plastic temporary crown.

9. Green Vector Titanium Implant Screw , Tapered implant stage variable -tion dvuspiral thread 2x1.4mm (2.8mm) . The intraosseous portion has a double surface treatment : sandblasting , etching , and aluminum oxide kislotoy. Verhnyaya Supraspinous portion in the form of an abutment with a machined surface polishing .

Diameter, mm : 2.4, 3.0 , 3.3, 3.75 , 4.2, 5.0

Length in mm: 10 , 11.5, 13, 16 .

10. Green Vector Plus Titanium Implant Screw , antiratotsionny tapered implant with external hexagon 2.7mm variable double stranded thread 2x1.4mm (2.8mm) . Internal thread implantata 1 - 72 " under the abutment screw . Double surface treatment: sand , aluminum oxide and etching.
Diameter, mm : 3.0, 3.3.
Length in mm: 10 , 11.5, 13, 16 .

11. Green Orthodontic Implant Kit $\varnothing$ 3.0, 3.3mm, length 10 , 11.5, 13, 16 mm

-Green Vector Plus Titanium Implant
and medical titanium abutment

three screw orthodontic titanium abutments :

- first (500140) with the shape of a ball with hexagon socket 1.25 mm hex wrench , which fastens the abutment to the implant
- the second in the form of orthodontic wide breketas hexagon 1.25 mm at the center of the hex key , which fastens the abutment to the implant
- third as orthodontic narrow breketas hexagon 1.25 mm at the center for a hex wrench , which fastens the abutment to the implant.

All implants are made of titanium alloy , G5- grade (TI 6AL 4V)

All implants were sterilized in a single package or bundled .

1. Odnocnaya pack contains:

- implant
- nositetel implant
- cover screw

2 . The kit contains:

- implant
- abutment
- support the abutment
- healing abutment
- cover screw

Titanium implant attachment to Prefix

This top module (PREFIX), titanium parts , which consists of a body in the form of a cylinder of height 3mm Hexagon 2.43mm with a special locking connection or 2.1mm and 2.42mm external hexagon . Green prefix for implant plus made of titanium alloy medical grey5 . The surface has a dual - processing sandblasting and acid obrabotka. Speralevidnaya changing mikrozbu 4x0 , 4mm . Differ Green prefix for implant plus in diameter. Diameters Green prefix for implant plus -3,75 ppm , 4.2 ppm , 5.0 ppm , 6.0 ppm . Is integral part of the upper module Green Classic Plus Titanium Implant and Green Universal Plus Titanium Implant. that its outer hexagon Allen enters the lower module plus implant and forms a fixed connection .

Abutments

Abutment	is attached to the implant part of the structure , provides support for the future crown or prosthesis
Superstructure	Superstructure on the implant , the collective name for all types of abutments and Formers

Titanium abutment Slim Platform: straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Standard titanium abutment : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Screw straight titanium abutment with the shoulder , no shoulder , orthodontic
 Aesthetic titanium abutment : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Aesthetic titanium abutment anatomical : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Wide straight titanium abutment with the shoulder without ledge.
 Plastic Temporary anatomical abutment : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Narrow straight titanium abutment ;
 Zirconia Abutment : straight without ledge , with a ledge , with a slope angle 15 °, 25 °;
 Titanium abutments for milling : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Anatomic titanium abutment : straight without ledge to ledge , angled slopes with 15 °, 25 °;
 Wide angle for milling titanium abutment 15 °, 25 °
 - Locator abutment ;
 - Abutments with spherical attachmentami ;
 - Direct aesthetic abutment screw -connector;
 - Screw the abutment connector;
 - Esthetic Abutment Screw ;
 - Plastic Abutment for casting;
 - Screw titanium abutment for simultaneous loading;
 - Temporary corner titanium abutment for simultaneous downloads 15 °, 25 °;

Method prosthesis depends on the design of the implant .

Inside the main types of implants osteointegriruemyh an internal thread for attachment thereto or prosthetic abutment designs using screws . For most implants , except thread, fittings are to prevent the rotational torque between the implant and abutment (or the prosthetic structure) , which is called " anti-rotation element " or " locking " the implant , in particular an internal or external hexagonal etc. .

Implant is inserted inside the item, which is intermediate between the intraosseous part and prosthetic design. This item is called " abutment . "

Abutments are divided into two groups:

1. Abutment tighten - one side abutment protruding from it - it's screw. Includes fixed abutment in the implant helped with this screw , which is completely immersed in the implant);
2. Abutment lockable (lock) - one side abutment protrudes from it in the form of the castle, which is connected to the implant antiratatsionnym element (abutment protruding from Allen enters the hexagon implant) .. Abutment hollow and is connected by a screw which passes through the entire thickness of the abutment and screwed the internal thread of the implant.

Prosthetics is divided into two large groups.

1. Removable prosthesis :

- Denture on Trimmer attachment (bardolder);
- Denture on spacers (attachments or magnets) .

2. Removable prosthetics :

- Crowns ;
- Bridges ;
- Dentures with screws (this type of prosthesis manufacturing method is similar to conventional denture but has such a broad base and almost does not rely on the oral mucosa) .

Fixed prosthesis divided into two groups by type of fixing : - cemented prostheses (crowns, bridges);

- Prostheses are fixed with screws (crowns, bridges , dentures screws. Prostheses for screw fixing can be called conditionally removable.

constant prosthetics divided into two groups and also by type of abutment connection with a permanent prosthesis :

abutment is final prosthesis

abutment connected to the implant , is the backbone of the final prosthesis , which is sometimes called superstructure ,

There are two types of fixing screws for the final prosthesis :

Screw in the implant screw passes through the thickness of the permanent prosthesis and through the abutment.

Screw in the abutment - the screw passes through the thickness of the permanent prosthesis and is screwed into the internal thread of the abutment , which has already been attached to the implant using existing abutment connector. Fixing permanent prosthesis for this type of "screw in abutment preferable.

Special accuracy required to prosthesis implant , depending on the chosen design and conditions of the oral cavity - can be achieved in different ways:

1 - way - for cement bridges and single crowns, cement-retained or screw ;

Method 2 and 3 using a standard or aesthetic abutment - shows prosthetic with bridges , single crowns and removable prosthesis with screw fixation .

In 2nd and 3rd methods there are several options :

variant using aesthetic abutment , which is used mainly for single crowns ;

- The option of using commonly used standard titanium abutments ;

- Option prosthetic implant denture with fastening on ball attachments ;

- The option of using universal aesthetic abutment , which is also suitable for implants other firms on the market ;

- Option prosthetics based implants with hexagon head .

All of the above implants are intended for use in surgical practice by reducing installation in the bone of the upper and lower jaws and create support for prosthetic devices , such as artificial teeth.

6. Contraindications

A. Full contraindication .

Imply the presence of certain diseases and conditions of the body , when surgery is an obvious risk to human health , as well as when there are untreatable diseases that make it impossible to achieve positive results implantation. These include :

a. Chronic disease decompensation .

b. Systemic disorders of coagulation and hemostasis

in . Immunodeficiencies , HIV seropositivity and any other infection .
of Mental Illness .

B. Relative contraindications or risk factors.

No contraindications , but only diseases that pose some difficulties to achieve the projected results , statistically reduce the effectiveness of implantation and may cause treatment failure . Risk factors should be attributed unfavorable anatomical conditions that require additional surgeries or non-standard approaches to treatment , as well as lifestyle, age of the patient , his intellectual level and emotional status , acute and chronic diseases , as well as pathology, which can be stabilized or homeostasis compensate for changes in the organs and systems of the body thanks to modern methods of treatment. These include :

1. Acute inflammatory diseases and acute viral infections;

2 . Chronic infectious diseases (tuberculosis, actinomycosis , etc.);

3. Exacerbation chronic diseases;

4. High risk of bacteremia (patients with prosthetic heart valves and have suffered bacterial endocarditis , rheumatic fever);

5. Recently transferred heart attack or stroke;

6. Pregnancy and lactation;

Treatment drugs impairing tissue regeneration (hormonal , radiation and chemotherapy , receiving immunosuppressive drugs , etc.);

Teens under the age of 18 years;

Osteopatii - diseases that adversely affect bone formation ;

10. Osteoporoz - overall decrease of bone tissue. Significant shrinkage of the jawbone and unfavorable anatomical conditions currently can not be considered as contraindications to implantation . Modern methods and tools used in surgical dentistry and dental implantology , allow for implantation in almost any anatomic conditions;

11. Osteomalacia - skeletal pathology is at which inadequate bone mineralization of the organic matrix while continuing normal skeletal mass and bone volume;

12. Diseases violating osteogenesis : Diseases of the thyroid gland, parathyroid disease , diabetes, pituitary disease , adrenal disease , blood disease;

13. Alcoholism and addiction not only cause mental changes, but also a variety of somatic disorders affecting bone formation;

14. Systemic connective tissue diseases : systemic lupus erythematosus , dermatomyositis, Segrena syndrome , a group of congenital systemic connective tissue diseases , inherited as an autosomal dominant trait : Knista syndromes , Hurler , Menkes , Gaucher's disease , Niemann -Pick, and various types of congenital dysplasia and dysostosis . They are accompanied by severe osteopathy .

15. Allergic disease. Implant and anufacture of nonbiological , biotolerable material in principle, should not cause an allergic reaction . However , a small number of people , certain metals can cause allergies and so-called metallosis;

16. Use by patiants drugs reduce blood clotting - in time (7 to 14 days) the abolition of medication that promotes liquefaction (decreased clotting) of blood , you can spend the implantation , and the need to work with friendly general patient's physician for additional medical correction somatic condition of the patient .

17. Patologic state maxillofacial and oral cavity : leukoplakia, stomatitis , xerostomia , dental caries , periodontal disease , gingivitis , periodontal disease , macroglossia , malocclusion , diseases of the TMJ , Bruxism and poor oral hygiene . Treatment of these diseases should be administered concurrently with implantation or implantation may be regarded as one of the treatments .

7. Basic precautions.

We can not guarantee the successful application of implants in all clinical situations . And it is recommended to reduce the risk and the surgical instruments used in orthopedic components "Implant System Technologie GmbH & Co. KG", since the combination of the other systems components can not correspond to each other, may cause mechanical damage to the tools , implants , prosthetic components , and also damage to body tissue or unsatisfactory aesthetic results .

8. Precautionary mesures in the preoperative period .

Before proceeding to implantation , it is necessary to draw up a program of inspection and treatment plan by answering the following questions:

- Is this clinical case suitable for implantation?
- There is a need for additional special physical examination and examination ?
- What are your plans for the restoration of the patient's dentition ?
- What bad habits are and what the patient may have an adverse effect these bad habits entering treatment using dental implants :
- effects of nicotine (mucosal disease);
- Muscle dysfunction and parafunksii (e.g., bruxism);
- Failure of oral hygiene and the like ?
- What violations related to the conditions of environmental impacts can affect the course of treatment?

In general, the choice of implant placement and prosthetic design must be carried out individually on the basis of general medical history and examination , and evaluation of mental and emotional state of the patient's dental status .

Precautions during surgical intervention.

In addition to comply with the mandatory requirements for all surgical procedures, dental implant with special requirements must be observed:

- Climatic conditions and the performance of the implants may be at a temperature of 32 ° C to -42 ° C;
- Mucosa incision should be made in the fixed mucosa instead of in the mobile mucosa of alveolar bone. Recommended cut slightly lingually or palatalno; thus when stitching the implant is fully lined with mucous membrane, and the seam is positioned over it. Thus, decreasing the risk of uncovering the implant and its infection. Necessary to ensure good blood flow delaminated areas mucosal;
- Minimal incision or perforation of the alveolar ridge mucosa while ensuring good blood flow implantable bone area, the postoperative period occurs with little collateral molasses;
- bone at preparing the implant, is necessary to preserve the maximum amount of bone in which the implant is installed (around the implant on all sides should be at least 1.5 mm of bone);
- At preparing osteotomy implant need irrigation cold sterile aqueous solution of drills and bone (do not overheat the bone) to save reproductive properties of bone tissue;
- before dissection of bone using a cylindrical or a round hole making carbide burs on top of the alveolar ridge to a depth of 2-3 mm, which is subsequently accurately mark the intended injection site of the implant;
- at preparing osteotomy for screw implant finish drilling the drill diameter 0.5 mm smaller than the diameter of the implant installed;

- when cases where bone density type D3-D4, then finish drilling drill 1.0mm to 3.0mm for less than the diameter of the implant installed,
- implant is packaged in a double sterile packaging. Opened the first package, we take for implant carrier and deduce from the second package. Further insertion of the implant in the bone prepared implant bed, twisting motion clockwise screwed implant, implant carrier or hex wrench 2.5 ppm, up to the top level of the alveolar crest. In this landmark final screwing is the vertical mark on the support implant, hex keys 2.5 ppm or abutment, which should have access to the vestibular side and serves as a reference for any of positioning the right abutment. When administered to an implant requires constant irrigation with cold sterile aqueous solution, to prevent overheating of the implant and the bone tissue;
- it is necessary to achieve the primary fixation, which is very important for successful implantation;
- Install the implant as long as possible and larger diameter;
- Install the implant between two layers of cortical bone which increases the chances of implantation;
- Install the implant angle, identical slope natural teeth;
- installing implants parallel with each other and to the remaining teeth;
- Set the maximum number of implants to achieve a reliable support for the final prosthesis and less pressure on each implant;
- with single defect must install the implant diameter not less than 3.75 mm;
- the volume of bone around the implant must not be less than 1.5 mm on all sides;
- the minimum distance between the implants should be 2-3 mm, to create conditions for the regeneration of bone tissue;
- After installing the implant to the bone level, release the screw connecting the implant with the carrier and remove the implant carrier. If the implant is not set to the level of the bone, use the hex wrench 2.5 mm (431825) for manual administration using key "ratchera" or (430025) for machine administration by fiziodispensera.

Further, the choice of methodology is dependent on the clinical situation, the patient's health status, the number and position of the implants, bone density, presence or absence of keratinized gingiva, hygiene, stability and wishes temporary prosthesis and the patient's mood.

10. STERILITY.

All implants are supplied sterile (sterilization by gamma radiation) and are intended for single use only

13

until expiry date . Do not use implants if the packaging is damaged or opened previously . Sterilization of all other surgical supplies held steam sterilized at a temperature of 135o F. S./274o

Abutments , healing abutments and other superstructures supplied sterile . Before use, they must be sterilized in an autoclave according to the manufacturer - at 121 ° C for 40 min. 30-minute dry cycle. Note that sets with implants , which is the abutment or healing cap , supplied sterile .

Tools are supplied non-sterile and require subsequent sterilization according to surgical protocol.

11. Side effects.

Treatment with implants can lead to inflammation and edemas mucosa around the implant , bone resorption , biological and mechanical damage bone and nerve , fatigue fracture of the implant or portions thereof .

12. Disclaimer.

Specialist using implants is responsible to determine the suitability of the product for a particular patient in a particular situation. Specialist is solely responsible for any direct or indirect damages arising in consequence of any errors in making professional decisions or professional practice , using implants.

13. Limited Warranty.

Company «Green Implant System Technologie GmbH & Co.» Committed to continuous improvement of its products and therefore reserves the right to modify design, product modification or engineering use , price changes , rules and conditions without notice.

Warranty : «Green Implant System Technologie GmbH & Co.» Provides no other warranty, except for guarantee by the absence of defects in the materials from which the products are manufactured , and defects due to production errors . In the event of a defect , pre- notify the company in writing before your return .

«Green Implant System Technologie GmbH & Co.» At its discretion will repair , replacement or refund for defective products.

In case of unsuccessful implantation company agrees to donate the other implant on the following terms :

- Within 6 months from the moment the first signs indicating the failure of implantation , the company sent a written notice of the application unsuccessful implantation follow-up act , drawn up by the company provided form. Corresponding X-ray and rejected by the implant. It is fully guaranteed implantation , provided the company and includes exceptional measures related to the process of implantation.

Buyer is responsible for the use of these products , both separately and in conjunction with other manufacturers' products . «Green Implant System Technologie GmbH & Co.» Strongly recommends additional training to work with implants and strictly adhere to the instructions for use, follow the rules of storage and transportation.

14.Warning.

Acceptable that a licensed dentists and doctors of their order. Sold only on prescription.

15. Necessity to use the notation .

Sterilization by gamma radiation;

Requires a doctor's prescription;

For single use;

Should read the instructions carefully.

Expiry date;

lot number;

date of manufacture;

For more information and indications for use of dental implants should contact the company released sources (website , catalog, user guide , research , etc.).

16. Packaging and Transportation.

Implants "Implant System Technologie GmbH & Co. KG" comes in a double pack to maintain sterility. The packaging consists of two cylindrical sections, wherein one (inner) compartment implant carrier and the covering screw and the second (outer) chamber with an inner compartment and an implant support. Implants for marking cover the outer packaging of various colors and on the outermost packaging is a label with the name of the implant length, diameter, date of sterilization, item, catalog numbers, etc. Color-coded cap outer diameter of the implant packaging shows:

purple markings - diameter2, 4mm
yellow markings - diameter3, 0mm
red markings - diameter3, 3mm;
blue markings - diameter3, 7mm;
green - labeling diameter4, 2mm;
- black marking - diameter5, 0mm;
white markings - diameter6, 0mm

Implant in double sterile packaging, packed in a cardboard box. Also inside the cardboard packaging has a sticker in duplicate showing all the data on the implant. One sticker pasted in history, and the second is labeled on the can extract. On the carton has a sticker with the name of the implant length, diameter, date of sterilization, item, catalog numbers, etc. Also there is an indication of the implant manufacturers, address, website, etc.

Carrying implants in double sterile packaging, packed in cardboard packaging at temperatures above 30 C is below -10 C.

17. Storage.

Storage implants should only be in the original packaging, in a dry dark place at a temperature not exceeding 35 C below 5 C.

Implants in an open or damaged packaging as well as the presence of defects implants non-refundable for free exchange.

Implants can not be re-sterilized.

Prosthetics should be kept in their original packaging. Before their introduction / fitting requires sterilization using steam sterilizer.

Maintenance and repair services - are not subject.

Shelf life: 5 years.

Disposal: Dispose of in accordance with the rules in the medical institutions of the country where the product sold.

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