Medtronic

Medtronic
Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com

ATTESTATION

STATE OF CALIFORNIA
COUNTY OF ORANGE

I, the document custodian, attest that the preceding or attached document(s) is(are) true, exact, complete, and unaltered photocopies made by me of the following document Appendix to IFU for ARES Antibiotic-Impregnated Ventricular and Peritoneal Catheters Kits, issued by Medtronic Inc.

Document Custodian: 
Gilbert Laim
Regulatory Affairs Specialist
Medtronic Inc.

Date: 17-Jun-2021


The document was presented to me and signed by the document custodian, Gilbert Laim

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California


County of Orange

Subscribed and sworn to (or affirmed) before me on this 17th day of June 2021 by Gilbert Laim, proved to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Signature

• (Seal)

ARES Antibiotic-Impregnated Ventricular and Peritoneal Catheters in kits

Appendix to IFU

Name of medical device	ARES Antibiotic-Impregnated Ventricular and Peritoneal Catheters in kits, variants: 1. Ventricular catheter ARES kit impregnated antibiotic, model 91101, contents: - antibiotic impregnated ventricular catheter outer diameter 2,5 mm, length 23 cm - 1 piece; - right angle clip - 1 piece; - stainless steel stylet - 1 piece. 2. ARES Kit of ventricular and peritoneal catheter, model 95001, contents: - antibiotic impregnated ventricular catheter outer diameter 2,5 mm, length 23 cm -1 piece; - antibiotic impregnated peritoneal catheter outer diameter 2,5 mm, length 120 cm - 1 piece; - right angle clip - 1 piece; - stainless steel stylet - 1 piece.																																						
Information on potential consumers of medical device, use conditions	Primary consumers are neurosurgeons and their clinical staff. Secondary consumers are patients with a neurological condition treated by the drainage of CSF from the lateral ventricles.																																						
Technical specifications of MD	All specifications are subject to tolerance of $\pm 5\%$ unless otherwise indicated <table border="1"> <thead> <tr> <th>№.</th><th>Parameter</th><th colspan="2">Values</th></tr> <tr> <th>1</th><th>Catheter</th><th>Ventricular</th><th>Peritoneal</th></tr> </thead> <tbody> <tr> <td>1.1</td><td>Mass/Weight</td><td>1.05 $\pm$ 0.01 g</td><td>5.27 $\pm$ 0.01 g</td></tr> <tr> <td>1.2</td><td>Outer diameter</td><td>0.25 cm</td><td>2.54 mm $\pm$ 0.127 mm</td></tr> <tr> <td>1.3</td><td>Inner diameter</td><td>1.27 mm $\pm$ 0.127 mm</td><td>1.27 mm $\pm$ 0.127 mm</td></tr> <tr> <td>1.4</td><td>Distance from the proximal end of the catheter to mark 3</td><td>3 cm$\pm$0.317 cm</td><td>-</td></tr> <tr> <td>1.5</td><td>Distance from the proximal end of the catheter to mark 5</td><td>5 cm$\pm$0.317 cm</td><td>-</td></tr> <tr> <td>1.6</td><td>Distance from the proximal end of the catheter to mark 10</td><td>10 cm$\pm$0.317 cm</td><td>-</td></tr> <tr> <td>1.7</td><td>Distance from the proximal end of the catheter to mark 15</td><td>15 cm$\pm$0.317 cm</td><td>-</td></tr> </tbody> </table>			№.	Parameter	Values		1	Catheter	Ventricular	Peritoneal	1.1	Mass/Weight	1.05 $\pm$ 0.01 g	5.27 $\pm$ 0.01 g	1.2	Outer diameter	0.25 cm	2.54 mm $\pm$ 0.127 mm	1.3	Inner diameter	1.27 mm $\pm$ 0.127 mm	1.27 mm $\pm$ 0.127 mm	1.4	Distance from the proximal end of the catheter to mark 3	3 cm $\pm$ 0.317 cm	-	1.5	Distance from the proximal end of the catheter to mark 5	5 cm $\pm$ 0.317 cm	-	1.6	Distance from the proximal end of the catheter to mark 10	10 cm $\pm$ 0.317 cm	-	1.7	Distance from the proximal end of the catheter to mark 15	15 cm $\pm$ 0.317 cm	-
№.	Parameter	Values																																					
1	Catheter	Ventricular	Peritoneal																																				
1.1	Mass/Weight	1.05 $\pm$ 0.01 g	5.27 $\pm$ 0.01 g																																				
1.2	Outer diameter	0.25 cm	2.54 mm $\pm$ 0.127 mm																																				
1.3	Inner diameter	1.27 mm $\pm$ 0.127 mm	1.27 mm $\pm$ 0.127 mm																																				
1.4	Distance from the proximal end of the catheter to mark 3	3 cm $\pm$ 0.317 cm	-																																				
1.5	Distance from the proximal end of the catheter to mark 5	5 cm $\pm$ 0.317 cm	-																																				
1.6	Distance from the proximal end of the catheter to mark 10	10 cm $\pm$ 0.317 cm	-																																				
1.7	Distance from the proximal end of the catheter to mark 15	15 cm $\pm$ 0.317 cm	-																																				

Available on request

1.8	Distance from the proximal end of the catheter to the screenprint location	20.0 cm±0.317 cm	29.87 mm
1.9	Catheter length	23 cm ±0,5 cm	120 cm ±1.27 cm
2.0	Radius	0.05 cm	-
2.1	16 through holes of diameter	0.046 cm±0.02 cm	-
2.2	Distance between the nearest holes	0.08 cm	-
2.3	Distance from catheter tip to first hole	0.53 cm	-
2.4	Catheter tip length	0.15 cm+0.05 cm-0	-
2.5	Distance between holes	0,15 cm	-
2.6	Tensile Strength	Minimum 700 PSI	Minimum 700 PSI
2.7	Elongation	Minimum 525%	Minimum 525%
2.8	Dynamic Breaking Strength	10% Elongation	10% Elongation
2.9	Crush Resistance	5,680 lbf	5,680 lbf
2	Stylet		
2.1	Surface roughness	≤ 0.32 μm	
2.2	(Rockwell) hardness	not less 30 HRC	
2.3	Weight	3 ± 0.01 g	
2.4	Total length	259.8 mm ±2.54 mm	
2.5	Outer diameter of the stylet	1.041 mm ±0.025 mm	
2.6	Effective stylet length	241.3 mm ±2.54 mm	
2.7	Circle diameter	17.8 mm	
2.8	Circle radius	8.9 mm	
2.9	Distance	1.397 mm ± 0.381 mm	
3	Right Angle Clip		
3.1	Weight	< 1.0 g	
3.2	Two holes with a radius	0.064 mm	
3.3	Radius on both sides	0.317 mm	
3.4	Diameter	3.86 mm	
3.5	Distance between the centers of the circles	6.096 mm	

Available on request

	3.6	Height	9.86 mm	
	4	Package		
	4.1	Package Type	The ARES catheter's packaging includes a polypropylene thermoform tray with a heat sealed Tyvek lid inside a labeled Tyvek-Poly peel pouch as the sterile barrier system	
	4.2	Package sizes with tolerance (±)		
	4.2.1.	Ventricular catheter:		
		Tray	12.25" (311.2 mm) ×3.94 (100.1 mm) ± 0.1 (L×W) 0.53" (13.5 mm) ± 3.0 mm (Depth)	
		Lid	(308 mm×100 mm) ± 3 mm (L×W)	
		Pouch	(406.4 mm×135 mm) ± 3 mm (L×W)	
		Shelf Carton	(438.28 mm×313.26 mm) ± 3 mm (L×W)	
		Shipping Carton (Box)	490 mm×364 mm×268 mm (L×W×H)	
	4.2.2.	Catheter Kit:		
		Tray	12.25" (311.2 mm) ×3.94 (100.1 mm) ± 0.1 (L×W) 0.53" (13.5 mm) ± 3.0 mm (Depth)	
		Lid	(308 mm×100 mm) ± 3 mm (L×W)	
		Pouch	(406.4 mm×135 mm) ± 3 mm (L×W)	
		Shelf Carton	(438.28 mm×313.26 mm) ± 3 mm (L×W)	
		Shipping Carton (Box)	490 mm×364 mm×268 mm (L×W×H)	
The list and description of medical device materials with direct or indirect contact with the patient's body (human body)	Medical Device name	Material	Body contact	Sterilization Method
	Ares Ventricular Catheter (91101)	Catheter Tubing:	Permanent contact with cerebrospinal fluid, bone, and tissue (Catheter, Drugs & Right Angle Clip) External communicating contact with bone, tissue / Less than 24hrs (Stylet - Used for catheter insertion)	Sterile Steam
		Polydimethylsiloxane		
		Polydimethylsiloxane		
		Barium Sulfate		
		Length Marker and Numerical Printing:		
		Silicone dispersion		
		Graphite		
		Crosslinker		
		Tip Fill:		
		Silicone Elastomer		
		Tantalum		
		Drug:		
		Clindamycin Hydrochloride		
		Rifampin (Rifampicin)		

Available on request

Medtronic

Medtronic Inc.
710 Medtronic Pkwy.
Minneapolis, MN 55432 USA
Tel 763-514-4000
www.medtronic.com

		<i>Stylet:</i>		
		Stainless Steel		
		<i>Right Angle Clip:</i>		
Ares Catheter Kit (95001)		Polypropylene homopolymer		
		<i>Ventricular Catheter, Catheter Tubing:</i>		
		Polydimethylsiloxane		
		Polydimethylsiloxane		
		Barium Sulfate		
		<i>Ventricular Catheter, Length Marker and Numerical Printing:</i>		
		Silicone dispersion		
		Graphite		
		Crosslinker		
		<i>Ventricular Catheter, Tip Fill:</i>		
		Silicone Elastomer		
		Tantalum		
		<i>Peritoneal Catheter, Catheter Tubing:</i>		
		Polydimethylsiloxane		
		Polydimethylsiloxane		
		Barium Sulfate		
		<i>Peritoneal Catheter, Length Marker and Numerical Printing:</i>		
		Silicone dispersion		
		Graphite		
		Crosslinker		
Primary Sterile Packaging for Ares Ventricular Catheter (91101), Ares Catheter Kit (95001)		<i>Drug:</i>		
		Clindamycin Hydrochloride		
		Rifampin (Rifampicin)		
		<i>Stylet:</i>		
		Stainless Steel		
		<i>Right Angle Clip:</i>		
		Polypropylene homopolymer		
		<i>Ares Tray:</i>		
		Polypropylene Homopolymer		
		<i>Ares Lid:</i>		
		Tyvek		
		<i>Ares Lid, Ink:</i>		
Packaging for terminally sterilized Ares Ventricular Catheter (91101), Ares Catheter Kit (95001)		Black Ink		
		<i>Ares Lid, Adhesive:</i>		
		Polymeric adhesive		
		<i>Ares Tray Insert:</i>		
		Tyvek		
		<i>Ares Tray Insert, Ink:</i>		
		Black ink		
		<i>Pouch, Top Web:</i>		
		PET/Polyolefin Film		
		<i>Pouch, Bottom Web:</i>		

Available on request

Medtronic

Medtronic Inc.

710 Medtronic Pkwy.
Minneapolis, MN 55432 USA
Tel 763-514-4000
www.medtronic.com

		Tyvek			
		<i>Pouch, Ink:</i>			
		Black Ink;			
Conditions and procedures for the removal of the medical device	ARES catheter is a permanently implantable medical device. The instructions for use do not provide information on the removal of the implant (explantation). The implant is removed at the discretion of the surgical operator based on the patient's needs.				
Possibility and methods of integration with other medical devices	The ARES catheters are intended to be used as a component of a hydrocephalus shunt system. The ARES ventricular and peritoneal catheters can be connected to any Medtronic Neurosurgery (MNS) shunt valve (Registration Certificate # FSZ 2008/03461).				
Methods and apparatus for disinfection and sterilization cleaning	The product is designed for single patient use only. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.				
Labelling and Packaging of medical device, package contents	Each ARES catheter is supplied separately, in individual sterile packaging, Right angle clip and stainless steel stylet included with ventricular catheter only. Each ARES catheter kit includes one ARES ventricular catheter complete with integrated rectangular clamp and stainless steel stylet and one peritoneal catheter, as well as instructions for use. This product is for single use only. Reuse, reprocessing or re-sterilization of this product is prohibited. Re-sterilization, re-treatment or re-use may compromise the structural integrity of the product and/or risk product contamination resulting in patient injury, illness or death. Do not use if the integrity of the packaging is compromised. Medtronic Neurosurgery is not responsible for the function of the product which has been re-sterilized ARES Ventricular Antibiotic Impregnated Catheter Kit, Model 91101: Each package of ARES Ventricular Antibiotic Impregnated Catheter kit, Model 91101, includes: - antibiotic-impregnated ventricular catheter, 2.5 mm outer diameter, 23 cm length - 1 pc, - right angle clip - 1 pc, - stainless steel stylet - 1 pc. and the instruction for use of the device. Set of Ares ventricular and peritoneal catheters impregnated with antibiotics, model 95001: Each package of the Ares Ventricular and Peritoneal Catheter Kit impregnated with Antibiotics, Model 95001, includes: - antibiotic-impregnated ventricular catheter, outer diameter 2.5 mm, length 23 cm - 1 pc, - Peritoneal catheter impregnated with antibiotics, outer diameter 2.5 mm, length 120 cm - 1 pc - right angle clip - 1 pc,				

Available on request

Medtronic

Medtronic Inc.
 710 Medtronic Pkwy.
 Minneapolis, MN 55432 USA
 Tel 763-514-4000
www.medtronic.com

- stainless steel stylet - 1 pc. and the instruction for use of the device.		
Description of the symbols	Symbol	Description
		European Conformity. This symbol means that the device fully complies with European Directive 93/42/EEC.
		Caution
		Sterilized using steam or dry heat
		Use by
		Package Contents
		Do not reuse
		Do not resterilize

Available on request

Medtronic


Medtronic Inc.

710 Medtronic Pkwy.

Minneapolis, MN 55432 USA

Tel 763-514-4000

www.medtronic.com

		Reference Number
		Lot Number
		Caution: U.S. Federal Law restricts this device for sale by or on the order of a physician
		Store away from direct sunlight
		Upper limit of temperature
		Magnetic resonance safe
		Open here
		Manufacturer
		Manufactured in

Available on request

	Consult Instructions for Use	
BARIUM-IMPREGNATED	Barium-Impregnated	
	Do not use if the package is damaged	
Non toxic	Non toxic MD	

Information in Russian sticker

The final packaging for the Russian market is accompanied with additional sticker.
The following information should be printed on the Russian-language sticker (word order and text font may vary):

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable).

Sample of the Russian sticker for Ventricular catheter ARES impregnated antibiotic kit, model 91101:

	Ventricular catheter ARES impregnated antibiotic kit, model 91101, contents: - antibiotic impregnated ventricular catheter outer diameter 2,5 mm, length 23 cm - 1 piece; - right angle clip - 1 piece; - Stainless steel stylet - 1 piece.	
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer's authorized representative:	
	Registration Certificate № _____ from _____	

Available on request

		Date of manuf./Lot/UBD: See on pack	See IFU See packaging symbols
	Non toxic		
	Sample of the Russian sticker for ARES Kit of ventricular and peritoneal catheter kit, model 95001:		
		ARES Kit of ventricular and peritoneal catheter kit, model 95001, contents: - antibiotic impregnated ventricular catheter outer diameter 2,5 mm, length 23 cm -1 piece; - antibiotic impregnated peritoneal catheter outer diameter 2,5 mm, length 120 cm - 1 piece; - right angle clip - 1 piece; - stainless steel stylet - 1 piece.	
		Manufacturer:	Code UPN
		Country of origin:	Code CFN
		Manufacturer's authorized representative:	
		Registration Certificate № from	
		Date of manuf./Lot/UBD: See on pack	See IFU See packaging symbols
	Non toxic		
Procedure and conditions for the disposal or destruction of a medical device	The hazard class of medical waste is B in accordance with local requirements in Russia. Disposal of the product must be carried out in accordance with legislation and disposed of appropriately.		
Shelf life	Guaranteed shelf life (shelf life): 18 months		
Expected Lifetime	Expected life (service): Expected lifespan of the product, baring complications or adverse events, is expected to be the lifespan of the patient.		
Operating conditions	Catheters	Normal Body Temperature 32°C–42°C External atmospheric pressure and humidity at the implant site must be 100 %.	
	Right Angle Clip	Normal Body Temperature 32°C–42°C External atmospheric pressure and humidity at the implant site must be 100 %.	
	Stylet	Normal body temperature 32°C–42°C (insulated by catheter, used for catheter placement only, not implanted), Humidity 100%	
Conditions of transportation and storage	Storage Conditions:		

Available on request

The list of national and international regulatory documents/standards, which matches medical device	Store this product below 15 – 27 °C (59 - 80 °F), away from direct light. Do not remove the product from the packaging until it will be used. Transportation Conditions: Air temperature: from -29 ° C to 60 ° C Humidity (relative): 30-80% Atmospheric pressure: from 84.0 kPa to 106.7 kPa		
	Number	Date	Title
Quality System			
EN ISO 13485 2016/AC:2016 Medtronic devices- Quality Management Systems- Requirements for regulatory purposes			
Design Verification and Validation			
EN ISO 7197 2009 Neurosurgical implants -- Sterile, single-use hydrocephalus shunts and components			
MEDDEV 2.1.3 Rev 3 Medical Devices Incorporating, as an Integral part, An Ancillary Medicinal Substance: Guidance Document			
General			
ISO 14630 2009 Non-active surgical implants - General requirements			
Risk Management			
EN ISO 14971 2012 Medical Devices — Application of Risk Management to medical devices			
Labeling			
EN 1041 2008 Information Supplied by the Manufacturer of Medical Devices			
EN ISO 15223-1 2016 Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied -- Part 1: General requirements			
Packaging			
ISO 11607-1 2009+A1:2014 Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems			
ISO 11607-2 2006+A1:2014 Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes			
Biocompatibility			
EN ISO 10993-1 2009/AC:2010 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process			
ISO 10993-3 2003 Biological evaluation of medical devices — Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity			

Available on request

	ISO 10993-5	2009	Biological evaluation of medical devices — Part 5: Tests for in vitro cytotoxicity
	ISO 10993-6	2007	Biological evaluation of medical devices — Part 6 Tests for local effects after implantation
	ISO 10993-10	2010	Biological evaluation of medical devices — Part 10: Tests for irritation and skin sensitization
	ISO 10993-11	2006	Biological evaluation of medical devices — Part 11 Tests for systemic toxicity
	Sterilization		
	ISO 17665-1	2006	Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
	ISO 17665-2:2009	2009	Sterilization of health care products - Moist heat - Part 2: Guidance on the application of ISO 17665-1
	EN 285	2006/A2:2009	Sterilization - Steam sterilizers - Large sterilizer
	EN 556-1	2001/AC:2006	Sterilization of medical devices — Requirements for terminally sterilized devices to be labeled “Sterile”
	ISO 11138-3	2006	Sterilization of health care products - Biological indicators - Part 3: Biological indicators for moist heat sterilization processes
	Clinical		
	MEDDEV 2.7.1	Rev 4	Clinical Evaluation: A Guide For Manufacturers and Notified Bodies
	Drug Specific		
	EP Monographs	Ph. Eur 9.0	For: Rifampin CAS 13292-46-1
	EP Monographs	Ph. Eur 9.0	For: Clindamycin CAS 58207-19-5
Customer Service	USP Monograph	USP 40-NF 35	For: Rifampin CAS 13292-46-1
	USP Monograph	USP 40-NF 35	For: Clindamycin CAS 58207-19-5
	Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia:		
	Medtronic, LLC, 123112, Moscow, Presnenskaya nab., 10, 9 floor, room III, com. 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78 www.medtronic.ru		

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Gilbert Laim, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

APES IFU M054364-1001

Description of documents

Is a true and exact copy of the correct and complete original document.

[Signature]

Signature of document custodian (affiant)

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State of California

County of Orange


Subscribed and sworn to (or affirmed) before me on this 14th day of July,
 2020 by Gilbert Laim, proved

to me on the basis of satisfactory evidence to be the person~~s~~ who appeared before me.

[Signature]

Signature

(Seal)

Medtronic

ARES™ ANTIBIOTIC-IMPREGNATED VENTRICULAR AND PERITONEAL CATHETERS

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Rx only

Description

The ARES Catheters are fabricated from radiopaque barium-impregnated silicone elastomer tubing and are kink and compression resistant. The catheters are subjected to a treatment process by which the silicone tubing is impregnated with rifampin and clindamycin hydrochloride.

The ARES silicone catheters have been shown in laboratory studies to reduce the colonization of gram positive bacteria on the tubing surface. Laboratory studies were conducted using staphylococcus aureus and staphylococcus epidermidis.

The quantities of rifampin and clindamycin hydrochloride used to impregnate the catheters are only a small fraction of a therapeutic dose of these two antibiotics, and have no potential for any systemic therapeutic effect.

ARES Antibiotic-Impregnated Ventricular Catheter

The ARES ventricular catheter measures 23 cm in length, 0.13 cm in inner diameter, and 0.25 cm in outer diameter. Lengths are marked in 1 cm intervals starting from 3 cm to 15 cm from the proximal catheter tip, thus enabling the surgeon to gauge the depth of penetration of the catheter into the lateral ventricle. The proximal end of the catheter has 32 flow holes—four lines of eight holes spaced at 90° intervals around the catheter circumference. Components supplied with the ARES Ventricular Catheter include a pre-loaded stainless steel stylet and a Right Angle Clip, which is included to facilitate placement and use of the ventricular catheter (Fig. 1).

ARES Antibiotic-Impregnated Peritoneal Catheter

The ARES peritoneal catheter measures 120 cm in length, 0.13 cm in inner diameter, and 0.25 cm in outer diameter. There are no length markers or wall slits on the catheter, and the tip is open ended. The catheter may be trimmed to the proper length.

Indications

For use in the treatment of hydrocephalus as a component of a shunt system when draining or shunting of cerebrospinal fluid (CSF) is indicated.

Contraindications

Shunting of CSF into the peritoneal cavity, right atrium, or other areas of the body should not be carried out if there is infection in any area of the body. Additionally, shunting into the atrium of patients with congenital heart disease or other serious cardiopulmonary abnormalities, is contraindicated.

Do not implant this device in patients with known hypersensitivity to rifampin or clindamycin hydrochloride.

Use of this device is contraindicated in patients receiving anticoagulants or known to have a bleeding diathesis.

Warnings and Precautions

- Do not immerse the ARES catheter in antibiotic solutions. When immersing the catheter in sterile water or saline, keep the time the catheter is immersed to a minimum. A pale orange color may be imparted to the immersion solution.
- Avoid contacting implantable products with lint, glove tail, oily residue from skin or based soaps, or other surface contaminants. Introduction of contaminants could result in improper performance of components, foreign body reactions, or allergic reactions.
- Improper use of instruments in handling or implanting shunt products may result in the cutting, shifting or crushing of components. Such damage may lead to loss of shunt integrity, and necessitate premature surgical revision of the shunt system.
- Care must be taken to ensure that particulate contaminants are not introduced into shunt components during implantation, testing, or handling. Introduction of contaminants could result in improper performance, over- or under-drainage of the shunt system. Particulate matter which enters the shunt system may result in shunt occlusion, or may hold pressure-flow controlling mechanisms open, resulting in overdrainage.
- The encircling ligatures used to fasten catheters to connectors should securely fasten the catheters to the connectors. However, over-tightening the ligatures could eventually result in cutting through the silicone tubing.
- Care must be taken in the routing of catheters to prevent kinking and needless abrasion along their course. Abrasion can result in premature catheter failure (fracture). The rim of the twist drill or burr hole may be trimmed to provide a beveled notch where the ventricular catheter emerges and is curved to be adjacent to the skull.
- Patients with hydrocephalus shunt systems must be kept under close observation in the postoperative period for signs and symptoms that suggest shunt malfunction. Clinical findings may indicate infection, shunt obstruction, and/or overdrainage of CSF.
- There have been rare reports of patients exhibiting signs of hypersensitivity associated with silicone shunts in the medical literature. Treatment consisted of steroids or replacement with modified-silicone shunt systems. The risks of the underlying condition should be weighed against the benefits of the therapy when treating patients with known or suspected hypersensitivity to silicone.
- Shunt obstruction may occur in any component of the shunt system. The ventricular catheter may become occluded by particulate matter (e.g. blood clots, brain fragments, and bacterial colonization), by investment of the catheter tip in choroid plexus, by embedding of the catheter in brain tissue, or by coaptation of the ventricular walls in the presence of overdrainage ("slit ventricles"). Shunt obstruction may occur due to growth of the child, or physical activities which result in the disconnection of the shunt components or withdrawal of a distal catheter from its intended drainage site.
- Separation of the components of a shunt system due to disconnection or catheter fracture has been reported. In some cases, the separated shunt components have migrated into the cerebral ventricles, the heart, or the peritoneal cavity. Shunt systems may fail due to mechanical malfunction, leading to under- or over-drainage.
- Malfunction or obstruction of the shunt system may lead to signs and symptoms of increased intracranial pressure if the hydrocephalus is not compensated. In the infant, common symptoms may include increased tension of the anterior fontanel, congestion of scalp veins, listlessness, drowsiness and irritability, vomiting and nuchal rigidity. In older children and adults, common symptoms may include headaches, vomiting, blurring of vision, double vision, deterioration of consciousness, and various abnormal neurological findings.
- Overdrainage of CSF may predispose development of a subdural hematoma or hygroma or collapse of the lateral ventricular walls leading to obstruction of the ventricular catheter.
- If the ventricular catheter becomes bound to the choroid plexus or adjacent brain tissue by fibrous adhesions, it is suggested that the catheter should not be removed forcibly. Gentle rotation of the catheter may help to free it and allow for non-forcible removal. It is advised that the catheter be left in place rather than risk intraventricular hemorrhage which may be caused by forcible removal.

Instructions for Use

- A variety of surgical techniques may be used in placing the catheters. Site of placement is at the discretion of the surgeon. The catheters should be checked for patency at the time of surgery.
- A stainless steel stylet, packaged with each ventricular catheter, is designed to facilitate introduction of the catheter into the lateral ventricle. The catheter is packaged with the stylet inserted in the lumen. Once the catheter is placed into the ventricle, the stylet is withdrawn from the catheter.

(Note: Due to the moist heat sterilization process some visible transfer of antibiotic from the catheters to the stylet and packaging materials may occur. This will not impact the safety or performance of the product.)

- The Right Angle Clip may be used to bend the ventricular catheter to an approximate 90° angle bend where it exits the twist drill or burr hole (Fig. 2).

(Note: The clip must be removed when the catheter is used with a Burr Hole Valve.)

The clip may be used as a marker for planned depth of catheter insertion by sliding it the appropriate distance from the proximal end of the catheter prior to inserting the catheter. In doing so, the stylet in the catheter. After the catheter is properly positioned in the ventricle, press the extracranial portion of the catheter into the split tubular segment of the clip to form the right angle bend (Fig. 4). Avoid stretching the catheter when it is pressed into the clip. If the catheter is to be placed directly into the ventricle through a tubular introducer, the clip must be removed prior to insertion of the catheter through the introducer.

- It is recommended that the right angle clip be secured to adjacent tissue by passing sutures through the two suture holes on the sides of the clip (Fig. 2). If the right angle clip is not used, it is recommended that the surgeon trim the rim of the twist drill or burr hole to provide a beveled notch where the catheter emerges and curves to be adjacent to the skull.

CAUTION: LOW TEAR STRENGTH IS A CHARACTERISTIC OF MOST UNREINFORCED SILICONE ELASTOMER MATERIALS. CARE MUST BE TAKEN WITH THE HANDLING AND PLACEMENT OF THE SILICONE ELASTOMER TUBING TO AVOID CUTS, NICKS, OR TEARS.

How Supplied

Each ARES catheter is available separately and is individually packaged sterile and non-pyrogenic. With the ventricular catheter ONLY, an integral Right Angle Clip and stainless steel stylet are included. Each ARES catheter kit includes an ARES ventricular catheter with an integral Right Angle Clip, stainless steel stylet and an ARES peritoneal catheter.

The enclosed product is designed for single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Do not use if package has been previously opened or damaged. Medtronic Neurosurgery is not responsible for the performance of any product that has been re-sterilized.

Patient Education

It is the responsibility of the physician to educate the patient and/or their representatives regarding CSF shunting. This should include a description of the complications associated with implantable medical systems and an explanation of potential alternative products and treatments.

Storage

Store this product below 27°C (80°F), away from direct light. Do not remove the product from the packaging until it will be used.

Complications

- Complications associated with ventriculoatrial and ventriculoperitoneal CSF shunting systems may be similar to those experienced in any surgical procedure carried out under local and/or general anesthesia. These include reactions to drugs and anesthetic agents, electrolyte imbalance, and excessive blood loss, particularly in infants. On rare occasions, the patient may exhibit a reaction due to sensitivity to the implant. Use of sharp instruments while handling these devices can nick or cut the silicone elastomer, resulting in leakage and necessitating valve revision. Care must also be taken when closing incisions to ensure that the valves are not cut or nicked by suture needles.
- In all CSF shunting procedures, the most common complications are due to obstruction of the system as described under Warnings and Precautions. Obstruction may occur in any component of the system due to plugging by brain fragments, blood clots, and tumor cell aggregates. Obstruction may also occur because of the separation of shunt components or kinking and/or coiling of the catheter. This may predispose migration of the ventricular catheter into the lateral ventricle and the distal catheter into the peritoneum, the heart and pulmonary arterial tree, or any other structure in which the catheter is implanted. Growth of the child may cause the distal catheter to be withdrawn from the peritoneum into tissue planes where the fluid cannot be absorbed, or from the atrium into the internal jugular vein. In addition to the types of shunt obstructions described above, there are other potentially serious complications. Local and systemic infections are not uncommon with shunting procedures. Usually, they are due to organisms inhabiting the skin, particularly *Staphylococcus epidermidis*. Other pathogens circulating in the bloodstream may colonize the shunt and, in the majority of patients, require its removal.
- In 1973, Robertson et al. summarized the incidence of infection in ventriculoatrial and ventriculoperitoneal shunts reported up to that time. The incidence of infection in ventriculoatrial shunting varied from 7 to 33%. Infection in ventriculoperitoneal shunting occurred in 5 to 10% of the patients in most of the reports. Because ventriculoatrial shunting predisposes the spread of bacteria into other organs of the body, ventriculoperitoneal shunting is considered less devastating.
- In 1969, Kestle et al. reported significant reductions in infection (less than 4%) with the use of antibiotics, short duration of surgery, surgical experience, and control of the operating room environment, e.g. designated operating room, limited personnel and traffic, covered skin surfaces. The article states that results can also be obtained without the use of antibiotics, but with rigorous perioperative control of the environment.
- The use of prophylactic antibiotics in shunted patients is somewhat controversial because their use may predispose infection by more resistant organisms. The decision to use antibiotics prophylactically rests with the attending physician and/or surgeon.
- Shunting into the peritoneum may fail because of investments of the catheter in loops of bowel or in the greater omentum. Perforation of the bowel by the peritoneal catheter with subsequent development of peritonitis has been described.
- Overdrainage of CSF may result in excessive reduction of CSF pressure and predispose the development of a subdural hematoma or hydroma, and excessive reduction of ventricular size, leading to obstruction because of impingement of the ventricular walls on the inlet holes in the catheter. In the infant, this excessive reduction of pressure will cause marked depression of the anterior fontanel, overriding of cranial bones, and may convert communicating hydrocephalus into obstructive hydrocephalus.
- The incidence of epilepsy after ventricular shunting procedures has been reported. The study also indicated that the incidence of seizures increased with multiple catheter revisions.

Returned Goods Policy





Products must be returned in unopened packages, with manufacturer's seals intact, to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling. Determination of a product defect or mislabeling will be made by Medtronic Neurosurgery which determination will be final. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days.

Warranty

A. Standard Limited Warranty. Medtronic Neurosurgery warrants to the original end user purchaser ("Purchaser") that the enclosed single use implantable product ("Product") purchased by Purchaser, at the time of delivery to Purchaser, shall be substantially free from defects in material and workmanship. Medtronic Neurosurgery makes no warranty (express, implied, or statutory) for Products that are modified (except as expressly contemplated herein) or subjected to unusual physical stress, misuse, improper operation, neglect, improper testing, use in combination with other products or components other than those for which the Products were designed, or use in any manner or medical procedure for which the Products are not indicated.

B. Remedy. Purchaser's exclusive remedy and Medtronic Neurosurgery's sole liability for breach of the foregoing warranty shall be, at Medtronic Neurosurgery's sole option and election, to replace the Product or credit Purchaser for the net amount actually paid for any such Product; provided that (i) Medtronic Neurosurgery is notified in writing within ninety (90) days after Purchaser's receipt of the Product that such Product failed to conform, including a detailed explanation in English of any alleged nonconformity; (ii) such Product is returned to Medtronic Neurosurgery within ninety (90) days after Purchaser's receipt of the Product F.O.B. 125 Cremona Drive, Goleta, California 93117, U.S.A. or as otherwise designated by Medtronic Neurosurgery; and (iii) Medtronic Neurosurgery is reasonably satisfied that the claimed nonconformities actually exist. Except as expressly provided in this paragraph, purchaser shall not have the right to return Products to Medtronic Neurosurgery without Medtronic Neurosurgery's prior written consent.

C. Exclusion of Other Warranties. EXCEPT FOR THE LIMITED WARRANTY PROVIDED IN (A) ABOVE, MEDTRONIC NEUROSURGERY GRANTS NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, AND MANUFACTURER SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC NEUROSURGERY NEITHER ASSUMES NOR AUTHORIZES ANY OTHER PERSON TO ASSUME ANY OTHER LIABILITIES ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF ANY PRODUCT.

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© 2016 Medtronic, Inc.