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CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

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State of California)

County of Orange)


On 20 Jan 2021 before me, Steve Oberlander, Notary Public
Date Here Insert Name and Title of the Officer

personally appeared Michelle Mirzaee
Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Steve Oberlander
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Appendix to IFU

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Michelle Ball 20 JAN 2021

Michelle Mirzaee
International Regulatory Affairs Specialist

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Michelle Mirzaee, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

React IFU

Description of documents

Is a true and exact copy of the correct and complete original document.

Michelle Hall

Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 27th day of March,
2020 by Michelle Mirzaee, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Joshua Michael Crawford
Signature

(Seal)

Medtronic

DWGS71009-001 Rev. C

WORD DRAWING

IFU, React Catheter OUS

General Notes Unless Otherwise Specified

1. Artwork Revision:
 - C (04/2019)
2. Inspect per:
 - QP00001
3. Paper:
 - 30lb Offset Opaque
 - Supplier to call out on cert/packing slip
4. Final Size/ Fold:
8 1/4" x 10 3/4" +/- 1/8"
5. Print:
 - Black
6. Text Translations: (Top to bottom and left to right in order)

English (en)	Portuguese (pt)	Polish (pl)	Croatian (hr)
French (fr)	Finnish (fi)	Turkish (tr)	Lithuanian (lt)
German (de)	Danish (da)	Norwegian (no)	Latvian (lv)
Italian (it)	Greek (el)	Slovak (sk)	Slovenian (sl)
Spanish (es)	Czech (cs)	Romanian (ro)	Brazilian Portuguese (pt-br)
Swedish (sv)	Hungarian (hu)	Bulgarian (bg)	Korean (ko)
Dutch (nl)	Russian (ru)	Estonian (et)	Arabic (ar)
7. Packaging:
 - Carton Pack
 - Supplier to call out bundle size on cert/packing slip
8. Suggested Vendor:
 - OCPC
2485 Da Vinci
Irvine, CA 92614
Phone: 949-475-1900

REV.	DESCRIPTION	APPROVED BY
C	Update to add React 71, Precaution and preparation and use section	See MAP Agile

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React™ Catheter

DEVICE DESCRIPTION

The catheter is a single lumen, flexible, variable stiffness composite catheter. The catheter shaft has a hydrophilic coating that spans the distal 40 cm to reduce friction during use. The catheter shaft is visible under fluoroscopy. The catheter dimensions are included on the individual device label. The catheter is introduced through a guide catheter or sheath and into the intracranial vasculature and guided over a neurovascular guidewire and/or microcatheter. The proximal end of the catheter has a luer fitting to allow attachment of accessories and infusion of liquids through the system. The catheter is provided sterile, non-pyrogenic, and is intended for single use only.

INDICATIONS FOR USE

The React™ 68 Catheter and React™ 71 Catheter are indicated for the introduction of interventional/diagnostic devices into the peripheral and neurovasculature.

The React™ 68 Catheter and React™ 71 Catheter are also indicated for the removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the neurovasculature.

CONTRAINDICATIONS

There are no known contraindications.

COMPATIBILITY

Refer to product label for device dimensions. Refer to labeling provided with other medical technologies to determine compatibility.

WARNINGS

- The catheter should only be used by physicians who have received appropriate training in interventional techniques.
- Do not reuse. The device is intended for single use only. Discard the catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- The catheter has not been tested for use with automated high-pressure contrast injection equipment, do not use this equipment with the device because it may damage the device.

PRECAUTIONS

- Do not use open, kinked, or damaged devices.
- Do not autoclave.
- Use the catheter in conjunction with fluoroscopic visualization.
- Inspect the catheter before use to verify that its size and condition are suitable for the specific procedure.
- Do not advance or withdraw the catheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of appropriate flush solution.
- If flow through the catheter becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Torquing the catheter may cause damage which could result in kinking and possible separation along the catheter shaft. Should the catheter become severely kinked, withdraw the catheter.
- An appropriate anticoagulation therapy should be applied per institutional guidelines.
- Operators should take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute occlusion
- allergic reaction and anaphylaxis from contrast media
- complication at puncture site
- complication of radiation exposure (e.g., alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia)
- embolism
- false aneurysm formation
- infection
- inflammation
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection or perforation

RECOMMENDED PREPARATION AND USE

1. Gently remove the catheter and packaging hoop from the pouch by grasping the packaging hoop and slowly pulling the catheter and hoop out of the pouch.
2. Flush package hoop with saline before removing catheter.
3. Remove the catheter from the packaging hoop.
4. Inspect the catheter for kinks or other damage. If any damage is observed, replace with new device.
5. Connect a rotating hemostasis valve (RHV) to the hub of the catheter and flush the inner lumen with heparinized saline.
6. Introduce the catheter into the vasculature through the sheath, and over a guide wire using a percutaneous entry technique of choice. A split sheath is provided in the packaging to provide support and facilitate the introduction of the tip of the catheter into the sheath.
NOTE: The maximum diameter recommended for a guide wire is 0.038" (0.97 mm). When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).
7. Under fluoroscopic guidance, advance the system over the guide wire until the desired position is attained.
8. Remove the guide wire.
9. If introducing interventional devices, refer to the instructions for use provided by the manufacturer.
10. If performing aspiration, apply aspiration to the catheter.
11. If necessary, additional attempts at flow restoration may be used.
12. Obtain a post-treatment angiogram in accordance with the physicians preferred technique.

STORAGE

Store in a cool, dry and dark place. See the product label for the device "Use-by date". Do not use the device beyond its labeled "Use-by date".

LIMITED WARRANTY

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2019 Medtronic.
Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.