State of Minnesota

SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Warming Blanket, Spinal Underbody

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Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Temperature Management System that includes Warming Blanket, Spinal Underbody (Model 57501) to the Bair Hugger Unit for patient warming and Bair Hugger Unit for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.

A handwritten signature in black ink, appearing to read 'Regina', is written over a horizontal line.

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A handwritten date '4/13/2017' is written in black ink over a horizontal line.

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
A handwritten signature in black ink, appearing to read 'Irma Vargas', is written over a horizontal line.

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INSTRUCTIONS FOR USE WARMING BLANKET, SPINAL UNDERBODY

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket, Spinal Underbody (Model 57501) to the Bair Hugger Unit for patient warming (**hereinafter** the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient. The blanket is a lightweight, single use, and is not made with natural rubber latex.

Indications for Use

Warming Blanket, Spinal Underbody (Model 57501) is intended to prevent unintended hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold.

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- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- The patient should face perforated side of the blanket. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 or 700 series warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.
- ⚠ **WARNING: To reduce the risk of patient injury or death due to altered drug delivery:**
- Do not use a warming blanket over transdermal medication patches.
- ⚠ **WARNING. To reduce the risk of injury due to interference with ventilation:**
- Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.
- ⚠ **WARNING: To reduce the risk of injury due to patient falls:**
- Do not use a warming blanket to transfer or move the patient.
- ⚠ **CAUTION: To reduce the risk of cross-contamination:**
- This warming blanket is not sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.
- ⚠ **CAUTION: To reduce the risk of fire:**
- This product is classified as Class I Normal Flammability as defined by the Consumer product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.
- ⚠ **CAUTION: To reduce the risk of thermal injury:**
- Do not use if primary packaging has been previously opened or is damaged.
- ⚠ **CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:**
- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

Warming Blanket, Spinal Underbody

Warming Blanket, Spinal Underbody is designed to be used as an underbody, convective warming blanket for warming patients during surgical procedures involving the spine while utilizing an open frame spinal surgery table.

The blanket is positioned on the procedure table before the patient arrives and provides full, unrestricted patient access.

This special design of the blanket does not interfere with the adjustment of patient support pads typical of spinal surgery tables and provides clinicians with full patient visualization. The blanket

Model 57501 is rectangular in shape and is sized for adults and intermediate-to-large pediatric patients. Features include patented fluid outlets (drain holes), which minimize pooling of fluids. A removable head section allows visualization of the patient's face throughout the procedure. Perforations on the side of the blanket allow the table frame to pass through the blanket and integrated tie strips secure the blanket to the table frame. A head drape is included, which helps maintain warm air around the head of the intubated patient. The port to connect hose for air supply is located on the head side of the mattress. The device for patients warming Bair Hugger with the hose for air supply is connected with the port.

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Mode of Action

A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to an inflatable blanket positioned over, under or around the patient. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin); warming the patient.

Materials

Component	Material
Blanket	Spunbond Polypropylene
Hose Card	Cardboard (24 Point Solid Bleach Sulfate)
Patient Adhesive	Acrylic 3M 1504XL HiTack Transfer Adhesive TIS

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

Instructions

1. Open and unfold warming blanket to locate the white product label on the blanket. Unfold the blanket so the product label faces the patient.
2. Position the warming blanket onto the surface of the patient footboard or onto the surface of the sling assembly. Using the adhesive tape strip on the underside of the blanket, secure the blanket to the table to prevent it from slipping (Figure A).


3. Pass the warming blanket between the patient footboard (◆) and the patient top thigh set (●) through the frame to the underside of the frame of the spinal surgery table (Figure B). Unfold the blanket under the top thigh pad set (●), top hip pad set (●), patient chest pad (▲), and under the patient head support plate (■).


4. Draw the end of the warming blanket to the end of the spinal surgery table and secure it to the end of the table frame using the tie strips located at the end of the blanket (Figure B).

5. Open the four perforations on both sides of the warming blanket and position the blanket around the edges of the spinal surgery frame (Figure C).


6. Before the patient has been positioned, draw tie strips across the top of the frame of the spinal surgery table and around each side of the patient chest pad (Figure D, Option 1) or directly to the frame (Figure D, Option 2).


7. Pillows and other positioning modalities use for the lower extremities should be positioned under the blanket.

8. Optional: Tear out the circular (face) section at the head of the warming blanket to visualize the patient's eyes and airway from under the table (Figure E).

⚠ Warning: If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.

9. Insert the end of the Bair Hugger warming unit hose into the hose port (Figures F and G). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

⚠ Warning: Do not treat patients with the Bair Hugger hose alone. Always attach the hose to a Bair Hugger blanket before providing therapy.

10. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific warming unit model)

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11. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the warming unit hose from the warming blanket and discard the blanket per Rules for Disposal.

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APOSTILLE

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1. Country: United States of America

This Public Document: Warming Blanket, Lithotomy Underbody

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Temperature Management System that includes Warming Blanket, Lithotomy Underbody (Model 58501) to the Bair Hugger Unit for patient warming and Bair Hugger Unit for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.

A handwritten signature in black ink, appearing to read 'Regina Feferman'.

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
Date

State of Minnesota)
County of Washington) SS


Subscribed and sworn to before me this 13th day of April 2017.

A handwritten signature in black ink, appearing to read 'Irma Vargas'.

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INSTRUCTIONS FOR USE WARMING BLANKET, LITHOTOMY UNDERBODY

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket, Lithotomy Underbody (Model 58501) to the Bair Hugger Unit for patient warming (hereinafter the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. . A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

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Indications for Use

Warming Blanket, Lithotomy Underbody (Model 58501) is intended to prevent unintended hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold. The temperature management system can be used with adult and pediatric patients.

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⚠ CONTRAINDICATION: To reduce the risk of thermal injury:

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⚠ WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.

- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- The patient should face perforated side of the blanket. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 or 700 series warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.
- ⚠ **WARNING: To reduce the risk of patient injury or death due to altered drug delivery:**
 - Do not use a warming blanket over transdermal medication patches.
- ⚠ **WARNING: To reduce the risk of injury due to interference with ventilation:**
 - Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.
- ⚠ **WARNING: To reduce the risk of injury due to patient falls:**
 - Do not use a warming blanket to transfer or move the patient.
- ⚠ **WARNING: To reduce the potential for decreased leg perfusion:**
 - Do not tie the tie strips tightly around the patient's legs as this could act as a tourniquet.
- ⚠ **CAUTION: To reduce the risk of cross-contamination:**
 - This warming blanket is not sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.
- ⚠ **CAUTION: To reduce the risk of fire:**
 - This product is classified as Class I Normal Flammability as defined by the Consumer product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.
- ⚠ **CAUTION: To reduce the risk of thermal injury:**
 - Do not use if primary packaging has been previously opened or is damaged.
- ⚠ **CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:**
 - 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

Warming Blanket, Lithotomy Underbody

Warming Blanket, Lithotomy Underbody is designed to be used as an underbody, convective warming blanket for warming patients during surgical procedures. The lithotomy blanket provides full, unrestricted patient access and flexibility during surgical procedures involving the abdomen, perineum, pelvic cavity and lower extremities. It is rectangular in shape and is sized for adults and intermediate- to-large pediatric patients. Features include a pair of pass-through slits designed for flexible patient positioning and accommodates the use of a draw sheet. Outer edges of the blanket draw up next to the patient when inflated to surround them in warmth. A non-inflated area at the head of the blanket accommodates different height positions, while tape and tuck flaps provide stability by securing the blanket to the OR table. Fluid outlets minimize pooling of fluids on the surface of the blanket. Generous perineal cut-out base of blanket provides patient access and does not interfere with the clinician. A clear head drape is attached, which keeps warm air around the intubated patient's head.

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Materials

Component	Material
Blanket	Spunbond Polypropylene
Hose Card	Cardboard (24 Point Solid Bleach Sulfate)
Patient Adhesive	Acylic 3M 1504XL HiTack Transfer Adhesive TIS

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
Temperature range (Internal requirements) 39°C $\pm 3^\circ\text{C}$


Instructions

Note: In cases where fluid runoff is expected, place an absorbent cloth beneath the model 58501 warming blanket.

1. Lay the non-perforated side of the warming blanket (side with adhesive) on the table. The perforated side (side with small holes) must be toward the patient in direct contact with the patient's skin. Align the lower edge of the warming blanket with the end of the procedure table. Using the adhesive tape strips on the underside of the blanket, secure the blanket to the table to prevent it from slipping (Figure A).


2. Tuck the warming blanket's tuck flaps under the mattress or operating table pad to provide stability (Figure B).

Continue with Step 3 or go to **ALTERNATE PATIENT POSITION.**

3. **For Lithotomy Position Using a Draw Sheet** If using a draw sheet or securement device (i.e. safety strap, tape), tear open the perforations along both sides of the warming blanket. NOTE: Do not open the side perforations if a draw sheet or a securement device are not used. Before you inflate the blanket, lay the draw sheet across the mid-section of the blanket and position the patient on the draw sheet.
4. **OPTIONAL:** Tie strips located at the base of the warming blanket may be opened and loosely tied around the patient's upper thigh. (Figure C).


⚠ Warning: Do not tie tightly as the tie strips around the patient's legs could act as a tourniquet.

5. **OPTIONAL:** The non-inflated center arc-shaped section may be removed by tearing along the perforation (Figure D).


6. Bring the draw sheet over the patient's arms and down through the side perforations on the warming blanket (Figure E). Tuck the ends of the draw sheet securely. (Figure F). If a securement device is being used, pass the securement device through the side perforation on the blanket.

⚠ **Warning:** If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.

7. A hose port is provided at the head of the warming blanket. Insert the end of the Bair Hugger warming unit hose into the hose port (Figures G and H). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

⚠ **Warning:** Do not treat patients with the Bair Hugger hose alone. Always attach the hose to a Bair Hugger blanket before providing therapy.

8. If the patient is intubated and ventilated, position the head drape over the patient's head and neck.

⚠ **Warning:** Do not allow the head drape to cover the patient's head or airway when the

patient is not mechanically ventilated.

9. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific warming unit model)

⚠ Caution: Patient Monitoring Recommendations:


- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
10. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per Rules for Disposal.

ALTERNATIVE OPTIONS FOR PATIENT POSITIONING

Supine Position Using Arm Boards

Before placing the patient on the warming blanket, open the side perforations. Extend the patient's arms through the open perforations and position on the arm board (Figure 1).

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Packaging configuration

Each blanket is sealed in a printed plastic pouch (non-sterile). After being sealed, blankets are packed in either quantities of five or ten blankets per carton depending on commercial offer.

Package size:

Length: 15.37 in / 39.04 cm $\pm$ 1

Width: 7.56 in / 19.2 cm $\pm$ 1

Height: 15.75 in / 40.01 cm $\pm$ 1

Environmental conditions**Operating conditions**

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Blanket is single use, shelf life is unlimited.

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APOSTILLE

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1. Country: United States of America

This Public Document: Warming Blanket, Sterile Cardiac

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Temperature Management System that includes Warming Blanket, Sterile Cardiac (Model 63000) to the Bair Hugger Unit for patient warming and Bair Hugger Unit for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.

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

INSTRUCTIONS FOR USE WARMING BLANKET, STERILE CARDIAC

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket, Sterile Cardiac (Model 63000) to the Bair Hugger Unit for patient warming (**hereinafter** the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. . A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

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Indications for Use

Warming Blanket, Sterile Cardiac (Model 63000) is intended to prevent unintended hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold.

- STERILE.
- Federal law (USA) restricts this device to sale by or on order of a licensed healthcare professional.

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- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave patients unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- The patient should face perforated side of the blanket. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 or 700 series warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.
- Do not place the warming blanket directly over a dispersive electrode pad.

⚠ WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

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⚠ WARNING: To reduce the risk of injury due to interference with ventilation:

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⚠ WARNING: To reduce the potential for injury due to patient falls:

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⚠ WARNING: To maintain proper aseptic technique:

- Any surgical incision in the lower extremities must be closed before applying the blanket to the lower part of the body.

⚠ CAUTION: To reduce the risk of cross-contamination:

- This warming blanket is sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.

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⚠ CAUTION: To reduce the risk of thermal injury:

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⚠ CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:

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Warming Blanket, Sterile Cardiac

Warming Blanket, Sterile Cardiac is designed to be used with warmed or cooled air during cardiac surgeries. The design is comprised of tubular inflatable channels fused to a perforated underside. Approximately 1/3 of the lower portion of the Blanket is non-inflatable and is comprised of a foot drape that is designed to protect the patient's feet during warming. An integral trunk extension is attached to an inlet located in the lower half of the blanket. This trunk extension is attached to the hose of the Bair Hugger warming unit and allows the unit's hose to remain outside the sterile field. A clear, polyethylene window, integral to the blanket, allows access to the groin area. The blanket has a waist belt design, i.e. after surgical field treatment, the blanket in the form of a wide belt is placed in the waist area of the patient, and after surgical wound suturing in the area of the saphenous vein harvesting, the blanket is completely rolled out and inflated after removal of the transverse clip from aorta.

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
Instructions


1. Leave the warming blanket folded up until ready to begin warming therapy.

⚠ Warning: Any surgical incision in the lower extremities must be closed before applying the sterile warming blanket to the lower part of the body.

2. Remove the backing from the adhesive tape strip and adhere the warming blanket to the patient's waist (Figure A).

3. Unfold the warming blanket and position the clear window to allow access to the groin. Extend the warming blanket so that the end is outside the sterile field (Figure B).

4. Attach the sterile hose of the warming blanket to the Bair Hugger warming unit hose (Figure C). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

⚠ Warning: Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.

5. Tuck the edges of the warming blanket under the patient.

Optional: Place a sterile drape over the Bair Hugger warming blanket to increase effectiveness.

6. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific warming unit model)

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- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
7. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per Rules for Disposal.

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All sterile product lots are labeled to include a) whether or not the product is sterile, b) the type of sterilization method used (ETO in this case), and c) the expiration date (date product is intended to be used by).

Packaging configuration

Each blanket is sealed in a Tyvek pouches (sterile). After being sealed in the pouch, blankets are packed in either quantities of five or ten blankets per carton depending on commercial order.

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Sterile blanket is single use, and shelf life constitutes 3 years from sterilization.

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APOSTILLE

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1. Country: United States of America

This Public Document: Warming Blanket, Full Surgery Access

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To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Temperature Management System that includes Warming Blanket, Full Surgery Access (model 63500) to the Bair Hugger Unit for patient warming and Bair Hugger Unit for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.

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

INSTRUCTIONS FOR USE WARMING BLANKET, FULL SURGERY ACCESS

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket, Full Surgery Access (model 63500) to the Bair Hugger Unit for patient warming (hereinafter the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. . A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

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Indications for Use

Warming Blanket, Full Access Underbody is intended to prevent unintended hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold.

- NOT STERILE.

- Federal law (USA) restricts this device to sale by or on order of a licensed healthcare professional.

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WARNING:

Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

CAUTION:

Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury.

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- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave patients unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- The patient should face perforated side of the blanket. Always place the perforated side (with the small holes) directly in contact with the patient's skin.

- In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 or 700 series warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.

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WARNING. To reduce the risk of injury due to interference with ventilation:

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WARNING: To reduce the risk of injury due to patient falls:

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Warming Blanket, Full Access Underbody

Warming Blanket, Full Access Underbody is designed for use on patients during cardiac, abdominal, hepatic, vascular, cosmetic, and trauma surgeries. The blanket is flat in the center (no air channels) and distributes airflow toward the patient through the peripheral channels. Drainage holes for fluid runoff are located in the center of the blanket. Adhesive strips are attached on the non-patient side to secure it to a bed or table to prevent the blanket from slipping. Side tuck flaps can also be used to secure the blanket. Continuous seam seals are located along both sides of the perforations so a draw sheet can be used with the blanket and to allow positioning the patient's arms in the supine, prone, or lateral position. Hose ports are located at both ends of the blanket; one hose port has a plug that can be moved to the other port if the clinician desires to use the port with the plug. The hose from a Bair Hugger warming unit connects to the blanket through one of the ports. The blanket also includes a drape that may be used if the patient is intubated and ventilated.

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Blanket	Spunbond Polypropylene
Hose Card	Cardboard (24 Point Soild Bleach Sulfate)
Patient Adhesive	Acylic 3M 1504XL HiTack Transfer Adhesive TIS

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

Instructions

Note: In cases where fluid runoff is expected, place an absorbent cloth beneath the model 63500 warming blanket.

1. Lay the non-perforated side of the warming blanket (side with adhesive) on the table. The perforated side (side with small holes) must be toward the patient in direct contact with the patient's skin. The head of the blanket should extend 6 inches (15 cm) beyond the table pad. Using the adhesive tape strips on the underside of the warming blanket, secure the blanket to the table to prevent it from slipping (Figure A).

2. Tuck the warming blanket's side flaps under the mattress or operating table pad to provide stability (Figure B).


Continue with Step 3 or go to **ALTERNATE PATIENT POSITIONS.**

3. Supine Position Using a Draw Sheet

If using a draw sheet or securement device (i.e. safety strap, tape), tear open the perforations along both sides of the warming blanket. Lay the draw sheet across the midsection of the warming blanket and position the patient on the draw sheet.

Note: Do not open the side perforations if a draw sheet or securement device are not used.

4. Bring the draw sheet over the patient's arms and down through the side perforations on the warming blanket (Figure C). Tuck the ends of the draw sheet securely under the table pad (Figure D).

5. If a securement device (i.e. safety strap, tape) is being used, pass the securement device through the side perforation on the blanket.


Warning: If a securement device is used, ensure the warming channels are not occluded.

6. Hose ports are provided at either end of the warming blanket for clinician preference. Remove the hose port card from the port (Figure E) you wish to use. Place the hose port card on the port, which is not used during warming. Insert the end of the Bair Hugger warming unit hose into the hose port (Figure F). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.


Warning: Do not treat patients with the Bair

Hugger warming unit hose alone. Always attach the hose to a Bair Hugger blanket before providing warming therapy.

E

F

7. If the patient is intubated and ventilated, position the head drape over the patient's head and neck.

Warning: Do not allow the head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

8. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific Warming Unit Model)

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
9. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per Rules for Disposal.

ALTERNATE PATIENT POSITIONS


Supine Position Using Arm Boards

Before placing the patient on the warming blanket, open the side perforations. Extend the patient's arms through the open perforations and position on the arm board (Figure G). Continue with Step 5.


Lateral Position

Before placing the patient on the warming blanket, open the applicable side perforation. Extend the patient's lower arm through the open perforation. Position the upper arm on top of the blanket's side channel (Figure H). Continue with Step 5.

Prone Position

Before placing the patient on the warming blanket, open the blanket's side perforations and tear out the circular (face) section at the head of the blanket. Extend the patient's arms through the side perforations (Figure I). Continue with Step 5.

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APOSTILLE

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1. Country: United States of America

This Public Document: Warming Blanket, Sterile Cardiac Access

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
To whom it may concern


In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Temperature Management System that includes Warming Blanket Sterile Cardiac Access (model 64500) to the Bair Hugger Unit for patient warming and Bair Hugger Unit for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.

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
INSTRUCTIONS FOR USE WARMING BLANKET, STERILE CARDIAC ACCESS

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket, Sterile Cardiac Access (model 64500) to the Bair Hugger Unit for patient warming (hereinafter the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

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Indications for Use:

Warming Blanket, Cardiac Access is intended to prevent unintentional hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold. The temperature management system can be used with adult and pediatric patients.

- **STERILE**
- Federal law (USA) restricts this device to sale by or on order of a licensed healthcare professional.

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* In the operating room, do not use this warming blanket with any device other than a Bair Hugger 500 or 700 series warming unit.

- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blanket.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds on the warming device Bair Hugger. Unplug the warming unit and contact a qualified service technician.
- If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket directly over a dispersive electrode pad.

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- Do not use if primary packaging has been previously opened or damaged.

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Warming Blanket, Cardiac Access

Warming Blanket, Cardiac Access is designed for use during cardiac surgeries requiring sterile field, or whenever clinically indicated. The design comprises of an inflated tube that is taped around the patient, a footpad that surrounds the feet, which is used to protect the patient from thermal injury, and a hose inlet port at the centermost point of the inflated tube. The hose from a Bair Hugger warming unit connects to the blanket through the hose port at the feet of the patient.

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
Temperature range ((Internal requirements) 39°C ± 3°C

Instructions


1. Position the foot pouch over the patient's feet with the hose port facing away from the patient (Figure A). Ensure that the foot pouch completely covers the feet before starting warming therapy.

2. Remove the backing from the adhesive tape strip (Figure B).


3. Adhere the warming blanket to the outer portion of the patient's legs, starting at the ankles and moving toward the waist (Figure C). Keep the perforated side (the side with small holes) toward the patient in direct contact with the patient's skin.

Warning: If a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.

4. Insert the end of the Bair Hugger warming unit hose in the warming blanket hose port (Figure D). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to insure secure attachment.

Warning: Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger blanket before providing therapy.

5. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator Manual for your specific Warming Unit Model)

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- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
 - Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
- Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the hose from the warming blanket and discard the blanket per Rules for Disposal.

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APOSTILLE

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1. Country: United States of America

This Public Document: Instruction for Use for Bair Hugger Warming Blanket

2. has been signed by Irma Vargas

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Irma Vargas, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: April 18, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9465552-1

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

3M Health Care Business

3M Center
2510 Conway Ave., Bldg. 275-5W-06
St. Paul, MN 55144 U.S.A.
651 733 1110

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M HealthCare Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for Bair Hugger Warming Blanket for patients warming with accessories, Series 500 or 700 in English is truthful and accurate.


Regina Feferman
International Regulatory Affairs Associate
3M Health Care

4/18/2017
Date

State of Minnesota)
County of Washington) SS

Subscribed and sworn to before me this 18th day of April 2017.

Irma Vargas, Notary Public
My Commission Expires: January 31, 2020

INSTRUCTIONS FOR USE


**Pre-Op
Outpatient
Blanket**


**and Full
Warming Blanket**


**Body Chest
Warming
Blanket**


**Access Pediatric
Underbody
Warming Blanket**

**Multi-Access
Warming Blanket**

**Upper Body Warming
Blanket (Adult)**

**XL Upper Body Warming
Blanket**

**Lower
Warming
Small**

**Body Blanket,
Surgical Access Full
Body Warming
Blanket**

Description of the product

Bair Hugger Temperature Management System includes: Warming Blanket to the Bair Hugger Unit for patient warming (**hereinafter** the blanket or product) and Bair Hugger Unit for patients warming with accessories, Series 500 or 700, made by 3M. The system was designed to prevent unintended hypothermia. A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to the mattress/blanket. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

The blanket is a lightweight, single use, and is not made with natural rubber latex.

Model	Description/Identifier
11101	111 Pre-Op and Outpatient Warming Blanket
30000	300 Full Body Warming Blanket
30500	305 Chest Access Warming Blanket
31000	310 Pediatric Underbody Warming Blanket
31500	315 Multi-Access Warming Blanket
62200	622 Upper Body Warming Blanket (Adult)
52301	523 XL Upper Body Warming Blanket
52500	525 Lower Body Warming Blanket, Small
57000	570 Surgical Access Full Body Warming Blanket

Indications for Use

Warming Blankets are intended to prevent unintended hypothermia. In addition, the blankets can be used to provide patient thermal comfort when conditions exist that may cause patients to feel too warm or too cold.

- NOT STERILE.
- Federal law (USA) restricts this device to sale by or on order of a licensed healthcare professional.

Intended use

Therapy is intended to prevent hypothermia. 3M mattresses can be used in all clinical conditions in which the patient can feel too warm or too cold. They can also be used in the operating room to control patient temperature.

Contraindications, Warnings and Cautions

Explanation of Signal Word Consequences

WARNING:

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION:

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

CONTRAINDICATION: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNING: To reduce the risk of thermal injury:

- Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.
- Do not allow the patient to lie on the warming unit hose.
- Do not allow the warming unit hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- The patient should face perforated side of the blanket. Always place the perforated side (with the small holes) directly in contact with the patient's skin.
- In the operating room, do not use these models warming blankets with any device other than a Bair Hugger 500 or 700 series warming unit.
- Do not use a Bair Hugger 200 series warming unit in the operating room.
- Do not use a Bair Hugger 800 series patient adjustable warming unit with any Bair Hugger warming blankets.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.
- Do not place the warming blanket directly over a dispersive electrode pad.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket over transdermal medication patches.

WARNING: To reduce the risk of injury due to interference with ventilation:

- Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

WARNING: To reduce the potential for injury due to patient falls:

- Do not use a warming blanket to transfer or move the patient.

CAUTION: To reduce the risk of cross-contamination:

- This warming blanket is not sterile and is intended for single patient use ONLY. Placing a sheet between the warming blanket and the patient does not prevent contamination of the product.

CAUTION: To reduce the risk of fire:

- This product is classified as Class I Normal Flammability as defined by the Consumer product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

CAUTION: To reduce the risk of thermal injury:

- Do not use if primary packaging has been previously opened or is damaged.

CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

Pre-Op and Outpatient Warming Blankets

The pre-op and outpatient warming blankets are typically used for prewarming and to enhance patient comfort during outpatient procedures or after surgery. This model is not intended for use in the operating room. The blanket is equipped with tubular channels and covers the patient from the shoulders to the feet, providing expanded unventilated part of the blanket in the shoulder area, which makes it possible to tuck the blanket under the patient in the seated or semi-reclined position. Unventilated part of the blanket is also an important feature in the blanket design in the area of the patient's feet. This blanket design prevents the heated air access from the lower part of the blanket in the foot area ensuring safety for the patient. The Bair Hugger connects to the blanket using the port that located directly above the unvented part in the foot area of the blanket.

Full Body Warming Blanket

Full Body Warming Blanket is intended for use outside the operating room such as the Post Anesthesia Care Unit (PACU). The blanket has a tubular, channeled design that covers the patient from the shoulders to the feet that provides even heat distribution. The design includes a non-inflated foot drape that is positioned over the patient's feet to guard against thermal injury in the heat-sensitive foot and lower leg area. The design includes a hose inlet port just above the foot drape and the hose from a Bair Hugger warming unit connects to the blanket through this port.

Chest Access Warming Blanket

Chest Access Warming Blanket is designed specifically for the thoracic surgery patient featuring a clear plastic access panel over the lateral thoracic region that provides clear, easy viewing to the chest area. The blanket has a tubular, channeled design that covers the patient from the shoulders to the feet. The design includes a non-inflated foot drape that is positioned over the patient's feet to guard against thermal injury in the heat-sensitive foot and lower leg area. The design includes a hose inlet port just above the foot drape and the hose from a Bair Hugger warming unit connects to the blanket through this port.

Pediatric Full Body Warming Blanket


Pediatric Full Body Warming Blanket is a small version of the blanket version 30000 (model 300) and is intended for use outside the operating room such as the Post Anesthesia Care Unit (PACU). The blanket has a tubular, channeled design that covers the patient from the shoulders to the feet that provides even heat distribution. The design includes a non-inflated foot drape that is positioned over the patient's feet to guard against thermal injury in the heat-sensitive foot and lower leg area. The design includes a hose inlet port just above the foot drape and the hose from a Bair Hugger warming unit connects to the blanket through this port.

Multi-Access Warming Blanket

Multi-Access Warming Blanket is designed to provide full body coverage while allowing access to multiple portions of the body in post-op. The blanket has a tubular, channeled design that covers the patient from the shoulders to the feet that provides even heat distribution. The design includes a non-inflated foot drape that is positioned over the patient's feet to guard against thermal injury in the heat-sensitive foot and lower leg area. The design includes a hose inlet port just above the foot drape and the hose from a Bair Hugger warming unit connects to the blanket through this port.

Upper Body Warming Blanket

Upper Body Warming Blanket is designed for use during surgeries on the lower half of the body. The blanket has a tubular, channeled design that covers the patient's upper torso, arms, neck, and head of the patient. Two hose inlet ports provide flexibility in positioning the warming unit on either side of the patient. The hose from a Bair Hugger warming unit connects to the blanket through one of these ports. The blanket also has optional tie strips on each arm that can be used to secure the blanket after inflation. A clear head drape is attached, which keeps warm air around the intubated patient's head and provides visualization by the clinician.

XL Upper Body Warming Blanket

XL Upper Body Warming Blanket is designed for use during surgeries on the lower half of the body and is designed for use on large, obese and bariatric patients. The blanket has a tubular, channeled design that covers the patient's upper torso, arms, neck, and head of the patient. The product is similar to the model 62200 but includes one additional air channel across the torso for a larger chest, is 10" (25.4 cm) longer in length and 12" (30.48 cm) wider in width to accommodate a higher rise through the chest and to maximize warming in the chest, shoulders and upper arms. Two hose inlet ports provide flexibility in positioning the warming unit on either side of the patient. The hose from a Bair Hugger warming unit connects to the blanket through one of these ports. The blanket also has optional tie strips on each arm that can be used to secure the blanket after inflation. A clear head drape is attached, which keeps warm air around the intubated patient's head and provides visualization by the clinician.

Lower Body Warming Blanket

Lower Body Warming Blanket is designed for use during surgeries on the upper half of the body and can be used in the supine, lateral or prone position. The blanket has a tubular, channeled design that covers the patient's lower torso, legs and feet. The design includes a non-inflated foot drape that is positioned over the patient's feet to guard against thermal injury in the heat-sensitive foot and lower leg area. Adhesive strips are designed to secure the blanket to the patient. The port for connection with the air hose from the Bair Hugger is located immediately above the unventilated area at the bottom of the blanket.

Surgical Access Warming Blanket

Surgical Access Warming Blanket is a tubular and channeled designed blanket that covers an adult patient from the shoulders to the feet and intended for use primarily during surgeries on the spine when the patient is prone. It can also be used for abdominal surgeries. Part of the design includes an un-inflated foot drape which is positioned over the patient's feet, since lower extremities such as feet can exhibit poor perfusion and are susceptible to thermal injury; this design prevents heat from being delivered to this area as a safety measure. Located approximately in the center of the blanket is a rectangular opening with adhesive tape along the perimeter. The blanket has two hose ports in the shoulder area on opposite sides of the blanket. A head drape is included, which helps maintain warm air around the head of the intubated and ventilated patient.

Accessories

No accessories are associated with this product line.

Mode of Action

A warming unit forms the base of the system by generating warmed air. The warmed air is transferred from the warming unit through a flexible hose. The hose is connected to an inflatable blanket positioned over, under or around the patient. The warm air is then passed through perforations on the underside of the blanket (that is, the side facing the patient's skin), warming the patient.

Materials

Component	Material
Blanket	Spunbond Polypropylene
Hose Card	Cardboard (24 Point Soild Bleach Sulfate)
Patient Adhesive	Acylic 3M 1504XL HiTack Transfer Adhesive TIS

Design Specification

Model	Name/Identifier	Size / Dimension / Weight	Components
11101	111 Pre-Op and Outpatient Blanket	84 x 36 in (213 x 91 cm) ±15% 5 oz (142 g) ±15%	Standard shoes covers: 1 pair Hose Port: 1 Diameter 55 mm ±15%
30000	300 Full Body Warming Blanket	84 x 36 in (213 x 91 cm) ±15% 5 oz (142 g) ±15%	Hose Port: 1шт. Diameter 55 mm ±15%
30500	305 Chest Access Warming Blanket	72 x 36 in (183 x 91 cm) ±15% 5 oz (133 g) ±15%	Drape (2) 24 x 24 in (61 x 61 cm) ±15% Hose Port: 1 Diameter 55 mm ±15%
31000	310 Pediatric Full Body Warming Blanket	60 x 36 in (152 x 91 cm) ±15% 5 oz (142 g) ±15%	Hose Port: 1 Diameter 55 mm ±15%

31500	315 Multi-Access Warming Blanket	84 x 36 in (213 x 91 cm) $\pm 15\%$ 5 oz (142 g) $\pm 15\%$	Hose Port: 1 Diameter 55 mm $\pm 15\%$
62200	622 Upper Body Warming Blanket	78 x 24 in (198 x 61 cm) $\pm 15\%$ 4 oz (113 g) $\pm 15\%$	Head Drape 24 x 24 in (61 x 61 cm) $\pm 15\%$ Hose Port: 2 Diameter 55 mm $\pm 15\%$
52301	523 XL Upper Body Warming Blanket	84 x 36 in (213 x 91 cm) $\pm 15\%$ 5 oz (153.1 g) $\pm 15\%$	Head Drape 24 x 24 in (61 x 61 cm) $\pm 15\%$ Hose Port: 2 Diameter 55 mm $\pm 15\%$
52500	525 Lower Body Warming Blanket	60 x 36 in (152 x 91 cm) $\pm 15\%$ 4.0 oz (113 g) $\pm 15\%$	Hose Port: 1 Diameter 55 mm $\pm 15\%$
57000	570 Surgical Access Warming Blanket	84 x 36 in (213 x 91 cm) $\pm 15\%$ 6 oz (170 g) $\pm 15\%$	Drape 24 x 16 (61 x 41 cm) $\pm 15\%$ Hose Port: 2 Diameter 55 mm $\pm 15\%$

In accordance with IEC 80601-2-35

Temperature restriction


Contact temperature no more than 48°C

Mean contact surface temperature no more than 46°C

Temperature range (Internal requirements) 39°C $\pm$ 3°C


Instructions

1. Place the perforated side of the 3M™ Bair Hugger™ warming blanket (the side with small holes) directly on top of the patient in contact with the patient's skin (Figure A).

2. Where applicable, remove the backing from the adhesive tape strip and adhere the warming blanket to the patient (Figure B). This prevents air from flowing toward the surgical site.

Optional: Place one cloth blanket or sheet on top of the warming blanket to increase effectiveness.


Warning: Do not place patient securement device (i.e. safety strap or tape) over the warming blanket.

Warning: Do not place the warming blanket directly over a dispersive electrode pad.

3. Insert the end of the Bair Hugger warming unit hose in the hose port (Figure C). Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. Support hose to ensure secure attachment.

Warning: Do not treat patients with the warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket before providing warming therapy.

NOTE: See special considerations for Bair Hugger blankets shown below.

4. Select the desired temperature setting on the warming unit to initiate warming therapy. (See the Operator's Manual for your specific warming unit model)

Caution: Patient Monitoring Recommendations:


- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

5. Based on the warming unit model utilized, turn the unit off or to standby mode to discontinue warming therapy. Disconnect the warming unit hose from the warming blanket and discard the blanket per Rules for Disposal

Special Considerations:

Dual Port for Models 62200, 52301, 57000

Two hose ports are provided for clinician preference. Place the removable hose port card in the hose port that is not being used during warming therapy (Figure D).


Head Drape for Models 62200, 52301 and 57000

If the patient is intubated and ventilated, lay the head drape over the patient's head and neck (Figure E); otherwise, tuck the drape between the channels of the warming blanket, away from the patient's head.

Warning: Do not allow the warming blanket or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.


E

Models 62200 and 52301 (Optional)


Pull the tie-strip tabs, centered along the upper and lower edges of the warming blanket. Tie these strips to prevent the inflated blanket from lifting away from the patient. To use the upper body blanket with one

armboard, one half of the warming blanket can be tied off by wrapping tape around it or by tucking it beneath the patient (Figure F).

Blanket, Model 31500 (Optional)


To use the access panels, tear the uncut tab at the edge of the warming blanket. Remove the backing from the tape strip on the center of the blanket. Fold the access panel back and press against the exposed tape. Pull the panel away from the tape to release (Figure G).

Blanket, Model 30500 (Optional)

Chest Access Blanket

Raise the clear plastic drape to provide patient care to the upper body area (Figure H).

Labeling and Packaging

Labeling

Labeling is provided in multiple languages. Per ISO 15223-1 and EN 1041 international and regional symbols for medical devices are used. Instructions for use are provided in the form of an IFU booklet which is contained within each carton.

Packaging configuration

Each blanket is sealed in a printed plastic pouch (non-sterile). After being sealed, blankets are packed in either quantities of five or ten blankets per carton depending on commercial order.

Package size:

Length: 15.37 in / 39.04 cm $\pm$ 1

Width: 7.56 in / 19.2 cm $\pm$ 1
Height: 15.75 in / 40.01 cm $\pm$ 1

Environmental conditions

Operating conditions

Recommended temperature: from 15°C to 25°C

Maximum relative humidity: 80% (up to 31°C) with linear decrease in relative humidity up to 50% at 40°C.

Maximum altitude elevation 2000 m

Transport and storage temperature: from -20°C to 45°C

Medical device is transported by all means of transport.

Disposal

Medical products should be disposed by the local public companies in accordance with SanPiN 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for the Management of Medical Waste", waste class B).

Warranty

3M warrants this product will be free from defects in material and manufacture. Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Blankets are single use, shelf life is unlimited.

Manufacturer

3M Company, 3M Health Care, USA

3M Center, 2510 Conway Avenue, Bldg. 275-5W-06, Saint Paul, Minnesota, 55144, USA

Authorized representative of the manufacturer on the territory of RF

ZAO "3M Russia"

Russia 121614, Moscow, Krylatskaya str. bldg.17/3

+7 (495) 784 7474 / +7 (495) 784 7475 (fax)

3mra@mmm.com

3M, BAIR HUGGER and the BAIR HUGGER logo are trademarks of 3M.

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© 2016, 3M. All rights reserved.