OPERATING MANUAL


for the medical product
Application cap "BougieCap"

Manufacturer: Ovesco Endoscopy AG

Meets standard international and Russian requirements of GOST 50444-92

The product was produced according to Directive 93/42/EEC from 14.06.1993, appendix 2, class 1.

October 23, 2020

Christoph Damm, Dipl.-Ing. (FH)

Head of Quality Management & Regulatory Affairs

Ovesco Endoscopy AG

Friedrich-Miescher-Strasse 9
72076 Tübingen
Germany

ooo | ovesco

Ovesco Endoscopy AG
Friedrich-Miescher-Straße 9
D-72076 Tübingen
www.ovesco.com

Table of contents

1. Medical device design
2. Technical requirements
3. Area of application
4. Operation information
5. Using the device
6. Maintenance information
7. Repair information
8. Labeling
9. Safety information
10. Disposal information
11. Expiration date and period of storage
12. Transportation and storing
13. Warrianties
14. List of international certificates and a declaration of conformity
15. List of the national and Russian standards
16. The main steps of the engineering and manufacturing process

1. MEDICAL DEVICE DESIGN

Application cap BougieCap with accessories (further – application cap) intended for using into the upper gastrointestinal tract during the endoscopic diagnostic and surgical interventions for expanding the target lesion.

Due to ergonomic features of the device, it efficiently helps to increase the visualization during the procedure.


Fig 1. Application cap

Table 1. Specification of application caps

Ref.No.	Endoscope diameter (a)	Outer diameter (b)	Lenght (c)	Mounting depth (d)
400.31.01	5,5 – 6 mm	7,1 mm	14,7 mm	5 mm
400.31.02	5,5 – 6 mm	8,1 mm	16,6 mm	5 mm
400.31.03	5,5 – 6 mm	10,1 mm	20,6 mm	5 mm
400.32.01	9,8 – 10,3 mm	12 mm	26,8 mm	9 mm
400.32.02	9,8 – 10,3 mm	14 mm	28,8 mm	9 mm
400.32.03	9,8 – 10,3 mm	16 mm	33,8 mm	9 mm

Table 2. Product range of the device:

400.31.01	- Application cap with outer diameter 7,1 mm, for the scopes with diemeter of a distal part 5,5 – 6 mm - Instruction for users
400.31.02	- Application cap with outer diameter 8,1 mm, for the scopes with diemeter of a distal part 5,5 – 6 mm. - Instruction for users
400.31.03	- Application cap with outer diameter 10,1 mm, for the scopes with diemeter of a distal part 5,5 – 6 mm. - Instruction for users
400.32.01	- Application cap with outer diameter 12 mm, for the scopes with diemeter of a distal part 9,8 – 10,3 mm. - Instruction for users
400.32.02	- Application cap with outer diameter 14 mm, for the scopes with diemeter of a distal part 9,8 – 10,3 mm.

	- Instruction for users
400.32.03	- Application cap with outer diameter 16 mm, for the scopes with diameter of a distal part 9,8 – 10,3 mm. - Instruction for users

Supply kit:

- 3 application caps with the same diameter
- Instruction for users

2. Technical requirements

Requirements:

1. Total weight of the device is not more than 100 g.
2. The device is a sterile one. The sterilization was produced by ethylene oxide.
3. The device is a single use.
4. The device is Sterile, Non-pyrogenic, Non-toxic.
5. Materials:

Components of the system	Material/Trademark	Standard	Manufacturer
Application cap	- Styrene Acrylonitrile SAN LURAN 388S	- EN 71-9	- Ineos Styrolution, Germany

6. Cracks, pockets, dints, scratches, chipped places, burrs, burns, slag, particles of the rubbing and polishing materials are not acceptable.
7. Fixation force of the application cap to the distal end of the scope should be less than 40H.
8. The force for the dismantling of the application cap from the distal end of the scope should be less than 40H.
9. Overall dimensions of the package that is subject to final sterilization (in mm) for BougieCap: 160±10 x 90±10 mm
10. The material of the sterile package: PET (Polyethylenterephthalat)/PE (Polyethylen) 12/50 um PEEL, high strength paper (Tyvek).
11. The seam strength of the sterile package is defined as $\geq 1,5N / 15mm$.

12. The surface weight of the sterilization materials. These are 70 g/m² (-5/+3) for the sterilization paper and 75 g/m² for Tyvek.

3. AREA OF APPLICATION

The product is used exclusively in stationary conditions in departments of hospitals and clinics. The endoscopist injecting the must be experienced in the administration technique.

4. OPERATION INFORMATION

Do not use if packaging has been damaged.

Do not use after the expiration time.

The condition of the using the device:

- Temperature – 10-40 °C
- Humidity 40 – 100 %

Application cap should meets the requirements of the biological body fluids and secretions, that it's contact during the procedure with a temperature 32 – 42 °C for the device category U6 according to the GOST R 50444.

Application cap is a sterile, disposable device.

The insertion of the application cap should only be performed under endoscopic observation.

The procedure should be stopped if the tissue resistance too great. The procedure should be stopped if the application cap was damaged or lost. Application cap should be removed from the body of the patient by using the endoscopic accessories.

Depends on the potential risk the device classified as 1 (according to the order №4 from 06.06.2012 of the Russian Ministry of Health).

OKPD2 – 32.50.13.110

The code of the device according to Russian medical device classification (according to the order №4 dated 06.06.2012 of the Russian Ministry of Health) – 271550.

INDICATION FOR USE

Application cap BougieCap with accessories (further – application cap) intended for using into the upper gastrointestinal tract during the endoscopic diagnostic and surgical interventions for expanding the target lesion.

CONTRADICTIONS

The application cap must not be used if flexible endoscopic procedures are contraindicated.

The application cap could be used only with flexible endoscope.

The application cap must not be used if the bougienage of a stenosis or stricture is contraindicated.

POSSIBLE COMPLICATIONS

- Damage or perforation of the tissue of the upper gastrointestinal tract
- Bleeding
- Inflammation or infection of the tissue
- Loss of the application cap during the procedure, which could trigger a disturbance in the passage.

5. USING THE DEVICE

5.1 Before use

Please inspect and check the packing carefully before using the device. Products with broken packing are forbidden to be used.

The physician must guarantee that the correct size of application cap is selected for the procedure and for the endoscope.

5.2 Application the device

1. Place the application cap to the distal end of the scope. Slide the endoscope as far as it goes into the cone-shaped tip of the cap.


Fig 3. Mounted application cap on the distal end of the endoscope

2. Fix application cap using the fixing tape that are included in the set. Tension the fixing tape slightly to attach them without any folds. The endoscope must be dry and free from any contamination before the tapes are attached.


Fig 4. Mounted fixing tape

3. Guide the endoscope with the attached application tape until before the target lesion.
4. The punctiform opening on the tip also allows to insert it using a guide wire (max. diameter of the guide wire is 1 mm).
5. Use gentle rotary motions to slide the endoscope with the application cap carefully. If there is too much resistance or the tissue visibly tears, cease the dilation via bougie or, if required, select a narrower diameter.
6. If you need to use an additional working channel, you should mount it on the scope.

6. MAINTENANCE INFORMATION

The product warranty period is 3 years from the day of manufacture.

It is guaranteed that the product has no defects connected with production and materials quality. The product would be replaced if there was no physical/mechanical intervention to the product, inconsistent with the requirements of the operation manual. Make sure that any malperformance appears not due to neglect or noncompliance with instructions. The storage period is 3 years. In case of malfunction under compliance with the requirements of the operation manual, the manufacturer shall replace the product.

The device cannot be re-sterilized. Using of repeated reprocessing device could harm the patient.

For all questions, please contact official distributor in Russian Federation: *Endo Stars LLC*, legal address: 197101, Russia, Saint-Petersburg, Bolshaya Monetnay, 27, lit. A, 12H ; tel./fax: +7-800-555-56-57, e-mail endo@endo-stars.com.

7. REPAIR INFORMATION

Single use. Unrepairable.

8. LABELING

- Marked:
- The description of the device
- REF
- Dimensions
- Diameter
- Word or symbol "non pyrogen" and "non toxic"
- The individual packaging should contain the following symbols:
 - CE mark and identification number of the notified body
 - See instruction for use
- LOT
- Ref number
- Manufacturer
- Date of manufacture
- Use until
- EO sterilized
- Do not reuse
- Do not resterilized
- Do not use if packaging damaged
- Latex-free
- Items per package

Transport labeling should be marked on the pouch, carton boxes or on the package. Labeling is stuck or connected in other reliable ways that guarantee the safety of the goods and labeling.

The signs "Fragile. Caution ", " Protect from moisture "and the inscription:" Storage conditions - 2 ", " Top should be glue to the package.

9. SAFETY INFORMATION

Do not use if packaging has been damaged.

Do not use after the expiration time

The product was produced according to Directive 93/42/EEC from 14.06.1993, appendix 2, class 1

Single use: yes

10. DISPOSAL INFORMATION

According to the Sanitary and Epidemiological Requirements for Disposal of Medical Waste: Sanitary and Epidemiological Regulations and Norms 2.1.7.2790-10. Class B

11. EXPIRATION DATE AND PERIOD OF STORING

Sterility period: 3 years

12. TRANSPORTATION AND STORAGING

The device is packaged in appropriate packaging with instruction for users. Any transportation meet the requirements of temperature mode to the condition 5 (OG4) GOST 15150 is allowed.

Keep the product in a dry place, free from the damaging effects of temperature, humidity, ventilation, sunlight, dust and air containing salt, etc.

After transportation, please pay your special attention to possible damages.

Avoid the places with gas and chemicals.

Take care from humidity.

The device is hold the temperature and humidity according to terms of transportation and storage group 5 (OG4) GOST 15150.

The tools in transport packaging should be sustainable in the process of transportation under GOST R 50444.

Storage period is – 3 years.

13. WARRANTIES

a. The manufacturer grants conformity of the product to this technical documentation subject to the condition of correct using, transporting and storing.

b. All question should be send to an official distributor in Russian Federation:

Endo Stars LLC, 197101, Saint-Petersburg, Bolshaya Monetnaya, 27, Lit. A, 12H,
Phone/fax 8-800-555-56-57, email: endo@endo-stars.com

14. LIST OF INTERNATIONAL CERTIFICATES AND A DECLARATION OF CONFORMITY


- Certificate ISO 13485:2016 № 50970-14-01 от 20.05.2019
- CE Certificate № 50970-16-05 от 20.05.2019
- Certificate of the company's registration dd. 20.01.2016

15. LIST OF THE NATIONAL AND RUSSIAN STANDARDS

The list of the applied RF national standards and international standards for the purpose of safety, effectiveness and quality provision, as well as standards with the control methods: GOST R 50444-92, Directive 93/42/EEC, EN ISO 13485:2016.

Manufacturing of the device is performed in accordance with the Directive 93/42/EEC dd. 08.11.2015, appendix 2 for the class 1 items.

16. THE MAIN STEPS OF THE ENGINEERING AND MANUFACTURING PROCESSES

Notary Daniel Westa

Schaffhausenstrasse 77 • D-72072 Tuebingen

Phone +49 (0)7071 / 13870-0 • Fax +49 (0)7071 13870-99

mail@notar-westa.de

www.notar-westa.de

Authentication of Signature

Herewith I certify that

Mr. **Christoph Damm**, born on May 5th, 1980,

with business address in 72076 Tübingen, Friedrich-Miescher-Straße 9

– identified by German ID-card –

has confirmed the preceding signature as his signature.

city of Tuebingen, October 28th, 2020

federal state of Baden-Wuerttemberg / Germany

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