The State of Texas

Secretary of State

Not for use within the United States of America

This Apostille only certifies the signature, the capacity of the signer and the seal or stamp it bears. It does not certify the content of the document for which it was issued.

Certificate Validation available at www.sos.state.tx.us

APOSTILLE

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- | | |
|--------------------------------|---|
| 1. Country | United States of America |
| This public document | |
| 2. has been signed by | JUDY ANN PEREZ |
| 3. acting in the capacity of | Notary Public, State of Texas |
| 4. and bears the seal/stamp of | JUDY ANN PEREZ,
Notary Public, State of Texas,
Commission Expires: 12-11-23 |

CERTIFIED

- | | |
|---------------------------------------|---------------------|
| 5. at Austin, Texas | 6. on June 29, 2020 |
| 7. by the Secretary of State of Texas | |
| 8. Certificate No. 11944856 | |
| 9. Seal | 10. Signature: |

A handwritten signature in black ink, appearing to read "Ruth R. Hughes", written over a horizontal line.

Ruth R. Hughes
Secretary of State

GF/dp

This certificate is attached to a **26** page document dealing with/entitled **Instruction for Users** and dated **June 17, 2020**.
of pages

ACKNOWLEDGEMENT CERTIFICATE

STATE OF Texas

COUNTY OF Dallas

On this 19th day of June, 2020.

Before me, the undersigned notary public, personally appeared

Keith Wiley

Name of Personal Acknowledging

proved to me through satisfactory evidence of identification, which was

Drivers License

Type of Identification

to be the person whose name is signed on the preceding or attached document and acknowledged to me that he/she signed it voluntarily for its stated purpose.

Judy Ann Perez

Signature of Notary Public


December 11, 2023

Commission Expiration Date of Notary Public

Set of Q-Optics binocular loupe in various options, Quality Aspirators, Inc
IFU for Medical Device

All-Russian product
classification code:
OKPD2 32.50.13.190

Approved by:

Signature/Seal

17.06.2020

**INSTRUCTION FOR USERS
for Medical Device:**

Set of Q-Optics binocular loupe in various options
Manufacturer: *Quality Aspirators, Inc., USA*

Meet with the requirements of international and Russian standards
GOST R 50444-92, GOST-R IEC 60601-1-2010

Contents

- 1. Product Designation**
- 2. Supply Package**
- 3. General Use and Operation Information**
- 4. Indications**
- 5. Contraindications**
- 6. Product Description**
- 7. Conditions and Rules of Product Operation**
- 8. Product Materials**
- 9. Product Specifications**
- 10. Labeling**
- 11. Packaging**
- 12. Transportation**
- 13. Storage**
- 14. Maintenance**
- 15. Minor Repairs and Warranties**
- 16. Disposal and Destruction Information**
- 17. List of National and International Standards**
- 18. List of Data Concerning International Certificates**
- 19. The main steps of the engineering and manufacturing processes**
- 20. Electromagnetic compatibility tables**

1. PRODUCT DESIGNATION

Set of binocular loupes are designed for scaling and high-definition image creation in course of diagnostics and surgeries.

2. SUPPLY PACKAGE

Set of Q-Optics binocular loupe in various options:

- i. Q-Optics TTL Binocular Loupes 2,5x consisiting of:
 1. Q-optics TTL Binocular Loupe 2,5x – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo case -1 pc
- ii. Q-Optics TTL Binocular Loupes 3,0x consisiting of:
 1. Q-optics TTL Binocular Loupe 3,0x – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- iii. Q-Optics TTL Binocular Loupes 3,5x consisiting of:
 1. Q-optics TTL Binocular Loupe 3,5x – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- iv. Q-Optics TTL Binocular Loupes 2,5x High Resolution consisiting of:
 1. Q-optics TTL Binocular Loupe 2,5x High Resolution – 1 pc

2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- v. Q-Optics TTL Binocular Loupes 3,5x Prism consisting of:
1. Q-optics TTL Binocular Loupe 3,5x Prism – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- vi. Q-Optics TTL Binocular Loupes 4,0x Prism consisting of:
1. Q-optics TTL Binocular Loupe 4,0x Prism – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- vii. Q-Optics TTL Binocular Loupes 4,5x Prism consisting of:
1. Q-optics TTL Binocular Loupe 4,5x Prism – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- viii. Q-Optics Flip-Up Binocular Loupes 2,5x consisting of:
1. Q-Optics Flip-Up Binocular Loupe 2,5x – 1 pc
 2. LED Headlight – 1 pc

3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- ix. Q-Optics Flip-Up Binocular Loupes 3,0x consisiting of:
1. Q-Optics Flip-Up Binocular Loupe 3,0x – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- x. Q-Optics Flip-Up Binocular Loupes 3,5x consisiting of:
1. Q-Optics Flip-Up Binocular Loupe 3,5x – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- xi. Q-Optics Flip-Up Binocular Loupes 2,5x High Resolution consisiting of:
1. Q-Optics Flip-Up Binocular Loupe 2,5x High Resolution – 1 pc
 2. LED Headlight – 1 pc
 3. Q-optics Battery Pack – 1 pc
 4. Battery Charger – 1 pc
 5. Anti-Curing Filter – 1 pc
 6. Shirt Clip – 1 pc
 7. Spiral Wrap -1 pc
 8. Sport Strap -1 pc
 9. Combo Case -1 pc
- (hereinafter – Binocular Loupes)

3. GENERAL USE AND OPERATION INFORMATION

Product use conditions: general, plastic and dental surgery units, ENT units of healthcare facilities.

Set of the binocular loupes represent a multiple-use product.

Set of the binocular loupes represent a non-sterile medical product.

In terms of potential application risk, the medical light grades as class 1 (according to the Decree of the Ministry of Health of the Russian Federation No. 4н dated 06.06.2012).

All-Russian product classification code by types of economic activities– OKPD2
32.50.13.190

Connection of any products not included in the binocular loupe supply package is strictly forbidden.

Set of the binocular loupes are used in a patient environment.

Patient touching of any parts of the binocular loupes under maintenance is strictly forbidden.

4. INDICATIONS

Set of the binocular loupes are used for scaling and high-definition image creation in course of diagnostics and surgeries.

5. CONTRAINDICATIONS

Set of the binocular loupes are customized for particular treating specialists; never use a product not customized for you.

Do not use the binocular loupes to light the eye area.

6. PRODUCT DESCRIPTION

Set of Q-Optics binocular loupe in various options:

1. Q-Optics Binocular Loupes

Figure 1

Q-Optics Binocular Loupes are an integral part of the following supply package:

Q-Optics TTL Binocular Loupes 2,5x, Q-Optics TTL Binocular Loupes 3,0x, Q-Optics TTL Binocular Loupes 3,5x, Q-Optics TTL Binocular Loupes 2,5x High Resolution, Q-Optics TTL Binocular Loupes 3,5x Prism, Q-Optics TTL Binocular Loupes 4,0x Prism, Q-Optics TTL Binocular Loupes 4,5x Prism, Q-Optics Flip-Up Binocular Loupes 2,5x, Q-Optics Flip-Up Binocular Loupes 3,0x, Q-Optics Flip-Up Binocular Loupes 3,5x, Q-Optics Flip-Up Binocular Loupes 2,5x High Resolution.

Q-Optics Binocular Loupes represent a product designed to create required work conditions for doctors and/or medical personnel.

Binocular Loupes are an optical device and are made with a precision accuracy and with account of individual parameters of a treating specialist. Working distance is set with a high precision, wherefore the binocular loupe ensures normal working position of a treating specialist with regard to a wide range of body position changes without the necessity to adjust focus distance to individual working position.

Production range for the binocular loupes varies depending on individual preferences of treating specialists:

1. Q-Optics Flip-Up Binocular Loupes (zooming: 2,5x, 3,0x, 3,5x; working distance: 25-56cm, field of view: 4-10 cm, weight $46-48 \pm 3$ g)


Figure 2 (Assembly of Lenses)

2. Q-Optics TTL Binocular Loupes (zooming: 2,5x, 3,0x, 3,5x; working distance: 25-66cm, field of view: 3-10cm, weight $36-38 \pm 3$ g) - a lower-weight model.


Figure 3 (Assembly of Lenses)

3. Q-Optics TTL Binocular Loupes, High Resolution (zooming: 2,5x; working distance: 25-66cm, field of view: 7-10cm, weight 43 ± 3 g) – a model with an ultra-dense lens system. Eyepieces of the loupes represent a combination of various types of lenses, which heightens image definition.


Figure 4 (Assembly of Lenses)

4. Q-Optics Flip-Up Binocular Loupes, High Resolution (zooming: 2,5x; working distance: 25-56cm, field of view: 7-10cm, weight $55 \pm 3g$) - a model with an ultra-dense lens system. Eyepieces of the loupes represent a combination of various types of lenses, which heightens image definition.

Figure 5 (Assembly of Lenses)

5. Q-Optics TTL Binocular Loupes Prism (zooming: 3,5x, 4,0x, 4,5x; working distance: 25-68cm, field of view: 5-10cm, weight $54 - 58 \pm 5g$) – an ergonomic prismatic loupe with a robust rigid construction.


Figure 6 (Assembly of Lenses)

2. LED headlight, Q-optics battery pack, battery charger, anti-curing filter, shirt clip, spiral wrap, sport strap, combo/loupe case.

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The LED light with its components represent a product used to create required work conditions for doctors and/or medical personnel. Technical properties of this product include creation of necessary light flow through the built-in LED light sources for confident execution of diagnostics and surgeries using the binocular loupes.

7. CONDITIONS AND RULES OF PRODUCT OPERATION

Prior to the very first work with the binocular loupes, the product must be adjusted for comfortable operation.

- In order to adjust the nose pad of the binocular loupes so that the rim would be situated higher along the bridge, carefully bend the lateral curves toward each other (see Figure 15).

Figure 15 (Nose pad)

- In order to adjust the nose pad of the binocular loupes so that the rim would be situated lower along the bridge, carefully bend the lateral curves from each other (see Figure 16).


Figure 16 (Nose pad)

- For safer and more comfortable product operation, we recommend to bend the temple pads (see Figure 17).

Figure 17 (Temple pad of the binocular loupes)

The binocular loupes ensure zooming of the field under examination of a treating specialist, but do not correct vision. For vision correction, use dioptric inserts made in accordance with your eye doctor's prescription.

Using dioptic inserts:

- Insertion of the dioptic inserts:

Put the slot located on the left side of the dioptic inset's opposite the mounting screws (two metal pins) located on the internal part of the rim as Figure 18 shows. The mounting screw cap must fit the slot on the dioptic insert's clamp.


Figure 18 (dioptic insert's clamp mounting)

Move the clamp from the mounting screws as Figure 19 shows. Move from right to left on the internal part of the rim.


Figure 19 (dioptic insert's clamp mounting)

To remove a dioptic insert, move the clamp of the mounting screws as Figure 20 shows. Move from left to right on the internal part of the rim.


Figure 20 (Removal of a dioptic insert)

Having inserted dioptic inserts (if necessary), the specialist must remove the nameplate with the manufacturer's name of the moving from top to bottom as Figure 21 shows. The nameplate must be removed in order to release the mounting screws to mount the LED light.


Figure 21 (Nameplate removal)

Upon removal of the nameplate, it is necessary to mount the LED light. In order to do that, match the slot of the LED light clamp with the mounting screw, move the light's clamp toward the mounting screws as Figure 22 shows.


Figure 22 (LED light mounting)

Note: Working with the binocular loupes, always use the LED light for improved illumination and visualization of the examined area.

Thereafter, put on the binocular loupes with the mounted light, fix the sport strap using the flexible spiral wrap and the shirt clip as it is comfortable for you.

Upon completing examination, put on the binocular loupes and the light. In order to charge the battery, connect charger to the Q-optics battery. A red indicator light will come on on the battery (charging process); when the indicator light turns green the charging is complete.

8. PRODUCT MATERIALS

Set of Q-Optics binocular loupe in various options

No.	Syste, Component Name	Color Marks and Materials	Material/ Trademark	Name	Human Contact Type
1	LED Headlight	N/A	Working part: 1. Optical flint Framework: 2.Polycarbonate Case: 3.ABS resin	1. BK-7/SF-10 2.SC1004A-EM59P 3.VH-0810 G 7643	Brief, with uninjured skin
2	Q-optics Battery Pack	N/A	Case: 1.ABS resin 2. Stainless steel	1. VH-0810 G 7643 2. STS301	Brief, with uninjured skin
3	Anti-Curing Filter	N/A	Optical flint ABS resin	1. VH-0810 G 7643 2. VH-0810 G 7643	Brief, with uninjured skin
4	Shirt Clip	N/A	1. SAN LURAN 388S 2. Stainless steel	1. EN 71-9 2. STS301	Brief, with uninjured skin
5	Spiral Wrap	N/A	SAN LURAN 388S	EN 71-9	Brief, with uninjured skin
6	Battery Charger	N/A	Case: 1.ABS resin 2. Stainless steel	1. VH-0810 G 7643 2. STS301	Brief, with uninjured skin
7	Sport Strap	N/A	PE Dyneema Purity UG dtex 110 TS100	ISO 10993 – 3,4,5,6,10,11,12	Brief, with uninjured skin
8	Binocular loupes	N/A	1.Optical flint Framework: 2.Polycarbonate 3. Titanium	1. BK-7/SF-10 2.SC1004A-EM59P 3.Grades 1,3,4,7,11	Brief, with uninjured skin
9	Combo case	N/A	ABS resin	VH-0810 G 7643	Brief, with uninjured skin

9. PRODUCT SPECIFICATIONS

9.1 Overall dimensions of the products:

9.1.1 Binocular Loupes:

- Length – 150 ± 30 mm
- Width – 145 ± 10 mm
- Height – 60 ± 15 mm

9.1.2 LED Headlight

- Length – 25 ± 5 mm
- Width – 20 ± 3 mm
- Height – 20 ± 5 mm

9.1.3 Q-optics Battery Pack

- Length – 77 ± 5 mm
- Width – 54 ± 5 mm
- Height – 35 ± 5 mm

9.1.4 Battery Charger

- Length – 70 ± 5 mm
- Width – 30 ± 5 mm
- Height – 44 ± 5 mm

9.1.5 Anti-Curing Filter

- Length – 28 ± 5 mm
- Width – 22 ± 5 mm
- Height – 9 ± 2 mm

9.1.6 Shirt Clip

- Length – 50 ± 5 mm
- Width – 13 ± 5 mm
- Height – 15 ± 5 mm

9.1.7 Spiral Wrap

- Length – 36 ± 5 mm
- Width – 4 ± 1 mm
- Height – 4 ± 1 mm

9.1.8 Sport Strap

- Length – 29 ± 3 mm
- Width – 18 ± 5 mm
- Height – 10 ± 2 mm

9.1.9 Combo Case

- Length – 270 ± 20 mm
- Width – 230 ± 20 mm
- Height – 85 ± 15 mm

9.2 The product weight comprises:

9.2.1 Binocular Loupes:

- Q-Optics TTL Binocular Loupes – not over 40g
- Q-Optics Flip-Up Binocular Loupes – not over 50g
- Q-Optics TTL Binocular Loupes High Resolution – not over 50g
- Q-Optics Flip-Up Binocular Loupes High Resolution – not over 60g
- Q-Optics TTL Binocular Loupes Prism – not over 65g

9.2.2 LED Headlight – not over 25 g

9.2.3 Q-optics Battery Pack – not over 150g

9.2.4 Battery Charger – not over 130g

9.2.5 Anti-Curing Filter – not over 5 g

9.2.6 Shirt Clip– not over 10g

9.2.7 Spiral Wrap – not over 10g

9.2.8 Sport Strap – not over 15g

9.2.9 Combo Case – not over 600g

9.3 Zooming of the binocular loupes:

- Q-Optics TTL Binocular Loupes – 2,5x, 3,0x, 3,5x
- Q-Optics Flip-Up Binocular Loupes – 2,5x, 3,0x, 3,5x
- Q-Optics TTL Binocular Loupes High Resolution – 2,5x
- Q-Optics Flip-Up Binocular Loupes High Resolution – 2,5x
- Q-Optics TTL Binocular Loupes Prism – 3,5x, 4,0x, 4,5x

9.4 Working distance of the binocular loupes:

- Q-Optics TTL Binocular Loupes – 250 – 660mm
- Q-Optics Flip-Up Binocular Loupes – 250 – 560mm
- Q-Optics TTL Binocular Loupes High Resolution – 250 – 660mm
- Q-Optics Flip-Up Binocular Loupes High Resolution – 250 – 560mm
- Q-Optics TTL Binocular Loupes Prism – 250 – 680mm

9.5 Field of view of the binocular loupes:

- Q-Optics TTL Binocular Loupes – 30 – 100mm
- Q-Optics Flip-Up Binocular Loupes – 40 – 100mm
- Q-Optics TTL Binocular Loupes High Resolution – 70 – 100mm

- Q-Optics Flip-Up Binocular Loupes High Resolution – 70 – 100mm
- Q-Optics TTL Binocular Loupes Prism – 50 – 100mm

9.6 The LED light must provide illumination of at least 12000lx at a working distance of 250mm

9.7 Parameters of the Q-optics Battery Pack:

- Number of lithium-ion batteries – 1 piece
- Voltage of lithium-ion batteries – 4.2V
- Capacity of lithium-ion batteries – 4500mAh
- Battery charging time – at least 4h
- Working time of the product with a battery – at least 8h

9.8 Parameters of the charger for the battery:

- Charger input voltage – 100-240V
- Charger output voltage – 9V
- Charger frequency – 50/60Hz

9.9 The LED light must provide a light spot of 75mm at a working distance of 460mm

9.10 Cracks, blisters, nicks, scratches, chips, заусенцы, burrs, drosses, grinding and polishing material particles are not admissible.

9.11 Products must withstand temperature and air humidity effects in the process of transportation and storage under group-5 conditions (OZH4 – particularly severe) according to GOST 15150.

9.12 during operation, products must be resistant to environment climatic aspects effect in accordance with the requirements to products designed for UKHL 4.2 (moderately cool climate) according to GOST 15150.

9.13 Products in a transportation package must be resistant to mechanical during transportation and comply with GOST-R 50444.

10. LABELING

10.1. The product labeling must comply with indications in the user manual and the requirements of GOST-R 50444-92.

10.2. The packaging must contain the following basic product information:

- name of the manufacturer,
- manufacturer's address and phone number,
- package contents description,

Additional product information inclusion is admissible.

Reference number


Lot number

Warning

Manufacturer

Serial number

11. PACKAGING

Individual packaging consists of a carton box with inserted dividers and foam fillers.

12. TRANSPORTATION

The product is packed into an appropriate packaging together with operational documentation. The product transportation conditions must correspond to storage conditions 5 according to GOST ГОСТ 15150 (at an ambient temperature of -50-50°C and relative humidity of up to 100%)

13. STORAGE

Conditions for storage of the product in a transportation package at the manufacturer's (user's) warehouses must correspond to storage conditions 2 according to GOST 15150 (at an ambient temperature of -40-50°C and relative humidity of up to 98%)

14. MAINTENANCE

Cleaning:

Binocular Loupes:

- The binocular loupes must be cleaned with a watered or alcohol-wetted soft cloth. Do not use any solvents or detergents.
- Sterilization of the loupes with any other agents is forbidden.

LED light:

- Before cleaning, make sure the device is off.
- The outer surface of the product may be cleaned with a damp cloth, 3-% penozone with 0.5-% solution of a detergent like Lothos or 1-% chloramine solution.
- To disinfect the outer surface, use a soft sterilizing agent in accordance with the manufacturer's instructions. Clean the light's lens with gentle circular motions. Never clean the light's lens moving 'back-and-forth' – this may result in scratches on the light's lens.
- Application of abrasives or rough detergents is also forbidden.
- Prevent the medical light from exposure to liquids.
- Sterilization of the light with any other agents is forbidden.

15. MINOR REPAIRS AND WARRANTIES

The user does not carry out minor repairs.

Warranty for materials and quality of the medical light is provided for one (1) year from the initial trade date. If a product breaks down due to defects in the material or improper quality within the one (1)-year warranty period, the product is repaired or replaced.

In case of warranty event, contact the manufacturer's agent:

Dentrum Ltd., Varshavskoye shosse, 152, 8, 199, Moscow 117405, phone: 8-495-740-56-20

16. DISPOSAL AND DESTRUCTION INFORMATION

Medical devices, packing material and accessories must be disposed in accordance with applicable national laws and standards. In accordance with the Sanitary Regulations and Norms (SanPiN) 2.1.7.2790, the cart grades as A-class – epidemiologic non-hazardous wastes.

17. LIST OF NATIONAL AND RUSSIAN STANDARDS


The list of applicable RF national standards and international standards for safety, efficiency and quality assurance, as well as standards with test methods is as follows: GOST-R 50444-92, GOST-R IEC 60601-1-2010, EN ISO 13485, FDA certificate No. 1781-11-2018.

The product is manufactured in accordance with FDA and ISO 13485 requirements and the list of devices No. 1781-11-2018 for class-1 products.

18. LIST OF DATA CONCERNING INTERNATIONAL CERTIFICATES

- Certificate EN ISO 13485 dated 27.11.2018
- Certificate FDA for foreign government bodies
- Certificate EC dated 27.11.2018
- Company registration certificate

19. THE MAIN STEPS OF THE ENGINEERING AND MANUFACTURING PROCESSES

20. ELECTROMAGNETIC COMPATIBILITY TABLES

20.1 Electromagnetic Compatibility Tables

Table 1 – MANUFACTURER'S Guidelines and Statement concerning ELECTROMAGNETIC RADIATION for MEDICAL EQUIPMENT and MEDICAL SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Radiation		
Radiation Value Check	Compatibility	Electromagnetic Environment – Guidelines
RF radiation CISPR11	Group 1	RF energy is used only for internal operation. Therefore, RF radiation is very insignificant and will probably not cause interferences for electronic installations nearby. Suitable for use in all institutions other than household facilities – thus, the device is connected directly to a low-voltage public power network providing power supply to buildings used for household purposes.
RF radiation CISPR11	Class A	
Harmonic IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Compatible	

20.2 Manufacturer's Guidelines and Statement Concerning Electromagnetic Immunity for Medical Equipment and Medical Systems

Table 2 – MANUFACTURER'S Guidelines and Statement concerning ELECTROMAGNETIC IMMUNITY for MEDICAL EQUIPMENT and MEDICAL SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Immunity			
Immunity Check	EN/IEC 60601 Test Level	Immunity Requirements Compliance Level	Electromagnetic Environment – Guidelines
Static electricity EN/IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	The floor must be wooden, concrete or tiled. For synthetic floors, relative humidity must be at least 30%.
Microsecond pulse disturbances EN/IEC 61000-4-4	±2kV network ±1kV inputs/outputs	±2kV network ±1kV inputs/outputs	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment.
Surge voltage/Surge current EN/IEC 61000-4-5	±1kV differential ±2kV general	±1kV differential ±2kV general	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment.
Voltage dips / Failures EN/IEC 61000-4-11	Dip>95% for 0.5 cycle Dip60% for 5 cycles Dip30% for 25 cycles Dip>95% for 5 seconds	Dip>100% for 0.5 cycle Dip60% for 5 cycles Dip30% for 25 cycles Note 1	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment. Whenever continuous operation is required during power interruptions, it is recommended to connect to an uninterruptible power supply or an accumulator battery.
Industrial frequency 50/60Hz Magnetic field EN/IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields must correspond to standard power frequency magnetic fields in an industrial and household or hospital environment.

Note 1: The equipment under test turns off, and the operator must manually restart it in 5 seconds after the AC power outage.

20.3 Manufacturer's Guidelines and Statement Concerning Electromagnetic Immunity for Medical Equipment and Medical Systems That Do Not Represent Life-Support Equipment or Systems

Table 3 – MANUFACTURER'S Guidelines and Statement concerning ELECTROMAGNETIC IMMUNITY for MEDICAL EQUIPMENT and MEDICAL SYSTEMS that do not represent LIFE-SUPPORT EQUIPMENT OR SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Immunity

Immunity Check	EN/IEC 60601 Test Level	Immunity Requirements Compliance Level	Electromagnetic Environment – Guidelines
Induces radio waves/interferences EN/IEC 61000-4-6	3V (r.m.s. voltage) from 150kHz to 80MHz	(3)V (r.m.s. voltage)	The below calculated/specified minimal distance must be observed between the portable and movable communications equipment: $D=(3,5/\sqrt{1})(P \text{ sq. root})$ $D=(3,5/E1)(P \text{ sq. root})$ from 80 to 800 MHz $D=(7/E1)(P \text{ sq. root})$ from 800 MHz to 2,5 GHz
Radiated radio waves EN/IEC 61000-4-3	3V (r.m.s. voltage) from 80MHz to 2.5GHz	(3)V/m	Where P – max. power in Watts, and D – recommended distance in meters. The field intensity created by transmitters and determined via site EM survey, must not exceed the immunity requirements compliance levels (V1 and E1). Interferences may arise near devices equipped with transmitters.

20.4 Recommended Distances between Portable and Movable Communications Equipment and Equipment that does represent not a Life-Support Equipment or System

Table 4 – Recommended Distances between Portable and Movable Communications Equipment and EQUIPMENT that does represent not a LIFE-SUPPORT EQUIPMENT or SYSTEM

Max output power (W)	Distance (m) from 150 kHz to 80 MHz $D = (3,5/\sqrt{1})(P \text{ sq. root})$	Distance (m) from 150 kHz to 80 MHz $D = (3,5/\sqrt{E1})(P \text{ sq. root})$	Distance (m) from 150 kHz to 80 MHz $D = (7/\sqrt{E1})(P \text{ sq. root})$
0.01	0.1166	0.1166	0.2333
0.1	0.3689	0.3689	0.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.333

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