Medtronic

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
Medtronic

Medtronic Xomed, Inc.
6743 Southpoint Drive North
Jacksonville, Florida 32216, USA
www.medtronic.ru

tel 800 874 5797

Custodian:

Gabriela B. Moreira
Regulatory Affairs Manager
Medtronic Xomed, Inc.

Dec. 14, 2018

Date

IFU


NuVent EM Sinus Dilation System

ATTESTATION


STATE OF Florida
COUNTY OF DUVAL

On this Friday, December 14, 2018, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following document, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Gabriela B. Moreira, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Notary:

Official Notary Signature

United States of America
State of Florida
County of Duval

Subscribed and sworn to before me this 20 day of December 2018

2018


NuVent™

EM Sinus Dilation System

For Use with Medtronic EM Computer-Assisted Surgery Systems

1830717FRT, 1830617FRT, 1830517FRT, 1830707MAX, 1830607MAX,
1830507MAX, 1830717SPH, 1830617SPH, 1830517SPH, 18INFKIT

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








EN

Indications for use

The NuVent EM Sinus Dilation System is intended for use in sinus procedures when surgical navigation or image-guided surgery may be necessary to locate and move tissue, bone or cartilaginous tissue surrounding the drainage pathways of the frontal, maxillary, and sphenoid sinuses.

The EM Sinus Dilation system is used in conjunction with the Medtronic computer-assisted surgery system.

The Medtronic computer-assisted surgery system and its associated applications are intended as an aid for precisely locating anatomical structures in either open or percutaneous procedures. Their use is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as the skull, can be identified relative to a CT- or MR-based model, or digitized landmarks of the anatomy.

The system and its associated applications should be used only as an adjunct for surgical guidance. They do not replace the surgeon's knowledge, expertise, or judgment.

Device description

The NuVent EM Sinus Dilation System comprises sterile, single-use instruments that combine electromagnetic (EM) "plug and play" tracking capability with the pathway expansion effects of balloon dilation technology and an inflator. Each of the three types of sinus seekers (frontal, maxillary and sphenoid) has a unique shape and angle that allows for entry into the sinus outflow tract. The inflator consists of a plunger, barrel and extension tube.

Each sinus seeker is intended for use in conjunction with the Fusion software on a Medtronic computer-assisted surgery system. Inside each sinus seeker is an EM tracker. The emitter on the EM Computer-Assisted Surgery System generates a low-energy magnetic field to locate the tracker mounted inside of the sinus seeker. The software displays the location of the sinus seeker's tip within multiple patient image planes and other anatomic renderings. After confirmation of placement, the sinus seeker's balloon can be inflated with saline solution by using the inflator to expand the outflow of the targeted sinus.

Contraindications

None known.

Warnings

- **The sinus seekers are fixed, rigid instruments that are not designed to bend. Do not bend or alter the instrument, as navigational inaccuracy will occur.**
- Do not open sterile barrier packages or containers until surgical use.
- Deflate balloon prior to removal from anatomy.
- Only use device in conjunction with an endoscope.
- Only use sterile saline as inflation solution.
- Do not reuse the seeker or inflator, as these are single patient use devices.
- Do not submerge the sinus seeker.
- Only use the sinus seeker in conjunction with Medtronic Inflator Part Number: 18INFKIT.
- Do not place metal objects in magnetic field.
- Do not use the sinus seeker to perforate thick bone.
- Use of the incorrect balloon size may cause damage to the sinus balloon and / or surrounding anatomy.
- Use of an undersized sinus balloon may result in failure to properly treat the target sinus.
- Prior to use, examine the product packaging for damage, deterioration, and expired shelf life. Replace it with an unused system and contact Medtronic Technical Service.
- Only use device in conjunction with Medtronic computer-assisted surgery system.
- You must re-register the patient if you replace the EM ENT Patient Tracker during any part of the procedure.
- Do not hold the seeker at the proximal end of the handle or apply excessive force to the handle, as this may cause navigational inaccuracy. Refer to Figure 7.

Instruction for use

The NuVent EM Sinus Dilation System includes an inflator and a sinus seeker (Figure 1). For proper operation, follow the proceeding instructions which cover inflator preparation, priming, inflator/sinus seeker connection, and EM Computer-Assisted Surgery System.

Preparing the inflator

1. Inspect product packaging prior to use. Do not use product if packaging is damaged.
2. Remove the inflator tray from the sterile pouch and sterile tray.
3. Connect the extension tube to the Luer fitting of the inflator (Figure 2).
4. Fill sterile container with saline solution. Place the extension tube in the container.

Priming the inflator

5. With the inflator pointed up and the extension tube immersed in the saline solution, prime the inflator by pulling and pushing on the inflator plunger, drawing saline solution into and out of the inflator (Figure 3). Repeat at least five times to remove all air bubbles from the extension tube and inflator.
6. Inspect the extension tube and inflator for air bubbles. If air bubbles are still present, prime the inflator again. Repeat until air bubbles are no longer present in the inflator and extension tube.
7. Remove the extension tube from the sterile saline solution and set the inflator aside while preparing the sinus seeker.

Connecting the sinus seeker to the EM computer-assisted surgery system

8. Surgeon should use patient CT scan (hard copies or on Medtronic computer-assisted surgery system) to determine the correct balloon size based on the bony margins of the outflow tract.
9. Remove the tray from the sterile pouch.
10. Place the tray on the sterile preparation surface with the connector and connecting wire on top.
11. Without separating trays, remove the connector and connecting wire from the top tray (Figure 4).
12. Prior to use, examine the connector and cable for damage and deterioration. Do not use product if damaged.
13. Insert the connector to the EM Computer-Assisted Surgery System and unwrap wire.
14. Remove the top tray.

Connecting the sinus seeker to the inflator

15. Grasp the sinus seeker handle in the recessed area indicated by the two arrows on the tray and remove sinus seeker from the bottom tray (Figure 5).
16. Connect the extension tube to the luer fitting of the sinus seeker (Figure 6).

Checking the sinus seeker/EM computer-assisted surgery system recognition

17. Confirm that the fusion software displays the correct sinus seeker.
18. Prior to use, confirm accuracy of the sinus seeker by touching the tip of the sinus seeker to several bony anatomical landmarks and observing the cross hairs on the monitor.
19. If accurate, proceed with navigation.
20. If accuracy is not acceptable, visually inspect the seeker. If the seeker is bent or damaged, replace the sinus seeker and contact Medtronic Technical Services. If accuracy is still not acceptable with a new seeker, re-register the patient.

For additional information about the Fusion software, refer to the Fusion ENT Pocket Guide (9734349).

Dilating the sinus ostia

21. Hold seeker at the distal end of the handle (Figure 8). Do not use excessive force.
22. Insert the prepared sinus seeker into the desired anatomical location using an endoscope and EM navigation. The probe description and balloon diameter will be displayed on the Fusion software.
23. Confirm the proper positioning of the tip of the sinus seeker using EM navigation. The cross hairs on the screen indicate the sinus seeker's tip.
24. Squeeze the inflator plunger into the inflator over the course of at least two seconds until the inflator indicator protrudes (Figure 9). Hold inflation for up to five seconds. This action inflates the balloon to approximately 10 atm (Figure 10). If the indicator does not protrude, prime the inflator again.
25. Upon completion of up to two dilations per sinus ostium, deflate the balloon before removing the sinus seeker from patient's anatomy by pulling back on the inflator handle until it stops.

Note: Only remove the sinus seeker when the balloon is completely deflated.

26. Balloon may be wrapped with a gentle twisting motion to reduce its profile before re-insertion (Figure 11).
27. Prime the inflator again when changing sinus seekers.

Sterilization

Discard the seeker, inflation device and extension tube after use in accordance with all applicable regulations.

The NuVent EM Sinus Dilation System and the NuVent Inflator with Extension tube are provided in a sterile condition and are single-use only. Do not reprocess, re-sterilize, or re-use this device. Attempts to reprocess or re-sterilize the device may be ineffective and may compromise its structural integrity. Any re-use creates a risk of contamination which could result in patient injury, illness, or death.

International customer service

International customers should contact their local Medtronic Xomed office.

Note: When contacting Customer Services and Technical Support, please have the appropriate product number, date of purchase, lot number, and nature of inquiry available.

Limited warranty

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
- (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
- (3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.

(4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.

(5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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ATTESTATION


STATE OF Florida
COUNTY OF DUVAL

On this Thursday, June 13, 2019, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following document, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Gabriela B. Moreira, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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Subscribed and sworn to before me this 13th day of June 2019

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Inflation pressure supplied to balloon: the lowest pressure – 8 atm.; nominal pressure – 10 atm.; the highest pressure – 12 atm. Tolerance ± 0.1 atm (10.1

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All unspecified tolerances for size in mm are assumed to be:

x: x.x ± 0.1

x.xx ± 0.01

x.xxx ± 0.005

x.xxxx ± 0.0005.

The volume element is a solid plastic rod, with a different length for each seeker configuration. Each volume element has approximately the same volume as the balloon in a given seeker. In this way the volumetric sum of the balloon plus the volume element will remain constant across all the seeker configurations. This results in a more consistent inflation pressure, and a simple inflator design.

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Each seeker is factory calibrated and the distance from the CTA to the distal tip of the device is stored in the seeker's connector, which the user will connect to the Navigation System. The emitter on the system generates a low energy EM field which allows the system to locate the probe in the 3-dimensional space based on the location of the CTA and its regulated the distance from the instrument tip.

System Accuracy:

95% confidence and 99.5% reliability of navigated tip accuracy of less than 2.82 mm:

Confidence level - 95%		
Percent of population in interval - 99.5%		
Statistics		
N	Mean	StDev
9	1.174 mm	0.353
95% Upper Tolerance Bound		
Normal Method	Nonparametric Method	Achieved Confidence
2.784 mm	x[9]	4.4%

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
Description of all accessories, medical devices or non-medical devices, intended for use with the medical device submitted for registration (if applicable)

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The balloon seeker and inflator are intended for use with to the Fusion Navigation System (Registration Certificate № FSZ 2009/03987). The navigation system displays images of patient anatomy derived from CT-scan data. The figure below shows the connections.

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© 2014 Medtronic, Inc.