Serial No. 1000
Form No. 6

CERTIFICATION OF COPY

I, the undersigned, Asher Afriat,adv, Notary holding license no. 2063591, hereby certify that I was presented with:

a copy of the original document

Drawn up in the English language in which I am fluent.

I confirm, with my signature and seal, that the attached document marked with the letter 'A' is an exact copy of the document or information presented to me.

Today 26/05/2020.

Asher Afriat,Adv

Notary's seal

Notary fee	91.00NIS
VATax	15.47NIS
Total	106.00NIS
Required	

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Instructions For Use (IFU)

MicronJet 600 Microneedle for intradermal injections

NanoPass Technologies Ltd.
3 Golda Meir St., Nes Ziona 7403648, Israel

2020

English

The MicronJet 600

For any further assistance/support/questions, please contact the manufacturer.

Manufacturer:

NanoPass Technologies Ltd.
3 Golda Meir St., Nes Ziona 7403648, Israel
Tel 972-8-9462905
Fax 972-50-8966787
info@nanopass.com

EU Representative:

MedNet GmbH
Borkstrasse 10-48163, Muenster, Germany
Phone 49-251-322660
Fax 49-251-3226622

Authorized Representative in Russia:

Limited Liability Company «RAIFARM» (RAIFARM LLC),
Ul. Malaya Dmitrovka 4, office 8,
Moscow, Russia 127006
Tel. +7 495 937-56-08, Fax +7 495 937-56-14

MicronJet 600 Instructions for Use, PK_3033 Rev.01, Jan 2020

Description:

The MicronJet 600 is a disposable micro-needle device that is mounted on a standard syringe in the same manner as a conventional needle.

The MicronJet 600 comprises an integral array of 3 micro-needles; each micro-needle 0.6mm in length.

The MicronJet 600 is used to inject liquid substances, allowing for controlled intradermal delivery in any procedure which requires administration of substances to the dermal compartment.

The MicronJet 600 is latex free

Field of use: General practice

Conditions of use: in hospital and clinical conditions.

Indications:

The MicronJet 600 is intended for intradermal injections of any substance or drug approved for delivery by this delivery route. For infants, child and adult patients.

Contraindications:

- The MicronJet 600 should not be used on skin abrasions, open wounds, cuts and scars.
- The MicronJet 600 should not be used on rashes, skin infections or any other area of damaged or diseased skin.

Warnings and precautions:

- The MicronJet 600 is intended for single use, by healthcare professionals
- Re-use may lead to infection or other illness/injury
- Used or un-used devices should not be re-capped or removed from syringe unless there is no alternative or such action is required by a specific medical procedure.
- The MicronJet 600 is not intended for aspiration of liquids
- If the device appears to be damaged or broken, discard and use a new device

Adverse events:

Local adverse reactions, including edema, erythema and discoloration of skin at the injection site have been observed, most likely related to the injected substance, and not to the MicronJet 600.

Adverse reactions related to use of the MicronJet 600 include small risk of micro bleeding (self-limiting blood dotting, ceasing within normal clotting time). Other potential device related adverse events may rarely include local infection at the injection site.

Storage and handling:

- MicronJet 600 should be stored at room temperature. In accordance with USN-NF room temperature is 15° C – 30°C. Guaranteed shelf-life – 5 years.
- The MicronJet 600 is individually packed in a soft blister, sterilized by ethylene oxide (EtO).
- Do not use the MicronJet 600, in the event that the soft blister is either open and/or damaged.
- Dispose of all devices in accordance with national/institutional guidelines, and in their absence, discard in a sharps container
- Device can be transported by any closed vehicles as general cargo, in accordance with applicable transportation rules.

Shelf-life of the medical device:

The shelf-life for the specified storage conditions is 5 years. During this period of time devices should be stored at room temperature (15°C – 30°C).

The device should be used immediately after opening and discarded immediately after use.

Interactions with other medical devices: The device should be mounted on a standard syringe in the same manner as a conventional needle. Harmful interactions with other medical devices are not known.

Requirements for the maintenance and repairs of the medical device:

MicronJet 600 contains no serviceable parts, and is not subject to maintenance and repairs.

Special measures for disposal of the medical device:

Used devices should be immediately disposed as potentially infected sharp medical products. Medical personnel must follow the national/institutional guidelines for medical waste management and universal precaution guidelines.

Warning: single-use device! Discard in a sharps container!

Non-used expired devices should be disposed as non-infected sharp medical products.

Supply method:

The device is supplied as a single sterile MicronJet 600 device, packed in a soft blister for single use. Each dispenser carton contains 100 units.

Transportation and storage conditions for the medical device:

Storage - MicronJet 600 and NANOSOFT products should be stored at room temperature (United States Pharmacopeia-National Formulary (USP-NF) between 15-30°C (59-86°F). Warranty storage period – 5 years.

Transportation - The product can be transported by any means of covered transport in accordance with the transportation rules applicable to this type of transport.

Vibration and drop test were performed as part of packaging design verification according to ASTM D-999-08 and ASTM 5276-98. The purpose of these tests was to verify that MicronJet 600 blister package is resistant to transportation vibration. Blister primary packages were visually inspected prior to secondary package and testing. Blister primary package and dispenser secondary package were inspected following testing. Results of visual test were found to meet the acceptance criteria.

Warranty:

This product was developed, manufactured, tested, and packaged in conformity with all applicable requirements.


The manufacturer's warranty applies to the health care product before the expiry date specified on the package, provided that it has not been opened or damaged and that it has been correctly transported and stored.

The shelf-life for the specified storage conditions is 5 years. During this period of time devices should be stored at room temperature (15° C – 30°C).

Principal technical characteristics of the "MicronJet 600" medical device Microneedle:

Characteristic	Specification
Height (length) of microneedle	600 ± 30 µm (0,6 mm)
Microchip length	2300 ± 20 µm
Microchip width	880 ± 20 µm
Spaces between microneedles	600 ± 20 µm
Pore diameter	61 (+5, - 10) µm
Adapter outer diameter (A)	6 ± 0,1 mm
Cap outer diameter (B)	7,5 ± 0,1 mm
Microneedle adapter length (C)	19 ± 0,1 mm
Protective cap length (D)	22,5 ± 0,1 mm
The length of the device assembly (the cap is put on the adapter)	26,5 ± 0,1 mm
Female Luer connector	As specified in ISO 80369-7:2016
Needles integrity	No broken needles & edges
Integrity of adhesive layer	No cracks
Holes opening	No needle blocking
Luer-chip separation pressure	≥ 2000 kPa (20 bar)

Assembled device flow rate	74.6±4.5 mm ³ /sec (no less than 31G ½ “ TSK flow rate (according to ISO 7864:2016)
Microneedle Chip material	Single Crystal Silicon coated with Silicon Dioxide, no trade name, (Manufacturer: Micro to Nano Co., Ltd. (M2N), Korea)
Female Luer (adapter) material	Granular medical grade polycarbonate, Makrolon RX 2530 (Manufacturer: Bayer Material Science, Germany)
Blue line	Transparent dark blue ink, BL1411 (Manufacturer: amcor, UK)
Adhesive	Loctite® 3211, (Manufacturer: Henkel Technologies, Germany)
Cover material	Polypropylene 13R9A (Manufacturer: Flint Hills Resources, USA)

Labeling symbols definitions

Name of the medical device	MicronJet 600
Quantity	100 pcs
Product type	Microneedles for Intradermal Injection
Name and address of the manufacturer	
Lot number	
Reference number	
Expiry date	
The word “Sterile” and the type of sterilization: Ethylene Oxide	
Do not use if the package has been damaged	
Refer to the Instruction for Use	
Do not re-use	
Storage conditions	Store at room temperature
Non-pyrogenic	Non-pyrogenic




Applicable Standards:

- Conforms to ISO 13485: 2016, which is certified by Certificate No. 2121377 EN ISO **13485: 2016**
- Conforms to Council Directive 93/42 / EEC relating to medical devices, as certified **by EC** Certificates 2120317CE01
- The risk assessment report for the medical device was prepared in accordance with ISO **14971: 2012**.
- The product partially complies with ISO 7864-2016 due to design features

Russian Federation national standards product complies with:

- GOST R 50444-92
- GOST ISO 11607-1-2018
- GOST R ISO 15223-1-2014
- Partially complies with GOST ISO 7864-2011

Injection Guide

	1 <u>Device preparation for injection</u> • Draw injection solution into the syringe without purging • Open the soft blister and firmly connect the MicronJet 600™ to the syringe. • Remove the plastic cover. • Purge air and excess fluid leaving the required volume in the syringe. • Make sure injection site is dry prior to injection
	2 • Hold syringe flanges between thumb and fingers • The blue line must face you at all times • Syringe flanges are not the orientation guide as they may not always be parallel to the blue line; use the blue line for orientation
	3 • Stretch the skin downwards throughout the entire injection • Insert the device at a 45° angle • Press the device into the skin until a shallow indentation is formed (refer to FAQ for further information)
	4 • Inject slowly, expect strong resistance • A white wheal (bleb) will form at the injection site • Upon completion, wait 3 seconds and remove syringe.

Instructions For Use (IFU)

NANOSOFT Microneedle for intradermal injections

NanoPass Technologies Ltd.
3 Golda Meir St., Nes Ziona 7403648, Israel

2020

English

NANOSOFT

For any further assistance/support/questions, please contact the manufacturer.

Manufacturer:

NanoPass Technologies Ltd.
3 Golda Meir St., Nes Ziona 7403648, Israel
Tel 972-8-9462905
Fax 972-50-8966787
info@nanopass.com

EU Representative:

MedNet GmbH
Borkstrasse 10-48163, Muenster, Germany
Phone 49-251-322660
Fax 49-251-3226622

Authorized Representative in Russia:

Limited Liability Company «RAIFARM» (RAIFARM LLC),
Ul. Malaya Dmitrovka 4, office 8,
Moscow, Russia 127006
Tel. +7 495 937-56-08, Fax +7 495 937-56-14

Description:

The NANOSOFT is a disposable micro-needle device that is mounted on a standard syringe in the same manner as a conventional needle.

The NANOSOFT comprises an integral array of 3 micro-needles; each micro-needle 0.6mm in length.

The NANOSOFT is used to inject liquid substances, allowing for controlled intradermal delivery in any procedure which requires administration of substances to the dermal compartment.

The NANOSOFT is latex free

Field of use: General practice

Conditions of use: in hospital and clinical conditions.

Indications:

The NANOSOFT is intended for intradermal injections of any substance or drug approved for delivery by this delivery route. For infants, child and adult patients.

Contraindications:

- The NANOSOFT should not be used on skin abrasions, open wounds, cuts and scars.
- The NANOSOFT should not be used on rashes, skin infections or any other area of damaged or diseased skin.

Warnings and precautions:

- The NANOSOFT is intended for single use, by healthcare professionals
- Re-use may lead to infection or other illness/injury
- Used or un-used devices should not be re-capped or removed from syringe unless there is no alternative or such action is required by a specific medical procedure.
- The NANOSOFT is not intended for aspiration of liquids
- If the device appears to be damaged or broken, discard and use a new device

Adverse events:

Local adverse reactions, including edema, erythema and discoloration of skin at the injection site have been observed, most likely related to the injected substance, and not to the NANOSOFT.

Adverse reactions related to use of the NANOSOFT include small risk of micro bleeding (self-limiting blood dotting, ceasing within normal clotting time). Other potential device related adverse events may rarely include local infection at the injection site.

Storage and handling:

- NANOSOFT should be stored at room temperature. In accordance with USN-NF room temperature is 15° C – 30°C. Guaranteed shelf-life – 5 years.
- The NANOSOFT is individually packed in a soft blister, sterilized by ethylene oxide (EtO).
- Do not use the NANOSOFT, in the event that the soft blister is either open and/or damaged.
- Dispose of all devices in accordance with national/institutional guidelines, and in their absence, discard in a sharps container
- Device can be transported by any closed vehicles as general cargo, in accordance with applicable transportation rules.

Shelf-life of the medical device:

The shelf-life for the specified storage conditions is 5 years. During this period of time devices should be stored at room temperature (15°C – 30°C).

The device should be used immediately after opening and discarded immediately after use.

Interactions with other medical devices: The device should be mounted on a standard syringe in the same manner as a conventional needle. Harmful interactions with other medical devices are not known.

Requirements for the maintenance and repairs of the medical device:

NANOSOFT contains no serviceable parts, and is not subject to maintenance and repairs.

Special measures for disposal of the medical device:

Used devices should be immediately disposed as potentially infected sharp medical products. Medical personnel must follow the national/institutional guidelines for medical waste management and universal precaution guidelines.

Warning: single-use device! Discard in a sharps container!

Non-used expired devices should be disposed as non-infected sharp medical products.

Supply method:

The device is supplied as a single sterile NANOSOFT device, packed in a soft blister for single use. Each dispenser carton contains 30 units.

Transportation and storage conditions for the medical device:

Storage - MicronJet 600 and NANOSOFT products should be stored at room temperature (United States Pharmacopeia-National Formulary (USP-NF) between 15-30°C (59-86°F). Warranty storage period – 5 years.

Transportation - The product can be transported by any means of covered transport in accordance with the transportation rules applicable to this type of transport.

Vibration and drop test were performed as part of packaging design verification according to ASTM D-999-08 and ASTM 5276-98. The purpose of these tests was to verify that MicronJet 600 blister package is resistant to transportation vibration. Blister primary packages were visually inspected prior to secondary package and testing. Blister primary package and dispenser secondary package were inspected following testing. Results of visual test were found to meet the acceptance criteria.

Warranty:

This product was developed, manufactured, tested, and packaged in conformity with all applicable requirements.

The manufacturer's warranty applies to the health care product before the expiry date specified on the package, provided that it has not been opened or damaged and that it has been correctly transported and stored.


The shelf-life for the specified storage conditions is 5 years. During this period of time devices should be stored at room temperature (15° C – 30°C).

Principal technical characteristics of the "NANOSOFT" medical device

Microneedle:

Characteristic	Specification
Height (length) of microneedle	600 ± 30 µm (0,6 mm)
Microchip length	2300 ± 20 µm
Microchip width	880 ± 20 µm
Spaces between microneedles	600 ± 20 µm
Pore diameter	61 (+5, - 10) µm
Adapter outer diameter (A)	6 ± 0,1 mm
Cap outer diameter (B)	7,5 ± 0,1 mm
Microneedle adapter length (C)	19 ± 0,1 mm
Protective cap length (D)	22,5 ± 0,1 mm
The length of the device assembly (the cap is put on the adapter)	26,5 ± 0,1 mm
Female Luer connector	As specified in ISO 80369-7:2016
Needles integrity	No broken needles & edges
Integrity of adhesive layer	No cracks

Holes opening	No needle blocking
Luer-chip separation pressure	≥ 2000 kPa (20 bar)
Assembled device flow rate	74.6 ± 4.5 mm ³ /sec (no less than 31G ½ “ TSK flow rate (according to ISO 7864:2016)
Microneedle Chip material	Single Crystal Silicon coated with Silicon Dioxide, no trade name, (Manufacturer: Micro to Nano Co., Ltd. (M2N), Korea)
Female Luer (adapter) material	Granular medical grade polycarbonate, Makrolon RX 2530 (Manufacturer: Bayer Material Science, Germany)
Blue line	Transparent dark blue ink, BL1411 (Manufacturer: amcor, UK)
Adhesive	Loctite® 3211, (Manufacturer: Henkel Technologies, Germany)
Cover material	Polypropylene 13R9A (Manufacturer: Flint Hills Resources, USA)

Symbols definitions

Name of the medical device	NANOSOFT
Quantity	30 pcs
Product type	Microneedles for Intradermal Injection
Name and address of the manufacturer	
Lot number	 LOT
Reference number	 REF
Expiry date	
The word “Sterile” and the type of sterilization: Ethylene Oxide	 STERILE EO
Do not use if the package has been damaged	
Refer to the Instruction for Use	
Do not re-use	
For healthcare professional use	Rx Only

Non-pyrogenic symbol	
Storage conditions	Store at room temperature

Applicable Standards:

- Conforms to ISO 13485: 2016, which is certified by Certificate No. 2121377 EN ISO 13485: 2016
- Conforms to Council Directive 93/42 / EEC relating to medical devices, as certified by EC Certificates 2120317CE01
- The risk assessment report for the medical device was prepared in accordance with ISO 14971: 2012.
- The product partially complies with ISO 7864-2016 due to design features

Russian Federation national standards product complies with:

- GOST R 50444-92
- GOST ISO 11607-1-2018
- GOST R ISO 15223-1-2014
- Partially complies with GOST ISO 7864-2011


Injection Guide

Device preparation for injection


- Draw injection solution into the syringe without purging
- Open the soft blister and firmly connect the NANOSOFT™ to the syringe.
- Remove the plastic cover.
- Purge air and excess fluid leaving the required volume in the syringe.
- Make sure injection site is dry prior to injection

1


- Hold syringe flanges between thumb and fingers.
- The blue line must face you at all times.
- Syringe flanges are not the orientation guide as they may not always be parallel to the blue line; use the blue line for orientation.

2

- Stretch the skin downwards throughout the entire injection.
- Insert the device at a 45° angle.
- Press the device into the skin until a shallow indentation is formed.

3

- Inject slowly, expect strong resistance.
- A white wheal (bleb) will form at the injection site.

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