Eurotrol

Instruction for use


Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations

Manufactured by Eurotrol B.V., the Netherlands

Approved by

Eurotrol B.V.


Date: 2020-10-16

(Signature)

Daniel Philippens

QA/RA Director

Stamp

Eurotrol B.V.
PO Box 722, 6710 BS Ede
Keplerlaan 20, 6716 BS Ede
The Netherlands
www.eurotrol.com
T: +31 318 695777

Eurotrol B.V. · P.O. Box 722, 6710 BS Ede · Keplerlaan 20, 6716 BS Ede · The Netherlands
+31 318 695777 · office@eurotrol.com · www.eurotrol.com

BIC ABNANL2A · IBAN NL92ABNA0561838909 · Trade register Arnhem 09067367 · VAT No. NL009139448B01

Eurotrol

Eurotrol GlucoTrol-NG

Intended purpose/Intended use

Eurotrol GlucoTrol-NG is an assayed glucose control intended for professional use in the verification of the precision and accuracy of the HemoCue Glucose systems.

IVD Medical Device

Eurotrol GlucoTrol-NG complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark. Eurotrol GlucoTrol-NG is for in vitro diagnostic use only.

Summary

Eurotrol GlucoTrol-NG is a quality control material with a known glucose concentration. Eurotrol GlucoTrol-NG may be used in the periodic verification of the precision and accuracy of the HemoCue Glucose systems. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol GlucoTrol-NG is produced at four physiologically relevant levels (Level 1, 2, 3, 4) for the Russian market. Each vial contains 1.0 mL of purified bovine plasma with glucose.

Storage and stability

Eurotrol GlucoTrol-NG must be stored in the refrigerator at 2-8 °C (35-46 °F). Stored unopened at this temperature it is guaranteed stable as indicated until the expiration date on the outer box. After opening the vial, it is stable for 30 days when properly recapped and stored at 2-8 °C (35-46 °F).

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15-30 °C (59-86 °F).
2. Gently mix the vial 8-10 times before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the cuvette in the analyzer for immediate analysis.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. Store in the refrigerator at 2-8 °C (35-46 °F).
6. If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol GlucoTrol-NG is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing of vials.

Eurotrol B.V. · P.O. Box 722, 6710 BS Ede · Keplerlaan 20, 6716 BS Ede · The Netherlands

+31 318 695777 · office@eurotrol.com · www.eurotrol.com

BIC ABNANL2A · IBAN NL92ABNA0561838909 · Trade register Arnhem 09067367 · VAT No. NL009139448B01

APOSTILLE

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1. Country: THE NETHERLANDS
This public document
2. has been signed by **W.W. Blankestijn**
3. acting in the capacity of official of the Chamber of
Commerce and Industry
4. bears the seal/stamp of aforesaid Chamber of
Commerce and Industry

Certified

5. in Arnhem
6. on 22-10-2020
7. by the registrar of the district court of Gelderland
8. no. 20-2305
9. Seal/stamp:
10. Signature:

T.J. van der Linden

Eurotrol

4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol GlucoTrol-NG is not to be used as a calibrator.

Expected results

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Glucose systems. The values are only valid on the HemoCue Glucose systems as referred to in the assay sheet. Values for "HemoCue whole blood systems" are applicable for HemoCue Glucose systems reporting whole blood values. Values for "HemoCue plasma equivalent systems" are applicable for HemoCue Glucose systems reporting plasma equivalent values. Results obtained should fall within the expected range specified for the control. Not all HemoCue Glucose systems are available world wide.

For further information or technical assistance, please contact Eurotrol B.V.:

Phone +31 318 695 777

Fax +31 318 695 770

E-mail office@eurotrol.com

www.eurotrol.com

Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands

Quality controls for IN VITRO diagnostics

Symbols used

Attention, see instructions for use

Use by

In Vitro Diagnostic Medical Device

Manufacturer

Batch code

Temperature limitation

Reference number

CE mark

AN01064A01 20180130

Eurotrol B.V. · P.O. Box 722, 6710 BS Ede · Keplerlaan 20, 6716 BS Ede · The Netherlands

+31 318 695777 · office@eurotrol.com · www.eurotrol.com

BIC ABNANL2A · IBAN NL92ABNA0561838909 · Trade register Arnhem 09067367 · VAT No. NL009139448B01

Eurotrol

Eurotrol GlucoTrol-NG

Level 1

Valuesheet (example)

146.001.002

907 60

2020-07

		HemoCue whole blood «B-Glu»		HemoCue plasma equivalent «P-Equiv»	
HemoCue Glucose systems	Measuring unit	Target	Expected range	Target	Expected range
HemoCue Glucose 201 systems	mmol/L	2.4	1.6-3.2	2.6	1.8-3.4
	mg/dL	43	29-57	47	33-61
HemoCue Glucose 201RT systems	mmol/L	2.6	1.8-3.4	2.9	2.1-3.7
	mg/dL	47	33-61	52	38-66

Barcode for HemoCue Glucose whole blood systems (B-Glu)

Barcode for HemoCue plasma equivalent systems (P-Equiv)

Barcode for HemoCue Glucose 201RT whole blood systems

Barcode for HemoCue Glucose 201RT plasma equivalent systems

AN01058A01 20180130

AN01407A01 20180330

Eurotrol B.V. · P.O. Box 722, 6710 BS Ede · Keplerlaan 20, 6716 BS Ede · The Netherlands

+31 318 695777 · office@eurotrol.com · www.eurotrol.com

BIC ABNANL2A · IBAN NL92ABNA0561838909 · Trade register Arnhem 09067367 · VAT No. NL009139448B01

Eurotrol

Addendum to the operational documentation*
of medical device

Instruction for use

Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations

2020

* Available on request

Eurotrol B.V. · P.O. Box 722, 6710 BS Ede · Keplerlaan 20, 6716 BS Ede · The Netherlands
+31 318 695777 · office@eurotrol.com · www.eurotrol.com
BIC ABNANL2A · IBAN NL92ABNA0561838909 · Trade register Arnhem 09067367 · VAT No. NL009139448B01

Eurotrol

Table of Contents

Intended Use	7
Composition	7
Key consumer features (not part of the manufacturer's review)	8
Service maintenance and repair	9
The device does not require maintenance and cannot be repaired.	9
Analytical and diagnostic performance	9
Precision performance	9
Metrological traceability	9
Value ranges	10
Mode of application	10
Procedures/Instruction for use	11
Shelf life	12
Storage and transportation	12
Utilization	13
Standards	13

Eurotrol

Intended Use

Eurotrol GlucoTrol-NG is an assayed glucose control intended for professional use in the verification of the precision and accuracy of the HemoCue Glucose systems.

In Vitro Diagnostic Medical Device (IVD)

Eurotrol GlucoTrol-NG complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Eurotrol GlucoTrol-NG is for in vitro diagnostic use only.

Overview

Eurotrol GlucoTrol-NG is a quality control material with a known glucose concentration. Eurotrol GlucoTrol-NG can be used to periodically monitor the precision and accuracy of HemoCue Glucose systems. Perform control measurements in accordance with local regulations.

Description of the analyte, data regarding its scientific substantiation (not part of the manufacturer's review)

Glucose is one of the most important and basic carbohydrates of peripheral blood, a source of energy supply for cells and tissues of the body. The concentration of glucose in the blood is regulated by a group of hormones, the most important of which are produced by the pancreas (insulin, glucagon and others). The level of glucose in the blood reflects the state of carbohydrate metabolism, and its changes (hyper- and hypoglycemia) are associated with impaired hormonal regulation, primarily with a decrease in insulin levels or the effectiveness of the manifestation of its regulatory effect. An increase in blood glucose is hyperglycemia, the main symptom of such a serious and widespread disease as diabetes. The glucose content in the blood of patients with diabetes mellitus varies significantly during the day, depends on food intake, stressful situations, physical exercises, changes with insulin therapy, glucose-lowering drugs. In addition, the concentration of glucose in the blood can change with pancreatitis, thyroid dysfunction, renal failure and some liver diseases. It is known that determining the level of glucose in the blood is the most required test in the menu of clinical diagnostic laboratories (CDL).

Reagents

Eurotrol GlucoTrol-NG is produced with four physiologically relevant glucose levels (Level 1, 2, 3, 4) for the Russian market. Each vial contains 1.0 ml of purified bovine plasma with glucose.

Composition

Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations

Variations:

1. Level 1, 2 x 1.0 ml in a package, includes:

- Vial with a red cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

Eurotrol

2. Level 2, 2 x 1.0 ml in a package, includes:

- Vial with a white cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

3. Level 3, 2 x 1.0 ml in a package, includes:

- Vial with a blue cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

4. Level 4, 2 x 1.0 ml in a package, includes:

- Vial with a green cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

Values of quality controls and acceptable ranges are indicated in the certificate of analysis for each lot of the product.

Note:

1. Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations are made with the use of blood products of animal origin. Bovine materials are not biohazardous to humans.

Among other things, it excludes the possibility of containing the surface antigen of hepatitis B virus (HBsAg), hepatitis C virus (HCV) and human immunodeficiency virus (HIV-1 / HIV-2). This material does not contain the causative agents of transmissible spongiform encephalopathy (TSE).

2. The product does not contain materials that come into direct or indirect contact with the body of the patient and personnel using the product when meeting the requirements of operational documentation (instructions for use).

3. The product contains no medicines or pharmaceutical substances.

Key consumer features (not part of the manufacturer's review)

Each level of QC is designed for 50 tests.

Eurotrol

Sterilization method of the device

The device does not need sterilization.

Software

N/A.

Service maintenance and repair

The device does not require maintenance and cannot be repaired.

Principal Method

The principle of operation (principle of the analytical method) of Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations is based on measuring the amount of glucose contained in control solutions on a HemoCue Glucose 201+ analyzer in the Quality Control mode. The obtained measurement results are compared with the certified values of the Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations of execution, as a result of which the functional characteristics of the HemoCue Glucose 201+ analyzer are determined.

This is measured using a dual wavelength spectrophotometric technology. Measurements are performed at 660 nm and 840 nm.

Analytical and diagnostic performance

Precision performance

The glucose concentration of all four levels of Eurotrol GlucoTrol-NG was measured frequently during a period of 20 (working) days on the HemoCue target instruments (HemoCue Glucose 201+ analyzer for Glucose measuring and HemoCue «Glucose 201» microcuvettes). The glucose concentration was measured in duplicate twice a day. For each day new dropper bottles were used from the same batch. Dropper bottles were stored according to the instructions for use of Eurotrol GlucoTrol-NG.

Based on the results of this evaluation, we conclude that the within-run precision and total precision is within the clinical acceptable performance criteria (<11% (CV%)).

Metrological traceability

At value assignment the mean assay values are determined and are derived from measurements performed on 6 HemoCue Glucose 201+ analyzers and 6 HemoCue Glucose 201RT analyzers. For value assignment 9 dropper bottles of each produced batch of GlucoTrol-NG are analyzed using 9 cuvettes (3 from each cuvette batch) per analyzer (1 cuvette per dropper bottle). HemoCue Glucose 201+ analyzer is traceable to reference quality assessment methodology - an Isotope Dilution Gas Chromatography - Mass Spectrometry (ID GC-MS) method. The technique was calibrated according to the standard of the National Institute of Standards and Technology (NIST), (hereinafter: NIST standard). The GC / MS method is monitored in an external IFCC management program for reference laboratories, thus representing an independent quality assessment procedure for reference laboratories, that is, it is a reference quality assessment technique. HemoCue calibrates reference analyzers and reference batches of microcuvettes several times a year using this reference technique. These analyzers and batches of microcuvettes are used as a reference for all manufactured analyzers and microcuvettes in HemoCue.

Eurotrol

HemoCue Glucose 201 RT analyzer is traceable to the YSI 2300 STAT Plus (manufactured by YSI Incorporated Life Sciences, USA (№ in the GRSI of the Russian Federation: 57156-14)) method and ID GC-MS method.

Value ranges

To cover the physiological relevant range of glucose and the measuring range of the HemoCue Glucose systems, Eurotrol GlucoTrol-NG is manufactured at four different levels for the Russian market.

Device name	Article number	Glucose, mg/dL (mmol/L)
Eurotrol GlucoTrol-NG QC for Glucose measuring Level 1	146.001.002	45 (2,5)
Eurotrol GlucoTrol-NG QC for Glucose measuring Level 2	146.002.002	108 (6,0)
Eurotrol GlucoTrol-NG QC for Glucose measuring Level 3	146.003.002	180 (10,0)
Eurotrol GlucoTrol-NG QC for Glucose measuring Level 4	146.004.002	315 (17,5)

The measuring range of HemoCue Glucose 201+ analyzer is 0-400 mg/dL (22,21 mmol/L) for whole blood and 0-440 mg/dL (24,43 mmol/L) for plasma. HemoCue Glucose 201+ analyzer is the system with a measuring method based on enzymatic glucose dehydrogenasis forming NADH, which in the presence of diaphorase produces a colored formazan with MTT (a tetrazolium salt). This is measured using a dual wavelength spectrophotometric technology. Measurements are performed at 660 nm and 840 nm.

Clinical Performance and Evidence

The clinical performance of a quality control such as Eurotrol GlucoTrol-NG cannot be evaluated based solely on the QC. Only the combination of the IVD analyzer, reagents and QC have a clinical performance as these need to be used, and can only be used together. Therefore, the clinical performance of the IVD analyzer applies also to the QC and no separate data of the clinical performance of the QC Eurotrol GlucoTrol-NG is available.

Test samples

Sampling is not provided.

Mode of application

Equipment and consumables (for manual testing)

- HemoCue Glucose 201+ analyzer, Registration Certificate № P3H 2016/3622 dated 03.02.2016.
- HemoCue «Glucose 201» microcuvettes for HemoCue Glucose 201+ analyzer (Registration Certificate №ФC3 2011/10948 dated 29.01.2016);
- general laboratory equipment.

Eurotrol

Procedures/Instruction for use

1. Allow the vial to stand for 15 minutes at room temperature 15-30 °C.
2. Gently mix the vial 8-10 times before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the cuvette in the analyzer for immediate analysis.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. Store QCs at 2-8 °C (35-46 °F).
6. If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

Calculated test

Concentrations of test substances are calculated automatically by the HemoCue Glucose 201+ analyzer.

Expected results

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Glucose systems. The values are only valid on the HemoCue Glucose systems as referred to in the assay sheet. Values for "HemoCue whole blood systems" are applicable for HemoCue Glucose systems reporting whole blood values. Values for "HemoCue plasma equivalent systems" are applicable for HemoCue Glucose systems reporting plasma equivalent values. Results obtained should fall within the expected range specified for the control. Not all HemoCue Glucose systems are available world wide.

Contradictions

No.

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol GlucoTrol-NG is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing of vials.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol GlucoTrol-NG is not to be used as a calibrator.

Eurotrol

Shelf life

The shelf life of Eurotrol GlucoTrol-NG QCs for Glucose measuring in variations is 18 months. Expired QCs are not to be used.

Storage and transportation

Storage

Samples of Eurotrol GlucoTrol-NG should be stored at 2-8 ° C (35-46 ° F). When these storage conditions are met, the samples remain stable until the expiration date indicated on the outer packaging. After opening the bottle, the QC remains stable for 30 days, provided that the bottle is properly closed and stored at a temperature of 2-8 ° C (35-46 ° F).

The storage of control solutions in the manufacturer's packaging should be carried out at a temperature of 2 to 8 ° C for the entire shelf life in refrigerators or refrigerators that provide a regulated temperature regime with daily temperature recording.

When stored closed at this temperature, the stability of the QC is guaranteed until the expiration date, as indicated on dropper bottles and a carton box.

Transportation

Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations, should be transported in the original packaging of the manufacturer, which protects the product from external influences, on all modes of transport in covered vehicles in accordance with the transport rules applicable to this type of transport with temperature from 2 to 8 ° C for the entire shelf life.

Manufacturer's warranty

Claims about the quality of the medical device are sent to the manufacturer or an authorized representative in the Russian Federation.

Manufacturer:

Eurotrol B.V.

Phone: +31 318 695 777

Fax: +31 318 695 770

Email: office@eurotrol.com

www.eurotrol.com

Keplerlaan 20, 6716 BS, Ede, The Netherlands

Eurotrol

Authorized Representative in Russia:

LLC «RADIOMETER»

Phone: +7 495 228 92 06

Email: moscow@radiometer.ru

Website: www.radiometer.ru, www.radiometer.com

Address: Russian Federation, 109004, Moscow, Stanislavskogo street, 21, building 3, room I

Utilization

Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations, must not be disposed of with normal waste.

In order to avoid possible environmental pollution, the waste of the used Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations, as well as QCs with expired shelf life, must be disposed of in accordance with local state and / or hospital regulations requirements.

(The following information is not part of the manufacturer's evaluation)

Used Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations and also with expired shelf life should be disposed of in accordance with SanPiN 2.1.7.2790-10 as potentially infected medical waste (class B). When disposing of Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations it is necessary to act in accordance with hygienic standards and follow the official instructions for disposal and destruction of potentially infectious and chemical material adopted in the Russian Federation and in this health facility. The key issues in the waste management strategy are the optimization of tactics and approaches to end products, as well as the minimization of possible environmental pollution.

When using Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations, class B wastes are generated, which are classified and disposed of in accordance with San PiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for handling medical waste."

Standards

The following is a list of standards in accordance with which medical device is produced: "Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations":

- CEN-EN-ISO 13485: 2016 Medical devices. Quality Management Systems. Regulatory requirements
- CEN-EN-ISO 14971: 2012 Medical devices. Application of risk management to medical devices
- CEN-EN-ISO 18113-1: 2011 Medical diagnostic products in vitro. Information supplied by the manufacturer (labeling). Part 1. Terms, definitions and general requirements
- CEN-EN-ISO 18113-2: 2011 Medical diagnostic products in vitro. Information supplied by the manufacturer (labeling) - Part 2: In vitro diagnostic professional reagents
- CEN-EN-ISO 23640: 2015 In vitro diagnostic medical equipment
- CEN-EN 13612: 2002 Laboratory evaluation of the technical characteristics of medical devices

Eurotrol

- CEN-EN-ISO 11737-2: 2009 Sterilization of medical devices. Microbiological methods. Part 2. Sterility tests carried out during validation of the sterilization process
- ISO 15223-1: 2016 Medical devices. Symbols used for marking on medical devices, labels and accompanying documentation. Part 1. Basic requirements.
- ISTA-7D: 2007 Thermally controlled shipping packaging for systematic delivery of goods
- CLSI-EP25-A: 2009 Assessment of the stability of in vitro diagnostic reagents
- CLSI-EP05-A3: 2014 Accuracy Assessment for Quantitative Measurement Methods



The following list was not part of the manufacturer's essential principles review, the conclusions regarding compliance were drawn by Radiometer LLC:

Eurotrol GlucoTrol-NG QCs for Glucose measuring, in variations also comply with the following standards:

- GOST R 52905-2007 "Medical laboratories. Safety requirements".
- SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the treatment of medical waste"
- GOST R 51088-2013 Medical products for in vitro diagnostics. Reagents, reagent kits, test systems, control materials, culture media. Product Requirements and Supporting Documentation
- GOST R ISO 15223-1-2014 Medical devices. Symbols used for labeling on medical devices, labels and accompanying documentation. Part 1. Basic requirements
- GOST 15150-69. Machines, devices and other technical products.
- Versions for different climatic regions. Categories, conditions of operation, storage and transportation regarding the impact of climatic environmental factors
- GOST R 51352-2013 Medical products for in vitro diagnostics. Test methods.
- GOST ISO 18153-2011 Medical devices for in vitro diagnostics. Measurement of quantities in biological samples. Metrological traceability of the values of the catalytic concentration of enzymes attributed to calibrators and control materials.
- GOST R ISO 23640-2015 Investigation of the stability of in vitro diagnostic reagents
- GOST R ISO 18113-2-2015 Medical products for in vitro diagnostics. Information provided by the manufacturer (labeling). Part 2. In vitro diagnostic reagents for professional use
- Order of the Ministry of Health of Russia dated 06.06.2012 No. 4n (as amended according to the Order of the Ministry of Health of the Russian Federation of September 25, 2014 No. 557n).
- Order of the Ministry of Health of Russia dated 01.09.2014, No. 2n.
- Order of the Ministry of Health of Russia dated January 19, 2017 No. 11n.

Eurotrol

Symbols used

Quality Controls for in vitro diagnostic

Eurotrol GlucoTrol-NG

Level 1

Valuesheet (example)

REF 146.001.002

LOT 907 60

2020-07

		HemoCue whole blood «B-Glu»		HemoCue plasma equivalent «P-Equiv»	
HemoCue Glucose systems	Measuring unit	Target	Expected range	Target	Expected range
HemoCue Glucose 201 systems	mmol/L	2.4	1.6-3.2	2.6	1.8-3.4
	mg/dL	43	29-57	47	33-61
HemoCue Glucose 201RT systems	mmol/L	2.6	1.8-3.4	2.9	2.1-3.7
	mg/dL	47	33-61	52	38-66

Barcode for HemoCue Glucose whole blood systems (B-Glu)

Barcode for HemoCue plasma equivalent systems (P-Equiv)

Barcode for HemoCue Glucose 201RT whole blood systems

Barcode for HemoCue Glucose 201RT plasma equivalent systems

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