State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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1. Country (Pays): **United States of America**

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2. has been signed by (a été signé par) **Christopher James Devet**

3. acting in the capacity of (agissant en qualité de) **Notary Public, State of Minnesota**

4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Christopher James Devet, Notary Public, State of Minnesota

Certified Attesté

5. at (à) **St. Paul, Minnesota**

6. Date (Date) **10/11/2021**

7. by (par) **Secretary of State, State of Minnesota**

8. File No (Sous n°) **1262835200156**

9. Seal / Stamp (Sceau / Timbre):

10. Signature (Signature):

Steve Simon

Steve Simon
Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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
2510 Conway Avenue, St. Paul,
Minnesota, 55144-1000, USA

Dear Sir/Madame,

In my capacity as Regulatory Affairs Specialist, 3M ESPE Oral Care Solutions Division,
St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Filtek™ One
Bulk Restorative in refills in English is truthful and accurate.


Regina Feferman
Regulatory Affairs Specialist
3M ESPE
Oral Care Solutions Division


October 6, 2021
Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 6 day of October 2021

CHRIS DEVET Notary Public
My Commission Expires: JANUARY 31, 2025

ctions for Use

PE Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative in ns.

General information

Filtek™ One Bulk Fill Restorative in versions, is a visible-light activated, restorative composite optimized to create fast and easy restorations. This material provides excellent strength and low wear for durability and improved esthetics. The material can be placed and cured up to 5 mm deep, enabled by a stress relieving resin system and optimized optical properties. Filtek™ One Bulk Fill Restorative in versions serves to enhance the 3M ESPE lineup of restorative materials by improving the esthetic properties of a bulk fill material to allow for broader use in both posterior and anterior restorations.

Filtek™ One Bulk Fill Restorative in versions is offered in a range of five tooth-Shade shades. All shades are radiopaque. The fillers are a combination of a non-agglomerated/non- aggregated 20 nm silica filler, a non-agglomerated/non-aggregated 4 to 11 nm zirconia filler, an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles), and an ytterbium trifluoride filler consisting of agglomerate 100 nm particles. The inorganic filler loading is about 76.5 % by weight (58.5 % by volume). Filtek™ One Bulk Fill Restorative in versions contains AUDMA, AFM, diurethane-DMA, and 1, 12-dodecane-DMA. Filtek™ One Bulk Fill Restorative in versions is applied to the tooth following use of a methacrylate-based dental adhesive, such as manufactured by 3M ESPE, which permanently bonds the restoration to the tooth structure.

Filtek™ One Bulk Fill Restorative in versions is a modification of predicate device, Filtek™ Bulk Fill Posterior Restorative, also manufactured by 3M ESPE Dental Products.

The formulation was modified by adjusting the ratios of fillers in the product to improve the esthetics of the cured restorative material. Two additional tooth-Shade shades of the restorative have also been added. No new ingredients were introduced.

Both Filtek™ One Bulk Fill Restorative in versions and predicate device Filtek™ Bulk Fill Posterior Restorative are packaged in traditional dispensers for dispensing restorative on a pad outside the mouth, and single-dose capsules for dispensing restorative intraorally. The currently marketed 3M ESPE 5707SD Restorative Dispenser is used for dispensing predicate device Filtek™ Bulk Fill Posterior Restorative and other 3M ESPE restoratives from capsules. The same 5707SD Restorative Dispenser will also be used to dispense Filtek™ One Bulk Fill Restorative capsules.

Variations of MD

Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative in versions:

1. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade A1 (20 capsules of 0.2 g), instructions for use - 2 pcs.
2. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade A2 (20 capsules of 0.2 g), instructions for use - 2 pcs.
3. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade A3 (20 capsules of 0.2 g), instructions for use - 2 pcs.
4. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade B1 (20 capsules of 0.2 g), instructions for use - 2 pcs.
5. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade C2 (20 capsules of 0.2 g), instructions for use - 2 pcs.

6. Material for restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative (trial sample) material for anterior and posterior teeth shade A2 (5 capsules 0.2 g each), instructions for use - 2 pcs.
7. Material for restoration Filtek™ One Bulk Fill Restorative material for anterior and posterior teeth shade A2 (1 pc., 4 g) in a dispenser, instructions for use - 2 pcs.
8. Material for the restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade A2 (1 pc., 4 g) in a dispenser, instructions for use - 2 pieces.
9. Material for restoration Filtek™ One Bulk Fill Restorative material for anterior and posterior teeth shade A3 (1 pc., 4 g) in a dispenser, instructions for use - 2 pcs.
10. Material for the restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade B1 (1 pc., 4 g) in a dispenser, instructions for use - 2 pieces.
11. Material for the restoration of anterior and posterior teeth Filtek™ One Bulk Fill Restorative shade C2 (1 pc., 4 g) in a dispenser, instructions for use - 2 pcs.
12. Material for restoration Filtek™ One Bulk Fill Restorative (trial sample) material for anterior and posterior teeth shade A2 (1 pc., 1 g) in a dispenser, instructions for use - 2 pcs.

Medical Device Classification (MD)

Russian Federation: 2a

Sterility: Non-sterile

Invasiveness: Non-invasive

Duration of contact: Direct contact— more than 30 days. Contact with tissues / bones / dentin.

Intended use

Dental composite material for anterior and posterior restorations.

Indications

Filtek™ One Bulk Fill Restorative in versions is indicated for use in:

- Direct anterior and posterior restorations (including occlusal surfaces)
- Base / liner under direct restorations
- Core build-ups
- Splinting
- Indirect restorations including inlays, onlays and veneers
- Restorations of deciduous teeth
- Extended fissure sealing in molars and premolars
- Repair of defects in porcelain restorations, enamel, and temporaries.

Contraindications

Avoid use of this product in patients with known acrylate allergies.

Adverse effects

There are no known side effects.

Potential consumer MD

Trained dental professionals.

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the

...trial sample
...tooth shade
...necessary and discontinue future use of the product.

Cautionary Information for Dental Personnel

Capsules may be warmed (Do not warm syringes). This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed. 3M SDS information can be obtained from 3M.com/dental or contact your local subsidiary.

Possibilities and ways of integration with other medical devices

1. Sof-Lex Diamond Polishing System is used for contouring, finishing and polishing and is registered in Russian Federation, RZN 2015/2556.
2. To bond Filtek™ One Bulk Fill Restorative in versions to tooth structure, use of Single Bond Universal Adhesive is recommended, RZN 2017/5981.
3. This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range.
4. Vitrebond Glass Ionomer can be used for pulp protection and to line areas of deep cavity excavation, FSZ 2010/07297).
5. The capsules are dispensed using the 5707SD Restorative Dispenser, which is sold separately and has been marketed globally since 2010 for dispensing 3M ESPE restorative materials from capsules:

Item	Catalogue / Order Number
Restorative Dispenser	5707SD

Principles of operation

When irradiated by light, the methacrylate functionalities of the restorative material, i.e., ERGP-DMA, DDDMA, UDMA, AFM1, and surface-treated Silica, Zirconia and Zirconia/Silica fillers undergo, in conjunction with the photoinitiator system, a light- induced polymerization to form a hard composite that is bonded to the tooth structure with a dental adhesive.

Mode of application¹

Preparation

1. **Prophy:** Teeth should be cleaned with pumice and water to remove surface stains.
2. **Shade Selection:** Prior to isolation of tooth, select the appropriate shade(s) of Filtek™ One Bulk Fill Restorative in versions using a standard VITAPAN® classical shade guide.
3. **Isolation:** A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

Direct Restorations

4. Cavity Preparation:

- 4.1 Anterior restorations: Use conventional cavity preparations for all Class III, IV and V restorations.
- 4.2 Posterior restorations: Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should be left in the internal form of the preparation that would interfere with light transmission and therefore, the hardening of the restorative material.

¹ Trial samples are intended to educational and research purposes and are not intended for use on patients.

5. Placement of Matrix:

5.1 Anterior restorations: Mylar strips and crown forms may be used to minimize the amount of material used.
5.2 Posterior restorations: Place a thin dead-soft metal, or a pre-contoured mylar or a pre-contoured-metal matrix band and insert wedges firmly. Burnish the matrix band to establish proximal contour and contact area. Adapt the band to seal the gingival area to avoid overhangs.

Note: The matrix may be placed following the enamel etching and adhesive application steps if preferred.

6. Pulp protection: If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ ESPE™ Vitrebond™ or Vitrebond™ Plus Light Cure Glass Ionomer. Vitrebond or Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation.

7. Adhesive System: To bond Filtek™ One Bulk Fill Restorative in versions to tooth structure, use of Single Bond Universal Adhesive is recommended. Refer to adhesive system product instructions for full instructions and precautions for the products. After curing the adhesive, continue to maintain isolation from blood, saliva and other fluids and proceed immediately to placement of Filtek™ One Bulk Fill Restorative in versions.

Note: Follow the adhesive system instructions for use for recommended silane treatment during repair of ceramic restorations, followed by the adhesive application.

8. Delivery:

Dispensing the Composite:

8.1 Syringe: Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of restorative when dispensing is completed, turn the handle counterclockwise a half to stop paste flow. Immediately replace syringe cap. If not used immediately, the material should be protected from light.

8.2 Single-Dose Capsule: Insert capsule into the 3M™ ESPE™ Restorative Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

9. Placement:

9.1 Avoid intense light in the working field. Exposure of paste to intense light may cause premature polymerization.

9.2 Capsule: Start dispensing in the deepest portion of the preparation, holding the tip close to the preparation surface. Withdraw the capsule tip slowly as the cavity is filled, and avoid lifting the tip out of dispensed material while dispensing, to reduce voids. When dispensing has been completed, drag the capsule tip against the cavity wall while withdrawing from the operative field. For proximal areas, hold the tip against the matrix to aid material flow into the proximal box.

9.3 Slightly overfill the cavity to permit extension of composite beyond cavity margins. Contour and shape with appropriate composite instruments.

10. Curing: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure. Use light cure chart to determine appropriate cure times and conditions for all shades.

Indication	Increment Depth	All halogen lights (with output of 550-1000 mW/cm ²)	LED lights (with output 1000-2000 mW/cm ²)
Core Build-up and Class II Direct Restorations	5 mm	20 sec occlusal, 20 sec buccal, 20 sec lingual	10 sec occlusal, 10 sec buccal, 10 sec lingual
All indications listed (except Core Build-up and deeper Class II Direct Restorations)	4 mm	40 sec	20 sec
Anterior or shallow Class I Restorations	≤ 3mm	20 sec	10 sec

For Class II restorations, remove the matrix band prior to the buccal and lingual curing steps.
Contouring: Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with 3M™ ESPE™ Sof-Lex™ Finishing Strips.

12. Adjust Occlusion: Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

13. Finish and Polishing: Polish with the 3M™ ESPE™ Sof-Lex™ Finishing and Polishing System.

Indirect Procedure for Inlays, Onlays or Veneers

1. Dental Operatory Procedure

1.1 Shade selection: Choose the appropriate shade(s) of Filtek™ One Bulk Fill Restorative in versions prior to isolation.

1.2 Preparation: Prepare the tooth.

1.3 Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used.

2. Laboratory Procedure

2.1 Pour the impression of the preparation with die stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used.

2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base, then cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model on an adequate articulator.

2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.

2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is required.

2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, and then add another thin layer.

2.6 Add the first increment of composite to the floor of the preparation, stay short of the margins, and follow the cure recommendations described in the Direct Restoration section (Step 10).

2.7 Place and cure additional increments of composite. Allow for the last increment (incisal) to include the contact areas.

2.8 Place the die back into the articulated arch. Add the last increment of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process following the cure times in the Direct Restoration section (Step 10).

2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.

2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, and then polish as noted above in Direct Restorative steps 11-13.

3. Dental Operatory Procedure

3.1 Roughen the interior surfaces of the indirect restoration.

3.2 Clean the prosthesis in a soap solution in an ultrasonic bath and rinse thoroughly.

3.3 Cementation: Cement the prosthesis using a 3M™ ESPE™ resin cement system.

Technical characteristics of MD

Capsule: Dimensions capsule, (L×Ø) (mm×mm): $27.813 \pm 0.381 \times 7.722 \pm 0.127$; Tip length (mm): $8.130 \pm 1\%$; Tip diameter (mm): $2.760 \pm 1\%$; Weight (g): $1.0 \pm 1\%$.

Capsule cap: Dimensions, (L×Ø) (mm×mm): $8.229 \pm 0.127 \times 5.588 \pm 0.127$; Weight (g): $0.1 \pm 1\%$.

Dispenser: Volume (g): ((1 and 4) $\pm 1\%$)*; Cylinder inner diameter (mm): $8.03 \pm 1\%$; Cylinder outer diameter (mm): $11.05 \pm 1\%$; Hole diameter (mm): $4.64 \pm 1\%$; Total length (excluding cap) (mm): $141.58 \pm 1\%$; Cap length (mm): 66.904 ± 0.381 ; Shaft handle width (mm): 15.723 ± 0.254 ; The length of the inner shaft (mm): 74.676 ± 0.381 ; Cap width (mm): 11.303 ± 0.127 ; Height (mm): $26.55 \pm 1\%$; Weight (excluding material) (g): $10.7 \pm 1\%$.

* Dispenser filling volume depends on MD version

L-Pop blisters with adhesive Single Bond Universal: Dimensions capsule, (L×W×H) (mm×mm×mm): $(109.83 \times 23.91 \times 3.81) \pm 1\%$; Weight (g): $0.7 \pm 1\%$.

Values of the main parameters, characteristics (properties) of the material: Depth of cure (mm): ≥ 4.15 ; Sensitivity to ambient light (s): 60; Flexural strength (MPa): ≥ 144 ; Diametral tensile strength (MPa): $\geq 74,2$; Water sorption ($\mu\text{g}/\text{mm}^3$): ≤ 27.8 ; Solubility ($\mu\text{g}/\text{mm}^3$): ≤ 0.3 ; Radiopacity (mm): 3.02; Density (g/cm^3): 1.9.

Shelf life

Shelf life: 36 months from the manufacturing date.

Exploitation time

Not applicable. Medical device is disposable.

Use/operation conditions of MD

Operating and Storage:

1. This product is designed to be used at room temperature. If desired, the product may be warmed in a commercial warmer prior to use (no higher than 70°C / 158°F , no longer than 1 hour); for capsules only.

2. The product is best stored at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months.

Ambient temperatures routinely higher than 27°C / 80°F may reduce shelf life. See outer package for expiration date.

3. Do not expose restorative materials to intense light.

4. Do not store materials in proximity to eugenol containing products.

	Operating condition	Storage conditions
Temperature	2°C to 27°C	2°C to 27°C
Humidity	Up to 85 %	Up to 85 %
Pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa

Transportation:

This product is designed to be transported at storage conditions. During transportation the products may be affected short term temperature variations with no impact to product safety and performance.

Transportation of medical device is carried out by all type of transport.

	Transportation conditions
Temperature	2°C to 27°C
Humidity	Up to 85 %
Pressure	700 hPa to 1060 hPa

Compliance with standards applicable to this MD

ISO 7405, ISO 4049, ISO 6874, EN 980, EN 1041, EN ISO 13485, EN ISO 14971, ISO 10993-1, ISO 10993-3, ISO 10993-5, ISO 10993-6, ISO 10993-10, ISO 10993-11, ISO 10993-12, ISO 10993-18, IEC 62366-1.

Information about medicines and (or) pharmaceutical substances included in MD

Not applicable for this type of product.

Cleaning and Disinfection

The multiple-use syringe dispenser is not intended for direct patient contact. Use new, uncontaminated gloves when handling the syringe dispenser. Directions for cleaning and low level disinfection of the syringe dispenser are provided below:

Step 1 (Cleaning):

Use a CaviWipes™, or equivalent cleaning wipe, and wipe the entire surface of the device thoroughly for at least 30 seconds and until no visible soil remains on the device.

Step 2 (Disinfecting):

Use a new CaviWipes™, or equivalent alcohol-quaternary ammonium disinfectant wipe, to disinfect the entire surface of the device by keeping wet for the contact time listed on the disinfectant label.

Requirements for technical support and repair of MD

This medical device is not serviceable or repairable.

List of animal and / or human origin materials

Human blood derivatives: Medical devices does not contain or deliver human blood derivatives.

Animal-derived substances: No materials of animal origin are intentionally being added to the products.

MD disposal rules and conditions

Disposal methods: Dispose of content/container in accordance with the local/regional/national/international/regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.3684-21 to the waste class B.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Name and legal address of manufacturer MD and authorized representative in Russian Federation

- Developer:

3M ESPE Dental Products

Address: 2510 Conway Avenue St. Paul, MN 55144-1000 USA.

- Manufacturer:

3M ESPE Dental Products

Address: 2510 Conway Avenue St. Paul, MN 55144-1000 USA.

- **Authorized representative of the manufacturer in the Russian Federation:**



JSC «3M Russia»

Address: 108811, Almaty, st. Moscow, p. Moskovsky, Kievskoye sh., 22nd km, dom. 6, p. 1

Phone: +7 (495) 784 7474;

Fax: +7 (495) 784 7475.

Symbol Glossary

Reference number & Symbol / Title	Symbol	Description of symbol
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community.
Use-by date		Indicates the date after which the medical device is not to be used.
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
Use-by date		Indicates the manufacturer's catalogue number so that the medical device can be identified.
Catalogue number		Indicates the temperature limits to which the medical device can be safely exposed.
Temperature limit		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
Caution		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
CE Mark 0123		Indicates that U.S. Federal Law restricts this device to sale by or on the order of.
Rx Only		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Green Dot		Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
Do not re-use		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC).
Low-density polyethylene		

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