Handwritten signature and date: 11/17/20

The State of Maryland

Office of the Secretary of State

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille does not certify the content of the document for which it is issued.

Apostille

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1. Country: United States of America
This public document
2. has been signed by **Barbara H. Meiklejohn**
3. acting in the capacity of **Clerk of the Circuit Court for Montgomery County**
4. bears the seal/stamp of the **Circuit Court for Montgomery County**

Certified

5. at Annapolis, Maryland
6. the **26th day of August, 2020**
7. by The Secretary of State of Maryland
8. **No. 525119**
9. Seal

10. Signature

Handwritten signature of John C. Woburn

Secretary of State

Suite 103, 2141 P Street NW, Washington, DC 20037

Affidavit

State of Maryland)
) SS:

State of Maryland, Montgomery County, Sct.

In the Office of the Clerk of the Circuit Court for Montgomery County

I, Barbara H. Meiklejohn, Clerk of the Circuit Court for Montgomery County, Maryland, a court of record, hereby certify that **VERA V CHAWLA** was a commissioned/ appointed and qualified Notary Public commencing on the 9th day of August, 2016.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the Circuit Court for Montgomery County this 18th day of August, 2020.

A handwritten signature in black ink, appearing to read 'Barbara H. Meiklejohn', written over a horizontal line.

Barbara H. Meiklejohn

Clerk of the Circuit Court for Montgomery County

Suite 103, 2141 P Street NW, Washington, DC 20037

Affidavit

State of Maryland)
) SS:
County of Montgomery)

Before me, the undersigned notary public, this day, personally, appeared to me known or proven, who being duly sworn according to law, deposes the following:

I, Vera Chawla, hereby certify that the attached document is

- ☒ an original document
- ☐ a true copy of an original document

to the best of my knowledge and belief.

V. Chawla
Signed

Subscribed and sworn to before me this

12 day of August, 2020.

VERA V. CHAWLA
NOTARY PUBLIC
MONTGOMERY COUNTY
MARYLAND

Notary Public
MY COMMISSION EXPIRES JULY 29, 2020

To whom it may concern,

To be submitted to the Competent Health Authorities of the Russian Federation.

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© 2010-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Avanta, MEDRAD Avanta, FluiDots, and MEDRAD FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the MEDRAD® Avanta sterile disposable products can be found on page 2 of these Instructions for Use.

 **Bayer Medical Care Inc.**
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

 **Bayer Medical Care B.V.**
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

Importer BMC B.V.:
Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht,
The Netherlands

English	5	Français	6
Deutsch	7	Italiano	8
Español (España)	9	Русский	10
Svenska	11	日本語	12
Português (Brasil)	13	Nederlands	14
Hrvatski	15	Türkçe	16
Dansk	17	Български	18
Suomi	19	Ελληνικά	20
Magyar	21	Srpski	22
Česky	23	Norsk	24
Polski	25	Română	26
Slovenský	27	Español (LatAm)	28
한국어	29	Português (Europeu)	30
Slovenščina	31		

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86615878 Rev. A
February 20, 2019

[illegible]

Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the sterile Disposable Sets properly.

Important Safety Notice: The MEDRAD Avanta Sterile Disposable Sets are intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD Avanta Fluid Management Injection System and angiographic procedures and techniques.

Indications for Use: The MEDRAD Avanta Sterile Disposable Sets are specifically intended for use in the x-ray angiography environment. They are designed to administer intravascular radio-opaque contrast compounds, and common flushing agents at various volumes and flow rates into humans for use in diagnostic and interventional angiographic procedures performed in cardiology, radiology, and vascular surgery.

Contraindications: For use only with the MEDRAD Avanta Fluid Management Injection System. The MEDRAD Avanta Sterile Disposable Sets are not intended to be used for chemotherapy and are not intended to administer fluids other than intravascular contrast agents and common flushing solutions.

Restricted Sale: Rx Only.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector, or attempt an injection until all air has been removed from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD Fluidots indicators (where applicable) to verify that fluid is present and to reduce the chance of air embolization during an injection.

Do not attempt to aspirate fluid from the waste port stopcock when open to the PIV. Doing so may increase the risk of inducing air into the fluid path.

Do not attempt to aspirate fluid from the hemodynamic port. Doing so may increase the risk of inducing air into the fluid path.

Before connecting the Single-Patient Sterile Disposable Set (SPAT) to the patient, ensure all stopcocks and open ports are closed to air and all air has been removed from the fluid path before injection. Improper manipulation of the waste port stopcock as well as any fluid port left open to air may increase the risk of introducing air into the fluid path.

Biological contamination can result from reusing disposable items or failure to follow sterile technique during setup or use. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or operator injury may result from mishandling of the sharp spike. Use care in handling and inserting the spike into the contrast bottle and the saline bag/bottle.

Storage of filled syringes can promote bacterial growth. Syringes from Bayer are designed to be filled just prior to the procedure. Discard unused filled syringes.

Patient injury could result if syringe is not properly engaged. Improper engagement may cause air embolization or under-volume delivery.

Patient or operator injury may result from contrast media leaks or ruptures to the syringe or connector tubing. Use of greater pressure or occlusions in the fluid path may result in leaks or ruptures to the syringe or connector tubing. Ensure that the fluid path is open; do not exceed appropriate pressure.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Inadvertent aspiration can occur if the Multi-Patient disposable contrast tubing set is kinked or blocked during a syringe fill process. Verify multi-patient contrast tubing is free of kinks before attempting to fill the syringe.

Patient injury can result if the patient is connected during any contrast or saline purge. Ensure the patient is disconnected before any purge.

Blood can contaminate the multi-patient portion of the disposable set. Ensure that blood does not appear past the pressure isolation valve in the single patient disposable set.

Ensure vent cap on spike(s) are open if using bottles. If closed, a vacuum may be created and air may be introduced.

Do not exceed the five patients per day usage requirement for the Multi-Patient Disposable Set (MPAT). This is the usage approved by appropriate regulatory bodies.

Use of non-Bayer supplied disposables, including administration sets and additions to administration sets, such as but not limited to bleed back control devices and pressure transducers may cause patient injury if not properly connected or flushed. These devices must be compatible with your system. Refer to manufacturer's instructions for proper use of these devices.

Use of extension tubing connected to the Single Patient Disposable Set may cause patient injury and will degrade system performance. Do not connect extension tubing to the Single Patient Disposable Set.

Cautions

Component damage may occur if disposables are not installed properly. Assure all connections are secure. This will help minimize leaks, disconnections, air introduction, and component damage. Do not use tools to over tighten connections or to assist in removal of disposables.


Do not use methyl alcohol on disposable sets. Use of methyl alcohol will degrade the integrity of the disposable sets.

All contrast and saline media should be used in accordance with the manufacturer's indication for use.

NOTE: To perform the following procedures, the operator must select the Case Setup screens on the DCU. Refer to the injector operation manual for additional information.

Multi-Patient Sterile Disposable Set (MPAT) Installation (usable for up to five patients). The Multi-Patient Disposable Set consists of the following items. (Also see MPAT Setup illustration):

Ventrad Contrast Spike	2nd HPCT w/ male luer
Contrast Drip Chamber	Saline Spike
Low Pressure Contrast Tube (LPCT)	Saline Drip Chamber
High Pressure Contrast Tube (HPCT)	Low Pressure Saline Tube (LPST)
Contrast Flow Control Valve (CFCV)	Post Drip Chamber Saline Tube
150 ml Multi-Patient Syringe	

General Setup Instructions

Be sure that the patient is disconnected and the injector is on.

1. Install syringe. (See Install Multi-Patient Syringe below).
2. Install saline tubing. (See Install Saline Tubing below).
3. Install contrast tubing. (See Install Contrast Tubing below).
4. Rotate head to upright position, fill syringe and remove dust caps.
5. Purge multi-patient disposables. Catch the discharged fluid in a suitable container and tap disposable connection points and pressure jacket to remove air.
6. Rotate head downward.
7. Attach MEDRAD Avanta Sterile Hand Controller, and maintain sterility, if applicable.

NOTE: Refer to Operation Manual for MEDRAD Avanta Fluid Management Injection System for further information.

After step 7 in the General Setup Instructions is complete, go to SPAT Installation.

Install Multi-Patient Syringe

NOTE: Ensure patient is not connected.

Select Multi-Patient Syringe from Case Setup.

1. Retract piston.
2. Insert syringe into pressure jacket and raise drop front arms.
3. Engage the piston to the plunger.

Install Saline Tubing

NOTE: Ensure patient is not connected.

Select Multi-Patient Disposable from Case Setup, then select Install Saline Tubing.

1. Install drip chamber.
2. Route tubing through pump and close door.
3. Route tubing around peg.
4. Insert tubing into gross air detector and close door.
5. Spike saline bag.
6. Fill saline drip chamber completely and snap into position before use.

NOTE: To differentiate the saline and contrast lines of the MPAT, the contrast line will be identified as the line connected to the stopcock that is placed in the CFCV.

Install Contrast Tubing

Select Install Contrast Tubing from Case Setup.

1. Install drip chamber.
2. Connect tubing to syringe.
3. Install stopcock.
4. Insert tubing into gross air detector and close door.
5. Spike contrast bottle.
6. Fill contrast drip chamber completely and snap into position before use.

Go to step 4 in General Setup Instructions.

Single-Patient Sterile Disposable Set (SPAT) Installation (must be replaced for each patient use). The Single-Patient Disposable Set consists of the following items (Also see SPAT Setup illustration):

Rotating luer (A)	High Pressure Saline Tube w/ male luer connector (D)
Stopcock for waste port (B)	High Pressure Contrast Tube w/ female luer connector (E)
Pressure Isolation Valve (PIV) Assembly / Hemodynamic pressure monitoring port (C)	

Single-Patient Sterile Disposable Set (AVA 500 SPAT Angio) is indicated for use when a user attachable transducer is unnecessary and does not contain item (C) pictured below.

Select Single-Patient Disposable from Case Setup.

Be sure that the patient is disconnected.

1. Remove all dustcaps.
2. Connect the SPAT saline tube to MPAT saline tube.
3. Connect the SPAT contrast tube to MPAT contrast tube.

To purge air out of SPAT:

NOTE: Ensure patient is not connected.

When using AVA 500 SPAT L continue with step 4. When using AVA 500 SPAT ANGIO proceed to step 5.

4. Attach transducer to hemodynamic port, verify the transducer stopcock is open, and the waste port stopcock is off to waste port.
 - a. Attach 20cc saline syringe to the transducer port.
 - b. Hold SPAT parallel to the patient with the PIV and 20cc syringe connection facing upward at a 45° angle.
 - c. Forcefully flush through the transducer and PIV assembly until all air is expelled from both components. Flick or tap PIV to facilitate air removal (Repeat as necessary).
 - d. Close transducer port stopcock. Remove 20cc syringe.
5. Initiate contrast and saline purge, tap SPAT to MPAT saline connection point and wasteport stopcock to remove air while purging is in progress.
6. Repeat steps above as needed, until all air is purged.
7. Carefully inspect the entire SPAT and MPAT to verify all air is removed before connecting to patient.

NOTE: Disconnect the SPAT from the MPAT immediately following each procedure.

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
OFFICE OF THE SECRETARY OF STATE


NOTICE

This notice is given by the Secretary of State acting pursuant to the Proclamation of Governor Lawrence J. Hogan, Jr. (the "Governor") dated March 5, 2020 proclaiming a state of emergency and a catastrophic health emergency related to COVID-19, as it may be amended or renewed from time to time, and the Order of the Governor dated March 12, 2020 entitled "Extending Certain Licenses, Permits, Registrations, and Other Governmental Authorizations, and Authorizing Suspension of Legal Time Requirements", as it may be amended from time to time, after finding that the action(s) described herein will not endanger the public health, welfare, or safety.

The four year statutory term of a notary public set forth in Section 18-103(c)(1) of the State Government Article, Maryland Annotated Code shall be extended for all notaries public with an active notary commission. Furthermore, the notice and timing requirements set forth in Section 18-103(d), (e)(6) of the State Government Article, Annotated Code of Maryland are suspended. Any effect of Title 18 of the State Government Article, Annotated Code of Maryland inconsistent with the foregoing is hereby suspended.

Effective immediately and signed on this 20th day of March, 2020, this Notice shall remain in effect until 30 days after the state of emergency has been terminated and the proclamation of the catastrophic health emergency has been rescinded.

John C. Wobensmith
Secretary of State

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Certified

5. at Annapolis, Maryland
6. the **26th day of August, 2020**
7. by The Secretary of State of Maryland
8. **No. 525127**
9. Seal

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A handwritten signature in black ink, reading "Barbara H. Meiklejohn", written over a horizontal line.

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

12 day of August, 2020.

To whom it may concern,

To be submitted to the Competent Health Authorities in Federal Republic of Russia

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola, Pennsylvania 15051-0780, USA) hereby state that attached document (the addendum to the instruction for use) for medical device “Шприц одноразовый для нескольких пациентов к устройствам для внутривенного введения MEDRAD Avanta” is true.

////////////////

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1 Bayer Drive
Indianola (PA), United States
15051

Best regards,

www.bayer.com

Emily Bartolomeo
Manager, Global Regulatory Strategy Devices
Global Regulatory Affairs, Radiology
Bayer Medical Care, Inc.

Medical device name:

Single-use syringe for multiple patients to intravenous devices MEDRAD Avanta

Medical device purpose

MD is designed for intravascular administration of X-ray contrast agents and normal flushing solutions in various volumes and rates into human body during diagnostic and interventional angiographic procedures used in cardiology, radiology and vascular surgery.

Possible undesirable effects of the medical device use

No possible undesirable effects identified by manufacturer. All possible warnings that are necessary for the safe use of the product are indicated in the instructions for use.

Technical characteristics

Parameter	Value
Syringe capacity	Not more than 150 mL
Syringe barrel overall length	218.8 mm
Syringe barrel outer diameter	45.9 mm
Maximum pressure	8273 kPa

Shelf life

Shelf life – 2 years from sterilization if storage and transportation requirements are met.

Transportation and storage conditions

Medical device can be transported by any type of enclosed transport subject to transport and storage requirements below:

Temperature	-20°C to +60°C
Relative humidity	5 % to 100 %
Ambient pressure	57 kPa to 106 kPa

Operation conditions

Temperature	+15°C to +30°C
Relative humidity	20 % to 75 %
Ambient pressure	70 kPa to 106 kPa

* Efficient and safe operation of the device is only possible in the specified environmental conditions limits (temperature, humidity, ambient pressure).

Warranty liabilities

The warranty lifetime is not provided for due to the specifics of the medical device.

The manufacturer warrants that the products are free from defects.

Manufacturer bear no responsibility nor liability in case of misuse or non-adherence to the Instructions For Use.

Reclamation

Issues related to the medical device quality must be referred to the manufacturer's representative in the Russian Federation at the address below:

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