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ATTESTATION

STATE OF FLORIDA
COUNTY OF DUVAL

On this 27th day of August 2020, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following Instructions for Use (IFU) of HydroCleanse™, issued by Medtronic Xomed, Inc., located at 6743 Southpoint Drive N, in Jacksonville, Florida USA 32216, and presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.


Alexandra Oliver

Alexandra Oliver
Regulatory Affairs
Medtronic Xomed, Inc.

Sworn to (or affirmed) and subscribed before me by means of ☒ physical presence OR ☐ online
notarization on this 27th day of AUGUST, 2020, by

ALEXANDRA OLIVER

Judith L. Fisher
(Official Notary Public signature)

JUDITH L. FISHER

Judith Fisher


Personally known ☒ OR produced identification ☐

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








2

Indicated use

HydroCleanse™ sinus wash delivery system is intended to irrigate and suction within the paranasal sinuses.

Indications for use

HydroCleanse™ sinus wash delivery system is indicated for use whenever irrigation of the paranasal sinuses is desired.

Device description

HydroCleanse™ sinus wash delivery system is a manual device that enables a user to simultaneously and continuously deliver pressurized saline irrigation to the paranasal sinuses and suction out excess fluid. The system consists of three major components: suction catheter tip, irrigator handle system, and syringe system.

Suction catheter tips include three different angles: 13°, 70°, 120°. The irrigator handle system includes an irrigation catheter with a fan-shaped spray nozzle and a rotating dial wheel for larger spray coverage, a spin lock for suction catheter tip attachment, and a suction connector that connects to a vacuum pump to provide suction. The syringe system includes a syringe with self-returning handgrip, a 3-way luer connector with check valves to prevent fluid backflow, and a spike for direct connection into the saline bag. Spray action is initiated by depressing the handgrip and the syringe refills when the handgrip self-returns.

Contraindications

None known.

Warnings

Do not irrigate if the sinus cavity is compromised either intraoperatively or prior to surgery (for example, dehiscence, erosion, or fistula), or an anatomical anomaly is known or suspected.

Do not use if the sterile package has been opened, or is damaged.

Precautions

Do not attempt to re-sterilize or reuse single use disposable device. Re-sterilization or reuse may result in device damage and patient injury.

Dispose of the used HydroCleanse™ sinus wash delivery system following healthcare facility guidelines on proper disposal of contaminated materials.

It is important to monitor the orbit for pressure changes and/or swelling during irrigation and to monitor for changes to vital signs that could be suggestive of increased ocular or cerebral pressure.

Verify the suction line is working properly to minimize ingestion of irrigation fluid during the procedure. The device may be configured to a suction only instrument when a higher suction rate is needed.

Setup

Remove the HydroCleanse™ sinus wash delivery system from its package.

Attach the irrigation tubing assembly to the syringe through the 3-way luer connector (figure 2).

Attach the light blue luer connector located at the open end of the irrigator tubing onto the light blue dial wheel located on the irrigator handle system (figure 3).

Insert the spike connector on the solution tubing into the saline bag (approximately 100 cm above the patient) and close the pinch clamp.

Grip the handgrip on the syringe system with one hand and pull off the clear wrap on the handle to open it.

Open the pinch clamp and prime the handgrip until the fluid path is filled with saline irrigation solution (figure 4).

Note: You must prime the hand piece until the fluid path is filled with saline irrigation solution before any saline will exit the spray nozzle.

To use suction, insert the suction connector on the irrigator handle system into the suction line (the vacuum range is -0.3 bar to -0.7 bar) (figure 5).

Use

Place the device into the desired location under endoscopic visualization.

Irrigate within the sinus spaces and suction out excess fluid.

Rotate the spray nozzle

You may use your index finger to turn the dial wheel clockwise or counterclockwise to rotate the spray nozzle. Use a combination of clockwise and counterclockwise motions to avoid twisting the irrigator tubing (figure 6).

Change the suction catheter tip

You will need to switch the suction catheter tip depending on which sinuses you are irrigating. The HydroCleanse™ sinus wash delivery system comes with the 13° suction catheter tip attached. Select the properly degreed suction catheter tip (13°, 70°, 120°) for the type of sinus irrigation you intend to perform.

- Frontal sinuses: When accessing the frontal ostia/recesses, use the 70° suction catheter tip.
- Sphenoid sinuses: When accessing the sphenoid sinus ostia, use the 13° suction catheter tip.
- Maxillary sinuses: When accessing the maxillary ostia/ethmoid infundibula, use the 120° suction catheter tip.

Use the following procedure to switch between suction catheter tips:

1. Grip the dark blue front spin lock in one hand and grip the wings at the base of the suction catheter tip in the other hand (figure 7).
2. Turn the front spin lock counterclockwise until loose.
3. Pull the suction catheter tip away from the spin lock to remove it from the irrigation catheter (figure 8).
4. Select the desired suction catheter tip.
5. Slide the suction catheter tip over the irrigation catheter at a desired irrigation angle until it is seated against the inside of the front spin lock.
6. While holding the wings of the suction catheter tip in one hand, turn the spin lock in a clockwise direction until it is firmly secure.

Change the direction of the suction catheter tip

While holding the wings of the suction catheter tip in one hand, turn the dark blue front spin lock in a counterclockwise direction until the suction catheter tip can twist between your fingers.

Twist the suction catheter tip until the curve is facing the desired direction.

While holding the wings of the suction catheter tip in one hand, turn the spin lock in a clockwise direction until it is firmly secure.

Change to suction only instrument

Grip the irrigator handle system in one hand and the rear light blue spin lock near the dial wheel in the other.

Turn the rear spin lock in a counterclockwise direction until the irrigation catheter is free.

Pull back on the dial wheel until the irrigation catheter is free from the irrigator handle system (figure 9).

Attach the white cap (attached at the base of the irrigator handle system) to the rear of the irrigation port (figure 10).

Sterility warning

The HydroCleanse™ sinus wash delivery system is provided STERILE and is intended for SINGLE PATIENT USE ONLY. DO NOT RESTERILIZE THE HydroCleanse™ sinus wash delivery system. MEDTRONIC ASSUMES NO LIABILITY FOR DEVICES THAT HAVE BEEN RESTERILIZED BY HEALTH CARE FACILITIES.

Storage

Store the HydroCleanse™ sinus wash delivery system at room temperature.

Troubleshooting

If the HydroCleanse™ sinus wash delivery system does not function properly, contact Medtronic Xomed, Inc. customer service.

Customer service information

For further information regarding the use of this product or to report any problems, please contact Medtronic Xomed, Inc., or contact your local distributor.

Limited warranty

- A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.
- B. To qualify for this LIMITED WARRANTY, the following conditions must be met:
 1. The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 2. The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 3. Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.
 4. The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.
 5. Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.
- C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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ATTESTATION

STATE OF FLORIDA
COUNTY OF DUVAL

On this 27th day of August 2020, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following Appendix to HydroCleanse™, issued by Medtronic Xomed, Inc., located at 6743 Southpoint Drive North, and presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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120° Suction catheter tip

- Hardness: 70 ± 5 SHORE A.
- Stiffness: deflection shall be 0.6 ± 0.02 mm under 1N at middle point between any simply support 70mm length.
- Detach force >70 N/15s (static).
- Leakage pressure >200 kpa (2 bar)/15S.
- Catheter tip angle is 120°
- Overall catheter length is $104,9 \pm 0,5$ mm.
- Tube length is 89.8 ± 0.5 mm.
- Tube diameter OD/ID is 3.20 mm / 2.60 mm.
- Luer length is 20 mm.
- Luer width is 14.15 mm.
- Luer inner diameters are 3.45 ± 0.10 mm/ 4.4 mm (distal/proximal)

70° Suction catheter tip

- Hardness: 70 ± 5 SHORE A.
- Stiffness: deflection shall be 1.3 ± 0.2 mm under 1N at middle point between any simply support 70 mm length.
- Detach force >70 N/15s (static).
- Leakage pressure >200 kpa (2 bar)/15S.
- Catheter tip angle is 70°
- Overall catheter length is 111.5 ± 0.5 mm.
- Tube length is 99 mm.
- Tube diameter OD/ID is 3.20 mm / 2.60 mm.
- Luer length is 20 mm.
- Luer width is 14.15 mm.
- Luer inner diameters are 3.45 ± 0.10 mm/ 4.4 mm (distal/proximal)

13° Suction catheter tip

- Hardness: 70 ± 5 SHORE A.
- Stiffness: deflection shall be 1.3 ± 0.2 mm under 1N at middle point between any simply support 70mm length.
- Detach force >70 N/15s (static).

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- Leakage pressure >200kpa (2 bar)/15S.
- Catheter tip angle is 13°
- Overall catheter length is 119.9 ± 0.5 mm.
- Tube length is 108.5 ± 0.5 mm.
- Tube diameter OD/ID is 3.20 mm / 2.60 mm.
- Luer length is 20 mm.
- Luer width is 14.15 mm.
- Luer inner diameters are 3.45 ± 0.10 mm/ 4.4 mm (distal/proximal)

Irrigation catheter

- Leakage pressure in air 400kpa(4bar)/15s.
- Irrigation catheter length is 153.6 ± 0.5 .
- Catheter tube OD/ID is 1.95 ± 0.05 mm/ 1.55 ± 0.05 mm.
- Tube end radius is 0.77.
- Nozzle area end width Outer/Inner is 0.42 mm/ 0.33 ± 0.05 mm.

Syringe

- Detach force >200N.
- 10 ml syringe.
- Leakage (Vacuum: 88kpa/0.88bar, 60s).
- Length of the syringe is 118.90 mm.
- Diameter of the syringe is 16.30 mm.
- Flange width is 31.80 mm
- Barrel flange diameter is 18.90 mm

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Information in Russian sticker

The final packaging for the Russian market is accompanied with additional sticker.

The following information should be printed on the Russian-language sticker (word order and text font may vary):

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable).

Sample of the sticker:

	Name of medical device	
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer's authorized representative:: Regestretion Certificate № _____ от _____	
	Date of manuf../Lot/UBD: See on pack	See IFU See packaging symbols

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© 2018 Medtronic, Inc.