Instructions for use of a medical device

Name of a medical device: OKVision BioTwin universal solution for contact lens care (see Attachment)

Manufacturer: Oy Finnsusp AB, Paaskykalliontie 5, FI-21420 Lieto, Finland

Composition: Disodium EDTA, Twin-EdaXyl, hyaluronic acid, Polyhexanide 0.0001% in sterile isotonic buffered solution.

Material of the component	Grade or specification	Percentage, %
Polyhexanide	Azelis	0,0001
Hyaluronic acid	Contipro	0,01
Twin-Edaxyl : Xylitol EDG (ethyl-6-O-dekanoyyli glucoside)	Sigma Aldrich Chiral Vision	0,0125
Disodium EDTA	VWR, Sigma Aldrich	0,09
Sodium Borate	VWR	0,1
NaCl	VWR, Sigma Aldrich, AkzoNobel	0,6
Boric acid	VWR	0,6
Purified water	-	>98

	Material of the component	Grade or specification	Percentage, %
Container	PEHD (high-density polyethylene)	LyondellBasell	>98
	Antimicrobial agent containing silver: Alphasen MB SLA0050122	Clariant	
Caps	PEHD (high-density polyethylene)	LyondellBasell	>98
Right cap and container	White pigment COLBN008	Clariant	Less than 100 ppm
Left cap	Blue pigment COLBL010	Clariant	Less than 100 ppm

Packing:

Bottle	HDPE	Promens	100%
Tip	LDPE	Coreplast	100%
Protective ring	PP	Muovitalli Lehtonen	100%
Cap	PP	Muovitalli Lehtonen	100%

BioTwin contains 0.01% hyaluronic acid, which is classified as a medicinal product in European Union. The molecular weight of the hyaluronic acid used is approximately 0.85 MDa and corresponds to the class of medical devices.

The product is intended for use as a care agent for contact lenses, for cleaning, disinfecting, storing, rinsing and moisturizing. The intended action of the hyaluronic acid in the product is to bind water on contact lenses and thus moisturize them and make contact lens wear more comfortable. The intended action of hyaluronic acid does not involve pharmacological, immunological, or metabolic mechanisms.

Field of application and purpose of a medical device: Field of application – Ophthalmology. Designed for cleaning, disinfecting, storing, rinsing and moisturizing of all types of contact lenses.

Properties and efficiency: Twin-EdaXyl, a double effect innovative active complex with lubricates the surface of contact lenses and removes deposits during daily care. Biotechnologically produced hyaluronic acid moisturizes contact lenses and eyes at the same time. Made in Finland and developed in collaboration with the Finnish Allergy, Skin and Asthma Federation

Indications for use: Designed for cleaning, disinfecting, storing, rinsing and moisturizing of all types of contact lenses.

Contraindications for use: Do not use the solution if you have a hypersensitivity to any of its components. Do not use container if you are sensitive to silver.

Instructions for use:

Clean, rinse and disinfect contact lenses daily to prevent deposition of foreign substance and microorganisms on their surface.

1. Thoroughly wash, rinse and dry your hands.
2. Remove one lens and place it in the palm of your hand. Apply a few drops and gently rub the lens with your other hand for 10-15 seconds.
3. Thoroughly rinse the lens with the solution and place it in the container. Repeat the same procedure for the other lens.
4. Fill the lens container with solution and seal tightly for disinfection for at least 4 hours or overnight. Always rinse your contact lenses with fresh solution before using them. After inserting the contact lens into the eye, pour out the used solution, rinse the container with fresh solution and let it dry.

Precautions:

Important!

Contact lenses can only be used safely after consultation with a contact vision specialist or ophthalmologist. Strictly follow these instructions and your healthcare professional's advice to ensure the safe and appropriate use of your contact lenses.

- If you experience irritation or discomfort in your eyes, remove your contact lenses immediately and seek the help of an ophthalmologist.
- Do not use the solution if the integrity of the control ring on the vial is broken.
- Do not touch the tip of the vial or pour the solution into other containers to avoid compromising sterility.
- Never reuse solution
- After use, always close the bottle carefully and store it at room temperature (15-30 °C).
- Never mix solution with other contact lens care products. It is not recommended to change the type of solution during the period of use of contact lenses of planned replacement.
- Never use tap water to rinse contact lenses and containers.
- If the lenses have not been used for more than two weeks, disinfect them again before use.
- Make sure the container is clean and change it regularly. Rinse container with fresh solution and allow to dry open.
- Keep solution out of the reach of children.

Possible side effects:

Allergic reactions are possible. If allergic reaction occurs, stop use and consult a doctor.

Interaction with other drugs:

The use of any medicinal products for the eyes simultaneously with the solution is allowed only under the supervision of a doctor.

Product form: 120 ml and 360 ml plastic bottles with a tamper evident cap. The antibacterial container only comes with a 360 ml bottle. The 120 ml solution is supplied without a lens container. The container material contains silver ions, which have antimicrobial properties. For ease of use of lenses, container lids have different markings and colors. The lid marked L is dark blue and the lid marked R is white.

Labeling:

Oy Finnsusp AB, Paaskykalliontie 5, FI-21420 Lieto, Finland

MD

0537

STERILE

A

The main technical and functional characteristics of the medical device:

The solution is produced in compliance with international and national standards. In terms of biological safety, the product meets the requirements of ISO 10993-1. The medical device is non-toxic. The solution is sterile. The solution is dispensed under aseptic conditions. The osmolality of the solution should amount to 270-360 mOsm/ kg, pH - 7.1-7.7, appearance- clean, colorless, homogeneous, without impurities, the concentration of hyaluronic acid should amount to 0.006-0.011%.

Parameter	Specification for the finished product
pH	7.1-7.7
Osmolality	270-360 mOsm/kg
Visual	Clear, colorless, homogenous, no particles
Sterility	No growth
Bottle size	120 ml 360 ml
HA concentration	0.006-0.011 %

Requirements for the use and operation of a medical device:

It is recommended for a private use. This medical device is intended for use at home, as well as for use in medical institutions.

Shelf life: 3 years from the manufacturing date. The solution should be used within 4 months from the date of opening the bottle. Do not use after the expiration date printed on the bottle. Do not store the solution in the refrigerator or out in the sun.

Conditions of transportation and storage: 15 – 30°C. Keep away from direct sunlight.

Warranty: This product is designed, manufactured, tested and packaged in accordance with all applicable requirements. The manufacturer guarantees high-quality manufacturing of products in accordance with the received documentation and instructions.

Disposal Procedures: Dispose of products properly according to established procedures. Upon completion of use, this product does not preset a potential biohazard. Any such devices should be handled and disposed of in accordance with good medical practice and applicable local, state and national laws and regulations. Solution: Dispose of by flushing with plenty of water. Plastic packaging: dispose of by incineration or recycling, in accordance with applicable regulations.

Authorized representative:

For all questions, including complaints about products, repairs and maintenance, contact the Authorized Representative in the Russian Federation:

OKVision Retail Ltd

63B, bld. 4, Mikhalkovskaya Str., Moscow, 125438

Tel: 8 (495) 602-05-51

E-mail: office@okvision.ru

Website: www.okvision.ru


Latest revision: 02/2021

OKVision BioTwin universal solution for contact lens care, product forms:

- 1) OKVision BioTwin universal solution for contact lens care in a 120 ml bottle;
- 2) OKVision BioTwin universal solution for contact lens care in a 360 ml bottle + antimicrobial lens container.

I hereby confirm that the information
above is accurate.


Lined 1st April, 2021


RIIKKA JÄRVINEN
Oy Finnsusp AB

This is to certify the own hand signature of
RIIKKA JÄRVINEN

and that ~~he~~ she is ~~legally~~ legally authorized to sign
for the company **Oy Finnsusp AB**
Tampere, 1st of April 2021

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
This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.
This Apostille does not certify the content of the document for which it was issued.

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