sebia

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
2018/12

NAME OF THE MEDICAL PRODUCT

Reagent kit for qualitative and quantitative determination of protein fractions by capillary electrophoresis for in vitro diagnostics (CAPI 3 PROTEIN(E) 6)

Note: Further on, the text may contain "reagent kit", "product", "CAPI 3 PROTEIN(E) 6 kit".

INTENDED USE

For *In Vitro* Diagnostic Use. It is designed for the separation and quantification of human serum protein fractions by capillary electrophoresis. The kit of reagents is used in conjunction with the CAPILLARYS 3 automatic capillary electrophoresis system. It is an auxiliary tool in diagnostics.

FUNCTIONAL PURPOSE.

An auxiliary tool in diagnostics.

Determination of quantitative and qualitative changes in the main protein fractions of blood serum is used to diagnose and control the treatment of patients at risk of protein abnormalities.

SAMPLES FOR ANALYSIS

Human blood serum is used for the analysis.

DESCRIPTION OF THE TARGET ANALYTE:

Normal serum proteins separate into six major fractions (gamma-globulins, beta-2-globulins, beta-1-globulins, alpha-2-globulins, alpha-1-globulins, and albumin) to determine quantitative and qualitative changes.

FIELD OF APPLICATION. INFORMATION ABOUT POTENTIAL CONSUMERS

For *In Vitro* Diagnostic Use. For professional use in clinical-diagnostic laboratories, by qualified personnel in the field of clinical laboratory diagnostics.

Note: Use of the reagents kit was shown in conjunction with the following instruments (including, but not limited to): CAPILLARYS 3 OCTA, CAPILLARYS 3 TERA systems. Please refer to instruction manuals or contact manufacturer's authorized representatives to clarify the possibility of using the reagents kit with other SEBIA capillary electrophoresis systems.

NOTE: In this instruction sheet, the name "CAPILLARYS 3" is used for the SEBIA CAPILLARYS 3 OCTA and CAPILLARYS 3 TERA automated instruments.

PRINCIPLE OF THE TEST

Protein electrophoresis is a well established technique routinely used in clinical laboratories for screening samples for protein abnormalities. Besides the electrophoresis techniques performed on different media, including agarose gel, the capillary electrophoresis has been developed to provide a complete automation of this testing with fast separation and good resolution. It is defined as a technique of electrokinetic separation carried out in a tube of internal diameter lower than 100 μm filled with a buffer composed of electrolytes. In many aspects, the methodology can be considered as an intermediary type of technique between classical zone electrophoresis and liquid chromatography.

The CAPILLARYS 3 instrument uses the principle of capillary electrophoresis in free solution that is the most common form of capillary electrophoresis. With this technique, charged molecules are separated by their electrophoretic mobility in an alkaline buffer with a specific pH. Separation also occurs according to the electrolyte pH and electroosmotic flow.

The CAPILLARYS 3 instrument performs all sequences automatically to obtain a protein profile for qualitative or quantitative analysis.

The CAPILLARYS 3 OCTA instrument has 8 capillaries functioning in parallel allowing 8 simultaneous analyses. The CAPILLARYS 3 TERA instrument have 12 capillaries functioning in parallel allowing 12 simultaneous analyses. A sample dilution with buffer is prepared and injected by aspiration at the anodic end of the capillary. A high voltage protein separation is then performed and direct detection of the proteins is made at 200 nm at the cathodic end of the capillary. The capillaries are immediately washed with a wash solution and prepared for the next analysis with buffer.

By using alkaline pH buffer, proteins are detected in the following order: gamma globulins, beta-2 globulins, beta-1 globulins, alpha-2 globulins, alpha-1 globulins and albumin with each zone containing one or more proteins.

The electrophoregrams can be interpreted visually to screen for any pattern abnormalities. Direct detection provides accurate relative quantification of individual protein fractions.

REAGENTS AND MATERIALS SUPPLIED IN THE CAPI 3 PROTEIN(E) 6 KIT

ITEMS	PN 2503
CAPILLARYS PROTEIN(E) 6 buffer (ready to use)	3 vials, 700 mL each
Filter CAPILLARYS	4 filters
Package insert	1

The shelf life of the reagent kit is determined by the reagent component of the set with the nearest expiration date. The maximum shelf life of the kit is 36 months from the date of manufacture.

FOR OPTIMAL MANAGEMENT OF TRACEABILITY : All reagents from the same kit must be

used together.

TO OBTAIN THE EXPECTED PERFORMANCES : The package insert instructions must be observed.

WARNING : Do not use marketed deionized water, such as water for ironing for example (risk of important capillaries damage). Use only water with ultrapure quality, such as injection grade water.

1. CAPILLARYS PROTEIN(E) 6 BUFFER

The vial with the buffer is marked with an RFID (radio frequency) tag for recognition by the instrument and a label with white inscriptions on a pink background.

Vial dimensions (LxWxH, mm): (125x65x145) $\pm 5\%$

Cap diameter, mm: 50 $\pm 5\%$

Cap height, mm: 19 $\pm 5\%$

Reagent volume, ml: 700 $\pm 5\%$

Maximum capacity of the vial (ml): 750 ml $\pm 5\%$

Weight of vial with reagent, g: 820 $\pm 5\%$

Empty vial weight, g: 60 $\pm 5\%$

Reagent pH: 9.9 ± 0.5

Composition

Water, boric acid, sodium hydroxide, anhydrous sodium sulfate, Proclin300.

Preparation

The buffer is ready to use. It contains : buffer solution pH 9.9 ± 0.5 .

Storage, stability and signs of deterioration

Store the buffer at room temperature (15 to 30 °C) or refrigerated (between 2 and 8 °C).

It is stable until the expiration date indicated on the kit package or buffer vial labels.

Avoid storage close to a window or to a heat source. Do not use the buffer after the expiration date.

NOTE : When analysis buffer is stored between 2 and 8 °C, it is recommended to allow reagent to come to room temperature prior to use.

DO NOT FREEZE.

Once the buffer vial has been opened and positioned on the CAPILLARYS 3 instrument, it is stable for a maximum of 2 months (accumulated). If the buffer vial is planned to be

used for more than 2 months, it must be removed from the instrument after each use and stored at room temperature (15 to 30 °C) or refrigerated (between 2 and 8 °C), it is then stable until the expiration date indicated on the buffer vial label.

Discard buffer if it changes its appearance, e.g., becomes cloudy due to microbial contamination.

2. FILTER CAPILLARYS

Use

Disposable filters for filtration of analysis buffer and distilled or deionized water (used for capillaries rinsing).

Filter dimensions (LxWxH, mm): (50x50x25) $\pm$ 5%

Packaged filter weight, g: 3,2 $\pm$ 5%

Filter diameter, mm: 26 $\pm$ 5%

Pore diameter: 5 μ m

The filter is made of polypropylene.

IMPORTANT : When kit replacement, change systematically all the filters. Wear clean gloves for handling and installation of filters.

Screw one filter at the connector situated at the extremity of each tube plunging in vials of buffer and distilled or deionized water. When setting filters on the instrument, rinse the connectors and the tubes with distilled or deionized water.

Storage

Before use, store the filters in their sealed package in a dry place at room temperature (15 to 30 °C) or refrigerated (between 2 and 8 °C).

The procedure and conditions for disposal

Dispose of all kit materials in a safe and acceptable manner that complies with SanPiN 2.1.7.2790-10 and SP 1.3.2322-08. The hazard class of waste when used is B (epidemiologically dangerous) and is subject to disposal (collection of item 4. SanPiN 2.1.7.2790-10), disinfection/neutralization (item 5. SanPiN 2.1.7.2790-10), temporary storage (item 6. SanPiN 2.1.7.2790-10), transportation (item 7. SanPiN 2.1.7.2790-10), accounting and control (item 8. SanPiN 2.1.7.2790-10).

Cardboard boxes and liners, as well as the entire set, after the expiration date or if the packaging is damaged, are disposed of as medical waste of class A (epidemic safe) or as waste of class 5 (practically non-dangerous) in accordance with the criteria of environmental regulation.

Waste generated as a result of operation

Equipment waste presents a mixed hazard, both chemical and biological. Biological risk can be reduced using Sodium hypochlorite solution.

Concentration and contact time:

Final chlorine %	Dilution	Action	Contact time
0,5 %	Dilution 1/5 from a 2.6 % solution	Inactivation of Conventional Transmissible Agents	For at least 15 minutes
	or Dilution 1/20 from a 9.6 % solution	Inactivation of Unconventional Transmissible Agents Group II	For at least 15 minutes
2 %	Dilution 4/5 from a 2.6 % solution or Dilution 1/5 from a 9.6 % solution	Inactivation of Unconventional Transmissible Agents Group III	For at least 1 hour

In practice:

Let a contact time of one night so that the chlorine ions are eliminated on contact with treated medium proteins. After removal of biological risk and without risk of chemical, effluents can be discharged into the sink.

In cases of chemical residue, the user must arrange for special steps to eliminate such waste.

Wastes may be discharged into sewers only if:

- beforehand they are decontaminated with Sodium hypochlorite solution (then chemical hazard only),
- The chlorine residual is below the thresholds of local rejection.

If these thresholds are exceeded, decontaminated waste must be kept in closed and resistant plastic containers hermetically closed and specifically identified, and then properly disposed by a specialized agency.

REAGENTS REQUIRED BUT NOT SUPPLIED

WARNING : See the safety data sheets.

1. DISTILLED OR DEIONIZED WATER

Use

For capillaries rinsing in automated instrument for capillary electrophoresis CAPILLARYS 3, SEBIA.

It is recommended to use filtered distilled or deionized water (on a filter with a porosity $\leq 0.45 \mu\text{m}$) and with a conductivity lower than $3 \mu\text{S/cm}$, which corresponds to a resistivity higher than $0.33 \text{ M}\Omega\cdot\text{cm}$.

To prevent microbial proliferation, change the water every day.

IMPORTANT : Before filling the rinse container, it is recommended to wash it with plenty of distilled or deionized water.

2. CAPICLEAN (registration certificate in Russia PY № ФС3 2009/05453 from 25.11.2009)

Composition

The vial of CAPICLEAN concentrated solution contains: proteolytic enzymes, surfactants and additives nonhazardous at concentrations used, necessary for optimum performances.

Use

For sample probe cleaning in automated instrument for capillary electrophoresis CAPILLARYS 3, SEBIA, during the CAPICLEAN cleaning sequence.

IMPORTANT :

When less than 500 samples are analyzed within a week, launch a CAPICLEAN cleaning sequence at least once a week.

When less than 500 samples are analyzed within a day but more than 500 samples are analyzed within a week, launch a CAPICLEAN cleaning sequence after every 500 analyses.

When more than 500 samples are analyzed within a day, launch a CAPICLEAN cleaning sequence once a day.

See the instruction sheets of CAPICLEAN CAPILLARYS 3 and the instruction manual of CAPILLARYS 3, SEBIA.

Storage, stability and signs of deterioration

Store CAPICLEAN refrigerated (between 2 and 8 °C). It is stable until the expiration date indicated on the vial label. DO NOT FREEZE. Precipitate or combined particles in suspension (flocules) may be observed in the CAPICLEAN vial without any adverse effects on its utilization. Do not dissolve this precipitate or these particles. It is recommended to collect only the supernatant.

3. SODIUM HYPOCHLORITE SOLUTION (for sample probe cleaning)

Preparation

Prepare a sodium hypochlorite solution (2 % to 3 % chloride) by diluting 250 mL 9.6 % chloride concentrated solution to 1 liter with cold distilled or deionized water.

Use

For the sample probe cleaning in the CAPILLARYS 3 instrument (weekly maintenance in order to eliminate adsorbed proteins from the probe).

See the instruction manual of CAPILLARYS 3, SEBIA.

Storage, stability and signs of deterioration

Store the working chlorinated solution at room temperature (15 to 30 °C) in a closed container, it is stable for 3 months. Avoid storage in sunlight, close to heat and ignition source, and to acids and ammonia.

4. WASH SOLUTION (registration certificate in Russia PY № ФС3 2009/05453 from 25.11.2009)

Preparation

The vial of the stock wash solution should be diluted up to 750 mL with distilled or deionized water.

After dilution, the wash solution contains an alkaline solution pH $\approx$ 12.

Use

For washing the capillaries after electrophoretic separation.

IMPORTANT :

When wash solution vial replacement, change systematically the filter. Wear clean gloves for handling and installation of the filter.

- Before placing the wash solution vial in the instrument, it is recommended to wash the connector and the tube with plenty of distilled or deionized water to avoid salts deposit.
- Screw the filter at the connector situated at the extremity of the tube plunging in the wash solution vial.

Storage, stability and signs of deterioration

Store the stock and working wash solutions in closed containers at room temperature (15 to 30 °C) or refrigerated (between 2 and 8 °C). The stock wash solution is stable until the expiration date indicated on the kit or wash solution vial label.

Working wash solution is stable for 3 months.

After dilution and immediate installation of the vial in the instrument, the solution is stable for 3 months (if the working wash solution is stored out of the instrument before use, this time of 3 month storage must take into account the time during which the solution is stored outside the instrument). Discard working wash solution if it changes its appearance, e.g., becomes cloudy due to microbial contamination.

NOTES :

The assays that were performed for the validation of reagents demonstrated that, for the different solutions and using an adapted equipment for the reconstitution volume, a variation of $\pm 5\%$ on the final volume has no adverse effect on the analysis. The distilled or deionized water used to reconstitute solutions, must be free of bacterial proliferation and mold (use a filter $\leq 0.45\ \mu\text{m}$) and have a conductivity lower than $3\ \mu\text{S/cm}$, which corresponds to a resistivity higher than $0.33\ \text{M}\Omega\cdot\text{cm}$.

5. Normal control serum for electrophoresis (NORMAL CONTROL) (registration certificate in Russia PY № ФСЗ 2009/05453 from 25.11.2009)

Composition

The control serum is obtained from a pool of normal human sera. The serum is in a stabilized lyophilized form.

NOTE: Depending on the vials, the lyophilizate may present a variable color and its aspect may be more or less compact, without any adverse effects on the functionality of the reconstituted control.

Intended use

This serum is designed for the quality control of electrophoretic quantification of human serum proteins, lipoproteins, cholesterol and apolipoproteins on HYDRAGEL and for serum proteins on CAPILLARYS and MINICAP instruments.

Procedure

Reconstitute the lyophilized serum with 1 mL of distilled or deionized water.

Allow to stand for 30 minutes and mix gently (avoid formation of foam).

The reconstituted Normal Control Serum should be used as any human serum sample.

Storage, stability and signs of deterioration

Before reconstitution, store the lyophilized Normal Control Serum refrigerated (at 2 to 8 °C). It is stable until the expiration date indicated on the vial labels or the package.

Store the reconstituted Normal Control Serum at 2 - 8 °C, it should be used as any human serum sample. Due to the risk of microbial contamination and denaturation, use the reconstituted serum within one week.

The reconstituted serum may also be frozen (in aliquots) and stored at - 18 / - 30 °C for 6 months maximum. Each aliquot of the control serum may be frozen / thawed 20 times ; the freezing must be performed immediately after use and the thawing (at 15 - 30 °C) immediately before use.

6. Hypergamma control serum for electrophoresis (HYPERGAMMA CONTROL)

(registration certificate in Russia PY № ФС3 2009/05453 from 25.11.2009)

Composition

The control serum is obtained from a pool of normal human sera. The serum is in a stabilized lyophilized form.

NOTE: Depending on the vials, the lyophilizate may present a variable color and its aspect may be more or less compact, without any adverse effects on the functionality of the reconstituted control.

Intended use

This serum is designed for the quality control of electrophoretic quantification of human serum proteins, lipoproteins, cholesterol and apolipoproteins on HYDRAGEL and for serum proteins on CAPILLARYS and MINICAP instruments.

Procedure

Reconstitute the lyophilized serum with 1 mL of distilled or deionized water.

Allow to stand for 30 minutes and mix gently (avoid formation of foam).

The reconstituted Hypergamma Control Serum should be used as any human serum sample.

Storage, stability and signs of deterioration

Before reconstitution, store the lyophilized Hypergamma Control Serum refrigerated (at 2 to 8 °C). It is stable until the expiration date indicated on the vial labels or the package. Store the reconstituted Hypergamma Control Serum at 2 - 8 °C, it should be used as any human serum sample. Due to the risk of microbial contamination and denaturation, use the reconstituted serum within one week.

The reconstituted serum may also be frozen (in aliquots) and stored at - 18 / - 30 °C for 6 months maximum. Each aliquot of the control serum may be frozen / thawed 20 times; the freezing must be performed immediately after use and the thawing (at 15 - 30 °C) immediately before use.

Allow to stand for 30 minutes and mix gently (avoid formation of foam).

The reconstituted Normal Control Serum should be used as any human serum sample.

Storage, stability and signs of deterioration

Before reconstitution, store the lyophilized Normal Control Serum refrigerated (at 2 to 8 °C). It is stable until the expiration date indicated on the vial labels or the package.

Store the reconstituted Normal Control Serum at 2 - 8 °C, it should be used as any human serum sample. Due to the risk of microbial contamination and denaturation, use the reconstituted serum within one week.

The reconstituted serum may also be frozen (in aliquots) and stored at - 18 / - 30 °C for 6 months maximum. Each aliquot of the control serum may be frozen / thawed 20 times ; the freezing must be performed immediately after use and the thawing (at 15 - 30 °C) immediately before use.

6. Hypergamma control serum for electrophoresis (HYPERGAMMA CONTROL)

(registration certificate in Russia PY № ФС3 2009/05453 from 25.11.2009)

Composition

The control serum is obtained from a pool of normal human sera. The serum is in a stabilized lyophilized form.

NOTE: Depending on the vials, the lyophilizate may present a variable color and its aspect may be more or less compact, without any adverse effects on the functionality of the reconstituted control.

Intended use

This serum is designed for the quality control of electrophoretic quantification of human serum proteins, lipoproteins, cholesterol and apolipoproteins on HYDRAGEL and for serum proteins on CAPILLARYS and MINICAP instruments.

Procedure

Reconstitute the lyophilized serum with 1 mL of distilled or deionized water.

Allow to stand for 30 minutes and mix gently (avoid formation of foam).

The reconstituted Hypergamma Control Serum should be used as any human serum sample.

Storage, stability and signs of deterioration

Before reconstitution, store the lyophilized Hypergamma Control Serum refrigerated (at 2 to 8 °C). It is stable until the expiration date indicated on the vial labels or the package.

Store the reconstituted Hypergamma Control Serum at 2 - 8 °C, it should be used as any human serum sample. Due to the risk of microbial contamination and denaturation, use the reconstituted serum within one week.

The reconstituted serum may also be frozen (in aliquots) and stored at - 18 / - 30 °C for 6 months maximum. Each aliquot of the control serum may be frozen / thawed 20 times; the freezing must be performed immediately after use and the thawing (at 15 - 30 °C) immediately before use.

EQUIPMENT AND ACCESSORIES REQUIRED

1. System of automatic capillary electrophoresis CAPILLARYS 3, in variants:

A. CAPILLARYS 3 OCTA, ref. 1245;

B. CAPILLARYS 3 TERA, ref. 1246;

connected to the computer with installed data processing software PHORESIS.

2. Sample racks supplied with automatic capillary electrophoresis instrument CAPILLARYS 3:

- CAPILLARYS 3 & MC SAMPLE RACKS

- CAPILLARYS 3 & MC CALIBRATORS / CONTROLS RACK (Sample rack No. 0)

- CAPILLARYS 3 & MC SHUTDOWN RACK.

3. Water Container

4. Wash Container

5. Waste Container internal

6. CAPILLARYS 3 & MC SWITCH RACK FOR PROTEIN(E) 6, to launch automatically a technique change to PROTEIN(E) 6 procedure on the CAPILLARYS 3 instrument.

7. CAPI 3 REAGENT CUPS, SEBIA, to prepare biological samples to analyze with the automated instrument. One reagent cup is intended for the analysis of 8 samples with CAPILLARYS 3 OCTA and 12 samples with CAPILLARYS 3 TERA

WARNING : After use, reagent cups with biological samples have to be handled with care. When the analysis is completed, reagent cups must be discarded with biological waste products and they must NEVER be reused.

Storage : Before use, store the reagent cups in their sealed package in a clean and dry place and at a temperature comprised between 2 and 30 °C.

8. CAPI 3 BINS FOR USED REAGENT CUPS, SEBIA, for automated collection of used reagent cups in CAPILLARYS 3 OCTA and CAPILLARYS 3 TERA. To place at the location intended for this purpose.

WARNING : Bins containing used reagent cups with biological samples have to be handled with care.

9. Collection tubes for samples without anticoagulant with a diameter of 13-16 mm and a height of 75-100 mm (placed in the device without a lid) or micro-tubes with a capacity of 1.5 ml, located on empty test tubes for samples with a diameter of 13 mm, used as holders.

WORK WITH SAMPLES FOR ANALYSIS

Sample collection and storage

Fresh serum samples are recommended for analysis. Serum must be collected following established procedures used in clinical laboratory testing. Samples can be stored up to 10 days between 2 and 8 °C.

For longer storage, samples should be frozen at - 18 / - 30 °C within 8 hours of collection. Frozen sera are stable for 2 months.

Proteins of the samples stored at 2 to 8 °C or between 15 and 30 °C, degrade, particularly the C3 complement for which the degradation kinetics is very rapid at 15 - 30 °C and is clearly visible beyond 3 days.

A serum stored between 2 and 8 °C or between 15 and 30 °C has a beta-2 fraction that

gradually decreases and may appear distorted (with small additional fractions appearing on the gamma side and / or beta-1 following the deterioration of C3 complement) and an alpha-2 fraction whose shape can be slightly changed.

Beyond 10 days between 2 and 8 °C or 3 days between 15 and 30 °C, the beta-1 fraction deforms by expanding, and the beta-2 fraction strongly decreases.

Depending to the samples, during their storage beyond 10 days at 2 to 8 °C or 3 days at 15 and 30 °C, the automated integration of fractions by the software for data processing may potentially be disturbed.

NOTE : *Each laboratory must ensure that the samples are transported in optimal conditions for their integrity (1).*

(1) ISO 15189 : Medical laboratories - Requirements for quality and competence.

Sample preparation

Use undiluted serum samples.

Upon storage at 2 to 8 °C or after freezing, some serum (particularly those containing cryoglobulin or cryogel) may become viscous or develop turbidity. At room temperature, these samples can be directly analyzed.

Samples containing a polymerized immunoglobulin may be used without any treatment. It is advised to observe the serum aspect before analysis (cases of hemolysis, cryoglobulins or turbidity).

Samples to avoid

- Avoid hemolysed serum samples. Hemolysis induces a double alpha-2 zone.
- Avoid aged, improperly stored serum samples, beta-2 fraction would be decreased.
- Avoid plasma samples. Fibrinogen migrates in beta-2 position (shoulder on beta-2 or superimposed with the beta-2 zone with possibly an increase of this fraction). When it is present in some samples (plasma, serum not totally defibrinated or patient with anticoagulant treatment), fibrinogen may interfere on the analysis and makes interpretation inaccurate (suspicion of monoclonal band or beta-2 fraction increase). When analyzing an aged plasma sample (not recommended), the C3 complement which is labile over the time is partially degraded, the beta-2 zone then corresponds essentially to fibrinogen.

NOTE : *The technical characteristics of biological sample collection tubes are described in the user manual of the CAPILARYS 3 system and in the available documentation on pre-analytical phase for bio-medical analysis (data provided by the tube manufacturers, guides and recommendations on biological sample collection...).*

PROCEDURE

The CAPILARYS 3 system is a multiparameter instrument for serum proteins analysis on 8 or 12 parallel capillaries in the following sequence :

- sample racks identification by RFID (Radio Frequency Identification),

- bar code reading of sample tubes (for up to 8 tubes),
- sample dilution from primary tubes into reagent cups,
- capillary washing,
- injection of diluted samples,
- protein analysis and direct detection on capillaries.

The manual steps include :

- placement of sample tubes in sample racks,
- placement of sample racks on the CAPILLARYS 3 instrument,
- removal of sample racks and sample tubes after analysis,
- removal of bins for used reagent cups.

PLEASE CAREFULLY READ THE CAPILLARYS 3 INSTRUCTION MANUAL.

I. PREPARATION OF ELECTROPHORETIC ANALYSIS

1. Switch on CAPILLARYS 3 instrument and computer.
2. Wait until the instrument is completely initialized.
3. Start the PHORESIS software installed on the computer for data processing.
4. The CAPI 3 PROTEIN(E) 6 kit is intended to run with "PROTEIN(E) 6" analysis program from the CAPILLARYS 3 instrument.

To select "PROTEIN(E) 6" analysis program and place the CAPILLARYS PROTEIN(E) 6 buffer vial in the instrument, please read carefully the CAPILLARYS 3 instruction manual. If necessary, place the vial with the reconstituted wash solution in the instrument.

5. The sample rack contains 8 positions for sample tubes. Position up to 8 sample tubes on each sample rack ; the bar code of each tube must be visible in the openings of the sample rack.
6. Take a pack of new reagent cups by holding the handle and place it on the automated loading system for cups of CAPILLARYS 3 ; then, remove the flange (a message will be displayed when reagent cups are missing).
7. Place a new bin for used reagent cups into the CAPILLARYS 3 instrument at the location intended for this purpose.
8. Slide the sample rack(s) into the CAPILLARYS 3 instrument through the opening in the right side of the instrument. Up to 15 sample racks can be introduced successively and continuously into the instrument.

NOTE : When analyzing a control serum, it is advised to use the sample rack No. 0 intended for control serum.

9. Remove analyzed sample racks from the plate on the left side of the instrument.
 10. If necessary, take off carefully the bin containing used reagent cups and discard it.
- WARNING : Bins containing used reagent cups with biological samples have to be handled with care.*

DILUTION - MIGRATION - DESCRIPTION OF THE AUTOMATED STEPS

1. Sample rack identification by RFID.
2. Bar codes are read on primary sample tubes.
3. Serum samples are diluted in analysis buffer and the sample probe is rinsed after each sample.
4. Capillaries are washed.
5. Diluted samples are injected into capillaries.
6. Migration is carried out under constant voltage, controlled by Peltier effect for about 4 minutes.
7. Proteins are detected directly by scanning at 200 nm and data of the obtained protein electrophoretic pattern are transmitted from the instrument to the computer equipped with the software for data processing.

II. RESULT ANALYSIS

At the end of each analysis, the corresponding data are transmitted by the instrument to the software for data processing and an electrophoretic profile appears on the screen of the computer. Relative quantification of individual zones is made automatically and profiles can be analyzed.

With the total protein concentration, the software will calculate each fraction concentration. The electrophoregrams are interpreted visually for pattern abnormalities.

Electrophoretic profiles are visualized by default using the re-drawn mode : then, the alpha-1 fraction is closer to albumin. Optionally, the standard mode allows to visualize the initial pattern obtained with raw data.

III. END OF ANALYSIS SEQUENCE

At the end of each analysis sequence, the operator must start the shutdown procedure of the CAPILLARYS 3 instrument in order to store capillaries in optimal conditions.

IV. FILLING OF REAGENT CONTAINERS AND MANAGEMENT OF DISPOSABLES

The CAPILLARYS 3 instrument has an automatic control for reagents (by using RFID labels) and for disposables (reagent cups and bins for used cups).

IMPORTANT : It is necessary to respect the designed position for wash solution, rinse and waste containers.

On the screen of the CAPILLARYS 3 instrument, the "Main compartment" menu for reagents management displays information when it is necessary to perform one of the following tasks :

- place a new buffer vial and / or,
- fill the container with working wash solution and / or,
- fill the container with filtered distilled or deionized water for rinsing capillaries and / or,
- empty the waste container.

WARNING : Do not use marketed deionized water, such as water for ironing for example (risk of important capillaries damage). Use only water with ultrapure quality, such as injection grade water.

IMPORTANT : Before filling the rinse container, it is recommended to wash it with plenty of distilled or deionized water.

QUALITY CONTROL

It is advised to include a control serum with each sequence of analysis.

RESULTS

Values

Direct detection at 200 nm in capillaries yields relative concentrations (percentages) of individual protein zones.

Reference values (mean $\pm$ 2 SD) for individual major electrophoretic serum protein zones have been established from a healthy population of 246 adults with normal triglycerides levels (men and women) :

	CAPI 3 PROTEIN(E) 6
Albumin	55.8 – 66.1 %
Alpha-1 globulins	2.9 – 4.9 %
Alpha-2 globulins	7.1 – 11.8 %
Beta globulins (cumulative)	8.4 – 13.1 %
Beta-1 globulins	4.7 – 7.2 %
Beta-2 globulins	3.2 – 6.5 %
Gamma globulins	11.1 – 18.8 %

It is recommended each laboratory establish its own reference values.

NOTE : Reference values have been established using the standard parameters of the PHORESIS software (smoothing 2 and automatic drift).

Interpretation

A distortion of the electrophoretic pattern compared to a normal one indicates an abnormality, especially the appearance of an additional thin peak in the gammaglobulins zone.

The C4 complement migrates between beta-1 and beta-2 zones ; CRP migrates in beta-2 position, see ELECTROPHORETIC PATTERNS.

A relative increase of the beta-2 zone compared to the beta-1 zone, without any clinical context of inflammatory disease, must be a warning signal for necessary complementary analyses.

In case of doubt concerning the interpretation of the pattern and / or the positioning of minima (particularly during the analysis of an external control), it is necessary to overlay the obtained pattern with that of the control serum for electrophoresis NORMAL CONTROL (registration certificate in Russia PY № ФC3 2009/05453 from 25.11.2009). A monoclonal component may be suspected in the serum sample when a single protein

electrophoretic pattern is delayed or distorted or in the case of impossibility for the software to redraw the albumin / alpha-1 zone. The following warning message is then displayed on the electrophoretic pattern "Warning: Migration time out of range" with a red warning signal. This red warning signal is also displayed on the curves mosaic and in the result table for the sample concerned. To confirm the presence of a monoclonal component in such sample, it is necessary to treat the sample with beta-mercaptoethanol and to repeat the analysis on the sample after reducing treatment. In this case (i) prepare 1 % beta-mercaptoethanol (BME, or 2-mercaptoethanol, 2 ME) in Fluidil (registration certificate in Russia PY № ФС3 2009/05453 from 25.11.2009), (ii) the CAPILLARYS 3 instrument ready waiting for rack, add 100 µL of this reducing solution to 300 µL neat serum, (iii) vortex and wait for 15 minutes maximum, then follow the standard procedure.

IMPORTANT : After reducing treatment with beta-mercaptoethanol, the sample must be analyzed without any delay ; no introduced sample rack must be waiting for analysis in the CAPILLARYS 3 instrument.

When many electrophoretic patterns show the same warning signal, call SEBIA Technical Service. As an aid in interpretation of serum protein electrophoregrams, see BIBLIOGRAPHY.

Alpha-2 zone :

In some samples and according to the haptoglobin phenotype, alpha-2 zone can be split, see ELECTROPHORETIC PATTERNS.

Interference and Limitations

See SAMPLES FOR ANALYSIS.

Lipoproteins / triglycerides or biliary pigments (with a characteristic yellow – green color of the serum) at high concentration in the sample may lead to the visual impression of a bisalbuminemia on the electrophoretic pattern.

In the case of suspected contamination between two samples (very rare), due to the presence of some monoclonal components (in high concentration for example), it is recommended to repeat the test on those investigated samples (contaminating and potentially contaminated) by reversing their analysis order.

Due to the resolution and sensitivity limits of zone electrophoresis, it is possible that some monoclonal components may not be detected with this method.

A monoclonal component may be not detected (i.e., polymerized immunoglobulin spread or hidden in the polyclonal background). Conversely, a slight distortion of the electrophoretic pattern may indicate the presence of a monoclonal immunoglobulin. In all cases, the clinical context must be analyzed and if a gammopathy is suspected, it is then recommended to perform an immunotyping analysis on the sample. If an uncertainty persists, confirm the result by an immunofixation technique on agarose gel.

Troubleshooting

In case the test fails to perform while the instruction for the preparation and storage of materials, and for the procedure were carefully followed, contact the technical support

PERFORMANCE DATA

	Parameter	Specification
1.	Coefficient of variation:	
1.1.	Repeatability, CV, %	
1.1.1.	Albumin	CV ≤ 2%
1.1.2.	Alpha-1 globulins	CV ≤ 10 %
1.1.3.	Alpha-2 globulins	CV ≤ 10 %
1.1.4.	Beta-1 globulins	CV ≤ 10 %
1.1.5.	Beta-2 globulins	CV ≤ 10 %
1.1.6.	Gamma globulins	CV ≤ 10 %
1.2.	Reproducibility between capillaries, CV, %	
1.2.1.	Albumin	CV ≤ 2%
1.2.2.	Alpha-1 globulins	CV ≤ 10 %
1.2.3.	Alpha-2 globulins	CV ≤ 10 %
1.2.4.	Beta-1 globulins	CV ≤ 10 %
1.2.5.	Beta-2 globulins	CV ≤ 10 %
1.2.6.	Gamma globulins	CV ≤ 10 %
1.3.	Acceptable variation of results for parallel definitions of the same sample by different products of the same lot, %	
1.3.1.	Albumin	Acceptable variation ≤ 2,0 %
1.3.2.	Alpha-1 globulins	Acceptable variation ≤ 2,0 %
1.3.3.	Alpha-2 globulins	Acceptable variation ≤ 2,0 %
1.3.4.	Beta-1 globulins	Acceptable variation ≤ 2,0 %
1.3.5.	Beta-2 globulins	Acceptable variation ≤ 2,0 %
1.3.6.	Gamma globulins	Acceptable variation ≤ 2,0 %
2.	Linearity	+/- 10% of the assigned value in the range of Albumin concentrations from 10 g/dl to 49.9 g/dl

Analytical specificity

Interference :

Lipoproteins / triglycerides or biliary pigments (with a characteristic yellow – green color of the serum) at high concentration in the sample may lead to the visual impression of a bisalbuminemia on the electrophoretic pattern.

Interfering factor	Specification
Bilirubin (≤ 342 mM/l or 20 mg/dl)	Acceptable relative deviation for the fraction: Albumin < 4% Alpha-1 < 10 % Alpha-2 < 10 % Beta-1 < 10 % Beta-2 < 10 % Gamma < 10 %
Triglycerides (≤ 37 mM/l or 33 g/l)	Acceptable relative deviation for the fraction: Albumin < 4% Alpha-1 < 10 % Alpha-2 < 10 % Beta-1 < 10 % Beta-2 < 10 % Gamma < 10 %

Sensitivity

Serial dilutions of one serum sample with an abnormal protein in gamma zone at 1.198 g/dL was electrophoresed using the CAPI 3 PROTEIN(E) 6 procedure.

The highest dilution with a discernible abnormality corresponded to a concentration of 19 mg/dL of the abnormal protein.

NOTE : According to the position of the abnormal protein and polyclonal background in the gamma zone, the detection limit may vary.

Accuracy – Internal correlation

The internal concordance study of the CAPI 3 PROTEIN(E) 6 procedure performed with the CAPILLARYS 3 instrument was evaluated in a study based on the Clinical Laboratory Standards Institute (CLSI - USA) EP09-A2 guideline "Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition (Interim Revision)".

The results for percentages (%) of each protein fraction were analyzed using statistical tools recommended by CLSI.

NOTE : The results presented below have been obtained from 1 internal accuracy study. The analyzed serum samples were provided by 2 laboratories in France.

All the samples were exactly treated the same way with both techniques and followed

the same guidelines in regards to sample integrity.

The levels of each protein fraction were measured in 118 serum samples, including 18 normal samples and 100 pathological samples, both by electrophoretic separations obtained with the CAPI 3 PROTEIN(E) 6 procedure on the CAPILLARYS 3 instrument and another commercially available capillary electrophoresis technique for proteins analysis. The measured values of percentages of the 6 protein fractions from both procedures were analyzed by a linear regression statistical procedure. The obtained results have demonstrated a perfect correlation between both procedures for proteins quantification.

The results of linear regression analysis are tabulated below ($y = \text{CAPI 3 PROTEIN(E) 6}$), the sensibility and specificity of the CAPI 3 PROTEIN(E) 6 procedure compared to the reference procedure have been calculated using the recommended method (Wendling, 1986).

Fraction	Correlation coefficient	y-Intercept	Slope	Range of % values CAPI 3	Sensibility (%)	Specificity (%)
Albumin	0.999	- 0.230	1.013	24.3 – 74.7	93.7	100.0
Alpha-1	0.998	- 0.056	0.981	2.0 – 16.1	90.7	98.4
Alpha-2	0.998	- 0.383	1.016	4.2 – 23.1	97.6	98.7
Beta-1	0.998	- 0.131	1.004	2.3 – 41.0	78.1	98.8
Beta-2	1.000	0.028	1.013	1.3 – 45.7	96.9	98.8
Gamma	1.000	- 0.329	1.014	1.7 – 64.2	97.1	100.0

Contraindications to use.

Not applicable to in vitro diagnostic products

Possible side effects.

Not applicable to in vitro diagnostic products

Transportation and Use conditions:

The product can be transported by all types of covered transport in accordance with the rules of cargo transportation applicable to this type of transport, providing a transport temperature of 2-30 °C in a place protected from light, at a normal pressure of 101.3 kPa (760 mm Hg), and relative humidity of up to 80%.

Use conditions:

- temperature: 15-30 °C
- pressure 101.3 kPa (760 mmHg)
- relative humidity up to 80%.

LIST OF APPLICABLE STANDARDS:

EN ISO 13485:2016 ISO 13485	<i>Medical devices - Quality management systems - Requirements for regulatory purposes.</i>
EN 13612:2002 EN 13612:2002/AC:2002	<i>Performance evaluation of in vitro diagnostic medical devices.</i>
EN ISO 23640:2015 ISO 23640:2011	<i>In vitro diagnostic medical devices - Stability testing of in vitro diagnostic reagents.</i>
EN 13641: 2002	<i>Elimination or reduction of risk of infection related to in vitro diagnostic reagents.</i>
EN 13975:2003	<i>Sampling procedures used for acceptance testing of in vitro diagnostic medical devices. Statistical aspects.</i>
EN 14136:2004	<i>Use of external quality assessment schemes in the assessment of the performance of in vitro diagnostic examination procedures.</i>
EN ISO 14971:2012 ISO 14971:2007	<i>Medical devices — Application of risk management to medical devices.</i>
EN ISO 15193:2009 ISO 15193:2009	<i>In vitro diagnostic medical devices — Measurement of quantities in samples of biological origin — Requirements for content and presentation of reference</i>
EN ISO 15194:2009 ISO 15193:2009	<i>In vitro diagnostic medical devices — Measurement of quantities in samples of biological origin — Requirements for certified reference materials and the content of supporting documentation.</i>
EN ISO 15223-1:2016 ISO 15193:2016	<i>Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General requirements</i>
EN ISO 18113-1:2011 ISO 18113-1:2011	<i>In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling). Part 1: Terms, definitions and general requirements</i>
EN ISO 18113-2:2011 ISO 18113-2: 2009	<i>In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling). Part 2: In vitro diagnostic reagents for professional use</i>

PRECAUTIONS FOR USE IN RUSSIA:

Using the kit follow the requirements of GOST R 52905-2007 " medical Laboratories. Safety requirement».

Only personnel who have been trained to work with the product and equipment should be allowed to work with the kit.

Before starting work, read the instructions and use the kit strictly for its intended purpose.

Since the product is necessary for the analysis of potentially infected materials, when working with the kit, use personal protective equipment that is allowed for use in the territory of the Russian Federation. It is recommended to use protective gloves, safety glasses, and work protective clothing.

It is forbidden to eat, use cosmetics and smoke in the premises intended for working with sets.

Clean and disinfect spilled samples or reagents using disinfectants in accordance with MU-287-113.

Remove unused reagents in accordance with the requirements of SanPiN 2.1.7.2790-10.

All used disposable materials and sample residues must be treated with disinfectants and then disposed of in accordance with the requirements of SanPiN 2.1.7.2790-10 and MU-287-113.

If liquid components of the kit come into contact with the skin: wash off with plenty of soap and water. In case of contact with eyes: rinse with plenty of water, then consult a doctor.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN, MEDICINES

The product does not contain materials of animal or human origin. The product does not contain medicinal products

DATA OF SERVICE LIFE OR SHELF LIFE.

The shelf life of the kit of reagents is determined by the reagent component of the kit with the nearest expiration date. The maximum shelf life of the kit is 3 years from the date of manufacture. To each vial is given in the form of a label "exp:".

After opening, when stored on board the analyzer, use the buffer for 60 days.

It is forbidden to use the product after the expiration date.

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