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The Hand Piece is designed to convert electrical energy from a compatible HARMONIC Generator to mechanical motion for the instrument blades. This Hand Piece is intended for use with a compatible HARMONIC Generator.

The Hand Piece is permanently attached to a cord which connects to the front of the Generator.

The Hand Piece is a re-useable instrument with a limited service life. The Hand Piece is programmed with a counter to limit the service life to 95 procedures. The generator will give a Hand Piece error after 95 procedures are completed. The number of activations during a procedure is not limited, and the counter will not log a procedure until the Hand Piece is unplugged from the generator or the generator is powered down. The number of remaining procedures is indicated on the display of the generator during power up with the generator software 1.09 and via Bio Medical Mode (BME) with the prior versions of generator software.

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The Hand Piece is packaged with a Blade Wrench and Test Tip. The Blade Wrench is used to secure instrument blades to the Hand Piece. The Hand Piece, Blade Wrench, and Test Tip are packaged non-sterile and must be sterilized per their insert instructions prior to use. Refer to a compatible HARMONIC Generator User Manual before using these instruments.

5.2. Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instru-

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ment (refer to Warnings and Precautions).

1. Prior to use, remove and discard the protective cover over the Blade Mount Surface.
2. Before initial use, clean and sterilize the Hand Piece as per recommended Instructions for Reprocessing included in this insert.
3. Prior to each use, visually inspect the Hand Piece for damage that may have occurred during handling and cleaning, and verify that the instrument functions properly and is not in need of maintenance or repair. Refer to Inspection section for details of the visual inspection.

5.3. Hand Piece Attachment

Refer to a compatible HARMONIC Generator User Manual for Hand Piece attachment and system operation instructions.

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– GEN11 Generator (Registration certificate № ФСЗ 2012/11599)

Handle used with HARMONIC tools:

HBC05	Ball Coagulator, 5mm x 31cm
HC325	Curved, 5mm, 5mm x 32cm
HDH05	Dissecting Hook, 5mm x 32cm
WAVE18S	HARMONIC Wave
ACE14S	Curved Shear, Scissor Handle, 5mm x 14cm
ACE23E, ACE23J	Curved Shears with Ergonomic Handle, 5mm x 23cm
ACE36E, ACE36J	Curved Shears with Ergonomic Handle, 5mm x 36cm
ACE23P	Curved Shears, 5mm x 23cm Rationalized!

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
HARMONIC Hand Piece

Main elements:









1. Test Tip
 2. Blade Wrench
 3. Hand Piece
 4. Blade Mount Surface
- The Blade Wrench is used to secure instrument blades to the Hand Piece.

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



Cable assembly:

Coating for Handpiece Cable Assembly: poly(para-xylylene) Parylene N
 Strain Relief Overmold for Handpiece: Polydimethyl/methylvinylsiloxane elastomer NuSil MED-4950
 End Cap: Aluminum 6061-T6 ASTM B209
 End Cap – Coating: Metallic Chromium Coating Metallic Chromium (99.0% Chromium and Oxygen)

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	Humidity Limitation
	Authorized Representative in the European Community
	Authorized Representative in the United States
	Electrical and electronic equipment. Return waste to a collection system or treatment and recycling facilities. Applicable in the EU. Follow decontamination instructions before returning waste.
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. *
	Unit Quantity
	Serial number
	«PCT» Compliance Symbol
	Indicates that the Electrical and Electronic Product (EEP) contains toxic and hazardous substances or elements over the defined maximum concentration values defined above. The names and content of the toxic and hazardous substances or elements shall be provided in the product instruction manual. The number "10" in the symbol refers to the environment friendly use period; after their environmentfriendly use period expires, they shall be put into the recycling system (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products). *
	Fragile. Caution!
	This way up.
	Caution! Refer to the instructions for use

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Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Audible high-pitched tones are an abnormal condition and an indicator that the blade or Hand Piece is not operating properly. The tones may be an indicator that the Hand Piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high temperatures and user or patient injury.
- Do not use the HARMONIC Hand Piece without the proper adaptor. Failure to use the proper adaptor may result in user or patient burn injury.
- To prevent burn injury, discontinue use if the Hand Piece temperature becomes uncomfortable to hold.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by the activation of the HARMONIC instruments in fatty tissue are potentially flammable.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blades should not be in contact with the patient, drapes or flammable materials while not in use.

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- During and following activation in tissue, the instrument blades may become hot. Avoid unintended blade contact with tissue, drapes, surgical gowns, or other unintended sites after activation.

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intended for patient contact. To prevent burn injury, avoid direct tissue contact with the Hand Piece and adaptor or take preventative measures to protect tissue that comes in contact with the Hand Piece and adaptor. Refer to the instrument inserts for applicable duty cycle (proportion of ON/OFF time) information.

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- Verify compatibility with generators prior to use. This handpiece is not compatible with the ULTRACISION® Generator (GEN01/GEN32).

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- Do not use this Hand Piece with the HS2 or DSH10 blade system. Use of these blades on the Hand Piece may result in excessive temperatures and risk of burns.

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11. INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

Not applicable. The devices do not contain medicinal substances and are not used in combination with medicinal substances.

12. DESCRIPTION OF STERILIZATION METHOD

The HARMONIC Hand Piece can withstand 95 sterilization cycles.

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Instructions for Reprocessing

- This device is supplied non-sterile and must be cleaned and sterilized prior to initial and each subsequent use.
- Thorough cleaning is necessary to ensure correct mechanical functioning of the device.
- Each device must be treated individually and matte surfaces must be protected against scratches.
- The design of this device permits effective cleaning and reliable sterilization. The user must ensure that cleaning and sterilization are performed in accordance with guidelines, standards and National Health Authority requirements.

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- Clean device as soon as possible after use.
Remove excessive solids using disposable wipes, enzymatic foam or per hospital procedure. Soiled devices should be separated from unsoiled devices. Soiled Devices should be covered with a towel dampened with tap, sterile or critical water to prevent blood and/or debris from drying. Encrusted soil can increase the time and effort required to complete processing.
- If the recommendations in the steps above are not possible, the device may be placed in a container and immersed in tap, sterile or critical water and covered for transport to the processing area.
Caution: To minimize the chance of corrosion, do not use saline for cleaning.
- Follow the Cleaning and Sterilization instructions below.

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The user must ensure that cleaning and sterilization are conducted in accordance with the appropriate guidelines, standards, or National Health Authority requirements. Hospital grade low foaming neutral pH to mildly alkaline enzymatic detergents, hospital grade low foaming neutral pH detergents (pH 7-9) and hospital grade low foaming mildly alkaline detergents with a pH of ≤ 11 (prepared according to manufacturer's instructions) should be used. An enzymatic detergent and a mildly alkaline detergent were used in cleaning efficacy studies and are types of products that can be used throughout the cleaning process. "Detergent(s)" will be used generically throughout the rest of these instructions. Ensure that detergent manufacturer instructions are followed for concentration, temperature and water quality throughout the cleaning process. When using mildly alkaline detergents, critical water should be used throughout the cleaning process.

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1. Immerse and soak the device for a minimum of 10 minutes up to 60 minutes in an enzymatic detergent or equivalent prepared according to manufacturer's instructions for concentration, water quality and temperature.
2. Remove the device from the detergent.
3. Immerse the device in a clean sink with cold to warm tap water – not to exceed 40°C. Using a toothbrush style soft bristle (nylon) brush or soft cloth, clean the exterior of the device paying particular attention to crevices, threads and seams.

Caution: Do not use scouring pads, abrasive cleaners, or a steel wire brush,

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which may cause damage to the surfaces of the device.

4. Rinse the exterior of the device thoroughly under cold to warm running tap water for a minimum one minute.

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Following the pre-cleaning step, the validated automated or manual processing method must be performed. Refer to standard hospital protocol and local/national practices to determine whether disinfection is required. Note that disinfection is mandatory in some areas of the world; e.g. EU. There are different A0 requirements in various countries; e.g. $A0 \geq 600$ (90°C for 1 minute) or $A0 \geq 3000$ (90°C for 5 minutes). The user must ensure that cleaning, disinfection and sterilization are conducted in accordance with the appropriate local guidelines, standards, or National Health Authority requirements.

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Phase	Recirculation Time (minutes:seconds)	Temperature	Detergent Type	Validated Parameters
Pre-Wash	02:00	Cold tap water	N/A	2 minutes cold tap water
Enzyme Wash	02:00	Tap water <40°C (104°F)	Neutral pH enzymatic detergent or mildly alkaline detergent	no
Wash	03:00	Tap water $\geq 40^\circ\text{C}$ (104°F)	Neutral pH detergent or mildly alkaline detergent	5 minutes tap water 55-60°C
Rinse	02:00	Critical water $> 40^\circ\text{C}$ (104°F)	N/A	2 minutes critical water
*Thermal Disinfection	01:00 – 10:00	Critical water $\geq 90^\circ\text{C}$ (194°F)	N/A	5 minutes at 90°C

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Drying	30:00	≥115°C (239°F) – not To exceed 134°C (273°F)	N/A	30 minutes at 115°C
--------	-------	--	-----	------------------------

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Note: The automated washer program may or may not include an enzymatic wash phase. Single or combined washing phases should be performed for a minimum of 5 minutes. Temperatures used during the enzymatic and/or wash phases should meet the requirements of the detergent manufacturer's recommendations.

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Manual Cleaning Method

1. Immerse the device for a minimum 5 minutes in a detergent prepared in accordance with manufacturer's recommendations for concentration, water quality and temperature.
2. Remove the device from the detergent.
3. Immerse the device in a clean sink with cold to warm critical water – not to exceed 40°C. While immersed, use a toothbrush style, soft bristle (nylon) brush or soft cloth to clean the exterior of the device. Pay particular attention to crevices, threads and seams.

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Manually clean the upper and lower contact rings using a cotton swab

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dipped in isopropyl alcohol until visibly clean.

Inspection

Devices do not have an indefinite functional life. All surgical devices are subject to a degree of wear and tear as a result of normal use. They must therefore undergo a regular and precise visual check before each use:

- Inspect the hand piece blade mount surface for pitting and scratches;
- Inspect the hand piece blade mount surface for cracks, bends, or other significant forms of physical damage;
- Inspect the hand piece threads for physical damage;
- Inspect the blade mount surface for excessive side play;
- Inspect the hand piece body for press fit joint separation;
- Inspect the hand piece body for cracks;
- Inspect the cable for damage;
- Ensure the complete removal of soil from all Hand Piece surfaces. If soil or blood is present, clean the device following the Pre-Cleaning and Automated or Manual Cleaning Instructions.

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Prepare Device for Sterilization

1. Place and wrap device in an approved sterile barrier wrap material per hospital protocol and in accordance with ISO11607-1 and any other local and/or National requirements using standardized wrapping techniques. Or place device in a sterilization tray wrapped using an approved sterile barrier material and standardized wrapping techniques.

2. Place wrapped device or wrapped tray into autoclave.

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3. Run the device through the autoclave using one of the settings below.

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STERRAD® Systems

1. Clean and thoroughly dry the HARMONIC Hand Piece and Test Tip following the steps listed in the Hand Piece and Test Tip Cleaning section.

2. Using standardized wrapping techniques, place and wrap the device in an approved STERRAD® System-compatible polypropylene sterilization wrap in accordance with hospital protocol, the STERRAD® Operator's Manual, and ISO11607-1. Or place device in a sterilization tray wrapped using an approved sterile barrier material and standardized wrapping techniques approved for use in a STERRAD® System. **Do not use paper pouches or sterilization wraps containing cellulose or cotton.**

3. Sterilize the Hand Piece and Test Tip within the STERRAD® System according to the instructions provided in the STERRAD® Operator's Manual.

Refer to the table below for STERRAD® Systems and corresponding sterilization cycle times.

STERRAD® System	Purpose	Cycle Time
STERRAD® 100NX	For general surgical instrumentation	Standard 47 minutes
STERRAD® 100S	Rigid & Tubing	Short 54 minutes
STERRAD NX®	(2x400 mm SS; 1x350 mm tubing)	Standard 28 minutes
STERRAD® 200	Rigid & tubing	Standard 75 minutes

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13. CONDITIONS OF TRANSPORTATION AND STORAGE

Transportation and Storage conditions:

Designation	Range
Temperature	-35 °C to +54 °C
Relative humidity	10 % to 80 %
Barometric pressure	700 hPa to 1060 hPa

Sterilized products should be stored in a dry, clean environment, protected from direct sunlight, dust, pests, and extreme temperature and humidity. Refer to sterilization wrap manufacturer's IFU for limits on sterile product storage time.

14. MD DISPOSAL RULES AND CONDITIONS

Some internal components of the Hand Piece contain Lead (PZT Lead-Zirconate-Titanate) Disposal should be performed according to local requirements and regulations.

The unused device is disposed in conformity with the rules of handling medical wastes of Class A.

The used device must be disposed in conformity with the rules of handling medical wastes of Class B.

15. WARRANTY

A storage warranty period - the manufacturer guarantees the reliability for the intended use, subject to the right transportation, and storage within the number of uses as identified within the product IFU. Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

16. NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION**DEVELOPER:**

Ethicon Endo-Surgery, LLC
4545 Creek Road

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AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

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4. SCOPE AND POTENTIAL CONSUMERS OF MD

Scope: surgery:

The HARMONIC System instruments (including Handpieces) can be used in general, plastic, pediatric, gynecologic, urologic, exposure to orthopedic structures (such as spine and joint space), sealing and transection of lymphatic vessels, ENT (Ears, Nose, Throat), including tissues of the soft palate, oral structures, and oropharyngeal airway, thoracic surgery, including mobilization of the Internal Mammary Artery (IMA), or other open or endoscopic procedures.

Potential consumers of MD:

Qualified Surgeons.

4.1. Use/operation conditions of MD

Operation conditions	
Humidity	10 % to 85 %
Temperature	+10 °C to +30 °C
Pressure	700 hPa to 1060 hPa

5. DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS AND PRINCIPLE OF OPERATION OF MD

5.1. General description of MD

The hand piece contains an acoustic transducer that converts the electrical energy supplied by the generator to mechanical motion. The transducer is connected to an ultrasonic wave guide/amplifier which amplifies the motion produced by the transducer and relays it to the end effector.

The mechanical motion from the acoustic amplifier advances to the end effector, transmitting high frequency ultrasonic energy which enables hemostatic cutting and/or coagulation of tissue.

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Device Description

The Harmonic Blue Hand Piece is designed to convert electrical energy from a compatible Harmonic Generator to mechanical motion for the instrument blade. The Blue Hand Piece is intended for use exclusively with a compatible Harmonic Generator. The Blue Hand Piece is permanently attached to a blue cord which connects to the front of the Generator.

The Blue Hand Piece is a re-useable instrument with a limited service life. The Blue Hand Piece is programmed with a counter to limit the service life to 100 procedures. The generator will give a Hand Piece error after 100 procedures are completed. The number of activations during a procedure is not limited, and the counter will not log a procedure until the hand piece is unplugged from the generator or the generator is powered down. The number of remaining procedures is indicated on the display of the generator during power up with generator software 1.09 and via Bio Medical Mode (BME) with prior versions of generator software.

Note: The generator is placed in BME mode by depressing the Standby and Power Level Down Arrow buttons simultaneously while powering up the generator. (Hold down both buttons until power-up sequence is completed).

The Harmonic Blue Test Tip is compatible with the torque wrench packaged with the single patient use instruments.

5.2. Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

- 1 Prior to use, remove and discard the protective cover over the Blade Mount Surface.
- 2 Before initial use, clean and sterilize the Hand Piece and test tip as per

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recommended Instructions for Reprocessing included in this insert.

3 Prior to each use, visually inspect the Hand Piece for damage that may have occurred during handling and cleaning, and verify that the instrument functions properly and is not in need of maintenance or repair. Refer to Inspection section for details of the visual inspection.

5.3. Hand Piece Attachment
Refer to a compatible Harmonic Generator User Manual for Hand Piece attachment and system operation instructions.

5.4. The possibility and ways of integrating with other medical devices

Handles are designed for use with a compatible HARMONIC generator:

– GEN11 Generator (Registration certificate № ФС3 2012/11599)

Handle used with HARMONIC tools:

SNGCB	HARMONIC SYNERGY Dissecting Hook, Curved Blade and Torque Wrench
SNGHK	HARMONIC SYNERGY Dissecting Hook, Curved Blade and Torque Wrench
SNGHK2	HARMONIC SYNERGY Combination Hook Blade with Torque Wrench
HAR9F	HARMONIC FOCUS® Shears + Adaptive Tissue Technology
HAR17F	HARMONIC FOCUS® Long Shears + Adaptive Tissue Technology

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6. LIST AND DESCRIPTION OF MD MATERIALS

Body:

Mid Housing: Aluminum 6061-T6 ASTM

Mid Housing: Metallic Chromium Coating Metallic Chromium (99.0% Chromium and Oxygen).

Alternative:

Mid Housing: Chromium Coating / HEEF 25 R (500)

Nose Cone Assembly:

Nose Cone: Aluminum 6061-T6 ASTM B209

Nose Cone: Metallic Chromium Coating Metallic Chromium (99.0% Chromium and Oxygen)

Alternative:

Nose Cone: Chromium Coating / HEEF 25 R (500)

Cable:

End Cap: Aluminum 6061-T6 ASTM B209

End Cap: Metallic Chromium Coating Metallic Chromium (99.0% Chromium and Oxygen)

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

Outlet Housing: Polyphenylsulfone RADEL R-5800 NT
Cable – Coating: poly(para-xylylene) Parylene N
LIM Connector: Polydimethyl / methylvinylsiloxane (Silicone) NuSil MED-4950
LIM Connector – Coating: poly(para-xylylene) Parylene N
Flex Cord Relief: Polydimethyl / methylvinylsiloxane (Silicone) NuSil MED-4950
Flex Cord Relief– Coating: poly(para-xylylene) Parylene N

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








7. DATA ABOUT LABELING OF MD AND ITS PACKAGING

Interpretation of symbols used on the package.

	Catalogue number
	Non-sterile
	European Conformity*
	Refer to Instruction Manual/Booklet for Information Related to Safety
	Manufacturer
	date of manufacture
	Do not use the HARMONIC Blue Hand Piece with Ultracision® Generator (GEN01/GEN32)
	For use with Blue Hand Piece only.
	Temperature Limit

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	Humidity Limitation
	Authorized Representative in the European Community
	Authorized Representative in the United States
	Electrical and electronic equipment. Return waste to a collection system or treatment and recycling facilities. Applicable in the EU. Follow decontamination instructions before returning waste.
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. *
	Unit Quantity
	Serial number
	«PCT» Compliance Symbol
	Indicates that the Electrical and Electronic Product (EEP) contains toxic and hazardous substances or elements over the defined maximum concentration values defined above. The names and content of the toxic and hazardous substances or elements shall be provided in the product instruction manual. The number “10” in the symbol refers to the environment friendly use period; after their environmentfriendly use period expires, they shall be put into the recycling system (People's Republic of China' symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products). *
	Fragile. Caution!
	This way up.
	Caution! Refer to the instructions for use

* - packaging and labelling of medical devices may not contain noted symbols

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Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.
- Audible high-pitched tones during activation are an abnormal condition and an indicator that the blade or Hand Piece is not operating properly. The tones may be an indicator that the Hand Piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high temperatures and user or patient injury.
- To prevent burn injury, discontinue use if the Hand Piece temperature becomes uncomfortable to hold.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by the activation of the HARMONIC instruments in fatty tissue are potentially flammable.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blades should not be in contact with the patient, drapes or flammable materials while not in use.
- During and following activation in tissue, the instrument blades may become hot. Avoid unintended blade contact with tissue, drapes, surgical gowns, or other unintended sites after activation.

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- The Hand Piece meets the international safety standard EN 60601-1 for user contact and is not

intended for patient contact. To prevent burn injury, avoid direct tissue contact with the Hand Piece or take preventative measures to protect tissue that comes in contact with the Hand Piece. Refer to the instrument inserts for applicable duty cycle (proportion of ON/OFF time) information.

- Handle the Hand Piece carefully, as damage may shift resonant frequency.

- Do not bang or drop the Hand Piece.

- Do not clean the Hand Piece electrical connector with alcohol.

- Verify compatibility with generators prior to use. This Hand Piece is not compatible with the UltraCision® Generator (GEN01/GEN32).

- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC System. Use of such products may lead to unanticipated results and possible injury to the user or patient.

- Use this Hand Piece only with a compatible HARMONIC Generator to avoid potential electric shock hazard.

- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

- In case of system failure, ensure the availability of the appropriate back up equipment relevant to the specific procedure.

- Reuse and improper reprocessing or resterilization may cause functional degradation of the device and/or lead to device failure that in turn may result in patient injury or illness.

- To minimize the chance of corrosion, do not use saline for cleaning.

- Follow the instructions and warnings issued by the suppliers of any cleaning agents and equipment used. Avoid exposing the instruments to hypochlorite solutions, as these will promote corrosion. Do not use cleaning agents with chlorine or chloride as they are corrosive.

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- The design of this device permits effective cleaning and reliable sterilization. The user must ensure that cleaning and sterilization are performed in accordance with guidelines, standards and National Health Authority requirements.

Point of Use Treatment (in the operating room)

- Clean device as soon as possible after use.
Remove excessive solids using disposable wipes, enzymatic foam or per hospital procedure. Soiled devices should be separated from unsoiled devices. Soiled Devices should be covered with a towel dampened with tap, sterile or critical water to prevent blood and/or debris from drying. Encrusted soil can increase the time and effort required to complete processing.
 - If the recommendations in the steps above are not possible, the device may be placed in a container and immersed in tap, sterile or critical water and covered for transport to the processing area.
- Caution: To minimize the chance of corrosion, do not use saline for cleaning.
- Follow the Cleaning and Sterilization instructions below.

Containment and Transportation

Surgically used devices may be considered bio-hazardous and should be safely transported to a designated processing area in accordance with hospital procedures and local policies.

Cleaning (after procedure)

Cleaning should be performed as soon as possible after use, preferably within one hour of use. The device has been designed for thorough cleaning and safe sterilization, without disassembly. It is the responsibility of the end user to ensure that the cleaning is performed using appropriate equipment, materials and personnel to achieve the desired result.

Caution: Follow the instructions and warnings issued by the suppliers of any cleaning agents and equipment used. Avoid exposing devices to hypochlorite solutions, as these will promote corrosion. Do not use cleaning agents with chlorine or chloride as they are corrosive.

The user must ensure that cleaning and sterilization are conducted in accordance

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with the appropriate guidelines, standards, or National Health Authority requirements. Hospital grade low foaming neutral pH to mildly alkaline enzymatic detergents, hospital grade low foaming neutral pH detergents (pH 7-9) and hospital grade low foaming mildly alkaline detergents with a pH of ≤ 11 (prepared according to manufacturer's instructions) should be used. An enzymatic detergent and a mildly alkaline detergent were used in cleaning efficacy studies and are types of products that can be used throughout the cleaning process. "Detergent(s)" will be used generically throughout the rest of these instructions. Ensure that detergent manufacturer instructions are followed for concentration, temperature and water quality throughout the cleaning process. When using mildly alkaline detergents, critical water should be used throughout the cleaning process.

The chemical quality of the water used during processing can impact device safety. Facilities should use the recommended water quality requirements for device processing in accordance with local guidance (such as AAMI TIR 34, Water for the reprocessing of medical devices), detergent manufacturers and these instructions for use.

Pre-cleaning Instructions (PRIOR to Automatic or Manual Cleaning)

1. Immerse and soak the device for a minimum of 10 minutes up to 60 minutes in an enzymatic detergent or equivalent prepared according to manufacturer's instructions for concentration, water quality and temperature.
 2. Remove the device from the detergent.
 3. Immerse the device in a clean sink with cold to warm tap water – not to exceed 40°C. Using a toothbrush style soft bristle (nylon) brush or soft cloth, clean the exterior of the device paying particular attention to crevices, threads and seams.
- Caution:** Do not use scouring pads, abrasive cleaners, or a steel wire brush, which may cause damage to the surfaces of the device.
4. Rinse the exterior of the device thoroughly under cold to warm running tap water for a minimum one minute.

Note: The use of ultrasonic cleaners is not recommended for the Hand Piece.

Following the pre-cleaning step, the validated automated or manual processing

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method must be performed. Refer to standard hospital protocol and local/national practices to determine whether disinfection is required. Note that disinfection is mandatory in some areas of the world; e.g. EU. There are different A0 requirements in various countries; e.g, $A0 \geq 600$ (90°C for 1 minute) or $A0 \geq 3000$ (90°C for 5 minutes). The user must ensure that cleaning, disinfection and sterilization are conducted in accordance with the appropriate local guidelines, standards, or National Health Authority requirements.

Automated Cleaning Method (Preferred Method)

Following the pre-cleaning step, the device must be cleaned/disinfected in a validated automated washer-disinfector in compliance to ISO 15883 or equivalent standard, using a typical non-lubricating Instrument, Vario TD Cycle or a program similar to the following cycle, providing the minimum parameters identified below are met:

Phase	Recirculation Time (minutes:seconds)	Temperature	Detergent Type	Validated Parameters
Pre-Wash	02:00	Cold tap water	N/A	2 minutes cold tap water
Enzyme Wash	02:00	Tap water <40°C (104°F)	Neutral pH enzymatic detergent or mildly alkaline detergent	no
Wash	03:00	Tap water $\geq 40^\circ\text{C}$ (104°F)	Neutral pH detergent or mildly alkaline detergent	5 minutes tap water 55-60°C
Rinse	02:00	Critical water $>40^\circ\text{C}$ (104°F)	N/A	2 minutes critical water
*Thermal Disinfection	01:00 – 10:00	Critical water $\geq 90^\circ\text{C}$ (194°F)	N/A	5 minutes at 90°C
Drying	30:00	$\geq 115^\circ\text{C}$ (239°F) – not To exceed 134°C (273°F)	N/A	30 minutes at 115°C

Note: Do not expose the device to a lubrication cycle.

Note: The automated washer program may or may not include and enzymatic wash phase. Single or combined washing phases should be performed for a

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minimum of 5 minutes. Temperatures used during the enzymatic and/or wash phases should meet the requirements of the detergent manufacturer's recommendations.

Note: Manually clean the upper and lower contact rings using a cotton swab dipped in isopropyl alcohol until visibly clean.

Disinfection

Thermal disinfection can be included as part of the automatic cleaning cycle. The device may be thermally disinfected at a temperature between 90°C (194°F) and 100°C (212°F) for one (1) minute up to ten (10) minute(s) to render the device safer for handling.

Manual Cleaning Method

1. Immerse the device for a minimum 5 minutes in a detergent prepared in accordance with manufacturer's recommendations for concentration, water quality and temperature.

2. Remove the device from the detergent.

3. Immerse the device in a clean sink with cold to warm critical water – not to exceed 40°C. While immersed, use a toothbrush style, soft bristle (nylon) brush or soft cloth to clean the exterior of the device. Pay particular attention to crevices, threads and seams.

Caution: Do not use scouring pads, abrasive cleaners, or a steel wire brush, which may cause damage to the surfaces of the device.

4. Rinse the exterior of the device thoroughly with critical water for a minimum one minute.

Drying

Dry the device with a clean, absorbent cloth. If a washer-disinfector without a drying cycle is used, a hot-air dryer may be used but the device must not be heated to temperatures above 134°C (273°F). The cavities of device parts must be blown out with clean, compressed air, up to 20 psi, in order to eliminate the possibility of any residual water remaining.

Note: Manually clean the upper and lower contact rings using a cotton swab dipped in isopropyl alcohol until visibly clean.

Inspection

Devices do not have an indefinite functional life. All surgical devices are subject to a degree of wear and tear as a result of normal use. They must therefore

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undergo a regular and precise visual check before each use:

- Inspect the hand piece blade mount surface for pitting and scratches;
- Inspect the hand piece blade mount surface for cracks, bends, or other significant forms of physical damage;
- Inspect the hand piece threads for physical damage;
- Inspect the blade mount surface for excessive side play;
- Inspect the hand piece body for press fit joint separation;
- Inspect the hand piece body for cracks;
- Inspect the cable for damage;

Ensure the complete removal of soil from all Hand Piece and test tip surfaces. If soil or blood is present, clean the device following the Pre-Cleaning and Automated or Manual Cleaning Instructions.

Discard and replace the instrument if damage or degradation is present. Defective instruments (malfunctioning or difficult to operate instruments and those with stress, corrosion, cracking or pitting, etc.) increase the risk to the patient and must therefore not be used.

Sterilization

Steam (moist heat) sterilization is recommended in an approved, pre-vacuum cycle listed below. The steam sterilizer must be FDA cleared, if performed in the United States, and validated to the requirements of any local standards and guidance such as EN285, ISO17665 or AAMI/ANSI ST8.

Prepare Device for Sterilization

1. Place and wrap device in an approved sterile barrier wrap material per hospital protocol and in accordance with ISO11607-1 and any other local and/or National requirements using standardized wrapping techniques. Or place device in a sterilization tray wrapped using an approved sterile barrier material and standardized wrapping techniques.
2. Place wrapped device or wrapped tray into autoclave.
3. Run the device through the autoclave using one of the settings below.

The following steam sterilization cycle parameters are approved for use.

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Sterilizer Type	Method	Cycle Time (minimum time at Temperature)	Temperature Set Point	Preconditioning Pulses (minimum)	Drying Time
Pre-vacuum	Wrapped	3 minutes	134°C (273°F)	3	15 minutes

Sterilization in US or other applicable countries

Sterilizer Type	Method	Cycle Time (minimum time at Temperature)	Temperature Set Point	Preconditioning Pulses	Minimum Drying Time
Pre-vacuum	Wrapped	4 minutes	132°C (270°F)	3	15 minutes
Pre-vacuum	Wrapped	3 minutes	134°C (273°F)	3	15 minutes
Gravity	Wrapped	15 minutes	132°C (270°F)	Not applicable	15 minutes
Gravity	Wrapped	30 minutes	121°C (250°F)	Not applicable	15 minutes

Notes:

- The above table includes the minimum temperature and cycle time validated to assure sterility.
- Based on steam autoclave tolerances, the actual autoclave temperature can exceed the set point temperature by a maximum of +3°C (+5°F).
- Health Authorities in some regulated regions do not accept unwrapped sterilization methods; e.g. in the EU. Please review the appropriate guidelines, standards and National Health Authority guidelines when determining acceptable steam sterilization process parameters for use in each respective country.
- Immediate-use steam sterilization is only intended for individual devices and should only be performed when approved by local policies. Immediate-use steam sterilization of devices is not supported. It is the sole responsibility of the user to validate immediate-use steam sterilization if performed.

STERRAD® Systems

- Clean and thoroughly dry the HARMONIC Hand Piece and Test Tip following the steps listed in the Hand Piece and Test Tip Cleaning section.
- Using standardized wrapping techniques, place and wrap the device in an approved STERRAD® System-compatible polypropylene sterilization wrap in

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