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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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Single Use Sterile Flexible Ureterorenoscope, model US31B-12

Operation Manual

CE 0197

Version: A/6


Effective Date: 31 august 2020 r.

Signature: _____

Position: CEO

Shanghai AnQing Medical Instrument Co., Ltd.

Flexible Ureterorenoscope

Operation Manual

Only to be used with Video Processor (Model: EOS-B-01)

approved by SHANGHAI ANQING

CE 0197

Version: A/6

Effective date: August 31, 2020

Shanghai AnQing Medical Instrument Co., Ltd

Contents

User Notice.....	3
Flexible Ureterorenoscope.....	4
Device Description.....	4
Intended Use.....	8
Applicability of Endoscopy and Endoscopic Treatment	9
Safety Information.....	9
How Supplied	12
Handling and Storage.....	12
Operational Instructions	12
Disposal.....	17
EMC Information	17

User Notice

Thank you for your expression of confidence in the SHANGHAI ANQING brand name. Flexible Ureterorenoscope is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from SHANGHAI ANQING.

The SHANGHAI ANQING Flexible Ureterorenoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the technique of Ureterorenoscope.

We ask that you read this manual carefully before proceeding to work with this product. Keep this manual available for reference when needed.

Please follow the instructions exactly, and make proper care and maintenance.

A WARNING indicates that the personal safety of the patient or physician may be involved.

Disregarding a Warning could result in injury to the patient or physician.

Please contact your local distributors or SHANGHAI ANQING directly if there is a problem with the product. We will provide you with excellence service and help.

We are committed to providing users with more detailed technical information.

flexible Ureterorenoscope

Only to be used with the designed Video Processor approved by SHANGHAI ANQING)

Devices supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call SHANGHAI ANQING.


The product is single use. Do not re-use, reprocess or re-sterilize the product. Re-use, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose the product and packaging in accordance with hospital, administrative and/or local government policy.

Device Description

flexible Ureterorenoscope is a flexible inserted portion intended for visual examination and treatment of the ureter (the upper urinary tract that connects the kidney to the bladder). It is inserted through the external urethral orifice and bladder during ureteroscopy. Anatomical images are transmitted to the user by a video processor produced by SHANGHAI ANQING Company with a CMOS chip at the distal end of the endoscope and the images showing on a monitor. This device is commonly used for biopsies, endoscopic lithotripsy, and treatment of upper urinary tract malignancies, incising strictures, and repairing ureteropelvic junction obstructions. This is sterile and single use disposable device with human body contact duration about 45min -2h.

flexible Ureterorenoscope is equipped with a sterile disposable plastic cable sheath as an accessory to protect cable from contamination.

1. Instrument Channel Outlet
2. LED Cold Light
3. Camera Lens
4. Distal End
5. Bending Section
6. Insertion Tube
7. Insertion Section
8. Irrigation Valve
9. Irrigation Luer Taper
10. Product Specification
11. Angulation Control Knob
12. Operating Handle
13. Instrument Channel
14. Instrument Channel Valve
15. Angulation Control Mark
16. Signal Port

Components and Functions

1. **Instrument Channel Outlet**
Endo therapy accessories extend from the instrument channel outlet.
2. **LED Cold Light**
Illuminate a surgical site.
3. **Camera Lens**
Collect endoscopic image, and transmit it to the user by a video processor produced by SHANGHAI ANQING and the images shows on a monitor.
4. **Distal End**
Distal End, integrated with LED cold light source and camera system.
5. **Bending Section**
The section moves the distal end of the endoscope when the UP/ DOWN angulation control knobs are operated.
6. **Insertion Tube**
Insertion tube may be inserted into the patient's body.
7. **Insertion Section**
It is inserted into patients' body with or without sheath to complete check and operation.
8. **Irrigation Valve**
Used to open and close irrigation valve.
9. **Irrigation Luer Taper**
Catheter is connected by luer taper to supply water.
10. **Product Specification**
US31A-12 or US31B-12
11. **Angulation Control Knob**
When this knob is turned in the "U" direction, the bending section moves UP; when the knob is turned in the "D" direction, the bending section moves DOWN; when the knob is turned in the "0" direction, the bending section is kept STRAIGHT.
12. **Operating Handle**
Synchronous movement between scope tip and operating levers on the handle, zero delay from the surgeon's finger to the scope tip.
13. **Instrument Channel**

The instrument channel functions as channel for the insertion of endo-therapy accessories.

4. Instrument Channel Valve

Rotate instrument channel valve to fix endo-therapy accessories.

15. Angulation Control Mark

D for downward, U for upward and 0 for straight.

16. Signal Port

Connect Signal cord to Signal port to supply power and transmit data.

Specification

Table 1 Flexible Ureterorenoscope Specification

Model	Distal End Outer diameter	Instrument Channel Inner diameter	Working length	Total length	Bending Section Angulation range	Resolution (pixel)
US31A-12	3.1mm (9.3Fr)	1.2mm (3.6 Fr)	650 mm	920 mm	275°up/275°down	160K
US31B-12						

When this knob is turned in the "U" direction, the bending section moves UP; when the knob is turned in the "D" direction, the bending section moves DOWN; when the knob is turned in the "0" direction, the bending section is kept STRAIGHT.

Table 2 Operating and Storage Environment

Working Environment	Ambient Temperature	10°C~40°C
	Relative Humidity	30%~80%
	Atmospheric Pressure	860hPa~1060hPa
	Supply Voltage	AC 100V~240V
	Frequency	50Hz~60Hz
	Ambient Condition	The device should be far away from the interference of strong magnetic field environment.
Storage and Transportation Environment	Ambient Temperature	-10°C ~+55°C
	Relative Humidity	10%~80%
	Atmospheric Pressure	500hPa~1060hPa

Table 3 Performance Index

Performance Index		
Camera System	Image Sensor	Color CMOS
	Resolution	400(H)x400(v)
	Clear Observation Range	3mm~100mm
	Field of view	110°±5%
	Direction of View	0°±5°
	Resolution power/ (lp/mm)	≥7.73lp/mm
	Light/lux	> 5000lux
Sterilization	Sterilization Method	EO
	Sterilization Validity	3 years
Classification regarding to protecting against electric shock	Electric Shock	II
	Protect Against Electric Shock	BF

The Flexible Ureterorenoscope should be used in combination with Video Processor manufactured by SHANGHAI ANQING. Video Processor provides light source and image transmitted for Flexible Ureterorenoscope.

Intended Use

Flexible Ureterorenoscope is inserted through the external urethral orifice and bladder and used for diagnosis and treatment of ureter and renal pelvic diseases.

Contraindications

Contraindications for this device are those specific to urinary tract endoscopy.

- Diagnostic or therapeutic ureterorenoscopy is contraindicated in people with an untreated urinary tract infection.
- Other contraindications to therapeutic ureterorenoscopy (e.g., lithotripsy, endopyelotomy, tumor therapy) are more numerous and can mirror those associated with the corresponding open surgical interventions. Patients on anticoagulants or with coagulopathies should be managed appropriately.

Potential Adverse Events

Possible complications include, but may not be limited to:

- Bleeding • Fever
- Avulsion • Sepsis
- Stenosis / Stricture • Renal failure
- Inflammation • Perforation (ureter, renal pelvis or bladder)
- Laceration • Hematuria

• Urteral Reflux
• Hematoma
• Urethral damage
• Infection

Applicability of Endoscopy and Endoscopic Treatment

There are official standards on the applicability of endoscopy and endoscopic treatment that are defined by the hospital's administrations or other official institutions, such as academic societies on endoscopy, follow those standards. Before starting endoscopy and endoscopic treatment, thoroughly evaluate its properties, purposes, effects, and possible risks (their nature, extent, and probability). Perform endoscopy and endoscopic treatment only when its potential benefits are greater than its risks.

Fully explain to the patient the potential benefits and risks of the endoscopy and endoscopic treatment as well as any examination/treatment methods that can be performed in its place, and perform the endoscopy and endoscopic treatment only after obtaining the consent of the patient.

Even after starting the endoscopy and endoscopic treatment, continue to evaluate the potential benefits and risks, and immediately stop the endoscopy/treatment and take proper measures if the risks to the patient become greater than the potential benefits.

Safety Information

Please read these safety instructions carefully before using the instrument. Before using the instrument on the patient, it is imperative that you be acquainted with how the instrument operates and is controlled.

The Flexible Ureterorenoscope must be used according to the recognized medical rules and procedures of endoscopy.

Warning

- Read the following safety instructions carefully before using the Flexible Ureterorenoscope to avoid putting your patients, personnel or yourself at risk.
- Inspect the Flexible Ureterorenoscope prior to each surgical procedure to ensure that it functions correctly and that its sterile package integrity.
- Every time before the Flexible Ureterorenoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged Flexible Ureterorenoscope or damaged accessories must not be used.
- Combinations of medical devices are only then assured to be safe if
 - a) They are identified as such in the respective instruction manual or
 - b) The intended purpose and the interface specifications of the devices used in combination permit this.

Pay careful attention to the instructions and interface specifications of medical devices used in combination.

Only use accessories (guide probes, cleaning brush, etc.) that compatible with this Flexible Ureterorenoscope and can be operated in line with generally accepted technical standards and occupational safety guidelines.

- Do not allow any electrical components or electrical instruments to come into contact with the Flexible Ureterorenoscope.
- The Flexible Ureterorenoscope must not be used during the discharge of a cardiac defibrillator. Remove the Flexible Ureterorenoscope from the operation area before discharge takes place.
- Danger of explosion if flammable gases are used in the immediate proximity of the Flexible Ureterorenoscope. It must not be used in the presence of flammable or combustible media.
- The temperature of the distal end of the endoscope may exceed 41°C due to intense endoscopic illumination. Surface temperatures over 41°C may cause mucosal burns. Always maintain a suitable distance necessary for adequate viewing while using the minimum level of illumination for the minimum amount of time. Do not use close stationary viewing or leave the distal end of the endoscope close to the mucous membrane for a long time without necessity.
- If a malfunction should occur during a procedure on a patient, stop the procedure immediately. Place the distal tip of the Flexible Ureterorenoscope in a neutral, nondeflected position and remove the video endoscope slowly and carefully from the patient.
- The device is not intended to be used within the MR environment.
- When the camera system is used to examine a patient, do not allow metal parts of the endoscope or its accessories to touch metal parts of other system components. Such contact may cause unintended current flow to the patient.
- Always set the minimum required brightness.
- Do not continue observation in the proximity to tissue or keep the distal end of the endoscope in contact with living tissue for a long time. It may cause patient burns.
- When discontinuing the use of the endoscope, be sure to turn the light source OFF so that the endoscope does not irradiate unnecessary light.
- If the endoscopic image dims during use, blood, mucus or debris may adhere to the light guide on the distal end of the endoscope. Carefully withdraw the endoscope from the patient and remove the blood or mucus to obtain optimum illumination and to ensure the safety of the examination. If you continue to use the endoscope in such a condition, the distal end temperature may rise and cause mucosal burns. It may also cause patient and/or operator injury.
- Do not touch the electrical contacts inside the device connectors.
- Do not apply excessive force to the device and/or other instruments connected. Otherwise, damage and/or malfunction can occur.
- Do not connect or disconnect the device connector while the camera system is turned ON. Connecting or disconnecting the device while this camera system is ON may destroy the CMOS. Turn the camera system OFF before connecting or disconnecting the device.
- In the event that the live endoscopic image is lost, do not advance or insert the Flexible Ureterorenoscope. Doing so can cause patient injury such as perforation, avulsion, hemorrhage or urothelial damage.
- Do not use excessive force while advancing or withdrawing the scope. Doing so can cause patient injury such as perforation, hemorrhage or urothelial damage or damage to the

ureterorenoscope.

Do not use excessive force when advancing or withdrawing an accessory within the ureterorenoscope. Doing so can cause patient injury such as perforation, avulsion, hemorrhage, urothelial damage or damage to the Flexible Ureterorenoscope.

- Do not withdraw a laser fiber back into the Flexible Ureterorenoscope while the laser is firing. Doing so may cause patient injury and/or scope damage.
- Do not open the handle of the Flexible Ureterorenoscope. Doing so can damage the waterproof seals and result in risk of electric shock.
- Do not excessively bend the flexible shaft or the articulating section of the Flexible Ureterorenoscope as this may break or kink the shaft.
- If damage to the Flexible Ureterorenoscope occurs or it stops functioning during a procedure, stop using the Flexible Ureterorenoscope immediately. Continue the procedure with a new Flexible Ureterorenoscope, as appropriate.
- Do not insert a wet, contaminated, or damaged connector cable plug into the Flexible Ureterorenoscope as poor video performance or damage to the scope may result.

Cautions

- Do not expose the Flexible Ureterorenoscope to direct sunlight or excessive heat, to avoid damage to the device.
- Do not steam sterilize (autoclave) the Flexible Ureterorenoscope.
- The maximum permissible temperature for decontamination, sterilization, transport and storage is 55°C.
- Do not hit the distal tip of the Flexible Ureterorenoscope against hard objects, and never use a sharp object to remove dirt.
- Never lay heavy objects on the Flexible Ureterorenoscope.
- Flexible Ureterorenoscope should not be used in combination with video processor from other manufacturers. The use of non-compatible units may cause injury to the patient or damage to the instrument. It must be ensured that parts of the endoscope that are not electrically insulated do not come into contact with conductive surfaces or live parts of other unit, at any time.
- The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases, actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.
Do not use force to turn the deflection lever, to avoid either injuring the patient or damaging the instrument
- Use accessories approved by the manufacturer only. Do not connect electronic devices other than approved by the manufacturer to the USB port.
- Be sure that this instrument is not used adjacent to or stacked with other equipment (other than the components of this instrument or system) to avoid electromagnetic interference.

Supplied

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. Do not use if the device is damaged.

Handling and Storage

Store the Flexible Ureterorenoscope in a cool, dry, dark place.

Storage/transport conditions: Temperature: -10°C ~55°C; Humidity: 10%~80%

Caution:

Do not expose the Flexible Ureterorenoscope to direct sunlight or excessive heat.

Caution:

The distal tip of the Flexible Ureterorenoscope must be in the neutral position during storage.

Caution:

The sheaths of the Flexible Ureterorenoscope must never be kinked or coiled up too tightly. It is best to suspend the Flexible Ureterorenoscope for storage or to keep it in a large container.

Combination with Other Devices

Only to be used with Video Processor (Model: EOS-B-01) approved by SHANGHAI ANQING.

Flexible sheath has a 3.6 Fr ID working channel.

Operational Instructions

Unpacking the device

Check for missing items and evidence of shipping damage.

File any complaints immediately with SHANGHAI ANQING or your local distributors.

Operating conditions

Ambient temperature: 10°C ~40°C

Humidity: 30%~80%

Electromagnetic compatibility

This instrument has been tested and complies with the EMC limits as per EN/ IEC 60601-1-2 CISPR 11 Class B.

Handling

Note: Damage to the Flexible Ureterorenoscope resulting from incorrect operation is not covered by

manufacturing's warranty.

The following instructions outline the correct operation of the Flexible Ureterorenoscope and the corresponding accessories. This is not an introduction to the techniques of ureteroscopy. Instruction in endoscopic technique can be found in the relevant medical literature.

Handling the Flexible Ureterorenoscope correctly

Grasp the handle of the Flexible Ureterorenoscope in your left or right hand, as shown in the picture. With this kind of control grip, the Flexible Ureterorenoscope can be operated with one hand.

Handling the sheath

 **Caution:** Do not carry the Flexible Ureterorenoscope by the sheath. Do not pull, clamp or twist the sheath.

 **Caution:** The sheaths of the Flexible Ureterorenoscope must never be kinked or coiled up too tightly. It is best to suspend the Flexible Ureterorenoscope for storage.

Warning:

Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

If you lose your orientation, pull the sheath back slightly in the neutral position until vision is restored.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, turn the UP/ DOWN Angulation Control Knobs in the "0" direction, and slowly pull the sheath of the Flexible Ureterorenoscope out of the patient's body.

Handling the Angulation Control Knob

Note:

The neutral "0" position of the distal tip is straight, without any deflection whatsoever.

Warning:

Suitable accessories, such as laser fibers, must only be introduced into the instrument channel of the Flexible Ureterorenoscope with the distal tip in the neutral, non-deflected "0" position.

Once the instrument has been introduced, actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

Warning:

Before removing the Flexible Ureterorenoscope from the patient, the distal tip of the Flexible Ureterorenoscope must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

no longer be properly controlled (sticking etc.) via the deflection lever, or if the instrument develops any other fault, the Flexible Ureterorenoscope must not be used.

If a functional defect occurs during use on the patient, discontinue use immediately. Put distal tip of the Flexible Ureterorenoscope into the neutral, non-deflected "0" position and remove the Flexible Ureterorenoscope slowly and carefully from the patient.

Caution:

The resistance that can be felt when the deflection lever deflects the tip increases as the deflection angle increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.



Caution:

Do not use force to turn the deflection lever, to avoid either injuring the patient or damaging the instrument.

The distal tip can be deflected 275° upward or downward by continuously actuating the deflection lever slowly and carefully.

Before use, the user must always have full control over the deflection and must develop a sure feel for the deflection mechanism and in particular for the neutral position.

Mounting plastic cable sheath

1. Plastic cable sheath ring
Buckle the shell of power port.
2. Plastic cable sheath
Protect the cable from contamination.
3. Plastic cable sheath box
Storage the plastic sheath.

Step 1 Buckle sheath ring to shell of signal port on the Ureterorenoscope.

Step 2 Pass signal cable through plastic cable sheath and connect to signal port on the

Ureterorenoscope.

Pull sheath box along the signal cable to cover the cable.

Inspection prior to use

Before using, inspect the package for any damage to ensure packaging of the sterile product intact. If the package is opened or damaged.

Warning:

First check the Flexible Ureterorenoscope and any accessories used in combination with it, are in good condition. Damaged Flexible Ureterorenoscope or damaged accessories must not be used.

Inspecting the surfaces

Before Flexible Ureterorenoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Flexible Ureterorenoscope with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

Inspect all the supply cables for damage (such as breaks, cracks, twisted or squashed sections...).

Warning:

The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

Caution:

Do not squeeze the distal tip of the Flexible Ureterorenoscope too tightly.

Checking the deflection

The distal tip is deflected by means of the deflection lever. Turn the deflection lever to "D" the distal tip should move downward. Turn the deflection lever in the opposite direction to produce an upward movement.

Warning:

Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device!

Move the deflection lever smoothly and slowly.


Warning:

Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective Flexible Ureterorenoscope must not be used.

g the LED function

video processor as described above and connect to the Flexible Ureterorenoscope. Emission from the distal end should be recognizable. Light intensity is manually regulated by the Video

Inspecting the Video Processor

Prepare the video processor as described above and connect to the Flexible Ureterorenoscope. Check whether the image is clear.

Warning:

If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the Flexible Ureterorenoscope from the patient's body gently and carefully. If the Flexible Ureterorenoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.





Caution:

The video processor must be switched off before connecting the Flexible Ureterorenoscope. Before each use or after a change of viewing modes/setting, check to ensure the view observed through the endoscope provides a live image and has the correct image orientation.

Inspecting the instruments and the instrument channel

Check the instrument channel by introducing the appropriate instrument (max. diameter of instruments: 3.6Fr.)

Caution:

Accessories, which are bent or twisted may no longer function properly, are difficult to pass through the channel, and may cause serious damage to the channel. They must therefore be replaced!

Warning:

Unsterilized instruments pose a potential risk of infection on arrival at the recipient's premises.

Disposal

Flexible Ureterorenoscope does not contain latex and phthalate.

After use, dispose the plastic cable sheath and Flexible Ureterorenoscope according to national laws and regulations.

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...cables other than those specified or provided by the MANUFACTURER of ... result in increased electromagnetic emissions or decreased electromagnetic ... and result in improper operation.

...communications equipment (including peripherals such as antenna cables and ... should be used no closer than 30 cm (12 inches) to any part of this model, ... specified by the MANUFACTURER. ... degradation of the performance of this equipment could result.

... and manufacturer's declaration – Electromagnetic emissions

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manufacturer's declaration – Electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The customer or owner of this model should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
Electromagnetic discharge IEC 61000-4-2	ENCLOSURE Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air PATIENT coupling Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air Signal input/output parts Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Same as left	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material that tends to produce static, the relative humidity should be at least 30%.
Electrical transient/burst IEC 61000-4-4	Input a.c. power Port ± 2 kV 100 kHz repetition frequency Signal input/output parts Port ± 1 kV 100 kHz repetition frequency	Same as left	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input a.c. power Port Line-to-line $\pm 0,5$ kV, ± 1 kV Line-to-ground $\pm 0,5$ kV, ± 1 kV, ± 2 kV Signal input/output parts Port Line-to-ground ± 2 kV	Same as left	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11 Voltage interruptions IEC 61000-4-11	Input a.c. power Port 0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles Single phase: at 0° 0 % U_T ; 250/300 cycle	Same as left	Mains power quality should be that of a typical commercial or hospital environment. If the user of instrument requires continued operation during power mains interruptions, it is recommended that the instrument be powered from an uninterruptible power supply or a battery.

IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
ENCLOSURE Port 30 A/m 50Hz or 60Hz	Same as left	It is recommended to use this instrument by maintaining enough distance from any equipment that operates with high current. If image distortion occurs, it may be necessary to position the Model further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

the AC mains voltage prior to application of the test level.

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




100 V RMS at 1 kHz


... in all situations.

... is affected by absorption and reflection from structures, objects and people.

... RF transmitters as determined by an electromagnetic site survey should be less than the ... in each frequency range.

... may occur in the vicinity of high-frequency electrosurgical equipment and/or other ... with the following symbol:

	Description	Symbol	Description
	Caution		IEC 60601 Type BF shock protection
	Consult instructions for use		EU Authorized Representative
	Do not re-use.		Batch code
	Do not use if package is open or damaged.		Serial number
	Do not re-sterilize		Date of manufacture
	Sterilized using ethylene oxide		Manufacturer
	Keep dry		Use-by date
	Upward		Waste management
	Keep away from sunlight		Temperature limit -10°~55°
	Humidity limitation 10%~80%		Atmospheric pressure limitation 500hPa~1060hPa
	Fragile		CE mark

**EU Authorized
Representative**

Emergo Europe
Prinsessegracht 20, 2514 AP The Hague,
The Netherlands
Tel: +31.70.345.8570

**Legal
Manufacturer**

Shanghai AnQing Medical Instrument Co., Ltd.
Room 202, 2nd Floor, No. 7 Building
No. 150 Cailun Road,
Zhangjiang High-Tech Park
Shanghai, China

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

2018-09-29 2021-09-28 QTY 1

(01)06970034210130
(17)210928
(10)180909US

INNOVEX
ANQING MEDICAL
Flexible
Ureterorenoscope
REF US31B-12
LOT 180909US
2021-09-28

INNOVEX
ANQING MEDICAL
Flexible
Ureterorenoscope
REF US31B-12
LOT 180909US
2021-09-28

INNOVEX
ANQING MEDICAL
Flexible
Ureterorenoscope
REF US31B-12
LOT 180909US
2021-09-28

INNOVEX
ANQING MEDICAL
Flexible
Ureterorenoscope
REF US31B-12
LOT 180909US
2021-09-28

Shanghai AnQing Medical Instrument Co., Ltd.
Room 202, 2nd Floor, No.7 Building, 150 Cailun Road,
Zhangjiang High-Tech Park, 201210 Shanghai, CHINA

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THIS IS TO CERTIFY THAT: the signature of Yan Hang, the CEO of Shanghai Anqing Medical Instrument Co., Ltd. and the seal of the said company on the annexed DOCUMENT are genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhou Pingfan

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(Date: Mar. 25, 2021)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Video Processor EOS-B-01

Operation Manual

CE 0197

Version: A/7

Effective Date: 31 august 2020 r.

Signature: _____

Position: CEO

Shanghai AnQing Medical Instrument Co., Ltd.

Video Processor (Model: EOS-B-01)

Operation Manual

CE 0197

Version: A/7

Effective date: August 31, 2020

Shanghai AnQing Medical Instrument Co., Ltd

Contents

1	User Notice	3
2	Safety Information	4
3	Device Purpose	8
4	Initial Device Startup	11
5	Operating the Device.....	14
6	Care/ Maintenance/Repair.....	19
7	Trouble Shooting.....	21
8	EMC Information	22
9	Disposal	25

1 User Notice

- Thank you for your expression of confidence in the SHANGHAI ANQING brand name. Video Processor (Model: EOS-B-01) is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from SHANGHAI ANQING.
- The Video Processor (Model: EOS-B-01) may only be used by persons with an appropriate medical qualification and who are acquainted with the technique of endoscopy.
- We ask that you read this manual carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead
 - To life-threatening injuries of the patient,
 - To severe injuries of the surgical team, nursing staff or service personnel, or
 - To damage or malfunction of device and / or accessories.
- Please follow the instructions exactly, and make proper care and maintenance.
- A WARNING indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.
- Please contact your local distributors or SHANGHAI ANQING directly if there is a problem with the product. We will provide you with excellence service and help.
- We are committed to providing users with more detailed technical information.
- The manufacturer reserves the right to modify the appearance, graphics, and technical data of the supplied product through continued product development.

2 Safety Information

Exclusion of liability

The manufacturer is not liable for direct or consequential damage and the warranty is null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the manual are not adhered to;
- non-authorized persons perform repairs, adjustments, or alterations on or to the device or accessories;
- non-authorized persons open the device;
- the prescribed inspection and maintenance schedules are not adhered to.

Receipt of technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on or to the device or accessories.

Authorized service technician

Only an authorized service technician may perform repairs, adjustments, or alterations on the device or accessories and use the service menu. Any violation will void the manufacturer's warranty. Authorized service technicians are only trained and certified by the manufacturer.

Intended Use

The Video Processor is used with the specified endoscope designed by ANQING during minimally invasive surgery. The Video Processor provides power and processes the images for medical electronic endoscope.

User qualification

The device is only intended to be used by trained physicians who are familiar with endoscopic operations.

Care and maintenance


The service and maintenance of the device and its accessories have to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check the integrity and function of the device before each use.

Waste management

This symbol indicates that the waste of electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately instead. The user must discard and dispose the product by local laws and hospital management specifications.


Hazard

Warning!

Technique and procedures

Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

Warning!

Check all factory settings.


Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

Warning!

Original accessories

For your own safety and that of your patient, use only original accessories.

Warning!

Not explosion-proof

Never install and operate the video system center in locations where:

- the concentration of oxygen is high;
- oxidizing agents (such as nitrous oxide (N₂O)) are present in the atmosphere;
- flammable gases are present in the atmosphere;
- flammable liquids are near.

Otherwise, explosion or fire may result because this video system center is not explosion-proof.

Warning!

Risk of electrical shock

To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

Warning!

Professional qualification

This manual does not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualification working under the direction and supervision of a physician.

Warning!

Function test

The function test must be performed prior to each surgery.

Warning!

Sterile mediums and accessories

Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

Warning!

Replacement device and accessories

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within easy reach to be able to finish the operation with the replacement components.

Warning!

Cleaning the device

Do not sterilize the device.

Warning!

Condensation / Water penetration

Protect device from moisture. Do not use if moisture has penetrated the device.

Warning!

Device-inherent dangers

Read the warnings specific to this device in chapter-3.4 "Device-inherent Dangers and Information".

Warning!

Device defect

If a device defect is suspected or confirmed, do not use it. Make sure the device can no longer be used until a qualified service technician conducts the appropriate tests and repairs.

 Warning!

Powered accessory

The residual current flowing through the patient could increase when using endoscopes with electrically powered accessories.

 Warning!

Obvious defects

Never use the device if it has obvious defects, especially if these involve the power plugs or the mains power supply connection cables. In this case have the device repaired by authorized service personnel.

 Caution!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

 Caution!

Endoscope

The device may only be connected with specified endoscope designed by SHANGHAI ANQING.

 Caution!

Electrical Interference

Chapter 8 "Guidelines and Manufacturer's Statement – Electromagnetic Compatibility")

Electrical interference with other devices or instruments was practically eliminated when developing this device and none was detected during testing. However, if you still detect or suspect such interference, please follow these suggestions:

Move this, the other or both devices to a different location;

Increase distance between used devices;

Consult an electro-medical expert.

 Caution!

Leakage current

If more power consuming devices are connected simultaneously to one socket by means of distribution boxes the sum of the individual leakage currents may exceed the tolerated limit values.

Device Purpose

3.1 Structure

The Video Processor is made by the control unit, foot switch, power adaptor, data transmission cable, USB connect cable, HDMI-DVI connector, HDMI cable and AV cable. The control unit is mainly composed of the buttons and main board's components.

3.2 Accessories list

No.	Accessories	Unit	Qty	Remark
1	Video Processor	ea	1	Standard config
2	Foot switch	ea	1	Standard config
3	Power adaptor	ea	1	Standard config
4	Data transmission cable	ea	1	Standard config
5	USB connect cable	ea	1	Standard config
6	HDMI-DVI connector	ea	1	Standard config
7	HDMI cable	ea	1	Standard config
8	AV cable	ea	1	Standard config

Warning!

Components and parts of Video Processor have no stored energy, but please collected and disposed them separately according to local regulations for disposal of electronic devices, or which maybe result to unacceptable risk disassemblers or others.

Warning!

Components and parts of Video Processor no radioactive sources and induced radioactive, there is no related safe risk during life management of Video Processor.

3.3 Specification

Specification		
Function	White balance	Adjust white balance manually
	Brightness	With 20%,40%,60%,80%,100% level

	Frozen image	With frozen image function
	Photo/Video	Take photo and video, store and play back
	Switch view frame	With round, octagonal and panoramic frames
	Display zoom	Zoom in and out, 1X~4X
	Memory	With 32G SD card
	Foot switch control	Take photo and frozen image by foot switch
Working condition	Temperature	10°C ~ 40°C
	Humidity	30% ~ 80%
	Air Pressure	860hPa ~ 1060hPa
	Power supply	AC 100V~240V
	Frequency	50Hz ~ 60Hz
	Environment condition	Away from the interference of magnetic field
Storage/transport condition	Temperature	-10°C ~ +55°C
	Humidity	10% ~ 80%
	Air Pressure	500hPa ~ 1060hPa
Product Life	5 years	
Electric shock	Class II	
Protect electric shock	BF	

3.4 Device-inherent Danger and Information

Warning!

ON/OFF switch

The device is only completely disconnected from the mains power supply if the power plug is unplugged from the shockproof safety socket.

Warning!

Installation

Electrical connections of the operating room where the equipment is used must comply with the corresponding national requirements.

Warning!

Danger of burns

Contact with the open end of an endoscope or a light cable connected to the light source can result in burns due to the high energy of the emitted light. Avoid long-lasting, motionless contact between the distal tip of the endoscope and the patient's tissue during endoscopic surgery. To avoid danger of burns, do not leave an active light-light cable, with or without a connected endoscope, on or near the patient.

Warning!

Danger of ignition

Always comply with national regulations concerning avoiding the danger of ignition caused by electrostatic charges.

Warning!

Device type BF

Degree of protection against electric shock. The endoscope must not have an electrical conductive connection to the casing of the Video Processor.

Caution!

Ambient temperature

If the Video Processor is operated within an enclosed device tower or rack system, the operating temperature within this tower may be higher than the ambient room temperature. Therefore, always install the equipment in an environment compatible with the manufacturer's rated ambient temperature.

Caution!

Peripheral devices

Additional peripheral equipment connected to interfaces of the medical monitor has to meet the requirements of the following specifications: EN 60601-2-18 for endoscopic devices and EN 60601-1 for electrical medical devices. All configurations have to comply with EN 60601-1-1 specifications. Whoever connects additional equipment to signal output or signal input is obliged to meet requirements of the standard EN 60601-1-1.

 Caution!

Device interfaces

The device operator may never touch the interfaces of the device and the patient at the same time.

 Caution!

Mobile telephones

The interference caused by electromagnetic waves (e.g. mobile telephones) affecting electronic devices is well known. These should be avoided.

 Caution!

Endoscopes and accessories

The outer surface of the endoscope and each endoscopically usable piece of additional equipment has to be checked to ensure that the area is free of rough spots and/or sharp edges that could harm the patient.

 Caution!

Keep fluids away from all electrical equipment. If fluids are spilled on or into the unit, stop operation of the video system center immediately and contact the manufacturer.

 Caution!

Do not prepare, inspect or use this video system with wet hands.

 Caution!

Use accessories approved by the manufacturer only, or it will interfere with other medical electronic equipment used in combination.

4 Initial Device Startup

Delivery Inspection

Always check all parts and accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

Scope of Delivery

Check the delivered equipment for completeness. Compare delivered parts with the enclosed packing list.

Return the Device

If it becomes necessary to return the device, use of the original packaging is required. The manufacturer does not take responsibility for damage that has occurred during transportation if the damage was caused by inadequate transport packaging.

Please make sure that all required information has been supplied:

- Name of owner
- Address of owner
- Device type
- Serial number of the equipment (see identification plate)
- Detailed description of defect

4.2 Preparing the Device

It is the responsibility of the operator to ensure the device is safe and functions as intended when being used.

Setting up the Device

Place the device on a level surface and install in a dry environment. The ambient temperature and humidity must meet the requirements mentioned.

Caution!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device.

Incorrect voltage can cause Errors and malfunctions and may destroy the device.

While using the device the patient must be treated and kept under observation with the usual medical care.

Establish sterile conditions, provided they are required.

Caution!

Do not use in an area where heating sources or an opening or vent for air conditioning and ventilation are present. Do not expose the unit to direct sunlight, excessive dust, vibrations, or mechanical shocks.

Warning!

Do not pull on the data transmission cable.

Never squeeze, compress, bend, or otherwise twist or spindle the data transmission cable. This can damage the wires inside of the cable and cause image loss or cable failure.

4.3 Connecting the Equipment

Caution!

The manufacturer is not liable for direct or consequential damages and the warranty becomes null and void due to improper or incorrect Video Processor assembly.

Caution!

Ensure that the data transmission cable plug is correctly inserted into the socket. Make sure the red dot mark is in the correct position of the Video Processor.

Caution!

The Video Processor may be plugged in or removed only if the Video Processor is switched off (no hot-plugging). If the Video Processor plug is connected while the Video Processor is switched on, the Video Processor must be switched off and on again.

4.3.1 Connecting Monitor and Endoscope

Use cable to connect Video Processor

- with a monitor;
- a specified endoscope;
- and any additional devices if required.

4.3.2 Connecting the Main Connection Cable

Caution!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

Warning!

ON/OFF switch

The device is only completely disconnected from the mains power supply if the power plug is unplugged from the shockproof safety socket.

4.3.3 Disconnect device from power supply

Disconnect the device from the mains power supply (pull cable out off the Disconnect device from power supply grounded safety wall socket) if the device is not being used for several days or longer.

Warning!

Always grasp the power plug when disconnecting the device from the power supply. Never pull on the cable itself.


4.3.4 Performing White Balance


A white balance has to be performed at the beginning of each surgical procedure using light source and endoscope. Hold a white area (e.g. a white sheet of paper) in front of the optical system to fill the screen completely with the white area.

- Press White Balance key.
- The white balance (WB OK) is displayed in the bottom right area of the connection monitor.
- After the auto-completion, the Video Processor is ready for use.

5 Operating the Device

System configuration


Interface

- ① Start-up button
- ② Foot switch interface
- ③ SD card interface
- ④ Data transmission line interface
- ⑤ TV interface
- ⑥ HDMI interface
- ⑦ HDMI interface
- ⑧ USB interface
- ⑨ Power adapter interface


ating

- Connect endoscope to Video Processor.
- Connect image processor to power adapter via ⑨ Power adapter interface.
- Plug in data transmission line and connect to the endoscope. The red dot on the plug of data transmission line is pointed at the red dot position of ④ Data transmission line interface.


- Connect to the monitor via ⑥ HDMI connector, ⑦ HDMI connector or ⑤ TV connector.
- The operator holds the endoscope with his left hand. Press 1-Start-up button of the image processor to turn the device on, and the switch button lights up in green.
- Users can make adjustments according to the Menu -User interface to facilitate use, and can also make relevant adjustments during use.

Menu-User Interface


Power On


Power up, detect Date/ Time.
If the time has not been saved back to the initial time, enter the time setting screen.
(The left and right keys are date/time options, the up and down keys are numeric settings; the Enter key is set to complete.)
If the power-on detection time does not return to the initial time, the time setting will not be skipped.

Enter folder name setting, input the patient's name in English letters or numbers.
After the name is set, the patient's information will be stored in this folder.
The folder name will be displayed on the UI interface (see Figure).

After the time and folder settings are completed, enter the Preview mode, as shown in the figure.

Take Picture

Press camera key to take a photo.
At the same time, the photo logo is displayed in the upper right corner of the UI interface for 2 seconds. The first photo taken is displayed as 001.

EDG (Vascular Enhancement)

Press EDG key to enhance the display of vascular enhancement.


MWB (Manual White Balance)

Press MWB key to adjust the white balance.
After the white balance is finished, display "WB OK" for 3 seconds in the lower right corner of the display.

Zoom+-

Zoom in, zoom out.
The magnification is displayed in the lower right corner, synchronized with the operation. The rate is 1~4 times.

LED Left/ Right

Adjust LED brightness.
The default value is 60%. 0%, 20%, 40%, 60%, 80%, 100% ranks are displayed in the lower right corner of the display.

OSD

Press the OSD key to switch the view frame.
Press again to switch to other frame (round, octagonal, panoramic).
The corresponding view frames indicator will emit blue light.

FRZ

Press FREEZE button to freeze screen.
You can also press Freeze after zoom+/- . When the picture is frozen, the word "Freeze" appears in the lower right corner.

Return

Press the Return key, the menu function exits and returns to the previous step.

Video

Press the video key to start the video.
During the video recording, you can press the photo key to take photos.
When video is recorded, the video logo and the time of video recording will appear on the upper right corner of the screen.
The blue indicator shows that video is on/off.

Replay

Press the replay key to start playback, watch photos and videos.
In this mode, the left and right keys can choose photos and videos.
Enter key is for video play / pause.
Press replay key to exit the replay mode.

Enter (Menu)

In preview mode, press Enter key to enter Menu settings screen.

- **Menu Function —Name**
Add or modify patient information.
- **Menu Function – Date / Time**
Set date / time
- **Menu Function – Language**

Set different languages

Menu Function – TV mode

Set TV output mode

Menu Function – Default Setting

Reset the default setting

Menu Function – Version

Show software version

Care/ Maintenance/Repair

Caution!

The Video Processor (or parts of the Video Processor) must not be sterilized (e.g. through autoclaving). The parts may be damaged!

Cleaning the Video Processor

The surfaces of the Video Processor including video processor and accessories can be cleaned with medical alcohol (75%). Use a soft, slightly moistened cloth to wipe off the surfaces.

Caution!

The device should be cleaned and disinfected before and after each use. Do not use disinfectant agents and cleaning methods except the cleaning instructions.

Caution!

Liquids or moisture may not penetrate the interior of the Video Processor.

Maintenance Service

For devices used in medical institutions, in addition to inspections before every use by the users or operators, we recommend that the person in charge of the maintenance of the devices to inspect the parts described in the Operation Manual quarterly. If you have a suspicion that the product is abnormal, please do not use the product and inspect it following the instructions in Chapter 7 "Troubleshooting." If you still have a suspicion that the product is abnormal after such inspection, please contact the local sales representative or distributors.

Shanghai AnQing Medical Instrument Co., Ltd. or its authorized professionals must maintain the products in use, and inspect or check whether the parts described in the Operation Manual are abnormal or not. If an abnormality is found or you have a suspicion that the product is abnormal, please record the information and contact local authorized agent or Shanghai AnQing Medical Instrument Co., Ltd.

Warranty Service

Shanghai AnQing Medical Instrument Co., Ltd. guarantees that: for design, materials or

workmanship defects (except for those listed in the exception clauses), we provide free warranty services for the owners registered upon purchase within one (1) year from the date of purchase (warranty period).

If any faults in or damages to a product due to normal use or non-human factors during the warranty period, please contact your local authorized representative or distributors for registration first. If Shanghai AnQing Medical Instrument Co., Ltd. or its authorized professionals judged that the faulty product is qualified for the free warranty services, they will repair or replace the product's parts and components, then return the product to the customer.

If Shanghai AnQing Medical Instrument Co., Ltd. or its authorized professionals judged that the faults are caused by human factors such as improper use or abnormal operations after inspection, the maintenance is not under the warranty and the customer should pay for the maintenance. Shanghai AnQing Medical Instrument Co., Ltd. or its authorized professionals will provide the definite fault causes and maintenance budget to the customer. After communication and acknowledgment with the customer, they will repair or replace the parts and components, then return the product to the customer.

We are responsible for maintaining products or parts returned to our company due to materials or workmanship defects within the warranty period at no charge. If the warranty period of the product or parts under warranty or maintained by Shanghai AnQing Medical Instrument Co., Ltd., has expired, we will provide a 30-day warranty when the faults caused by design and workmanship defects.

Repair Service

Shanghai AnQing Medical Instrument Co., Ltd. guarantees that: except for the listed exception clauses, we provide non-free warranty services for products or parts purchased beyond one (1) year but within five (5) years since the date of purchase.

Please contact your local representative or distributors for registration first. Shanghai AnQing Medical Instrument Co., Ltd. or its authorized professionals will provide the definite fault causes and repair budget to the customer. After communication and acknowledgment with the customer, they will repair or replace the parts, then return the product to the customer.

We are responsible for repairing products or parts returned to our company due to materials or workmanship defects within the warranty period at no charge. If the warranty period of a product or parts under warranty or maintained by Shanghai AnQing Medical Instrument Co., Ltd., has expired, we will provide a 30-day warranty when the faults caused by design and workmanship defects.

Warning

- No part of this product should be maintained /repaired by the customer. Do not disassemble, modify or attempt to maintain this product on your own; doing so may cause injuries to patients or operators, or damage to the devices.

Unauthorized disassembly, modification or maintenance attempts may cause severe safety hazards! The customer is responsible for accidents caused by unauthorized disassembly, modification or maintenance attempts! Shanghai AnQing Medical Instrument Co., Ltd. will not provide any service for such customers.

- **This product belongs to precision instruments. Products sent for repair should be transported in their original package. Our company will not be liable for products damaged or lost during transport.**

Notes:

- The warranty and maintenance/repair services are not applicable to damage caused by human factors such as misuse, improper use, negligence, abnormal operations, unauthorized modifications or adjustments.
- The warranty and maintenance /repair services are not applicable to products whose original serial number or other identification marks have been removed or damaged.
- The warranty and maintenance/repair services are not applicable to products maintained by external persons without authorization or written permission from Shanghai AnQing Medical Instrument Co., Ltd.
- The warranty and maintenance/repair services are not applicable to products that are missing parts or components, or products maintained with parts or components that are not authorized by Shanghai AnQing Medical Instrument Co., Ltd.

7 Trouble Shooting

Error symptoms	Troubleshooting
Device does not function. Device-on indicator at Video Processor is not lit.	• Check power connection • Return device to manufacturer
No picture on monitor	• Check device-on indicator on Video Processor • Check device-on indicator on monitor • Are all connections established properly and are all cables plugged in firmly? • Is the connected signal type enabled in the monitor? • As a control measure, connect Video Processor to different monitor • Connect Video Processor to different endoscope • Return device to manufacturer
Display color not satisfactory	• Check all cable connections • Performing white balance • Use control/reference monitor if necessary • Return devices to manufacturer
Distorted image	• Return devices to manufacturer

EMC Information

This model should be used in the electromagnetic environment specified. The customer or user shall ensure that it is installed and used in such environments.

Warning!

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning!

Use of accessories and cables other than those specified or provided by the MANUFACTURER of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning!

PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this model, including cables specified by the MANUFACTURER.

Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration – Electromagnetic emissions

This model is intended for use in the electromagnetic environments specified below. The customer or the user of this model should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-Guidance
Conducted Emissions	CISPR 11 Group 1 Class B	This instrument uses RF (Radio Frequency) energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated Emissions	CISPR 11 Group 1 Class B	This instrument is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Current Emissions	IEC 61000-3-2	This instrument's harmonic emissions are low and are not likely to cause any problem in the typical commercial power supply connected to this instrument.
Voltage Fluctuation and Flicker	IEC 61000-3-3	This instrument stabilizes its own radio variability and has no effect, such as flicker in lighting apparatus.

ance and manufacturer's declaration – Electromagnetic immunity

This model is intended for use in the electromagnetic environment specified below. The customer or user of this model should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	ENCLOSURE Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air PATIENT coupling Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air Signal input/output parts Port ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Same as left	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material that tends to produce static, the relative humidity should be at least 30%
Electrical transient/burst IEC 61000-4-4	Input a.c. power Port ± 2 kV 100 kHz repetition frequency Signal input/output parts Port ± 1 kV 100 kHz repetition frequency	Same as left	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input a.c. power Port Line-to-line $\pm 0,5$ kV, ± 1 kV Line-to-ground $\pm 0,5$ kV, ± 1 kV, ± 2 kV Signal input/output parts Port Line-to-ground ± 2 kV	Same as left	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11 Voltage interruptions IEC 61000-4-11	Input a.c. power Port 0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles Single phase: at 0° 0 % U_T ; 250/300 cycle	Same as left	Mains power quality should be that of a typical commercial or hospital environment. If the user of instrument requires continued operation during power mains interruptions, it is recommended that the instrument be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
RATED power frequency magnetic fields IEC 61000-4-8	ENCLOSURE Port 30 A/m 50Hz or 60Hz	Same as left	It is recommended to use this instrument by maintaining enough distance from any equipment that operates with high current. If image distortion occurs, it may be necessary to position the Model further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
NOTE U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – Electromagnetic immunity

This model intended for use in the electromagnetic environment specified below. The customer or the user of this model should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	Input a.c. power Port 3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz PATIENT coupling Port 3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz Signal input/output parts Port 3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz	Recommended separation distance: 30cm
Radiated RF EM fields IEC 61000-4-3	ENCLOSURE Port 3 V/m 80 MHz – 2,7 GHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	

80 % AM at 1 kHz

These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strength from fixed RF transmitters as determined by an electromagnetic site survey should be less than the compliance level in each frequency range

Electromagnetic interference may occur in the vicinity of high-frequency electrosurgical equipment and/or other equipment marked with the following symbol:

9 Disposal

This product belongs to electronic devices. It cannot be disposed as unclassified municipal waste. It must be collected and disposed separately.

When the lifetime of the product ends, the user must discard and dispose the product by local laws and hospital management specifications.

Symbol	Description	Symbol	Description
	Caution		Keep dry
	Consult instructions for use		Upward
	IEC 60601 Type BF shock protection		Keep away from sunlight
	IEC 60601 Class II device		Fragile
	Batch code		Waste management
	Serial number		Temperature limit -10°~55°
	Humidity limitation 10%~80%		Atmospheric pressure limitation 500hPa~1060hPa
	Date of manufacture		Use-by date

Manufacturer

CE 0197

CE mark

EC REP

EU Authorized Representative

EC REP

**EU Authorized
Representative**

Immergo Europe
Prinsessegracht 20, 2514 AP The Hague,
The Netherlands
Tel: +31.70.345.8570

1
2

**Legal
Manufacturer**

Shanghai AnQing Medical Instrument Co., Ltd.
Room 202, 2nd Floor, No. 7 Building
No. 150 Cailun Road,
Shanghai High-Tech Park
Shanghai, China
Tel: +86 21 6111 7375
Fax: +86 21 6111 7374

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