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

Подпись/Signature
Joseph M. Vitullo 08/10/2020

JOSEPH M. VITULLO
Notary Public
COMMONWEALTH OF MASSACHUSETTS
My Commission Expires On
March 13, 2026

Instruction for use

Reusable surgical instruments for orthopedic operations for anterior or posterior cruciate ligament reconstructive surgery in the knee joint with accessories: Instrumentation Tray, Instrumentation Tray lid, Guide Carriage, Obturator for Retrograde Reamer Bullet, Retrograde Reamer Bullet Inserts for use with 2.4mm Beath Pin, Restore Drill point pin with eyelet, Outside In Suture Passer

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procedures.

INSTRUCTIONS FOR USE

Refer to the implant specific package insert for use of an instrument with the implant.

CONTRAINDICATIONS

None known.

WARNINGS

Appropriate personal protective equipment and attire should be worn when cleaning and sterilizing soiled surgical instruments.

Instruments may have sharp cutting edges that will puncture skin.

PRECAUTIONS

It is imperative that the surgeon and operating theater staff are fully conversant with the appropriate surgical technique for the surgical instruments.

Inspect for damage prior to use. Replace a damaged, worn or bent instrument. Do not attempt to straighten, sharpen or repair.

Do not use a reusable instrument for any purpose other than its intended use as that may result in damage to the instrument.

Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used.

Do not exceed 137°C during sterilization.

Highly alkaline conditions can damage products with aluminium parts.

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Complex devices, such as those with tubes, hinges, retractable features, matted surfaces, and textured surface finishes, require special attention during cleaning. Manual pre-cleaning of such device features is required before automated cleaning processing.

Avoid exposing instruments to hypochlorite solutions, as these will promote corrosion.

Clean delicate instruments separately from other instruments.

Repeated processing has minimal effects on instrument life and function.

End of useful instrument life is generally determined by wear or damage from handling or surgical use. Inspect instruments between uses to verify proper functioning. (See "Inspection and Functional Testing" instructions below.)

CLEANING AND STERILIZATION INSTRUCTIONS

These cleaning and sterilization instructions have been validated for preparing reusable instruments for reuse. It is the responsibility of the end user to ensure that the cleaning and sterilization is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. This normally requires validation and routine monitoring of the process.

Any deviation from these instructions should be evaluated for effectiveness and potential adverse consequences.

PREPARATION

- Clean instruments as soon as possible after use. If cleaning must

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be delayed, wet the instruments in a compatible liquid solution to prevent drying and encrustation of surgical soil.

- Avoid prolonged exposure to saline to minimize the chance of corrosion.
- Remove excessive soil with a disposable wipe.

Prior to sterilization, instruments should be cleaned by either automated or manual means described below.

AUTOMATED CLEANING

Instrument types with such complex design features as cannulations, lumens, (i.e. tubes), holes, hinged joints, box locks, spring-loaded retractable features, or flexible shafts; must be manually pre-cleaned prior to automated cleaning to improve the removal of adherent soil.

Manual Pre-Cleaning Instructions

- Flush internal areas with a warm detergent solution. Pay close attention to lumens, cannulations, holes, threads, crevices, seams, and any hard to reach areas. While flushing, actuate any moveable mechanisms, (hinged joints, box locks, or spring-loaded features), to free trapped blood and debris.
If the components of the instrument can be retracted, retract or open the part while flushing the area.
- For instruments with flexible shafts, flex the instrument while flushing.
- For instruments with lumens, brush the internal surfaces of the lumen.
- Rinse the instrument thoroughly with warm tap water.
- Ultrasonically clean the instrument for 10 minutes in neutral pH

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detergent, prepared in accordance with the manufacturer's instructions.

Loading Instructions

- Load instruments into washer disinfectant so that hinges are open and cannulations and holes can drain.
- Place heavier instruments on the bottom of containers. Do not place heavy instruments on top of delicate instruments.
- For instruments with concave surfaces, such as curettes, place instrument with the concave surface facing downward to facilitate draining.

Automated Cleaning Instructions

- Clean in a validated washer disinfectant using the "INSTRUMENTS" cycle and a pH neutral cleaning agent intended for use in automated cleaning. Effective automated cleaning can be achieved using the following parameters: two cold water pre-washes (minimum hold of two minutes each); a hot water enzyme wash (minimum hold of 4 minutes); followed by a hot water pH neutral detergent wash (minimum hold 3 minutes at 60°C); a hot water rinse (minimum time 20 seconds); a thermal rinse (minimum hold of 1 minute at 82.2°C) and dry (minimum hold 5 minutes at 95°C).
- After automated cleaning, proceed to inspection, instructions below.

MANUAL CLEANING

Cleaning Instructions

- Prepare an enzymatic cleaning solution in accordance to the manufacturer's instructions.
- Soak soiled instruments for a minimum of 10 minutes. Do not exceed 30°C.

When cleaning, fully immerse the instrument in the cleaning solution to avoid aerosol generation. Use a soft bristle brush to remove all traces

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of blood and debris; pay close attention to any hard to reach areas, textured surfaces, or crevices. When cleaning Instruments with cannulations or lumens (i.e. tubes), or holes, use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub the cannula, lumen, or hole. Push in and out, using a twisting motion to remove debris. Use a syringe filled with enzymatic cleaning solution to flush hard to reach internal areas. When cleaning articulating instruments, (those with moveable parts), brush with a soft non-metallic bristle brush to remove all traces of blood and debris. Pay close attention to threads, crevices, seams, and any hard to reach areas.

Actuate any moveable mechanisms, such as hinged joints, box locks, or spring-loaded features, to free trapped blood and debris. If the components of the instrument can be retracted, retract or open the part while cleaning the area.

For instruments with flexible shafts, flex the instrument under the cleaning solution. Rinse the instrument thoroughly with warm water.

- Ultrasonically clean the instrument for 10 minutes in neutral pH detergent, prepared in accordance with the manufacturer's instructions. Rinse the instrument thoroughly with warm water.

Rinsing Instructions

Rinse the instrument thoroughly with filtered water at 38–49°C, agitate for a minimum of 1 minute; repeat rinse at least two additional times. When rinsing, pay particular attention to flush the cannulations, lumens, or holes with warm tap water. Also pay particular attention to internal areas and moveable parts. Actuate moveable parts while rinsing. If the components of the instrument can be retracted, retract or open the part while rinsing the area. For instruments with flexible shafts, flex the instrument under the rinse solution.

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Drying Instructions

Dry the instrument immediately after final rinse. Use filtered compressed air to dry internal areas. After drying proceed to Inspection instructions.

INSPECTION AFTER CLEANING

- Inspect all instruments before sterilization or storage to ensure the complete removal of soil from surfaces, tubes and holes, moveable parts.
- If areas are difficult to inspect visually, check for blood by immersing or flushing the instrument in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. Rinse instruments thoroughly after using hydrogen peroxide solution.
- If soil is still present, reclean the instrument.

MAINTENANCE

Between uses, lubricate moving parts with a water-soluble lubricant in accordance with the manufacturer's instructions.

INSPECTION AND FUNCTIONAL TESTING

- Visually inspect the instrument and check for damage and wear.
- Cutting edges should be free of nicks and have a continuous edge.
- Jaws and teeth should align properly.
- Moveable parts should have smooth movement without excessive play.
- Locking mechanisms should fasten securely and close easily.

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- Long, thin instruments should be free of bending and distortion.

PACKAGING

- For instruments that are provided in sets, instrument trays may be used for containment during sterilization.
- Package trays/instruments with a barrier wrap material in accordance with local procedures, using standardized wrapping techniques such as those described in ANSI/AAMI ST46-1993.

STERILIZATION

Trays

Place instrument in the appropriate location within the sterilization tray. If a specific location is not identified for an instrument, it may be placed in the general purpose area (pin-mat) ensuring that the sterilant has adequate access to all surfaces including difficult to reach areas, lumens, etc.

- Use a validated, properly maintained and calibrated steam sterilizer.

Effective steam sterilization can be achieved using the following cycles:

Cycle Type	Temperature	Minimum Exposure Time
Prevacuum	132°C (270°F)	4 min
Prevacuum (outside the U.S.)	134–137°C	3 min

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STORAGE

Store sterile packaged instruments in a manner that provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.

Storage conditions:

Temperature: 15-25 °C

Relative humidity: No special conditions. The device is made of non-moisture sensitive material.

Barometric pressure: No special conditions

Transportation conditions:

Temperature: From -10° C to + 50°C

Relative humidity: No special conditions. The device is made of non-moisture sensitive material.

Barometric pressure: No special conditions

Warranty and reclamation

Repeated processing has minimal effects on instrument life and function. End of useful instrument life is generally determined by wear or damage from handling or surgical use. Inspect instruments between uses to verify proper functioning.

Instructions for inspection and functional testing are specified in the IFU.

The manufacturer provides warranty regarding disadvantages of medical device within the whole device usage period which is defined by achievement of natural wear.

The manufacturer does not assume any additional warranty obligations,

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passing larger 4.8mm TWISTR Reamer
Useful to create pilot hole in patients with hard bone

Drill point pin with eyelet

Enables user to pre-drill with drill pin prior to using the retroreamer, as well as used to pass suture during the procedure in to the joint for facilitation of graft passage.

Outside In Suture Passer

Fork device used to pass utility suture into joint to facilitate graft passage
Fits inside TWISTR Reamer Sleeve w/suture.
Forked tip holds suture.

Instrumentation Tray

The Instrument Sterilization Tray is intended to hold DePuy Mitek Instruments during the sterilization process and for storage.

Instrumentation Tray lid

Instrumentation Tray lid and Instrumentation Tray are intended to hold Instruments.

Symbols of marking

Table «The deciphering of the symbols, which using in the marking MD»

	Catalogue number
 YYYY-MM-DD	Date of Manufacture (YYYY-MM-DD)
	Caution! Consult instructions for use

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




	Number of Units
	Nonsterile
	Federal Law (United States of America) restricts this Device to sale by, or on the order of a physician.
	Consult instructions for use
	Product certification code in the European Union
	Batch code
	Manufacturer
	Authorized representative in the European community
	Recycle

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List and description of md materials

Name	Material of MD
Guide Carriage	SS
Obturator for Retrograde Reamer Bullet	SS
Retrograde Reamer Bullet Inserts for use with 2.4mm Beath Pin	SS

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**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER
IN RUSSIA:**

Johnson & Johnson LLC,
121614, Moscow, Krylatskaya str., 17/2.
Tel. +7 (495) 580-7777, Fax: +7 (495) 580-7878
DL-JNJRUMedicalQA@ITS.JNJ.com

ПРОИЗВОДИТЕЛЬ: Medos International SARL
Chemin-Blanc 38, CH-2400
Le Locle, Switzerland, Швейцария

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