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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of DONGGUAN ZSR BIOMEDICAL TECHNOLOGY CO., LTD on the annexed DOCUMENT is genuine. The attached Russian version of the said DOCUMENT is true to the English text.

China Council for the Promotion
of International Trade

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Authorized
Signature: Lyu Cuilian

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(Date: Mar. 24, 2021)

证明书查询 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Disposable Circumcision Anastomat

ZSR-DCA-34, ZSR-DCA-30A

ZSR-DCA-26, ZSR-DCA-22

ZSR-DCA-18, ZSR-DCA-14,

ZSR-DCA-10

This instruction for use is used for the model of:

ZSR-DCA-34, ZSR-DCA-30A, ZSR-DCA-26, ZSR-DCA-22, ZSR-DCA-18, ZSR-DCA-14, ZSR-DCA-10.

PLEASE READ ALL WARNINGS AND INSTRUCTIONS BEFORE USE. USE THE PRODUCT ONLY FOR ITS INTENDED USE AS DESCRIBED IN THESE INSTRUCTIONS. IMPROPER USE CAN RESULT IN SERIOUS OF FATAL ILLNESS/INJURY, IMPROPER TREATMENT OF PROPERLY DAMAGE

ALL THE CIRCUMCISION ANASTOMAT ARE STERILIZED AND SINGLE USE ONLY FOR HYG

Dongguan ZSR Biomedical Technology Company Limited

Intended Use

This device is used for clinical male genital circumcision suture surgery and is suitable for patients with phimosis, paraphimosis and redundant prepuce.

Contraindications:

1. Man with malformed prepuce or glans (including the person whose prepuce and glans are widely and tightly adhered together) should avoid circumcised use this product.
2. Man whose prepuce and glans are infected and the person with edema around glans should avoid circumcised use this product.
3. Man who is suspected to have penile carcinoma should avoid circumcised use this product.
4. Man who suffers from systematic disease (such as bleeding tendency, hypopro-teinemia, serious cardiovascular disease, etc.) and mental disease should avoid circumcised use this product.

Side effect

Complications (there is no side effect with disposable circumcision anastomat) :

1. Bleeding (patients with infection or serious phimosis need to be separated before surgery , or strenuous exercise after the surgery may course bleeding)
2. Infection (The urine is not cleaned up and causes infection)
3. Edema (patients with infection or choose a smaller size will have the chance to course edema)

Potential consumer

Medical professionals: urologists, surgeons

Application area:

Urology: male, child, general

The principle of operation

The operating principle of anastomat is similar with stapler. Disposable circumcision anastomat through the mechanical driving device to push the circular blade and the suturing nail which are placed in the product in advance. After pass through silicone pad and two layer of tissues, the suturing nail will be blocked by the front pin-resisting seat and they will be bent inwards, in the form of the letter B to bent inward to make two layers of tissues and pad together, and the circular blade cut the extra prepuce at the same time. Silicone pad performs the function of high-quality hemostasis and formation of an aesthetic seam. The suturing nail is made of 316L stainless steel, this material has good biocompatibility. As suturing nails are lined up with equal distance and suturing tightness is appropriate .Thus, the defects of the manual suture which are the seaming is too sparse or too dense and the ligation is too tight or too loose can be avoided, which not only guarantees the good healing of issues, but also largely shorten the operation time.

The silicone pad, together with the suturing nail, remains on the proximal part of the circumcised foreskin of the penis. Their removal (pad and suturing nail) occurs spontaneously after the healing of the suture line (7-10 days). On the third day after the operation, you can cut the silicone ring (gasket) with scissors in the gap between the two brackets (this process will help the suturing nail fall out).

To determine the size of the anastomat, each product is supplied with a Slizer- a cardboard ruler with 7 holes cut out, corresponding to the diameter of the anastomats with an indication of their sizes. The slizer is subjected to gas sterilization, but does not remain sterile, due to the lack of individual packaging. Before use, it must be treated with any skin antiseptic. Size of the slicer (length 19.0 cm x width 5.5 cm $\pm$ 5%).

Anastomat complies with the requirements of Directive 93/42 / EEC on medical instruments.

Classification: class of potential risk of use 2A.

Type of medical device according to the nomenclature classification: 351600. OKPD2 code 32.50.13.190

Description of all parts of the product

The anastomat consists of an internal concave seat and an outer host shell


The internal concave seat - this is the internal conical removable part of the anastomat. The inner surface of the riser is used to protect the glans penis from possible damage during the application of the anastomosis; on its outer surface, the

foreskin is fixed for circumcision. The concave seat in its proximal part has the shape of a bell from the outside with an emphasis for suturing nail for their automatic bending along the seam line. The base of the concave seat is a metal pin that has a threaded connection in its distal part, through which it is attached to the adjustment knob.

The outer shell has two flexible handgrip, an adjustment knob, a circular blade, suturing nail, silicone pad, a protective cap and a safety pin. The circular blade, suturing nail and silicone pad are located inside the outer shell and are protected by a protective cap and safety pin.

The silicone pad is an important integral part of the anastomat and is integrated directly at the location of the stitching staples (upper section of the staples compartment). It provides stabilization of the suture in the event of an early erection, performs the function of high-quality hemostasis and the formation of an aesthetic suture.

The structure and name of the parts are shown below:

No.	Name of the part	No.	Name of the part
1	concave seat	10	connect part
2	nail seat	11	adjusting knob
3	circular anvil	12	safety pin
4	nail-cabin protective cap	13	fixed frame
5	nail-cabin	14	mover
6	suturing nail	15	flexible handgrip support
7	host shell	16	spring
8	screw	17	the end cutting body
9	flexible handgrip	18	circular blade

Technical specification

Disposable circumcision anastomat is available in the following versions:

1. Disposable Circumcision Anastomat, 10 mm (ZSR-DCA-10) in the composition of: - Disposable Circumcision Anastomat, 10 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
2. Disposable Circumcision Anastomat, 14 mm (ZSR-DCA-14) in the composition of: - Disposable Circumcision Anastomat, 14 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
3. Disposable Circumcision Anastomat, 18 mm (ZSR-DCA-18) in the composition of: - Disposable Circumcision Anastomat, 18 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
4. Disposable Circumcision Anastomat, 22 mm (ZSR-DCA-22) in the composition of: - Disposable Circumcision Anastomat, 22 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
5. Disposable Circumcision Anastomat, 26 mm (ZSR-DCA-26) in the composition of: - Disposable Circumcision Anastomat, 26 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
6. Disposable Circumcision Anastomat, 30 mm (ZSR-DCA-30A) in the composition of: - Disposable Circumcision Anastomat, 30 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.
7. Disposable Circumcision Anastomat, 34 mm (ZSR-DCA-34) in the composition of: - Disposable Circumcision Anastomat, 34 mm- 1 piece; - Slizer – 1 piece; - instructions for use – 1 piece.

The anastomat is a single-stage stapler, with a knife and a lever drive of the sewing mechanism and a screw connection of the concave seat with the main body of the anastomat.

Material of manufacture: stainless steel, polycarbonate, polypropylene, ABS plastic, silicone, zinc alloy No. 3, brass, aluminum alloy. The product does not contain any medicinal products in its composition.

The materials from which the suturing nail parts in contact with the tissues are made must be resistant to aggressive biological fluids (blood). Metal parts of the anastomat must be made of corrosion-resistant steels and alloys.

- There should be no cracks, scratches, burrs, grinding marks or other mechanical defects visible to the naked eye on the outer metal surfaces of the blade:

- The knife blade should be sharp with no chipping, cracks, spalled areas and dents. The width of the cutting edge-no more than 8 microns, the height of irregularities - no more than 10 microns

- on the surface of anastomat there should be no cracks, nicks, scratches, burrs, shells, pores, the thrown out places, stratifications, burns, dross, particles of materials of grinding, polishing and traces of greasing

- Move the moving parts of the machine should be free, without jerks and jams.

- The force of pressing the handles necessary for living should be no more than 200 N.

- nails should not fall out of the bracket slots under the action of its own weight. When stitching paper nails should be completely pushed out of the slots.

The principal dimensions of anastomatic are listed below:

Nº	Feature name	ZSR-DCA-10	ZSR-DCA-14	ZSR-DCA-18	ZSR-DCA-22	ZSR-DCA-26	ZSR-DCA-30A	ZSR-DCA-34
1	External diameter of the anastomat, mm	ø21.0±1.0	ø26.0±1.0	ø30.0±1.0	ø34.0±1.0	ø38.0±1.0	ø42.0±1.0	ø46.0±1.0
2	Diameter of a seam, mm	ø14±0,3	ø18±0,3	ø22±0,3	ø26±0,3	ø30±0,3	ø34±0,3	ø38±0,3
3	The diameter of the blade, mm	ø10±0,3	ø14±0,3	ø18±0,3	ø22±0,3	ø26±0,3	ø30±0,3	ø34±0,3
4	Height of the blade, mm	20,0±0,5	20,0±0,5	21,0±0,5	24,0±0,5	26,0±0,5	26,0±0,5	26,0±0,5
5	The angle of sharpening the blade, no more than, deg.	30						
6	Inner diameter of the concave seat, mm	ø5,0±0,5	ø9,0±0,5	ø13,0±0,3	ø17,0±0,5	ø21,0±0,3	ø24,0±0,3	ø28,0±0,5
7	Outer diameter of the concave seat, mm	ø14,2±0,3	ø18,2±0,3	ø22,2±0,3	ø26,2±0,3	ø30,2±0,3	ø34,2±0,3	ø38,2±0,3
8	Depth of the concave part of the concave seat, mm	16,0±2,0	16,0±2,0	18,0±2,0	18,0±2,0	18,0±2,0	22,0±2,0	22,0±2,0
9	Type of seam	single row, multi-row, longitudinal, circular						
10	The number of the nails, pcs	8	10	12	15	18	20	28
11	The weight of the tool, no more than, g	185,0	190,0	195,0	205,0	215,0	220,0	230,0
12	Size of the U-shaped nails (diameter x height of the leg x width of the back), not more than, mm	ø0,3 x 2,7 x 3,5						
13	Height and shape of the closed nail, not more than, mm	1,8 B-shaped						

The product is operated inside of organs and cavities of human body in the presence of biological fluids, tissue, and skin secretions (blood, urine, etc.) and must perform its functions and to be serviceable in the process of operation at a nominal temperature of 32 to 42°C, should be resistant to the effects of biological fluids and secretions of the tissues with which they come into contact in the process of operation. - relative humidity 30-85%; - atmospheric pressure-from 84.0 to 106.7 kPa (from 630 to 800 mm Hg. art.).

- Depending on the mechanical effects, the anastomat belongs to wearable, portable products, and is not intended for use when carrying and moving within a medical facility.

Operation Process

The operation shall be prepared according to the hospital rules.

Before each use, prepare and check the product in accordance with the instructions below. If there are any minor deviations from normal operation, stop using the products and use a replacement product. Damage or malfunction of the product may reduce its safety for the patient or the operator, that is, create a risk of infection, tissue irritation, or injury to the patient.

1. You should always have in stock are ready for immediate use replacement products.

2. Verification

- check the sterile packaging for the presence of breaks, spills, or violations of integrity. If the above-mentioned packaging defects are detected, the sterility of the instrument may be violated. In this case, do not use the product.
- an external inspection of the product: open the package, applying standard rules of reimbursement, make sure that no damage anastomatic.

3. Operation.

1. Do various inspections to the patient and measure the size of glans with a measuring scale. The slicer must be treated with any skin antiseptic before use. The place between the glans and 2/3 of coronary sulcus will be measured (the experienced doctor can instinctively choose the correct size of product without use measuring scale). Then a appropriate disposable circumcision anastomat was selected. Then disinfect the surgical spot and do local anaesthesia, or choose other effective anesthesia ways. For the patient who has adherent prepuce, over short frenum or other congenital anomalies, carry on prepuce cutting and suturing operation after the medical treatment is done and blood is stopped.

2. Anticlockwise adjust the location adjusting knob at the end of the product, and gently take out the concave seat at the front of the product, as shown in figure 1;

3. Putting the concave seat inside the prepuce (if the prepuce is adherent, the adhesion shall be tripped. If it is a phimosis, the preputial orifice shall be expanded or cut open with scissors. Several small preputial orifices can be cut open to put the concave seat on the glans. But one preputial orifice can not be cut too deep and close to coronary sulcus so as to avoid that it exceeds the circumcision part, and thus the wound can't be sutured), and make the prepuce turn out and wrap the outside of concave seat and cover the glans. Make the concave seat tilt and maintain the inclination similar with that of the coronary sulcus, and appropriately leave a bit more frenum so as not to cut the frenum too short. The prepuce naturally fits the curve top of concave seat, as shown in figure 2;

4. Gently remove the nail-cabin protective cap which is used to protect the suturing nails to not fall off, and slowly insert the screw of will-fixed concave seat into the center hole at the main body of anastomat. Then, clockwise screw the location adjusting knob to make the plane at the end of screw parallel with the plane at the back of location adjusting knob. The tightness is subject to the specific circumstances of patients. The skilled and experiences doctor can appropriately screw it by feeling, as shown in figure 3;

5. After adjust the knob in place, make sure that the place where the prepuce will be cut is right and the prepuce is not loose. After everything are confirmed without any problem, remove the safety pin, as shown in figure 4;

6. Slowly and horizontally hold the flexible handgrip with strength and trigger the anastomat for about 5 to 10 seconds (at this time, use dry gauze to wrap the operation part to prevent blood leakage), as shown in figure 5;

7. Loosen the flexible handgrip and assist to reset it. Anticlockwise and slowly unscrew the location adjusting knob. After the knob is back out for 5mm to 7mm, observe the cutting situation. Press the knob forward and separate the concave seat from the main body of anastomat. Slowly unscrew and back out the anastomat as well as the cut prepuce, and observe if the prepuce is stuck or not. If the prepuce is stuck, slightly press it with fingers or scissors, and then the prepuce will be separated naturally. Don't aggressively move arms or excessively screw the concave seat to avoid the incision from being torn, as shown in figure 6;

8. After remove the instrument, check the wound. If the wound is not completely sutured due to carelessness, the wound can be appropriately sutured. The wound should be flat and the suturing should be uniform, as shown in figure 7;

9. Disinfect the wound and wrap it with dry gauze, and then dress it with an elastic bandage for 3 to 5 hours, as shown in figure 8.

Operation Chart

Postoperative Cares

1. After the operation, the patient shall follow doctor's advice, coordinate with postoperative care and use drugs. Before postoperative recovery, avoid over-activity and try best to avoid erection. Do not drink alcohol, eat seafood and spicy food. Within one month after operation, it is forbidden to have sexual intercourse;
2. The wound shall be dressed with elastic self-adhesion bandage after operation. The bandaged surgical dressing can also be stitched, which can prevent the bandage from being loose and falling apart. The suggestion for the time of binding up wound is 7 to 10 days;
3. About the postoperative pain, the wound will be slightly painful after the anesthetic effect disappears. But several hours later, the pain will gradually be eased and relieved. If the pain is very serious, the patient can take analgesic drugs to relieve the symptom;
4. The postoperative edema usually doesn't need any measures to be taken. The edema will fade away by itself. If it is serious, the patient shall go to hospital for treatment.
5. During the wound healing period after the operation, the liquid will ooze from the wound. When there are circumstances like that the wound becomes red, swollen, painful, heating, etc, the patient shall go to hospital for treatment to avoid the symptom from becoming more serious;
6. After operation, if the wound is torn due to accident or strenuous exercise, the patient shall immediately go to hospital for treatment;
7. After operation, the patient shall pay attention to avoiding strenuous exercise. It will not be affected if the work doesn't include lots of activities. One week later, the patient can take a shower. After shower and later everyday, iodophor can be used to disinfect the wound; and then the wound can be simply dressed by patient himself;
8. On the third day after operation, the patient shall go back to the hospital for inspection and let the doctor to observe the recovery condition to master the situation of patient. At this moment, cut off white soft oblate ring by scissors at the gap between each two nails (this process can help nails fall out);
9. After operation, the patient shall maintain the cleanliness of wound and use dilute iodophor to clean up the wound once in the morning and once in the evening, and wear breathable loose trousers and keep the wound dry, and do routine postoperative care to prevent infection.

Notes

1. Do not use if package is damaged or out of validity date;
2. His product is single use only, please don't repeatedly disinfect it for use;
3. This product shall be used by qualified doctors to do operation;
4. Must select the product of proper size to do operation. Do not to use the product of improper size to do operation;
5. After operation, this product shall be professionally handled as the medical waste; avoiding circular blade hurts others;
6. Please take appropriate rest under the doctor's suggestion after operation. Pay attention to the local hygiene and the anti-inflammatory treatment. Please follow the doctor's advice to use drugs. It is forbidden to take strenuous exercises to avoid hematoma;
7. In one week after the wound healing, the suturing nail will slowly fall off. Due to the different physiques, if the suturing nail doesn't fall off in 3 weeks, the patient shall go back to hospital to remove it manually;
8. Around the prepuce of high-sensitivity physique, if there are slight inflammations, local edema may appear. At this time, anti-inflammatory and anti-allergic drugs can be used to control that;
9. Pay attention to hygiene after operation; control and eliminate the factors that may cause infections.

Marking

The primary packaging label contains: - trademark of the manufacturer; - address of the manufacturer; - the name of the product; - batch number; - catalogue number; - the number of products in the package; - shelf life; - appointment (in international English); - contraindication (in international English); - the main conditions for storage and transportation (international English language); - color indicator of sterility of EO; m- QR-code; - manipulation signs corresponding to: «Caution: Consult accompanying documents», «Single use only», «Batch code», «Catalogue number», «The product complies with current European Union directives», «Use-by date», «Sterilized using ethylene oxide», «Do not use if package is damaged», «Authorized Representative in the European Community».

Individual labelling of the outer packaging contains:

The individual outer packaging must bear the following markings: - trademark of the manufacturer; - address of the manufacturer; - the name of the product; - batch number; - catalogue number; - item code; - number of products in the package; - production date; - shelf life; - product description (in international English); - compliance with standards; - manipulation signs corresponding to: «Caution: Consult accompanying documents», «Single use only», «Batch code», «Catalogue number», «The product complies with current European Union directives», «Fragile», «This way up», «Use-by date», «Use-by date», «Sterilized using ethylene oxide», «Do not use if package is damaged», «Authorized Representative in the European Community», «Keep Away from rain», «Recycle».

The group package label contains:

trademark of the manufacturer; - address of the manufacturer; - the name of the product; - batch number; - catalogue number; - item code; - number of products in the package; - production date; - shelf life; - product description (in international English); - compliance with standards; - manipulation signs corresponding to: «Caution: Consult accompanying documents», «Single use only», «Batch code», «Catalogue number», «The product complies with current European Union directives», «Use-by date», «Sterilized using ethylene oxide», «Do not use if package is damaged», «Authorized Representative in the European Community».

Layout of additional marking for delivery to the Russian Federation

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Packaging

Packaging requirements-must comply with the requirements of EN ISO 11607-1: 2009 and must ensure the safety of the product during transportation and storage. Individual consumer packaging anastomat should not have damages and should be sealed, must not allow the penetration of microorganisms.

Anastomat in consumer packaging, Packed in transport containers, must be resistant to mechanical stress.

Individual primary packaging: 1 piece of anastomat packed in a blister pack (individual sterile packaging). Package size 20.4 x 13.4 x 59 cm. Packing material: manufacturer Shenzhen Ganglin Co., Ltd front part: TYVEK PETG, back part: dialyzed paper. Anastomat in individual primary packaging is placed in an individual outer packaging-a cardboard box.

Complete set (individual outer packaging):

- Anastomat in individual primary packaging – 1 piece;
- instructions for use-1 piece,
- slizer (to select the required size of the anastomat) - 1 piece.

Package size: 21, 5x14, 1x65 cm.


Group packaging: a cardboard box containing 30 pieces of anastomates in individual outer packaging. Package size: 69, 8 x 44, 8 x 22.8 cm

Transportation

And we the anastomats are transported by all modes of transport

During transportation of the anastomosis must be ensured:

- ambient temperature– - 50 to 70 °C;
- the relative humidity of the surrounding air is not more than 80% at a temperature of plus 25 ° C;

- atmospheric pressure - from 84.0 to 106.7 kPa (from 630 to 800 mm Hg).

Safety requirement

The anastomat must be sterile.

The materials from which the parts of the anastomat having contact with the patient are made must have toxicological safety.

Storage and shelf life

The item must be stored in the warehouses of the supplier and the consumer.

Sterile, packaged products in transport or group or consumer packaging must be resistant to climatic influences in the temperature range from minus 50 to plus 70 ° C, out of the reach of outsiders, in ventilated warehouses located in any macroclimatic areas, avoid direct sunlight.

- the relative humidity of the surrounding air is not more than 80% at a temperature of plus 25 ° C;

- atmospheric pressure - from 84.0 to 106.7 kPa (from 630 to 800 mm Hg).

The warranty period of use of machines, which are in consumer packaging, is 2 years from the date of manufacture, subject to tightness of the package.

Utilization

1. after the application of the anastomat, it must be destroyed in accordance with the generally accepted method;

2 Disposing of anastomat should be carried out in accordance with the developed in the institutions by the rules of waste disposal and in the following sequence:

- processing of anastomatic after use in dezinficiruyusch solution in accordance with the manufacturer des.tools;

- heat treatment -incineration, pyrolysis.

Hazard class of medical waste based on the characteristics of the morphological composition: class B.

Proper disposal does not harm the environment.

Warranty obligation

The manufacturer guarantees the compliance of the product with all the parameters and properties specified in this document, in compliance with the conditions of transportation, storage, use in accordance with the instructions for use during two years.

Authorized representative for the Russian Federation. Complaints Department for the Russian Federation:

Limited liability company "OMEGAMEDSERVIS",

190020, St. Petersburg, St. Paper, d. 17, office 125.

(812) 495-97-74 (multi)

el. address: info@omegamed.ru

Dongguan ZSR Biomedical Technology Company Limited

Unit 448, Qingfeng Road, Taihu Village, Sanzhong Village Committee, Qingxi Town, 5236 51

Dongguan City, Guangdong Province, R.P.China

Fax: + 86 769 - 87890862 Tel: +86 769-87328976

Email: sales@zsrbiomedical.com Website : www.ZSRbiomedical.com

Renault-Petersen Limited

Адрес: 5 Bankside, Hanborough Business Park, Witney OX29 8 LJ UK

Tel: +44 1993882779 Fax: +44 1993880110

Character descriptions

	Caution: Consult accompanying documents		Use-by date
	Do not use if package is damaged		Batch code
	Sterilized using ethylene oxide		Catalogue number
	Single use only		Fragile
	Authorized Representative in the European Community		This way up
	The product complies with current European Union directives		Keep Away from rain
	Manufacturer		Recycle

Applied Standards

1	MDD 93/42/EEC
2	MEDDEV 2.7.1 Evaluation of clinical Data
3	MEDDEV 2.12-1 rev. 8 Vigilance
4	MEDDEV 2.12/2 rev2 Post market clinical follow-up studies A Guide for Manufactures and Notified Bodies.
5	2013 473 EU on the audits and assessments performed by notified bodies in the field of medical devices
6	EN ISO 13485:2012 Medical devices-Quality management systems-Requirement for regulatory purposes
7	EN ISO 14971:2012 Medical Devices-Application of risk management to medical devices
8	EN ISO 15223-1:2012 Medical devices.Symbols to be used with medical devices labels,labeling and information to be supplied. Part 1: General requirements
9	EN 1041:2008 Information supplied by the manufacture of medical device
10	EN 980:2008 Symbols for use in the labeling of medical devices
11	EN ISO 10993-1:2009/AC:2010 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
12	ISO 10993-5:2009 Biological evaluation of medical devices -Part 5:Tests for in vitro cytotoxicity
13	ISO 10993-7:2009 Biological evaluation of medical devices -Part 7 Ethylene oxide sterilization residuals
14	EN ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
15	ISO 11135-1:2007 Sterilization of health care products -Ethylene oxide-Part 1:Requirement for development ,validation and routine control of a sterilization process for medical devices.
16	ISO 11737-1:2006 Sterilization of medical devices-Microbiological methods-Part1: Determination of a population of microorganisms on products.
17	ISO 11737-2:2013 Sterilization of medical devices-Microbiological methods-Part2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
18	ISO 11138-1:2006 Sterilization of health care products-Biological indicators-Part 1:General requirements
19	ISO 11138-1:2006 Sterilization of health care products- Biological indicators-Part 2 Biological indicators for ethylene oxide sterilization processes
20	ISO 11607-1:2006 Packaging for terminally sterilization medical devices-Part 1: Requirements for materials, sterile barrier systems and packaging systems
21	EN 556-1: 2001/AC: 2006 Sterilization of medical devices-Requirements for medical devices to be designed "STERILE" - Part1:Requirements for terminally sterilized medical device
22	EN ISO 11140-1:2009 Sterilization of health care products - Chemical indicators - Part 1: General requirements (ISO 11140-1:2005)
23	ISO 5832-1:2007 Implants for surgery-Metallic materials-Part 1: Wrought stainless steel

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