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China Council for the Promotion of International Trade
China Chamber of International Commerce

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China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhang Bo

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(Date: Aug. 26, 2020)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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

184541

OPERATIONS MANUAL

DKT200-2A Laboratory Dry Bath Incubator

PREFACE

Thank you for choosing a laboratory dry bath incubator (hereinafter referred to as device). This operations manual describes functions and operation of the device. Please read this manual carefully to use the product correctly.

Opening Check

Check the product and components according to the packing list when opening the package for the first time. If there are any discrepancies, contact the seller or the manufacturer.

Manufacturer: Hangzhou Miu instruments Co., Ltd.

Address: Room 201, building 1, No. 1418-3 Moganshan Road, Hangzhou, China

Telephone: +86-(0)571-87653906

Fax: +86-(0)571-87653907

E-mail: sales@miulab.com

Website: <http://www.miulab.com>

Authorized representative of the manufacturer in the Russian Federation:

PlazmoVita LLC

Address: Russia, 115280, Moscow Leninskaya Sloboda Str., 9, prem. 2 room 37-3, office 2

Telephone: 8 (968) 4339241

E-mail: plasmovita@plasmovita.com.

Purpose: The DKT200-2A Laboratory Dry Bath Incubator is used to cultivate a fibrin clot in clinical laboratory diagnostics.

Potential consumer: Laboratory Technician of any clinical diagnostics laboratory.

The device is not contraindicated.

No side effects related to the device use have been identified.

Medicines, animal and human materials are not included in the device.

Applicable standards: EN 61010-1:2010, IEC 61010-2-010:2014, EN 61326-1-2013, EN 15223-1:2012.

The device is delivered non-sterile and is not sterilizable by the user.

Device life: 5 years.

In the device covered by the European Directive 2012/19/EU (WEEE Directive on Waste Electrical and Electronic Equipment), so it must be disposed of separately from household waste.

At the end of its useful life, the device should be disposed of in accordance with national regulations in the country of use of the medical device for this type of waste.

Warnings and Safety Instructions

Important safety information in the operation

Users should be fully acquainted with the principle of operation of the device before use. Please read this operations manual carefully prior to use the device.

Do not operate without reading the operations manual. Read the recommendations and instructions below and take action accordingly.

2 Safety

Observe the basic instructions and warnings below when operating, servicing and repairing the device. This affects the expected lifetime of the device and also the protective equipment.

This product is standard for indoor use.

Read this manual carefully before working with the product. This device is intended for use in a laboratory environment by persons qualified in good safe laboratory practice.

The operator must not open or repair the device himself. Failure to do so will result in a breach of the warranty repair or an accident may occur. The manufacturer warrants the repair of the device in accordance with the warranty description.

The grounding of the AC source must be a good protection against electric shock. The 3-pin connector supplied with the incubator power cable is a protection device that must match the parameters of a specific grounded outlet.

The temperature of the metal block is very high during normal operation and may cause burns or boiling of liquids. Never touch the product with any part of the body.

Cap tubes before placing in the block. Liquid may spill in the block or on the device when the tube lid is opened, resulting in damage to the block or device.

Verify that the voltage is correct. Verify that the rated load of the electrical outlet is not below the specified value. The power cord must be replaced with the same type of cable if damaged. Make sure that no objects are placed on the power cable. Hold the connector while removing the plug. Do not pull the power cable. Do not place the cable in the aisles.

Place the device in a location with low temperature, low dust, no moisture, no sunlight or bright light, good ventilation, no corrosive gases or strong magnetic fields, away from central heating systems, portable stoves, and other heating sources. Do not place the product in damp or dusty areas. A vent on the block is for air supply. Do not close or cover the vent. The distance between devices must be at least 100 cm when

multiple devices are present.

The main switch is on the rear of the device. Turn the switch to “I” to power on the device and to “O” to power off the device.

Turn off the power when you finish the operation. If the product is not used for an extended period of time, disconnect the plug and cover the product with a cloth to prevent dust.

In the following cases, remove the plug from the connector immediately and contact the supplier:

- Liquid ingress into the product.
- Water or fire.
- Abnormal operating conditions, such as abnormal sound or odor.
- Drop of the device or damage to the outer shell.
- An obvious change in functioning.

3 Maintenance of the device

The block wells should be cleaned with alcohol-soaked tissue to ensure good heat conduction between the block and tube and absence of contamination. If stains occur on the device, wipe them off with a cloth.

Turn off the power while cleaning the product.

Do not allow cleaning fluid to drip into the well during cleaning.

The use of aggressive cleaning agents is strictly prohibited.

The product does not contain any user-serviceable parts and should not be opened, or no independent attempts of its repair must be undertaken. For all types of faults, please describe their symptoms in detail and contact the manufacturer's authorized representative for advice and repair or qualified service personnel.

Contents

Chapter 1: Introduction..... - 6 -

Chapter 2: Specifications..... - 6 -

1. Standard operating conditions..... - 6 -

2. Basic parameters and functional characteristics - 6 -

Chapter 3: Basic operating instructions..... - 8 -

1. Construction - 8 -

2. Operation panel - 9 -

3. Buttons..... - 9 -

Chapter 4: Operation guide - 9 -

1. Temperature and time setting..... - 9 -

2. Start and stop..... - 10 -

3. Temperature calibration - 10 -

4. Replacement of metal block..... - 12 -

Chapter 5: Error analysis and troubleshooting- 13 -

Chapter 6: Additional information for use in the territory
of the Russian Federation.....- 14 -

Annex 1: Electrical wiring diagram of the DKT200-2A

Laboratory Dry Bath Incubator - 15 -

Chapter 1: Introduction

The Laboratory Dry Bath Incubator is a microprocessor-controlled thermostat with high temperature control accuracy. Suitable for parallel sample preparation. It is used as an alternative to traditional water bath devices. Widely used for cultivation of specimens in clinical laboratory diagnostics.

Characteristics

1. LCD display. Excellent compact design.
2. High heat speed, even heat, precise temperature control, high stability, low power consumption and no noise.
3. Built-in temperature calibration, automatic fault detection and audible alarm.
4. Built-in, safe and reliable thermal protection increases the life of the incubator.
5. Compact, impermeable design, takes up little space, moves easily and freely.
6. Various blocks for easy replacement, easy cleaning and disinfection.

Chapter 2: Specifications

1. Standard operating conditions:

Air Temperature: 5-30°C

Relative air humidity: ≤ 70%

Voltage and frequency used: 110/220 V, 50/60 Hz

2. Key parameters and features:

Laboratory Dry Bath Incubator	DKT200-2A
Temperature Control Range	5°C~150°C
Temperature Setting Range	5°C~150°C
Temperature Stability @ 40°C~100°C	≤ ± 0,5°C
Temperature Stability @ >100°C	≤ ± 1°C
Block Temperature Uniformity @ 40°C	≤ ± 0,3°C
Block Temperature Uniformity @ >100°C	≤ ± 0,5°C
Temperature Display Accuracy	0.1°C
Heating Speed	≤ 30min (20°C to 150°C)
Time Range	1 min~99 h 59 min
Sample Capacity	2 individually controlled blocks
Voltage	AC 220/110 V
Frequency	50/60Hz
Power	400 W
Fuse	250 V, 3 A/6 A
Dimensions (width x depth x height)	290 x 260 x 95 mm
(Net) weight	4.1 kg




omposition of the DKT200-2A Laboratory Dry Bath Incubator:

- Incubator - 1 piece
- Power cable - 1 piece
- Screwdriver - 1 piece
- Operations Manual - 1 piece
- Heating block for 10 ml tubes - 1 piece
- Heating block for 5 ml tubes - 1 piece
- Heating press-block with two notches - 1 piece
- Heating press-block with three notches - 1 piece
- Press for heating press-block with two notches- 2 pieces
- Press for heating press-block with three notches - 1 piece
- Raiser - 1 piece


Chapter 3: Basic Operating Instructions

This chapter provides a brief description of the mechanical design of the device, the buttons and functions of each button, and instructions for pre-power-up preparation. Read this chapter before using the incubator for the first time.

1. Construction

2. Operation panel

3. Buttons

Down arrow button: Decreases the setting value

Up arrow button: Increases the setting value

Start/Stop Start/stop button: To start, press the button after setting the temperature and time. Press this button for 2 seconds to stop the device.

Chapter 4: Operation Guide

1. Setting the temperature and time

a) Turn on the power, the screen will look like the picture on the right. The device returns to its original state and you hear a beep.

b) After two seconds, the values change to the current temperature of the block and the time setting of the last operation. In the figure on the right, "28.5" is the current temperature of the block and "00:35" is the last time setting, which means 35 minutes.

c) Press   or  , the temperature value changes to the setting value of the last operation, the last digit of the setting value flashes. Press   or  to change the temperature setting. Pressing and holding the button changes the value ten, one hundred or one thousand times faster. When the desired value is reached, release  or  , the device will automatically confirm and save the value.

d) Press   or , the last digit of the time setting flashes.

Press   or  to change the time setting. Pressing and holding the button changes the value ten, one hundred or one thousand times faster. When the desired value is reached, release

 or , the device will automatically confirm and save the value.

Notes: The time value "00:00" means there is no time setting for the operation, the block is running continuously at the set temperature.

2. Start and stop


a) when the temperature and time settings have been reached, press the Start/Stop button to start the device. An audible alarm will sound when the temperature is increased.

The "." dot in the temperature value flashes intermittently when the device is heating up. When the temperature reaches the setting value, the "." stops flashing and the ":" colon symbol in the time value flashes and the time counts down.

When the time ends, the process stops and a beep sounds. The LCD display shows the current temperature of the block and the time "oUEr", indicating that the operation is complete.

b) the product enters standby mode when it is finished. Press  or  to reset the temperature and time, then press the Start/Stop button to start a new operation. Press the Start/Stop button without changing the settings to restart the operation.

c) Press and hold the Start/Stop button for 2 seconds during operation to stop. Press the Start/Stop button again to continue.

3. Temperature calibration

The temperature of the device is calibrated at the factory before release. If there is a difference between the actual temperature and the displayed temperature, you can calibrate the temperature as follows:

NOTE

1) the device has 2 temperature calibration points for reliability. Linear calibration values are 40°C and 100°C. Temperature accuracy is within $\pm 0,5^{\circ}\text{C}$ after calibration.

2) The ambient and block temperatures must not exceed 35°C during calibration.

Calibration methods:

a) turn on the device and it goes into standby mode. Check that the current temperature on the display does not exceed 35°C. If not, wait until it drops below 35°C.

b) place the olefin oil in one of the wells in the block, and then place the thermometer in that well (the thermometer accuracy should be 0.1°C and the thermometer ball should be fully immersed in the olefin oil in the block well). Insulation materials should be provided to allow adiabatic process to be carried out in the block (see Figure A below).


Figure A

NOTE: To ensure accurate calibration, read the actual temperature no earlier than 20 minutes after the temperature reaches the calibration points.


c) Press and simultaneously hold  and  to go to the temperature calibration screen, see figure on the right. The time display shows "Ad Jt", the temperature value represents the current temperature, program the automatic temperature control to 40°C .

When the temperature reaches 40°C , the decimal separator starts to flash. Wait for at least 20 minutes, the actual temperature of the thermometer is 39.6°C , press  or  to change the temperature value to 39.6. Press the Start/Stop button to confirm.

The program saves the value.

The temperature rises to 100°C automatically.

d) When the temperature reaches 100°C , the decimal separator flashes. Wait for at least 20 minutes, the actual temperature of the thermometer is 101.5°C , press  or  to change the temperature value to 101.5. Press the Start/Stop button to confirm.

e) After temperature calibration, press and simultaneously to   exit temperature calibration and return to standby mode.


Note

After the temperature calibration, the actual temperature of the block is displayed.


Simultaneously press and hold and during   calibration to exit the calibration procedure, the calibration is invalid.

4. Replacing the metal block


a) Remove the two screws securing the block to the heating panel using a screwdriver.


b) Attach the raiser to the center well of the block.

c) Pull out the raiser together with the block.

- d) Disconnect the raiser from the block, then attach the raiser to another block if necessary. Place the raiser and secure it with screws using a screwdriver.

Note: The procedure for replacing the block requires caution and attention.

Chapter 5: Error Analysis and Troubleshooting

No.	Failure	Possible Causes	Measures of elimination
1	No signal on display when switched on	No power	Check the power cable connection
		Fuse triggers	Replace the fuse (250 V, 3 A, $\varnothing 5 \times 20$)
		Faulty switch	Replace the switch
		Other	Contact the seller
2	The actual temperature and the displayed temperature are significantly different	Faulty sensor or loose module connections	Contact the seller
3	The display shows ERR and the alarm sounds	Faulty sensor or room temperature below zero	Contact the seller
4	The block does not heat up	Sensor failure	Contact the seller
		Damaged solid state relay	
		Faulty heater	
5	The button is not functioning	Button malfunction	Contact the seller

Chapter 6: Additional information for use in the territory of the Russian Federation

To dispose of the medical device:


In the Russian Federation, the centrifuge is to be disposed of at the end of its useful life in accordance with Sanpin 2.1.7.2790-10 "Sanitary and Epidemiological requirements for the Management of Medical waste".

Hazard class of medical wastes based on the characteristics of the morphological composition: Class B (epidemiologically hazardous wastes).

Annex 1:

Electrical wiring diagram of the DKT200-2A Laboratory Dry Bath Incubator

(The diagram below is for reference. All information is subject to change without prior notice).

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