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Instructions for Use
Series DWSBM / DWSRM / DWSRD

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Contents

For Your Safety and that of Your Patients.	3
Intended Use	4
Before starting.	4
Connect potentialfree contacts	5
Set time control.	5
Set the pressure switching points	6
Set pressure switch for tank pressure alarm	6
Operation	7
If the system does not work.	7
Möller Control System	8
Adjustment of base load	9
System pressure time value.	10
Advice of status in automatic running	11
Pressure.	12
Operation hour counter	13
Switch over from double to triple running (option).	14
Interchange over time control or base load.	15
Adjustment of time control.	16
Fault – Cause – Remedy	18
Overview Switch cabinet in general.	19
Maintenance.	19
Disposal	19
Technical Data	20
Performance characteristics	20
Operating characteristics	20
Order List	22

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Dräger Medical ANSY GmbH

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the pressure control panel requires full understanding and strict observation of these instructions for use.

The apparatus is only to be used for the purposes specified here.

Maintenance

The cabinet is maintenance-free.

A regular operational test should be carried out together with the maintenance of the plant.

Maintenance of the cabinet only by authorized personnel.

Only original parts are to be used for maintenance.

CAUTION!

Not for use in areas of explosion hazard.

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Safe connection with other electrical equipment

No electrical connections to equipment which is not listed in these instructions for use.

Liability for proper function or damage

The liability for proper function of the cabinet is irrevocably transferred to the owner or operator if the cabinet is improperly serviced or repaired by personnel not authorized by the manufacturer or if the cabinet is used in a manner not conforming to the intended use.

The manufacturer is not liable for damage which will be caused by non-compliance with the recommendations given above.

The warranty and liability provisions of the commercial terms of sale and delivery of the manufacturer are not modified by the recommendations given above.

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Intended Use

Compressed air control panel for fully-automatic operation and monitoring of compressed air systems with two or three compressors of the type DWSBM, DWSRM or DWSRD.

Operation is assured even if one compressor fails – or if two compressors fail in a system with three compressors.

CAUTION!

The equipment must not be used in the presence of inflammable gases – risk of fire!

Before starting

Starting the electrical part of the system

The following points must be noted when starting the electrical part of the system:

- Check all screws on devices, terminals and rails and retighten if necessary.
- Check that main voltage matches the operating voltage specified on the control panel; check primary voltage of control transformer.
- Set motor protection switch to rated motor current. (at ambient temperatures $> 20^{\circ}\text{C}$: current $+10\%$)
- Motors with direct starting:
Check motor voltage and switching. Fit jumpers on motor terminal board accordingly. Interchange two phases to change direction of rotation.
- Motors with star-delta starting:
Remove jumpers on motor terminal board, connect leads to control panel in phase sequence. ($U1/V1/W1-U2/V2/W2$ - motor to $U1/V1/W1-U2/V2/W2$ -control panel)
- To change direction of rotation: interchange phases after fuse or at $U1/V1/W1$ and $U2/V2/W2$.
- Set star-delta time-lag relay in accordance with motor starting characteristic, but not to less than 5 seconds.
- Set motor protection relay to 0,58 times the rated motor current.
- Check safety measures and observe the regulations of the local utilities.

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
Connect potentialfree contacts

Connect contacts to outgoing terminal strip X7 (1), see circuit diagram.

The contacts can be used to transmit malfunction messages.

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Set time control

For controlling of change-over for the compressors.

- Set running period at 24 hours for each compressor.
- Note the corresponding instructions for use.

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Set the pressure switching points

The pressure gets measured by an analogue pressure absorber and is processed in the Möller-PLC.

The switching points have to be set at the control panel. For a detailed description see point "Möller control system".

Set pressure switch for tank pressure alarm

Required tools: wrench (wrench size 10 mm)


- 1 Close ball valve (2) = turn handle to horizontal position.
- 2 Set required tank pressure (see schedule page 20) = turn drain valve (2) to the left to bleed measurement line.
- 3 Check pressure on pressure gauge. Set switching point:
- 4 Loosen nut on adjusting screw.
- 5 Set switching point = turn adjusting screw.
Set higher switching point = turn adjusting screw to the right.
Set lower switching point = turn adjusting screw to the left.
Check switching point by measuring continuity on potential-free contact.
- 6 Tighten nut on adjusting screw.
- 7 Open ball valve = turn handle to vertical position.

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- 1 Switch on power supply =
Set main switch (1) to "ON"

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
Function of the buttons:

- 1 = DEL: delete (not required)
- 2 = ALT: advice of status
- 3 = push buttons: select points of the menu, set pressure value
- 4 = OK: switch, store
- 5 = ESC: switch back, cancel

The Möller control system has the following functions for:


- Adjustment of push buttons
- Advice of status in automatic running
- Switch from double to triple running (option)

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Adjustment of base load

To set the corresponding base load, the service button has to be changed in the switch cabinet in order to indicate the switch points.

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
Basic load (switch on pressure compressor 1)
 Peak load (switch on pressure compressor 2)
 Reserve load (switch on pressure compressor 3 – only active in triple running)
 Switch-off pressure (switch off pressure of all compressor)

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System pressure Time value tank:

(time until system-pressure is reached 0.3 bar over the basic load value otherwise fault message indicated)

System pressure time value compressor 1-3:

(time until system-pressure of the selected compressor must be achieved otherwise fault message indicated. Signal over the external systempressure button.)

The corresponding pressure will be selected by the buttons > and < and activated by OK. The value is selected and must be activated by OK. The required value will be accepted by the button OK.

Further on the separate system-pressure overview can be switched on and off.


In order to leave the menus you have to push the button ESC and reset the service button in the switch cabinet.

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Advice of status in automatic running

In order to see the advice of status you have to switch by the buttons < or >.

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The following information are to read:

C1: off (compressor is switched off)

C1: Automatic off (compressor is in automatic function and has not received the switch-on point)

C1: Automatic on (compressor is in automatic function and is running)

C1: malfunction (compressor is not operable)

C3: inactive (only in double operation)

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
Pressure

Pressure: the actual pressure in bar


Switch off pressure: adjusted switch off pressure.

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You can read the following information:

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
Switch over from double to triple running (option)

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In order to adjust the corresponding operation, the service button has to be changed in the switch cabinet.

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The corresponding menu will be selected by the buttons > and <. The required value will be adjusted by the **arrow buttons**.

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
Interchange over time control or base load:

In order to adjust the corresponding function, you have to switch the service button in the switch cabinet.

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The corresponding menu will be selected by the buttons > and <.

The required value will be adjusted by the **arrow buttons**.

Arrow up: permutation over base load

Arrow down: permutation over time control

To leave the menu you have to switch the service button in the switch cabinet.

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Adjustment of time control

The internal weekly time control provides 4 disk spaces (A/B/C/D). In order to get to the parameter level you first have to activate the *-button and after that the **OK**-button. The line **parameter** will be selected by the arrow buttons and confirmed by **OK**.

The week-time component **HW01** will be addressed by the arrow buttons and confirmed by **OK**.

The following Parameter can be changed by the **arrow buttons** and by **OK**.

A/B/C/D	disk space
DY1	day-begin
DY2	day-end
ON	turn on time
OFF	turn-off

The input will be confirmed by **OK**. You get to the previous picture by **ESC**.

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Fault – Cause – Remedy

The following table is intended as a guide for locating and remedying the causes of any faults which may arise.
Only qualified specialists may carry out work for which the equipment must be opened.

Fault	Cause	Remedy
Motor does not start (manual mode)	Power supply not connected to the motor Motor damaged	Main switch "ON" Check main voltage Check circuit-breaker Check motor protection relay Check mode switch
Motor does not start (automatic mode)	Pressure wrongly set on PLC	Set pressostat switching points correctly
Pressure does not build up	2 Phases wrongly connected	Restore correct phase sequence Call for Service
Motor fuses are tripped	Motor current exceeds rated amperage of fuse	Check rated amperage of fuse in accordance with circuit diagram Check motor data in accordance with circuit. Check motor Call for Service
Compressors do not switch between basic load / peak load / reserve load	Failure of control system	Check control panel Call for Service
Red indicator lights (control voltage)	Failure of control voltage	Install replacement fuse Call for Service

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Overview Switch cabinet in general

See circuit-diagram, enclosed to the cabinet.

Maintenance

The following operational test should be made together with the maintenance of the plant:

- Check cutting-in / cutting-off points of the Möller control system.
- Check setting of pressure switch for tank pressure alarm.
- Check correct functioning of indicators (triple plants).

The corresponding instructions for use must be observed when servicing the compressor.

Disposal

At the end of the equipment's useful service life:

- Dispose of the equipment in accordance with national waste disposal regulations
- or
- return the equipment to an authorized disposal company for correct disposal.

Further information can be obtained from the local public authorities.

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90 39 412 - GA 6965.115 en © Dräger
© Dräger Medizintechnik GmbH
Dräger Medizintechnik GmbH