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
Ludwig-Kegel-Straße 16
79853 Lenzkirch, Germany

2019, rev. 1

1. MEDICAL DEVICE NAME AND COMPOSITION

ATMOS Chair Comfort with accessories, versions:

- ATMOS Chair Comfort basic;
- ATMOS Chair Comfort sync.
- Operating instructions

Accessories:

- Headrest.

2. INTENDED USE

- As a treatment chair for treatment and examination on patients who are seated or lying down in clinical areas, medical practices, company doctors, dentists (prophylaxis).
- As an appropriate and comfortable surface for supporting the patient during treatment/examination. As an aid for the persons performing the procedure (in order to improve ergonomics and quality during the work).

2.1 Intended users

Exclusively for personnel in the appropriate department who have received the appropriate training and instruction, i.e. doctors, assistants, nurses, etc.

3. CONTRAINDICATIONS

Not applicable.

ATMOS Chair Comfort is not intended for the following use:

- As a work chair.
- As an operating chair/table.
- As a substitute for a bed.

4. CLASSIFICATION OF MEDICAL DEVICE

	<i>ATMOS Chair Comfort</i>
Medical device type (according to Russian Nomenclature Classification of Medical Devices by Type)	139100
Classification acc. to Appendix IX EC Directive 93/42/EEC	I
Marking CE	CE


5. USAGE

The ATMOS Chair Comfort can be adjusted to suit the individual patient. The seat height can be adjusted hydraulically using a pump lever, and the height and angle of the headrest can be altered manually. The backrest can be tilted manually up to a fully horizontal position. The top section of the chair can be rotated, and the arm rests can be folded up individually.


5.1 Device overview:

- Front view

BASIC


SYNC

1. Headrest
2. Backrest
3. Arm rest
4. Backrest release lever
5. Equipotential bonding connection
6. Base
7. Foot support
8. Leg rest (Sync version)
9. Seating surface

- Rear view

1. Backrest release lever
2. Seat height pump lever
3. Headrest locking screw (height)
4. Headrest locking screw (angle)


	Basic	Sync
Hydraulically adjustable seat height	x	x
Manually adjustable backrest	x	x
Synchronised leg rest		x

5.2 Connecting the device

1. Clean and disinfect the product prior to first use.
2. Connect the equipotential bonding if required.
3. The equipotential bonding must be connected if required by the spatial environment in which the product is going to be used. Refer to the specifications of standard IEC EN 60601-1.
4. Perform a function check.

5.3 Positioning the patient

Ensure that the patient sits in the middle of the seating surface. The patient must not be allowed to sit on the arm rest, backrest or leg rest. Do not allow patients to place their full weight on the foot support.

Risk of tipping due to overloading on one side.

1. Move the backrest {1} to an upright position.
2. Move the chair {2} to its bottom position.
3. Do not let the patient sit down or stand up until this is done.

5.4 Adjusting the Chair


CAUTION

Risk of crushing by moving parts.

There is a risk of the patient's fingers being crushed:

- The patient's arms must remain on the arm rest at all times while you are adjusting the patient chair.

5.4.1 Adjusting the arm rests

Folding the arm rests up:


- Fold the arm rest {1} up in a controlled manner.

Folding the arm rests down:

The arm rests can be damaged if allowed to drop down under their own weight.

- Fold the arm rest {1} down in a controlled manner.


5.4.2 Adjusting the foot support (Sync version)

Unfolding the foot support:

- Pull the handle and fold the foot support down.

Folding the foot support away:

- Fold the foot support up.


5.4.3 Adjusting the backrest

Mechanical:

- Push the backrest release lever down, pushing or pulling on the backrest at the same time.
- Let go of the backrest release lever as soon as the backrest reaches the required position.

» The backrest is fixed in place.


5.4.4 Adjusting the headrest

Angle:

- Undo the headrest locking screw (angle) {1} and move the headrest to the required position.
- Tighten the headrest locking screw (angle) {1}.


» The headrest is fixed in place at the current angle.

Height:

- Undo the headrest locking screw (height) {1} and move the headrest to the required position.
- Tighten the headrest locking screw (height) {1}.

» The headrest is fixed in place at the current height.

5.4.5 Adjusting the seat height

Upwards

- Pump the pump lever {1} to raise the seat height.

Downwards


- Push the pump lever {1} as far as it will go.

5.4.6 Turning the Chair

- Turn the top section of the chair.
- To fix the patient chair in place, push the pump lever {1} up using your foot.

5.5 Adjusting the shock position

- Adjust the angle of the backrest up to a fully horizontal position.

5.6 Troubleshooting

The device has been subjected to a thorough quality control in the factory. If a fault should occur despite this care, you may be able to resolve it yourself. See the table below.

6. TECHNICAL SPECIFICATION

Seat height	56 cm (rigid seat upholstery) $\pm$ 20 mm
Max. lifting capacity	150 kg (safe working load)
Height of backrest	Approx. 600 mm (without headrest) $\pm$ 50 mm Approx. 780 mm (headrest in bottom position) $\pm$ 50 mm Approx. 980 mm (headrest in upper position) $\pm$ 50 mm
Rotation	360°
Backrest inclination	+7° up to -90°
Ambient conditions for transport/storage: • Temperature • Humidity without condensation • Air pressure	-10...+50°C 20...80% without condensation 70,0...106,0 kPa (525...800 mm Hg)
Ambient conditions for operation: • Temperature • Humidity without condensation • Air pressure	+10...+35°C 20...80% without condensation 70,0...106,0 kPa (525...800 mm Hg)
Maximum operational altitude	$\leq$ 2000 m
Dimensions (LxWxH): • Basic version • Sync version	820 x 670 x 1280 mm $\pm$ 30 mm 700 x 670 x 1280 mm $\pm$ 30 mm
Weight: • Basic version • Sync version	63 kg $\pm$ 5kg 59 kg $\pm$ 5kg
Period tests	Recommended: Test every 24 months.

7. SAFETY REQUIREMENTS

Please read and pay attention to the safety instructions prior to using the product.

General safety information

- Installation and start up may need to be carried out by the operating company's technical personnel, or those of an authorized specialist dealer.
- Never leave elderly persons, children, injured persons, persons with disabilities or other patients who require assistance alone on the patient chair.
- The Declaration of Conformity shall be rendered null and void if any changes are made to the product by the customer or any third parties. This includes, but is not limited to: conversions of all types, the use of third party accessories, and the removal of warning and information signs.
- If any signs are missing from your device or no longer legible, contact ATMOS and replace the signs in question.

Danger for users, patients and third parties

- Keep the device fully functional at all times.
- Malfunctions could cause injury to you and your patients:
 - Prior to each use, check whether the device is damaged. Do not operate the device if you notice any damage.
 - Perform a function check at least once a week, and also at any other time if you have any concerns with regard to the safety of the device.
 - Observe the specifications regarding period tests.
 - Only use original accessories and original spare parts from ATMOS.
 - Assembly, new settings, modifications, extensions and repairs may only be carried out by persons who are authorized by ATMOS.
- Risk of tipping due to overloading on one side:
 - The patient chair may tip over, injuring you or your patient.
 - Ensure that the patient sits in the middle of the seating surface. The patient must not be allowed to sit on the arm rest, backrest or leg rest.
 - Move the patient chair to the bottom position and the backrest to an upright position before the patient sits down or stands up.
- Risk of infection due to patient secretion on the device!
- Risk of diseases being transmitted:
 - Clean and disinfect the device after every use.
 - Clean and disinfect the device in accordance with the operating instructions.
- Only a fully functional product meets the safety requirements of users, patients and third parties. As such, observe the following instructions regarding the product.

Avoid damage to the device

- Damage caused by rolling the patient chair in contravention of intended use.
- Do not allow patients to place their full weight on the foot support.

8. CLEANING AND DISINFECTION

The device is provided not sterile. Does not require sterilization before use.

The device is reusable.

It is recommended to document always all maintenance work and part replacements in writing.

It is the responsibility of the user to ensure that the required results for cleaning and disinfection are adhered to. Validation and routine monitoring of the procedure will generally be necessary.

CAUTION

Risk of infection due to secretion on the product and accessories.

Risk of diseases being transmitted.

- Always wear disposable gloves during all cleaning operations.
- Clean the product after each use.
- Clean and disinfect the device in accordance with the operating instructions.

8.1 *Cleaning and disinfecting the surfaces*

Clean the device after each use. Disinfect the patient chair if it becomes contaminated by secretion or other substances, at least once per day is necessary.

1. Move the chair to its bottom position.
2. Moisten a clean, non linting cloth with lukewarm, soapy water or a conventional plastic cleaning agent.
3. Clean the surfaces with the damp cloth.
4. Disinfect the surfaces with a clean, non-linting cloth and a recommended disinfectant.
5. Leave the device to dry before using it again.

8.2 *Cleaning and disinfecting the upholstery*

Always remove any oil, grease or sweat from the upholstery immediately. Clean the patient chair after each use. Disinfect the patient chair if it becomes contaminated by secretion or other substances, at least once per day is necessary.

1. Move the chair to its bottom position.
2. Use warm, mild, soapy water.
3. Clean the upholstery using a damp micro-fibre cloth or a soft brush.
4. Disinfect the upholstery with a clean, non-linting cloth and a recommended disinfectant. Do not allow the disinfectant to soak in for longer than necessary.

8.3 *Recommended disinfectants*

ATTENTION

Unsuitable cleaning agents and disinfectants.

Risk of damage to the device surface, corrosion damage and stress cracks.

Do not use cleaning agents or disinfectants that contain the following ingredients:

- Irritants,
- Substances that may damage plastics or steel components or dissolve lubricants.

ATTENTION

Unsuitable cleaning agents and disinfectants.

There is a risk of damage to the upholstery.

Do not use cleaning agents or disinfectants that contain the following ingredients for the upholstery:

- Solvents or chlorides,
- Polish or wax polish,
- Chemical cleaning agents,

- Oils, greases or alcohol.

Use cleaning agents with a pH value of between 6 and 8.

Observe the manufacturer's specifications with regard to concentration and soaking-in time.

Use a suitable surface disinfectant.

9. TRANSPORT CONDITIONS

The device should be transported only in a shipping carton that is padded and offers sufficient protection.

If damage occurs during transport:

1. Document and report the transport damage.
2. Send the device to ATMOS. Check whether all the fixing screws are fully tightened.


Ambient conditions for transportation:

- Temperature: -10...+50°C
- Relative humidity: 20...95% without condensation
- Air pressure: 70,0...106,0 kPa (525...800 mm Hg)

10. PACKAGING

Packaging: Carton.

Chair is covered with foil and there is a box with headrest inside.

11. LABELING

ATMOS Chair Comfort models labeling example:

	 Greiner GmbH Wettstrasse 1 D-74385 Pleidelsheim	
REF	6612085 HYD,MAN,KO,FU 535.5000.0 ATMOS Chair Comfort basic	
	2000-00-00 For: ATMOS MedizinTechnik GmbH & Co. KG	
V in:	V ~, Hz I in: max.	A IP
OP:		
 : 150kg		Polster / Cushion:
UDI:		Q-Control:
SN		
20194539700050001		

ATMOS Chair Comfort models labeling example (version for the Russian Federation):

	 Greiner GmbH Wettstrasse 1 D-74385 Pleidelsheim	
REF	6612085 HYD,MAN,KO,FU 535.5000.0 ATMOS Chair Comfort basic	
	2000-00-00 For: ATMOS MedizinTechnik GmbH & Co. KG	
V in:	V ~, Hz I in: max.	A IP
OP:		
   : 150kg		Polster / Cushion:
UDI:		Q-Control:
SN		
20194539700050001		

Shipping labeling contains following handling symbols

- Keep dry;
- Top;
- Maximum load constraint.

Shipping packaging contains the following information:

- Manufacturer name;
- Manufacturer contact information (address, phone, fax, web-site, e-mail);
- Device number according to manufacturer's numbering system;
- Production date (year, month);
- "Please refer to instructions prior to use" symbol.

On the device and type plate

	Follow operating instructions (blue)
	Warning; take extra care to observe
	This product complies with the relevant requirements of the EU Directives
	Manufacturer
	Manufacturing date
SN	Serial number
REF	Order number
	Professional disposal
	Maximum load

12. STORAGE CONDITIONS

Ambient conditions for storage:

- Temperature: -10...+50°C
- Relative humidity: 20...95% without condensation
- Air pressure: 70,0...106,0 kPa (525...800 mm Hg)

13. DISPOSAL

13.1 Packaging

Please recycle the packing.

13.2 ATMOS Chair Comfort

The device does not adversely impact the environment.

The ATMOS Chair Comfort must be disposed of in compliance with the applicable regulations.

The bin symbol on the product or in the manual indicates that this product must not be disposed of with household waste. ATMOS MedizinTechnik GmbH & Co. KG has a waste management system. This also includes free return of ATMOS Chair Comfort.

Arrangements for the return of all unusable equipment to ATMOS for correct removal and environmentally friendly disposal should be discussed with the local distributor or dealer of the country in which the equipment was purchased.

Do not dispose of the product together with household waste.

The product contains oil. Oil must be disposed of professionally, and must not be allowed to penetrate into the soil or water.

1. Clean and disinfect the device.
2. Dispose of the product professionally and according to the country specific laws and regulations.

The housing is fully recyclable. Observe to the country-specific laws and regulations.

13.3 Hydraulic unit and gas spring

Observe the following instructions if you are unable to return the product to ATMOS or your specialist dealer for disposal.

CAUTION

Risk of injury and crushing due to high internal pressure in the gas springs and sudden movements of the piston rods.

There is a risk of injury to the eyes, face and hands:

- Wear suitable protective equipment.
 - Cover the vent.
 - Take extra care when operating the release mechanism after the gas springs have been removed.
1. Depressurize the gas spring. Always observe the gas spring manufacturer's specifications. The manufacturer is named on the type plate.
 2. Drain the oil used for the gas spring.
 3. Drain the oil out of the hydraulic unit.
 4. Dispose of the oil properly and professionally.
 5. Recycle the hydraulic unit and gas spring.

14. WARRANTY POLICY

Operational period is not less than 5 years.

ATMOS MedizinTechnik GmbH & Co. KG provides a warranty of one year from the date of purchase on all ATMOS Chair Comfort versions, provided that no other warranty terms have been agreed in the contract of sale. Device and material defects detected during this period will be repaired by ATMOS MedizinTechnik GmbH & Co. KG.

Notwithstanding any alternative terms of warranty that may have been specified in the contract of sale, the conditions of warranty specified here exclude and replace all other verbal warranties or warranties in writing that were explicitly agreed or tacitly implied.

The conditions of warranty specified here apply exclusively to defective devices. In no case will ATMOS MedizinTechnik GmbH & Co. KG accept liability for accidental or indirect losses (including loss of profit), bodily injury, damages or costs directly or indirectly resulting from use of the ATMOS Chair Comfort.

Complaints in respect of device operation can be submitted to the manufacturer ATMOS MedizinTechnik GmbH & Co. KG.

15. MAINTENANCE AND SERVICE

Maintenance, repairs and period tests may only be carried out by persons who have the appropriate technical knowledge and are familiar with the product. The person must possess the necessary test devices and original spare parts required to carry out these measures.

ATMOS recommends: Work should be carried out by an authorized ATMOS service partner. This ensures that the repairs and testing are carried out professionally, that original spare parts are used and that warranty claims remain unaffected.

15.1 Period tests

Please comply with the country-specific guidelines regarding regular testing.

ATMOS recommends performing a test every 24 months.

15.2 Function check

Perform a function check at least once a week, and also at any other time if you have any concerns with regard to the safety of the device:

1. Check whether the chair is damaged.
2. Move the chair to its bottom position, paying attention to any unusual noises.
3. Check that all the moving parts can be adjusted without any problems:
 - Backrest
 - Arm rests
 - Headrest
 - Foot support
4. Check whether the chair can be turned without any problems.

15.3 Lubricating moving parts

Lubricate the moving parts at regular intervals to keep them fully functional:

1. Remove any dirt from the moving parts.
2. Lubricate the axes and joints using commercially available oil or grease.
3. Check whether all the fixing screws are fully tightened.

15.4 Sending in the device

1. Clean and disinfect the product and accessories in accordance with the instructions.
2. Place used accessories with the product.
3. Fill in the "Delivery complaint/return shipment" form and the corresponding decontamination certificate.
4. This form is enclosed with each delivery and can be found at www.atmosmed.com
5. The device must be well padded and packed in suitable packaging.
6. Place the "Delivery complaint/return shipment" form and the corresponding decontamination certificate in an envelope.
7. Affix the envelope to the outside of the package.
8. Send the product to ATMOS or your dealer.


16. COMPLAINTS

ATMOS MedizinTechnik GmbH & Co. KG
Ludwig- Kegel-Straße 16
79853 Lenzkirch, Germany
Phone: 0049 7653 689- 0
Fax: 0049 7653 689-190
e-mail: atmos@atmosmed.de

17. AUTHORIZED REPRESENTATIVE IN RUSSIA

ATMOS Medical LLC
105005, Russia, Moscow, Poslannikov pereulok, 5, bld. 8
Phone: 8-800-707-08-94
e-mail: atmosmed@atmosmed.ru
Web: www.atmos-med.ru

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