	INSTRUCTIONS FOR USE Spike Infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

Versions:

1. **A-6425/118:** Infusion-aspiration cannula with air filter RoweSpike II
2. **A-6226/118:** Infusion-aspiration cannula with air and hydrophilic filter RoweSpike II 5 µm

Area of application:

Spike Infusion-aspiration in versions are applied for withdrawal of liquids from flasks and bags, drugs' mixing and solutions' preparation in process of infusion therapy.

INDICATIONS:

Spike infusion-aspiration cannulas of different versions are indicated for preparation, mixing and multiple subsequent withdrawal of sterile solutions and medicines, with possibility of filtering and binding of particles and inclusions, aerosols, air bubbles from solutions, prevention of infections and reducing the risk of in-hospital contamination.

CONTRAINDICATION:

Relative contraindications are:

- Allergic reaction to any materials and components that make up the device.
- Active disagreement of patient to use the device.

There are no absolute contraindications.

SIDE EFFECTS AND COMPLICATIONS:

- Allergic reaction to materials and components of the device
- No other has been identified.

APPLICATION:

- can be used with flasks 3 – 1000 ml.
- availability of attached cap and air filter with small pore size provides maximum protection from contamination and toxic aerosols, extracted during solution of dry substances
- medical device should be used by nursing medical personnel for preparation and multiple withdrawal of solutions in special procedural premises and pharmacies of medical hospitals.

INSTRUCTION FOR USE:

Spike infusion-aspiration cannulas of different versions:

- Check the integrity of the package and the expiry date of the device;
- Carefully open the package, remove the device, preserving its sterility;
- Connect to a container with solution or container for solution preparation;
- Make sure that the cap is tightly closed, the container is ready to be used;
- Connect a syringe or an extension set, infusion system to the Luer-Lock connector;
- Start preparing the required medicinal solution, or collect the ready solution for use from the patient;
- After use, treat and dispose of the device in the manner prescribed.

Safety requirements and precautions:

- Must be used only by qualified and trained personnel
- Before using the device, it is necessary to read the attached instruction
- Single-use devices must not be cleaned and re-sterilized
- Must not be used if package integrity is violated
- Devices must be used immediately after opening sterile packaging
- After use, the device should be properly disposed of and handled

Operation conditions

Requirements for environmental conditions during use	Value
Air temperature, °C	0 to +25
Air humidity (relative), without condensate, %	50 to 60
Atmospheric pressure, kPa	60 to 106

Transportation and storage condition.

- The following requirements must be observed during transportation and storage:
- Protection from direct sunlight
- Transportation and storage in shipping multi-container, which protects the device from dust, mechanical damage

«RoweMed AG – Medical 4 Life» - Germany
 Тел. +49 (0)3871 451 280 Факс +49 (0)3871 451 282 - <http://www.rowemed.de/>
 Importer: M.Schilling Medical Products RUS, Ltd., Russia, Moscow, 115280, Leninskaya sloboda, d. 26,
 premise XXXII, office 3 Тел.: +7(495) 662-72-95

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CE 0482

Rev. 1.0 / 08.05.2017

RoweMed Medical 4 Life	INSTRUCTIONS FOR USE Spike infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

- Isolation from corrosive chemicals
- Protection from mechanical impact (compression, bending, twisting, stretching, etc.)
- Distance to heating appliances at least 1m
- Absence of open fire

Transportation

Requirements for environmental conditions during transportation	Value
Air temperature, °C	-10 to +25
Air humidity (relative), without condensate, %	<60
Atmospheric pressure, kPa	86 to 106

Storage

Requirements for environmental conditions during storage	Value
Air temperature, °C	0 to +25
Air humidity (relative), without condensate, %	<60
Atmospheric pressure, kPa	86 to 106

Warranty

LLC M.Schilling Medical Products Rus, as an official and authorized representative of RoweMed AG – Medical 4 Life, Germany, guarantees proper resolution of all the issues associated with possible claims as regards physical damage, assortment, incompleteness and quality of the products manufactured by the company in accordance with the legislation of the Russian Federation and the Law "On protection of consumer rights" (No. 2300-1 of 07/02/1992). The company is not responsible for the use of manufactured medical devices outside of intended purpose and not in accordance with the instructions for use.

Warranty period - 5 years from the date of sterilization (which is considered to be manufacturing date)

ATTENTION:

- alcohol and lipid resistant
- latex - free
- PVC, DEHP-free
- shelf-life – 5 years from the date of sterilization (which is considered to be manufacturing date)

Procedure and conditions for disposal

When disposing of the medical device, act in accordance with all the requirements of local authorities and the policy of this healthcare institution for the prevention of environmental contamination (including infectious).

If the packaging is damaged, the device must be disposed of in accordance with the rules for handling Class A medical waste.

Used device must be disposed of in accordance with the rules for handling Class B medical waste.

"SINGLE USE"

The product may be used in sterile condition, one-time. Cleaning or re-sterilization are excluded. Re-using single-use products bears potential risks for both, patient and user. Impairment of material properties and inadequate sterility associated with such re-uses may lead to illness, damage or death of the patient.

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APPENDIX TO INSTRUCTION FOR USE

Information about manufacturer of medical device

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http://www.rowemed.de

Adders of the manufacturing site

RoweMed AG – Medical 4 Life, Juri-Gagarin-Ring 4 – 19370 Parchim, Germany

Description of the principles that underline the operation of the medical device and their aspects

Infusion therapy is a treatment method which is based on intravenous or subcutaneous injection of various medicinal solutions and products for the purpose of normalizing water-electrolyte, acid-alkaline balance of the body and correcting (or preventing) pathological losses of the body. All kinds of single-use consumables are used for infusion therapy, including:

Spike infusion-aspiration cannulas of different versions are used for the preparation, multipoint sampling or mixing solutions and medications from small- and large-volume, multi-dose vials, polymer bags for the purpose of using ready-made solutions and medications in injections and infusions in patients.

Ability and ways of integration with other medical devices

Spike infusion-aspiration cannulas of different versions can be used both separately and form various types of closed infusion systems that work with great efficiency and solve all the main issues of infusion therapy. The devices can be used together with various devices for intravenous injection, with intermittent or prolonged infusions, for connecting or disconnecting with syringes, infusion sets, infusion plugs, infusion systems.

Materials

Infusion-aspiration cannula with air filter RoweSpike II:

Cannula body – ABS plastic; white dye

Spike needle – ABS plastic; white dye

Membrane – polytetrafluoroethylene

Cover – polypropylene, green dye

Protective cap – high-density polyethylene

Infusion-aspiration cannula with air and hydrophilic filter RoweSpike II 5 µm:

Cannula body – ABS plastic; white dye

Spike needle – ABS plastic; white dye

Membrane – polytetrafluoroethylene

Hydrophilic filter – polyethylene terephthalate

Cover – polypropylene, blue dye

Protective cap – high-density polyethylene

List of materials of animal and (or) human origin

Not applicable for this MD

List of medicinal products for human use contained in the MD

This MD does not contain medicinal products

Possibility and special aspects of the use of the medical device in people with medical devices implanted in the body, pregnant women, women during breastfeeding, children, adults with chronic diseases

Spike infusion-aspiration cannulas of different versions can be used in people with medical devices implanted in the body, pregnant women, women during breastfeeding, children, adults with chronic diseases. No special aspects of the use of the medical devices have been identified.

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CE 0482

Rev. 1.0 / 08.05.2017

RoweMed Medical 4 Life	INSTRUCTIONS FOR USE Spike infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

Information on the potential effect of the use of the medical device on the ability to operate vehicles, mechanisms
Not applicable for this MD

Symbols
Symbols used on the label of sterile packaging of the device

Symbol	Explanation	Symbol	Explanation
	Catalog number		European Union certification mark
	Lot code		Caution! Refer to the instruction for use
	Use before		Refer to the instruction for use
	Re-use prohibited		Does not contain natural latex
	Do not use if the packaging is damaged		Does not contain phthalates
	Keep away from sunlight. The cargo should be protected from sunlight		Nonpyrogenic
	Keep away from moisture. It is necessary to protect the cargo from moisture		Non-toxic
	Do not re-sterilize		Mark of conformity in Russia
	Sterilization with ethylene oxide		Address of manufacture

RoweMed Medical 4 Life	INSTRUCTIONS FOR USE Spike infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

Symbols used on the label of shipping packaging

Symbol	Explanation
	Catalog number
	Lot code
	Use before
	Re-use prohibited
	Do not use if the packaging is damaged
	Keep away from sunlight. The cargo should be protected from sunlight
	Keep away from moisture. It is necessary to protect the cargo from moisture
	Do not re-sterilize
	Sterilization with ethylene oxide
	Address of manufacture

Symbol	Explanation
	Number of devices per package
	European Union certification mark
	Caution! Refer to the instruction for use
	Refer to the instruction for use
	Does not contain natural latex
	Does not contain phthalates
nonpyrogenic	Nonpyrogenic
non-toxic	Non-toxic
	Mark of conformity in Russia

«RoweMed AG – Medical 4 Life» - Germany
Тел. +49 (0)3871 451 280 Факс +49 (0)3871 451 282 - <http://www.rowemed.de/>
Importer: M.Schilling Medical Products RUS, Ltd., Russia, Moscow, 115280, Leninskaya sloboda, d. 26,
premise XXXII, office 3 Тел.: +7(495) 662-72-95

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Rev.1.0 / 08.05.2017

20110707

RoweMed Medical 4 Life	INSTRUCTIONS FOR USE Spike infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

Parameters and characteristics*

*Deviations of characteristics within $\pm 5\%$ range are allowed, unless otherwise specified

Spike infusion-aspiration cannulas of different versions:

1. Infusion-aspiration cannula with air filter RoweSpike II

Parameters	
Length, mm	62
Radius of the wide part, mm	28
Weight, g	4.3
Cannula (material, type)	Plastic
Needle length, mm	21
Needle's working length, mm	21
Needle's internal diameter:	
- diameter of the hole for liquid, mm	1.3
- diameter of the hole for air, mm	1.25
Outer diameter, mm	4
Needle sharpening angle	26°
Luer-Lock connectors	"Female" type
Working pressure, MPa	≤ 0.1
Filled volume, mL	0.25
Pore size of antibacterial hydrophobic ventilation filter, μm	0.1
Pore size of filter from soluble particles (hydrophilic filter), μm	-
Length of the protective cap, mm	22
Outer diameter of the protective cap, mm	6
Internal diameter of the protective cap, mm	4
Nominal density of the material, ABS plastic kg/m^3	1040

2. Infusion-aspiration cannula with air and hydrophilic filter RoweSpike II 5 μm

Parameters	
Length, mm	62
Radius of the wide part, mm	28
Weight, g	4.3
Cannula (material, type)	Plastic
Needle length, mm	21
Needle's working length, mm	21
Needle's internal diameter:	
- diameter of the hole for liquid, mm	1.3
- diameter of the hole for air, mm	1.25
Outer diameter, mm	4
Needle sharpening angle	26°
Luer-Lock connectors	"Female" type
Working pressure, MPa	≤ 0.1
Filled volume, mL	0.25
Pore size of antibacterial hydrophobic ventilation filter, μm	0.1
Pore size of filter from soluble particles (hydrophilic filter), μm	5
Length of the protective cap, mm	22
Outer diameter of the protective cap, mm	6
Internal diameter of the protective cap, mm	4
Nominal density of the material, ABS plastic, kg/m^3	1040

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Rev. 1.0 / 08.05.2017

20110707

RoweMed Medical 4 Life	INSTRUCTIONS FOR USE Spike infusion-aspiration cannulas	Product sterile
		Revision: 1.0 / 13.09.2018

Applicable national standards

The list of applicable national standards is published on the web-site of importer (Empowered representative). Also it can be provided by importer to the customer in paper or electronic version by request.

Requirements for installation

Spike infusion-aspiration cannulas of different versions

- Connect to a container with solution or container for solution preparation.
- Make sure that the cap is tightly closed, the container is ready to be used.

For more information, see clause "Usage".

Basic composition of package

Spike infusion-aspiration cannulas of different versions:

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One "Instruction for use" in paper format is inserted into multiple (consumer) packaging.

Information on maintenance, information on verification

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I, the undersigned notary,

Regina N i e m a n n

with the official residence in Schwerin,
in the business premises 19061 Schwerin, Weinbergstr. 16

authenticated publicly as real and as today before me personally executed the
Name signature of

Mr. Dr. Dirk **F o r b e r g e r**,
born 03. September 1971,
business resident in 19370 Parchim, Juri-Gagarin-Ring 4,
me, the notary, known by person,

acting here not in his own name, but as a sole representative entitled member of the
board for the in the commercial register of the district court Schwerin under
HR B 7789 registered corporation in company

RoweMed AG - Medical 4 Life


based in Parchim

business address: 19370 Parchim, Yuri Gagarin Ring 4.

Due to this on 03 April 2019, the retrieval of data from the electronically controlled
commercial register of the district court Schwerin, HR B 7789, I, the attorney notary,
hereby certify that Dr. Dirk **F o r b e r g e r** a single member of Executive Board of
RoweMed AG - Medical 4 Life with registered office in Parchim.

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