NINDEX

The software for monitoring of anesthesia depth

User Guide for Neurosoft Devices

Controles S.A y Dr. D. Cibils
Av. Rivera 3314
Montevideo URUGUAY

Tel. +5982 6220651
Fax +5982 6222048
www.controles.com

The authorized representative of the manufacturer in Russia: Neurosoft Ltd., Ivanovo

NINDEX software to monitor the anesthesia depth includes:

- NINDEX software to monitor the anesthesia depth;
- NINDEXMN software to review the monitoring data;
- Rockey4ND hardware dongle;
- User Manual.

NINDEX Software version: 3.4.x
Document version: April 3, 2012

LIABILITY LIMITATIONS / DISCLAIMER

The Software is provided in an “as is” condition. The User is entirely responsible for any risk regarding the quality and performance of the Software itself. The Author makes no other warranty regarding the Software and specifically (a) waives any warranty related to the fitness for any particular purpose, (b) denies that the Software will operate without interruptions or errors and (c) denies that errors of the Software (if there are any) will be corrected.

The Author assumes no liability nor grants any other express or implicit warranty in relation with any harmful consequences (damages, loss of profit, etc.) that might arise against the User or third parties as a result of the Software use.

CONSIDERATIONS FOR NINDEX MONITORING

Clinical judgment should always be used when interpreting the NINDEX values. ***Reliance on the NINDEX alone for anesthetic drug delivery is not recommended.*** As at any monitoring, artifacts or low signal quality may cause improper NINDEX values. Potential artifact sources are poor skin-to-electrode contact, head and body motion, eye blinking, muscle activity, improper electrode placement or electrical interference. NINDEX values should be interpreted with caution in patients with revealed neurological disorders or patients taking psychoactive drugs, due to limited clinical experience in those applications.

TABLE OF CONTENTS

TABLE OF CONTENTS.....	4
INTRODUCING THE NINDEX SOFTWARE	6
Indications for use	6
General information.....	6
SAFETY PRECAUTIONS.....	9
Cautions and notes	9
CAUTIONS	9
SOFTWARE SETUP	10
Setup checklist.....	10
Software setup	10
System requirements.....	10
Setup	11
Starting and ending the software	11
Screen size setup.....	12
External program communication setup.....	12
Setup of data exchange with EEG system using Neuron-Spectrum.NET software	12
Setting of connection with Neuron-Spectrum.NET external program	13
Setup of data exchange with EEG system using Neuro-IOM.NET software	15
Neuro-IOM.NET connection setup.....	17
USING THE SOFTWARE.....	18
Electrode placement	18
Electrode to EEG acquisition device connections	21
Neuro-MEP-Micro (Neuro-EMG-MS).....	21
Neuron-Spectrum-1	21
Monitoring start	22
Monitoring display	22
The numeric region.....	23
Monitoring general data region	24

Graphics region	25
Buttons region	26
The Menu Dialog	26
The Settings Dialog	28
The Event Dialog	29
The Event Edition Dialog	30
The Patient Data Dialog	32
Accelerator keys	32
Monitoring end	33
SOFTWARE FEATURES	34
Anesthesia depth monitoring	34
Checking of skin-to-electrode contact	35
Alarms	35
Data storage to disk	36
Monitoring data	36
Software settings	36
Disk storage managing	36
SOFTWARE UPDATE	37
MONITORING DATA REVIEW	38
Introduction	38
Opening the monitoring data for review	38
The monitoring data review display	39

INTRODUCING THE NINDEX SOFTWARE

Indications for use

NINDEX software is meant for a continuous evaluation of the brain state basing on EEG processing, intended for use in the operating room, intensive care unit and clinical research centers.

The result of the EEG processing is a number called NINDEX, which may be used as an aid in monitoring the effects of anesthetic drugs in male and female patients above the age of 2 years, having no revealed neurological disorder and not taking psychoactive drugs.

General information

NINDEX software is intended to monitor the anesthesia depth basing on EEG processing. It continuously calculates a number called NINDEX, which correlates with the depth of patient anesthesia.

NINDEX software runs under control of Windows XP SP2, Windows XP SP3, Windows Vista or Windows 7 operating systems.

NINDEX software provides a simple user friendly Graphical User Interface (GUI) which consists of:

- Current NINDEX number
- Current BSR number
- NINDEX trend plot
- Raw EEG waveform plot
- Signal Quality (SQ) indicator
- Current time and date
- Monitoring start time

- Patient identifier
- Electrode impedance levels (optional)
- BSR trend plot (optional).

NINDEX software GUI has two buttons for starting and ending a case, and for displaying the program Menu Dialog. The Menu Dialog allows to mark user defined events, to configure the GUI appearance and to configure sound alarms.

NINDEX software generates sound alarms if current NINDEX value is greater than a high threshold or lower than a low threshold. Alarms may be silenced by the user.

NINDEX software automatically stores monitoring data to disk. NINDEXMN software for Windows, which can be downloaded from www.controles.com website, allows reviewing recorded trend data and event markers.

NINDEX software may be used with one of the following EEG acquisition devices, manufactured by Neurosoft Ltd (www.neurosoft.ru):

- Neurophysiological systems for EEG, EP and EMG on the following variants (Neurosoft Ltd. (www.neurosoft.ru)):
 - Neuron-Spectrum-1 8-channel digital EEG system with 1 polychannel for ECG recording and breath channel;
 - Neuron-Spectrum-2 16-channel digital EEG system with 1 polychannel for ECG/EMG/EOG/GSR recording, breath channel, auditory and visual stimulation channels, pattern stimulation channel;
 - Neuron-Spectrum-3 19-channel digital EEG system with 1 polychannel for ECG/EMG/EOG/GSR recording, breath channel, auditory and visual stimulation channels, pattern stimulation channel;
 - Neuron-Spectrum-4 21-channel digital EEG system with 1 polychannel for ECG/EMG/EOG/GSR recording, breath channel, auditory and visual stimulation channels, pattern stimulation channel;

- Neuron-Spectrum-4/P 21-channel digital EEG system with 1 polychannel for ECG/EMG/EOG/GSR recording, breath channel, auditory and visual stimulation channels, pattern stimulation channel;
 - Neuron-Spectrum-4/EP 21-channel digital EEG system with 1 polychannel for ECG/EMG/EOG/GSR recording, breath channel, auditory and visual stimulation channels, pattern stimulation channel;
- 2-channel portable digital EMG and EP system with a built-in miniature dedicated keyboard Neuro-MEP-Micro (Neuro-EMG-MS delivery set), (Neurosoft Ltd. (www.neurosoft.ru)).

NINDEX software is intended for use with external programs. In this mode, the external program transmits the EEG signal to NINDEX software, which transmits calculated NINDEX, BS and SQ values back to the external program. It may be used with one of the following external programs by Neurosoft Ltd. (www.neurosoft.ru):

- Neuron-Spectrum.NET;
- Neuro-IOM.NET.

SAFETY PRECAUTIONS

Cautions and notes

The terms **caution** and **note** have specific meanings in this manual:

- A CAUTION is a piece of advice allowing to prevent actions or situations that could lead to data inaccuracy.
- A NOTE provides useful information regarding a function or procedure.

CAUTIONS

The use of accessories other than those specified may result in increased electrical interference on the EEG, decreasing the reliability of data provided by the software.

Proper electrode placement is essential for the reliability of the data provided by the software.

Good skin-to-electrode contact is essential for the reliability of the data provided by the software.

Strictly follow the electrode placement procedure described in this manual.

To proper patch to skin adhesion, avoid wetting the electrode self-adhesive part when placing electrodes.

The data provided by the software is not considered reliable if the software indicates a bad electrode contact warning.

The software requires the c:\eegdata folder for proper operation. If this folder is accidentally deleted then the monitoring and configuration data will be lost.

In order not to lose monitoring data by means of automatic deletion, back up the data after maximum 100 monitoring hours.

SOFTWARE SETUP

Caution:

Carefully read this section before software setup.

This section provides NINDEX software setup instructions. It includes:

- Setup checklist
- Software setup
- Starting and ending the software
- Screen size setup
- External program communication setup

Setup checklist

One NINDEX software license consists of:

- One CD including the software installer InsNINDEX-mm-x.y.z.exe, the software installer InsNindexMN.exe and NINDEX software user manual
- One Rockey4ND hardware dongle.

Software setup

System requirements

Recommended system requirements:

- Windows XP SP2, Windows XP SP3, Windows Vista or Windows 7 operating systems
- Display: 8 inches 15:9 and higher, resolution 1024 x 768 pixels and higher
- USB port to attach the hardware dongle
- USB port to attach the EEG acquisition device

- Free space on disk: 100 MB to install the program and 150 MB to store monitoring data
- Main memory: 1 GB and higher
- Processor: Intel Pentium with 1.6 GHz clock rate and higher.

Setup

To start the NINDEX program setup, as local administrator, double click InsNINDEX-mm-x.y.z.exe installer and follow the instructions.

The installer software does the following:

- Creates a c:\eegdata folder for NINDEX monitoring data and configuration data storage

Caution:

The software requires the c:\eegdata folder for proper operation. If this folder is accidentally deleted then the monitoring and configuration data will be lost.

- Creates a folder for NINDEX program installation
- Creates a “NINDEX” shortcut on the Desktop
- Creates a “NINDEX” group on Windows’ Start Menu All Programs list.

Starting and ending the software

To start NINDEX software follow these steps:

1. Connect the hardware dongle to an USB port
2. Double click “NINDEX” shortcut on the Desktop.

To end NINDEX software follow these steps:

1. Click the Menu button 
2. Click the “Close NINDEX” button.

Screen size setup

Upon executing NINDEX software for the first time, it shows up a “Screen Dimensions” dialog (Figure 1). Enter the screen dimensions (in mm) and then click the OK button.


Figure 1

External program communication setup

Note:

When used with an external program, the NINDEX software will not start the EEG processing until it receives the electrode impedance values.

Setup of data exchange with EEG system using Neuron-Spectrum.NET software

Working with Neuron-Spectrum devices Neuron-Spectrum.NET software can be used as anesthesia depth monitor. To organize the correct interaction between NINDEX and EEG system using Neuron-Spectrum.NET software it is necessary to setup the hardware port of Neuron-Spectrum.NET.

Connect Neuron-Spectrum-1 system to PC and run Neuron-Spectrum.NET software.

Choose “Settings” toolbar button in the program main window. In the drop down list that appears choose “Change” item. At that the setting window appears (see Figure 2).

Choose the “Anesthesia depth monitor” tab. Activate the “Allow data transfer for anesthesia depth calculation” by ticking the corresponding check box. Choose the required channel for data transmission (polychannel E1 by default). TCP data port is also specified on this tab (502 by default). At that the NINDEX software works as a client and external software as a server. Press OK.

Setting of connection with Neuron-Spectrum.NET external program

To set connection between NINDEX and Neuron-Spectrum.NET open “Settings” menu in Neuron-Spectrum.NET main window, choose “External programs” item from the drop down menu and click on “List...” item. In the window that appears click on “Add” button. In the “heading” field enter NINDEX and specify the filepath in the command line. By default the program is run by C:\Program Files\NINDEX\linhsbc.exe file.

Run NINDEX software using “External programs” menu. In the NINDEX main window that appears click on the  button. At that the program menu appears. Click on the “Connection” button. Enter the port address (502 by default). You can check the address in the “Anesthesia depth monitor” menu of Neuron-Spectrum.NET program. To start working with NINDEX press on the monitoring start button () In the program menu.


Figure 2

Setting window of Neuron-Spectrum.NET program

Herewith, if the current electrode montage contain the chosen channel, then data during the exam will be displayed in the Neuron-Spectrum.NET software and will be transmitted to the NINDEX software (it should be launched) (Figure 3). NINDEX program calculate the current anesthesia depth index (NINDEX value 100 indicates that patient is awake, 0 indicates patient's unconsciousness, and NINDEX values within 35 to 65 corresponds to insensibility to pain stimulus), displays it, and transmits to Neuron-Spectrum.NET software. Current NINDEX and burst suppression values are displayed in the Neuron-Spectrum.NET status line, and NINDEX trend can be observed in the "EEG trend" window. NINDEX trend allows to observe its dynamics during the exam. Besides, Neuron-Spectrum.NET software can inform you when current NINDEX value is out of specified normal values.


Figure 3

NINDEX software for anesthesia depth index calculation

Setup of data exchange with EEG system using Neuro-IOM.NET software

Neuro-IOM.NET software working with device of Neuro-MEP series can be used as anesthesia depth monitor.

To organize the correct interaction between NINDEX software and digital EEG system using Neuro-IOM.NET software it is necessary to setup the hardware port of Neuro-IOM.NET system.

Connect the Neuro-MEP-4 device to PC and run Neuro-IOM.NET software.

Choose the "Test templates" item in the "Settings/Change" menu. Start new exam. Press OK. Choose NINDEX template from the list. Change amplifier type and other settings if it is necessary. Open "Test/Edit test" drop down list. Choose the required amplifier (for example, Neuro-MEP-4) from the drop down list. To change the amplifier channel from which the signal will be acquired, use "Leads" tab. NINDEX software should be launched on the PC before monitoring start. In a few moments after monitoring start values of calculated parameters will appear in the anesthesia depth window.


Figure 4

The window of anesthesia depth monitoring if you use Neuro-IOM.NET software.

For correct calculation of parameters it is necessary to perform impedance check.

The indicator located in the top left corner of the screen shows the signal quality. If it is green – the signal quality is good, yellow – satisfactory, red – bad signal, it is impossible to calculate parameters.

To setup the window use “Setup..” local menu command. In the dialog box that appears (Figure 5) you can change some parameters.


Figure 5

The window of anesthesia depth setup using Neuro-IOM.NET software

“Enable” — anesthesia depth monitor activity flag.

“Lead” — choose the lead by which parameter calculation will be performed from the drop down list. The list of leads is determined by used test template.

“Calculation parameters” — choose parameters to be calculated. Values of chosen parameters can be displayed in trend window or connected with alarms.

Neuro-IOM.NET connection setup

The “Anesthesia depth monitor” (Figure 6) contains settings for work with NINDEX software.


Figure 6

Figure 6 – Window of Neuro-IOM.NET connection setup

“TCP port” — the port for data exchange with NINDEX software.

USING THE SOFTWARE

Caution:

Carefully read this section before using the software in a clinical setting.

This section provides instructions for using NINDEX software in a clinical environment. It includes the following sections:

- Electrode placement
- Electrode to EEG acquisition device connections
- Monitoring start
- Monitoring display
- The Menu Dialog
- Accelerator keys
- Monitoring end.

Electrode placement

NINDEX measurement uses one EEG channel recording acquired by 2 electrodes at locations 1 and 2 of Figure 7, and a third reference electrode at location 3 of Figure 7. Electrode locations 1, 2 and 3 are defined as follows:

- Location 1: intersection of the horizontal line across the middle of the forehead and the vertical line across the middle of the eyebrow
- Location 2: in front of the ear on the same side of the head as location 1, with a minimum distance of 10 cm to location 1
- Location 3: intersection of the horizontal line across the middle of the forehead and the vertical line across the nose.


Figure 7

Electrodes for NINDEX monitoring should be pre-gelled silver/silver chloride electrodes, e.g. Swaromed ECG electrodes.

Caution:

The use of accessories other than those specified may result in increased electrical interference on the EEG, decreasing the reliability of data provided by the software.

For electrode placement follow steps in Figure 8.

A	Rub the skin 10 times in a circular motion using a gauze and 70% isopropanol or alcohol	
B	Remove the electrode protective foil	
C	Wet the electrode with a drop of saline, with special care to not wet the self-adhesive part	
D	Attach the electrode to skin, ensuring complete adherence of self-adhesive part	

Figure 8

Caution:

Proper electrode placement is essential for the reliability of the data provided by the software.

Good skin to electrode contact is essential for the reliability of the data provided by the software.

Strictly follow the electrode placement procedure described in this manual.

Electrode to EEG acquisition device connections

Neuro-MEP-Micro (Neuro-EMG-MS)

Connect the electrode at location 1 to the “1+” input, the electrode at location 2 to the “1-” input and the electrode at location 3 to the connector for common wire (Figure 9).


Figure 9

In case you use Neuro-MEP-Micro (Neuro-EMG-MS) choose the CFM (aEEG) montage in Neuron-Spectrum.NET software.

Neuron-Spectrum-1

Connect the electrode at location 1 to the “+” input, the electrode at location 2 to the “-” input and the electrode at location 3 to the connector for common wire (Figure 10).


Figure 10

Monitoring start

Once the electrodes are correctly placed and connected to the EEG acquisition device, to start the monitoring click the [Monitoring Start / Stop] button .

Monitoring display

Figure 11 shows the monitoring display.


Figure 11

The monitoring display size should be approximately 16 cm x 10 cm. The monitoring display consists of four regions: the numeric region, the monitoring general data region, the graphics region and the buttons region.

The numeric region

The numeric region is in the upper left corner of the monitoring display. By default it consists of the rectangles "NINDEX", "SQ" and "BS", and optionally it adds the rectangle "Z". It is updated once every 3 s.

"NINDEX" rectangle in the upper left corner displays current NINDEX value in large, easily readable font. The number displayed disappears when SQ is red (low signal quality).

The "BS" rectangle displays current BSR (Burst Suppression Ratio) value. BSR is a calculated parameter to give the user an indication when an isoelectric EEG condition may exist. It is the percentage of time over the last minute period that the signal is considered to be in the suppressed state. For example, BSR = 20 means isoelectric EEG over 20% of the last minute review. The number displayed disappears when SQ is red (low signal quality).

The "SQ" (Signal Quality) rectangle displays current EEG signal quality value, which may be green, yellow or red. SQ relates to the percentage of discarded data (because of artifact detection) over the last minute of EEG processing. SQ is green if the percentage of discarded data is less than 25%; SQ is yellow if the percentage of discarded data is greater than 25% and less than 50%, or if the percentage of discarded data is greater than 50% and 15% or more of the data is considered to have ocular artifacts; SQ is red if the percentage of discarded data is greater than 50% and less than 15% of the data is considered to have ocular artifacts.

Figure 12

The “Z” rectangle may be added to the numeric region by selecting “Yes” in the “Show impedance” option in the “Settings” dialog (see “The Settings Dialog”). This rectangle shows impedance of electrode 1 to electrode 3 and impedance of electrode 2 to electrode 3 (see “Electrode placement”), in units of 4 kΩ. For example “0-1” means that electrode 1 to electrode 3 impedance is less than 4 kΩ y and that electrode 2 to electrode 3 impedance is greater than 4 kΩ and less than 8 kΩ.

Monitoring general data region

The monitoring general data region is located in the upper center of the monitoring display. It consists of the rectangles “Patient”, “Date”, “Time” and “Start time”.

“Patient” rectangle displays the two numbers that identify the monitoring. The first number is a 4-digit serial number unique to each NINDEX hardware dongle. The second is a 5-digit monitoring number which automatically increments upon each monitoring start.

“Date” rectangle displays current date in DD-MM-YYYY format.

“Time” rectangle displays current time in 24 hour format.

“Start Time” rectangle displays monitoring starting time in 24 hour format.

Graphics region

The graphics region consists of the EEG Waveform Display, in the upper right of the monitoring display, and the NINDEX Trend Display on the bottom left of the monitoring display.

The EEG Waveform Display shows a plot of raw EEG waveform over a 5-second time period. It is updated once every 1.5 seconds. Horizontal resolution is 1 second/div and vertical resolution can be set to 50 $\mu\text{V}/\text{div}$, 25 $\mu\text{V}/\text{div}$ or 12,5 $\mu\text{V}/\text{div}$ (see “The Settings Dialog”). Discarded data because of artifact detection is crossed by yellow lines (Figure 13).


Figure 13

NINDEX Trend Display shows a red plot of NINDEX value over a 90 minute, 180 minute or 270 minute time period. It is updated once every 15 s. Current NINDEX trend value is shown with a gray cursor, and is calculated as the median of NINDEX current value over last 15 second time period. The cursor disappears and the trend will not be displayed if SQ is red (low signal quality).

Due to relationship with clinical anesthesia and the emergence phase, NINDEX range from 45 to 99 is considered the most interesting for the user. NINDEX Trend Display vertical scale resolution is transformed by a nonlinear function for optimal visual resolution in this range.

The BSR trend plot may be added to the NINDEX Trend Display by selecting “Yes” in the “Show BS trend” option in the “Settings” dialog (see “The Settings Dialog”). BSR trend is plotted in blue. Current BSR trend value is calculated as the median of BSR current value over last 15 second time period. The trend will not be displayed if SQ is red (low signal quality).

Horizontal NINDEX Trend Display resolution can be set to 10 min/div, 20 min/div or 30 min/div (see “The Settings Dialog”).

If the monitoring time period exceeds NINDEX Trend Display time period, NINDEX trend plot is shifted one horizontal division towards left and the time origin is incremented.

Buttons region

The buttons region is located in the bottom right of the monitoring display. It consists of the [Menu] button  and the [Monitoring Start / Stop] button .

The Menu Dialog

To open the Menu Dialog click the [Menu] button .

The Menu Dialog consists of buttons with specific functions.

If monitoring is not activated, Menu Dialog functions are (Figure 14):


Figure 14

- Settings: shows the Settings Dialog
- Event Edition: shows the Event Edition Dialog
- Communications: shows the Communications Dialog
- Cancel: closes the Menu Dialog
- About ...: displays the software release, the NINDEX hardware dongle serial number, the software building (Neuro-EMG-MS, NS1 or Net-NSx) and general information about the software manufacturer Controles S.A.
- Close NINDEX: ends NINDEX software.

If monitoring is activated, Menu Dialog functions are (Figure 15):


Figure 15

- Settings: shows the Settings Dialog
- Event: shows the Event Dialog
- Cancel: closes the Menu Dialog
- About ...: displays the software version, the NINDEX hardware dongle serial number, the software building (Neuro-EMG-MS, NSx or Net-NSx) and general information about the software manufacturer Controles S.A.

- Close NINDEX: ends NINDEX software.
- Patient Data: shows the Patient Data Dialog
- Electrodes Check: enabled only if the EEG acquisition device is Neuro-EMG-MS or Neuron-Spectrum-1. Checks skin to electrodes contact (see "Checking of skin-to-electrode contact").

The Settings Dialog

Figure 16

Settings Dialog (Figure 16) consists of the following options:

- Trend Plot Period (min): NINDEX Trend Display time period. May be 90 min (10 min/div), 180 min (20 min/div) o 270 min (30 min/div)
- EEG Amplitude: EEG Waveform Display Amplitude. May be x1 (50 uV/div), x2 (25 uV/div) or x4 (12,5 uV/div)
- Show BS Trend: If "Yes" adds BSR trend to NINDEX Trend Display (see "Monitoring display")
- Show Impedance: If "Yes" adds the "Z" rectangle to the numeric region of the monitoring display (see "Monitoring display")
- High Alarm: if "Yes" enables the High Alarm (see "Alarms")
- Low Alarm: if "Yes" enables the Low Alarm (see "Alarms").

The Event Dialog


Figure 17

The Event Dialog (Figure 17) allows the user to add event markers in monitoring time. The user can create her/his own event types (see “The Event Edition Dialog”). The event markers are saved to disk and added to NINDEX Trend Display. Each event marker consists of a label identifying the event type and a time stamp relative to monitoring start time. The default event type list is:

- Induction
- Intubation
- Light Hypnosis (Isoflurane < 1%)
- Standard Hypnosis (1% < Isoflurane < 2%)
- Deep Hypnosis (Isoflurane > 2%)
- EEG Suppression
- Agent off
- Awakening.

To add an event marker, follow these steps in the Event Dialog:

1. Expand the Event Description drop down list by clicking the expand arrow 

2. Select the corresponding Event Description
3. Click the “OK” button.

The Event Edition Dialog


Figure 18

The Event Edition Dialog (Figure 18) displays the Event Type Edition List, which allows the user to replace the default event types by custom event types. The Event Type List may include up to 14 event types.

Each event type consists of:

- Identifier (Id column): 2-digit number that identifies the event type, from 01 to 14
- Label: character string that identifies the event type in NINDEX Trend Display and when the event marker is saved to disk. Up to 12 characters
- Description: character string that identifies the event type in the Event Dialog (see “The Event Dialog”). Up to 50 characters.

To edit an event type Label follow these steps in the Event Edition Dialog:

- 1) Click on the specific label string
- 2) Enter the new Label string in the edit box
- 3) On the keyboard press the “Enter” key.

Note:

Changes are not saved to disk until the “Save” button is clicked.

To edit an event type description follow these steps in the Event Edition Dialog:

- 1) Click on the specific Description string
- 2) Enter the new Description string in the edit box
- 3) On the keyboard press the “Enter” key.

Note:

Changes are not saved to disk until the “Save” button is clicked.

To add an event type, follow these steps in the Event Edition Dialog:

- 1) Click the “Add” button
- 2) Edit the event type Label and Description as described above.

Note:

Changes are not saved to disk until the “Save” button is clicked.

To delete an event type, follow these steps in the Event Edition Dialog:

- 1) Click on the specific event type “Id” column
- 2) Click the “Delete” button.

Note:

Changes are not saved to disk until the “Save” button is clicked.

The Patient Data Dialog

The screenshot shows a dialog box titled "Patient Data". It contains the following fields and controls:

- Name:** A text input field.
- Gender:** Two radio buttons labeled "Male" and "Female".
- Age (years):** A text input field with the value "0".
- Weight (kg):** A text input field with the value "0".
- Hospital:** A text input field.
- Comments:** A text input field.
- Buttons:** "Save" and "Cancel" buttons at the bottom right.

Figure 19

The Patient Data Dialog (Figure 19) allows editing patient data (name, gender, age and weight), the Hospital name and general comments. When the “Save” button is clicked the edited data is saved to disk.

Accelerator keys

NINDEX software allows adding specific parameters in monitoring time by pressing some accelerator keys. These functions are available:

- “D” key: allows adding current dose of Propofol and Remifentanyl. The data is added to the NINDEX Trend Display and saved to disk along with a time stamp
- “H” key: allows adding current value of hemodynamic variables. The data is added to the NINDEX Trend Display and saved to disk along with a time stamp
- “C” key: allows editing a comment. The edited text is saved to disk along with a time stamp
- “T” key: allows adding patient’s temperature. The data is added to the NINDEX Trend Display and saved to disk along with a time stamp
- “A” key: shows a help screen with accelerator key functions.

Monitoring end

To end the monitoring click the [Monitoring Start / Stop] button .

SOFTWARE FEATURES

Caution:

Carefully read this section before using the software in a clinical setting.

This section provides a description of NINDEX software main features. It includes the following chapters:

- Anesthesia depth monitoring
- Checking of skin-to-electrode contact
- Alarms
- Data storage to disk.

Anesthesia depth monitoring

The EEG acquisition device amplifies and the digitizes patient's EEG signal, then sends it directly to NINDEX software or to an external program which in turn sends it to NINDEX software.

NINDEX software filters the EEG digital data and then analyzes it for artifact. Artifact detection is necessary to correctly interpret the EEG. As NINDEX software is intended to provide only highly reliable data, it is therefore programmed not to display data whenever excessive artifact corruption occurs.

Data considered artifact-free is processed using digital signal processing techniques to compute a large number of EEG features, which are input into a multivariate statistical algorithm to do an automatic EEG classification according to a refinement of Kugler EEG scale. The result is expressed as the NINDEX number, which correlates to the patient's anesthetic depth such as NINDEX = 99 indicates awake state and NINDEX = 0 indicates flat line EEG.

Checking of skin-to-electrode contact

NINDEX software is intended for use with EEG acquisition devices that are able to check skin to electrode impedance.

If the EEG acquisition device is Neuro-EMG-MS or Neuron-Spectrum-1, then the impedance checking is done by the NINDEX software at the monitoring start, or when the “Electrodes check” button is clicked (see “The Menu Dialog”).

When used with an external program, the NINDEX software receives the impedance levels from the external program.

Note:

When used with an external program, the NINDEX software will not start the EEG processing until it receives the electrode impedance values.

If the impedance levels of the electrodes are not within acceptable limits, then the software displays a bad electrode contact warning.

Caution:

The data provided by the software is not considered reliable if the software shows up a bad electrode contact warning.

Alarms

NINDEX software generates two sound alarms: a high alarm and a low alarm. Any alarm may be silenced by the user (see “The Menu Dialog”).

A high alarm condition is triggered when NINDEX current value exceeds 75, and remains active if NINDEX current value is greater than 73.

A low alarm condition is triggered when NINDEX current value is below 36, and remains active if NINDEX current value is less than 38.

When an alarm is active, the software generates a 6-second alarm sound every one minute.

Data storage to disk

Monitoring data

Recorded data for the XXXX-YYYYY monitoring (see “Monitoring general data region”) is stored in the “c:\eegdata\patientXXXX-YYYYY” folder.

Features of monitoring data storage to disk:

- Data storage: NINDEX and BSR trend, electrode impedance, raw EEG waveform, event markers, patient data, date and time
- Storage time period: 15 s
- Storage capability: greater than 100 hours.

Software settings

Software settings are stored in the “c:\eegdata” folder.

Disk storage managing

The maximum amount of storage for monitoring data allowed by the software is enough for more than 100 monitoring hours.

The software verifies total monitoring data storage size upon startup. If the amount of storage exceeds the acceptable maximum, 25% of monitoring data is deleted using FIFO algorithm.

Caution:

In order not to lose monitoring data by means of automatic deletion, back up the data after maximum 100 monitoring hours.

SOFTWARE UPDATE

To update NINDEX software, follow these steps:

1. Contact the software manufacturer (Controles S.A.) to get the installer software for the update (see contact information on the foot of the first page of this manual).
2. Execute the installer software for the update and follow the instructions.

MONITORING DATA REVIEW

This section provides instructions to review NINDEX software's monitoring data using NINDEXMN software. It includes:

- Introduction
- Opening the monitoring data for review
- The monitoring data review display.

Introduction

NINDEXMN software for Windows 2000 or greater allows the user to review NINDEX software's monitoring data, including the review of trend data, event markers and patient data.

Opening the monitoring data for review

To review stored data for a specific monitoring follow these steps:

1. Click the "File Open"  button in the program Toolbar
2. In the opening dialog box browse to the monitoring data folder, then select the *.htr file and then click the Open button. In the example below (Figure 20), the monitoring data folder is "c:\eegdata\patient020-00040" and the *.htr file is "p00040m.htr".


Figure 20

The monitoring data review display

The monitoring data review display is shown in Figure 21. The NINDEX Trend Display is on the top and the BSR Trend Display is on the bottom. Horizontal resolution for NINDEX and BSR Trend Display is 10 min/div. All marked events are displayed in the NINDEX Trend Display. The NINDEXMN window title shows the monitoring number, patient’s name and age (if available), the monitoring date and the monitoring time. In the example below (Figure 21), the monitoring identifier is “paciente020-00040”, patient’s name is “Juan”, patient’s age is 13 years, the monitoring date is 17/03/2011 and the monitoring time is 09:25.


Figure 21

To review monitoring comment list click View on the program menu option and then click Comments.

To review monitoring patient's data click View on the program menu and then click Patient data.

To view the NINDEX Trend Display with constant (i.e. linear) vertical resolution follow these steps:

- 1) Click Options on the program menu and then click Display settings
- 2) Enable the check Linear NINDEX plot
- 3) Click the OK button.

In order to not display trend values if the impedance levels of the electrodes are not within acceptable limits (see below) follow these steps:

- 1) Click Options on the program menu and then click Display settings

- 2) Enable the check Do not display if bad electrodes
- 3) Click the OK button.

To set the acceptable limit of the electrode impedance levels for the review follow these steps:

- 1) Click Options on the program menu and then click Thresholds
- 2) Set the desired limit value on "Maximum electrode impedance"
- 3) Click the OK button.

Total broach, enumerate, fasten
and stamped 41 pages

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CONTROLES S.A.