Robert C. Meagher  
Senior Director  
World Wide Regulatory Affairs  
Carestream Health, Inc.  
150 Verona Street  
Rochester, New York 14608  
Telephone 585-627-6528

« 20 » January 2016

**Carestream**

**INSTRUCTION FOR USE**


Dental X-ray Unit Digital 3D tomographic CS 8100, versions CS 8100 3D, CS 8100 3D Access, with accessories

**Photo, diagram, drawing medical device**

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
## SAFETY INSTRUCTIONS

- ♣ You are responsible for the operation and maintenance of this unit. Only legally qualified persons can operate this unit. They **MUST** have training to use the radiological equipment. Do **NOT** open the cover of the unit. When necessary, have a trained authorized service technician carry out inspection and maintenance operations.
- ♣ Install this unit in an X-ray room that complies with current installation standards. From this location, you must be able to maintain visual or audio communication with the patient and be able to access the Acquisition interface module during exposure.
- ♣ This unit must be permanently connected to the ground with a fixed power supply cable. To avoid the risk of electric shock, this equipment must **ONLY** be connected to a mains supply with protective earth.
- ♣ Do **NOT** operate the unit if there is the threat of an earthquake. Following an earthquake, ensure that the unit is operating satisfactorily before using it again. Failure to observe this precaution may expose patients to hazards.
- ♣ X-ray equipment is hazardous to patients and the operator if you do not observe the exposure safety factors and operating instructions.
- ♣ Considering radiation safety of pediatric population, protocol for Acquisition on Pediatric patients must be followed. For more information on imaging pediatric patients more safely and effectively, refer to FDA Pediatric X-ray Imaging webpage:  
<http://www.fda.gov/radiationemittingproducts/radiationemittingproductsandprocedures/medicalimaging/ucm298899.htm>
- ♣ Do **NOT** place objects within the field of operation of the unit.
- ♣ The patient should wear a protective lead-lined shoulder apron with a thyroid collar, unless other Radiation Protection Protocols apply locally.
- ♣ While adjusting the height of the unit, ensure that the patient is kept clear of the mechanism.
- ♣ When the unit is not in use, ensure that the ON/OFF switch is set to OFF (O).
- ♣ If the unit develops a fault, switch it to off (O), display an “Unserviceable” notice and contact a service technician.
- ♣ To dispose of the unit or its components, contact a service technician.
- ♣ Ask the patient to refrain from moving during the entire period of exposure.

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- ♣ Ask the patient to remain still until the unit arm has stopped moving and the RESET movement has completed.
- ♣ Do NOT use this unit in conjunction with oxygen-rich environments. This unit is not intended for use with flammable anesthetics or flammable agents.
- ♣ It is not recommended to use accessories other than those specified in this document and sold by Carestream Health.

#### Computer

- ♣ Do NOT place the computer and the peripheral equipment connected to it in the immediate vicinity of the patient in the unit. Leave at least 1.83m distance between the patient and the unit. The computer and the peripheral equipment must conform to the IEC 60950 standard.
- ♣ See your computer installation guide for details of the data processing system and screen. Leave a sufficient amount of clear space around the CPU to ensure that it is properly ventilated.
- ♣ To obtain maximum image quality and visual comfort, position the screen to avoid direct light reflections from internal or external lighting.
- ♣ Always use Microsoft Windows Update to make sure that the latest security patches are correctly installed.

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### CS 8100 3D Overview

The CS 8100 3D Family comprises:

- CS 8100 3D: panoramic modality and dental volumetric reconstruction modality, obtaining 3D of the selected area and receiving 3D of the upper and lower jaws fully.
- CS 8100 3D Access: panoramic modality and dental volumetric reconstruction modality limited to 3D focused teeth acquisition. The system can be expanded field of view and receive 3D upper and lower jaws completely by purchasing additional licenses.

### Mobile Components

Figure illustrates the:

- Up and down movement of the unit
- Rotation and translation movement of the rotative arm

**Important: The patient can enter either the left side or the right side of the unit.**

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


### General Functional Components

Figure illustrates the general functional components of the unit.

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

- 1 ON/OFF button
- 2 Emergency stop knob
- 3 Unit rotative arm
- 4 X-ray source assembly
- 5 Unit column
- 6 Head and chin rest
- 7 Temple support
- 8 Digital sensor
- 9 Unit head
- 10 X-ray remote control
- 11 Computer hosting the imaging and the acquisition software

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## Head and Chin Rest

Figure illustrates the functional components of the head and chin rest.


- 1 Positioning Panel
- 2 Temple support adjustment knob
- 3 Temple support
- 4 Bite block
- 5 Chin rest
- 6 Hand grips

### Positioning Accessories

The following accessories are used when positioning a patient. They are delivered with the CS 8100 3D.

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
Table lists the positioning accessories.

| Accessory                                                                           | Description                          |
|-------------------------------------------------------------------------------------|--------------------------------------|
|    | Panoramic chin rest                  |
|    | Sinus chin rest                      |
|    | TMJ nose rest                        |
|    | Standard bite block                  |
|    | Bite block for edentulous patients   |
|                                                                                     | Single use sheaths for bite block    |
|  | Single use sheaths for 3D bite block |
|  | 3D bite block support                |

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
|                                                                                   |                     |
|-----------------------------------------------------------------------------------|---------------------|
|   | 3D bite block       |
|  | 3D molar bite block |
|  | 3D child bite block |

### Positioning Panel

The positioning panel is a console on the head and chin rest that enables you to correctly position and align a patient before you acquire an image.

#### 1 Height Adjustment buttons:

- Adjusts the height of the unit to the height of the patient.
- When the unit is not in use, places the rotative arm in a parallel position to the unit head, leaving more free space around the unit.

To do so, press and hold both buttons until the Ready Indicator LED switches off. When you release both buttons, the Ready Indicator LED will flash before the rotative arm turns to the parallel position.


2 Ready Indicator LED: Green indicates that the unit is ready for acquisition.

3 Temple support adjustment knob: Opens and closes the temple supports when you turn the knob.

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
## X-ray Remote Control Overview

The X-ray remote control enables you to launch a radiological image acquisition via the exposure button from outside the X-ray room. You must press and hold the exposure button until the end of acquisition. Premature release of the exposure button interrupts the acquisition.

1. Exposure button: launches image acquisition.

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
For the minimum computer system requirements and configuration for the Imaging Software, see the CS 8100 3D Safety, Regulatory and Technical Specifications User Guide.

**Important: It is MANDATORY to check that the computer system configuration is compatible with the computer system requirements for the CS 8100 3D software. If necessary you must UPDATE your computer system configuration. The unit MUST be connected to the computer via a point-to-point Ethernet link and not via a LAN.**

### General Software Overview

The CS 8100 3D operates with the following software:

- Imaging software

The Carestream's Dental Imaging Software or CS Imaging Software are user-friendly working interfaces that were designed and developed specifically to support image diagnosis.

- Acquisition interface

The Acquisition interface is a user-friendly interface that was designed and developed specifically for the CS 8100 3D.

The Acquisition interface has the following features:

- The Panoramic Acquisition interface.
- The 3D Acquisition interface.

### Acquisition Interface

#### Panoramic Acquisition Interface Overview

The Panoramic Acquisition interface is the main interface of the CS 8100 3D that provides you with imaging acquisition functions.

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
#### 1. Information button:

- Show versions: Identifies Software and Firmware versions.
- Additional features: Contains information on additional activated options and date of expiration.
- Resetting Anatomy setting: Resets to the manufacturing parameter settings.
- Memorize Anatomy setting: Memorizes the user preference settings for each patient type (kV, mA and seconds).

#### 2. Preview Screen: Displays the acquired acquisition frames.

#### 3. Selected Parameter Display: Displays the current acquisition parameter settings.


#### 4. System Status Screen: Displays various alert or warning messages originating from the unit.

#### 5. Generator Cooling indicator: Indicates the automatic cooling time (mm:ss) required before the next acquisition.

#### 6. Reset button: Resets the unit to the patient-entry position to enable you to position a patient in the unit.

#### 7. Parallel position button: Leaves more free space around the unit when it is not in use by placing the rotative arm in a parallel position to the unit head.

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8. X-ray On/Off button: Enables or disables X-ray emission.
9. Ready Indicator LED.
  - Green indicates the unit is ready to start acquisition.
  - Black indicates the unit is not ready to start acquisition.
10. Exit button: Closes the Acquisition interface.
11. X-ray Emission indicator: Yellow, indicates the X-ray emission status.
12. Selector Button: Selects different acquisition setting options.
  - Click Program to select examination type options.
  - Click Patient to select patient type parameters.
  - Click Parameters if you need to manually adjust the exposure parameter options.

### Program Pane

The Program pane enables you to choose the different radiological exams. The diagram in the Program pane represents a jaw, with the R side of the diagram corresponding to the right side of the patient.

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Radiological exam options  
Full panoramic acquisition:

Click on all the areas of interest. Example:

#### Segmented panoramic acquisition:

Click on the following areas of interest that you want to examine:

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The selected area of interest will highlight. Example (anterior exam):

#### Maxillary Sinus acquisition:

Click: 

#### TMJ acquisition:

Click  for a TMJ acquisition.

The selected are of interest will highlight:

Double click   appears for a TMJ x4 acquisition.

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
Select one or both area of interest: 

### Patient Pane

The Patient pane enables you to choose and select the different patient parameters. The selected parameters for a patient influences the image quality because:

- The default kV and mA exposure settings are set for the selected patient mode.
- The image is reconstructed according to the patient dental arch morphology. The selected parameters must be based on the age and morphology (jaw size and bone density) of the patient.

Note: A selected patient mode sets the:

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
Click  if the patient is large.

### 2 Patient dental arch morphology:

Click  for a **normal** shape dental arch.

Click  for a **square** dental arch.

Click  for a **sharp** dental arch.

### Parameter Pane

The Parameter Pane enables you to choose exposure parameters for the radiological image acquisition. If the default parameter setting is not adapted to your patient type, you can manually adapt the parameter setting to the patient type and save this setting as the default setting.

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3D full upper and lower jaw acquisition:

-  Full upper and lower jaw exam.
-  Upper only jaw exam.
-  Lower only jaw exam.
-  Fast scan: low dose exam.

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To examine a required block of the upper teeth area of interest.

To examine a required block of the lower teeth area of interest

To enable an acquisition at a higher resolution. Select if you want to have an image with more precision.

Fast scan: low dose exam.

### 3D Patient Pane

The **Patient** pane enables you to choose and select the different patient parameters. The selected parameters for a patient influences the image quality because the default kV and mA exposure settings are set for the selected patient mode. The selected parameters must be based on the age and morphology of the patient.

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

Note: A selected patient mode sets the kV and mA:

Patient type parameters.  
Pediatric mode.

Click  if the patient is a child.

Adult mode:

Click  if the patient is small.

Click  if the patient is medium.

Click  if the patient is large.

### 3D Parameter Pane

The **Parameter** Pane enables you to choose exposure parameters for the radiological image acquisition. If the default parameter setting is not adapted to your patient or program type, you can manually adapt the parameter setting to the patient or program type and save this setting as the default setting.

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
To save a parameter setting for a specific patient type, click  and select Memorize Anatomy setting.

1 Exposure setting options:

 kV (kilovolt)  
 mA (milliampere)

2 Fine-tuning buttons:


Click  or  to fine-tune the kV and the mA.

### Switching on the Unit

Before switching on the unit, check that:

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- The installation of the unit is complete.
- The computer is switched on.
- The imaging software is not opened.

**Important: You must switch on the computer and wait for it to be ready for the connection before switching on the unit.**

To switch on the unit, follow these steps:

1. On the unit column, press the ON button.
2. You must wait for a minute for the connection between the unit and the computer to be established. If you start the imaging software before the connection is established an error message is displayed. Click OK, close the imaging software and wait for the connection to be established.
3. You can now proceed to start the imaging software.

#### Increasing the operating life of the X-ray tube

**Important: You must carry out the following procedures at the first loading of the unit or if the unit has not been used for a month to increase the operating life of the X-ray tube:**

1. In the **Panoramic Acquisition Window**, select the **Parameter** pane.
2. Select the following parameter setting: 70 kV - 6.3 mA.
3. Leave the X-ray room and close the door.
4. Launch an X-ray for the parameter setting by pressing and holding the X-ray remote control exposure button.
5. Repeat **Steps 2 to 4** for the following parameters:
  - 80 kV - 10 mA
  - 85 kV - 10 mA

The unit is now ready to be used for acquisition.

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### Firewall or Network Settings

You will not be able to access the **Acquisition** interface if you did not configured the firewall or network settings.

To configured the firewall or network settings, follow these steps:

1. On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Imaging Software) to open the imaging software. The **Windows Security Alert** dialog box displays. Select **Private networks**, **Public networks**, and click **Allow access**.

2. Click **Quit** on the following dialog box:

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
Data Communication Blocked. Please check your firewall or network settings.

Err. 5, Pos. 37

OK

3. Close the imaging software.

4. On your desktop, double-click  or  to re-open the imaging software.

Note: See your firewall documentation for information on how to configure the firewall or network settings.

### Accessing the Cephalometric Acquisition Interface

To access the **Acquisition** interface, follow these steps:

1. On your desktop, double-click  or  to open the imaging software.

**Important:** If the Windows Security Alert dialog box opens to block the imaging software, you must configure the firewall or network settings.

2 Find or create the patient record.

3 Double-click the patient record to access the imaging window or tab.

4 In the imaging window or tab:

Data Communication Blocked. Please check your firewall or network settings.

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• Click  or  to access the **Panoramic Acquisition** interface.

• Click  or  to access the **3D Acquisition** interface.  
The acquisition interface is displayed.

### Acquiring a Panoramic, a TMJ x2, or a Sinus Image for the Adult and Pediatric Patient

Before acquiring an image, check that you have:

- Reset the unit rotative arm to the patient entry position for patient to enter the unit.
- Selected the patient record.
- Accessed the imaging window.
- Accessed the **Panoramic Acquisition** interface.

### Preparing the Unit and Setting the Acquisition Parameters

To set the acquisition parameters, follow these steps:

1. In the **Panoramic Acquisition** interface, click the **Patient** button to access the **Patient** pane.

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
Select the patient type:

- Child
- Adult: Small, Medium, Large Jaw size
- Dental arch morphology

2. Click the **Program** button to access the **Program** pane.

Click on the area of study that you want to examine (for example, a panoramic) .

3. **(Optional)** If the default parameter setting is not adapted to your patient type, click the **Parameter** button:

- Select the appropriate parameters.

- Click  and select **Memorize Anatomy setting** if you want to save the new parameter settings as future default settings.

4. Position and fix the panoramic chin rest (A). Fit the bite block to the panoramic chin rest (Ba). Make sure that you click it well into place. If needed, use the edentulous bite block or the sinus chin rest (Bb).

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


Note: Use the panoramic chin rest (A) for a panoramic or TMJ x2 image acquisition. Use the sinus chin rest (Bb) for a sinus image acquisition. Use the TMJ nose rest for a TMJ x4 acquisition.


5. Insert a hygienic barrier over the bite block (C). Make sure that you cover the bite block completely with the hygienic barrier (D).

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
## Preparing and Positioning the Adult and Pediatric Patient

To prepare and position the patient, follow these steps:

1. Ask the patient to remove all metal objects.

2. Ask the patient to wear a lead apron with a thyroid collar. Ensure that the apron lays flat across the patient's shoulders.


3. Turn the adjustment knob to open the temple supports laterally.

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

4. Ask the patient to enter the unit. On the Positioning panel, press  and hold to raise the chin rest to the height of the patient.

Note: If the patient is too tall, tell the patient to sit on a stool.

5. Ask the patient to:

- Stand up straight (**do not position the feet slightly forward**).
- Grip both the lower handles of the head and chin rest.
- Relax and lower the shoulders for full motion of the unit rotative arm

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
6. Ask the patient to place the chin on the chin rest support (A).

7. Ask the patient to bite into the groove of the bite block (B).

8. Turn the adjustment knob to close the temple supports (C) and use the horizontal positioning indicators as visual aid to adjust the up or down inclination of the patient head (D1) for a Frankfort plane horizontal alignment.


9. Use both hands to align the vertical position of the head (D).

10. Ask the patient to:

- Close the eyes (E).
- Swallow.
- Remain still.

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- Breathe through the nose.
- Place the tongue on the palate.

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1. Leave the X-ray room and close the door. You must keep visual contact with the patient during acquisition.

**Important:** To stop the acquisition, if there is any problem, release the exposure button of the remote control or press the red emergency stop button.

2. Launch the X-ray with the remote control. Press and hold the exposure button until the end of the acquisition is indicated by the on-screen “Release

Handswitch” message. The  turns yellow and a warning sound is heard, indicating X-ray emission. When the acquisition ends, the **Acquisition** interface disappears and the acquired image is automatically transferred to the imaging Window.

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3. Check the image quality.
4. Do the following when the acquisition is finished:
  - Release the patient.
  - Remove the hygiene barrier from the bite block.

### Acquiring a TMJ x4 Image for the Adult and Pediatric Patient

Before acquiring an image, check that you have:


- Reset the unit rotative arm to the patient entry position for patient to enter the unit;
- Selected the patient record;
- Accessed the imaging window;
- Accessed the **Panoramic Acquisition** interface.

### Preparing the Unit and Setting the Acquisition Parameters

To set the acquisition parameters, follow these steps:

- 1 In the **Panoramic Acquisition** interface, click the **Patient** button to access the **Patient** pane.

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Select the patient type:

- Child
- Adult: Small, Medium, Large Jaw size
- Dental arch morphology

2. Click the **Program** button to access the **Program** pane.

3. Click **x2**. The selected are of interest will highlight:

4. Double click **x2**. **x4** appears for a TMJ x4 acquisition.

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- Select the appropriate parameters.


- Click  and select **Memorize Anatomy setting** if you want to save the new parameter settings as future default settings.

6. On the head and chin rest, remove  or , place  and cover it with a hygienic barrier.

### Preparing and Positioning the Adult and Pediatric Patient

To prepare and position the patient, follow these steps:

1. Ask the patient to remove all metal objects.

2. Ask the patient to wear a lead apron with a thyroid collar. Ensure that the apron lays flat across the patient's shoulders.

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

3. Turn the adjustment knob to open the temple supports laterally.

Note: You can adjust the temple supports manually to the back or front for a better head fit.


4. Ask the patient to enter the unit. **On the Positioning panel**, press  and hold to raise the chin rest to the height of the patient.

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Note: If the patient is too tall, tell the patient to sit on a stool.

5. Ask the patient to:

- Stand up straight (**do not position the feet slightly forward**).
- Grip both the lower handles of the head and chin rest.
- Relax and lower the shoulders for full motion of the unit rotative arm

6. Turn the adjustment knob to close the temple supports.

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1. Ask the patient to remain still, eyes closed, mouth closed, and tongue placed on the palate. Leave the X-ray room and close the door. You must keep visual contact with the patient during acquisition.

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**Important:** To stop the acquisition, if there is any problem, release the exposure button of the remote control or press the red emergency stop button.

2. Launch the X-ray with the remote control. Press and hold the exposure button until the end of the acquisition is indicated by the on-screen “Release

Handswitch” message. The  turns yellow and a warning sound is heard, indicating X-ray emission. When the acquisition ends, the **Acquisition** interface disappears and the acquired image is automatically transferred to the **imaging Window**.

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3. Ask the patient to stay in the same position but with mouth opened. Repeat Step 2 to acquire another image with the mouth opened.

4. Check the image quality.

5. Do the following when the acquisition is finished:

- Release the patient.

- Remove the hygiene barrier and the  from the head and chin rest.

### X-ray Dose Emission Information

#### Compliance with EURATOM 97/43 Directive

You can right-click on each image to display the estimated emitted dose received by the patient. You can use this information to calculate the effective dose received by the patient for the image.

The radiation emission dose is expressed in  $\text{mGy}/\text{cm}^2$ . This dose is measured at the primary collimator outlet. The dose is accurate to  $\pm 30\%$ .

### Acquiring a 3D Image of the Full Upper and Lower Jaw for the Adult and Pediatric Patient

Before acquiring an image, check that you have:

- Reset the unit rotative arm to the patient entry position for patient to enter the unit;

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- Selected the patient record;
- Accessed the imaging window;
- Accessed the **3D Acquisition** interface.

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1 In the **3D Acquisition** interface, click the **Patient** button to access the **Patient** pane.

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3. Click on one of the following programs and select the area of interest that you intend to examine:

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5. Position and fix the 3D bite block support (A). Fit either one of the three 3D bite blocks to the 3D bite block support (B). Make sure that you click it well into place.

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


See the table below for information on the bite block that you must use for the corresponding area of interest and the acquisition type.

| Full upper and lower jaw acquisition*                                               |                                                                                                         |                                                                                                               |                                                                                                               |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Acquisition option                                                                  | 1<br>3D bite block<br> | 2<br>3D molar bite block<br> | 4<br>3D child bite block<br> |
|    | 8x9<br> **             | 8x9<br> **                   | 8x8<br>                     |
|   | 8x5<br>               | 8x5<br>                     | 8x5<br>                    |
|  | 8x5<br>              | 8x5<br>                    | 8x5<br>                   |











\*NOT available for CS 8100 3D Access.

\*\*In Canada, the full upper and lower jaw acquisition is limited to 8x8 FOV instead of 8x9 FOV.

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
6. Remove the paper backing from the hygienic barrier (C). Cover the bite block with the hygienic barrier (D). Make sure that you cover the bite block completely with the hygienic barrier (E).

**Important: You must cover the bite block with FDA-cleared protective sheaths that are available from distributors.**

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


4. Ask the patient to enter the unit. On the **Positioning panel**, press  and hold to raise the chin rest to the height of the patient.

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
6. Ask the patient to bite into the 3D bite block (A).

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
7. Make sure that the patient places the upper incisors before the incisive stopper (B).

**Important: The spinal column and the nose of the patient must be aligned in a straight vertical line.**

8. Turn the adjustment knob to close the temple supports (C).

9. Use both hands to align the vertical position of the head (D). Ensure that the patient occlusion is fully horizontal (camper plane) (D1).

10 Ask the patient to:

- Close the eyes (E).
- Swallow.
- Remain still.
- Breathe through the nose.
- Place the tongue on the palate.

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**Important: To stop the acquisition, if there is any problem, release the exposure button of the remote control or press the red emergency stop button.**

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- Select the appropriate parameters.
- Click  and select **Memorize Anatomy setting** if you want to save the new parameter settings as future default settings.



5. Position and fix the bite block support (A). Fit the appropriate bite block to the bite block support (B). Make sure that you click it well into place.

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See the table below for information on the bite block that you must use for the corresponding are of interest and the acquisition type.

| 3D focused teeth acquisition                                                       |                                                                                            |                                                                                            |                                                                                            |
|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
| Acquisition option                                                                 | 3D bite block (1)                                                                          | 3D molar bite block (2)                                                                    | 3D child bite block (4)                                                                    |
|   | <br>5 x 5 | <br>5 x 5 | <br>4 x 4 |
|  | <br>5 x 5 | <br>5 x 5 | <br>4 x 4 |

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Handswitch” message. The  turns yellow and a warning sound is heard, indicating X-ray emission. When the acquisition ends, the **Acquisition** interface disappears and the volume is set under reconstruction. During reconstruction, release the patient and remove the hygiene barrier from the bite block.

3. Open and review the volume from you imaging software.

## X-ray Dose Emission Information

### Compliance with EURATOM 97/43 Directive

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
## 3D Object Acquisition

The CS 3D Object Acquisition comprises the following acquisition modes:

- **CS Restore:** Acquires 3D representations of digital X-ray images that are used for dental restorative workflows.
- **CS Model:** Acquires 3D representations of digital X-ray images that are used for dental orthodontic workflows.
- **Appliance:** Acquires 3D representations of guided surgical images for double-scan protocol producing acquisition files that are compatible with competitor solutions.

## General Acquisition Interface Overview

1. **Show versions button:** Identifies Software and Firmware versions.
2. **Animated display screen:** Guides you on the acquisition positioning

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4. **System Status Screen:** Displays various alert or warning messages originating from the unit.

5. **Generator Cooling indicator:** Indicates the automatic cooling time (mm:ss) required for the generator to reach 0 for a new acquisition.

6. **Reset button:** Resets the unit to the patient-entry position.

7. **Parallel position button:** Leaves more free space around the unit when it is not in use by placing the rotative arm in a parallel position to the unit head.

8. **X-ray On/Off button:** Enables or disables X-ray emission.

9. **Ready Indicator LED:**

- Green indicates the unit is ready to start acquisition.
- Black indicates the unit is not ready to start acquisition.

10. **Exit button:** Closes the Acquisition interface.

11. **X-ray Emission indicator:** Yellow indicates that X-ray emission is in progress.

12. **Acquisition Options:** Provides the various acquisition scenarios.

13. **Acquisition Export Mode:** Provides one of the following acquisition export modes to view the acquired images:

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
- Click
- Select

Click  for the triple tray impression mode.

Click  for the lower jaw mode.

Click  for the upper jaw mode.

### CS Model: Plaster Pane Overview

The CS Model plaster program pane becomes active when you:

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- Click
- Select

**Important:** To have a true occlusion alignment of the full jaw mode, make sure that the rear ends of both the occlusion models are cut flat. If they are not, you must select the full jaw mode with wax bite.

Click  for the full jaw mode.

Click  for the full jaw mode with wax bite.

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## CS Model: Impression Material Pane Overview

The CS Model impression material program pane becomes active when you:

- Click
- Select

Click  for the full jaw mode.

Click  for the full jaw mode with wax bite.

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Appliance: Radiological Guide Pane Overview

The  radiological guide program pane becomes active when you .

Click  for the lower jaw mode.

Click  for the upper jaw mode.

3D Object Acquisition Accessories

Table: 3D Object Acquisition Accessories

| Item                                                                                | Description                         | Quantity |
|-------------------------------------------------------------------------------------|-------------------------------------|----------|
|  | 3D dental impression holder support | 1        |

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
2. Place the dental impression holder (B) onto the 3D bite block support (A).
3. Place the black foam (C) on the dental impression holder (B).
4. Click the **Program** button to access the **Program** pane.
5. Click  to access the impression material acquisition mode interface.
6. Select the acquisition option that you require. See “Launching the Acquisition for the Impression Material” or “Launching the Triple Tray Impression Acquisition”.

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
#### Launching the Acquisition for the Impression Material

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
- **ALWAYS** use a metal-free tray for the Impression Material.
- **Make sure that Impression Material and the plastic tray are dry.**

To launch the acquisition, follow these steps:

1. Select either the lower jaw or upper jaw acquisition mode.

2. Follow the on-screen animated instruction on how to correctly position the impression at the center of the black foam.

3. Leave the X-ray room and close the door.
4. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.

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
The indicator in the acquisition turns yellow and a warning sound is emitted to indicate that X-ray emission is in progress.

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5. Once all items have been scanned, the **Acquisition** interface closes. Wait for the 3D object reconstruction.

The completely reconstructed object will display in the image browser.

### Launching the Triple Tray Impression Acquisition

To launch the triple tray impression acquisition, follow these steps:

**Important:**

- Make sure that the triple tray impression is dry.
- The triple tray impression must be level when you position it on the black foam. If it is not level, use a scalpel to cut and level it.

1. Select the triple tray impression acquisition mode.

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
## Preparing the Unit and Setting the Acquisition Options

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
2. Place the dental impression holder (B) onto the 3D bite block support (A).
3. Place the black foam (C) on the dental impression holder (B).
4. Click the **Program** button to access the **Program** pane.

5. Click  to access the plaster acquisition mode interface.
6. Select the acquisition option that you require. See “Launching the Acquisition for the Impression Material” or “Launching the Triple Tray Impression Acquisition”.

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


### Launching the Plaster Acquisition

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2. Follow the on-screen animated instruction on how to correctly position the upper jaw or lower jaw plaster at the center of the black foam.

3. Leave the X-ray room and close the door.
4. On the **Acquisition** interface, when  is green, press and hold the

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exposure button of the X-ray remote control until the end of acquisition.

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The reconstructed object will display in the image browser.

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
2. Follow the on-screen animated instruction on how to correctly position the triple tray impression on the black foam.

A

B

**Important:** You must place the triple tray impression on the black foam in the same orientation as when you take it out of the patient mouth.

3. Leave the X-ray room and close the door.

4. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.


**Important:** If you have a problem that requires you to stop the

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
acquisition, release the exposure button of the remote control or press the red emergency stop button.

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### Acquiring Images of CS Model Impression Material Objects

Before acquiring an image of the object, do the following:

- On your desktop, double-click  (Carestream's Dental Imaging Software)

or  (CS Imaging Software) to open the imaging software.

- Select the patient record.
- Access the imaging window.

• If you use Carestream's Dental Imaging Software, click  to access the **Acquisition** interface.

• If you use CS Imaging Software, click  then  to access the **Acquisition** interface.

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
To prepare the unit for acquisition, follow these steps:  
1. Insert the 3D bite block support (A) in the chin rest base.

2. Place the dental impression holder (B) onto the 3D bite block support (A).
3. Place the black foam (C) on the dental impression holder (B).
4. Click the **Program** button to access the **Program** pane.
5. Click  to access the plaster acquisition mode interface.
6. Select the acquisition option that you require. See “Launching the Acquisition for the Impression Material”.

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
**Launching the Acquisition for the Impression Material**

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### Important:

- **ALWAYS** use a metal-free tray for the Impression Material.
- Make sure that Impression Material and the plastic tray are dry.

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1. Select either one of the acquisition options..

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
The acquired image appears in the **Preview Screen**. If you select the option that requires the acquisition of the wax bite, see the following section for information: “Launching the Wax Bite Acquisition”.

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### Launching the Wax Bite Acquisition

To launch the wax bite acquisition, follow these steps:

**Important:** Make sure that the wax bite is dry.

2. Follow the on-screen animated instruction on how to correctly position the wax bite at the center of the black foam.

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
2. Leave the X-ray room and close the door.
3. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.

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4. Once all items have been scanned, the **Acquisition** interface closes. Wait for the 3D object reconstruction.

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**Acquiring Images of CS Model Plaster Objects**

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- On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Imaging Software) to open the imaging software.
- Select the patient record.
- Access the imaging window.

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
2. Place the dental impression holder (B) onto the 3D bite block support (A).
3. Place the black foam (C) on the dental impression holder (B).

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
## Launching the Wax Bite Acquisition

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**Important: Make sure that the wax bite is dry.**

1. Follow the on-screen animated instruction on how to correctly position the wax bite at the center of the black foam.

2. Leave the X-ray room and close the door.

3. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.

**Important: If you have a problem that requires you to stop the acquisition, release the exposure button of the remote control or press the**

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
### Acquiring Images of the Radiological Guide in the Patient Mouth

The Appliance acquisition export mode enables you to carry out the acquisition procedures for dual scan protocol that is needed to perform dental implant guided surgery.

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- On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Imaging Software) to open the imaging software.
- Select the patient record.
- Access the imaging window.
- Click  (Carestream's Dental Imaging Software) or  (CS Imaging Software) to access the Acquisition interface.

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

2. Ask the patient to wear a lead apron. Ensure that the apron lays flat across the patient's shoulders.

3. Ask the patient to do the following:

- Stand up straight or sit on a stool.
- Grip the lower handle on each side.
- Position the feet slightly forward.
- Relax and lower the shoulders to allow for the full motion of the unit rotative arm.

4 Place the radiological guide (A) and the bite index (B) in the patient's mouth. Ask the patient to remain in occlusion during the scanning operation.

5 Ask the patient to bite into the 3D bite block (A).

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

6. Make sure that the patient places the upper incisors before the incisive stopper (B).

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7. Turn the adjustment knob to close the temple supports (C).

8 Use both hands to align the vertical position of the head (D). Ensure that the patient occlusion is fully horizontal (camper plane) (D1).

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
2. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.

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**Important: If you have a problem that requires you to stop the acquisition, release the exposure button of the remote control or press the red emergency stop button.**

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3. The animated display in the **Acquisition** interface changes to demonstrate how to correctly position the next item to be scanned. Repeat steps 1 to 2.

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### Acquiring Images of the Radiological Guide Alone

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- or  (CS Imaging Software) to open the imaging software.

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- If you use Carestream's Dental Imaging Software, click  to access the **Acquisition** interface.

- If you use CS Imaging Software, click  then  to access the

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
## Acquisition interface.

### Preparing the Unit and Setting the Acquisition Options

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2. Place the dental impression holder (B) onto the 3D bite block support (A).
3. Place the black foam (C) on the dental impression holder (B).
4. Place the Radiological Guide on the black foam.

**Important:** You must place the Radiological Guide on the black foam in the same orientation as when you take it out of the patient mouth.

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
5. Click the **Program** button to access the **Program** pane.

6. Click  to access the Radiological Guide acquisition mode interface.

 appears in the Acquisition Export Mode field.

7. Select either one of the appliance acquisition options.

 **Important:** The animated display screen will guide you on how to position the acquisition object that you have selected.

 **WARNING:** Make sure that you position the acquisition object correctly to obtain quality data.

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3. Once all items have been scanned, the **Acquisition** interface closes. Wait for the 3D object reconstruction.

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## Acquiring Appliance-Based Images

**Important:** Make sure that you have a calibration phantom that is provided by Nobel Biocare before you perform a NobelGuide™ or NobelClinician™ calibration procedure.

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
Before acquiring images of the NobelGuide™ or NobelClinician™ calibration phantom, do the following:

- On your desktop, double-click  (Carestream's Dental Imaging Software) or  (CS Imaging Software) to open the imaging software.
- Select the patient record.
- Access the imaging window.
- If you use Carestream's Dental Imaging Software, click  to access the **Acquisition** interface.
- If you use CS Imaging Software, click  then  to access the **Acquisition** interface.

### Preparing the Unit and Setting the Acquisition Options of the NobelGuide™ or NobelClinician™ Calibration Phantom

To prepare the unit for acquisition, follow these steps:

1. Insert the 3D bite block support (A) in the chin rest base.

2. Place the dental impression holder (B) onto the 3D bite block support (A).

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

3. Place the black foam (C) on the dental impression holder (B).
4. Place the NobelGuide™ or NobelClinician™ calibration phantom on the black foam.


5. Click the **Program** button to access the **Program** pane.
6. Click  to access the Radiological Guide acquisition mode interface.
-  appears in the Acquisition Export Mode field.
7. Select either one of the appliance acquisition options.

**Important:** The animated display screen will guide you on how to position the acquisition object that you have selected.

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**WARNING:** Make sure that you position the acquisition object correctly to obtain quality data.

### Launching the X-Ray To Acquire an Image of the NobelGuide™ or NobelClinician™ Calibration Phantom

To launch the X-ray, follow these steps:

1. Leave the X-ray room and close the door.
2. On the **Acquisition** interface, when  is green, press and hold the exposure button of the X-ray remote control until the end of acquisition.

**Important:** If you have a problem that requires you to stop the acquisition, release the exposure button of the remote control or press the red emergency stop button.

The  indicator in the acquisition turns yellow and a warning sound is emitted to indicate that X-ray emission is in progress.

The acquired image appears in the **Preview Screen**.

3. The animated display in the **Acquisition** interface changes to demonstrate how to correctly position the next item to be scanned. Repeat steps 1 to 2.
4. Once all items have been scanned, wait for the 3D object reconstruction. The reconstructed object will display in the image browser.

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## LEGEND OF LABEL

| Symbols                                                                             | Means                                                                                                                                                                                                                                                                                       |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Type BF device symbol complying with the IEC 60601-1 standard.                                                                                                                                                                                                                              |
|    | In the European Union, this symbol indicates: DO NOT discard this product in a trash receptacle; use an appropriate recovery and recycling facility. Contact your local sales representative for additional information on the collection and recovery programs available for this product. |
|    | Warning                                                                                                                                                                                                                                                                                     |
|    | Refer to instruction manual/booklet                                                                                                                                                                                                                                                         |
|   | Manufactured date.                                                                                                                                                                                                                                                                          |
|  | Manufacturer's address                                                                                                                                                                                                                                                                      |
|  | Note                                                                                                                                                                                                                                                                                        |
|  | European conformity mark                                                                                                                                                                                                                                                                    |
| SN                                                                                  | Serial number                                                                                                                                                                                                                                                                               |

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**THE MATERIALS FROM WHICH IT IS MADE OR ITS MEDICAL DEVICE BASIC PARTS**

| Product Name                         | Material                                         |
|--------------------------------------|--------------------------------------------------|
| Temple support                       | Polycarbonate (Makrolon RX 1805)                 |
| Bite block                           | 104-701 polycarbonate polymer + carbon fiber 20% |
| Standard bite block                  |                                                  |
| Bite block for edentulous patients   |                                                  |
| 3D bite block                        | Polycarbonate (Makrolon RX 1805)                 |
| 3D molar bite block                  |                                                  |
| 3D child bite block                  |                                                  |
| Single use sheaths for bite block    | Ethylene methyl acrylate copolymer               |
| Single use sheaths for 3D bite block |                                                  |
| Panoramic chin rest                  | Polycarbonate (Makrolon RX 1805)                 |
| 3D bite block support                |                                                  |
| Sinus chin rest                      | Polycarbonate (Makrolon DR 1452)                 |

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## STERILIZATION

### Hygiene and Disinfection

#### Cleaning and disinfecting the Accessories

Cleaning and disinfecting the accessories that have contact with the mucous membranes

**CAUTION:** You **MUST** cover the standard bite block and the bite block for edentulous patients with FDA-cleared or CE mark protective sheaths that are available from distributors to use them between each patient. We recommend that you cover the TMJ nose rest and the 3D bite block with FDA-cleared or CE Mark protective sheaths that are available from distributors to use them between each patient.

The following accessories must first be cleaned and then steam-sterilized between each patient use:

- TMJ nose rest;
- Standard bite block;
- Bite block for edentulous patient;
- 3D bite block.

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Note: It is recommended that the accessory be reprocessed as soon as reasonably practical following use.

### Cleaning

To clean the accessories that have contact with the mucous membranes, follow these steps:

1. Remove and discard the protective sheath from the accessory.
2. Remove all visible soil by with disposable cloth or paper wipe.
3. Rinse at least 1 minute under running water to thoroughly clean the accessory from any excess soil.
4. Using a soft brush, apply medical enzyme detergent solutions (basically with a multi-enzymatic formula) to all surfaces of the accessory. Detergent manufacturer's directions must be strictly adhered to.
5. Rinse thoroughly under running water for at least 1 minute to remove detergent residue.
6. Dry the accessory with compressed air or hygiene disposable cloth.
7. Visually inspect the accessory for residual soil. If soil is visible, either repeat steps 2 to 5, or safely dispose of the accessory.

### Sterilizing with Steam Autoclave

To steam autoclave the accessory, once cleaning is complete, follow these steps:

**CAUTION:** You must use a medical autoclaving equipment cleared by the FDA in the USA or that is recognized by your Local Authority. You must always follow the operating parameters recommended by the manufacturer of the autoclaving equipment. Use

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### FDA cleared or CE mark standard packaging material.

1. Wrap the cleaned accessory using a standard packaging material for autoclaving.
2. Steam autoclave at 132°C (270°F) for 4 minutes in the USA or depending on your local regulation you can steam autoclave at 134°C (273°F) for 18 minutes.
3. Visually inspect the accessory for signs of deterioration. If any damage is noted, do not use the accessory and contact your representative.
4. Once sterilized, the accessory can be used immediately or stored dry and dust-free in its sterilization wrapping under temperature specified in section "CS 8100 3D Environmental Requirements" of the present guide.

### Cleaning and disinfecting the accessories that have skin contact

The following accessories must first be cleaned and then disinfected between each patient use:

- Panoramic chin rest;
- Sinus chin rest;
- Occlusion chin rest.

Note: It is recommended that the accessory be reprocessed as soon as reasonably practical following use.

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To manually clean the accessories that have skin contact, follow these steps:

1. Remove all visible soil by with disposable cloth or paper wipe.
2. Rinse at least 1 minute under running water to thoroughly clean the accessory from any excess soil.
3. Using a soft brush, apply medical enzyme detergent solutions (basically with a multi-enzymatic formula) to all surfaces of the accessory. Detergent manufacturer's directions must be strictly adhered to.
4. Rinse thoroughly under running water for at least 1 minute to remove detergent residue.
5. Dry the accessory with compressed air or hygiene disposable cloth.
6. Visually inspect the accessory for residual soil. If soil is visible, either repeat steps 1 to 4, or safely dispose of the accessory.

### Disinfecting

To disinfect the accessory, once the cleaning is complete, follow these steps:

1. Disinfect the accessory by using an EPA-registered hospital disinfectant for low-level activity or low-level disinfectant that is recognized by your Local Authority (for example, quaternary ammonium compounds and some phenolics). You must follow the disinfectant's manufacturer instructions for use, especially with respect to contact time.

**CAUTION:** If there is visible contamination of blood components should be cleaned with a hospital disinfectant that is registered with the US Environmental Protection Agency (EPA) as a mid-level disinfectant or disinfectant.

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