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Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the operational documentation. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: [Signature]
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
Erg.Nr. 2417/19
Geb. nach GebO

bezahl

Read Highlighted Changes: Revised July 2018.

Package insert Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT Syphilis TP

INTENDED USE

The ARCHITECT Syphilis TP assay is a chemiluminescent microparticle immunoassay (CMIA) for the qualitative detection of antibodies to *Treponema pallidum* (TP) in human serum and plasma, including specimens collected post-mortem (non-heart-beating). The ARCHITECT Syphilis TP assay is intended to be used as an aid in the diagnosis of Syphilis infection and as a screening test to prevent transmission of *Treponema pallidum* to recipients of blood, blood components, cells, tissue and organs.

SUMMARY AND EXPLANATION OF THE TEST

Syphilis is caused by infection with the bacterium TP¹ which can be transmitted congenitally or by sexual contact. The disease can evolve into a latent phase in which syphilis is clinically inapparent. Serological tests (nontreponemal and treponemal specific), in addition to patients' clinical history, are currently the primary methods for the diagnosis and management of syphilis.

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The ARCHITECT Syphilis TP assay is a two-step immunoassay for the qualitative detection of antibody to TP in human serum or plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

- Sample, microparticles coated with recombinant TP antigens (TpN15, TpN17 and TpN47) and Assay Diluent are combined. Anti-TP antibodies present in the sample bind to the TP coated microparticles.
- After washing, acridinium-labeled anti-human IgG and IgM conjugate is added to create a reaction mixture.
- Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
- The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-TP antibodies in the sample and the RLUs detected by the ARCHITECT iSystem optics.

The presence or absence of anti-TP antibodies in the specimen is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from a previous ARCHITECT Syphilis TP calibration. If the chemiluminescent signal in the specimen is greater than or equal to the cutoff signal, the specimen is considered reactive for anti-TP.

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Kit Contents

ARCHITECT Syphilis TP 8D06

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	8D06-32	8D06-42
	100	500
MICROPARTICLES	1 x 4.1 mL	1 x 16.3 mL
CONJUGATE	1 x 5.9 mL	1 x 26.3 mL
ASSAY DILUENT	1 x 6.3 mL	1 x 34.1 mL

MICROPARTICLES TP (*E.coli*, recombinant) antigen coated microparticles in HEPES buffer with detergent. Minimum concentration: 0.08% solids. Preservatives: sodium azide and other antimicrobial agents.

CONJUGATE Murine anti-IgG/anti-IgM acridinium-labeled conjugate in MES buffer with protein (bovine) stabilizer. Minimum concentration: (anti-IgG) 26.6 ng/mL / (anti-IgM) 1.34 ng/mL. Preservatives: sodium azide and other antimicrobial agents.

ASSAY DILUENT Syphilis TP Assay Diluent containing MES buffer with detergent. Preservatives: ProCln 950 and other antimicrobial agent.

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

Warnings and Precautions

- IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.²⁻⁵

The following warnings and precautions apply to: MICROPARTICLES / CONJUGATE	
	
WARNING:	Contains Polyethylene glycol octylphenyl ether (Triton X-405) and Sodium azide.
H318	Causes serious eye irritation.
EUH032	Contact with acids liberates very toxic gas.

Prevention	
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: **ARCHITECT DILUENT**

WARNING:	Contains Polyethylene glycol octylphenyl ether (Triton X-100) and Methylisothiazolone.
H319	Causes serious eye irritation.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between kits.
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.

- Septums **MUST** be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.
- When handling conjugate vials, change gloves that have contacted human serum or plasma, since introduction of human IgG/IgM will result in a neutralized conjugate.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage. Store in upright position.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

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When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT Syphilis TP assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

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Specimen Types

Verified specimen types to be used with this assay:

Specimen Types	Collection Tubes
Human serum	Serum
	Serum separator tubes
	Potassium EDTA
	Lithium heparin
Human plasma	Sodium heparin
	Sodium citrate
	CPD

- Other specimen collection tube types have not been tested with this assay.
- Performance has not been established for the use of body fluids other than human serum or plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating), for details refer to section **Testing of Cadaveric Blood Specimens**.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - grossly hemolyzed (> 500 mg/dL hemoglobin)
 - obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at ≥ 10,000 RCF (Relative Centrifugal Force) for 10 minutes before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.

Transfer clarified specimen to a sample cup or secondary tube for testing.

- Centrifuged specimens with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified specimen without the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum	Room temperature	≤ 72 hours
	2-8°C	≤ 7 days
Plasma	Room temperature	≤ 72 hours
	2-8°C	≤ 30 days

If testing will be delayed, serum or plasma should be removed from the clot, red blood cells, or separator gel.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

Testing of Cadaveric Blood Specimens

- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating) that have been collected up to 21.5 hours after death. Performance was established using 50 spiked and 50 non-spiked cadaveric blood specimens.⁷
- Testing of cadaveric blood specimens from patients with plasma dilution due to transfusions of > 2000 mL of blood or colloids within 48 hours, or > 2000 mL of crystalloids within 1 hour (or any combination thereof) prior to collection of the specimens have not been validated.
- Follow general standards and/or regulations for collection, storage and handling.
- Follow the tube manufacturer's processing instructions for serum or plasma collection tubes. After initial centrifugation, transfer the supernatant to a centrifuge tube and centrifuge at 10,000 RCF (Relative Centrifugal Force) for 10 minutes. If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.
- Cadaveric blood specimens can be stored for up to 7 days at 2-8°C or up to 1 day at 15-30°C following collection.
- No qualitative differences were observed for cadaveric blood specimens (nonreactive or spiked reactive) when subjected to up to 3 freeze/thaw cycles. However, multiple freeze/thaw cycles should be avoided.

PROCEDURE

Materials Provided

8D06 ARCHITECT Syphilis TP Reagent Kit

Materials Required but not Provided

- ARCHITECT Syphilis TP Assay file obtained from the ARCHITECT iSystem s-Assay CD-ROM found on www.abbottdiagnostics.com.
- 8D06-04 ARCHITECT Syphilis TP Calibrator
- 8D08-13 ARCHITECT Syphilis TP Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

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Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - **Invert the microparticle bottle 30 times.**
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
 - Once the microparticles have been resuspended, discard the cap and place a septum on the bottle. For instructions on placing septums on bottles refer to the Reagent Handling section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

 - Priority:
 - Sample volume for first test: 80 µL
 - Sample volume for each additional test from same sample cup: 30 µL
 - ≤ 3 hours on board:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 30 µL
 - > 3 hours on board: Additional sample volume required. For information on sample evaporation and volumes, refer to the ARCHITECT System Operations Manual, Section 5.
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Syphilis TP Calibrator and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles vertically and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each calibrator: 5 drops
 - for each control: 5 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures

Specimen Dilution Procedures

Specimens cannot be diluted for the ARCHITECT Syphilis TP assay.

Calibration

- **Test Calibrator 1** in replicates of three. The calibrator should be priority loaded.

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Once an ARCHITECT Syphilis TP calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT Syphilis TP assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures. Ensure that assay control values are within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

■ RESULTS

The ARCHITECT iSystem calculates the cutoff (CO) using the mean chemiluminescent signal (RLU) from three replicates of the Calibrator 1 and stores the result.

Calculation

The ARCHITECT Syphilis TP assay calculates a result based on a cutoff determined by the following calculation.

- $\text{Cutoff (CO)} = \text{Calibrator 1 Mean RLU} \times 0.20$
- $\text{S/CO} = \text{Sample RLU} / \text{Cutoff RLU}$
- The cutoff RLU is stored for each reagent lot calibration.

Interpretation of Results

- Specimens with S/CO values < 1.00 are considered nonreactive by the ARCHITECT Syphilis TP assay.
- Specimens with S/CO values ≥ 1.00 are considered reactive by the ARCHITECT Syphilis TP assay.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

■ LIMITATIONS OF THE PROCEDURE

- False positive results can be expected with any test kit. The proportion of these falsely reactive specimens is dependent upon the specificity of the test kit, specimen integrity, and the characteristics of the local population being screened.
- If the ARCHITECT Syphilis TP results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- No diagnostic test provides absolute assurance that a sample does not contain low levels of antibodies to TP, such as those present at a very early stage of infection. Therefore, a negative result at any time does not preclude the possibility of exposure to infection with syphilis. Additional information may be required for diagnosis.

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Precision

The ARCHITECT Syphilis TP assay precision is $\leq 15\%$ for the positive control. Precision was determined as described in the Clinical and Laboratory Standards Institute (CLSI), Protocol EP5-A2.⁶ Six samples consisting of four plasma based panels and both of Syphilis TP Controls were assayed in replicates of two at two separate times per day for twenty days ($n=80$ for each sample), using three lots of reagents. Data from this study are summarized in the following table.*

ARCHITECT Syphilis TP Precision

Sample	Lot	Mean S/CO	Repeatability ^a		Within-laboratory ^b	
			SD	%CV	SD	%CV
Negative Control	1	0.04	0.002	5.4	0.002	5.4
	2	0.04	0.005	14.3	0.005	14.7
	3	0.03	0.004	13.1	0.004	13.6
Positive Control	1	2.66	0.048	1.8	0.060	2.3
	2	2.65	0.066	2.5	0.094	3.6
	3	2.66	0.048	1.8	0.069	2.6
Panel D1	1	6.49	0.108	1.7	0.133	2.0
	2	6.57	0.099	1.5	0.144	2.2
	3	6.62	0.116	1.8	0.143	2.2
Panel D2	1	3.74	0.074	2.0	0.092	2.5
	2	3.73	0.059	1.6	0.091	2.4
	3	3.78	0.068	1.8	0.089	2.3
Panel D3	1	1.89	0.034	1.8	0.043	2.3
	2	1.88	0.035	1.9	0.046	2.6
	3	1.89	0.034	1.8	0.048	2.5
Panel D4	1	0.75	0.014	1.8	0.019	2.5
	2	0.74	0.017	2.3	0.020	2.8
	3	0.75	0.015	2.0	0.018	2.5

* Repeatability corresponds to within-run variance.

^b Within-laboratory variance is an accumulation of within-run, between-run and between-day variance.

* Representative data; results in individual laboratories may vary from these data.

Specificity

A study to evaluate performance was conducted at one internal site. Of 5174 BD specimens, 3 specimens were confirmed falsely reactive by ARCHITECT Syphilis TP assay, after resolution. Of 1262 hospitalized/diagnostic specimens, three specimen were falsely reactive by ARCHITECT Syphilis TP assay. The data from this study are summarized in the following table.*

Specificity Results Using Random Blood Donors and Hospitalized Patients

Category	n	IR [%]	RR [%]	Specificity	95% Confidence Interval
Overall Blood Donors	5174	3 [0.06]	3 [0.06]	99.94% (5171/5174)	99.83 - 99.99%
Blood Donor Plasma	2587	3 [0.12]	3 [0.12]	99.88% (2584/2587)	99.66 - 99.98%
Blood Donor Serum	2587	0 [0.00]	0 [0.00]	100.00% (2587/2587)	99.86 - 100%
Hospitalized/Diagnostic Specimens	1262	3 [0.24]	3 [0.24]	99.76% (1259/1262)	99.31 - 99.95%

* Representative data; results in individual laboratories may vary from these data.

Sensitivity

The ARCHITECT Syphilis TP assay demonstrated a sensitivity of $\geq 99.0\%$ in a study testing samples that were confirmed as true positive. The data from this study are summarized in the following table.*

Sensitivity Result

n	Lot	Reactive	Nonreactive	Observed Sensitivity
409	1	409	0	100%
409	2	409	0	100%
409	3	409	0	100%
409	4	409	0	100%

* Representative data; results in individual laboratories may vary from these data.

Analytical Specificity

This study was performed on the ARCHITECT i2000SR System.

The ARCHITECT Syphilis TP assay was evaluated for potential cross reactivity from specimens from individuals with medical conditions, and other disease states unrelated to a syphilis infection.

Six specimens reactive in ARCHITECT Syphilis TP and a treponemal chemiluminescent immunoassay (TP-CLIA) were also either positive (marked by * in table below) or indeterminate (marked by ° in table below) by confirmation testing. Confirmed reactive results were not unexpected because the specimens had been obtained from vendors based on the required disease condition and documentation of prior syphilis laboratory test results was not provided by the vendor.

One of the 10 CMV IgG specimens, 1 of the 10 HIV specimens, 1 of the 7 monoclonal hyper IgG, and 1 of the 10 anti-*E. coli* specimens (marked by ° in table below) were reactive in ARCHITECT Syphilis TP and were not confirmed with other treponemal or nontreponemal tests applied in this study.

One of the 11 gonorrhea specimens and 2 of the 6 HTLV-II specimens (marked by ° in table below) were reactive in ARCHITECT Syphilis TP but could not be tested by the other tests (TP-CLIA, RPR, TP-PA, Fluorescent Treponemal Antibody Absorption [FTA-ABS]) because specimens from these disease states could only be sourced as plasma specimens.

Clinical Category	N	Number of ARCHITECT Syphilis TP Reactive [n]
Chlamydia	15	1*
Cytomegalovirus (CMV) IgG	10	1*
Cytomegalovirus (CMV) IgM	10	1*
Anti-dsDNA Antinuclear	3	0
Epstein-Barr Virus (EBV) IgG	10	1*
Epstein-Barr Virus (EBV) IgM	24	0
Anti- <i>Escherichia coli</i> (E. coli)	10	1°
Gonorrhea	11	1*
HAVAb IgG	10	0
HIV IgM	4	0
Hemodialysis	10	0
Hepatitis A Virus (HAV)	10	1*
Hepatitis B Virus (HBV)	10	0
Hepatitis C Virus (HCV)	10	1*
Herpes Simplex Virus (HSV)	10	1*
Human Anti-Mouse Antibodies (HAMA)	10	0
Human Immunodeficiency Virus (HIV)	10	1*
Human T-Lymphotropic Virus-I (HTLV-I)	10	0
Human T-Lymphotropic Virus-II (HTLV-II)	6	2°

Clinical Category	N	Number of ARCHITECT Syphilis TP Reactive Results
Influenza Vaccine Recipient	20	0
Leptospirosis	6	0
Leptospirosis IgM	5	0
Lyme Disease	10	0
Monoclonal Hyper IgG	7	1b
Anti-Nuclear Antibody (ANA)	10	0
Polyclonal Hyper IgG	3	0
Pregnant	90	0
Rheumatoid Factor	10	0
Rubella IgG	10	0
Systemic Lupus Erythematosus	10	0
Toxoplasma gondii IgG	12	0
Toxoplasma gondii IgM	3	0
Transplant Recipient	10	0
Varicella Zoster Virus	10	0
Total	409	1b

^a Specimen was confirmed to be positive by other tests (TP-CLIA, RPR, TP-PA, FTA-ABS).

^b Specimen was not confirmed positive by other tests (TP-CLIA, RPR, TP-PA, FTA-ABS).

^c Specimen was reactive by ARCHITECT Syphilis TP but was not confirmed due to tube type limitations of the other tests (TP-CLIA, RPR, TP-PA, FTA-ABS).

^d Specimen was positive by TP-CLIA and indeterminate by TP-PA and FTA-ABS.

Interference

Potential interference from elevated levels of triglycerides, bilirubin, protein and hemoglobin in the ARCHITECT Syphilis TP assay is < 0.40 S/CO difference on negative specimens and < 20% S/CO difference on positive samples with the following concentrations.

Test Compound	Test Concentration
Triglycerides	3000 mg/dL
Bilirubin	20 mg/dL
Protein	12 g/dL
Hemoglobin	500 mg/dL

BIBLIOGRAPHY

- Meyer JC. Laboratory Diagnosis of Syphilis. *Curr Probl Dermatol*. 1996; 24: 1-11.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
- U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. *Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PTs)*. November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> Accessed October 24, 2017.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
ASSAY DILUENT	Assay Diluent
CONJUGATE	Conjugate
CONTAINS AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL NO.	Control Number
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF GERMANY	Product of Germany
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WARNING: EYE IRRITANT	Warning: Causes serious eye irritation.
WARNING: SENSITIZER	Warning: May cause an allergic reaction.
WASH BUFFER	Wash Buffer

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Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

Customer Service: Contact your local representative or find country-specific contact information on www.abbott diagnostics.com

Revised July 2018.

©2017, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITEC Syphilis TP Reagent (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8 in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in autorefrigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 12 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP) as amended, 1907/2006/EC (REACH) as amended, Regulation (EC) 1069/2009, Regulation (EC) 142/2011, CLASS 2013 (as amended), GHS (as amended), WHSR Part 7.1 (as amended), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980, EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23640, EN 13975, EN ISO 17511

POTENTIAL USER

The ARCHITEC Syphilis TP Reagent is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals and laboratories, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRA-INDICATIONS FOR USE

Not applicable.

STORAGE AND USE CONDITION

Storage of Medical device at a temperature of 2 °C to 8 °C should be carried out in cold rooms or in refrigerators providing a regulated temperature regime with a daily temperature record.

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1. Meyer JC. Laboratory Diagnosis of Syphilis. *Curr Probl Dermatol*. 1996; 24: 1-11.
2. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
6. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
7. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. *Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PTs)*. November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> По состоянию на: 24 октября 2017 г.

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Signed: ZK

Zaman Khan.

Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Revised June 2016.

INTENDED USE

The ARCHITECT Syphilis TP Calibrator is for the calibration of the ARCHITECT iSystem when used for the qualitative detection of antibody to *Treponema pallidum* (TP) in human serum and plasma. Refer to the ARCHITECT Syphilis TP reagent package insert for additional information.

CONTENTS

1 Bottle (4.0 mL) of ARCHITECT Syphilis TP Calibrator is prepared in recalcified human plasma (inactivated); reactive for anti-TP. Preservatives: sodium azide and other antimicrobial agents.

Calibrator	Color
CAL 1	Green*

* Dye: Green (Acid Yellow No.23 and Acid Blue No.9)

STANDARDIZATION

The calibrator is referenced to an Abbott internal reference standard. This internal reference standard is manufactured by dilution of high titer positive plasma with HEPES buffer with nonreactive human plasma.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

CAUTION: This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.^{1,4}

- The human plasma used in the Calibrator is reactive for anti-TP and nonreactive for anti-HCV, HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CAL 1	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

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STORAGE

- Calibrator is stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Calibrator must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

PREPARATION FOR USE

Calibrator may be used immediately after removal from 2-8°C storage.

Prior to use, mix by gentle inversion (5-10 times).

After each use, tightly close the cap and return the calibrator to 2-8°C storage.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult Instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Calibrator 1
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Germany
	List Number
	For use by or on the order of a physician only (applicable to USA classification only).

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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Revised June 2016.

INTENDED USE

The ARCHITECT Syphilis TP Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT ISystem when used for the qualitative detection of antibody to *Treponema pallidum* (TP) in human serum and plasma. Refer to the ARCHITECT Syphilis TP reagent package insert for additional information.

CONTENTS

2 Bottles (8.0 mL each) of ARCHITECT Syphilis TP Controls:
Negative Control and Positive Control prepared in recalcified human plasma. The Positive Control (inactivated) is reactive for anti-TP.
Preservatives: sodium azide and other antimicrobial agents.

Control	Color	Control Range S/CO
CONTROL -	Natural	≤ 0.40
CONTROL +	Blue*	1.25 - 3.75

* Dye: Acid Blue No. 9

PRECAUTIONS

- **BIOHAZARD**
- For In Vitro Diagnostic Use
- **FOR IN VITRO USE ONLY**

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in the Negative Control is nonreactive for anti-TP, anti-HCV, HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HIV-1/HIV-2.
- The human plasma used in the Positive Control is reactive for anti-TP and nonreactive for anti-HCV, HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CONTROL - / CONTROL +	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbott diagnostics.com or contact your local representative.

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STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Controls must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

PREPARATION FOR USE

Controls may be used immediately after removal from 2-8°C storage. Prior to use, mix by gentle inversion (5-10 times). After each use, tightly close the caps and return the controls to 2-8°C storage.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult Instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL +	Positive Control
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF GERMANY	Product of Germany
RANGE	Range
REF	List Number
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

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Signed: Zaman Khan

Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
B4J2S0

G45772R06

Read Highlighted Changes
Revised May 2019

HIV Ag/Ab Combo


Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Package Insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Key to symbols used

REF	List Number	CONTROL NO.	Control Number
IVD	<i>In Vitro</i> Diagnostic Medical Device	REAGENT LOT	Reagent Lot
LOT	Lot Number	REACTION VESSELS	Reaction Vessels
SN	Serial Number	SAMPLE CUPS	Sample Cups
	Expiration Date	SEPTUM	Septum
	Store at 2-8°C	REPLACEMENT CAPS	Replacement Caps
	Consult instructions for use	CONTAINS: AZIDE	Contains sodium azide. Contact with acids liberates very toxic gas.
	Manufacturer	GTIN	Global Trade Item Number
		PRODUCT OF GERMANY	Product of Germany
		WARNING: SEVERE IRRITANT	Warning: Severe Irritant

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

NAME

ARCHITECT HIV Ag/Ab Combo

INTENDED USE

The ARCHITECT HIV Ag/Ab Combo assay is a chemiluminescent microparticle immunoassay (CMIA) for the simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum or plasma including specimens collected post-mortem (non-heart-beating). The ARCHITECT HIV Ag/Ab Combo assay is intended to be used as an aid in the diagnosis of HIV-1/HIV-2 infection and as a screening test to prevent transmission of HIV-1/HIV-2 to recipients of blood, blood components, cells, tissue and organs. An ARCHITECT HIV Ag/Ab Combo result does not distinguish between the detection of HIV p24 antigen, HIV-1 antibody, or HIV-2 antibody reactivity.

SUMMARY AND EXPLANATION OF TEST

Acquired immunodeficiency syndrome (AIDS) is caused by two types of human immunodeficiency viruses, collectively designated HIV.¹⁻⁷

HIV is the etiologic agent of AIDS.^{1,3,9,7} HIV is transmitted by sexual contact, exposure to blood or blood products, and prenatal infection of a fetus or perinatal infection of a newborn.⁸ Antibodies against HIV are nearly always detected in AIDS patients and HIV infected asymptomatic individuals,^{8,9} and HIV infection is always detected in AIDS patients and seropositive individuals by culture or amplification of viral RNA and/or proviral DNA.^{8,10}

Phylogenetic analysis classifies HIV-1 into groups M (major), N (non-M, non-O), and O (outlier).^{4,5} Group M viruses have spread throughout the world to cause the global AIDS pandemic. In contrast, groups N and O are relatively rare and endemic to west central Africa.¹¹⁻¹⁷ However, group O infections have been identified in Europe and the USA.¹⁸⁻²² HIV-1 group M is composed of genetic subtypes (A, B, C, D, F, G, H, J, and K) and circulating recombinant forms (CRFs).^{5,23} The geographic distribution and regional predominance of HIV-1 subtypes and CRFs vary. All subtypes and many recombinant strains exist in Africa with CRF02_AG the predominant strain in west and west central Africa, subtypes A, C, and D predominant in east central Africa, and subtype C predominant in southern Africa.²³⁻²⁸ HIV-1 subtype B is the predominant subtype in the USA, Europe, Japan, and Australia. However, a significant percentage of new HIV-1 infections in Europe are caused by non-B subtypes.^{29,30} In Asia, subtype C is found in India, and CRF01_AE (formerly called subtype E) and subtype B are in Thailand and southeast Asia.³¹ South America predominantly has subtypes B and F.^{32,33}

Human immunodeficiency virus type 2 (HIV-2) is similar to HIV-1 in its structural morphology, genomic organization, cell tropism, *in vitro* cytopathogenicity, transmission routes, and ability to cause AIDS.^{8,8} However, HIV-2 is less pathogenic than HIV-1, and HIV-2 infections have a longer latency period with slower progression to disease, lower viral titers, and lower rates of vertical and horizontal transmission.³⁴⁻³⁷ HIV-2 is endemic to west Africa but HIV-2 infections, at a low frequency compared to HIV-1, have been identified in the USA, Europe, Asia, and other regions of Africa.^{31,37} HIV-2 is classified into genetic subtypes A-G with most infections caused by subtypes A and B.^{38,39}

The key immunogenic protein and antigenic target for serodetection of HIV infection is the viral (HIV) transmembrane protein (TMP). Antibodies against the TMP (anti-TMP) consistently are among the first to appear at seroconversion of HIV infected individuals.^{9,40-44} The anti-TMP response remains relatively strong throughout the course of the disease, as evidenced by the near universal presence of antibodies against the TMP in asymptomatic and symptomatic stages of HIV infection.^{9,40-44} TMPs from HIV-1 groups M and O and HIV-2 are represented in ARCHITECT HIV Ag/Ab Combo reagents by five recombinant antigens and two synthetic peptides derived from native TMP sequences. The rationale for including three pairs of TMPs is derived from the genetic diversity within HIV-1 and between HIV-1 and HIV-2.^{4,5,45-48} Serologic studies indicate that although HIV-1 and HIV-2 share multiple common epitopes in their core antigens, the envelope glycoproteins are much less cross-reactive.^{7,47-51} Antibodies elicited against the TMP (or portions of the TMP) of a viral strain within one group or type may react well, poorly, or not at all with the TMP (or portions of the TMP) from a viral strain of a different group or type.^{15,52-57} An exception may be antibodies elicited against HIV-1 group N.^{11,12}

Early after infection with HIV, but prior to seroconversion, HIV antigen(s) may be detected in serum or plasma specimens.⁵⁷⁻⁶⁵ The HIV structural protein most often used as the marker of antigenemia is the core protein, p24. The ARCHITECT HIV Ag/Ab Combo uses anti-HIV p24 in the reagents to detect HIV p24 antigen prior to seroconversion, thereby decreasing the seroconversion window and improving early detection of HIV infection.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT HIV Ag/Ab Combo assay is a two-step immunoassay to determine the presence of HIV p24 antigen and antibodies to HIV-1 (Group M and Group O) and HIV-2 in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemilux.

In the first step, sample, ARCHITECT i Wash Buffer, assay diluent, and paramagnetic microparticles are combined. HIV p24 antigen and HIV-1/HIV-2 antibodies present in the sample bind to the HIV-1/HIV-2 antigen and HIV p24 monoclonal (mouse) antibody coated microparticles. After washing, the HIV p24 antigen and HIV-1/HIV-2 antibodies bind to the acridinium-labeled conjugates (HIV-1/HIV-2 antigens [recombinant], synthetic peptides, and HIV p24 antibody [mouse, monoclonal]). Following another wash cycle, pre-trigger and trigger solutions are added to the reaction mixture. The resulting chemiluminescent reaction is measured as relative light units (RLUs). A direct relationship exists between the amount of HIV antigen and antibodies in the sample and the RLUs detected by the ARCHITECT i System optics. The presence or absence of HIV p24 antigen or HIV-1/HIV-2 antibodies in the specimen is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from an ARCHITECT HIV Ag/Ab Combo calibration. Specimens with signal to cutoff (S/CO) values greater than or equal to 1.00 are considered reactive for HIV p24 antigen or HIV-1/HIV-2 antibodies. Specimens with S/CO values less than 1.00 are considered nonreactive for HIV p24 antigen or HIV-1/HIV-2 antibodies.

Specimens that are initially reactive in the ARCHITECT HIV Ag/Ab Combo assay should be retested in duplicate. Repeat reactivity is highly predictive of the presence of HIV p24 antigen and HIV-1/HIV-2 antibodies. However, as with all immunoassays, the ARCHITECT HIV Ag/Ab Combo assay may yield nonspecific reactions due to other causes, particularly when testing in low prevalence populations. A repeatedly reactive specimen should be investigated further with sensitive, supplemental HIV-specific tests, such as immunoblots, antigen tests, and HIV nucleic acid tests. Supplemental testing of repeat-reactive specimens obtained from individuals at risk for HIV infection usually confirms the presence of HIV antibodies or HIV antigen, and HIV nucleic acid. A full differential diagnostic work-up for the diagnosis of AIDS and AIDS-related conditions includes an examination of the patient's immune status and a clinical history.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Reagent Kit, 100 Tests, 500 Tests, 4 x 100 Tests, 4 x 500 Tests

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT i Systems. Please contact your local distributor.

ARCHITECT HIV Ag/Ab Combo Reagent Kit (4J27)

- **[MICROPARTICLES]** 1 or 4 Bottle(s) (6.8 mL per 100 test bottle/27.0 mL per 500 test bottle) Microparticles: HIV-1/HIV-2 antigen (recombinant) and HIV p24 antibody (mouse, monoclonal) coated microparticles in TRIS buffered saline. Minimum concentration: 0.07% solids. Preservative: Sodium Azide.
- **[CONJUGATE]** 1 or 4 Bottle(s) (5.9 mL per 100 test bottle/26.3 mL per 500 test bottle) Conjugate: Acridinium-labeled HIV-1 antigens (recombinant), acridinium-labeled HIV-1/HIV-2 synthetic peptides, and acridinium-labeled HIV p24 antibody (mouse, monoclonal) conjugates in phosphate buffer with protein (bovine) and surfactant stabilizers. Minimum concentration: 0.05 µg/mL. Preservative: Sodium Azide.
- **[ASSAY DILUENT]** 1 or 4 Bottle(s) (5.9 mL per 100 test bottle/26.3 mL per 500 test bottle) Assay Diluent: HIV Ag/Ab Combo assay diluent containing TRIS buffer. Preservative: Sodium Azide.

Other Reagents

ARCHITECT i Pre-Trigger Solution

- **[PRE-TRIGGER SOLUTION]** Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

ARCHITECT i Trigger Solution

- **[TRIGGER SOLUTION]** Trigger Solution containing 0.35 N sodium hydroxide.

ARCHITECT i Wash Buffer

- **[WASH BUFFER]** Wash Buffer containing phosphate buffered saline solution. Preservatives: Antimicrobial Agents.

WARNINGS AND PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Safety Precautions

- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials are considered potentially infectious and be handled in accordance with the OSHA Standard on Bloodborne Pathogens,⁵⁶ Biosafety Level 2⁵⁷ or other appropriate biosafety practices^{58,59} should be used for materials that contain or are suspected of containing infectious agents.

The following warnings and precautions apply to this component:

- Assay Diluent

WARNING: Contains Polyethylene glycol octylphenyl ether (Triton X-100) and Sodium Azide
H318 Causes serious eye irritation.
EU032 Contact with acids liberates very toxic gas

Prevention

P264 Wash hands thoroughly after handling.
P280 Wear protective gloves / protective clothing / eye protection.

Response

P305+P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
+P338
P337+P313 If eye irritation persists: Get medical advice / attention.

This material and its container must be disposed of in a safe way.

- This product contains sodium azide. For a specific listing, refer to the **REAGENTS** section. Contact with acids liberates very toxic gas. This material and its container must be disposed of in a safe way.
- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.
- For information on the safe disposal of sodium azide and a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Handling Precautions

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between reagent kits.
- Prior to loading the ARCHITECT HIV Ag/Ab Combo Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE**, Assay Procedure section of this package insert.
- Septums **MUST** be used to prevent reagent evaporation and contamination, and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.
- To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
- Once a septum has been placed on an open reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts, which have no effect on assay efficacy.
- For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Storage Instructions

-8°C

- **2°C** The ARCHITECT HIV Ag/Ab Combo Reagent Kit must be stored at 2-8°C and may be used immediately after removal from 2-8°C storage. The reagent kit must be stored in an upright position.
- The ARCHITECT HIV Ag/Ab Combo Reagent Kit may be stored on board the ARCHITECT i System for a maximum of 30 days. After 30 days, the reagent kit must be discarded. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

- When stored and handled as directed, reagents are stable until the expiration date.
- Reagents may be stored on or off the ARCHITECT i System. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded. For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

- Before performing the assay, the ARCHITECT HIV Ag/Ab Combo assay file must be installed on the ARCHITECT i System obtained from an ARCHITECT i System Assay CD-ROM. For detailed information on assay file installation and on viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.
- For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.
- For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Verified Specimen Types

- Human serum (including serum collected in serum separator tubes)
- Plasma collected in:
 - potassium EDTA
 - sodium citrate
 - sodium heparin
 - ACD, CPDA-I, CPD
 - lithium heparin
 - potassium oxalate
 - plasma separator tubes
- The ARCHITECT i System does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the ARCHITECT HIV Ag/Ab Combo assay.
- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating), for details refer to section **TESTING OF CADAVERIC BLOOD SPECIMENS**.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed (> 500 mg/dL)
 - obvious microbial contamination
 - body fluids other than human serum and plasma
- Ensure complete clot formation in serum specimens before centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting time.
- Specimens from heparinized patients may be partially coagulated and contain fibrin. Draw the specimen prior to heparin therapy.
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.
- Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
- For optimal results, inspect all specimens for bubbles. Remove bubbles with an applicator stick prior to analysis. Use a new applicator stick for each specimen to prevent cross contamination.
- No qualitative performance differences were observed between experimental controls and more than 20 nonreactive or more than 20 spiked reactive specimens tested with elevated levels of bilirubin (20 mg/dL), triglycerides (3000 mg/dL), protein (4 - 12 g/dL), red blood cells (0.4% v/v), or hemoglobin (500 mg/dL).

Preparation for Analysis

- Specimens must be separated from clots or red blood cells using centrifugation as recommended by the tube manufacturer. Gravity separation is not sufficient for specimen preparation.
- To ensure consistency in results, specimens containing particulate matter or red blood cells, specimens that have been thawed, and specimens that require retesting must be transferred to a centrifuge tube and centrifuged at $\geq 10,000$ RCF (Relative Centrifugal Force) for 10 minutes prior to testing.
- Mix thawed specimens by inverting 10 times. Visually inspect the specimens for the absence of stratification. If layering or stratification is observed repeat inversion cycles until specimens are visibly homogeneous. Centrifuge prior to testing.
- Centrifuged specimens with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified specimen without the lipemic material.

Storage

- Serum or plasma specimens should be stored for no longer than 3 days at room temperature (15 to 30°C) or 14 days at 2 to 8°C following specimen collection. If a storage period greater than 3 days at room temperature or 14 days at 2 to 8°C is anticipated, the specimens should be removed from the clot, red blood cells, or separator gel, and the serum or plasma should be stored frozen at -20°C or colder.
- No qualitative performance differences were observed between experimental controls and 25 nonreactive or 25 spiked reactive specimens subjected to 6 freeze/thaw cycles; however, multiple freeze/thaw cycles should be avoided.

Shipping

- Before shipping specimens, it is recommended that specimens be removed from the clot, separator, or red blood cells.
- When shipped, specimens must be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances. Specimens may be shipped at 2-8°C (wet ice), or at -20°C or colder (dry ice). Do not exceed the 14-day storage time for specimens shipped at 2-8°C (wet ice).

TESTING OF CADAVERIC BLOOD SPECIMENS

- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating) that have been collected up to 17.5 hours after death. Performance was established using 50 spiked and 50 non-spiked cadaveric blood specimens²⁰.
- Testing of cadaveric blood specimens from patients with plasma dilution due to transfusions of > 2000 mL of blood or colloids within 48 hours, or > 2000 mL of crystalloids within 1 hour (or any combination thereof) prior to collection of the specimens have not been validated.
- Follow general standards and/or regulations for collection, storage and handling.
- Follow the tube manufacturer's processing instructions for serum or plasma collection tubes. After initial centrifugation, transfer the supernatant to a centrifuge tube and centrifuge at 10,000 RCF (Relative Centrifugal Force) for 10 minutes. If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.
- Cadaveric blood specimens can be stored for up to 14 days at 2-8°C or up to 3 days at 15-30°C following collection.
- No qualitative differences were observed for cadaveric blood specimens (nonreactive or spiked reactive) when subjected to up to 3 freeze/thaw cycles. However, multiple freeze/thaw cycles should be avoided.

PROCEDURE

Materials Provided:

- 4J27 ARCHITECT HIV Ag/Ab Combo Reagent Kit

Materials Required but not Provided:

- ARCHITECT i System
- ARCHITECT HIV Ag/Ab Combo Assay file, may be obtained from:
 - ARCHITECT i System e-Assay CD-ROM found on www.abbottdiagnostics.com
 - ARCHITECT i System Assay CD-ROM

- 4J27-03 ARCHITECT HIV Ag/Ab Combo Calibrator
- 4J27-12 ARCHITECT HIV Ag/Ab Combo Controls
- ARCHITECT i **PRE-TRIGGER SOLUTION**
- ARCHITECT i **TRIGGER SOLUTION**
- ARCHITECT i **WASH BUFFER**
- ARCHITECT i **REACTION VESSELS**
- ARCHITECT i **SAMPLE CUPS**
- ARCHITECT i **SEPTUM**
- ARCHITECT i **REPLACEMENT CAPS**
- Pipettes or pipette tips (optional) to deliver the volumes specified on the Patient or Control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the ARCHITECT HIV Ag/Ab Combo Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment:
 - Invert the microparticle bottle 30 times.
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles remain adhered to the bottle, continue inverting the bottle until the microparticles have been completely resuspended.
 - If the microparticles do not resuspend, DO NOT USE.
 - Once the microparticles have been resuspended, remove and discard the cap. Wearing clean gloves, remove a septum from the bag. Carefully snap the septum onto the top of the bottle.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Load the ARCHITECT HIV Ag/Ab Combo Reagent Kit on the ARCHITECT i System.
 - Verify that all necessary reagents are present. Ensure that septums are present on all reagent bottles.
- The minimum sample cup volume is calculated by the system and is printed on the Orderlist report. No more than 10 replicates may be sampled from the same sample cup. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
 - Priority: 150 μ L for the first ARCHITECT HIV Ag/Ab Combo test plus 100 μ L for each additional HIV Ag/Ab Combo test from the same sample cup.
 - ≤ 3 hours on board: 150 μ L for the first ARCHITECT HIV Ag/Ab Combo test plus 100 μ L for each additional ARCHITECT HIV Ag/Ab Combo test from the same sample cup.
 - > 3 hours on board: additional sample volume is required. For additional information on sample evaporation and volumes, refer to the ARCHITECT System Operations Manual, Section 5.
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient sample volume is present.
- Prepare calibrator and controls.
 - ARCHITECT HIV Ag/Ab Combo Calibrator 1 and Controls should be mixed by gentle inversion prior to use.
 - To obtain the recommended volume requirements for the ARCHITECT HIV Ag/Ab Combo Calibrator 1 and Controls, hold the bottles vertically and dispense 20 drops of calibrator or 10 drops of each control into each respective sample cup.
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For optimal performance, it is important to follow the routine maintenance procedures defined in the ARCHITECT System Operations Manual, Section 9.

Specimen Dilution Procedures

Specimens cannot be diluted for the ARCHITECT HIV Ag/Ab Combo assay.

Calibration

- To perform an ARCHITECT HIV Ag/Ab Combo calibration, test Calibrator 1 in replicates of three. Calibrator 1 should be priority loaded. A single sample of each ARCHITECT HIV Ag/Ab Combo Control must be tested to evaluate the assay calibration. Ensure that assay control values are within the S/CO ranges specified in the control package insert.
- Once an ARCHITECT HIV Ag/Ab Combo calibration is accepted and stored, all subsequent samples may be tested without further calibration unless one or both of the following occur:
 - A reagent kit with a new lot number is used
 - Controls are out of range
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 5.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT HIV Ag/Ab Combo assay is that a single sample of each control be tested once every 24 hours each day of use.

If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

The ARCHITECT HIV Ag/Ab Combo Control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

It is recommended that each laboratory establish a control range for the ARCHITECT HIV Ag/Ab Combo Positive Control 1 when a new lot of ARCHITECT HIV Ag/Ab Combo reagents is used.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B. The ARCHITECT HIV Ag/Ab Combo assay belongs to method group 5.

RESULTS

The ARCHITECT i System calculates the cutoff (CO) using the mean chemiluminescent signal (RLU) from three replicates of the Calibrator 1 and stores the result.

Calculation

The ARCHITECT i System calculates for the ARCHITECT HIV Ag/Ab Combo assay a result based on the ratio of the sample RLU(s) to the cutoff RLU for each specimen and control.

- Cutoff (CO) = Calibrator 1 mean RLU Value $\times$ 0.40
- S/CO = Sample RLU/Cutoff RLU
- The cutoff RLU is stored for each reagent lot calibration.

Interpretation of Results

- Specimens with S/CO values < 1.00 are considered nonreactive (NR).
- Specimens with S/CO values $\geq$ 1.00 are considered reactive (R).

NOTE: All specimens that are initially reactive must be centrifuged and retested in duplicate. Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section in this package insert.

ARCHITECT HIV Ag/Ab Combo Results

Initial Results (S/CO)	Retest Results	Final Result	Interpretation
R	Both tests are NR	NR	HIV p24 Ag and/or HIV-1/HIV-2 Ab not detected
R	One or both tests are reactive	R	Presumptive evidence of HIV p24 Ag and/or HIV-1/HIV-2 Ab; perform a supplemental assay
NR	No retest required	NR	HIV p24 Ag and/or HIV-1/HIV-2 Ab not detected.

- The Interpretation of Results for specimens with a final result of reactive by the ARCHITECT HIV Ag/Ab Combo assay and indeterminate by supplemental testing is unclear; further clarification may be obtained by testing another specimen taken three to six weeks later.
- ARCHITECT HIV Ag/Ab Combo and supplemental assay results should be interpreted in conjunction with the patient's clinical presentation, history, and other laboratory results.

NOTE: For details on configuring the ARCHITECT i System to use grayzone interpretations, refer to the ARCHITECT System Operations Manual, Section 2.

Flags

- Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

- If the assay results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits that employ mouse monoclonal antibodies.^{71,72} ARCHITECT HIV Ag/Ab Combo reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.⁷³ Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed. Additional information may be required for diagnosis.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT HIV Ag/Ab Combo demonstrated an imprecision of $\leq$ 14% for samples that were three times the cutoff value in a study where three calibrator lots, three control lots, and a panel consisting of four reactive samples were tested. The study was performed at four external sites (France, Italy, Switzerland, Germany) each using one instrument and one internal site using two instruments. Panel members were tested in replicates of three across two reagent lots at the external site and across three reagent lots at the internal site. Each combination of instruments, panel members, and reagent lots was tested in four runs, with the exception of one reagent lot at the internal site that was tested in six runs on one instrument. The intra-run and inter-run standard deviation (SD) and percent coefficient of variation (%CV) were analyzed with a variance components analysis⁷⁴ using a mixed analysis of variance model.⁷⁵ The data from the study are summarized in Table 1.^a

Table 1
ARCHITECT HIV Ag/Ab Combo Precision

Panel Member (units)	n	Mean	Intra-assay		Inter-assay ^a	
			SD	%CV	SD	%CV
Calibrator 1 (RLUs)	432	5629	201.8	3.6	235.7	4.2
NC ^b (S/CO)	1224	0.09	0.019	20.23	0.022	24.11
PC ^c 1 (S/CO)	1224	3.96	0.155	3.91	0.155	3.91
PC ^c 2 (S/CO)	1224	3.57	0.115	3.22	0.161	4.51
PC ^c 3 (S/CO)	1224	3.19	0.107	3.34	0.141	4.41
PM ^d 1 (S/CO)	432	2.57	0.098	3.82	0.137	5.33
PM ^d 2 (S/CO)	432	2.44	0.116	4.76	0.148	6.01
PM ^d 3 (S/CO)	432	2.36	0.075	3.16	0.087	4.12
PM ^d 4 (S/CO)	432	1.44	0.054	3.78	0.054	3.78

^a Inter-assay variability contains intra-assay variability.

^b Negative Control

^c Positive Control

^d Panel Member

^a Representative performance data are shown. Results obtained at individual laboratories may vary.

Specificity

The ARCHITECT HIV Ag/Ab Combo assay demonstrated a specificity of $\geq 99.5\%$ in a study where specimens from a blood donor population with an assumed HIV infection prevalence of zero were tested. The testing was performed at two external sites and one internal site on a total of 6385 serum and plasma specimens collected from four blood-donation centers. Data from this study are summarized in Table 2.*

Table 2
ARCHITECT HIV Ag/Ab Combo Specificity of Blood Donor Specimens

Specimen Type	n	Initially Reactive	Repeat Reactive	Specificity (%)	95% CI ^a
Plasma	2747	3	3	99.89	99.68 - 99.98
Serum	3618	6	4	99.89	99.72 - 99.97
Total	6385	9	7	99.89	99.77 - 99.96

^a Confidence Interval

* Representative performance data are shown. Results obtained at individual laboratories may vary.

Specimens from randomly selected hospitalized patients and specimens containing potentially interfering substances, which includes those from individuals with medical conditions unrelated to HIV infection, were tested at four different sites with the ARCHITECT HIV Ag/Ab Combo assay. Of the 2870 hospitalized patient specimens (HP) and the 322 specimens containing potentially interfering substances (IS), 29 HP and 12 IS specimens were confirmed as having HIV infection by confirmation testing. These specimens were excluded from the study. The data from the remaining 2841 HP specimens and the 310 IS specimens are summarized in Table 3.*

Table 3
ARCHITECT HIV Ag/Ab Combo Specificity of Hospitalized Patients and Interfering Substances

Specimen Type	n	Initial Reactive	Repeat Reactive	Specificity (%)
HP	2841	11	11	99.61
IS ^a	310	1	1	99.68
Total	3151	12	12	99.62

^a The IS specimens belonged to the following categories: Viral Infection (HBV, HSV, CMV, Rubella, HAV, HCV, EBV, HTLV-I, HTLV-II); fungal/yeast/protozoal/bacterial infection (*C. albicans*, *T. pallidum*, *T. gondii*, *E. coli*, *C. trachomatis*, *N. gonorrhoea*); autoimmune (rheumatoid factor [RF], antinuclear antibodies [ANA]); other conditions (pregnant females all trimesters, multiparous females, elevated IgG, elevated IgM, monoclonal gammopathy, flu vaccine recipients, HAMA, hemodialysis patients, hemophiliacs, multiple transfusion recipients).

* Representative performance data are shown. Results obtained at individual laboratories may vary.

One hundred specimens from patients at increased risk for HIV infection were tested by ARCHITECT HIV Ag/Ab Combo. Of these 100 specimens, 70 were repeat reactive by ARCHITECT HIV Ag/Ab Combo. Sixty-nine of these 70 repeat reactive specimens were positive by confirmation testing.* The specimens from patients at increased risk for HIV infection belonged to the following categories: Intravenous drug users, homosexual males, and sexually transmitted diseases.

* These are representative performance data. Results obtained at individual laboratories may vary.

Sensitivity

The ARCHITECT HIV Ag/Ab Combo assay demonstrated the sensitivity values provided in Table 4.* These values were determined in a study where specimens from individuals clinically diagnosed with HIV infection and classified disease status were tested.

Table 4
ARCHITECT HIV Ag/Ab Combo Sensitivity

Antibody Type	Total n	Reactive n	Sensitivity (%)
Anti-HIV-1	520	520	100.0
Anti-HIV-2	111	111	100.0
Anti-HIV gO ^a	6	6	100.0

^a An additional 29 diluted anti-HIV gO specimens were found to be reactive by the ARCHITECT HIV Ag/Ab Combo assay.

* Representative performance data are shown. Results obtained at individual laboratories may vary.

Sensitivity of the ARCHITECT HIV Ag/Ab Combo assay to detect p24 antigen and anti-HIV-1 was evaluated by testing sequential specimens from 31 seroconversion panels. These panels are commercially available and precharacterized for HIV antigen and antibodies. Table 5 shows data from 5 seroconversion panel members.*

Table 5
ARCHITECT HIV Ag/Ab Combo Seroconversion Sensitivity

Panel	Days Since 1 st Bleed	ARCHITECT		Western Blot ^a	HIV Ag ^a (S/CO)	PCR ^a Copies/mL
		HIV Ag/Ab Combo (S/CO)				
BBI	0	0.12		Neg ^b	0.1	4,000
PRB953	3	1.04		Neg	0.6	80,000
	7	11.71		Neg	8.8	700,000
	10	31.79		Neg	19.3	1,000,000
BBI	0	0.17		IND ^c (f24 ^d)	0.3	1,000
PRB955	3	1.86		IND (f24)	1.8	70,000
	7	17.48		IND (f24)	14.0	400,000
	12	24.62		IND (f24)	22.4	> 800,000
	14	34.17		IND (24)	21.3	700,000
BBI	0	1.09		Neg	2.3	80,000
PRB959	7	323.59		Neg	> 40.5	> 500,000
	9	192.13		Neg	> 40.5	> 500,000
	14	75.41		p24, gp160	14.4	> 500,000
	19	88.38		p24, gp160	0.7	200,000
	21	68.29		p24, gp160	0.7	200,000
	26	68.46		p24, gp160	1.1	300,000
BCP	0	0.10		p24	0.255	< 400
6247	2	0.09		p24	0.373	< 400
	7	0.09		p24	0.343	< 400
	9	0.09		p24	0.343	< 400
	14	0.09		p24	0.294	< 400
	16	0.09		p24	0.275	< 400
	21	1.11		p24	1.265	23,500
	23	7.65		p24	8.598	124,000
	28	155.50		p24	37.618	675,000
	30	129.37		Neg	39.941	> 750,000
NABI	0	0.12		IND ^c (p24)	0.45	2,328
SV404	8	0.13		IND (p24)	0.48	2,672
	11	0.13		IND (p24)	0.56	3,079
	15	0.23		IND (p24)	0.52	7,190
	18	1.33		IND (p24)	1.48	54,380
	22	6.51		IND (p24)	4.84	279,800

^a Data from the vendor.

^b Neg = no band or bands being observed.

^c IND = Indeterminant

^d f24 = faint p24 antigen band

* Representative performance data are shown. Results obtained at individual laboratories may vary.

The ARCHITECT HIV Ag/Ab Combo assay has an analytical sensitivity of < 50 pg/mL to HIV-1 p24 Ag. Antigen sensitivity was evaluated with the HIV Ag 2003 AFSSAPS panel across three lots of ARCHITECT HIV Ag/Ab Combo reagents. The results demonstrated an average sensitivity to HIV-1 p24 Ag of 18.06 pg/mL.* In addition, sensitivity was assessed using the HIV-1 p24 Antigen, First International Reference Reagent, NIBSC code: 90/536. The results demonstrated an average sensitivity of 0.87 IU/mL.* The WRAIR (Walter Reed Army Institute of Clinical Research, Rockville, MD) Clade Panel consisting of 32 members with the following subtypes was tested: A, B, B/A/B, C, D, G, F, Group O, CRF02_AG, CRF01_AE. For all panel members, the ARCHITECT HIV Ag/Ab Combo antigen sensitivity was better than antigen sensitivity of a comparative HIV antigen/antibody combo assay.*

* These are representative performance data. Results obtained at individual laboratories may vary.

BIBLIOGRAPHY

1. Barré-Sinoussi F, Chermann JC, Rey F, et al. Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). *Science* 1983;220:868-71.
2. Popovic M, Sarngadharan MG, Read E, et al. Detection, isolation, and continuous production of cytopathic retroviruses (HTLV-III) from patients with AIDS and pre-AIDS. *Science* 1984;224:497-500.
3. Gallo RC, Salahuddin SZ, Popovic M, et al. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 1984;224:500-3.
4. Letner T. Genetic subtypes of HIV-1. In: Myers G, Korber BT, Foley BT et al., editors. *Human Retroviruses and AIDS* 1996. Los Alamos, NM: Los Alamos National Laboratory, 1996:111-28-111-40. (Available on-line at <http://hiv-web.lanl.gov>)
5. Robertson DL, Anderson JP, Bradac JA, et al. HIV-1 nomenclature proposal: a reference guide to HIV-1 classification. In: Kuiken CL, Foley B, Hahn B, et al., editors. *Human Retroviruses and AIDS* 1999. Los Alamos, NM: Los Alamos National Laboratory, 1999:492-505. (Available on-line at <http://hiv-web.lanl.gov>)
6. Clavel F, Guétard D, Brun-Vézinet F et al. Isolation of a new human retrovirus from West African patients with AIDS. *Science* 1986;233:343-6.
7. Clavel F. HIV-2, the West African AIDS virus. *AIDS* 1987;1:135-40.
8. Schochetman G, George JR, editors. *AIDS Testing: A Comprehensive Guide to Technical, Medical, Social, Legal, and Management Issues*. 2nd ed. New York: NY Springer-Verlag; 1994.
9. Sarngadharan MG, Popovic M, Bruch L, et al. Antibodies reactive with human T-lymphotropic retroviruses (HTLV-III) in the serum of patients with AIDS. *Science* 1984;224:506-8.
10. Jackson JB, Kwok SY, Sninsky JJ, et al. Human immunodeficiency virus type 1 detected in all seropositive symptomatic and asymptomatic individuals. *J Clin Microbiol* 1990;28:16-9.
11. Simon F, Mauclet P, Roques P, et al. Identification of a new human immunodeficiency virus type 1 distinct from group M and group O. *Nature Medicine* 1998;4:1032-7.
12. Ayoub A, Souquères S, Njinku B, et al. HIV-1 group N among HIV-1-seropositive individuals in Cameroon. *AIDS* 2000;14:2623-5.
13. Günter LG, Hauser PH, Eberle J, et al. A new subtype of human immunodeficiency virus type 1 (MVP-5180) from Cameroon. *J Virology* 1994;68(3):1581-5.
14. Peeters M, Gueye A, Mboup S, et al. Geographical distribution of HIV-1 group O viruses in Africa. *AIDS* 1997;11:493-8.
15. Hunt JC, Golden AM, Lund JK, et al. Envelope sequence variability and serologic characterization of HIV type 1 group O isolates from Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(12):995-1005.
16. Hackett J Jr, Zekeng L, Brennan CA, et al. Genetic analysis of HIV type 1 group O p24⁹⁹ sequences from Cameroon and Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(13):1155-8.
17. Vanden Haesevelde M, Decourt J, De Leys RJ, et al. Genomic cloning and complete sequence analysis of a highly divergent African human immunodeficiency virus isolate. *J Virology* 1994;68(3):1586-96.
18. Hampel H, Sawitzky D, Stöfler-Mellicke M, et al. First case of HIV-1 subtype O infection in Germany. *Infection* 1995;23:369-70.
19. Loussert-Ajaka I, Chaix M-L, Korber B, et al. Variability of human immunodeficiency virus type 1 group O strains isolated from Cameroonian patients living in France. *J Virology* 1995;69:5640-9.
20. Soriano V, Gutiérrez M, García-Lerma G et al. First case of HIV-1 group O infection in Spain. *Vox Sang* 1996;71:66.
21. Rayfield MA, Sullivan P, Bondea CI, et al. HIV-1 group O virus identified for the first time in the United States. *Emerging Infectious Diseases* 1996;2:209-12.
22. Sullivan PS, Do AN, Robbins K, et al. Surveillance for variant strains of HIV: subtype G and group O HIV-1. *JAMA* 1997;278:292.
23. Peeters M. Recombinant HIV sequences: their role in the global epidemic. In: Kuiken C, Foley B, Hahn B, et al., editors. *HIV Sequence Compendium 2000*. Los Alamos, NM: Los Alamos National Laboratory, 2000:54-73. (Available on-line at <http://hiv-web.lanl.gov>)
24. Montavon C, Tours-Kane C, Liegeois F, et al. Most env and gag subtype A HIV-1 viruses circulating in West and West Central Africa are similar to the prototype AG recombinant virus IBNG. *J Acquir Immune Defic Syndr* 2000;23:363-74.
25. Brennan CA, Lund JK, Golden A, et al. Serologic and phylogenetic characterization of HIV-1 subtypes in Uganda. *AIDS* 1997;11:1823-32.
26. Vidal N, Peeters M, Mulanga-Kabeya C, et al. Unprecedented degree of human immunodeficiency virus type 1 (HIV-1) group M genetic diversity in the Democratic Republic of Congo suggests that the HIV-1 pandemic originated in Central Africa. *J Virology* 2000;74:10498-507.
27. Husseln M, Abebe A, Pollakis G, et al. HIV-1 subtype C in commercial sex workers in Addis Ababa, Ethiopia. *J Acquir Immune Defic Syndr* 2000;23:120-7.
28. Van Harmelen JH, Van Der Ryst E, Loubser AS, et al. A predominately HIV type 1 subtype C-restricted epidemic in South African urban populations. *AIDS Res Hum Retroviruses* 1999;15:395-8.
29. Couturier E, Damond F, Roques P, et al. HIV-1 diversity in France, 1996-1998. *AIDS* 2000;14:269-98.
30. Fransen K, Buve A, Nkengasong JN, et al. Longstanding presence in Belgians of multiple non-B HIV-1 subtypes. *Lancet* 1996;347:1403.
31. Weniger BG, Takebe Y, Ou C, et al. The molecular epidemiology of HIV in Asia. *AIDS* 1994;8(S2):S13-28.
32. Russell KL, Carcamo C, Watts DM, et al. Emerging genetic diversity of HIV-1 in South America. *AIDS* 2000;14:1785-91.
33. Tanuri A, Swanson P, Devare S, et al. HIV-1 subtypes among blood donors from Rio de Janeiro, Brazil. *J Acquir Immune Defic Syndr* 1999;20:60-8.
34. Marlink RG, Ricard D, Mboup S, et al. Clinical, hematologic, and immunologic cross-sectional evaluation of individuals exposed to human immunodeficiency virus type-2 (HIV-2). *AIDS Res Hum Retroviruses* 1988;4:137-48.
35. Simon F, Matheron S, Tamalet C, et al. Cellular and plasma viral load in patients infected with HIV-2. *AIDS* 1993;7:1411-7.
36. Kanki PJ, Travers KU, Mboup S, et al. Slower heterosexual spread of HIV-2 than HIV-1. *Lancet* 1994;343:943-6.
37. De Cock KM, Adjuoro G, Ekpin E, et al. Epidemiology and transmission of HIV-2: why there is no HIV-2 pandemic. *JAMA* 1993;270:2083-5.
38. Gao F, Yue L, Robertson DL, et al. Genetic diversity of human immunodeficiency virus type 2: evidence for distinct sequence subtypes with differences in virus biology. *J Virology* 1994;68:7433-47.
39. Yamaguchi J, Devare SQ, Brennan CA. Identification of a new HIV-2 subtype based on phylogenetic analysis of full-length genomic sequence. *AIDS Res Hum Retroviruses* 2000;16:925-30.
40. Dawson GJ, Heller JS, Wood CA, et al. Reliable detection of individuals seropositive for the human immunodeficiency virus (HIV) by Competitive Immunoassays using *Escherichia coli*-expressed HIV structural proteins. *J Infect Dis* 1988;157:149-55.
41. Montagnier L, Clavel F, Krust B, et al. Identification and antigenicity of the major envelope glycoprotein of lymphadenopathy-associated virus. *Virology* 1985;144:283-9.
42. Barin F, McLane MF, Allan JS, et al. Virus envelope protein of HTLV-III represents major target antigen for antibodies in AIDS patients. *Science* 1985; 228:1094-6.

43. Schulz TF, Aschauer JM, Hengster P, et al. Envelope gene-derived recombinant peptide in the serodiagnosis of human immunodeficiency virus infection. *Lancet* 1986;2:111-2.
44. Allan JS, Colgan JE, Barin F, et al. Major glycoprotein antigens that induce antibodies in AIDS patients are encoded by HTLV-III. *Science* 1985;228:1091-4.
45. HIV-1/HIV-2/SIV Complete genomes. In: Kuiken CL, Foley B, Freed E, et al., editors. *HIV Sequence Compendium 2002*. Los Alamos, NM: Los Alamos National Library; 2002:377-480.
46. Yamaguchi J, Bodelle P, Kaptué L, et al. Near full-length genomes of 15 HIV type 1 group O isolates. *AIDS Res Hum Retroviruses* 2003;19:979-88.
47. Hunt JC, Johnson-Paepke J, Boardway K, et al. Discrimination between HIV-1 and HIV-2 seropositive individuals using mouse monoclonal antibodies directed to HIV transmembrane proteins. *AIDS Res Hum Retroviruses* 1990;6:883-98.
48. Barin F, M'Boup S, Denis F, et al. Serological evidence for virus related to simian T-lymphotropic retrovirus III in residents of West Africa. *Lancet* 1985;2:1387-9.
49. Kanki PJ, Barin F, M'Boup S, et al. New human T-lymphotropic retrovirus related to simian T-lymphotropic virus type III (STLV-III_{mac}). *Science* 1986;232:238-43.
50. Kanki PJ, M'Boup S, Ricard D, et al. Human T-lymphotropic virus type 4 and the human immunodeficiency virus in West Africa. *Science* 1987;236:827-31.
51. Clavel F, Guyader M, Guétard D, et al. Molecular cloning and polymorphism of the human immune deficiency virus type 2. *Nature* 1988;324:691-5.
52. Cot M, Poulain M, Delagneau JF et al. Dual HIV-1 and HIV-2 infection in West Africa supported by synthetic peptide analysis. *AIDS Res Hum Retroviruses* 1988;4:239-41.
53. Simon F, Cot M, Lesager C, et al. Differentiation between HIV-1 and HIV-2 infection by radioimmunoassay and synthetic peptides in double reactive sera. *AIDS* 1989;3:401-4.
54. Norby E, Biberfeld G, Chiodi F, et al. Discrimination between antibodies to HIV and to related retroviruses using site-directed serology. *Nature* 1987;328:248-50.
55. Gürtler LG, Zekeng L, Simon F, et al. Reactivity of five anti-HIV-1 subtype O specimens with six different anti-HIV screening ELISAs and three immunoblots. *J Virol Methods* 1995;51:177-84.
56. Loussert-Ajaka I, Ly TD, Chaix ML, et al. HIV-1/HIV-2 seronegativity in HIV-1 subtype O infected patients. *Lancet* 1994;343:1393-4.
57. Devare SG, Desai SM, Dawson GJ, et al. Diagnosis and monitoring of HIV-1 and HIV-2 infection. In: Khan NC, Melnick JL, editors. *Human Immunodeficiency Virus: Innovative Techniques for Isolation and Identification (Monographs in Virology, vol. 18)*. Basel: S Karger; 1990; 18:105-21.
58. Kessler HA, Blaauw B, Spear J, et al. Diagnosis of human immunodeficiency virus infection in seronegative homosexuals presenting with an acute viral syndrome. *JAMA* 1987;258:1196-9.
59. Phair JP. Human immunodeficiency virus antigenemia. *JAMA* 1987;258:1218.
60. Allain J-P, Laurian Y, Paul DA, et al. Serological markers in early stages of human immunodeficiency virus infection in haemophiliacs. *Lancet* 1986;2(8518):1233-6.
61. Kenny C, Parkin J, Underhill G, et al. HIV antigen testing. *Lancet* 1987;1(8532):565-6.
62. Wall RA, Denning DW, Amos A. HIV antigenaemia in acute HIV infection. *Lancet* 1987;1(8532):566.
63. Stute R. HIV antigen detection in routine blood donor screening. *Lancet* 1987;1(8532):568.
64. Goudsmit J, De Wolf F, Paul DA, et al. Expression of human immunodeficiency virus antigen (HIV-Ag) in serum and cerebrospinal fluid during acute and chronic infection. *Lancet* 1986;2(8500):177-80.
65. von Sydow M, Gaines H, Sonnerberg A, et al. Antigen detection in primary HIV infection. *Brit Med J* 1988;296:238-40.
66. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
67. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2008.
68. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
69. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
70. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PS), November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> Accessed July 01, 2013.
71. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibody therapy. *Clin Chem* 1988;34:261-4.
72. Schreff RW, Faon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-85.
73. Boscasto LM and Stuart MC. Heterophilic antibodies, a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
74. Box GEP, Hunter WG, Hunter JS. *Statistics for experimenters; an introduction to design, data analysis, and model building*. New York, NY: John Wiley & Sons, Inc, 1978:510-39, 571-83.
75. SAS Institute, Inc. SAS Technical Report P-229, *SAS/STAT Software: Changes and enhancements, Release 6.07*. Cary, NC: SAS Institute, 1992:288-366.

The following US Patents are relevant to the ARCHITECT / System or its components. There are other such patents and patent application in the United States and worldwide.

5 468 646	5 543 524	5 545 739
5 565 570	5 669 819	5 783 699

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions.

SAS and SAS/STAT are property of their respective owner.

ABBOTT
Max-Planck-Ring 2
85205 Wiesbaden
Germany
+49-6122-580

May 2019

© 2004, 2019 Abbott Laboratories

■ ADDENDUM

to the Instruction for Use (for Russia)

SHIPPING

Recommended shipping conditions for device correspond to storage conditions provided in the Instructions for Use (2-8°C). Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated (GOST R 51088-2013). Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the warranty shelf life of device which corresponds to period between the date of manufacture and expiration date that are indicated on the label and/or in accompanying documentation. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

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- ARCHITECT i System
- Файл теста ARCHITECT HIV Ag/Ab Combo, доступный на:
 - ARCHITECT i System e-Assay CD-ROM на www.abbottdiagnostics.com
 - ARCHITECT i System Assay CD-ROM
- 4J27-03 Калибратор ARCHITECT HIV Ag/Ab Combo Calibrator
- 4J27-12 Контроли ARCHITECT HIV Ag/Ab Combo Controls
- ARCHITECT i **PRE-TRIGGER SOLUTION**
- ARCHITECT i **TRIGGER SOLUTION**
- ARCHITECT i **WASH BUFFER**
- ARCHITECT i **REACTION VESSELS**
- ARCHITECT i **SAMPLE CUPS**
- ARCHITECT i **SEPTUM**
- ARCHITECT i **REPLACEMENT CAPS**
- Пипетки или наконечники для пипеток (опционально) для объемов, указанных на экране заказа образца пациента или контроля.

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1. Barré-Sinoussi F, Chermann JC, Rey F, et al. Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). *Science* 1983;220:868-71.
2. Popovic M, Sarngadharan MG, Read E, et al. Detection, isolation, and continuous production of cytopathic retroviruses (HTLV-III) from patients with AIDS and pre-AIDS. *Science* 1984;224:497-500.
3. Gallo RC, Salahuddin SZ, Popovic M, et al. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 1984;224:500-3.
4. Leitner T. Genetic subtypes of HIV-1. In: Myers G, Korber BT, Foley BT et al., editors. *Human Retroviruses and AIDS 1996*. Los Alamos, NM: Los Alamos National Laboratory, 1996:III-28-III-40. (Available on-line at <http://hiv-web.lanl.gov>)
5. Robertson DL, Anderson JP, Bradac JA, et al. HIV-1 nomenclature proposal: a reference guide to HIV-1 classification. In: Kuiken CL, Foley B, Hahn B, et al., editors. *Human Retroviruses and AIDS 1999*. Los Alamos, NM: Los Alamos National Laboratory, 1999:492-505. (Available on-line at <http://hiv-web.lanl.gov>)
6. Clavel F, Guétiard D, Brun-Vézinet F et al. Isolation of a new human retrovirus from West African patients with AIDS. *Science* 1986;233:343-6.
7. Clavel F. HIV-2, the West African AIDS virus. *AIDS* 1987;1:135-40.
8. Schochetman G, George JR, editors. *AIDS Testing: A Comprehensive Guide to Technical, Medical, Social, Legal, and Management Issues*. 2nd ed. New York: NY Springer-Verlag; 1994.
9. Sarngadharan MG, Popovic M, Bruch L, et al. Antibodies reactive with human T-lymphotropic retroviruses (HTLV-III) in the serum of patients with AIDS. *Science* 1984;224:506-8.
10. Jackson JB, Kwok SY, Sninsky JJ, et al. Human immunodeficiency virus type 1 detected in all seropositive symptomatic and asymptomatic individuals. *J Clin Microbiol* 1990;28:16-9.
11. Simon F, Mauclet P, Roques P, et al. Identification of a new human immunodeficiency virus type 1 distinct from group M and group O. *Nature Medicine* 1998;4:1032-7.
12. Aypuba A, Souquiere S, Njinku B, et al. HIV-1 group N among HIV-1-seropositive individuals in Cameroon. *AIDS* 2000;14:2623-5.
13. Güntler IG, Hauser PH, Eberle J, et al. A new subtype of human immunodeficiency virus type 1 (MVP-5180) from Cameroon. *J Virology* 1994;68(3):1581-5.
14. Peeters M, Guéye A, Mboup S, et al. Geographical distribution of HIV-1 group O viruses in Africa. *AIDS* 1997;11:493-8.
15. Hunt JC, Golden AM, Lund JK, et al. Envelope sequence variability and serologic characterization of HIV type 1 group O isolates from Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(12):995-1005.
16. Hackett J Jr, Zekeng L, Brennan CA, et al. Genetic analysis of HIV type 1 group O p24⁹⁸⁻⁹⁹ sequences from Cameroon and Equatorial Guinea. *AIDS Res Hum Retroviruses* 1997;13(13):1155-8.
17. Vanden Haesevelde M, Decourt J, De Lays RJ, et al. Genomic cloning and complete sequence analysis of a highly divergent African human immunodeficiency virus isolate. *J Virology* 1994;68(3):1586-96.
18. Hampl H, Sawitzky D, Stöffler-Meilicke M, et al. First case of HIV-1 subtype O infection in Germany. *Infection* 1995;23:369-70.
19. Loussert-Ajaka I, Chaix M-L, Korber B, et al. Variability of human immunodeficiency virus type 1 group O strains isolated from Cameroonian patients living in France. *J Virology* 1995;69:5640-9.
20. Soriano V, Gutiérrez M, García-Lerma G et al. First case of HIV-1 group O infection in Spain. *Vox Sang* 1996;71:66.
21. Rayfield MA, Sullivan P, Bondia CI, et al. HIV-1 group O virus identified for the first time in the United States. *Emerging Infectious Diseases* 1996;2:209-12.
22. Sullivan PS, Do AN, Robbins K, et al. Surveillance for variant strains of HIV: subtype G and group O HIV-1. *JAMA* 1997;278:292.
23. Peeters M. Recombinant HIV sequences: their role in the global epidemic. In: Kuiken C, Foley B, Hahn B, et al., editors. *HIV Sequence Compendium 2000*. Los Alamos, NM: Los Alamos National Laboratory, 2000:54-73. (Available on-line at <http://hiv-web.lanl.gov>)
24. Montavon C, Toure-Kane C, Liegeois F, et al. Most env and gag subtype A HIV-1 viruses circulating in West and West Central Africa are similar to the prototype AG recombinant virus IBNG. *J Acquir Immune Defic Syndr* 2000;23:363-74.
25. Brennan CA, Lund JK, Golden A, et al. Serologic and phylogenetic characterization of HIV-1 subtypes in Uganda. *AIDS* 1997;11:1823-32.
26. Vidal N, Peeters M, Mulanga-Kabeya C, et al. Unprecedented degree of human immunodeficiency virus type 1 (HIV-1) group M genetic diversity in the Democratic Republic of Congo suggests that the HIV-1 pandemic originated in Central Africa. *J Virology* 2000;74:10498-507.
27. Hussein M, Abebe A, Pollekis G, et al. HIV-1 subtype C in commercial sex workers in Addis Ababa, Ethiopia. *J Acquir Immune Defic Syndr* 2000;23:120-7.
28. Van Harmelen JH, Van Der Ryst E, Loubser AS, et al. A predominately HIV type 1 subtype C-restricted epidemic in South African urban populations. *AIDS Res Hum Retroviruses* 1999;15:395-8.
29. Couturier E, Damond F, Roques P, et al. HIV-1 diversity in France, 1986-1998. *AIDS* 2000;14:289-96.
30. Franssen K, Buve A, Nkengasong JN, et al. Longstanding presence in Belgians of multiple non-B HIV-1 subtypes. *Lancet* 1996;347:1403.
31. Weniger BG, Takebe Y, Ou C, et al. The molecular epidemiology of HIV in Asia. *AIDS* 1994;8(S2):S13-28.
32. Russell KL, Carcamo C, Watts DM, et al. Emerging genetic diversity of HIV-1 in South America. *AIDS* 2000;14:1785-91.
33. Tanuri A, Swanson P, Devare S, et al. HIV-1 subtypes among blood donors from Rio de Janeiro, Brazil. *J Acquir Immune Defic Syndr* 1999;20:60-6.
34. Marlink RG, Ricard D, M'Boup S, et al. Clinical, hematologic, and immunologic cross-sectional evaluation of individuals exposed to human immunodeficiency virus type-2 (HIV-2). *AIDS Res Hum Retroviruses* 1988;4:137-48.
35. Simon F, Matheron S, Tamalet C, et al. Cellular and plasma viral load in patients infected with HIV-2. *AIDS* 1993;7:1411-7.
36. Kanki PJ, Travers KU, Mboup S, et al. Slower heterosexual spread of HIV-2 than HIV-1. *Lancet* 1994;343:943-6.
37. DeCock KM, Adjorlolo G, Ekpini E, et al. Epidemiology and transmission of HIV-2: why there is no HIV-2 pandemic. *JAMA* 1993;270:2083-6.
38. Gao F, Yue L, Robertson DL, et al. Genetic diversity of human immunodeficiency virus type 2: evidence for distinct sequence subtypes with differences in virus biology. *J Virology* 1994;68:7433-47.
39. Yamaguchi J, Dovare SG, Brennan CA. Identification of a new HIV-2 subtype based on phylogenetic analysis of full-length genomic sequence. *AIDS Res Hum Retroviruses* 2000;16:925-30.
40. Dawson GJ, Heller JS, Wood CA, et al. Reliable detection of individuals seropositive for the human immunodeficiency virus (HIV) by Competitive Immunoassays using *Escherichia coli*-expressed HIV structural proteins. *J Infect Dis* 1988;157:149-55.
41. Montagnier L, Clavel F, Krust B, et al. Identification and antigenicity of the major envelope glycoprotein of lymphadenopathy-associated virus. *Virology* 1985;144:283-9.
42. Barin F, McLane MF, Allan JS, et al. Virus envelope protein of HTLV-III represents major target antigen for antibodies in AIDS patients. *Science* 1985; 228:1094-6.
43. Schulz TF, Aschauer JM, Hengster P, et al. Envelope gene-derived recombinant peptide in the serodiagnosis of human immunodeficiency virus infection. *Lancet* 1986;2:111-2.
44. Allan JS, Colligan JE, Barin F, et al. Major glycoprotein antigens that induce antibodies in AIDS patients are encoded by HTLV-III. *Science* 1985;228:1091-4.
45. HIV-1/HIV-2/SIV Complete genomes. In: Kuiken CL, Foley B, Freed E, et al., editors. *HIV Sequence Compendium 2002*. Los Alamos, NM: Los Alamos National Library; 2002:377-480.
46. Yamaguchi J, Bodelle P, Kaptué L, et al. Near full-length genomes of 15 HIV type 1 group O isolates. *AIDS Res Hum Retroviruses* 2003;19:979-88.
47. Hunt JC, Johnson-Paepke J, Boardway K, et al. Discrimination between HIV-1 and HIV-2 seropositive individuals using mouse monoclonal antibodies directed to HIV transmembrane proteins. *AIDS Res Hum Retroviruses* 1990;6:883-98.
48. Barin F, M'Boup S, Denis F, et al. Serological evidence for virus related to simian T-lymphotropic retrovirus III in residents of West Africa. *Lancet* 1985;2:1387-9.
49. Kanki PJ, Barin F, M'Boup S, et al. New human T-lymphotropic retrovirus related to simian T-lymphotropic virus type III (STLV-III_{AGM}). *Science* 1986;232:238-43.
50. Kanki PJ, M'Boup S, Ricard D, et al. Human T-lymphotropic virus type 4 and the human immunodeficiency virus in West Africa. *Science* 1987;236:827-31.

51. Clavel F, Guyader M, Guétard D, et al. Molecular cloning and polymorphism of the human immune deficiency virus type 2. *Nature* 1986;324:691-5.
52. Cot M, Poulain M, Delagrange JF et al. Dual HIV-1 and HIV-2 infection in West Africa supported by synthetic peptide analysis. *AIDS Res Hum Retroviruses* 1988;4:239-41.
53. Simon F, Cot M, Lesager C, et al. Differentiation between HIV-1 and HIV-2 infection by radioimmunoprecipitation and synthetic peptides in double reactive sera. *AIDS* 1989;3:401-4.
54. Norby E, Biberfeld G, Chioldi F, et al. Discrimination between antibodies to HIV and to related retroviruses using site-directed serology. *Nature* 1987;329:248-50.
55. Gurtler LG, Zekeng L, Simon F, et al. Reactivity of five anti-HIV-1 subtype O specimens with six different anti-HIV screening ELISAs and three immunoblots. *J Virol Methods* 1995;51:177-84.
56. Loussert-Ajaka I, Ly TD, Chaix ML, et al. HIV-1/HIV-2 seronegativity in HIV-1 subtype O infected patients. *Lancet* 1994;343:1393-4.
57. Devare SQ, Desai SM, Dawson GJ, et al. Diagnosis and monitoring of HIV-1 and HIV-2 infection. In: Khan NC, Melnick JL, editors. *Human Immunodeficiency Virus: Innovative Techniques for Isolation and Identification (Monographs in Virology, vol.18)*. Basel: S Karger; 1990; 18:105-21.
58. Kessler HA, Blaauw B, Spear J, et al. Diagnosis of human immunodeficiency virus infection in seronegative homosexuals presenting with an acute viral syndrome. *JAMA* 1987;258:1196-9.
59. Phair JP. Human immunodeficiency virus antigenemia. *JAMA* 1987;258:1218.
60. Allain J-P, Laurian Y, Paul DA, et al. Serological markers in early stages of human immunodeficiency virus infection in haemophiliacs. *Lancet* 1986;2(8518):1233-6.
61. Kenny C, Parkin J, Underhill G, et al. HIV antigen testing. *Lancet* 1987;1(8532):565-6.
62. Wall RA, Denning DW, Amos A. HIV antigenemia in acute HIV infection. *Lancet* 1987;1(8532):566.
63. Stute R. HIV antigen detection in routine blood donor screening. *Lancet* 1987;1(8532):568.
64. Goudsmit J, De Wolf F, Paul DA, et al. Expression of human immunodeficiency virus antigen (HIV-Ag) in serum and cerebrospinal fluid during acute and chronic infection. *Lancet* 1986;2(8500):177-80.
65. von Sydow M, Gaines H, Sonnerberg A, et al. Antigen detection in primary HIV infection. *Brit Med J* 1988;296:238-40.
66. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
67. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
68. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
69. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
70. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining & Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PTs), November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> Accessed July 01, 2013.
71. Primus FJ, Kolley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibody therapy. *Clin Chem* 1988;34:261-4.
72. Schrott RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-85.
73. Boscatto LM and Stuart MC. Heterophilic antibodies; a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.

74. Box GEP, Hunter WG, Hunter JS. *Statistics for experimenters; an introduction to design, data analysis, and model building*. New York, NY: John Wiley & Sons, Inc, 1978:510-39, 571-83.
75. SAS Institute, Inc. SAS Technical Report P-229, *SAS/STAT Software: Changes and enhancements, Release 6.07*. Cary, NC: SAS Institute, 1992:289-366.

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HIV Ag/Ab Combo
Calibrator

INTENDED USE

The ARCHITECT HIV Ag/Ab Combo Calibrator is for the calibration of the ARCHITECT / System when used for the simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum or plasma. Refer to the ARCHITECT HIV Ag/Ab Combo reagent package insert for additional information.

CONTENTS

1 Bottle (4 mL) of ARCHITECT HIV Ag/Ab Combo Calibrator 1 (CAL1): Purified HIV viral lysate prepared in TRIS buffered saline with protein (bovine) stabilizer. Preservative: Sodium Azide.

Calibrator	Color
CAL1	Red ^a

^a Dye: Red D&C No. 33

STANDARDIZATION

To maintain a consistent HIV-1 p24 antigen sensitivity of < 50 pg/mL, the ARCHITECT HIV Ag/Ab Combo Calibrator has a HIV-1 p24 antigen concentration referenced to the international standard Agence française de sécurité sanitaire des produits de santé (AFSSAPS) panel.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents
- **CONTAINS: AZIDE** This product contains sodium azide. Contact with acids liberates very toxic gas.
- This material and its container must be disposed of in a safe way.
- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

STORAGE

- ARCHITECT HIV Ag/Ab Combo Calibrator is stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

PREPARATION FOR ANALYSIS

- Calibrator may be used immediately after removal from 2-8°C storage. Prior to use, mix by gentle inversion (5-10 times). After each use, lightly close the caps and return the calibrator to 2-8°C storage.

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BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization, *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute, *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety In Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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HIV Ag/Ab Combo Controls

INTENDED USE

The ARCHITECT HIV Ag/Ab Combo Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT i System when used for the simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum or plasma. Refer to the ARCHITECT HIV Ag/Ab Combo reagent package insert for additional information.

CONTENTS

4 Bottles (8 mL each) of ARCHITECT HIV Ag/Ab Combo Controls: Negative Control ([CONTROL]-), Positive Control 1 ([CONTROL+1]), and Positive Control 2 ([CONTROL+2]) prepared in recalcified human plasma. The Positive Control 1 (inactivated) is reactive for anti-HIV-1. The Positive Control 2 (inactivated) is reactive for anti-HIV-2. The Positive Control 3 ([CONTROL+3]) is purified HIV viral lysate prepared in TRIS buffered saline with protein (bovine) stabilizer, Preservatives for Negative Control, Positive Control 1 and Positive Control 2: sodium azide and antimicrobial agent. Preservative for Positive Control 3: sodium azide.

The controls are at the following S/CO ranges:

Control	Color	Control Range (S/CO)	RANGE
[CONTROL-]	Natural	0.00 - 0.50	
[CONTROL+1]	Blue ^a	1.20 - 11.50	
[CONTROL+2]	Yellow ^b	1.52 - 8.30	
[CONTROL+3]	Purple ^c	1.87 - 4.59	

^a Dye: Acid Blue No. 9

^b Dye: Acid Yellow No. 23

^c Dye: Acid Blue No. 9 and Red D&C No. 33

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

- **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.
- The human plasma used in the Positive Control 1 is reactive for anti-HIV-1 and nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HCV.
- The human plasma used in the Positive Control 2 is reactive for anti-HIV-2 and nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HCV.
- The human plasma used in the Negative Control is nonreactive for HBsAg, HIV-1 Ag or HIV-1 RNA, anti-HCV, and anti-HIV-1/HIV-2.
- **CONTAINS: AZIDE** This product contains sodium azide. Contact with acids liberates very toxic gas.
- This material and its container must be disposed of in a safe way.
- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

STORAGE

- ARCHITECT HIV Ag/Ab Combo Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

2°C/4°C

PREPARATION FOR ANALYSIS

- Controls may be used immediately after removal from 2-8°C storage. Prior to use, mix by gentle inversion (5-10 times). After each use, tightly close the caps and return the controls to 2-8°C storage.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

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