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
DECLARATION

I, the representative director Song Hee, Kim of Kims Med Co., Ltd
do hereby solemnly and declare:

1. That the attached document(s):

EZ Fan, Instruction for use

2. ☒ is (are) true and correct.

☐ is (are) true translation for the text(s) originally written in
the Korean language conscientiously believing the same to be
true and correct.

13 . May. 2021

Kims Med Co., Ltd

Representative director Song Hee, Kim

481-87-01096

KIMSMED Songhee Kim

#112&214, Bio-medical Material & Component Service
Center, Gwangju TechnoPark, 249, Chuam-ro, Buk-gu,
Gwangju, 61003, Republic of Korea
Tel +82-70-5208-5257 Fax +82-70-8611-3022

EZ Fan

Instruction for use

Model Name: REF 01KR

481-87-01096

KIMSMED Songhee Kim

#112&214, Bio-medical Material & Component Service
Center, Gwangju TechnoPark, 249, Chuam-ro, Buk-gu,
Gwangju, 61003, Republic of Korea
Tel +82-70-5208-5257 Fax +82-70-8611-3022

KIMS MED

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Description of EZ Fan Funnel

EZ Fan

The EZ Fan is a sterile funnel for the delivery of silicone gel implants during breast reconstruction surgeries.

It is a transparent, cone-shaped sleeve constructed of a flexible, polymeric film with a hydrophilic coating on the inside surface.

Note: HA is not absorbed into the human body during the surgical operation procedure.

Principle of operation

It delivers a sterile breast implant, minimizing the contact of the breast implant with surrounding tissue and making insertion easier.


Advantages of using EZ Fan funnel in comparison to the traditional manual insertion technique, without using any additional instruments, include the following: contactless delivery technique, possibility to visualize implant position during delivery, shortening duration of the surgery, reduction of surgeon fatigue, limitation of local pressure to the implant shell, less incision size and tissue injury degree during implant delivery.


Risk class – class 2a

Dimensional characteristics of EZ Fan funnel

Dimension (Unit: mm)

Tolerances : +/-1mm

EzFan funnel weight - 26,3 g $\pm$ 10%.

Technical characteristics of device raw materials

Raw materials	Specification	
Polyvinyl Chloride Film	Physical State	Solid
	Colour	clear
	Odor	none
	Density	1.30 ~ 1.40 gr/cm ³ (ASTM D792-2013)
	Heat deflection	62-65°C (ASTM D648-2016)
	Solubility in Water	$\leq 0.1\text{g}/1000\text{ml}$ at 23°C (ASTM E1148-2022)
Hydrophilic Coating Solution Top Coat*	Tensile strength	>16 MPa (ASTM D882 Method A)
	pH	7.00 ~ 9.00
Hydrophilic Coating Solution Base Coat*	Viscosity	500 ~ 900 cps(25°C)
	pH	6.00 ~ 9.00
	Viscosity	5 $\pm$ 2 cps(25°C)

*pH and Viscosity will be changed after the coating solution, this characteristics are relevant only for the solution.

Sterility of Product

Each EZ Fan funnel is sterilized by ethylene oxide in a sealed, double pouched configuration. Residual amount of EO is less than 4mg/unit. Sterility of this product can only be maintained if the package seals are intact.

IFU-01 (Date.2019-12-19)

Intended Use of EZ Fan funnel

The sterile funnel EZ Fan for delivery of silicone gel implants during reconstruction surgeries is intended to assist with the delivery of silicone gel Breast implants by providing a shell-tissue interface with less friction during insertion of the implant.

Intended user - Plastic surgeons qualified in breast augmentation

Indications

Installation of silicone gel breast implants during breast reconstruction.

Contraindication

Pregnant women

Side effects:

The risk of Cross infection when the instruction for use is violated.

1) Opening the package and extraction

1 Remove the EZ Fan funnel from its outer package in an **aseptic environment** using talc-free gloved hands. Peel open outer package.

2 Invert package over sterile field, allowing inner pouch to gently fall into sterile field, or by presenting the inner peel pouch to the scrubbed person using sterile technique.

3 Using sterile gloves, open inner pouch and remove EZ Fan funnel placing it into sterile field.

Note: the manufacturer KimsMed recommends that the funnel be prepared immediately prior to use. Hydration of the funnel should take place after the surgical pockets have been fully prepared. Prolonged soaking of the funnel will diminish the lubricity of the hydrophilic coating and affect the ease of placing the implant.

2) Trimming the funnel :

It is recommended to use sterile straight blade scissors to trim the distal end of the funnel for the implants with volumes greater than 315cc./4.5cm², to accommodate the size, texture, type and shape of the silicon gel implant being used to achieve your desired outcome, following the marks on the proximal tip of the funnel.

The sizing guide on the funnel refers to use with smooth round implants.

Information to trim the EZ Fan funnel

Smooth Round	Textured Round
Moderate & Moderate Plus: as shown	Moderate & Moderate Plus: +1.5 cm larger
High Profile: +0.5cm larger	High Profile: +2cm larger
Ultra High Profile: +1cm larger	Ultra High Profile: +2cm larger
Anatomical / Shaped Implant +2cm larger	
Form Stable Implants – 1 to 2 cm larger than indicated	

Note: The length of the distal end of the funnel does not dictate the incision length required. Follow the implant manufacturers' guidelines for determining and making the optimum incision based on the specific implant being used

3) HYDRATE:

Pour any sterile solution into EZ Fan funnel and close the upper open entrance of EZ Fan funnel with a hand. Then move the funnel back & forth for the hydration

019-12-19)

with holding the bottom and the upper of EZ Fan funnel, Pour out the saline out of EZ Fan funnel.

4) Confirm Sizing:


Insert the implant into the large proximal end of the hydrated funnel while holding over the basin of sterile saline/ sterile field. Use one hand to support the distal end of the funnel, and the other hand to manipulate behind the implant. Use this hand movement to push the implant forward through the tapered funnel and ejecting the implant into the basin of sterile saline/ on a sterile field. If the implant does not exit the small end of the funnel using gentle force, then remove the implant, re-trim, and repeat the test, gradually trimming larger opening sizes until the desired result is achieved.

5) Ensure Funnel Remains Moist:

Ensure the inner surface of the funnel remains moist. Load implant back into proximal end of funnel. Using one hand, apply slight pressure to the funnel to remove any slack, verify positioning of the implant orientation within the funnel and advance the implant to approximately 1cm from the distal end.

Note: The orientation of the implant relative to its resting position in the tissue pocket can be visually confirmed through the translucent fabric and manipulated as required prior to insertion into the surgical pocket.

When not in use, the Ez Fan should be kept moist to prevent the coating on the inner surface from drying completely.

6) Loading the funnel

The silicon gel implant is loaded into the trimmed and hydrated funnel. Round implants are placed directly into the funnel from the packaging, for what is called a "No-Bowl and No-Touch Technique". If you are receiving anatomical or shaved implants, the surgeon will load them into the funnel by hand, to ensure proper direction and orientation of the implant prior to placing it in the body. This is called Minimal-Touch technique.

7) PROPEL THE FUNNEL DELIVERY:

With the incision opening retracted, continue applying slight pressure on the implant within the funnel to keep the implant approximately 1 cm from the distal end. Insert the distal end of the funnel approximately 1 cm through the incision. Use one hand to support the funnel and use the opposite hand to employ a moderately forceful pushing/squeezing technique behind the implant to propel the implant forward into the center of the surgical pocket.

Note: Lubricity of the hydrophilic coating is diminished when allowed to dry and rewetted.

Caution!

- Check the expiry date before use. Do not use the expired devices.
- Check the integrity of the package before use. Do not use the product if the package has been opened or damaged.
- There are various manufacturers and types of implants, each with multiple variations of geometries, volumes, sizes, shapes, surface textures, cohesiveness, etc., as well as routine introductions of new models. The trim guidelines are only suggested trim locations. It is the responsibility of the user to assure the length of the distal funnel opening is adequate to insure the implant can pass through without any damage,

Breast Implants: The following trim guidelines are printed on the funnel, to accommodate

- various shapes and styles of implants based on their volume (cc's).
- Before delivering the implant into the surgical pocket, it is essential to confirm the distal end of the funnel is large enough for the implant to pass through. If the distal end is too small at this confirmation step, the funnel is designed to tear rather than damage the implant. However, when the end of the funnel is in the tissue pocket, the funnel may not tear prior to causing damage to the implant.
 - Using excessive force may damage the implant. *If the implant does not gently advance through the EZ Fan into the surgical pocket, similar to the amount of force required during the sizing confirmation step, STOP! Ensure the distal end is not folded, pinched, wrinkled or overly constricted. Verify that both the incision and **surgical pocket** are large enough to accommodate the implant. Verify the end of the funnel is not inserted more than 1 cm, is directed toward the centre of the surgical pocket, and the pathway is not obstructed by muscle tissue, or surgical instruments.*

Warnings:

- The device is single use and one patient only.
- The device is sterile, do not re-sterilize.
- Reuse or re-sterilization should not be performed under any conditions as contamination or infection may occur. Attempts to re-use or re-sterilize may cause diminished product performance, including loss of lubricity, potentially causing implant damage, including rupture.
- **DO NOT** expose the EZ Fan to lint, talc, sponges, towels, or other surface contaminants.

Materials of EZ Fan

No	Part	Raw Materials		Standards (CAS No.)	Composition (Wt%)
1	Body	Polyvinyl Chloride Films (formula: $(CH_3CH_2Cl)_n$)	Polyvinyl Chloride	93050-82-9	65%
			Dioctyl terephthalate	6422-86-2	35%
2	Coating Solution	Hydrophilic Coating Solution Top-coat	Ethanol 95%	64-17-5	≤50%
			D.I. water	7732-18-5	≤50%
			Hyaluronic Acid	9004-61-9	≤10%
		Hydrophilic Coating Solution Base-coat	Toluene	108-88-3	≤80%
			PBIN(Binder)	68891-50-9	≤10%
			1-Methoxy-2-propanol acetate	108-65-6	≤10%
3	Tyvek (1073B) Packages	PE/PET Film (35754-E, PerfecFlex)	PE (Polyethylene)		98%
			PET (Polyester)		2%

			Polyurethane adhesive		
		Tyvek film	Polyethylene	9002-88-4	100%
	3 rd package	Paper			
	4 th package	Double wall corrugated carton			

No materials of animal and (or) human origin were used for production of this medical device.

Package of EZ Fan

Type	Raw material	Size, mm (+/-1mm)	Q'ty (ea)	Remark
1 st package	PE/PET Film + Tyvek	200*125	1	Double packaging
2 nd package	PE/PET Film + Tyvek	275*200		
3 rd package	Paper	202 x 277 x 48	5	
4 th package	Double wall corrugated	515 * 420 * 300	100	

Sealing strength of the 1st and 2nd packages is not less than 3N/15mm

Storage and transportation conditions: Store at 1°C to 25°C, in a dry place. Keep away from sunlight.

Use at room temperature 15°C – 25°C

Shelf life: 2 Years

IFU-01 (Date.2019-12-19)

Disposal

EZ Fan relates to class B (epidemiologically hazardous waste) based on hazard classification, and is potentially infected waste. Used medical device should be disposed according to the medical waste hazard class and in accordance with local regulations. Medical devices with expired shelf life or in a damaged package should be disposed of as medical waste of Class A and in accordance with local regulations.

TROUBLE SHOOTING GUIDELINES FOR THE EZ Fan

ISSUE: Funnel is sized properly, but the implant does not move easily through the funnel.

Question: Is the inner surface of the EZ Fan slippery or lubricious?

Considerations:

- 1 Ensure the inner surface is hydrated.
- 2 Deliver the implants into the surgical pocket as soon as possible after hydration.
3. Prolonged soaking of the funnel will diminish the lubricity of the hydrophilic coating and affect the ease of placing the implant.

ISSUE: The funnel tore during the sizing step or during implant delivery.

Question: Was the end trimmed appropriately for the implant volume?

Considerations:

- 1 Open and trim a new funnel according to instructions.
- 2 The sizing template provides recommendations for trimming the funnel when delivering most smooth round implants.
- 3 Sizing the funnel for textured, shaped or high-profile implants may require larger trim size.
- 4 Confirm funnel is trimmed correctly and the implant moves easily through the funnel before inserting the funnel into the incision

Issue: Implant does not move through the funnel into the surgical pocket with minimal to moderate force similar to the amount of force used during the sizing confirmation step.

Questions: Is the end of the funnel open so the implant can move through funnel into the surgical pocket? Is the funnel oriented to allow the implant to pass into the center of the surgical pocket?

Considerations:

1. Ensure that the distal end of the funnel is not folded, pinched, wrinkled, or overly constricted in any way.
2. Verify that both the **incision** and **surgical pocket** are large enough to accommodate the implant.
3. Be sure the end of the funnel is not inserted too deep through the incision and into the pocket.
4. The funnel should be aligned to allow the implant to be directed to the center of the surgical pocket.
5. Be sure the pathway is not obstructed by muscle, tissue, or surgical instruments.

ISSUE: Surgeon delivers the implant on one side of the patient easier than the other.

Question: Is the surgeon right or left hand dominant?

Considerations:

- 1 For a Transaxillary Approach
 - a Right Hand Dominant surgeons can stand above the arm-board when delivering the right implant and below the arm-board when delivering the left implant.
 - b Left Hand Dominant surgeons can stand below the arm-board when delivering the right implant and above the arm-board when delivering the left implant.

2 For Periareolar, Inframammary

- a Right Hand Dominant surgeons can stand on patient's right side when delivering both left and right implants.
- b Left Hand Dominant surgeons can stand on patient's left side when delivering both left and right implants.

Labeling

On the funnel there is the following marking:

On the 2nd package (Back side), paper box and corrugated box there is the following label:

i) Main Label

Propriety name	EZ Fan (Sterile single use instrument for the delivery of silicone gel implants)	CE 2195	CE Mark
Model Name	REF REF_01KR		
Lot number	LOT XXXXXXXXX		
Date of Manufacture	YYMM YYYY-MM	Use By YYYY-MM	
Method of Sterilization	STERILE EO	Consult instructions for use	
Manufacturer	KEEP AWAY FROM SUNLIGHT Kims Med Co., Ltd. #112 #214 Bio Material and Component Service Center, Gwangju Techno Park 249 Chuam-ro Buk-gu Gwangju Korea Tel +82-70-5208-5257 Fax +82-70-8611-3022		
EU representative	EC REP Certification Experts B.V. Amerlandseweg 7, 3621 ZC Breukelen, The Netherlands Cautions Sterile unless package is open or damaged		

Marking of the third package also includes the number of the funnels inside – 5ea/package.

Labels on the 2nd package (Back side), paper box and corrugated box contain the following symbols

Symbols

	Date of manufacture		Manufacturer		Keep Away From Sunlight
	Lot number		EU representative		Temperature Limitation
	Consult instruction for use		Do not reuse		Cautions
	Article Name		Do not Resterilize		Use By
	Ethylene Oxide Sterile		CE Marking		Do not use if package is damaged

The list of applicable standards

No	Standard	Application
1	EN ISO 13485:2016	Medical devices - Quality management systems - Requirements regulatory purposes(ISO 13485:2016)

2	EN ISO 14971:2012	Medical devices - Application of risk management to medical devices(ISO14971:2012)
3	EN ISO 11135:2014	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
4	EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices
5	EN ISO 14644-1:2015	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
6	EN 62366:2008	Medical devices-part1: Application of usability engineering to medical devices
7	ASTM F1980 - 16	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
8	EN ISO 15223-1:2012	Medical devices: Symbols to be used with medical devices labels, labeling and information to be supplied
9	EN1041:2008+A1:2013	Information supplied by the manufacturer with medical devices
10	ISO 10993-1:2009	Biological evaluation of medical devices-Part1: Evaluation and testing within a risk management process
11	ISO 10993-13:2010	Biological evaluation of medical devices. Identification and quantification of degradation products from polymeric medical devices
12	ISO 17050-1:2010	Conformity assessment-supplier's declaration of conformity - part 1:General requirements

Maintenance and repair of a medical device


The medical device is not subject to maintenance and repair. It is single use only.

Warranty

The warranty ensures that the devices are free from defects during the stated shelf life.

Kims Med Co., Ltd.

#112, #214 Bio Material and Component
Service Center, Gwangju Techno Park,
249, Chuam-ro, Buk-gu, Gwangju,
Korea
E-mail : songkim@kimsmed.com www.
kimsmed.com Tel: +82-62-571-4645,
Fax: +82-62-571-4644

Certification Experts B.V.

Amerlandseweg 7, 3621 ZC Breukelen, The Netherlands

IFU-01 (Date.2019-12-19)

Claims Department/ Official Representative of Manufacturer in Russia

LLC CLOVERMED - Russia, 121099, Moscow, Noviy Arbat 34, bld 1, placement I, floor 5, room
1, Tel. +7 499 918 17 72; www.clovermed.ru

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attorney-in-fact of

KIMSMED CO., LTD. _____

Representative Director Song Hee, Kim

appeared before me and admitted
said principal's subscription to
the attached _____
DECLARATION _____

This is hereby attested on this
20th day of May 2021 at this office.

SHINHAN LAW & NOTARY OFFICE

Seoul Central

District Prosecutor's Office

#403-1, 19, Jong-ro, Jongno-gu,

Seoul, Korea

Yoon S. J.

Signature of the Notary Public
YOON SUK JUNG

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
3, Jan. 2005 Under Law No.325.

APOSTILLE

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1. Country : Republic of Korea

This public document

2. has been signed by YOON SUK JUNG

3. acting in the capacity of Notary Public

4. bears the seal/stamp of SHINHAN LAW AND NOTARY OFFICE

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