Apostille

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1. Country: United States of America

This public document:

2. has been signed by Dawn Renea Kenton

3. acting in the capacity of Notary Public

4. bears the seal/stamp of Dawn Renea Kenton, Notary Public,
Delaware

Certified


5. at Dover, Delaware

6. seventeenth day of August, A.D. 2020

7. by Secretary of State, Delaware Department of State

8. No.203483837

9. Seal/Stamp:

10. Signature:

Jeffrey W. Bullock, Secretary of State

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М.П. / Stamp

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Sworn to and subscribed before me this 13th
day of August, 2020.

SK

Dawn Renea Kenton, Notary Public

Instruction for Use

Revision 1.0

SDR[®] Plus Bulk Fill Flowable material in refills, with accessories

Contents

1. Name of the medical device	3
2. Classification of the medical device and accessories	3
3. Information regarding the developer of the medical device.....	3
4. Information regarding the manufacturer of the medical device.....	4
5. Manufacturing site	4
6. Intended use	4
7. Area of intended use.....	4
8. Indications, contraindications, adverse effects, precautions	4
9. Potential users.....	5
10. Technical description and characteristics	11
11. List and description of the materials brought in contact with the human body.....	16
12. Delivery forms.....	19
13. Labeling of the dental material	20
14. Packaging of the medical device.....	23
15. Information regarding materials of animal or human origin	32
16. Information regarding compatibility with other devices, materials, and etc.....	32
18. Information regarding drug products used in the medical device.....	32
19. Information regarding disinfection, cleaning and sterilization of the medical device	32
20. Requirements for maintenance and repair of the medical device.....	32
21. List of international and national regulations/standards applicable to the medical device	33
22. Conditions of use, storage and transportation.....	34
23. Shelf life.....	35
24. Procedure and conditions for disposal of the medical device.....	36
25. Warranty.....	36
26. Official representative.....	36

1. Name of the medical device

SDR® Plus Bulk Fill Flowable material in refills, with accessories:

1. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade, A1, A2, A3 – 15 pcs./pack;
2. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade, A1, A2, A3 – 50 pcs./pack;
3. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 4 pcs./pack, IFU – 1 pc., illustrated technical manual – 1 pc.;
4. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 10 pcs./pack, IFU – 1 pc., illustrated technical manual – 1 pc.;
5. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 2 pcs., shade A1 – 1 pc., shade A2 – 1 pc., shade A3 – 1 pc., IFU – 1 pc., illustrated technical manual – 1 pc.;
6. SDR® Plus Bulk Fill Flowable in 1 g syringes each – 2 pcs. (universal shade), syringe tips – 15 pcs., IFU – 1 pc., illustrated technical manual – 1 pc.
7. SDR® Plus Bulk Fill Flowable in 1 g syringes – 10 pcs (universal shade), syringe tips – 60 pcs.; IFU – 1 pc., illustrated technical manual – 1 pc.
8. SDR® Plus Bulk Fill Flowable in 1 g syringe – 1 pc. (universal shade), syringe tips – 3 pcs., IFU – 1 pc., illustrated technical manual – 1 pc., ID card – 1 pc., repeat order card – 1 pc.

Accessories:

- 1) Syringe tips – 60 pcs./pack.
- 2) SDR® storage box – 1 pc.

2. Classification of the medical device and accessories

All-Russian classifier of products by type of economic activity, OKPD 2: 32.50.11.000

Type of the medical device according to the Nomenclature Classification of Medical Devices:
281770

Class of potential risk related to the usage of the medical device in accordance with the Nomenclature Classification of Medical Devices: 2a.

3. Information regarding the developer of the medical device

Dentsply LLC also dba Dentsply Caulk

Address: 38 West Clarke Avenue, Milford, Delaware, 19963, USA

info@dentsplysirona.com

4. Information regarding the manufacturer of the medical device

Dentsply LLC also dba Dentsply Caulk

Address: 38 West Clarke Avenue, Milford, Delaware, 19963, USA

info@dentsplysirona.com

5. Manufacturing site

1. Dentsply LLC also dba Dentsply Caulk, 38 West Clarke Avenue, Milford, Delaware, 19963, USA
2. Dentsply De Trey GmbH, De-Trey-Street 1, 78467 Konstanz, Germany

6. Intended use

SDR® Plus Bulk Fill Flowable material in refills is intended for use as a base liner. It is designed to be used as a base in Class I and II restorations. It is also suitable as a stand-alone restorative material in non-occlusal-contact applications.

7. Area of intended use

Area of intended use – dentistry.

The medical devices are intended for use by dentists in dental clinics and dental laboratories

8. Indications, contraindications, adverse effects, precautions

8.1. Indications

SDR® Plus material direct restorative:

- Base in cavity Class I & II direct restorations
- Liner under direct restorative materials – Class II box liner
- Pit & Fissure Sealant
- Conservative Class I restorations
- Core build-up
- Class III and V restorations

8.2. Contraindications

SDR® Plus material is contraindicated:

- for use with patients who have a known hypersensitivity to methacrylate resins.

8.3. Adverse reactions:

The product may cause eye and skin irritation.

- Eye contact: irritation, possible corneal damage.
- Skin contact: irritation, possible allergic reaction. Red rash may appear on the skin.
- Mucous membrane contact: inflammation.

8.4. Warnings

SDR® Plus material contains polymerizable methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

The medical device must be handled only in gloves.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.

- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and alcohol and wash well thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.

- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

9. Potential users

The medical devices are intended for use by dentists in dental clinics and dental laboratories.

1. Description

SDR® Plus Bulk Fill Flowable is a one-component, fluoride-containing, visible light cured, radiopaque resin composite restorative material. It is designed to be used as a base in Class I and II restorations. It is also suitable as a stand-alone restorative material in non-occlusal-contact applications.

SDR® Plus material has handling characteristics typical of a “flowable” composite, but can be placed in 4 mm increments with minimal polymerization stress. SDR® Plus material has a self-leveling feature that allows intimate adaptation to the prepared cavity walls. When used as a base/liner, it is designed to be overlaid with a methacrylate based universal/posterior composite for replacing missing occlusal/facial enamel.

2. Step-by-step instructions

Cavity preparation

1. Prepare the cavity so that no residual amalgam or restorative material is left.

- Rinse surface with water spray and carefully dry it with air spray. Do not desiccate the tooth structure.
- Use a dental dam or cotton rolls to isolate the cavity from contamination.

Placement of matrix

For optimal proximal contacts proceed as follows:

1. Place a matrix (e.g. AutoMatrix® matrix system or Palodent® Plus/Palodent® V3 Sectional matrix system) and wedge. Burnishing of the matrix band will improve contact and contour. Prewedging/ ring placement is advocated.

Pulp protection, tooth conditioning/dentin pretreatment, adhesive application

Refer to adhesive manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated.

Proceed immediately to placement of SDR® Plus material.

Danger of injury due to excessive force.

1. Apply slow and steady pressure on the syringe plunger Compules® Tips Gun.
2. Do not use excessive force – Syringe or Compula® Tip rupture.

Compula® Tips

The pre-dosed Compula® Tip provides the combination of a unit dose Compules® Tip with a metal cannula Applicator Tip.

1. Load Compules® Tips Gun with pre-dosed Compula® Tip. Insert Compula® Tip into the notched opening of the Compules® Tips Gun barrel. Be certain that the collar on the Compula® Tip is inserted first.
2. Remove the colored cap from the Compula® Tip. The Compula® Tip may be rotated 360° to gain the proper angle of entrance into the cavity.
3. Dispense the material into the cavity preparation using a slow, steady pressure. **DO NOT USE EXCESSIVE FORCE.**

To remove the used Compula® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compula® Tip and remove.

Syringes

1. Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto pad, away from the patient field.
2. Attach disposable, black dispensing tip to end of the syringe. Turn tip clockwise ¼ to ½ turn to

assure that it is fully seated. Tug on tip to be sure that it is locked into the collar of the syringe.

3. Material should flow freely with gentle pressure. DO NOT USE EXCESSIVE FORCE. If more than gentle pressure is required, remove from patient field and check for obstruction.

4. DISCARD AND PROPERLY DISPOSE DISPENSING TIP IMMEDIATELY AFTER USE.

REPLACE

ORIGINAL CAP. Do not store syringe with black dispensing tip in place. STORE ONLY WITH ORIGINAL CAP.

Note: It is recommended to pull back slightly on the syringe plunger after use to prevent excessive flow of material.

Placement of SDR® Plus material

Note: As SDR® Plus material is designed to conform to cavity shape, proper contact areas must be established with the matrix system selected. Material will not distend matrix band. If needed, manually distend matrix band with a suitable hand instrument such as a condenser or ball burnisher prior to and/or during light curing.

1. Dispense SDR® Plus material directly into preparation site/tooth surface using slow, steady pressure. Begin dispensing at the deepest portion of the cavity, keeping tip close to cavity floor.

Gradually withdraw tip as cavity is filled. Avoid lifting the tip out of dispensed material while dispensing to minimize air entrapment. At the completion of dispensing, wipe tip against cavity wall while withdrawing from the operative field.

2. Within a few seconds, dispensed SDR® Plus material will self-level, eliminating the need for further manipulation with hand instruments. In case of overfill, use a flocked applicator tip wetted with residual adhesive to remove excess. Any visible air bubbles should be pierced with a clean, sharp explorer prior to curing.

3. When used as a Bulk Fill Base material, most cavities may be filled in one bulk increment (up to 4 mm) as needed to fill the cavity 2 mm short of the occlusal cavosurface¹. In deeper preparations, place material in 4 mm increments, thoroughly light curing each increment.

4. Alternatively, SDR® Plus material may be placed in a thin layer as a traditional “flowable” liner on exposed dentin.

Curing

SDR® Plus material is designed to be cured in increments up to a 4 mm depth/thickness.

¹ For fissure sealing, class I, III, V conservative restorations or stump build-up the material can be fully inserted inside the contour, in particular up to the edge of the hollow surface.

1. Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 550 mW/cm².

Pulp Injury due to excessive heat – curing output > 2,000 mW/cm².

1. Check Curing Light manufacturer's literature for stated output.
2. Consult Curing Light Manufacturer's Directions for compatibility curing recommendations.
3. Do not use table below for recommended curing times with output > 2,000 mW/cm².

Curing Recommendations ²

Shade	Light Output	Shade Light output	Shade Light output
Universal	Halogen and LED Lights 550-1,000 mW/cm ²	20 sec.	20 sec.
	High Power LED Lights 1,000-2,000 mW/cm ²	10 sec.	10 sec.
A1, A2, A3	Halogen and LED Lights 550-1,000 mW/cm ²	20 sec.	40 sec.
	High Power LED Lights 1,000-2,000 mW/cm ²	10 sec.	25 sec.

Inadequate polymerization due to insufficient curing.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output before each procedure.

Completion

When used as a base/liner, proceed immediately to complete the restoration using a methacrylate based universal/posterior restorative material following manufacturer's Directions for Use.

Technique Tip: In most applications, a body shade of occlusal restorative material provides the best esthetic result. In case of heavily stained dentin, use of an opaque shade may be required. Do not disturb or allow contamination of cured surface or exposed adhesive. If contamination occurs, follow adhesive manufacturer's Directions for Use for re-application of adhesive.

² For instructions on polymerization and compatibility information see the recommendations provided by the manufacturer of the light curing unit.

If surfaces remain undisturbed, additional application of adhesive between increments is neither necessary nor recommended.

Finishing and polishing

1. If contouring, finishing, and/or polishing of cured SDR® Plus material is necessary, use traditional rotary devices following manufacturer's Directions for Use.
2. Removal of resin flash is best accomplished with the Enhance® Finishing System. The Enhance® System will remove flash and finish restoration margins.
3. Polish final restoration using Enhance® PoGo® Polishing system or Prisma® Gloss™ Polishing Paste. (See complete Directions for Use supplied with polishing product chosen).

All shades of SDR® Plus material are radio-opaque, with a 1 mm radio-opacity equivalent to 2.6 mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin.

3. Hygiene

	Cross-contamination. Infection. 1. Do not reuse disposable materials. Dispose the materials according to the local regulations. 2. Reusable devices should be disposed accordingly.
---	---

3.1. Syringes –cross-contamination

	Cross-contamination. Infection. 1. The syringes cannot be reprocessed. 2. To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves. Do not reuse syringes if contaminated. 3. Dispose of contaminated syringes in accordance with local regulations.
---	--

For syringe, remove used dispensing tip and discard appropriately. Replace original syringe cap prior to storing. To prevent syringes from exposure to spatter or spray of body fluids or contaminated

hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination. Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument. Repeated liquid contact may damage label. Dry the syringe with a lint-free single use cloth.*

*** Note: Destruction of the label by vigorous wiping.**

Wipe syringe gently.


3.2. Compula® Tips – cross-contamination

	Cross-contamination. Infection. 1. Do not reuse Compula® Tips. Dispose Compula® Tips in accordance with local regulations.
--	---

10. Technical description and characteristics

13.1. Technical description

Table 1. Description of the medical device.

Device, description	Picture
SDR® Plus Bulk Fill Flowable in compula tips 0.25 g, (SDR® Plus Compula® Tips), universal shade Needle OD, mm: $0,9 \pm 1\%$ Needle ID, mm: $0,6 \pm 1\%$ Flange diameter, mm: $8,0 \pm 1\%$ Compula tip length, mm: $31,5 \pm 1\%$ Compula plug OD, mm: $4,25 \pm 1\%$ Weight, g: $1,0 \pm 1\%$	
SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shade A1 Needle OD, mm: $0,9 \pm 1\%$ Needle ID, mm: $0,6 \pm 1\%$ Flange diameter, mm: $8,0 \pm 1\%$ Compula tip length, mm: $31,5 \pm 1\%$ Compula plug OD, mm: $4,25 \pm 1\%$ Weight, g: $1,0 \pm 1\%$	
SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shade A2 Needle OD, mm: $0,9 \pm 1\%$ Needle ID, mm: $0,6 \pm 1\%$ Flange diameter, mm: $8,0 \pm 1\%$ Compula tip length, mm: $31,5 \pm 1\%$ Compula plug OD, mm: $4,25 \pm 1\%$ Weight, g: $1,0 \pm 1\%$	

**SDR® Plus Bulk Fill Flowable in compula
tips 0.25 g (SDR® Plus Compula® Tips),
shade A3**

Needle OD, mm: $0,9 \pm 1\%$

Needle ID, mm: $0,6 \pm 1\%$

Flange diameter, mm: $8,0 \pm 1\%$

Compula tip length, mm: $31,5 \pm 1\%$

Compula plug OD, mm: $4,25 \pm 1\%$

Weight, g: $1,0 \pm 1\%$

**SDR® Plus in 1 g syringes each (universal
shade)**

Syringe length, mm: $74,65 \pm 1\%$

Diameter, mm: $5 \pm 1\%$

Weight, g: $6,5 \pm 1$

Syringe tips

Needle OD, mm: $1,25 \pm 1\%$

Needle ID, mm: $0,97 \pm 1\%$

Length, mm: $29 \pm 1\%$

Diameter, mm: $2,54 \pm 1\%$

Weight, g: $0,26 \pm 1\%$

SDR® storage boxLength, mm: $89,0 \pm 1\%$ Width, mm: $66,0 \pm 1\%$ Height, mm: $35,5 \pm 1\%$ Weight, g: $45,0 \pm 1\%$ **13.2. Technical characteristics**

Table 2. Technical characteristics: universal shade

Characteristics	Specification
Physical mechanical properties	
Bending strength, MPa	Not less than 100
Barcol hardness	Not less than 85
N value characterizing bending flexibility, MPa	Not more than the strength value given in the GOST 31574-2012
Diametric strength or strength upon diametrical, MPa	Not less than 34 according to the GOST 31574-2012
Water absorption, $\mu\text{g}/\text{mm}^3$	Less or equal to 40 according to the GOST R 56924-2016
Water solubility, $\mu\text{g}/\text{mm}^3$	Less or equal to 7.5 according to the GOST R 56924-2016
Sensitivity to ambient light, c	60 according to the GOST 31574-2012
Curing depth, mm	Less or equal to 4 mm
Radiopacity, %	Not less than that of a radio opaque aluminum plate of the same thickness
Neck area polymerization time, s	Between 20 and 30
Digital value for color difference, designated as ΔE , between irradiated sample of cured material or water saturated sample and reference sample	Not more than 2.5 according to the GOST R 56924-2016
Aesthetic properties	
Transparency factor, %	Not less than 37.5
Appearance	Smooth paste, almost completely without foreign substances and excessive porosity
Discoloration	Discoloration should be barely

	noticeable
Adhesion properties	
Adhesion strength with dentine, MPa	Not less than 7 according to the GOST R 56924-2016

Table 3: Technical characteristics: shade A1

Characteristics	Specification
Physical mechanical properties	
Bending strength, MPa	Not less than 100
Barcol hardness	Not less than 85
N value characterizing bending flexibility, MPa	Not more than the strength value given in the GOST 31574-2012
Diametric strength or strength upon diametrical, MPa	Not less than 34 according to the GOST 31574-2012
Water absorption, $\mu\text{g}/\text{mm}^3$	Less or equal to 40 according to the GOST R 56924-2016
Water solubility, $\mu\text{g}/\text{mm}^3$	Less or equal to 7.5 according to the GOST R 56924-2016
Sensitivity to ambient light, c	60 according to the GOST 31574-2012
Curing depth, mm	Less or equal to 4 mm
Radiopacity, %	Not less than that of a radio opaque aluminum plate of the same thickness
Neck area polymerization time, s	Between 20 and 30
Digital value for color difference, designated as ΔE , between irradiated sample of cured material or water saturated sample and reference sample	Not more than 2.0 according to the GOST R 56924-2016
Aesthetic properties	
Transparency factor, %	Not less than 37.5
Appearance	Smooth paste, almost completely without foreign substances and excessive porosity
Discoloration	Discoloration should be barely noticeable
Adhesion properties	
Adhesion strength with dentine, MPa	Not less than 7 according to the GOST R 56924-2016

Table 4: Technical characteristics: shade A2

Characteristics	Specification
-----------------	---------------

Physical mechanical properties	
Bending strength, MPa	Not less than 100
Barcol hardness	Not less than 85
N value characterizing bending flexibility, MPa	Not more than the strength value given in the GOST 31574-2012
Diametric strength or strength upon diametrical, MPa	Not less than 34 according to the GOST 31574-2012
Water absorption, $\mu\text{g}/\text{mm}^3$	Less or equal to 40 according to the GOST R 56924-2016
Water solubility, $\mu\text{g}/\text{mm}^3$	Less or equal to 7.5 according to the GOST R 56924-2016
Sensitivity to ambient light, c	60 according to the GOST 31574-2012
Curing depth, mm	Less or equal to 4 mm
Radiopacity, %	Not less than that of a radio opaque aluminum plate of the same thickness
Neck area polymerization time, s	From 20 to 30
Digital value for color difference, designated as ΔE , between irradiated sample of cured material or water saturated sample and reference sample	Not more than 2.0 according to the GOST R 56924-2016
Aesthetic properties	
Transparency factor, %	Not less than 37.5
Appearance	Smooth paste, almost completely without foreign substances and excessive porosity
Discoloration	Discoloration should be barely noticeable
Adhesion properties	
Adhesion strength with dentine, MPa	Not less than 7 according to the GOST R 56924-2016

Table 5: Technical characteristics: shade A3.

Characteristics	Specification
Physical mechanical properties	
Bending strength, MPa	Not less than 100
Barcol hardness	Not less than 85
N value characterizing bending flexibility, MPa	Not more than the strength value given in the GOST 31574-2012
Diametric strength or strength upon diametrical, MPa	Not less than 34 according to the GOST 31574-2012

Water absorption, $\mu\text{g}/\text{mm}^3$	Less or equal to 40 according to the GOST R 56924-2016
Water solubility, $\mu\text{g}/\text{mm}^3$	Less or equal to 7.5 according to the GOST R 56924-2016
Sensitivity to ambient light, c	60 according to the GOST 31574-2012
Curing depth, mm	Less or equal to 4 mm
Radiopacity, %	Not less than that of a radio opaque aluminum plate of the same thickness
Neck area polymerization time, s	From 20 to 30
Digital value for color difference, designated as ΔE , between irradiated sample of cured material or water saturated sample and reference sample	Not more than 2.0 according to the GOST R 56924-2016
Aesthetic properties	
Transparency factor, %	Not less than 37.5
Appearance	Smooth paste, almost completely without foreign substances and excessive porosity
Discoloration	Discoloration should be barely noticeable
Adhesion properties	
Adhesion strength with dentine, MPa	Not less than 7 according to the GOST R 56924-2016

11. List and description of the materials brought in contact with the human body

Table 6. List of materials used

No.	Component	Materials
1.	SDR [®] Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR [®] Plus Compula [®] Tips) shades: universal shade	SureFil SDR resin
		Ethoxylated bisphenol A dimethacrylate
		Triethylene glycol dimethacrylate
		Trimethyl propane trimethacrylate
		Lumilux Blau
		Camphorquinone
		Ethyl-4-(dimethylamino)benzoate
		Butylated hydroxytoluene
		Uvinul H36206
		Aerosil
		Glass
		Aerosol OX-50
		Ytterbium trifluoride coated with SiO ₂ , silanized
		CAB-O-SIL TS-720
		Titanium dioxide (3328 CTFA BC)

		Silanized glass QuixFil
		Dye, iron oxide, yellow
		Dye, iron oxide, red (7080)
	Cap for compula tips	High density polyethylene
		Dye, white
	Plug for compula tips	Polyethylene
2.		Dye, white
	Tube	Stainless steel 1.4301
	Compula with a steel tube	Polyamide (nylon)
		Dye, grey
	SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: A1	SureFil SDR resin
		Ethoxylated bisphenol A dimethacrylate
		Triethylene glycol dimethacrylate
		Trimethyl propane trimethacrylate
		Lumilux Blau
		Camphorquinone
		Ethyl-4-(dimethylamino)benzoate
		Butylated hydroxytoluene
		Uvinul H36206
		Aerosil
		Glass
		Aerosol OX-50
		Ytterbium trifluoride coated with SiO ₂ , silanized
		CAB-O-SIL TS-720
		Titanium dioxide (3328 CTFA BC)
		Silanized glass QuixFil
		Dye, iron oxide, yellow
		Dye, iron oxide, red (7080)
	Cap for compula tips	High density polyethylene
		Dye, pink
	Plug for compula tips	Polyethylene
		Dye, white
	Tube	Stainless steel 1.4301
	Compula with a steel tube	Polyamide (nylon)
		Dye, grey
3.	SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: A2	SureFil SDR resin
		Ethoxylated bisphenol A dimethacrylate
		Triethylene glycol dimethacrylate
		Trimethyl propane trimethacrylate
		Lumilux Blau
		Camphorquinone
		Ethyl-4-(dimethylamino)benzoate
		Butylated hydroxytoluene
		Uvinul H36206
		Aerosil
		Glass

		Aerosol OX-50
		Ytterbium trifluoride coated with SiO ₂ , silanized
		CAB-O-SIL TS-720
		Titanium dioxide (3328 CTFA BC)
		Silanized glass QuixFil
		Dye, iron oxide, yellow
		Dye, iron oxide, red (7080)
	Cap for compula tips	High density polyethylene
		Dye, black
	Plug for compula tips	Polyethylene
		Dye, white
	Tube	Stainless steel 1.4301
	Compula with a steel tube	Polyamide (nylon)
		Dye, grey
4.	SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: A3	SureFil SDR resin
		Ethoxylated bisphenol A dimethacrylate
		Triethylene glycol dimethacrylate
		Trimethyl propane trimethacrylate
		Lumilux Blau
		Camphorquinone
		Ethyl-4-(dimethylamino)benzoate
		Butylated hydroxytoluene
		Uvinul H36206
		Aerosil
		Glass
		Aerosol OX-50
		Ytterbium trifluoride coated with SiO ₂ , silanized
		CAB-O-SIL TS-720
		Titanium dioxide (3328 CTFA BC)
		Silanized glass QuixFil
		Dye, iron oxide, yellow
		Dye, iron oxide, red (7080)
	Cap for compula tips	High density polyethylene
		Dye, red
	Plug for compula tips	Polyethylene
		Dye, white
	Tube	Stainless steel 1.4301
	Compula with a steel tube	Polyamide (nylon)
		Dye, grey
5.	SDR® Plus Bulk Fill Flowable material in syringes 1 g each (universal shade)	SureFil SDR # 999307 resin
		Ethoxylated bisphenol A dimethacrylate
		Triethylene glycol dimethacrylate
		Trimethyl propane trimethacrylate
		Lumilux Blau
		Camphorquinone

		Ethyl-4-(dimethylamino)benzoate
		Butylated hydroxytoluene
		Uvinul H36206
		Aerosil
		Glass
		Aerosol OX-50
		Ytterbium trifluoride coated with SiO ₂ , silanized
		CAB-O-SIL TS-720
		Titanium dioxide (3328 CTFA BC)
		Silanized glass QuixFil
		Dye, iron oxide, yellow
		Dye, iron oxide, red (7080)
	Syringe	Polypropylene
	Finger flange	Dye, black
		Polyethylene
6.	Syringe plug	Dye, black
		High density polyethylene
	Tube	Dye, white
		Stainless steel 1.4301
	Shaft	Polypropylene
		Dye, black
7.	SDR® storage box	Stainless steel 1.4301
		Dye, black
		Polyethylene
		Dye, white
		Tin
		Varnish

12. Delivery forms

1. SDR® Plus Bulk Fill Flowable material in refills, with accessories:

1. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade, A1, A2, A3 – 15 pcs./pack;

2. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade, A1, A2, A3 – 50 pcs./pack;

3. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 4 pcs./pack, IFU – 1 pc., illustrated technical manual – 1 pc.;

4. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 10 pcs./pack, IFU – 1 pc., illustrated technical manual – 1 pc.;

5. SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 2 pcs., shade A1 – 1 pc., shade A2 – 1 pc., shade A3 – 1 pc., IFU – 1 pc., illustrated technical manual – 1 pc.;

6. SDR® Plus Bulk Fill Flowable in 1 g syringes – 2 pcs each. (universal shade), syringe tips – 15 pcs., IFU – 1 pc., illustrated technical manual – 1 pc.

7. SDR® Plus Bulk Fill Flowable in 1 g syringes – 10 pcs (universal shade), syringe tips – 60 pcs.; IFU – 1 pc., illustrated technical manual – 1 pc.

8. SDR® Plus Bulk Fill Flowable in 1 g syringe – 1 pc. (universal shade), syringe tips – 3 pcs., IFU – 1 pc., illustrated technical manual – 1 pc., ID card – 1 pc., repeat order card – 1 pc.

13. Labeling of the dental material

16.1. Syringe labeling.

The labels are applied on the syringe packaging.

Dimensions: 16.25x12.5 mm ±5% (Width x Length).

The following information is indicated on the label according to the GOST R ISO 15223-1-2014:

- Name of the manufacturer;
- Name of the medical device;
- CE marking;
- Weight of the medical device;
- Curing time;
- Shade indication;
- Reference number of the medical device;
- Comment "Do not varnish"

Fig 1. Label contents

Table 7. Symbol description

Symbol	Description
--------	-------------

	CE marking
NO VARNISH	"Do not varnish"

16.2. Labeling of the consumer packaging

The labels are applied on the consumer packages according to the GOST R ISO 15223-1-2014.

- Name and address of the manufacturer;
- Manufacturer's trademark;
- Name of the medical device;
- Date of manufacture;
- Reference number of the medical device;
- Lot number;
- shelf life;
- "See instructions for use" label;
- "Caution! See instructions for use" label;
- "Do not reuse" label;
- "Protect from sunlight" label;
- "Temperature range" label;
- "Protect from moisture" label;
- Green dot;
- "Caution! This device can be purchased or ordered only by a doctor";
- EC marking;
- "Official representative in EC" label;
- MD (medical device) label

SDR® PLUS

SDR® Plus Bulk Fill Flowable in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade, A1, A2, A3 – 15 pcs./pack
See the shade on the label

Certificate of Registration № RZN #####/##### dated ##.##.####

Only for use by dentists.

Information regarding safety, indications and contraindications is given in the instructions for use.

Information regarding the terms and conditions of storage is given on the package

and in the instructions for use.

Keep the package till the end of the product use.

Manufacturer:

Dentsply LLC also dba Dentsply
Caulk
Address: 38 West Clarke Avenue,
Milford, Delaware, 19963, USA

Official representative:

Dentsply Sirona LLC,
Prospect Andropova 18, build. 6,
offices 1 39 40, Moscow, Russia
115432,
Tel.: +7 (495) 725 10 87

Place for label number

Country of origin:

USA


Fig. 2. Example of consumer packaging label

The number of Registration Certificate will be indicated on the Russian-language label on the consumer packaging

16.3. Symbols present on the consumer packaging

Table 8. Symbol description

Symbol	Description
	Name and address of the manufacturer
	Date of manufacture
	Lot number
	Shelf life
	See Instructions for Use
	Caution! See instructions for use
	Do not reuse
	Protect from sunlight

	Temperature range
	Protect from moisture
	Green dot label
	"Caution! This device can be purchased or ordered only by a doctor";
	CE marking
	MD (Medical Device) label

14. Packaging of the medical device


Fig. 3. Consumer packaging for SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g,
shades: universal shade, A1, A2, A3 – 15 pcs.

Packaging dimensions: (61 × 45 × 34) mm ± 5 %


Fig. 4 Consumer packaging for SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g,
shades: universal shade, A1, A2, A3 – 50 pcs./pack;

Packaging dimensions: (61 × 45 × 34) mm ± 5 %

Compula tips are placed inside the packaging according to the kit specification

Fig. 5. Compula tips packaging


Fig. 6. Consumer packaging for SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), shades: universal shade 4 pcs./pack, instruction for use – 1 pc., illustrated technical manual – 1 pc.

Packaging dimensions: 155 x 85 x 11 mm $\pm 5\%$


Fig. 7. Consumer packaging for SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 10 pcs./pack, IFU – 1 pc., illustrated manual – 1 pc.

Packaging dimensions: 155 x 85 x 11 mm $\pm 5\%$


Fig. 8. Consumer packaging for SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g (SDR® Plus Compula® Tips), universal shade – 2 pcs., shade A1 – 1 pc., shade A2 – 1 pc., shade A3 – 1 pc., IFU – 1 pc., illustrated technical manual – 1 pc.

Packaging dimensions: 155 x 85 x 11 mm ±5%


Fig. 9. SDR® Plus Bulk Fill Flowable material in syringes 1 g each – 2 pcs. (universal shade), syringe tips – 15 pcs., IFU – 1 pc., illustrated technical manual – 1 pc.

Packaging dimensions: (165 × 89 × 26) mm ±5 %

Fig. 10. SDR® Plus Bulk Fill Flowable material in syringes 1 g each – 10 pcs (universal shade), syringe tips – 60 pcs.; IFU – 1 pc., illustrated technical manual – 1 pc.

Packaging dimensions: (156 × 103 × 49) mm ±5 %

Syringes and tips are placed inside the packaging according to the kit specification

Fig. 11. Syringe and syringe tip packaging.

Fig. 12. SDR® Plus Bulk Fill Flowable material in syringes 1 g – 1 pc. (universal shade), syringe tips – 3 pcs., IFU – 1 pc., illustrative technical manual– 1 pc., ID card – 1 pc., repeat order card – 1 pc.
Packaging dimensions with the plastic bag: (200,0 × 125,4 mm) ±5 %


Fig. 13. Syringe tips – 60 pcs./pack.

Packaging dimensions: (90 × 43 × 43) mm ±5 %

Syringe tips are placed inside the polyethylene bag in quantity required by the kit specification

Fig. 14. Syringe tips in a polyethylene bag

15. Information regarding materials of animal or human origin

The dental material does not contain any substances of animal or human origin.

16. Information regarding compatibility with other devices, materials, and etc.

17. Compatible adhesives

The dental material is used after application of appropriate dentin-enamel adhesive. It is compatible with traditional dentin-enamel adhesives based on methacrylate, including adhesives for Dentsply Sirona light-cured composites.

Compatible restorative materials for enamel restoration

The dental material is used with appropriate restorative material or restorative material for filling third group of posterior teeth to replace occlusal/vestibular enamel and is compatible with traditional composite restoration materials based on methacrylate, including light-cured universal/dental filling composites by Dental Sirona intended for the manufacture of occlusal restorations of class I and II.

18. Information regarding drug products used in the medical device

The dental material does not contain any drug products.

19. Information regarding disinfection, cleaning and sterilization of the medical device

The medical device is not intended to be cleaned or disinfected. The medical device is not sterile and is not intended to be sterilized.

20. Requirements for maintenance and repair of the medical device

The medical device is not intended to be repaired or maintained.

21. List of international and national regulations/standards applicable to the medical device

ISO 14971 [2012] Medical devices – Application of risk management to medical devices;
EN 1641 [2009] Dentistry – Medical devices for dentistry - Materials;
ISO EN 13485 [2016] Medical devices – Quality management systems – Requirements for regulatory purposes;
EN ISO 7405:2011 Evaluation of biocompatibility of medical devices used in dentistry;
ISO 7491:2000(E) Dental materials – Determination of colour stability
ISO 4049:2000(E)- Dentistry – Polymer-based filling, restorative and luting materials;
ISO 10993-1:2003 Biological evaluation of medical devices. Part 1. Testing for skin irritation and sensitization;
ISO 10993-11 [2006] Biological evaluation of medical devices - Part 11 Tests for systemic toxicity;
ISO 11405: 2015 Dental materials – Testing of adhesion to tooth structure;
EN ISO 15223-1:2012 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements;
GOST R 56924-2016 (ISO 4049:2009) Dentistry. Polymer-based restorative materials;
GOST 31574-2012 Polymer-based dental restorative materials. Technical specifications. Methods of testing;
GOST ISO 7491-2012. Dental materials. Determination of colour stability of dental polymeric materials.
GOST ISO 14971-2011 Medical devices. Application of risk management to medical devices;
GOST ISO 13485-2017 Medical devices. Quality management systems. Requirements for regulation purposes;
GOST R ISO 7405-2011. Dentistry. Evaluation of biocompatibility of medical devices used in dentistry;
GOST ISO 10993-1-2011 Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and studies;
GOST R ISO 10993-2-2009 Medical devices. Biological evaluation of medical devices. Part 2. Animal handling requirements
GOST ISO 10993-5-2011 Medical devices. Biological evaluation of medical devices. Part 5. Cytotoxicity studies: in vitro methods
GOST ISO 10993-10-2011 Medical devices. Biological evaluation of medical devices. Part 10. Irritation and sensitization studies
GOST ISO 10993-11-2011 Medical devices. Biological evaluation of medical devices. Part 11. General toxicity studies;
GOST ISO 10993-12-2015 Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference samples

GOST R ISO 15223-1-2014 Medical devices. Symbols to be used with medical device labels, labeling, and information to be supplied. Part 1. General requirements;

GOST R 52770-2016 "Medical devices. Safety requirements. Methods of sanitary chemical and toxicological tests"

GOST CT 31214-2016 "Medical devices. Requirements for samples and documentation supplied for toxicological, sanitary chemical testing, sterility and pyrogenicity testing"

GOST 31209-2003 "Containers for blood and blood components. Chemical and biological safety requirements and testing methods"

p. 5.3 GOST 31576-2012 "Biological evaluation of medical devices. Classification and sample preparation"

GOST 14192-96 Labelling of goods;

GOST R 50444-92. Medical instruments, apparatus and equipment. General specifications.

22. Conditions of use, storage and transportation

Table 9. Conditions of use

Parameter	Value
Temperature	5 °C - 42 °C
Humidity	Up to 98% RH at 25 °C
Air pressure	70 kPa - 106 kPa

The devices are used inside the oral cavity containing biological liquids at nominal temperatures: 32 °C - 42 °C, for accessories to be used outside the oral cavity: 10 °C - 35 °C.

The devices are not intended to be repaired. No specific maintenance is required. The devices must be used only by qualified medical workers in clinical facilities.

Table 10. Storage conditions.

Parameter	Value
Temperature	2°C - 24 °C
Humidity	Not more than 70%
Air pressure	70 kPa - 106 kPa

After being stored in the refrigerator, as well as after transportation at sub-zero temperatures, the device should be warmed to room temperature before use.

The devices in original packaging should be stored in closed, dry and heated storage facilities protected from light. **Do not freeze.**

Inadequate storage conditions will shorten the shelf life and may lead to malfunction of the product.

- Keep out of direct sunlight and store in a well-ventilated place at temperatures between 2-24 °C/35-75 °F.
- Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.

The device can be transported by any type of closed vehicle in accordance with the transport regulations valid for particular type of vehicle.

Placement and fastening of boxes containing the devices should ensure their stable position, excluding the possibility of displacement and collision with each other or with the walls of the vehicle.

Table 11. Transport conditions

Parameter	Value
Temperature	2°C - 24 °C
Humidity	Not more than 70%
Air pressure	70 kPa - 106 kPa

The devices should be transported in accordance with the regulations valid for particular type of vehicles at temperature in the range between 2 °C and 24 °C and RH < 70%.

When transporting to the countries with unacceptable ambient temperature, the devices should be placed inside special containers which maintain the temperature and relative humidity in the specified range.

23. Shelf life

- 1) SDR® Plus Bulk Fill Flowable material in compula tips 0.25 g, (SDR® Plus Compula® Tips), shades: universal shade, shades A1, A2, A3 – 3 years;
- 2) SDR® Plus Bulk Fill Flowable material in syringes 1 g each (universal shade) – 3 years;
- 3) Syringe tips – unlimited;
- 4) SDR® storage box – 5 years.

24. Procedure and conditions for disposal of the medical device

The disposal of the medical device should be performed according to the RF SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the management of medical waste" sanitary rules and regulations. The medical waste hazard class is B.

The devices should be used only by qualified medical workers in clinical facilities. Upon disposal of apparatus adhere to the local regulations for disposal and processing of waste.

25. Warranty

The manufacturer bears no responsibility and pays no compensation for possible damage and accidents caused by:

- Use of components not included in the kit and which may interfere with the normal operation.
- Non-compliance with the instructions for use.

The user bears responsibility for testing the material for fitness for particular purpose not explicitly indicated in the instructions for use!

Please contact the official representative of the manufacturer in case of questions pertaining to the quality of the product.

26. Official representative

Dentsply Sirona LLC, Prospect Andropova 18, build. 6, offices 1 39 40, Moscow, Russia 115432, Tel.: +7 (495) 725 10 87

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