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Switzerland

INSTRUCTION FOR USE.

Ophthalmic surgical instruments: Handles:

1. Sutherland Handle NG, red
2. Sutherland Handle NG, blue
3. RENAISSANCE Handle NG, red
4. RENAISSANCE Handle NG, blue
5. Revolution Handle, red
6. Revolution Handle, blue
7. RENAISSANCE Advanced Handle

Manufacturer:

 Alcon Grieshaber AG

Winkelriedstrasse 52

CH-8203 Schaffhausen

Switzerland

INSTRUCTION FOR USE 1.

Ophthalmic surgical instruments: Handles:

- 1.Sutherland Handle NG, red
- 2.Sutherland Handle NG, blue
- 3.RENAISSANCE Handle NG, red
- 4.RENAISSANCE Handle NG, blue
- 5.Revolution Handle, red
- 6.Revolution Handle, blue

1. Description and Intended Use.

The Sutherland Handle NG, red; Sutherland Handle NG, blue; RENAISSANCE Handle NG, red; RENAISSANCE Handle NG, blue; Revolution Handle, red and Revolution Handle, blue (further – GRIESHABER Reusable Handles) are designed to be used together with GRIESHABER DSP Forceps Tips or Scissors Tips. These instruments are used in ophthalmic surgery. Different tip designs are available for cutting, grasping, ILM handling, etc. The GRIESHABER Reusable Handles are available in different colors to simplify recognition in the OR theatre.

These instruments are produced according to the latest methods and standards and on the highest level of function, safety and design.

The use of these instruments is obvious to specialized personnel.

2. Precautions and Warnings.

- The instrument may be used only by well-trained physicians and personnel.
- Prevent instrument from touching any hard surface, it could get damaged.
- Prior to each use, the instrument must be cleaned and sterilized.
- The instrument tip must be properly attached to the handle: Check that the dot on the tip does align with the mark on the handle and that the tip fits tightly.
- After attaching the instrument tip to the handle, but prior to each use, it must be thoroughly inspected for correct function, wear and damage.
- Sterilization instructions provided and validated by Alcon Grieshaber (see Chapter 4.4 Sterilization) are capable of sterilizing the instrument(s) for re-use. It remains the responsibility of the processor to insure that the processing as actually performed using equipment, materials and personnel in the facility achieve the desired result. This requires verification and routine monitoring of the process.

Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Please refer to nationally recognized standards, such as AAMI standards (see References 2, 3 and 4 in Chapter 4.4) or to your facility's standard procedures.

NOTE: Due to the potential for the accumulation of particulate and bio burden residues in the sterilizer water reservoirs, it is the surgical facility's responsibility to properly maintain the equipment and their associated filters to ensure the introduction of contaminant-free steam onto the device.

- To prevent from TASS (Toxic Anterior Segment Syndrome) thoroughly follow the steps of rinsing and flushing with water to remove any remnants or the reprocessing solutions.
- Refer to the safety instructions for the use of cleaning and disinfecting solutions and to the hospital's sterilization Instructions.
- The end or product life is recognized by decreased functionality and appearance of wear and tear (damage). The number of reprocessing cycles cannot be specified by number as it depends on many factors.
- The instrument has to be returned to the Customer Support Department of Alcon, if any damage or changes in its characteristics are observed. For this it's necessary to contact with authorized representative of manufacturer in Russian Federation.
- If any Instrument is received in defective condition, Alcon Grieshaber or authorized representative of manufacturer in Russian Federation has to be notified. Never use damaged or defective devices.

3. Directions for Use.

Set-up, prior to surgery:

- Attach tip to handle: Align dot at the tip with the groove in the handle. Slide the tip longitudinally onto the handle. While doing so, do not compress the levers; hold the handle at its rear part.
- Twist tip clockwise until the dot of the tip aligns with the marking on the handle. While doing this, two clicks will be noticeable.
- Check tip for tight fit.

The device is used at medical institutions of the ophthalmological profile.

4. Cleaning, disinfection, prion inactivation and sterilization.

4.1. General information for reprocessing:

- If in the medical opinion of the physician a patient with a prion related disease undergoes a high risk procedure, the instrument should be destroyed or be processed according to local requirements.
- To provide a long lifetime and smooth function, the following care should be taken:
 - The instrument should be cleaned and dried immediately after surgery. Any dried off remnants especially saline solution, left on/in the instrument could trigger irreversible corrosion and other damage.
 - All applied reprocessing procedures must be compatible with stainless steel and PEEK.
 - The instrument should not get exposed to oxidative or low pH cleaning solutions and especially not to NaOCl. Irreparable damage will result.
 - The instrument should not be exposed to temperatures higher than 140°C (284°F).
 - The instrument can withstand 1 molar NaOH-solution at room temperature 25°C (77°F).
 - The hot instrument should not get chilled in water or other liquid.
 - Allow the instrument to cool by air.
- Definitions: "water" use only deionized or distilled water. "Air" use only air from a syringe or dried and filtered compressed air. For each time water or air is named in the DFU these conditions must be kept.

4.2. Cleaning and disinfection.

Perform the following steps first with cleaning solution and then with disinfection solution.

- Soak the instrument in an ultrasonic bath filled with the relevant solution, for the minimum time recommended by the manufacturer of the solution.
- Rinse the instrument in an ultrasonic bath with water for at least one minute. Repeat this procedure at least twice.

Repeat the cleaning steps if the rinsing water is not clear or residue is present.

4.3. Additional care:

- Drying is recommended, if the instrument is not sterilized immediately. A flush with fresh isopropyl alcohol or ethanol solution (min. 70%) will facilitate drying. Dry the instrument thoroughly with air.
- The periodic use of an instrument milk (lubricant) will increase the lifetime and function.

4.4. Sterilization

- The Instrument must be cleaned and disinfected prior to sterilization.

- If any kind of packaging during sterilization is used, it must be suitable for steam sterilization and comply with EN 868 or AAMI/ISO 11607.
- Do not use packaging material as sterilization container.
- Sterilize the instrument using a steam sterilization cycle validated for this device by the hospital.

Sterilizer Type	Sample Configuration	Temperature	Minimum Exposure Time (Minutes)
Gravity Displacement	Wrapped	132°C (270°F)	20
	Unwrapped		10
Prevacuum pulsing ¹⁾	Wrapped	132°C (270°F)	4
	Unwrapped		
	Wrapped	134°C (273°F)	18 ²⁾
¹⁾ Prevacuum pulsing is based on cycle with four pulses.			
²⁾ European Standards.			

References:

- 1) ISO 17664: Sterilization of medical devices – information to be provided by the manufacturer for the processing of resterilizable medical devices.
- 2) AAMI TIR12: Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers.
- 3) AAMI TIR30: A compendium of processes, materials, test methods and acceptance criteria for cleaning reusable medical devices.
- 4) ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

5. Validation.

Alcon Grieshaber has performed validations to demonstrate the efficacy of the herein recommended procedures with the following setup:

- Enzymatic cleaner: CIDEZYME® by Johnson & Johnson, in 0.8% v/v
- Disinfectant: CIDEX® OPA Johnson & Johnson, undiluted, activated solution
- Autoclave: Amsco Eagle 300 Series

6. Conditions of storage, transportation and utilization. Shelf life.

Storage conditions:

Temperature	+5 °C to +30 °C
Relative humidity	30% to 85%

The end of life for the handles is determined by their wear and damage, but must not be longer than 5 years.

Transport conditions:

Temperature	+5 °C to +60 °C
Relative humidity	30% to 85%

Conditions of utilization

GRIESHABER Reusable Handles should be properly disposed as medical waste.

In case of clinic, warehouse, etc. medical device is utilized in accordance with medical waste circulation rules of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation».

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

7. Definition of symbols.

Consult accompanying documents

REF

Catalog number

LOT

Batch code

SN

Serial number

REUSABLE

Reusable

NON STERILE

Non sterile

Manufacturer

Date of Manufacture

8. Technical parameters.

Dimensions with tolerances of the GRIESHABER Reusable Handles:

Name	L* (mm)	I** (mm)	H (mm)
Revolution Handle, red Revolution Handle, blue	103,2±1,5	101,5±1,5	11±0,5
Sutherland Handle NG, red Sutherland Handle NG, blue	108,35±1,5	106,65±1,5	14±0,5
RENAISSANCE Handle NG, red RENAISSANCE Handle NG, blue	104,1±1,5	102,5±1,5	10±0,5

Notes

*L - indicates handle length in activated position, the pin is protruding.

**I - indicates handle length in not activated position.

Figures 1-3 are schematic only, representing the mechanism in a non-activated state (I) and the activation pin in an activated state (L).


Figure 1: Dimensions of Revolution Handle red and blue


Figure 2: Dimensions of Sutherland Handle NG red and blue


Figure 3: Dimensions of RENAISSANCE Handle NG red and blue

Weight of the GRIESHABER Reusable Handles

Device Name	Weight (g)
Revolution Handle, red Revolution Handle, blue	16
Sutherland Handle NG, red Sutherland Handle NG, blue	10.5
RENAISSANCE Handle NG, red RENAISSANCE Handle NG, blue	10

Technical characteristics of the GRIESHABER Reusable Handles:

Revolution Handle, red and Revolution Handle, blue					
Surface Roughness*		Hardness**		Travel	Actuation ratio***
Activation Pin	Front Part	Stainless Steel 1.4435	Stainless Steel 1.4310	Activation Pin	Handle / Tip
Ra = 0.8 µm	Ra = 0.8 µm	HB ≤ 215	HB ≤ 230	Not more 2 N	4.5/1
Sutherland Handle NG, red and Sutherland Handle NG, blue					
Surface Roughness*		Hardness**		Travel	Actuation ratio***
Activation Pin	Bayonet connector	Stainless Steel 1.4435		Activation Pin	Handle / Tip
Ra = 0.8 µm	Ra = 0.8 µm	HB ≤ 215		Not more 5 N	3.4/1
RENAISSANCE Handle NG, red and RENAISSANCE Handle NG, blue					
Surface Roughness*		Hardness**		Travel	Actuation ratio***
Activation Pin	Core	Stainless Steel 1.4435	Stainless Steel 1.4310	Activation Pin	Handle / Tip
Ra = 0.8 µm	Ra = 0.8 µm	HB ≤ 215	HB ≤ 230	Not more 5 N	3.3/1

* Parts for which the surface roughness has been indicated are selected based on the relevance of surface finish for instrument function.

**** Steel hardness is based on the EN 10088-3:2014 standard. No specific hardness specifications have been defined for the finished subcomponents, because it is not a required characteristic based on their intended use.**

***** The Handle / tip actuation ratio is an value indicating how far the surgeon must move his fingertip to actuate the forceps or scissors tip. For example 4.5/1 means 4.5 mm movement with the fingertip that translates to 1mm actuation of the tip and complete closure of the forceps or scissors tip**

9. Repair service. Warranties

If any instrument is received in defective condition, Alcon Grieshaber or the Authorized representative of manufacturer in Russian Federation has to be notified.

Do not use damaged or defective instruments.

The device must be used in accordance with indications and DFU. All questions regarding with the device's function can be addressed to Authorized representative of manufacturer. All repairs work for GRIESHABER Reusable Handles is carried out exclusively by Alcon Grieshaber. Contact Authorized representative in Russian Federation for assistance on returning devices.

Only cleaned, disinfected and sterilized instruments may be sent for repair.

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

Alcon Farmaceutika LLC, Russia,
109004, Russian Federation,
Moscow, Nikoloyamskaya str., 54
Telephone +7 495 961-13-33

Manufacturer:

Alcon Grieshaber AG

Winkelriedstrasse 52

CH-8203 Schaffhausen

Switzerland

INSTRUCTION FOR USE 2.

Ophthalmic surgical instruments: Handles:

RENAISSANCE Advanced Handle.

1. Description and Intended Use.

The RENAISSANCE Advanced Handle has been designed to be used with the GRIESHABER Disposable Forceps Tips or Scissors Tips. These instruments are used in ophthalmic surgery. Different tip designs are available for cutting, grasping, ILM handling, etc.

These instruments are produced according to the latest methods and standards and on the highest level of function, safety and design.

The use of these instruments is obvious to specialized personnel.

2. Precautions and Warnings.

- The instrument may be used only by well-trained physicians and personnel.
- The instrument tip must be properly attached to the handle: Check that the dot on the tip does align with the mark on the handle and that the tip fits tightly.
- After attaching the instrument tip to the handle, but prior to each use, it must be thoroughly inspected for correct function, wear and damage.
- Prevent instrument from touching any hard surface, it could get damaged.
- Prior to each use, the instrument must be cleaned and sterilized.
- Prior to each use, the instrument must be thoroughly inspected for correct function, wear and damage.
- Sterilization instructions provided and validated by Alcon Grieshaber (see Chapter 4.4 Sterilization) are capable of sterilizing the instrument(s) for re-use. It remains the responsibility of the processor to insure that the processing as actually performed using equipment, materials and personnel in the facility achieve the desired result. This requires verification and routine monitoring of the process.

Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Please refer to nationally recognized standards, such as AAMI standards (see References 2, 3 and 4 in Chapter 4.4) or to your facility's standard procedures.

NOTE: Due to the potential for the accumulation of particulate and bio burden residues in the sterilizer water reservoirs, it is the surgical facility's responsibility to properly maintain the equipment and their associated filters to ensure the introduction of contaminant-free steam onto the device.

- To prevent from TASS (Toxic Anterior Segment Syndrome) thoroughly follow the steps of rinsing and flushing with water to remove any remnants or the reprocessing solutions.
- Refer to the safety instructions for the use of cleaning and disinfecting solutions and to the hospital's sterilization Instructions.
- The end or product life is recognized by decreased functionality and appearance of wear and tear (damage). The number or reprocessing cycles cannot be specified by number as it depends on many factors.
- The instrument must be returned to the Customer Support Department of Alcon, if any damage or changes in its characteristics are observed. For this it's necessary to contact with authorized representative of manufacturer in Russian Federation
- If any Instrument is received in defective condition, Alcon Grieshaber or authorized representative of manufacturer in Russian Federation has to be notified. Never use damaged or defective devices.

3. Directions for Use.

Set-up, prior to surgery:

- Attach tip to handle: Align dot at the tip with the groove in the handle. Slide the tip longitudinally onto the handle. While doing so, do not compress the levers, hold the handle at its rear part.
- Twist tip clockwise until the dot of the tip aligns with the marking on the handle. While doing this, two clicks will be noticeable.
- Check tip for tight fit.

The device is used at medical institutions of the ophthalmological profile.

4. Cleaning, disinfection, prion inactivation and sterilization.

4.1. General information for reprocessing:

- If in the medical opinion of the physician a patient with a prion related disease undergoes a high risk procedure, the instrument should be destroyed or be processed according to local requirements.
- To provide a long lifetime and smooth function, the following care should be taken:
 - The instrument should be cleaned and dried immediately after surgery. Any dried off remnants especially saline solution, left on/in the instrument could trigger irreversible corrosion and other damage.
 - All applied reprocessing procedures must be compatible with stainless steel, PPSU and PEEK Polymer.

- The instrument should not get exposed to oxidative or low pH cleaning solutions and especially not to NaOCl. Irreparable damage will result.
- The instrument should not be exposed to temperatures higher than 140°C (284°F).
- The instrument can withstand 1 molar NaOH-solution at room temperature 25°C (77°F).
- The hot instrument should not get chilled in water or other liquid, - allow cooling by air.
- Definitions: "water" use only deionized or distilled water. "Air" use only air from a syringe or dried and filtered compressed air. For each time water or air is named in the DFU these conditions must be kept.

4.2. Cleaning and disinfection.

Perform the following steps first with cleaning solution and then with disinfection solution.

- Soak the instrument preferably in an ultrasonic bath filled with the relevant solution, for the minimum time recommended by the manufacturer of the solution.
- Rinse the instrument preferably in an ultrasonic bath with freshly filled water for at least one minute. Repeat this procedure at least twice.

Repeat the cleaning steps if the rinsing water is not clear or residue is present.

4.3. Additional care:

- Drying is recommended, if the instrument is not sterilized immediately. A flush with fresh isopropyl alcohol or ethanol solution (min. 70%) will facilitate drying. Dry the instrument thoroughly with filtered compressed air.
- The periodic use of an instrument milk (lubricant) will increase the lifetime and function.

4.4. Sterilization.

- The Instrument must be cleaned and disinfected prior to sterilization.
- If any kind of packaging during sterilization is used, it must be suitable for steam sterilization and comply with EN 868 or AAMI/ISO 11607.
- Do not use packaging material as sterilization container.
- Sterilize the instrument using a steam sterilization cycle validated for this device by the hospital.

Sterilizer Type	Sample Configuration	Temperature	Minimum Exposure Time (Minutes)
Gravity Displacement	Wrapped	132°C (270°F)	20
	Unwrapped		10
Prevacuum pulsing ¹⁾	Wrapped	132°C (270°F)	4
	Unwrapped		
	Wrapped	134°C (273°F)	18 ²⁾
¹⁾ Prevacuum pulsing is based on cycle with four pulses.			
²⁾ European Standards.			

References:

- 1) ISO 17664: Sterilization of medical devices – information to be provided by the manufacturer for the processing of resterilizable medical devices.
- 2) AAMI TIR12: Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers.
- 3) AAMI TIR30: A compendium of processes, materials, test methods and acceptance criteria for cleaning reusable medical devices.
- 4) ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

5. Validation.

Alcon Grieshaber has performed validations to demonstrate the efficacy of the herein recommended procedures with the following setup:

- Enzymatic cleaner: CIDEZYME® by Johnson & Johnson, in 0.8% v/v
- Disinfectant: CIDEX® OPA Johnson & Johnson, undiluted, activated solution
- Autoclave: Amsco Eagle 300 Series

6. Conditions of storage, transportation and utilization. Shelf-life.

Storage conditions:

Temperature	+5 °C to +30 °C
Relative humidity	30% to 85%

The end of life for the handles is determined by their wear and damage, but must not be longer than 5 years.

Transport conditions:

Temperature	+5 °C to +60 °C
Relative humidity	30% to 85%

Conditions of utilization.

RENAISSANCE Advanced Handle should be properly disposed as medical waste.

In case of clinic, warehouse, etc. medical device is utilized in accordance with medical waste circulation rules of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation».

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

7. Definition of symbols.

 Consult accompanying documents

REF Catalog number

LOT Batch code

SN Serial number

REUSABLE Reusable

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 Manufacturer

 Date of Manufacture

8. Technical parameters.

Dimensions with tolerances of the RENAISSANCE Advanced Handle

Name	L* (mm)	I** (mm)	H (mm)
RENAISSANCE Advanced Handle	107,1±1,5	106,1±1,5	9±0,5

Notes

*L - indicates handle length in activated position, the pin is protruding.

**l - indicates handle length in not activated position.

Figure 1 is schematic only, representing the mechanism in a non-activated state (l) and the activation pin in an activated state (L).


Figure 1: Dimensions of RENAISSANCE Advanced Handle

Weight of the RENAISSANCE Advanced Handle

Device Name	Weight (g)
RENAISSANCE Advanced Handle	10

Technical characteristics of the RENAISSANCE Advanced Handle:

Surface Roughness*		Hardness**		Travel	Actuation ratio***
Activation Pin	Core	Stainless Steel 1.4435	Stainless Steel 1.4310	Activation Pin	Handle / Tip
Ra = 0.8 µm	Ra = 0.8 µm	HB ≤ 215	HB ≤ 230	Not more 2 N	3/1

* Parts for which the surface roughness has been indicated are selected based on the relevance of surface finish for instrument function.

** Steel hardness is based on the EN 10088-3:2014 standard. No specific hardness specifications have been defined for the finished subcomponents, because it is not a required characteristic based on their intended use.

*** The Handle / tip actuation ratio is an value indicating how far the surgeon must move his fingertip to actuate the forceps or scissors tip. 3/1 means 3 mm movement with the fingertip that translates to 1mm actuation of the tip and complete closure of the forceps or scissors tip.

9. Repair service. Warranties.

If any instrument is received in defective condition, Alcon Grieshaber or the Authorized representative of manufacturer in Russian Federation has to be notified.

Do not use damaged or defective instruments.

The device must be used in accordance with indications and DFU. All questions regarding with the device's function can be addressed to Authorized representative of manufacturer.

All repairs work for RENAISSANCE Advanced Handle is carried out exclusively by Alcon Grieshaber. Contact Authorized representative in Russian Federation for assistance on returning devices.

Only cleaned, disinfected and sterilized instruments may be sent for repair.

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

Alcon Farmaceutika LLC, Russia,
109004, Russian Federation,
Moscow, Nikoloyamskaya str., 54
Telephone +7 495 961-13-33

Manufacturer:

 **Alcon Grieshaber AG**
Winkelriedstrasse 52
CH-8203 Schaffhausen
Switzerland

CE

Alcon Grieshaber AG
Winkelriedstrasse 52
CH-8203 Schaffhausen

Alcon
a Novartis company

TOTAL NUMBERED, BOUND
AND SEALED 35 PAGES

Alcon Grieshaber AG
Winkelriedstrasse 52
8203 Schaffhausen

APOSTILLE

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1. Country: Swiss Confederation

This public document

2. has been signed by Elisabeth Holzer

3. acting in the capacity of Urkundenbeamtin

4. bears the seal/stamp of Stadtkanzlei

Schaffhausen

Certified

5. at Schaffhausen

6. the 18.05.2015

7. by Staatskanzlei Kanton Schaffhausen

8. N°: 569/2015

9. Seal / stamp:

10. Signature:

A. Bär

Legalization

Seen for legalization of the signature

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