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Approved by
Technical Article Manager Heavy Incontinence & Skin Care
Health and Medical Solution
Essity Hygiene and Health AB
Susanne WEDEL

Essity Hygiene and Health AB
405 03 GÖTEBORG
SWEDEN

Susanne Wedel

Medical device

Breathable diapers with triple protection for incontinent patients TENA® Slip Plus
(Medium, Large), TENA® Slip Super (Medium, Large)

Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden

Instruction for use

2020

I, the undersigned, Gabriel Rudbeck, Notary Public
in the City of Stockholm, Sweden, certify that

Lena Susanne Wedel---

has/have issued and signed the foregoing document.

Fee

Stockholm, 2020-12-10

Crowns

Ex officio

Gabriel Rudbeck
Notary Public

Advokat-Notarjus Gabriel Rudbeck

Tel: 0046-8-20 59 90 e-mail: info@notarjus.gabrielrudbeck.se

Sveavägen 24-26, Box 3098, 103 61 Stockholm

[Handwritten signature]

1. Name of MD.

Breathable diapers with triple protection for incontinent patients TENA® Slip Plus (Medium, Large), TENA® Slip Super (Medium, Large)

2. Manufacturer and developer of MD: Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden

Production site: Essity Operations Poland Sp. z o.o., Ul. 3 Maja 30 A, 55-200 Olawa, Poland

3. Intended use of MD

Breathable diapers with triple protection for incontinent patients TENA® Slip Plus (Medium, Large), TENA® Slip Super (Medium, Large) are single-used multi-layer absorbent products with a cellulose absorbent layer containing gel-forming moisture-absorbing substances (superabsorbents) used to absorb and retain human urine. Purpose: designed to absorb and retain urine, are used to care for patients with urinary incontinence of heavy degree, in-patients, as well as in other cases.

Group of breathable diapers TENA® Slip Plus (Medium, Large) refer to size range for waist Medium 80-122 cm ; Large 96-144 cm. Absorption level for Plus subgroup is 6 "drops", indicated on consumer pack.

Group of breathable diapers TENA® Slip Super (Medium, Large) refer to size range for waist Medium 80-122 cm ; Large 96-144 cm. Absorption level for Super subgroup is 7 "drops", indicated on consumer pack.

The product is non-sterile, this is a medical device with a low degree of risk - class 1.

This medical device briefly contacts with intact human skin. Change the product no later than after 8 hours of continuous use (in case of polyuria syndrome - no more than 5 hours) or with full filling of the product. It is recommended to change the product immediately after defecation.

4. Indications for use of MD

Incontinence.

5. Contraindications for use of MD

Contraindications: individual intolerance to the materials of the product.

6. Side effects

Allergic reaction due to patient intolerance to the materials of the product.

7. Precautions

N/A

8. Information about potential consumers of MD

This medical device is intended for people suffering from incontinence of heavy degree. The product can be used in a home-care environment or at institutions.

9. Functional characteristics of MD and indications for use

Incontinence of heavy degree.

Name of MD	Total moisture absorbency (g) Rothwell not less
TENA® Slip Plus Medium	2100
TENA® Slip Plus Large	2450
TENA® Slip Super Medium	2350
TENA® Slip Super Large	2650

Name of MD	Sorption capacity after centrifugation (g)	Resorption, (rewetting) (g)	Absorption rate, (inlet time) (cm ³ /s)
TENA® Slip Plus Medium	Average 670 StDev 20,6	Average 0,7 StDev 0,37	Average 10,6 StDev 0,49
TENA® Slip Plus Large	Average 860 StDev 31,0	Average 0,5 StDev 0,13	Average 4,5 StDev 0,23
TENA® Slip Super Medium	Average 890 StDev 38,0	Average 0,2 StDev 0,03	Average 4,8 StDev 0,42
TENA® Slip Super Large	Average 1030 StDev 32,9	Average 0,3 StDev 0,07	Average 4,6 StDev 0,46

10. Risks of using a medical product, contraindications, expected and predictable side effects associated with the use of a medical product for the intended purpose.

None

11. Technical characteristics of MD


Measurement point	TENA® Slip Plus (Medium, Large), TENA® Slip Super (Medium, Large)	
	Medium (mm) ±mm	Large (mm) ± mm
Length of diaper (A)	800 ±15	940 ±15
Width of diaper (B)	690 ±15	850 ±15
Length of 1 st /2 nd absorption layer (C/D)	630/400 ±10	715/400 ±10
Width of 1st/2nd absorption layer (E/F)	240/150 ±10	255/150 ±10
Length of transfer layer (G)	270 ±20	270 ±20

Width of transfer layer (H)	90 ±10	90 ±10
Length of back panel (I)	200 ±10	250 ±10
Width of back panel (J)	210 ±5	290 ±5
Length of front panel (K)	170 ±10	200 ±10
Width of front panel (L)	170 ±5	230 ±5

Weight of MD.

Name of MD	Incontinence degree	Weight (g) not less
TENA® Slip Plus Medium	heavy	87 ±10%
TENA® Slip Plus Large	heavy	100 ±10%
TENA® Slip Super Medium	heavy	97,5 ±10%
TENA® Slip Super Large	heavy	110 ±10%

Scheme of MD

12. How to use MD

• Non - mobile

Remove the product from the bag and open the front and back side panel. With the user lying on their side, centre and smooth out the pad in the back and slide the pad between the legs from back to front. Make sure good fit in the groin. Centre and smooth out the pad in front. Fit the tapes, one at the time. The top tapes angled slightly downwards, and the bottom tapes angled slightly upwards. Check for comfort and fit. The pad allows for multiple refastening and adjustments of the tapes. Gently press and smooth the tapes on to the textile - like back -sheet. Give a gentle pull to anchor the tapes into the fibres of the back- sheet. This ensures the secure fixation.

• Semi-mobile

Remove the product from the bag and open the front and back side panel. With the user standing place the pad between the legs and ensure a good fit in the groin. Centre and smooth out the pad in front and back. Fit the tapes, one at the time. The top tapes angled slightly downwards, and the bottom tapes angled slightly upwards. Check for comfort and fit. The pad allows for multiple refastening and adjustments of the tapes. Gently press and smooth the tapes on to the textile-like back- sheet. Give a gentle pull to anchor the tapes into the fibres of the back -sheet. This ensures the secure fixation

13. Handling conditions

Requirements for environmental conditions during handling	
Air temperature, °C	+ 15°C - +35°C
Air relative humidity	not more 40%

14 . Materials used.

№	Component	Material/Brand	Supplier	MSDS (№; date)
1	Top sheet	Nonwoven	BERRY Spain	Oct 17, 2018
2	Standing gather	Nonwoven PC5FW-111/PC5FW113	Fitesa Germany	March 2, 2017
3	Transfer layer	Nonwoven highloft TER 50-12050W50A	Fiberweb Terno d'Isola S.r.l.	02 February 2016
4	Absorption core, Fluff	Pulp, BSK	GP Cellulose GmbH	May 07-2015
5	SAP	K-SAM, Sodium Polyacrylate, cross-linked	LG Chem, Ltd.	2013. 10. 22.
6	Front and back panels	Nonwoven	ATEX SpA	3/2020
7	Backsheet film	PP/CaCo3 film Type 12000 Aptra B 140	RKW	1/31/2018
8	Back-sheet nw	Spunbond nw	Fibertex	July 2nd 2018
9	Fixing elements	Y5965	Avery Dennison	28/10/14
10	Elastic threads	Spandex/elastane fiber	A&AT UK Limited	21 March 2018
11	Waist elastic film	CEX 836 WR/CEX-836WRT	Tredeggar Corporation	14 May 2019
12	Wetness indicator	KMELT PERFECTCARE W 3488	Colquimica S.A.	12/02/2016
13	Glue	KMELT PERFECTCARE W 2780	Colquimica S.A.	16/02/2018
14	Glue	NW 1002 ZP	HB Fuller	05.03.2015
15	Glue	Full-Care 8515 -NZP	HB Fuller	04/11/2016
16	Glue	Technomelt DM 4500 UV DISPOMELT	Henkel AG&Co.KGaA	17.09.2018

15.Labeling of consumer packs.

Labelling of the product includes model and absorption capacity.

Consumer pack is made of a polymer film.

Information on consumer pack:

Consumer pack contains the following info:

- product name
- inteded use
- size of MD
- size range for waist

- pack count;
- name and address of manufacturer (seller, supplier), trademark
- production site address;
- date of production;
- shelf life;
- storage conditions
- instructions for disposing of diapers
- rules how to use diapers (in the form of drawings);
- article
- Lot number;
- date and number of RC;
- bar code .

Graphic symbols on consumer pack

Symbol	Meaning
	The possibility of processing low density polyethylene (LDPE).
	Symbol of voluntary certification
	Throw in the ballot box
	Prohibition of re-use
	Do not dispose into drains
	Absorbency index
	CE Marking
	Standard icons for EU

Proskin	The name of a marketing concept that unites a group of products for people with incontinence
---------	--

16. Information about transportation of MD

The consumer packaging with the diapers enclosed in it should be complete. On the surface tears and other damage are not allowed. Temperature -30°C to +35°C; Air relative humidity no more than 80%. Diapers are transported only in a dry, clean transport that does not have extraneous odors.

17. Storage and Shelf life

The product in the packaging of the enterprise should be stored in a cool place - clean, dry, well ventilated areas; It is strictly forbidden to store products near sources of sharp odors (for example: paintwork materials, household chemicals, burns, solvents, etc., etc.); When storing and operating the product, direct sunlight and sudden changes in temperature should be avoided.

Environmental conditions during storage

Requirements for environmental conditions during storage	
Air temperature, °C	+ 15°C - +35°C
Air relative humidity	no more 40%

Shelf life of the diapers is 3 years from the production date indicated on the package.

18.UTILIZATION

The product must be disposed of as medical waste of hazard class "B"
Do not throw into drains.

19.Manufacturer's warranty

The manufacturer guarantees the quality of diapers under the conditions of operation, storage and transportation within 3 years from the date of manufacture date indicated on the package.

20.Authorized representative in Russia

During the warranty period, claims are made to the manufacturer at the address:
Russia, 117218 Moscow, ul.Krzhizhanovskogo, 14, building 3, phone: +7 (495) 967-33-67.

21. List of applicable national standards

Harmonized and other International Standards that are relevant to prove conformity to Medical Devices Regulation (EU) 2017/745 by Essity Hygiene and Health AB applied in full or in part.

ISO 13485:2016	Medical devices - Quality management systems – Systems requirements for regulatory purposes This international standard specifies requirements for a quality management system that can be used by an organization for the design and development, production, installation and servicing of medical devices, and the design, development, and provision of related services.
ISO 9001:2015	Quality management systems This international standard specifies requirements for a quality management system where an organization needs to demonstrate its ability to consistently provide product that meets customer and applicable regulatory requirements, and aims to enhance customer satisfaction through the effective application of the system, including processes for continual improvement of the system and the assurance of conformity to customer and applicable regulatory requirements.
ISO 14971:2012	Risk management This international standard is a framework for effective management by the manufacturer of the risks associated with the use of medical devices. It includes a procedure to investigate, using available information, the safety of a medical device, including accessories, by identifying hazards and estimating the risks associated with the device.
EN ISO 10993:2009	Biological evaluation a) the fundamental principles governing the biological evaluation of medical devices the definition of categories of devices based on the nature and duration of contact with the body; the selection of appropriate tests.
EN ISO 17025:2005	General Requirements for the competence of testing and calibrating laboratories This international standard specifies the general requirements for the competence to carry out tests and/or calibrations, including sampling. It covers testing and calibration performed using standard methods, non-standard methods and laboratory-developed methods.
EN 1041:2008	Information supplied by the manufacturer of medical devices This standard specifies the information to be supplied by a manufacturer for different categories of medical devices, as required by the relevant EEC Directives. It does not specify the language to be used for such information.
ISO 15223:2016	Medical devices – Symbols to be used with medical device labels, labeling and information supplied.

	This international standard identifies symbols conventionally used to convey information essential for proper use to the user and others for safe and effective use of medical devices.
ISO 7000:2012	Graphical symbols for use on equipment – Index and synopsis. This international standard provides a synopsis of graphical symbols, which are placed on equipment or parts of equipment of any kind in order to instruct the persons handling the equipment as to its use and operation.
EN ISO 14155:2011	Clinical investigation of medical devices for human subjects This international standard defines procedures for the conduct and performance of clinical investigations of medical devices. It specifies general requirements intended to --protect human subjects, --ensure the scientific conduct of the clinical investigation, --assist sponsors, monitors, investigators, ethics committees, regulatory authorities and bodies involved in the conformity assessment of medical devices. Part 1 of ISO 14155 a) specifies requirements for the conduct of a clinical investigation such that it establishes the performance of the medical device during the clinical investigation intended to mimic normal clinical use, reveals adverse events under normal conditions of use, and permits assessment of the acceptable risks having regard to the intended performance of the medical device, b) specifies requirements for the organization, conduct, monitoring, data collection and documentation of the clinical investigation of a medical device, c) is applicable to all clinical investigation(s) of medical devices whose clinical performance and safety is being assessed in human subjects. Part 2 of this standard provides requirements for the preparation of a Clinical Investigation Plan (CIP) for the clinical investigation of medical devices. The compilation of a CIP in accordance with the requirements of this standard and adherence to it will help in optimising the scientific validity and reproducibility of the results of a clinical investigation.
MEDDEV 2.4/1 rev. 9	Guidelines to the Classification of Medical Devices
MEDDEV 2.7.1 rev. 4	Guidelines on Medical Devices, Clinical Evaluation
MEDDEV 2.12/1 rev.8	Guidelines on a Medical Device Vigilance System
Commission Decision 2000/532/EC	European List of Waste

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